
Reference Manual for

METTLER TOLEDO Standard Interface Command Set (MT-SICS) MT-SICS 0 version 2.30 MT-SICS 1 version 2.20 MT-SICS 2 for Moisture Analyzer MJ33 version 2.30 MT-SICS 3 for Moisture Analyzer MJ33 version 1.30



Contents

	Overview of all commands	4
1	Introduction	6
2	Basic information on data interchange	9
2.1	Command formats	9
2.2	Response formats	10
2.2.1	Format of the response with weight value	10
2.2.2	Format of the response without weight value	11
2.2.3	Error messages	12
2.2.4	Tips for the programmer	12
3	Commands and responses	13
3.1	Commands and responses MT-SICS level 0	13
3.2	Commands and responses MT-SICS level 1	21
3.3	Commands and responses MT-SICS level 2	
	for Moisture Analyzer MJ33	23
3.4	Commands and responses MT-SICS level 3	
	for Moisture Analyzer MJ33	25
4	System configuration (MJ33 – computer)	37
5	What if?	38

Overview of all commands

Commands and responses MT-SICS level 0

10	Inquiry of all implemented MT-SICS commands	14
11	Inquiry of MT-SICS level and MT-SICS versions	15
12	Inquiry of instrument data	16
13	Inquiry of SW version and type definition number	16
14	Inquiry of serial number	16
15	SW-Identification number	17
S	Send stable weight value	17
SI	Send weight value immediately	18
SIR	Send weight value immediately and repeat	18
Z	Zero	19
ZI	Zero immediately	20
@	Reset	20

Page

25

Commands and responses MT-SICS level 1 for Moisture Analyzer MJ33

D	Display	21
DW	Weight display (Display show weight)	22

Commands and responses MT-SICS level 2 for Moisture Analyzer MJ33

DAT	Date	23
PWR	Power on/off	24
TIM	Time	24

Commands and responses MT-SICS level 3 for Moisture Analyzer MJ33

Control commandsHA01Reset application / escapeHA02Set factory settings

HAO2	Set factory settings	25
HAO3	Switch keypad on/off	25
HAO5	Start / end drying	26
HAO6	Trigger audio signal	26
HAO7	Report instrument status change	27
HAO8	Request printer records	28

Status inquiries

HA20	Inquiry of instrument status	
HA21	Inquiry of heating module position	29
HA22	Inquiry of last balance adjustment	29
HA23	Inquiry of last heating module adjustment	30
HA24	Inquiry of temperature	30
HA25	Inquiry of drying weights	31
HA26	Inquiry of drying data	32
HA27	Inquiry of drying result	33
Instrument se	ettings	
HA40	Inquiry / setting of language	34
HA40X	Inquiry / setting of menu parameters	35

1 Introduction

To enable you to integrate your Moisture Analyzer in your electronic data system in a simple manner and utilize their capabilities to the full, most Moisture Analyzer functions are available as appropriate commands via the data interface.

Standardization of the commands

All new METTLER TOLEDO balances and Moisture Analyzers support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS), which is divided into 4 levels, depending on the functionality of the equipment:

- MT-SICS level 0 Command set for the simplest balance, e.g. weighing cell
- MT-SICS level 1 Extension of the command set for standard balances, i.e. balances without integrated applications
- MT-SICS level 2 Extension of the command set family
- MT-SICS level 3 Application-specific commands as independent command set, e.g. MT-SICS level 3 for Moisture Analyzers MJ33

A particular distinguishing feature of this concept is that the commands combined in MT-SICS level 0 and 1 are identical for all balances and instruments. Both the simplest weighing balance and a Moisture Analyzer recognize the commands of MT-SICS level 0.

What do the commands of MT-SICS level 0 and 1 offer?

You can use the commands of MT-SICS level 0 and 1 to perform the following operations of the Moisture Analyzer via the interface:

- request weighing results (raw data),
- zero the balance,
- identify MT-SICS implementation (version number),
- identify the Moisture Analyzer (serial number),
- reset the Moisture Analyzer,
- control the display.

The commands of MT-SICS level 2 and 3 for the Moisture Analyzer

All additional higher level functions for the Moisture Analyzer are collected in the commands of MT-SICS level 2 and 3.

When creating your software application, please note that the commands of MT-SICS level 3 apply to your application and can not be supported by every METTLER TOLEDO balance or instrument.

Most commands in this manual are identical for the MJ33 and the other Moisture Analyzers of METTLER TOLEDO.

Additional documentation on data interface

Settings of the interface such as baud rate, number of data bits, parity, handshake protocols and connector pin assignment are described in the operating instructions of the peripheral instrument or cable in question.

How the Moisture Analyzer operates

Your Moisture Analyzer is used to determine the moisture content of virtually any substance. The instrument operates according to the thermogravimetric principle: At the start of the measurement the Moisture Analyzer determines the weight of the sample, the sample is then rapidly heated with the built-in dryer unit and the moisture evaporates. During the drying, the instrument continuously records the weight of the sample and shows the decrease in the moisture. On completion of the drying, the final result of the moisture or dry substance content of your sample is displayed.

Version number of the MT-SICS

Each level of the MT-SICS has its own version number which can be requested with the command 11 from level 0.

This manual describes

MT-SICS level 0, version 2.30 MT-SICS level 1, version 2.20 MT-SICS level 2 for Moisture Analyzer version 2.30 MT-SICS level 3 for Moisture Analyzer version 1.30

You can use the command 11 via the interface to request the MT-SICS level and MT-SICS versions implemented on your moisture analyzer.

Please make sure that the versions implemented on your moisture analyzer agree with those listed above.

During drying the Moisture Analyzer passes through the following instrument statuses:



The following instrument statuses also exist:

- Entry status
- Startup
- Taring
- Weight adjustment
- Temperature adjustment
- Error status

Some functions or commands can be executed only in particular instrument statuses. In the instrument status "Drying", for example, the value in the display can not be overwritten. If a command can not be executed for this reason, the Moisture Analyzer sends an appropriate message.

You will find detailed information on the functions of the Moisture Analyzer in the operating instructions.

2 Basic information on data interchange

Each command received by the Moisture Analyzer via the data interface is acknowledged by a response of the Moisture Analyzer to the transmitter (e.g. computer).

Commands and responses are data strings with a fixed format, and will be described in detail in chapter 3.

2.1 Command formats

Commands sent to the Moisture Analyzer comprise one or more characters of the ASCII character set. Here, the following must be noted:

- Commands must be entered in upper case letters (case sensitive!).
- The possible parameters of the command must be separated from one another and from the command name by a space (ASCII 32 dec., in this description represented as _).
- Each command must be closed by $C_{RL_{F}}$ (ASCII 13 dec., 10 dec.).

The characters C_RL_F , which can be inputted using the Enter or Return key of most entry keypads, are not listed in this description, but it is essential they be included for communication with the instrument.

Example

Command to Moisture Analyzer which writes Hallo into the display:

 D_{\sqcup} "Hallo" The command terminator $C_{R}L_{F}$ is not shown

Comment

The quotation marks "" must be inserted in the entry.


2.2 Response formats

All responses sent by the Moisture Analyzer to the transmitter (e.g. computer) to acknowledge the received command have one of the following formats:

- Response with weight value
- Response without weight value
- Error message

2.2.1 Format of the response with weight value

A general description of the response with weight value is the following.



Comment

 $C_{\text{R}}L_{\text{F}}$ will not be shown in this description.

Example

Response with stable weight value of 0.256 g: SuSuuuuuu 0.256ug

2.2.2 Format of the response without weight value

A general description of the response without weight value is the following. **ID** Response identification

ID 🔟 Statu	s ц Parameters CR LF	
1 characte 1 – x characters	9 r	
ц	Space (ASCII, 32 dec.)	
Status	Interface status of the Moisture Analyzer, see description of the commands and responses	
Parameters	Command-dependent response code	
C _R	Carriage Return (ASCII 13 dec.)	
$L_{\rm F}$	Line Feed (ASCII 10 dec.)	

Comment

 $C_R L_F$ will not be shown in this description.

Example

Response to Du"HALLO" when HALLO appears unabridged in the display: DuA.

2.2.3 Error messages

There are three different error messages. The identification always comprises two characters.

ID CR	ID CR LF			
ID	Error ider Possible	ntification error messages are		
	ES	Syntax error The Moisture Analyzer has not recognized the received command.		
	EL	Logical error The Moisture Analyzer can not execute the received command.		
C _R	Carriage Return (ASCII 13 dec.)			
L _F	Line Feed (ASCII 10 dec.)			

Comment

 C_RL_F will not be shown in this description.

2.2.4 Tips for the programmer

Command and response

You can improve the dependability of your application software by having your program evaluate the response of the Moisture Analyzer to a command. The response is the acknowledgement that the Moisture Analyzer has received the command.

Reset

To be able to start from a definite condition when establishing the communication between Moisture Analyzer and system, you should send a reset command to the Moisture Analyzer.

When the Moisture Analyzer is switched off, meaningless characters may be received or sent.

Quotation marks ""

Quotation marks included in the command must always be entered.

3 Commands and responses

The Moisture Analyzer receives commands from the system (e.g. computer) and acknowledges the command with an appropriate response.

The following sections contain a detailed description of all commands of the command set in alphabetical order with the associated responses. Commands and responses are always closed with C_RL_F . These termination characters are not shown in the following description, but they must always be entered with commands or sent with responses.

3.1 Commands and responses MT-SICS level 0

The commands of MT-SICS level 0 are available with even the simplest balances which support the METTLER TOLEDO Standard Interface Command Set.

l0 In	quiry of all imp	plemented MT-SICS commands
Command	IO	Send list of all implemented MT-SICS commands
Response	I0⊔B⊔x1⊔"1	• Command" x1 = ident of the MT-SICS level where the 1. Command belongs to.
	I0uBux1u"2 : :	• Command" 2nd (next) command implemented
	I0⊔A⊔x1⊔"1	ast Command" Last command implemented. End of the list
	I0UI	The list cannot be sent at present as another operation is taking place
Example		
Command	IO	Send list of commands
Response	I0⊔B⊔0⊔"I0	Level O command "IO" implemented
	I0UBU0U"I1 : I0UBU0U"S" : : I0UAU3U"HA	Level O command "11" implemented : Level O command "S" implemented :

- The 10 command lists all commands implemented in the present software.
- All level 0 commands are listed in alphabetical order before all commands of level 1 etc. This order corresponds the order how the commands are described in this manual.

l1 l	nquiry of MT-	SICS level an	d MT-SICS versions
Command	11	Inquiry of M	T-SICS level and MT-SICS versions
Response	I1⊔A⊔"x1	"ப"x2"ப"x3	"ப"x4"ப"x5"
		x1 = 0	Balance with MT-SICS level 0 (simplest balance)
		x1 = 01	Balance with MT-SICS level 0 and 1 (standard balance)
		x1 = 012	Balance with MT-SICS level 0, 1 and 2 (standard balance with extensions)
		x1 = 03	Balance with MT-SICS level 0 and 3 (simplest balance with a special application)
		x1 = 013	Balance with MT-SICS level 0, 1 and 3 (standard balance with a special application)
		x1 = 0123	Balance with MT-SICS level 0, 1, 2 and 3 (standard balance with extensions and a special application)
		x1 = 3	Application device with MT-SICS level 3 (not necessarily a balance)
		x2	Version of the implemented MT-SICSO commands
		xЗ	Version of the implemented MT-SICS1 commands
		x4	Version of the implemented MT-SICS2 commands
		x5	Version of the implemented MT-SICS3 commands
Example			
Command	I1	Inquiry of M	T-SICS level and versions
Response	I1⊔A⊔"3"	ப"2.30"ப"2	.20"⊔"2.30"⊔"1.30"
		3 Ar	oplication device with MT-SICS level 3

3	Application device with MT-SICS level 3
2.30	Level 0, version V2.30
2.20	Level 1, version V2.20
2.30	Level 2, version V2.30
1.30	Level 3, version V1.30

12 Inquiry of instrument data		
Command	12	Inquiry of instrument data
Response	I2⊔A⊔"text"	Instrument data as "text"
Example		
Command	12	Inquiry of instrument type
Responses	I2⊔A⊔"MJ33⊔Moisture-Analyzer⊔35.010⊔g"	

13 I	nquiry of SW versi	on and type de	finition number
Command	13	Inquiry of Moistu and type definitio	re Analyzer SW version n number
Response	I3⊔A⊔"TEXT"	Moisture Analyze definition numbe	r SW version and type r as TEXT
Example			
Command	13	Inquiry of SW ver number	rsion number(s) and type definition
Response	I3uAu"1.00u4	.10.5.93.43 " 1.00 4.10.5.93.43	Software version Type definition number

14 I	14 Inquiry of serial number		
Command	14	Inquiry of serial number	
Response	I4⊔A⊔"text"	Serial number as "text"	
Example			
Command	I4	Inquiry of serial number	
Response	I4⊔A⊔"012345	6789"	

15 SW-Identification number

Command Responses	I5 I5⊔A⊔"x"	Inquiry of SW-Identification number. SW-Identification number as Text. x: SW-Identification number.
	I5uI	Command understood, not executable at present.
Example Command	15	Inquiry of SW-Identification number.
Response	 I5⊔A⊔"12345	. ,

Comment

• The SW-Identification number is unique for every Software.

S S	end stable wei	ight value
Command	S	Send the current stable weight value
Response	oonse SuSuWeightValueuUnit Current stable weight value	
	SuI	Command not executable (Moisture Analyzer is currently executing another command)
	Sப+	Balance in overload range
	Sப-	Balance in underload range
Example		
Command	S	Send a stable weight value
Response	SuSuuuuu	ப 1.000பத The current, stable weight value is 1.000 g

Comment

• Timeout approx. 30 s.

		una ina manadimta ha
SI S	ena weight val	ue immediately
Command	SI	Send the current weight value, irrespective of balance stability
Response	SபSபWeight	ValueuUnit Stable weight value
	SபDபWeight	ValueuUnit Nonstable (dynamic) weight value
	SuI	Command not executable (Moisture Analyzer is currently executing another command)
	Sப+	Balance in overload range
	Sப-	Balance in underload range
Example		
Command	SI	Send current weight value
Response	ՏաDասասաս	2.907⊔g The current weight value is unstable (dynamic) and is 2.907 g
SIR S	end weight val	ue immediately and repeat
Command	SIR	Send the weight values repeatedly, irrespective of balance stability
Response	SபSபWeight	ValueuUnit Stable weight value
	SபDபWeight	ValueuUnit Nonstable (dynamic) weight value
	SuI	Command not executable (Moisture Analyzer is currently executing another command)
	Sப+	Balance in overload range

Example

Command	SIR	Send current weight values at intervals
Response	SuDuuuuu	12.907ug
	SuDuuuuu	12.850ug
	ՏսՏսսսսս	12.797⊔g
	ՏսՏսսսսս	12.775⊔g
	SuDuuuuuu	12.770⊔g
	••••	The Moisture Analyzer sends stable or nonstable weight values at intervals of 150 ms

Comment

• SIR is overwritten by the commands S, SI and @.

Z Zero			
Command	Z	Zero the Moisture Analyzer	
Response	Z⊔A	Zero setting performed, i.e. stability criterion and zero setting range complied with	
	ZuI	Command not executable as the Moisture Analyzer is not in the relevant instrument status (e.g. drying unit open)	
	Z⊔+	Upper limit of zero setting range exceeded	
	Z ப –	Lower limit of zero setting range exceeded	
Example			
Command	Z	Zero	
Response	ZபA	Zero setting performed	

Comment

• This command is equivalent to pressing the $\rightarrow 0/T \leftarrow$ key.

ZI	Zero immediately	
Command	ZI	Zero immedeately, i.e. stores immediately the current weight value, which can be stable or non stable (dynamic), as zero value.
Response	ZIUS	Zero setting performed, stable weight value
	ZI⊔D	Zero setting performed, non-stable (dynamic) weight value
	ZILI	Zero setting not performed (Moisture Analyzer is currently executing another command)
	ZILL	Command understood but not executable (e.g. certified version of balance)
	ZIu+	Upper limit of zero setting range exceeded
	ZIL-	Lower limit of zero setting range exceeded
Example 1 Command		Zero immediately
Response	ZIUS	Zero setting performed, weight value was stable
Example 2	2	
Command	ZI	Zero immediately
Response	ZILD	Zero setting performed, weight value was dynamic (non-stable)
@	Reset	
Command	e	Resets the interface to the condition found after switching on, but without a zero setting being performed.
Response	I4⊔A⊔"text"	Serial number of the Moisture Analyzer, the Moisture Analyzer is ready for operation.
Example Command	@	
Response	I4⊔A∟"111435	0697" Moisture Analyzer is reset, its serial number is 1114350697.
Commont		

Comment

All current commands are terminated on @ command

3.2 Commands and responses MT-SICS level 1 for Moisture Analyzer MJ33

The commands of MT-SICS level 1 are available with all standard balances which support the METTLER TOLEDO Standard Interface Command Set. With the MJ33 Moisture Analyzer, only the commands D and DW are supported.

D	Display	
Write into d	lisplay	
Command	Du"TEXT"	Write TEXT into Moisture Analyzer display
Response	D⊔A	TEXT appears unabridged left-aligned in the Moisture Analyzer display.
	D⊔R	The end of the text appears in the Moisture Analyzer display, the start is cut off
	DцI	Command not executable
	DuL	Command understood, parameter wrong
Example		
Command	Du"HALLO"	Write HALLO into the Moisture Analyzer display
Response	D⊔A	The full text HALLO appears in the Moisture Analyzer display
Clear displ	ay	
Command	Du" "	Clear Moisture Analyzer display
Response	DபA	Display cleared

Comments

• A display command can be cleared with the DW or Reset command.

• This command can be executed only in the instrument statuses "basic mode", "ready for taring", "weighing-in" and "ready for start".

DW Weight display (Display show Weight)

Command	DW	Switch display to weight mode
Response	DWLA	Display shows the current weight value
	DWuI	Command not executable

Comment

This command can be executed only in the instrument statuses "basic mode", "ready for taring", "weighing-in" and "ready for start".

3.3 Commands and responses MT-SICS level 2 for Moisture Analyzer MJ33

DAT	Date	
Inquiry of d	late	
Command	DAT	Inquiry of current date of the Moisture Analyzer
Response	DAT⊔A⊔dd∟	ப ளைப்ரூர்ரை "ddபாளப்yyyy" represents the date in the format dayபmonthபyear
Set date		
Command	DATuddumm	வ ்ரல்லை Set date in the format "ddபாmபyyyy"
Response	DATLA	Date has been set
	DATLL	Command not executed as the date format was not correct Inquiry of date of the Moisture Analyzer
Example		
Command	DAT	Current date of the Moisture Analyzer is 2 April 2000
Response	DATLAL02L	ــــــــــــــــــــــــــــــــــــــ

- The set date is retained after the reset command "@".
- Admissible years: 1.1.1901 31.12.2099.

PWR P	ower On/Off	
Command	PWRux	Switch Moisture Analyzer On or Off
		x = 0 Set Moisture Analyzer to standby mode
		x = 1 Switch Moisture Analyzer on
Response	PWRLA	Moisture Analyzer has been switched off successfully
	PWR⊔A I4⊔A⊔"text"	Moisture Analyzer with the serial number according to text has been switched on succesfully (see also 14 command)
	PWRLL	Command understood, parameter wrong

- In the standby mode, the interface remains active; but all commands except **PWR**, **HA07**, **HA20** and **@** are answered with EL.
- On switching on, the Moisture Analyzer also sends the serial number (see also 14 command).
- On switching off, all current commands are terminated.

TIM T	ime	
Inquiry of ti	me	
Command	TIM	Send current time of the Moisture Analyzer
Response	TIMuAuhhummu	SS
		"hhummuss" represents the time in the 24-hour format (hours/minutes/seconds)
Set time		
Command	TIMuhhummuss	Set time in 24-hour format (hours_minutes_seconds)
Response	TIMLA	Time has been set, clock running
	TIMLL	Command not executed as the time format is not correct (e.g. 22_67_25)
Example		
Command	TIM	Inquiry of time
Response	ТІМ⊔А⊔22⊔56⊔	11The current time of the Moisture Analyzer is22 hours, 56 minutes and 11 seconds

3.4 Commands and responses MT-SICS level 3 for Moisture Analyzer MJ33

All Moisture Analyzer specified commands are combined in MT-SICS level 3 for Moisture Analyzers HR73, HG53, HB43-S and MJ33.

HAO1 R	Reset application / escape				
Command	HA01	Reset application / escape			
Response	HA01⊔A	Application reset			

Comment

This command has the same effect as the Reset key, see operating instructions of the Moisture Analyzer. It terminates all current commands and activities.

HA02 S	et factory setti	ngs
Command	HA02	Set factory setting of the menu and method parameters
Response	НА02⊔А	Menu and and method parameters set to factory setting

Comment

All menu parameters are reset to factory settings except RS interface settings and language. This command terminates a drying.

HA03 S	Switch keypa	d on/off		
Command	НА03⊔х	x = 0	Keypad of Moisture Analyzer switched off	
		x = 1	Keypad of Moisture Analyzer switched on (factory setting)	
Response	наозца	Commar	Command executed	
	HA03⊔L	Commar	Command understood, parameter wrong	

HA05 S	tart / end dryi	ing
Command	нао5цх	 x = 0 End drying, possible only in instrument status "drying" x = 1 Start drying, possible only in instrument status "ready for start"
Response	НА05⊔А	Command executed
	HA05山I	Command not executable as the Moisture Analyzer is not in the relevant instrument status
	HA05⊔L	Command understood, parameter wrong

Comment

To abort current drying use **HA01** (Reset application).

HA06	Trigger audio signal				
Command	HA06	Trigger audio signal, e.g. at end of drying			
Response	на06ца	Command executed			

HA07 Report instrument status change

Command	HA07ux1	Report each internal status change x1 = 0 Switch off x1 = 1 Switch on
Response	HA07山A	Command executed
	HA07⊔A⊔x1	Status change (see HA20)
		$ \begin{array}{llllllllllllllllllllllllllllllllllll$
		x1 = 10 "Startup" x1 = 11 "Taring" x1 = 12 "Weight adjustment" x1 = 13 "Temperature adjustment" x1 = 101 "Error 1" x1 = 102 "Error 2" x1 = 10n "Error n"
	HA07LL	Parameter wrong (number, value range,)
	HA07山I	Response always available, hence not possible

- aborted with the HA01 command
- see also HA20 command
- also active in standby

HA08 Request printer records

Command	HA08⊔x1	Request printer records: $x1 = 0$ Do not send printer records $x2 = 1$ Send printer records	
Response	на08ца	Command executed	
наовыь Parameter wrong (number, value ran		Parameter wrong (number, value range,)	
	HA08⊔I	Response always available, hence not possible	

Comments

- The printer records use the 8-bit ASCII IBM table 4.
- Regardless of menu setting (see HA403).
- This setting is not stored.
- To reactivate the menu settings, use the "@" command.

HA20 Inquiry of instrument status

Command	HA20	Inquiry of instrument status		
Response	НА20∟А∟х	x = 0 x = 1 x = 2 x = 3 x = 4 x = 5 x = 6 x = 7 x = 10 x = 11 x = 12 x = 13 x = 101	Status: Status: Status: Status: Status: Status: Status: Status: Status: Status: Status:	"Standby" "Basic mode" "Ready for taring" "Weighing in" "Ready for start" "Drying" "End of drying" "Entry" "Startup" "Taring" "Weight adjustment" "Temperature adjustment"
		x = 10n	Status	"Error n", see operating instructions of the Moisture Analyzer

- With the message HA20_A_6 instrument status "End of drying", it is not apparent whether drying was ended correctly or terminated. This is possible only via the command HA25 – Inquiry of drying weight.
- Also active in standby.

HA21 Inquiry of heating module position

Command	HA21	Inquiry of	heating module position
Response	НА21⊔А⊔х	x = 0 x = 1	Heating module closed Heating module open

HA22 li	nquiry of las	st balance	adjustment
Command	HA22	Inquiry	of last successful balance adjustment
Response	HA22LAL	x1⊔x2⊔x3	பx4பx5பx6
·		xl	Number of the successful adjustments
		x2	Day of the last successful adjustment
		xЗ	Month of the last successful adjustment
		x4	Year of the last successful adjustment
		x5	Hour of the last successful adjustment
		x6	Minute of the last successful adjustment
Example			
Command	HA22		
Response	HA22LAL	15山02山04	ப2000ப09ப34
·		been p	of 15 successful balance adjustments have erformed.
		The las	st took place on April 02, 2000 at 9.34.

- The time of the last successful balance adjustment is specified in the 24-hour format.
- Possible years are 1999 ... 2099.
- The counter for the balance adjustments runs to 65535.

HA23 Inquiry of last heating module adjustment

Command	HA23	Inquiry of the last successful heating module adjustments				
Response	НА23⊔А⊔х1⊔	x2:x3:x4:x5:x6x1Number of successful adjustmentsx2Day of the last successful adjustmentx3Month of the last successful adjustmentx4Year of the last successful adjustmentx5Hour of the last successful adjustmentx6Minute of the last successful adjustment				
Example Command	HA23					
Response	НА23∟А∟15∟	HA23UAU15U02U04U2000U09U34 A total of 15 successful heating module adjustments have been performed. The last took place on April 02, 2000 at 9.34.				

- The time of the last successful heating module adjustment is specified in the 24-hour format.
- Possible years are 1999 ... 2099.
- The counter for the heating module adjustments runs to 65535.

HA24 Inquiry of temperature			
Command	HA24	Inquiry of current temperature	
Response	HA24 LALX	Current temperature in °C	
Example			
Command	HA24	Inquiry of current temperature	
Response	HA24⊔A⊔105	The temperature is 105 °C.	

HA25 Inquiry of drying weights

Command	HA25	Inquiry of drying weight of the last or current drying		
Response	HA25uAux1ux2ux3ux4			
		x1 x1 = 0 x1 = 1 x1 = 2 x1 = 3 x2	Drying running Drying ended	
		х3	Current weight or dry weight in grams	
		x4	Drying time (seconds)	
Example 1				
Command	HA25	Inquiry o	of drying weights	
Response	HA25⊔A⊔2	HA25LAL2L12.345L7.890L180 Drying has been ended regularly, wet weight 12.345 g, dry weight 7.890 g, drying time 180 seconds		
Example 2				
Command	HA25	Inquiry o	of drying weights	
Response	НА25⊔А⊔0	5பAப0ப0.000ப0.000ப0 No drying exists, e.g. as the battery was discharged		

Comment

Together with the command $\mathbf{HA07}$ – Report instrument status change – dryings can be shown in parallel on the host.

HA26 Inquiry of drying data

Command	HA26ux1	Inquiry o	f drying data in configurable display mode	
		x1 = 0	currently set display mode	
		x] =]	Grams	
		$x_1 = 2$		
		x1 = 3 x1 = 4	MC (moisture content), (factory setting) AM (ATRO moisture content)	
Response	HA26uAux1ux2ux3ux4ux5ux6			
		x1	Drying status	
		x1 = 0	No drying exists	
		x1 = 1	Drying running	
		x1 = 2	1 0	
		x1 = 3	Drying terminated	
		x2	Display mode	
		$x^2 = 1$	Grams	
		$x^2 = 2$		
		$x^2 = 3$ $x^2 = 4$	MC (moisture content), (factory setting) AM (ATRO moisture content)	
		x2 = 4 x3	Wet weight in grams	
		-		
		x4	Current weight or dry weight in grams	
		x5	Actual result in requested display mode	
		X6	Drying time (seconds)	
	HA26山L	Command understood, parameter wrong		
Example 1				
Command	HA26山3	Inquiry of drying data		
Response	HA26LAL2L3L4.762L3.066L35.61L497			
		Drying has been ended regularly, result requested in % moisture content, wet weight 4.762 g, dry weight 3.066 g, 35.61 % moisture content, drying ended		
		at 497 s		
		u +07 3		

Example 2

Command	HA26山2	Inquiry of drying data
Response	HA26UAU1U2	L2.672L2.467L92.33L143
		Drying is running, result requested in % dry content, wet weight 2.672 g, dry weight 2.467 g, 92.33 % dry content, drying for 143 seconds in progress

Comment

- If a drying is inexistent (e.g. after a RAM LOST), the parameters x3..x6 are set to 0.
- If the measuring results exceed the tolerances for ATRO result display (L-999.99 % AM or >999.99 % AD) the selected results in $x^2 = 4$ AM or $x^2 = 5$ AD will automatically be transferred in $x^2 = 3$ MC or $x^2 = 2$ DC respectively.

HA27 In	quiry of dryin	g result		
Command	HA27ux1	Inquiry o	Inquiry of drying data in configurable display mode	
		x1 = 0	currently set display mode	
		x] =]	Grams	
		x1 = 2	DC (dry content)	
		x1 = 3	MC (moisture content), (factory setting)	
		x1 = 4	AM (ATRO moisture content)	
Response	HA27⊔A⊔x1ı	⊥ x2		
·		x1 x2	Drying status (always 7 digit number) Display mode (g, %DC, %MC, %AM, %AD)	
	HA27山I	Respons	e not available (drying in progress)	
Example				
Command	HA27山3	Inquiry o	of drying result	
Response	HA27uAuu-73.25%MC		-	
		Drying re	esult -73.25 % MC	

Comment

If the measuring results exceed the tolerances for ATRO result display (L-999.99 % AM or > 999.99 % AD) the selected results in $x_1 = 4$ AM or $x_1 = 5$ AD will automatically be transferred in $x_1 = 3$ MC or $x_1 = 2$ DC respectively.

HA40 Inquiry / setting of language

Inquiry of lan Command Response	nguage HA40 HA40uAux	Inquiry of language currently set x Language (see below)
Setting lange	Jage	
Command	HA40 LX	Set language
		x = 1English $x = 2$ German $x = 3$ French $x = 4$ Italian $x = 5$ Spanish $x = 6$ Russian $x = 7$ Japanese (Nihongo) $x = 8$ Portuguese
Response	HA40山A	Language set
	HA40山L	Command understood, parameter wrong

HA40X Inquiry / setting of menu parameters

All commands on the inquiry / setting of menu parameters work similiarly. As an example the inquiry for the startmode setting and the actual setting of the startmode is shown.

Inquiry of menu parameters

-	HA401 HA401uAux1	Inquiry of startmode		
			Startmode automatic (factory setting) Startmode manual	
Example				
Command	HA401	Inquiry of s	startmode	
Response	HA401_A_1	Startmode manual		
Setting menu	Setting menu parameters			
Command	HA401ux1	Setting of startmode		
		x1 = 0	Startmode automatic (factory setting)	
		x] = 1	Startmode manual	
Response	HA401⊔A	Startmode	set	
	HA401uL	Command	l understood, parameter wrong	
Example				
Command	HA401u0	Setting startmode to automatic		
	HA401⊔A	Startmode	set	

List of menu parameter inquiries / settings


HA401	Startmode			
	x1 = 0	Startmode automatic (factory setting)		
	x] =]	Startmode manual		
HA402	Protection against change in the settings (menu protection)			
	x1 = 0	Menu protection off, changes possible (factory setting)		
	x1 = 1	Menu protection on, changes are not possible		
HA403	Printer			
	x1 = 0	Printout off		
	x1 = 1	Printout on (factory setting)		

Comment

After the usage of HA08 the command HA403 has no effect on the printout but only on the setting of printout in the menu

4 System Configuration (MJ33 – Computer)

The MJ33 Moisture Analyzer is equipped with a 9 pin female RS232C Interface connector. It can be connected to a computer using a cable with order number 11101051 (9 pin) or order number 11101052 (25 pin).



The standard configuration of MJ33 is: 2400 baud, 7 bits, even parity, no handshake. These settings may be adjusted in the menu of MJ33. Additionally the printout setting in the menu should be "off". This avoids that printout strings are sent to the computer. Please refer to the Operating instruction.

5 What if...?

Tips from actual practice when the communication between the system (e.g. computer) and the Moisture Analyzer does not function.

Establishing the communication

Test whether the unidirectional operation is working:

Connection

For bidirectional communication, at least three connecting lines are needed:

- Data line from the Moisture Analyzer (TxD signal with RS232 interface).
- Data line to the Moisture Analyzer (RxD signal with RS232 interface).
- Signal ground line (SG with RS232 interface).

Make sure that all these connections are in order. Check the connector pin assignment of the connection cables.

Interface parameters

For the transmission to function properly, the settings of the following parameters must match at both the computer and the Moisture Analyzer:

- Baud rate (send receive rate)
- Number of data bits
- Parity bit

Check the settings at both devices.

Handshake

For control of the transmission, in part separate connection lines are used (CTS/DTR). If these lines are missing or wrongly connected, the computer or Moisture Analyzer can not send or receive data.

Check whether the Moisture Analyzer is prevented from transmitting by handshake lines (CTS or DTR).

Set the parameter "handshake" for the Moisture Analyzer and the peripheral device to "No Handshake" or "none". The handshake lines now have no influence on the communication.

To protect your METTLER TOLEDO product's future: METTLER TOLEDO service assures the quality, measuring accuracy and preservation of value of all METTLER TOLEDO products for years to come.

Please send for full details about our attractive terms of service.

Thank you.



Subject to technical changes.

© Mettler-Toledo AG 2007 11781025 Printed in Switzerland 0703/2.12

Mettler-Toledo AG, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland Phone +41-1-944 22 11, Fax +41-1-944 30 60, Internet: http://www.mt.com