
Operating instructions

METTLER TOLEDO PB-L line of balances







Operating instructions in a nutshell



Press key briefly


Press and hold key down until the desired display appears



automatic sequence

Switching on





Switching off





Simple weighing



	0.00 g
o	1 182.03 g
	1250.00 g

Adjusting (calibration) external










Taring











MMV

0.00 g










95.97 g

Piece counting*
























Dynamic weighing* Automatic start (Dyn A)





Manual start (Dyn M)





Unit switching*



22.00 g	
0.78 oz	
0.10 %	_

20

147,25 g

 $0.00 \, g$

12.07 g

1000 %

10 15 %

244 PCS

^{*}These functions must be activated in the menu (section 4.3.2)


Contents


1	Getting to know PB-L balances	4
1.1	Introducing the PB-L line of balances	4
1.2	Layout of PB-L balances	5
1.3	Overview of key functions	6
2	Startup	7
2.1	Unpacking / standard equipment	7
2.2	Cautionary notes	
2.3	Setting up, leveling, preparations for weighing below the balance, connecting to power supply	8
2.4	Adjusting (calibration)	10
3	Weighing	11
3.1	On/off switching	11
3.2	Simple weighing	11
3.3	Taring	12
3.4	METTLER TOLEDO DeltaRange balances	13
4	Menu	14
4.1	Overview	14
4.2	Menu operation	15
4.3	Description of menu options	16

Functions	22
Piece counting	22
Dynamic weighing	24
Switching weight units	26
Technical data, optional equipment	27
Technical data	27
Interface	28
MT-SICS Interface commands and functions	29
Optional equipment	3
Dimensional drawings	32
Appendix	33
Typical printouts from METTLER TOLEDO RS-P26	33
	and LC-P45 printers

Getting to know PB-L balances

1.1 Introducing the PB-L line of balances






Balance features

- The PB-L balance line ranges from models with a readability of 0.001g through to models with a readability of 1g. The weighing ranges extend from 151g to 8.1kg.
- In addition to basic weighing operations such as weighing, taring and adjusting (calibration) the
 functions "Piece counting", "Percent weighing" or "Dynamic weighing" (automatic or manual start) can
 be activated
- METTLER TOLEDO DeltaRange balances also have a movable fine range, with 10 times smaller display increments, over the the entire weighing range.
- Several PB-L balances are fitted with a glass draft shield in the factory; with other models a **draft shield** is available as an **optional extra**.
- All balances in the PB-L range are fitted with an **RS232C interface** as standard.

Note

All PB-L balances are available as certified versions. Please ask your METTLER TOLEDO dealer for details.

1.2 Layout of PB-L balances



- 1 Keys
- 2 Display with backlight
- 3 Model plate with the following data:

"Max": maximum capacity


"d": readability

"Min": minimum capacity (recommended minimum load; only relevant for certified balances)

"e": verification scale interval (smallest display increment tested during certification; only relevant for certified balances)

- 4 Draft shield element
- **5** Weighing pan
- **6** Draft shield (standard supply with models PBxx3-L)
- 7 Leveling feet
- **8** Hanger opening for weighing below the balance (underside of balance)
- 9 AC adapter socket
- 10 RS232C interface
- 11 Lug for optional antitheft device
- 12 Leveling control

1.3 Overview of key functions




The balances have two operator control levels: the **weighing mode** and the **menu**. The function of each individual key depends on the operator control level and how long the key is pressed.

Key funct	ions in weighing mode			
Press briefly		Press and hold down		
On →0/T← C	Switch onZero/tareCancel function	Off	Switch off balance (standby mode)	
5	SwitchChange settings	F	 Call function A function must be activated in the menu, otherwise "F nonE" appears in the display 	
⊖	Transfer weighing data via interfaceConfirm settings	Cal/Menu	 Adjust (calibrate) Show menu (hold key down until MENU appears) 	

Key functions in menu mode				
Press	briefly 🛅	Press and hold down		
С	 Close menu (without saving changes) 			
5	Change settings			
ightharpoons	Select menu options	Cal/Menu • Save changes and close menu		

2 Startup

2.1 Unpacking / standard equipment



All PB-L balances are supplied in environmentally compatible packaging. The standard equipment for every balance comprises:


- AC adapter, to national codes
- Weighing pan with pan support or draft shield element
- Transparent plastic protective cover, fitted, to keep your balance clean at all times
- Retaining ring for protective cover, fitted, only on balances without draft shield
- **Draft shield**, standard supply for models PBxx3-L (for other models a draft shield is available as an optional extra) (see Section 6.4 Optional equipment)
- Operating instructions, to ensure optimum utilization of your balance's capabilities
- **CE declaration of conformity** (in separate brochure 11780294)

2.2 Cautionary notes



- PB-L balances must **not** be operated in **hazardous areas** with the standard-supply AC adapter.
- Before connecting the AC adapter, verify that the voltage printed on it corresponds to the local mains voltage.
 If this is not the case, please contact your local METTLER TOLEDO dealer.
- PB-L balances may only be used indoors in a dry environment.
- Use only with a tested AC adapter with limited and SELV output current.

2.3 Setting up, leveling, preparations for weighing below the balance, connecting to power supply



The optimum location

The correct location makes an important contribution to the accuracy of the weighing results of high-resolution analytical and precision balances.


Hence, ensure a

• stable, vibration-free position as horizontal as possible

Avoid

- direct sunlight
- excessive temperature fluctuations
- drafts


The best location is on a stable bench in a corner protected against drafts, as far away as possible from doors, windows, radiators or the louvers of air conditioners.




Antitheft device

PB-L balances are equipped with a lug for the optional antitheft device.

The antitheft device (cable with lock) is suitable for all models. It is available from METTLER TOLEDO unter the order number 590101.





Levelina

PB-L balances have a spirit level and two adjustable leveling feet to compensate for slight irregularities in the surface of the weighing bench. The balance is exactly horizontal when the air bubble is in the middle of the level glass.

Procedure

Adjust the two leveling feet appropriately until the air bubble comes to rest exactly in the middle of the glass:

Air bubble at "12 o'clock" turn both feet counterclockwise

Air bubble at "3 o'clock" turn left foot clockwise, right foot counterclockwise

Air bubble at "6 o'clock" turn both feet clockwise

Air bubble at "9 o'clock" turn left foot counterclockwise, right foot clockwise

Note

The balance should be leveled and adjusted (Section 2.4) each time it is moved to a new location.

Preparations for weighing below the balance

To carry out weighing operations below the balance, the special cover on the underside of the balance must be slackened (care: when doing this the balance may only be laid on its side, not turned upside down!), turned through 180° and retightened. This exposes the opening for the hanger, making weighing below the balance possible.


Connecting the power supply

- → Before connecting the AC adapter, check that the voltage printed on it corresponds to the local mains voltage. If this is not the case, please contact your local METTLER TOLEDO dealer.
- → Plug the AC adapter into the AC adapter socket on the balance, and connect to the power supply.
- → The balance performs a self-test. This test is finished when "OFF" appears.
- → Press the **«On»** key briefly: the balance is in operational readiness. Before any work is performed with the balance, it must be adjusted (Section 2.4).

Note

An optional AccuPac B-S (rechargeable external battery) can be used to operate all PB-L balances independently of the mains power supply

2.4 Adjusting (calibration)



To obtain accurate weighing results, the balance must be adjusted to match the gravitational acceleration at its location.

Adjusting is necessary

- before the balance is used for the first time
- at regular intervals during weighing service
- after a change of location

Procedure

To obtain accurate results, the balance must be connected to the power supply for 30 minutes in order to reach operating temperature before adjusting.


- → Have required adjustment weight ready
- → Unload weighing pan
- → Press and hold the «Cal/Menu» key down until "CAL" appears in the display, then release key. The required adjustment weight value flashes in the display.
- → Place adjustment weight in center of pan. The balance adjusts itself automatically.
- → When "0.00 g" flashes, remove adjustment weight. The adjusting is finished when the message "CAL done" appears briefly in the display, followed by "0.00 g". The balance is again in weighing mode and ready for operation.

Notes

- Certfied PB-L models cannot be adjusted by the user, depending on weights and measures legislation.
- The adjustment procedure can be terminated at any time with the «C» key. The message "Abort" appears briefly to confirm that adjustment has been canceled, and the balance reverts to weighing mode.

3 Weighing

3.1 On/off switching



Switching on

→ Remove any load from weighing pan and press «**On**» key briefly.


The balance performs a display test (all segments in the display light up briefly).

When zero is displayed, the balance is ready for operation.

Switching off


→ Press and hold the «**Off**» key down until "OFF" appears in the display. Release the key.

3.2 Simple weighing



- → Place weighing sample on the weighing pan.
- → Wait until the stability detector "o" disappears.
- → Read the result.

3.3 Taring



- → Place empty container on the balance.
- → The weight is displayed.
- \rightarrow Tare: press the « \rightarrow 0/T \leftarrow » key briefly.
- → Add weighing sample to container. The net weight is now displayed.


If the container is removed from the balance, the tare weight will be shown as a negative value.

The tare weight remains stored until the $\leftarrow 0.71$ key is pressed again or the balance is switched off.

Note

With METTLER TOLEDO DeltaRange balances (next Section), the fine range with its 10 times smaller display increments is available again after every taring operation.

3.4 METTLER TOLEDO DeltaRange balances



METTLER TOLEDO **DeltaRange balances** have a movable fine range with 10 times smaller display increments over their entire weighing range. In this fine range an additional decimal place always appears in the display.

The balance operates in the fine range

- after switching on
- after every taring operation

If the fine range is exceeded, the balance display automatically switches to coarser display increments.


4 Menu

4.1 Overview

In the menu you can change the weighing unit (for certified balances, only if the national weights and measures legislation allows), select additional functions and carry out various settings. A description of the individual menu options is given in Section 4.3.

Notes


- With certified balances, this menu option has a fixed setting and cannot be changed.
- 2) With certified balances, only those weighing units allowed by the appropriate national weights and measures legislation may be selected.
- 3) This menu option is only shown if "Host" has been selected in menu option 7 (Peripheral unit).
- 4) These menu options are only shown if "Host" or "Printer" has been selected in menu option 7 (Peripheral unit).



Menu option

Factory setting

4.2 Menu operation



Opening the menu

In weighing mode, press and hold down the **«Cal/Menu»** key until "MENU" appears in the display. Release the key: the 1st menu option is displayed.



Select menu options

The « > » key is used to select individual menu options with their current settings one after the other.



Change settings

The «Sa» key is used to change the setting at the selected menu option. Every time the key is pressed, the next setting is displayed. Once the desired setting appears in the display, the next menu option can be selected (see above) or you can close the menu (see following sections).



Saving settings and closing the menu

Hold the **«Cal/Menu»** key down until "StorEd" appears in the display. Release the key and the balance reverts to weighing mode. All changes are saved.

c


Abort


Press the «C» key briefly. The balance reverts to weighing mode. Changes are not saved.

Note

If no entry is made within 45 seconds, the balance reverts to weighing mode. Changes are not saved.

4.3 Description of menu options





4.3.1 Reset or recording of balance settings (1st menu option "RESET")


Reset balance settings

→ Select "Reset", press and hold down the «Cal/Menu» key until the message "r donE" confirms that all menu settings have been reset. The balance then reverts to weighing mode and works with the factory settings (Section 4.1).

Recording balance settings

→ Select "List" and hold down the «**Cal/Menu**» key until the message "StorEd" is displayed.

The current balance settings are transmitted to the peripheral device connected to the interface. To do this the setting "Printer" must always be selected at the 7th menu option (Peripheral unit). The current balance settings are saved at the same time.




4.3.2 Functions (2nd menu option / Operating Section 5)

In addition to simple weighing, the following functions can be selected:

F count Piece counting
F 100 % Percent weighing

F dYn A Dynamic weighing with automatic start F dYn M Dynamic weighing with manual start


F nonE No function, simple weighing (factory setting)



4.3.3 Weighing mode (3rd menu option)

This setting allows you to adapt the balance to the weighing mode. Select "Std" (standard) for all normal weighing processes or "doS" (dispensing) for dispensing liquid or powdery weighing samples. In this setting, the balance responds very quickly to smallest weight changes.

When set to "robuSt" (absolute weighing), the balance responds only to greater weight changes, and the weighing result is very stable.




4.3.4 Weighing unit 1 (4th menu option "UNIT 1")


Depending on requirements, the balance can operate with the following units (possible with certified balances only if permitted by national weights and measures legislation):

Unit g	gram	Conversion factor	Comments factory setting
kg	kilogram	1 kg = 1000 g	not with 0.1 mg and 1 mg balances
mg	milligram	1 mg = 0.001 g	with 0.1 mg and 1 mg balances
ct	carat	1 cf = 0.2 g	0 0
lb	pound	1 lb = 453.59237 g	
OZ	ounce	1 oz = 28.349523125 g	
ozt	troy ounce	1 ozt = 31.1034768 g	
GN	grain	1 GN = 0.06479891 g	
dwt	pennyweight	1 dwt ≈ 1.555173843 g	
mo	momme	1 mom ≈ 3.749999953 g	
m	Mesghal	1 msg ≈ 4.6083162 g	
H fl	Hong Kong taels	1 tlh ≈ 37.42900 g	
S tl	Singapore taels	1 tls ≈ 37.799366256 g	The Malaysian tael has the same value
† †I	Taiwan taels	1 tlt ≈ 37.499995313 g	•
cl	tical	1 tical ≈ 16.3293 g	

4.3.5 Weighing unit 2 (5th menu option "UNIT 2")

If it is required to show the weighing results in weighing mode in an additional unit by pressing the «S» key, the desired second weighing unit can be selected in this menu option. The same weighing units are available as under "UNIT 1", with the exception of the tael units ("H tl", "S tl" und "t tl").





This menu option allows you to switch the automatic zero correction on or off. When it is switched on, the zero point is automatically corrected for drift or contamination of the weighing pan.

The following settings are available:

Autozero switched on

The zero point is automatically corrected.

Autozero switched off

The zero point is not automatically corrected. This setting is advantageous for special applications (e.g. evaporation measurements).

Note

With certified balances, this setting is possible only with a resolution of e = 10d.



4.3.7 Peripheral unit (7th menu option)

At this menu option you can select the peripheral device connected to the optional RS232C interface. The balance automatically saves the appropriate settings (Section 4.3.8 - 4.3.12) for every peripheral device.

Printer Connected to a printer.

Host Connection to any desired peripheral device.

Aux. display Connection of an optional auxiliary display unit (communications parameters cannot be selec-

ted).



4.3.8 Send format (8th menu option)

Note: This menu option is only available if you have selected the setting "Host" in the 7th menu option (Peripheral unit)!

At this menu option you specify the data transfer format.

"SICS": The MT-SICS data transfer formats are used. Please refer to the "MT-SICS Reference Manual Basic-S balances 11780447" available from your METTLER TOLEDO dealer or download from the

Internet (www.mt.com/sics-classic). More Information please find in the Section 6.3.

"PM"*: The following PM balance data transfer formats are used:

S. Stb: பபபபப1.67890ப

S. Cont: Suuuu1.67890ug

SDLLL1.39110Lg

Additional information can be found on the Internet at www.mt.com/classic.

* unidirectional, no MT-SICS commands are accepted.



4.3.9 Send mode (9th menu option)


Note: This menu option is only available if you have selected the setting "Host" in the 7th menu option (Peripheral unit)!


At this menu option you specify how data is transmitted to a peripheral device.

S. oFF Send mode switched off

S. Stb The next possible stable value will be transferred after the «□→» key has been triggered.

S. Cont All values are transferred automatically.






The baud rate (data transfer rate) determines the speed of transmission via the serial interface. The unit is the baud (bd) = 1 bit/second.

The following settings are available: 600 bd, 1200 bd, 2400 bd, 4800 bd, 9600 bd and 19200 bd.

For problem-free data transmission the sending and receiving devices must be set at the same value.




4.3.11 Bit/Parity (11th menu option)

Note: This menu option is only available if you have selected the setting "Host" in the 7th menu option (Peripheral unit)!

At this menu option you can set the character format for the attached peripheral device.

7b–E 7 data bits/even parity 7h-no 7 data bits/no parity 8 data bits/no parity 8b-no 7 data bits/odd parity 7b-odd



4.3.12 Handshake (12th menu option)

Note: This menu option is only available if you have selected the setting "Host" in the 7th menu option (Peripheral unit)!


This menu option allows you to match the data transmission to different serial receivers.

HS off No handshake

HS SoFt Software handshake (XON/XOFF) HS HArd Hardware handshake (RTS/CTS)

5 Functions

5.1 Piece counting



Requirement

The function "F count" must be activated in the menu (Section 4).

 \rightarrow Place empty container on the balance and tare by briefly pressing the « \rightarrow 0/T \leftarrow » key.

Setting the reference: a reference weight must first be entered for piece counting:

→ Add a number of reference pieces to container. Possible numbers* are 5, 10, 20, 50, 100 and "no" (this setting deactivates the piece counting function). *approved balances min 10
Caution: Take into account minimum values: min. reference weight = 10d (10 digits), min. piece

 $weight^* = 1d (1 digit)! \ \ \, \hbox{\bf *approved balances} \ \ \, min \ \, 3e$


- **Note:** 1 digit corresponds to 1 display increment.
- → Hold the «F» key down until "SEt ... PCS" is displayed.
- → Repeatedly press the «与» key until the display equals the number of reference pieces entered.
- → Confirm the number of reference pieces with the «□→» key (or wait 7 seconds, in which case the number is adopted automatically). The current number of pieces (PCS = pieces) is displayed.

Note: The current reference weight remains stored until the reference setting is changed or the power supply is interrupted.

Switching between piece count and weight display

- → Add weighing sample to the container and read off number of pieces.
- → Press the «S» key. The weight is displayed.
- → Return to the piece count display by pressing the «S» key again.

5.2 Percent weighing




Requirement

The function "F 100 %" must be activated in the menu (Section 4).

Set target weight

- → Target weight (Reference weight, which corresponds to 100 %) in center of pan.
- → Hold the «F» key down until "SEt 100 %" is displayed.
- → Press the «S» key to select "SEt 100 %" or "SEt no %" (Percent weighing deactivated).
- → The «□→» key can be used briefly to confirm or automatic acceptance after 7 seconds. The target weight is specified.

Note: The current target weight remains stored until a new target weight is set or the power supply is interrupted.



Percent weighing / switching

- Place weighing sample in center of pan.
 The weight of the sample is displayed as a percentage of the target weight.
- → Press the «Si» key, The weight is displayed. (Unit 1 and and Unit 2, if activated).
- → Return to display in percent: pressing the «与» key again.

5.3 Dynamic weighing

Dynamic weighing is suitable for the weighing of unstable weighing samples. The mean value of the weighing results is determined over a specified time period (weighing time). The more unstable the weighing sample, the longer the selected weighing time.

Requirement

"F dYn A" for automatic start or "F dYn M" for manual start must be activated in the menu (Section 4). Factory setting is a weighing time of 3 seconds (t = 3").

Tare container


→ Tare: Press the $\ll \rightarrow 0/T \leftarrow \gg \text{ key}$.


Dynamic weighing with automatic start (F dYn A)

- ightarrow The «S» key can be used select the dynamic weighing. The display shows the symbol ightarrow.
- → Load weighing sample. As soon as the balance is relatively stable, weighing starts automatically. During the weighing time, a "count down" runs in the display.

→ Read off result.

The result of the dynamic weighing is displayed with * (= calculated value) and remains in the display until the weighing sample is removed from the weighing pan or the container.





Dynamic weighing with manual start (F dYn M)

- → The «S» key can be used select the dynamic weighing. The display shows the symbol ...
- → Load weighing sample.
- → Start weighing with the «□→» key. During the weighing time, a "count down" runs in the display.
- → Read off result.

The result of the dynamic weighing is displayed with * (= calculated value) and remains in the display until the weighing sample is removed from the weighing pan or the container.

Notes


- The weighing cycle with the same weighing sample can be restarted with the «□→» key.
- The «Sa» key can be used to switch between dynamic weighing and normal weighing.
- For weighing goods below 5 g the weighing must be started manually with the «□→» key, even for dynamic weighing with automatic start.

Changing the weighing time


- → Press and hold the «**F**» key, until "t = 3"" appears in the display.
- → Repeatedly press the «S» key, until the desired weighing time appears. Possible values are 3", 5", 10", 20", 1", 2".
- → Selection with the «□→» key briefly to confirm or by automatic acceptance after 3 seconds.

Note

The set weighing time remains stored until it is reset or the power supply fails.



5.4 Switching weight units



Requirement

Different weight units must be activated in the menu for unit 1 and unit 2 (Section 4).

→ The «与» key can be used at any time to toggle between the two weighing units selected in the menu ("UNIT 1" and "UNIT 2").

Notes:

- Switching between weight units may be blocked with certified balances, depending on national weights and measures legislation.
- This function is not available with dynamic weighing.

6 Technical data, optional equipment

6.1 Technical data

Standard equipment

- Protective cover, transparent, made from Barex
- AC adapter to national codes 100–240 VAC/50–60 Hz, 0.3 A 12 VDC, 0.84 A

Balance power input: 8–14.5 VAC, 50/60Hz, 6 VA or 9.5–20 VDC, 6 W

- Built-in RS232C interface
- · Draft shield with PBxx3-L models
- All models can weigh below balance.
- Display with backlight

Materials

- Housing: die-cast aluminum, painted
- Weighing pan: Chromium-nickel steel. X2CrNiMo 17 13 2 (1.4404)

Protection

- Protected against dust and water
- Pollution degree: 2
- Overvoltage category: class II
- EMC: see declaration of conformity (separate brochure 11780294)

Ambient conditions

The technical data are valid unter the following ambient conditions:

- Ambient temperature
 Relative humidity
 10 °C ... 30 °C
 15 % ... 80 % (
 - 15 % ... 80 % at 31 °C, linear decreasing to 50% at 40 °C noncondensing

Operability is assured at ambient temperatures between 5 and 40 $^{\circ}\text{C}.$


Technical data	PB153-L	PB303-L	PB303-L DeltaRange	PB602-L	PB1502-L	PB3002-L	PB3002-L DeltaRange
Readability	0.001 g	0.001 g	0.001 g*/0.01 g	0.01 g	0.01 g	0.01 g	0.01 g*/0.1 g
Max. capacity	151 g	310 g	60 g*/310 g	610 g	1510 g	3100 g	600 g*/3100 g
Repeatability (sd)	0.001 g	0.001 g	0.001 g*/0.008 g	0.01 g	0.01 g	0.01 g	0.01 g*/0.08 g
Linearity	0.002 g	0.002 g	0.01 g	0.02 g	0.02 g	0.02 g	0.1 g
Sensitivity temperature drift (10 °C 30 °C)	6 ppm/ °C	6 ppm/ °C	6 ppm/ °C	6 ppm/ °C	6 ppm/ °C	6 ppm/ °C	6 ppm/ °C
Settling time, typical	2 s	2 s	2 s	2 s	2 s	2 s	2 s
Adjustment weight (optional)	external 100 g ¹⁾	external 200 g ¹⁾	external 200 g ¹⁾	external 500 g 1)	external 1000 g ¹⁾	external 2000 g ¹⁾	external 2000 g ¹⁾
Backlight	yes	yes	yes	yes	yes	yes	yes
External dimensions of balance (W/D/H)	245x321x236	3 mm		245x321x89	mm		
External dimensions of packaging (W/D/H)	381x436x495	381x436x495 mm (0.082 m³)		381x436x27	3 mm (0.045 m³))	
Weighing pan	ø 100 mm	ø 100 mm		ø 180 mm			
Usable heigh of draft shield	141 mm	141 mm		-			
Net weight (with packaging)	4.9 kg (6.5 kg)			3.6 kg (4.7 k	(g)		

^{*} Fine range (DeltaRange) 1) Optional equipment

Technical data	PB1501-L	PB3001-L	PB5001-L	PB8001-L	PB8000-L
Readability	0.1 g	0.1 g	0.1 g	0.1 g	1 g
Max. capacity	1510 g	3100 g	5100 g	8100 g	8100 g
Repeatability (sd)	0.08 g	0.08 g	0.08 g	0.08 g	0.8 g
Linearity	0.1 g	0.1 g	0.1 g	0.1 g	1 g
Sensitivity temperature drift (10 °C 30 °C)	10 ppm/ °C	10 ppm/°C	10 ppm/ °C	10 ppm/ °C	10 ppm/°C
Typical stabilization time	1.5 s	1.5 s	1.5 s	2 s	1 s
Adjustment weight (optional)	external 1000 g ¹⁾	external 2000 g ¹⁾	external 2000 g 1)	external 4000 g ¹⁾	external 4000 g ¹⁾
Backlight	yes	yes	yes	yes	yes
External dimensions of balance (W/D/H)	245x321x89 mm				
External dimensions of packaging (W/D/H)	381x436x273 mm (0.045 m³)				
Weighing pan	ø 180 mm				
Net weight (with packaging)	3.6 kg (4.7 kg)				

¹⁾ Optional equipment

6.2 Interface



RS232C interface and interface accessories

Every PB-L balance is fitted with an RS232C interface for attachment to a peripheral device (e.g. printer or PC with a 9-pin male connector). Matching to a different device can be carried out in the menu (Sections 4.3.7-4.3.12).

You will find a detailed description of the available interface commands in the brochure "Reference Manual MT-SICS Basic-S balances 11780447" available from your METTLER TOLEDO dealer or download from the Internet (www.mt.com/sics-classic).

The wide range of features of the PB-L balances regarding documentation of the results can not be exploited to the full until a printer, e.g. the RS-P26 or LC-P45 from METTLER TOLEDO is attached. The printed results make a decisive contribution to a simple way of working in compliance with GLP/GMP.

6.3 MT-SICS Interface commands and functions

Many of the balances and scales used have to be capable of integration in a complex computer or data acauisition system.

To enable you to integrate balances in your system in a simple manner and utilize their capabilities to the full, most balance functions are also available as appropriate commands via the data interface.

All new METTLER TOLEDO balances launched on the market support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS). The commands available depend on the functionality of the balance.

Basic information on data interchange with the balance

The balance receives commands from the system and acknowledges the command with an appropriate response.

Command formats

Commands sent to the balance comprise one or more characters of the ASCII character set. Here, the following must be noted:

- Enter commands only in uppercase.
- The possible parameters of the command must be separated from one another and from the command name by a space (ASCII 32 dec., in this description represented as

).
- The possible input for "text" is a sequence of characters of the 8-bit ASCII character set from 32 dec to 255 dec
- Each command must be closed by C_pL_e (ASCII 13 dec., 10 dec.).

The characters $C_{\rm p}L_{\rm p}$, which can be inputted using the Enter or Return key of most entry keypads, are not listed in this description, but it is essential they be included for communication with the balance.

Example

S - Send stable weight value

Command	s	Send the current stable net weight value.				
Response	SuSuWeightVa	S⊔S⊔WeightValue⊔Unit				
		Current stable weight value in unit actually set under unit 1.				
	S⊔I	Command not executable (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).				
	S⊔+	Balance in overload range.				
	Su-	Balance in underload range.				
Example						
Command	S	Send a stable weight value.				
Response	SUSUUUUU10	0.00⊔g				
		The current, stable weight value is 100.00 g.				

29

The MT-SICS commands listed below is a selected list of available commands. For additional commands and further information please refer to the Reference Manual "MT-SICS Basic-S balances 11780447" downloadable from the Internet under www.mt.com/sics-classic.

S - Send stable weight value

Command s Send the current stable net weight value.

SI - Send value immediately

Command SI Send the current net weight value, irrespective

of balance stability.

SIR - Send weight value immediately and repeat

Command SIR Send the net weight values repeatedly,

irrespective of balance stability.

Z – Zero

Command **z** Zero the balance.

@ - Reset

Command @ Resets the balance to the condition found after

switching on, but without a zero setting being

performed.

SR - Send weight value on weight change (Send and Repeat)

Command SR Send the current stable weight value and then

send continuously the stable weight value after

every weight change.

The weight change must be at least 12.5 % of the last stable weight value, minimum = 30d.

ST – Send stable weight after pressing \Longrightarrow (transfer) key

Command ST Inquiry of actual status of the ST function.

SU - Send stable weight value with currently displayed unit

Command sv As the "s" command, but with the currently


displayed unit.


6.4 Optional equipment

6.4 Optional equipment					
AC adapter		Auxiliary display (RS232)		Printer, Report printer (RS-P26)	
AC adapter universal (EU, USA, AU, UK) 100–240 VAC/50–60 Hz, 0.3 A 12 VDC, 0.84 A	11120270	Auxiliary display for connection to the RS232C interface	12120057	Plain-paper printer, 24 characters, with additional functions (time, date)	12120788
		Draft shields		Protective cover	
AccuPac B-S Rechargeable external power source for 19	8 houre	Glass zylinder draft shield (for 0.001 g balances)	11137305	1 piece	11103681
weighing operation with no mains	o nours	Draft shield with sliding doors		Software	
connection	21254691	"mg" (165 mm) *	11137468	LabX direct balance	11120340
		Draft shield "mg" (141 mm) *	11103683	(software for easy transfer to PC)	
Adjustment weights		* 5			
Available as OIML weights (E1, E2, F1, with calibration certificate) For further details see METTLER TOLEDO Weights brochure	11795461	* Balances with a weighing pan Ø 180 mm require the weighing pan Ø 175 mm	11103680	Transport case For all PB-L models, with room for balance, draff shield	
or see www.mt.com/weights	11793401	Interface cable		165 mm / 141 mm and weight	11101050
		• RS9-RS25: (m/f), length 2 m	11101052		
Antitheft device		• RS9–RS9: (m/f), length 1 m	11101051		
Cable with lock (for all models)	00590101	RS9-RS9: (m/m), length 1 m RS232-USB converter cable	21250066 11103691		
Auxiliary display (RS/LC-BLD) Auxiliary display including RS cable for connection to the RS232C interface and separate AC adapter	00224200	Printer, Application printer (LC-P45) Plain-paper printer, 24 characters, with additional functions (time, date, statistic,			
		multiplier etc.)	00229119		

6.5 Dimensional drawings

All dimensions in millimeters (mm)





1) Opening for hanger

7 Appendix

7.1 Typical printouts from METTLER TOLEDO RS-P26 and LC-P45 printers

Function: Adjusting

, ,
- BALANCE CALIBRATION - 12.02.2007 09:52:12
METTLER TOLEDO
Type: PB3002-L
SNR: 1118015657
SW: 1.20
Weight ID: Weight: 2000.00 g
External Cal. done
Signature:
END

Function: **Piece counting**Printout with reference weight

PIECE APW:	COUNTING 0.99460	
Out of:	10	PCS
	27.000 27	g PCS

Function: Percent weighing

Ref.	olo	-	WEIGHING 10.008 g 100.00 %
			60.01 g 599.59 %

Function: Dynamic weighing

	DYNAMIC gh Time:	WEIGHING 2 s
DW	49.999	g

Function: **List**Printout of the current balance settings

LIST OF S 12.02.2007	09:52:12
METTLER TOLED Type: SNR: SW:	PB3002-L 1118015657 1.20
Application: Dynamic A	
Weighing Para Weighing Mod Unit 1 Unit 2 A.Zero	
Peripheral De P.Device Baud Bit/Parity Handshake	Printer 2400 7b-even Off
P.Device Sendmode Baud Bit/Parity Handshake	Host Off 9600 8b-no Soft

Function: Verification of the calibration (adjustment) with external weight.

Only possible with LC-P45. Function is triggered via the printer.

BALANCE TEST 12.02.2007 09:52:12
METTLER TOLEDO Type: PB3002-L SNR: 1118015657 SW: 1.20
Weight ID:
Target :
External test done
Signature:
END

Function: **Statistics**Only possible with LC-P45. Function is triagered via the

Function is triggered printer.

12.02.2007 1

12.02.200	7 10:44:07
ID	666
SNR:	1118015657
1	1100.15 g
2	1600.10 g
3	1699.95 g
n	3
x	1466.733 g
s	321.372 g
srel	21.91 %
min.	1100.15 g
max.	1699.95 g
dif.	599.80 g
	END

Function: **Multiplier**Only possible with LC-P45. Function is triggered via the

printer.

12.02.2007	08:23:22
ID	242
SNR:	1118015657
Factor	1.65
	588.43 g
*	970.9095

Notes

The operating instructions for the LC-P45 include a description of the functions that are triggered via that printer.

The RS-P26 prints all reports in English. This applies also to the LC-P45 reports that originate in the balance. In the case of reports triggered by the LC-P45, the following languages may be selected: German, English, French, Spanish or Italian.

7.2 What if ...?

Error/Error message	Cause	Rectification
Г	Overload	→ Remove sample from weighing pan, zero again (tare).
L J	Underload	→ Check whether weighing pan is positioned properly.
Error 1	No stability in taring or adjusting (calibration) when reference weight for piece counting is placed on pan	 → Wait for stability before pressing key. → Ensure more stable ambient conditions. → Remove weighing pan and clean if necessary
Error 2	Wrong adjustment weight on pan or none at all	→ Place required adjustment weight in center of pan.
Error 3	Reference number for piece counting too small	→ Increase number of reference pieces.
Error 4	Internal fault	→ Contact METTLER TOLEDO customer service.
Error 8	No standard calibration	→ Contact METTLER TOLEDO customer service.
20000 g	Wrong weighing pan or pan missing	→ Mount correct weighing pan.
Rbort	Adjustment aborted with the «C» key	

7.3 Maintenance and cleaning

Service

Regular servicing of your balance by a service technician prolongs its working life. Ask your METTLER TOLEDO dealer for details of servicing options.

Cleaning

Every now and then, clean the weighing pan, draftshield element, draftshield (depending on the model) and housing of your balance using a damp cloth. Your balance is made of high-quality, durable materials and can therefore be cleaned with a standard, mild cleaning agent.



Please observe the following notes

- On no account use cleaning agents, which contain solvents or abrasive ingredients, as this can result in damage to the terminal overlay.
- After working with chemicals, it is advisable to wash or clean the weighing pan and the bottom plate (if draft shield fitted).
- Although all materials are of high quality, corrosion may occur
 if corrosive substances are deposited on chrome steel for an
 extended period of time (and if air is excluded, for example by
 a coating of grease).
- Ensure that no liquid comes into contact with the balance or the AC adapter!
- Never open the balance or AC adapter they contain no components, which can be cleaned, repaired or replaced by the user.
- Soiled protective covers can be replaced on all balance types (see Optional equipment).

Disposal



In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste.

This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

To protect your METTLER TOLEDO product's future:
METTLER TOLEDO Service assures the quality, measuring accuracy and preservation
of value of all METTLER TOLEDO products for years to come.
Please send for full details about our attractive terms of service.
Thank you.



Subject to technical changes and to the availability of the accessories supplied with the instruments.

© Mettler-Toledo AG 2007 11780755B Printed in Switzerland 0710/2.12