Contents:

4 Membrane bodies
2 x 5 ml O₂ electrolyte
1 O-ring set
1 Watering cap

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www.mt.com/com

Changing the electrolyte of the InPro® 6950/6900 O₂ sensors

ENGLISH

 **Note:** The InPro 6950/6900 uses a special electrolyte which contains an oxygen scavenger. This electrolyte ensures fast response time and together with the guard ring increases the signal stability of the sensor. The electrolyte must be exchanged on a regular basis or when the sensor has been exposed to the air for more than 24 hours without a watering cap or if the sensor has been stored for several months.

If the membrane exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body has to be replaced.

 **Warning!** The O₂ electrolyte has an alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.**


If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.

When changing the electrolyte, please observe the following instructions (see also the following illustration):

 **Attention!** Make sure that this maintenance step is carried out in **clean place**.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.
4. Examine the O-rings visually for mechanical defects, and replace if necessary.
5. Half-fill the new membrane body with O₂ electrolyte.


 **Note:** The electrolyte bottle is equipped with a special pouring system. To ensure proper functioning, hold the bottle vertically, upside-down.

 **Note:** Make sure that all air bubbles are removed from the membrane body. Air bubbles can be removed by carefully tapping on the membrane body.

6. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue.
7. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
8. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated.

Equipment features 12 mm Sensor

T-6950 Type Kit with 4 membr., Order no. 52 206 106

T-6950 Gas Type Kit with 4 membr., Order no. 52 206 109

T-6900 Type Kit with 4 membr., Order No. 52 201 003

T-6900 R Type Kit with 4 reinf. membr., Order No. 52 201 109

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03/14 © Mettler-Toledo AG, Switzerland