Membrane Kit for InPro® 5000 CO₂ Sensors

Membran-Kit für InPro® 5000 CO₂ Sensoren

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A decorative graphic consisting of numerous thin, parallel green lines that form a large, stylized arrow pointing downwards and to the right, located in the bottom right corner of the slide.

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Instructions for electrolyte and membrane exchange of the InPro® 5000 CO₂ sensors


ENGLISH

 **Electrolyte exchange:**
After each calibration of the interior body (without membrane body) new electrolyte should be used.

 **Membrane body exchange:**
We recommend to use a new membrane body for each batch (sterilization). Prior to membrane body exchange the interior body (pH electrode) must be calibrated.

When replacing the membrane body and electrolyte, strictly observe the following instructions:

- 1.**
Unscrew the cap sleeve from the shaft and carefully pull it from the sensor.
- 2.**
Usually the membrane body remains in the cap sleeve and should be ejected by pushing it with the flat finger tip. Otherwise pull the membrane body from the interior body (pH electrode).
Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
- 3.**
Clean the interior body with deionized water and tip with a paper handkerchief or paper tissue.
- 4.**
Check all O-rings for mechanical defects and replace if necessary.
- 5.**
Half-fill the membrane body with CO₂ electrolyte and make sure that all bubbles are removed. Air bubbles can be removed by gently knocking on the membrane body.
- 6.**
Slip the membrane body over the interior body while holding the sensor in a horizontal position. The excess electrolyte will be displaced and must be removed with a paper tissue.
- 7.**
Carefully slip the cap sleeve over the fitted membrane body and screw it down.
- 8.**
Remove expelled electrolyte with a paper tissue.
- 9.**
Sensor is now ready for installation, sterilization and process calibration.

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08/03 © Mettler-Toledo (Schweiz) AG, Switzerland