EasyMax™ 102

Short Installation Instructions

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## Short Installation Instructions

### First Steps to Setting Up Your Instrument

**METTLER TOLEDO**

# Short Installation Instructions

Please follow these instructions carefully to ensure a secure and successful installation of your EasyMax™ 102 Instrument. Please check the EasyMax™ Product Catalog for more detailed information such as a fully illustrated part list for standard and optional equipment.

- Cool the instrument with salt-free tap water, ethanol, ethylene glycol, silicone oil, or a cryostat and purge it with dry air, nitrogen or argon before switching on the EasyMax™ 102.
- Check the minimum flow rates of purge gas and cooling liquid. If the given flow rates lie below these values, the instrument can sustain damage. See "Pre-Installation Guide" or the Operating Instructions.
- Attach both Pt100 temperature (Tr) sensors before you start the EasyMax™ 102! If you forget it, an error message is shown after the start procedure. Further information can be found in the Operating Instructions.


## Purge gas connection

Use the red PVC tubes for purge gas connection. Fit purge gas tubes to connectors with tube couplings included in the standard delivery. Both reactors can be purged simultaneously with gas.

## PVC tubes


Use uncolored PVC tubes to connect cooling liquid. Hose clamps will improve grip.

## Step 1 – Connect tubes, cables, and devices according to the scheme

### Rear side connections

- Cooling liquid inlet
- RS232 interface <sup>1)</sup>
- Power supply
- Network <sup>1)</sup>
- Touchpad
- Purge gas inlet
- Cooling liquid outlet

### Front side connections

- Purge gas outlet
- Overhead stirrer <sup>1)</sup>
- Pt100 temperature (Tr) sensor <sup>1)</sup>
- Thermostat
- USB port
- Power push button and status LED

<sup>1)</sup> Optional equipment

**Step 2 – Start up your instrument**


- Press the power push button on the EasyMax™ 102 – the start-up procedure shows this page on the touchpad followed by the **System Overview** page.



System Overview		Status: OK
4/17/2008 2:15 pm		
No experiment	No experiment	
Tr	56.8 °C	Tr
Tr-Tj	-6.4 K	Tr-Tj
Tj	63.2 °C	Tj
R	---- rpm	R

**Step 3 – Set the language, date, and time**

- Touch the **Settings** icon.
- Touch the **English** button.
- Touch “your” language button to close this page and transfer the input to the **Settings** page.

- Touch the **Date and Time** button.
- Touch the currently displayed date to change the date.
- Confirm with **OK** to transfer the input to the **Date and Time** page.
- Confirm with **OK** to go back to the **Settings** page.

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©06/2008 Mettler-Toledo AG
## **Mettler-Toledo AutoChem Inc.**
## **Mettler-Toledo AG, AutoChem**