Automatic feeding

Reliable item separation



LV11 automatic feeder – When individuality counts



Automatically controlled feeding

A photocell gate monitors the items passing through and interrupts the flow after each piece. If the weight value determined by the balance is valid, feeding is restarted. If required, at the end of a defined sample the automatic feeder empties itself.



Adjustable flow speed

The LV11 is started and controlled from the keyboard or via the interface. Depending on the size, quality or shape of the items, the vibration intensity, and therefore the flow speed, can be optimized at one of four levels. Easily visible LEDs indicate the selected level.

Constant automatic monitoring with the LV11 automatic feeder — item by item. Capsules, tablets or small electronic and mechanical components are fed to the balance quickly by the LV11 automatic feeder. Large weighing series are radically simplified. The LV11 is therefore especially suitable for statistical quality control (SQC) in the pharmaceutical industry.




Reliable data communication

Bidirectional RS232 or LocalCAN interfaces ensure troublefree data transfer to your computer or weighing system.



LV11











A computer or SQC system, a METTLER TOLEDO analytical or precision balance and an LV11 form a fully automatic control unit in:

- SQC compact systems (SQC14, SQC-AX, FreeWeigh.Net Compact)
- SQC network systems (FreeWeigh.Net)
- LabX balance
- Externally controlled systems (PC, SPS)

Technical data

230V (200...250V) or Power supply (switchable) 115V (100...125V) 50...60Hz Frequency Power consumption 15VA

Item size Ø 1...33mm Vibration levels

225 x 225 x 260mm (WxDxH) Dimensions

Weight 11.5kg RS232 interface parameters:

9600/8/No parity/Xon/Xoff (factory setting)

9600/8/No parity/DTR/DSR 2400/7/Even/Xon/Xoff 2400/7/Even/DTR/DSR

Accessories

LocalCan cable for connection of a balance with LocalCan interface is included as standard.

- RS232 cable for connection of a computer or SQC compact system 21900576
- RS232-Y cable for connection of an SQC14 and a balance 21900921
- Door with recess for XS/XP analytical balances 11106715
- Door with recess for XS/XP precision balances with high draft shield 11132711





