Operating Instructions

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
Mode d'emploi

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Istruzioni d'uso

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Education Line

pH meter EL20 / EL2

METTLER TOLEDO

日本語	Italiano	Español	Français	Deutsch	English
-----	----------	---------	----------	---------	---------

## Contents

<b>1. Introduction</b>	<b>3</b>
<b>2. Safety measures</b>	<b>3</b>
<b>3. Installation</b>	<b>4</b>
3.1    Installing the electrode arm (optional for EL20)	4
3.2    Installing the batteries (EL2)	4
3.3    FiveGo™ clip (EL2)	5
3.4    Installing the IP54 cover (EL2)	5
<b>4. Operating the Education Line pH meter</b>	<b>6</b>
4.1    Display and key controls	6
4.2    Calibration	8
4.2.1    Buffer groups	8
4.2.2    Selecting a predefined buffer group	8
4.2.3    Performing a 1-point calibration	8
4.2.4    Performing a 2-point calibration	8
4.3    Sample measurement	9
4.3.1    Performing a pH measurement	9
4.3.2    Performing a mV / Redox measurement	9
4.4    Temperature measurement	9
4.4.1    Manual temperature compensation	10
4.5    Analog output (EL20)	10
4.6    Using the memory (EL2)	10
4.6.1    Storing a reading	10
4.6.2    Recalling from memory	10
4.6.3    Clearing the memory	10
4.7    Self-diagnosis	11
4.8    Error message	11
<b>5. Maintenance</b>	<b>12</b>
5.1    Meter maintenance	12
5.2    Electrode maintenance	12
5.3    Disposal	13
<b>6. Accessories</b>	<b>14</b>
<b>7. Specifications</b>	<b>15</b>
<b>8. Appendix</b>	<b>16</b>
8.1    Buffer tables	16
8.2    Error Limits	17

## 1. Introduction

Thank you for purchasing this high quality METTLER TOLEDO meter. Ease of use, reliable results and ergonomic design – these are our goals for products.

The Education Line meters have an excellent price/performance ratio and many useful accessories (e.g. electrode clip and wrist strap for the portable meter) are already included in the standard delivery of the instruments. Other accessories such as electrodes, solutions, electrode arm for the bench meter or a carry bag for the portable meter are of course also available from METTLER TOLEDO.

You will find much helpful functionality in these meters – the Education Line is the best start into the pH world. A few of many exiting features are:

- Automatic endpoint format which freezes the measurement automatically when the value is stable;
- Three predefined buffer groups to select from;
- Automatic buffer recognition which allows you to calibrate with your buffers in any order;
- IP54 design (EL2).

## 2. Safety measures

### Measures for your protection

- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

### Measures for operational safety



- Do not unscrew the two halves of the housing!
- Dry off any liquid spills immediately! The instrument is not watertight (EL2 is IP54)!
- For EL2 instruments use batteries of the specified type only. Otherwise, proper operation cannot be guaranteed.
- Exclude the following environmental influences:
  - powerful vibrations,
  - direct sunlight,
  - atmospheric humidity greater than 80%,
  - corrosive gases present,
  - temperatures below 5 °C and above 40 °C,
  - powerful electric or magnetic fields!


### 3. Installation

Carefully unpack the meter. Keep the calibration certificate in a safe place.  
Insert the right adapter clip into the power adapter slot:


#### 3.1 Installing the electrode arm (optional for EL20)

The electrode arm can be fastened to EL20 on the left or the right side.


1. Remove the electrode stand rubber cover on that side of the bench meter where you prefer the stand.  
Keep the cover in a safe place;
2. Insert the electrode stand firmly into the opening as shown;
3. Fix the opening of the arm upon the stand from above, adjust it to a moderate height.

#### 3.2 Installing the batteries (EL2)

1. Pull the fixing point of the battery cover in the direction of the arrow, hold the lid with two fingers and remove it;
2. Insert the batteries in the battery compartment as shown;
3. Replace the battery cover.


### 3.3 FiveGo™ clip (EL2)

The FiveGo™ clip (Part no.51302960) is an electrode holder that can be placed next to the display on either side of the housing. To mount the clip, remove the cover from the clip's fixing point. Attach the clip by pressing it into the recess. Slide the shaft of the sensor into the clip from above.


### 3.4 Installing the IP54 cover (EL2)

Use the small plastic tool as shown:

## 4. Operating the Education Line pH meter

### 4.1 Display and key controls

1 Electrode condition (For electrode maintenance see section 5.2)

Slope: 95 - 105 %  
and offset:  $\pm$  (0 - 15) mV  
Electrode is in  
good condition

Slope: 90 - 94 %  
or offset:  $\pm$  (15 - 35) mV  
Electrode needs  
cleaning

Slope: 85 - 89 %  
or offset:  $\pm$  (>35) mV  
Electrode is defective

2 Calibration icon - calibration in progress

3 Measurement icon - measurement or calibration is running

4 Setup icon - instrument is in the setup mode

5 Battery status icon (EL2 only) shows the condition of the batteries - fully charged, half-charged or  
fully discharged (To replace batteries, see section 3.2)

6 pH/mV reading or slope in calibration process

7 Auto / Manual temperature compensation (see section 4.4)

8 Endpoint stability / Auto endpoint

9 Temperature during measurement or offset in calibration process

10 Calibration point Cal / Buffer group number MR / Memory number MR (EL2 only) /  
Error index Err

**Keys for EL20**

	<b>Press &amp; release</b>		<b>Press &amp; hold for 3 seconds</b>	
	- Start or endpoint measurement - Confirm setting, store entered value		- Turn auto endpoint on / off	
	- Start calibration		- Review the latest calibration data	
	- Meter on - Back to measurement screen		- Meter off	
	- Enter setup mode - Increase value during setting			
	- Switch between pH and mV measuring modes - Decrease value during setting			
	- Start self-diagnosis			

**Keys for EL2**

	<b>Press &amp; release</b>		<b>Press &amp; hold for 3 seconds</b>	
	- Start or endpoint measurement - Confirm setting, store entered value		- Turn auto endpoint on / off	
	- Start calibration		- Review the latest calibration data	
	- Meter on - Back to measurement screen		- Meter off	
	- Store current reading to memory - Increase value during setting - Scroll up through the memory		- Recall stored data	
	- Switch between pH and mV measuring modes - Decrease value during setting - Scroll down through the memory		- Enter setup mode	
	- Start self-diagnosis			

## 4.2 Calibration

### 4.2.1 Buffer groups

The Education Line pH meters allow you to perform 1- or 2- point calibrations. After selecting your calibration buffer group from one of the three predefined buffer groups in the meter, the buffers are automatically recognized and displayed during calibration (auto buffer recognition).

The three predefined groups are:

B1	1.68	4.01	7.00	10.01		(at 25°C)
B2	2.00	4.01	7.00	9.21	11.00	(at 25°C)
B3	1.68	4.00	6.86	9.18	12.46	(at 25°C)

Tables (B1...B3) for automatic temperature compensation are programmed in the meter for each buffer group (see also Appendix).

### 4.2.2 Selecting a predefined buffer group

Press the **Setup** key (EL20) or press and hold the **Mode/Setup** key (EL2) until the setup icon appears on the display and the MTC temperature blinks. Press **Read** to ignore – the current temperature unit blinks – press **Read** to ignore again. When the current buffer group blinks, the buffers belonging to this buffer group appear alternating on the display. Use **▲** or **▼** to select another buffer group and press **Read** to confirm your selection. The meter will exit automatically to the measurement screen.

### 4.2.3 Performing a 1-point calibration

Place the electrode in a calibration buffer and press **Cal**. The calibration icon and the measurement icon appear on the display. The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The relevant buffer value is displayed and stored; the measurement icon disappears from the display.

To finish the calibration and return to the sample measurement, press **Read**. The offset value and the slope are then shown on the display for 3 seconds. To reject the calibration, press **Exit** before the meter returns automatically after 3 seconds to the measurement screen.

#### Note:

With the 1-point calibration only the offset is adjusted. If the sensor was previously calibrated with multi-point calibration the previously stored slope will remain. Otherwise theoretical slope (-59.16 mV / pH) will be used.

### 4.2.4 Performing a 2-point calibration

- Step 1: Perform the first calibration as described above in "Performing a 1-point calibration".
- Step 2: Rinse the electrode with deionized water.
- Step 3: Place the electrode in the next calibration buffer and press **Cal**

The measurement icon appears. The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The relevant buffer value is displayed and stored; the measurement icon disappears from the display.

The offset value and the slope are then shown on the display for 3 seconds. To reject the calibration, press **Exit** before the meter returns automatically after 3 seconds to the measurement screen.

**Note:**

The use of a temperature sensor or an electrode with a built-in temperature sensor is recommended. If you use the MTC mode, you should enter the correct temperature value and keep all buffer and sample solutions at the set temperature.

To ensure the most accurate pH readings, you should perform a calibration once a day.

## 4.3 Sample measurement

### 4.3.1 Performing a pH measurement

Place the electrode in the sample and press **Read** to start the measurement: the measurement icon appears on the display and the decimal point blinks. The display shows the pH value of the sample. The automatic endpoint **A** is the default setting of the meter. When the signal has stabilized, the display freezes automatically, **/A** appears and the measurement icon disappears.

By pressing and holding **Read**, you can switch between the auto and manual endpoint modes. To manually endpoint a measurement, press **Read**: the display freezes and **/** appears.

**Stability criterion for pH and mV measurement** – The signal of the sensor input may not change by more than 0.1 mV in 6 seconds.

### 4.3.2 Performing a mV / Redox measurement

To switch to mV mode, press **Mode**.

To perform a mV measurement, follow the same procedure as for pH measurement.  
(For Redox measurements, connect a Redox sensor first.)

## 4.4 Temperature measurement

For better accuracy, we recommend the use of either a built-in or a separate temperature probe. If a temperature probe is used, **ATC** and the sample temperature are displayed.

**Note:**

The meter accepts an NTC 30 k $\Omega$  temperature sensor.

## Operating the pH meter

### 4.4.1 Manual temperature compensation

If the meter does not detect a temperature probe, it automatically switches to the manual temperature compensation mode and MTC appears.

To set the MTC temperature and the temperature unit, press **Setup** (EL20) or press and hold **Mode/Setup** (EL2) until the setup icon appears on the display and the MTC temperature blinks. Use **▲** or **▼** to increase or decrease the value to the temperature of your sample. Press **Read** to confirm the setting. The meter proceeds automatically to the temperature unit setting. Press **Exit** to leave the setup without changing the temperature unit or use **▲** or **▼** to select between °C (Celsius) or °F (Fahrenheit). Press **Read** to confirm the setting. The meter proceeds automatically to the calibration setting; press **Exit** to leave to setup. The default setting is 25 °C.

### 4.5 Analog output (EL20)

The analog output supplies a voltage proportional to the displayed measured value (pH, mV). The polarity of the output corresponds to the polarity at the measurement electrode. The voltage at the analog output is not temperature-compensated.

### 4.6 Using the memory (EL2)

#### 4.6.1 Storing a reading

The EL2 pH meter can store up to 30 endpointed results. Press **STO** when the measurement has endpointed. **M01** indicates that one result has been stored.

If you press **STO** when **M30** is displayed, **FUL** indicates the memory is full. To store further data you will have to clear the memory (see below).

#### 4.6.2 Recalling from memory

Press and hold **RCL** to recall the stored values from memory when the current measurement has endpointed.

Press **▲** or **▼** to scroll through the stored results. **R01** to **R30** indicates which result is being displayed. Press **Read** to exit.

#### 4.6.3 Clearing the memory

Continue pressing **▲** or **▼** to scroll through the stored results until **MRCL** appears. Then press **Read**, **CLR** blinks. Press **Read** again to confirm the deletion or press **Exit** to return to the measurement mode without deleting the data.

## 4.7 Self-diagnosis

Press and hold **Read** and **Cal** simultaneously until the meter displays the full screen. Each icon blinks one after the other. This way you may check whether all icons are correctly shown. The next step is to check that the keys are functioning correctly. This requires user interaction.

When **b** blinks, five icons are displayed. Press the five keys in any order. Each time you press a key an icon disappears from the screen, continue to press the other keys until all the icons have disappeared.

When the self-diagnosis has been completed successfully, **PAS** appears. If self-diagnosis fails, error message **Err 1** appears (see section 4.8).

**Note:**

You have to finish pressing all five keys within two minutes, otherwise **Err 1** appears and you will have to repeat the procedure.

## 4.8 Error message

Error 0	Memory access error	Reset to factory settings
Error 1	Self-diagnosis failed	Repeat the self-diagnosis procedure and make sure that you finish pressing all five keys within two minutes.
Error 2	Measured values out of range	Make sure that the electrode wetting cap has been removed and the electrode is properly connected and placed in the sample solution. If no electrode is connected, put the shorting plug in the socket.
Error 3	Measured buffer temperature out of range (5 ... 40 °C)	Keep the buffer temperature within the range for calibration
Error 4	Offset out of range	Make sure you have the correct buffer and that it is fresh. Clean or replace the electrode.
Error 5	Slope out of range	Make sure you have the correct buffer and that it is fresh. Clean or replace the electrode.
Error 6	Meter cannot recognize the buffer	Make sure you have the correct buffer and that it is fresh. Check that the buffer has not been used more than once during the calibration.
Error 9	The current data set has already been stored once	A measurement can only be stored once. Perform a new measurement to store a new data set.

## 5. Maintenance

### 5.1 Meter maintenance

Never unscrew the two halves of the housing!

The Education Line instruments do not require any maintenance other than occasional wipe with a damp cloth and the replacement of depleted batteries (EL2).

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is attacked by some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

### 5.2 Electrode maintenance

Make sure the electrode is filled with electrolyte solution if applicable.

For maximum accuracy, any filling solution that may have "crept" and encrusted the outside of the electrode should be removed with deionized water.

Always store the electrode according to the manufacturer's instructions and do not allow it to dry out. If the electrode slope falls rapidly, or if the response becomes sluggish, the following procedures may help. Try one of the following, depending the samples usually measured.

1. For fat or oil build-up, degrease the membrane with cotton wool soaked in either acetone or a soap solution.
2. If the sensor membrane has dried out, soak the tip of the electrode in 0.1 M HCl overnight.
3. If a protein build-up has occurred in the diaphragm, remove deposits by soaking the electrode in an HCl /pepsin solution (51340068).
4. If a silver sulfide contamination has occurred, remove deposits by soaking the electrode in a solution of thiourea (51340070).

After treatment a new calibration should be performed.

**Note:**

Cleaning and filling solutions should be handled with the same care as that given to toxic or corrosive substances.

### 5.3 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environment protection.

## Accessories

**6. Accessories**

Item	Order No.
Education Line pH bench meter	51302930
FiveEasy™ Power supply	51302950
FiveEasy™ electrode arm	51302951
Rubber covers (for electrode arm hole)	51302952
Shorting plug BNC	51302859
Education Line pH potable meter	51302920
FiveGo™ electrode clip	51302960
Rubber clip covers (for electrode clip slot)	51302961
Education Line Battery cover	51302963
Wrist strap	51302964
FiveGo™ carry bag	51302966
IP54 electrode caps	51302965
Sealing kit (O-rings for battery case)	51302967
Quick guide EL20	51710492
Quick guide EL2	51710494
pH 4.01 buffer sachets, 30 x 20mL	51302069
pH 4.01 buffer solution, 6 x 250mL	51340058
pH 7.00 buffer sachets, 30 x 20mL	51302047
pH 7.00 buffer solution, 6 x 250mL	51340060
pH 9.21 buffer sachets, 30 x 20mL	51302070
pH 9.21 buffer solution, 6 x 250mL	51300194
pH 10.01 buffer sachets, 30 x 20mL	51302079
pH 10.01 buffer solution, 6 x 250mL	51340231
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)	51302068
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)	51302080
LE407, combination pH sensor, plastic, gel electrolyte	51340330
LE409, combination pH sensor, glass, refillable	51340331
LE420, special pH sensor with PTFE junction, refillable	51340332
LE438, 3-in-1 pH sensor, plastic, ATC	51340242
LE501, Redox-sensor, glass	51340338
LE510, Redox-sensor, plastic	51340339
ATC probe, temperature sensor	51300164
HCl/Pepsin solution(removes protein contamination)	51340068
Thiourea solution (remove silver sulfide contamination)	51340070
Reactivation solution for pH electrodes	51340073

## 7. Specifications

	<b>Education Line pH meter EL20</b>	<b>Education Line pH meter EL2</b>
<b>Measurement range</b>	pH 0.00...14.00	pH 0.00...14.00
	-1999...1999 mV	-1999...1999 mV
	0 °C to 100 °C	0 °C to 100 °C
<b>Resolution</b>	0.01 pH	0.01 pH
	1 mV	1 mV
	0.1 °C	0.1 °C
<b>Limits of error</b>	± 0.01 pH	± 0.01 pH
	± 1 mV	± 1 mV
	± 0.5°C	± 0.5°C
<b>pH calibration</b>	Up to 2 points	Up to 2 points
<b>Isopotential point</b>	pH 7.00	pH 7.00
<b>Calibration buffer</b>	3 predefined groups	3 predefined groups
<b>Power requirements</b>	9 V DC, 0.4 W ± 10% For use with CSA Certified (or equivalent approved) power source, which must have a limited circuit output.	Rating : 6VDC, 5mA Powered with batteries : 4xAAA/LR3 1.5V or NiMH 1.2 rechargeable
<b>Size/weight</b>	200 x 175 x 52 mm / 0.6 kg	169 x 82 x 36 mm / 0.18 kg (without batteries)
<b>Display</b>	Liquid crystal	Liquid crystal
<b>pH input</b>	BNC, impedance > 10e+12 Ω	BNC, impedance > 10e+12 Ω
<b>T input</b>	Cinch, NTC 30 kΩ	Cinch, NTC 30 kΩ
<b>Reference input</b>	2mm banana socket	-
<b>Analog output</b>	2mm banana socket	-
<b>IP rating</b>	For indoor use only.	IP54 for system
<b>Battery life</b>	-	>500 working hours
<b>Ambient conditions</b>	Temperature: 5...40 °C	Temperature: 5...40 °C
	Relative humidity: 5 %...80 % (non-condensing)	Relative humidity: 5 %...80 % (non-condensing)
	Height up to 2000 m	Height up to 2000 m
	Installation category: II	Installation category: II
	Pollution degree: 2	Pollution degree: 2
<b>Materials</b>	Housing: ABS/PC reinforced	Housing: ABS/PC reinforced
	Window: Membrane	Window: polymethylmethacrylate(PMMA)
	Keypad: Membrane	Keypad: Membrane

## Appendix

### 8. Appendix

#### 8.1 Buffer tables

Education Line pH meters automatically correct for the temperature dependence of the buffer pH using the values given in the tables.

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

**Quality certificate.** Development, production and testing according to ISO9001.  
Environmental management system according to ISO14001.

**Worldwide service.** Our extensive service network is among the best in the world and ensures maximum availability and service life of your product.

**European conformity.** The CE conformity mark provides you with the assurance that our products comply with the most recent EU directives.

**On the Internet.** You will quickly find lots of essential information about our products, our services, and our company at  
<http://www.mt.com>

Subject to technical changes and to the availability of the accessories supplied with the instruments.

\*P51710490\*

© Mettler-Toledo AG 2007 ME-51710490 Printed in China 0702/2.52

**Mettler-Toledo AG, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach,**  
Tel. +41 44 806 77 11, Fax +41 44 806 73 50, Internet:<http://www.mt.com>