Enhanced InTrac[®] Series – When process safety is crucial!

Technical data



InTrac retractable housings from METTLER TOLEDO with TRI-LOCK[™] safety system enable the installation of electrodes/sensors in the chemical, pharmaceutical and food & beverage industry. Sensors can be calibrated and cleaned without interrupting the process which dramatically increases sensor life and reduces intensive maintenance.

Features

- TRI-LOCK™ safety system
- Cleaning and calibration without interrupting the process
- $-\ensuremath{\,\text{Various}}$ immersion depths to suit reactors, tanks or pipes sizes
- Large selection of materials of construction according to demanding application requirements
- Flexible process adaptations for the chemical, pharmaceutical and food & beverage industry
- Sterile double rinsing chamber (InTrac79Xe)
- Easy installation and withdrawal while the process is running
- ATEX (94/9/EC)

InTrac	7	7	7	е

Content

Operation and System	2
General operation	2
TRI-LOCK™ safety system	3
Complete measuring system	4
Pneumatic and other connections	5

Dimensional drawings	6
Dimensional drawings InTrac776e	6
Dimensional drawings InTrac777e/InTrac779e	7
Dimensional drawings InTrac797e/InTrac799e	8
Process connections	9

Technical specifications	12
Technical specifications InTrac776e	12
Suitable electrodes	13
Technical specifications InTrac77Xe/79Xe	14
Suitable electrodes/sensors	15

Ordering information	16
Product structure InTrac776e/InTrac777e/ InTrac779e	16
Product structure InTrac797e/InTrac799e	17
Accessories	19
Spare parts	19
Additional documentation	19



General operation

With InTrac retractable housings reliable measurements of pH/ORP, dissolved oxygen, CO₂, conductivity, and turbidity can be achieved in process applications. The large selection of materials in contact with the media and flexible process adaptation makes the housings suitable for applications in the chemical and the pharmaceutical industry and also meets the demands of engineering companies. The sterilizable InTrac79X e housing is designed to meet the highest demands of the food & beverage and pharmaceutical industry where sterile conditions are required. The housing can be operated manually or pneumatically and allows electrodes/sensors to be:

- separated from the process and retracted into the flushing chamber
- rinsed with water or cleaning solution
- calibrated
- completely un-installed
- sterilized (InTrac79Xe).


Maintenance position

InTrac housings are available in stainless steel DIN 1.4404/AISI 316L, DIN 2.4602/AISI Alloy C22 (i.e. Hastelloy C22) or titanium as well as in Polypropylene (PP), Polyvinylidenfluoride (PVDF) or Polyetheretherketone (PEEK) as wetted parts. Housings with short insertion length (70 mm or 100 mm) for pipe installation and longer insertion length (200 mm) for in-stallation in a reactor can be selected. A wide range of process fittings are available to meet your requirements (see section «process connections»). The InTrac79Xe is only available in stainless steel and with a short insertion length, for pipe installation > DN50 (2") or installation in a reactor.

Position indicator «Measurement» Position indicator «Maintenance» Additional O-ring

Measuring position


Electrode



prevents contamination

TRI-LOCK™ safety system


With TRI-LOCK[™] safety system from METTLER TOLEDO the risk of unintentional insertion of the electrode/sensor into the process or unintentional replacement of the electrode in the measuring position is minimized. This substantially increases the process safety of your measuring loop. The use of optional position indicators further improves operation safety. The pneumatic or inductive position indicators reliably signal that the electrode is in the measuring or in the maintenance position.



The positive position rod prevents the insertion tube from turning and therefore keeps it in an optimized position for flushing the sensor but also keeps the sensor/electrode in the best measuring position possible. In addition it reduces the risk of unintentional remove of the sensor/electrode while in the measuring position. The positive position rod is optional with all InTrac777e/779e housings as well as with InTrac797e and 799e housings.

The sensor detector prevents the housing from being inserted in the process without an installed sensor. The risk of unintentional release of process media is therefore eliminated. The sensor detector is standard with all InTrac housings.


The patented design of the immersion tube H = 100/200 guarantees in any position, safe containment of the process medium during sensor insertion and retraction.




The patented immersion tube comes as standard with all housings with insertion length of 100 mm and 200 mm. Housings with an insertion length of 70 mm do not offer this feature.

Complete measuring system

Basic measuring system



Measuring system combined with fully automatic EasyClean cleaning and calibration system



In combination with EasyClean cleaning and calibration system you can fully automate your pH measurement. The measuring system consists of an InTrac housing with postion indicators, a sensor, a transmitter and the EasyClean cleaning and calibration system. For complete installation in hazardous areas an EasyClean 400X is required in combination with a fully approved housing.


Pneumatic connections (continuously adjustable)

For pneumatic operation of the housing, connections 1 and 2 are to be linked to the air lines 1 and 2 as shown below.

Attention! Only the original fittings supplied with the housing may be used. If the position indicator system is to be used, the two connections **3** and **4** are to be linked to the position indicator valves **3** and **4**.

5/2-way valve

3/2-way valve



The operating pressure of the pneumatic system is 4...8 bars. The control air must be water and oil-free, and filtered (40 μ m). The hoses must have a minimum size of 6/4 mm. The speed of the movement can be continuously adjusted.

Note!

If the air pressure can go above 8 bars (also short spikes), the use of a pressure reducer is recommended.

Connections for position indicators:






left: pneumatic 3/2-way position indicator

inductive position indicator (Ex: 2-wire, non-Ex: 3-wire)

The lower connections are for the function «Measurement» the upper connections for «Maintenance».

right:




Manually operated InTrac776e with 70 mm (left), 100 mm (middle), and 200 mm (right) insertion length for pH/ORP electrodes with liquid electrolyte.

	Required clearance for installation	Overall length when inserted	Insertion length (DN25 type)	Insertion length (flange type)
InTrac776e/70	710 mm/27.95"	545 mm/21.46"	71 mm/2.78"	max. 46 mm*/1.81"
InTrac776e/100	710 mm/27.95"	545 mm/21.46"	100 mm/3.95"	max. 76 mm*/2.97"
InTrac776e/200	1110 mm/43.7"	645 mm/25.4"	200 mm/7.89"	max. 176 mm*/6.91"

The dimensions may slightly vary depending on the material and the process adaptation.

* Valid for DIN DN50-PN16 standard flange. With larger flanges the insertion length will become shorter.




Manually operated InTrac777 e with 70 mm (left), 100 mm (middle) and 200 mm (right) insertion length.

	Required clearance for installation	Overall length when inserted	Insertion length (DN25 type)	Insertion length (flange type)
InTrac777e/70	515 mm/20.3"	360 mm/14.2"	71 mm/2.78"	max. 44 mm*/1.73"
InTrac777e/100	515 mm/20.3"	360 mm/14.2"	100 mm/3.95"	max. 76 mm*/2.97"
InTrac777e/200	915 mm/36.0"	460 mm/18.11"	200 mm/7.89"	max. 176 mm*/6.91"

The dimensions may slightly vary depending on the material and the process adaptation.

* Valid for DIN DN50-PN16 standard flange. With larger flanges the insertion length will become shorter.



Manually operated InTrac797 e (left) and InTrac799 e (right).

	Required clearance for installation	Overall lengt when inserted	Insertion length (DN25 type)
InTrac797 e	715 mm/28.15"	460 mm/18.10"	92 mm/3.62"
InTrac799e	715 mm/28.15"	460 mm/18.10"	92 mm/3.62"

The dimensions may slightly vary depending on the material and the process adaptation.

Process connections

InTrac housings consist of a standardized construction with variable process connections.

Material		1.4404/	2.4602	Ti	PP	PVDF	PEEK
		316L	Alloy C22				
Process connection	Dimension						
For Ingold							
weld-in socket DN25	DN25	•	b)	b)	•	•	•
Flange DIN	DN32 PN16	•	a)	a)			
	DN40 PN16	•	a)	a)	•	٠	•
	DN50 PN16	•	a)	a)	•	•	•
	DN80 PN16	•	a)	a)	•	•	•
Flange ANSI	A150-1.5"	•	a)	a)	•	٠	•
	A150-2"	•	a)	a)	•	٠	•
	A150-3"	•	a)	a)	٠	•	•
NPT thread	NPT 1"	•	•	•	•	•	•
	NPT 1 ¹ /4"	•					
	NPT 1 ¹ /2"	•					
Tri-Clamp	TC 1.5"	•					
	TC 2"	•					
Tuchenhagen	Varivent DN50	•					

a) Flanges for housings made of DIN 2.4602/AISI Alloy C22 or titanium are supplied in DIN 1.4404/AISI 316L stainless steel together with a fitting plate made of PTFE (wetted part).


b) Cap nut for Ingold sockets are made of DIN 1.4404/AISI 316L stainless steel.

Other process connections on request.

Process connection for DN25 sockets

SS-type for InTrac776e/777e/779e*

SS-type for InTrac797e/799e* Polymer-type for InTrac776e/777e




* For housings used with turbidity sensors (InTrac 779e, InTrac 799e) please refer to the drawing page 8 for the different immersion tube.

Dimensional drawings

Process connection through DIN flange

SS-type for InTrac776e/777e/779e* Alloy C22/Ti-type for InTrac776e/777e/779e* Polymer-type for InTrac776e/777e



Process connection through ANSI flange


max. 75.5 / 175.5

SS-type for InTrac776e/777e/779e*

Alloy C22/Ti-type for InTrac776e/777e/779e*

max. 72.5 / 172.5

Polymer-type for InTrac776e/777e



* For housings used with turbidity sensors (InTrac 779e, InTrac 799e) please refer to the drawing page 8 for the different immersion tube.

Dimensional drawings


InTrac Series


Process connection through NPT thread


SS-type NPT 1"/1¹/4" for InTrac776e/777e/779e*

71

SS-type NPT 1¹/2" for InTrac776e/777e/779e*








Process connection through Tri-Clamp flange


SS-type for InTrac776e/777e/779e*

SS-type for InTrac797 e/799 e*



Process connection through Tuchenhagen Varivent DN50


SS-type for InTrac776e/777e/779e* SS-type for InTrac797 e/799e*



* For housings used with turbidity sensors (InTrac 779e, InTrac 799e) please refer to the drawing page 8 for the different immersion tube.

Technical specifications InTrac776 e

Ambient conditions						
Temperature	PP-body: 07	0 °C/32158 °	F			
	SS-body: -10	.70 °C/14158	₿ °F			
Process conditions						
Functional pressure range	manual: 05 bar / 073 psig					
	pneumatic: 0	.8 bar / 0116 p	osig			
Max. permissible pressure [PS]/[TS]*	PP	6 bar/20°C	87 psig/68 °F			
(linear decreasing)		0 bar/80°C	0 psig/176 °F			
	PVDF	6 bar/20°C	87 psig/68 °F			
		1 bar/110°C	15 psig/230 °F			
	PEEK	6 bar/20°C	87 psig/68 °F			
		1 bar/110°C	15 psig/230 °F			
	1.4404/316L/					
	2.4602/					
	Alloy C22, Ti	6 bar/130°C	87 psig/266 °F			
Operation	manual or pneu	umatic				
Insertion length	70 mm, 100 m	nm, 200 mm				
Wetted parts	DIN 1.4404/AIS	SI 316L, DIN 2.46	602/Alloy C22, titanium			
	PP, PVDF, PEEK					
Wetted O-rings		DM-FDA, Kalrez®-				
Not wetted parts	body: Polypropylene (PP) conductive or DIN 1.4404/AISI 316L					
	protective sleeve: Polypropylene (PP) conductive					
Weight	approx. 4.5 kg					
Outer dimensions						
Short version	length: approx. 545 mm in measuring position					
(70/100 mm)	length: approx. 710 mm for electrode removal (minimum)					
Long version	length: approx. 645 mm in measuring position					
(200 mm)			ectrode removal (minimum)			
Pneumatic connections	48 bar / 58.					
Air quality to ISO 8573-1	– air moisture	content class 4 (a	dew point +3°C)			
	- solids class					
		tent class 2 (0.1				
		ns for hoses 6/4	mm			
Flushing connections	26 bar / 29.					
(Water, steam or cleaning solutions)	1x connection I	N: thread G ¹ /8"				
		OUT: thread G 1/4				
	1x connection	TEMP: thread G 1	/8 ["]			
Position monitoring	pneumatic chea	ck-back (3/2 way	-valve); G ¹ /8"			
	inductive check	-back, non-Ex, M	12x1			
	inductive check	-back, Ex, M12x	1			



*Pressure rating according to PED 97/23/EC

Technical specifications

Pressure information	According to PED-Article 1, Section 2.2: «Pressure» is referenced to atmospheric pressure, e.g. an overpressure. Accordingly, a pressure in the vacuum region will be expressed as a negative pressure.
Explosion protection	According to ATEX guidelines (94/9/EG): 🐵 II 1/2 GD IIC T6/T5/T4/T3
(Valid for all housings with medium-wetted	IP6X T 69 °C/T 81 °C/T 109 °C/T 161 °C
parts made of metal)	SNCH 03 ATEX 3590 X
	According to FM guidelines: IS CL I, II, III, Div 1, GR ABCDEFG/T6
Certificates/approvals	Declaration of conformity CE
	Pressure equipment directive (PED) 97/23/EC
	Certificate of compliance with the order EN10204-2.1
	Inspection certificate 3.1B
	ATEX (94/9/EC), FM certificate, MaxCert™


Suitable electrodes

Insertion length	Length of electrode
H = 70/100 mm	250 mm
H = 200 mm	450 mm
pH electrodes	InPro2000
	465-50
ORP electrodes	Pt4865-50

InTrac Series

Technical specifications InTrac77Xe/79Xe

Ambient conditions			
Temperature	PP-body: 070 °C / 32158 °F		
	SS-body: -1070 °C / 14158 °F		
Ambient conditions	<u>.</u>		
Functional pressure range	manual: 05 bar / 073 psig		
	pneumatic: 08 bar / 0116 psig		
Max. permissible pressure [PS]/[TS]*	PP	6 bar/20°C	87 psig/68 °F
(linear decreasing)		0 bar/80°C	0 psig/176 °F
	PVDF	6 bar/20°C	87 psig/68 °F
		1 bar/110°C	15 psig/230 °F
	PEEK	6 bar/20°C	87 psig/68 °F
		1 bar/110°C	15 psig/230 °F
	1.4404/316L		
	2.4602/		
	Alloy C22, Ti	16 bar/140°C	232 psig/284 °F
Operation	manual or pneu	umatic	
Insertion length			
InTrac77Xe:	70 mm/100 mm/200 mm		
InTrac79Xe:	100 mm		
Wetted parts			
InTrac77Xe:	DIN 1.4404/AISI 316L, DIN 2.4602/AISI Alloy C22,		
	titanium, PP, PVDF, PEEK		
InTrac79Xe:	DIN 1.4404/AISI 316L		
Wetted O-rings	Viton®-FDA, EPDM-FDA, Kalrez®-FDA		
Non- wetted parts	body: Polypropylene (PP) conductive or DIN 1.4404/AISI 316L		
	protective sleeve: Polypropylene (PP) conductive		
Weight	approx 4.5 kg		
Outer dimensions			
Short version 77Xe	• • • •	360 mm in meas	
(70/100 mm)	length: approx. 515 mm for electrode removal (minimum)		
Long version 77Xe	length: approx. 460 mm in measuring position		
(200 mm)	length: approx. 915 mm for electrode removal (minimum)		
Version 79Xe	length: approx. 460 mm in measuring position		
(100 mm)	length: approx. 715 mm for electrode removal (minimum)		
Pneumatic connections	48 bar / 58116 psig		
Air quality to ISO 8573-1		content class 4 (d	lew point +3°C)
	– solids class 5 (filter 40 μm)		
	– max. oil cont	tent class 2 (0.1 i	mg/m ³)
	- air connectio	ns for hoses 6/4 i	mm



*pressure rating according to PED 97/23/EC

Technical specifications

Flushing connections		
InTrac77Xe	26 bar / 2987 psig	
(Water, cleaning solutions)	1 x connection IN: thread G ¹ /8"	
	1 x connection OUT: thread $G^{1/4}$ "	
	1 x connection TEMP: thread G ¹ /8"	
InTrac79Xe	26 bar / 2987 psig	
(Water, cleaning solutions, steam)	5 x connections: thread G ¹ /8"	
Position monitoring	pneumatic check-back (3/2 way-valve); G ¹ /8"	
·	inductive check-back, non-Ex, M12x1	
	inductive check-back, Ex, M12x1	
Pressure information	According to PED-Article 1, Section 2.2:	
	«Pressure» is referenced to atmospheric pressure, e.g. an overpressure.	
	Accordingly, a pressure in the vacuum region will be expressed	
	as a negative pressure.	
Explosion protection	According to ATEX guidelines (94/9/EG): 🐼 II 1/2 GD IIC T6/T5/T4/T3	
(Valid for all housings with medium-wetted	IP6X T 69 °C/T 81 °C/T 109 °C/T 161 °C	
parts made of metal)	SNCH 03 ATEX 3590 X	
	According to FM guidelines: IS CL I,II,III, Div 1, GR ABCDEFG/T6	
Certificates/approvals	Declaration of conformity CE	
	Pressure equipment directive (PED) 97/23/EC	
	Certificate of compliance with the order EN10204-2.1	
	Inspection certificate 3.1B	
	ATEX (94/9/EC), FM certificate, MaxCert™	

Suitable electrodes/sensors

InTrac777 e	
Insertion length	Length of sensor/electrode
H = 70/100 mm	220 mm or 225 mm
H = 200 mm	420 mm or 425 mm
рН	InPro3200, InPro4250, InPro4800
ORP	Pt4805
Dissolved oxygen	InPro6800/InPro6900 (220 mm only)
Conductivity	InPro7001-VP (225 mm only)

InTrac779e


Insertion length	Length of sensor
H = 70/100 mm	205 mm
H = 200 mm	407 mm
Turbidity	InPro8100, InPro8200

InTrac797 e

Insertion length	Length of sensor/electrode	
H = 100 mm	320 mm or 325 mm	
рН	InPro3030, InPro3100, InPro3300	
ORP	Pt4805	
CO ₂	InPro5000	
Dissolved oxygen	InPro6800/InPro6900 (320 mm only)	

InTrac799e Insertion length	Length of sensor
H = 100 mm	297 mm
Turbidity	InPro8100, InPro8200


Product structure InTrac776e/777e/779e



Housings made of Alloy C22 or titanium are supplied with cap nut made of stainless steel as standard.
 Flanges are also made of stainless steel and supplied with a sealing plate made of PTFE.

- Varivent and Tri-Clamp flanges are available in stainless steel only.

Product structure InTrac797e/799e



Housings with double flushing chamber are only available with DN25, Tri-Clamp or Varivent process connections.

Ordering information

InTrac Series

InTrac776e	Designation	Order no.
	InTrac776M/070/4404/D00/Vi/A00	52 403 220
	InTrac776M/100/4404/D00/Vi/A00	52 403 221
	InTrac776M/200/4404/D00/Vi/A00	52 403 222
	InTrac776P/070/4404/D00/Vi/A00	52 403 223
	InTrac776P/100/4404/D00/Vi/A00	52 403 224
	InTrac776P/200/4404/D00/Vi/A00	52 403 225
	InTrac776M/070/PVDF/D00/Vi/A00	52 403 226
	InTrac776M/100/PVDF/D00/Vi/A00	52 403 227
	InTrac776P/070/PVDF/D00/Vi/A00	52 403 228
	InTrac776P/100/PVDF/D00/Vi/A00	52 403 229
	InTrac776R/100/4404/D00/Vi/A00	52 403 230
	InTrac776I/100/4404/D00/Vi/A00	52 403 231
	InTrac776M/070/4404/D00/Ka/A00	52 403 232
InTrac777 e		
	InTrac777M/070/4404/D00/Vi/A00	52 403 205
	InTrac777M/100/4404/D00/Vi/A00	52 403 206
	InTrac777M/200/4404/D00/Vi/A00	52 403 207
	InTrac777P/070/4404/D00/Vi/A00	52 403 208
	InTrac777P/100/4404/D00/Vi/A00	52 403 209
	InTrac777P/200/4404/D00/Vi/A00	52 403 210
	InTrac777M/070/PVDF/D00/Vi/A00	52 403 211
	InTrac777M/100/PVDF/D00/Vi/A00	52 403 212
	InTrac777P/070/PVDF/D00/Vi/A00	52 403 213
	InTrac777P/100/PVDF/D00/Vi/A00	52 403 214
	InTrac777R/100/4404/D00/Vi/A00	52 403 215
	InTrac777I/100/4404/D00/Vi/A00	52 403 216
	InTrac777M/070/4404/N01/Vi/A00	52 403 217
	InTrac777M/100/4404/N01/Vi/A00	52 403 218
	InTrac777M/200/4404/N01/Vi/A00	52 403 219
InTrac779e		
	InTrac779M/070/4404/D00/Vi/A00	52 403 236
InTrac797 e		
	InTrac797M/100/4404/D00/Vi/A00	52 403 233
	InTrac797P/100/4404/D00/Vi/A00	52 403 235
	InTrac797M/100/4404/768/Vi/A00*	52 403 234
InTrac799 e		
	InTrac799M/100/4404/D00/Vi/A00	52 403 237

* with 768 adapter

All other configurations can be ordered according to the product key from pages 16 and 17.

Ordering information

InTrac Series

Accessories	Designation	Order no.
	Pt100-764/5 m	10 1003 1046
	Pt1000-764/5 m	10 1003 1166
	3/2-way reversion valve set (2 pcs), pneumatic	52 401 324
	Inductive check-back, non-Ex (2 pcs)	52 402 398
	Inductive check-back, Ex (2 pcs)	52 402 399
	Fitting Set SERTO	52 400 842
	Fitting Set SWAGELOK	52 400 843
	Fitting Set GYROLOK	52 400 844
	Fitting Set PVDF	52 400 768
	Fitting Set PVC	52 400 769
	pneumatic connection (throttle valve 1 pc*)	20 102 1049
	Positive position rod, long 7XXe	52 403 121
	* 2 pcs required per housing	
spare parts		
	O-ring set Viton [®] -FDA	52 403 172
	O-ring set Kalrez [®] -FDA	52 403 173
	O-ring set EPDM-FDA	52 403 174
	Set cable grommet 5 mm/7 mm	52 403 463
aaimonal accessories		
	Weld-in sockets	
	Weld-in sockets Weld-in socket L=40/DN25/inclined	00 764 1014
		00 764 1014 00 764 1302
	Weld-in socket L=40/DN25/inclined	
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined	00 764 1302
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined	00 764 1302 00 764 1015
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=40/DN25/straight	00 764 1302 00 764 1015 00 764 1016
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined	00 764 1302 00 764 1015 00 764 1016 00 764 1017
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=40/DN25/straight Weld-in socket L=50/DN25/straight	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=40/DN25/straight Weld-in socket L=50/DN25/straight Weld-in socket L=60/DN25/straight Blind plugs	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018 00 764 1019
Additional accessories	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=40/DN25/straight Weld-in socket L=50/DN25/straight Weld-in socket L=60/DN25/straight	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=50/DN25/straight Weld-in socket L=50/DN25/straight Weld-in socket L=60/DN25/straight Blind plugs Blind plug DN25 L=40/inclined Blind plug DN25 L=40/straight	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018 00 764 1019 00 764 2002
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=50/DN25/straight Weld-in socket L=50/DN25/straight Weld-in socket L=60/DN25/straight Blind plugs Blind plug DN25 L=40/inclined Blind plug DN25 L=40/straight	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018 00 764 1019 00 764 2002 00 764 2026
	Weld-in socket L=40/DN25/inclined Weld-in socket L=44/DN25/inclined Weld-in socket L=48/DN25/inclined Weld-in socket L=55/DN25/inclined Weld-in socket L=50/DN25/straight Weld-in socket L=50/DN25/straight Weld-in socket L=60/DN25/straight Blind plugs Blind plug DN25 L=40/inclined Blind plug DN25 L=40/straight	00 764 1302 00 764 1015 00 764 1016 00 764 1017 00 764 1018 00 764 1019 00 764 2002

For further accessories please contact your local METTLER TOLEDO representative.

Together with housings and electrodes/sensors METTLER TOLEDO offers a wide range of transmitters and EasyClean automated cleaning and calibration systems. Highly automated and application-specific measuring systems to meet your needs can be provided. For further information see www.mt.com

Additional documentation

Instruction manual InTrac Series (G/E/F)	52 403 124

METTLER TOLEDO Market Organizations

Sales and Service

Australia

Mettler-Toledo Ltd. 220 Turner Street Port Melhourne AUS-3207 Victoria Phone +61 31300 659 761 Fax +61 3 9645 3935 e-mail mtausprocess@mt.com

Austria

Mettler-Toledo GmbH Südrandstrasse 17 AT-1230 Wien +43 1 604 19 80 Phone Fax +43 1 604 28 80 infoprocess.mtat@mt.com e-mail

Brazil

Mettler-Toledo Ind. e Com. Ltda. Alameda Araguaia 451 - Alphaville BR-06455-000 Barueri/SP Phone +55 11 4166 74 00 Fax +55 11 4166 74 01 e-mail sales@mettler.com.br service@mettler.com.br

China

Mettler-Toledo Instruments (Shanghai) Co. Ltd. 589 Gui Ping Road Cao He Jing CN-200233 Shanghai Phone +86 21 64 85 04 35 +86 21 64 85 33 51 Fax e-mail mtcs@public.sta.net.cn

Croatia

Mettler-Toledo d.o.o. Mandlova 3 HR-10000 Zagreb +385 1 292 06 33 +385 1 295 81 40 Phone Fax e-mail mt-zagreb@mt.com

Czech Republic

Mettler-Toledo spol s.r.o. Trebohosticka 2283/2 CZ-100 00 Praha 10 Phone +420 2 72 123 150 +420 2 72 123 170 Fax e-mail sales mtcz@mt.com

Denmark

Mettler-Toledo A/S Naverland 8 DK-2600 Glostrup +45 43 27 08 00 Phone +45 43 27 08 28 Fax info.mtdk@mt.com e-mail



Management System certified according to ISO 9001 / ISO 14001

France

Mettler-Toledo Analyse Industrielle Sarl 30, Boulevard de Douaumont RP 949 F-75829 Paris Cedex 17 Phone +33 1 47 37 06 00 +33 1 47 37 46 26 Fax e-mail mtpro-f@mt.com

Germany

Mettler-Toledo GmbH ProzeBanalvtik Ockerweg 3 D-35396 Gießen Phone +49 641 507 333 +49 641 507 397 Fax e-mail prozess@mt.com

Great Britain

Mettler-Toledo LTD 64 Boston Road, Beaumont Levs GB-LE4 1AW Leicester Phone +44 116 235 7070 +44 116 236 5500 Fax e-mail enquire.mtuk@mt.com

Hungary

Mettler-Toledo Kereskedelmi KFT Teve u. 41 HU-1139 Budapest Phone +36 1 288 40 40 +36 1 288 40 50 Fax e-mail mthu@axelero.hu

India

Mettler-Toledo India Private Limited Amar Hill, Saki Vihar Road Powai IN-400 072 Mumbai Phone +91 22 2857 0808 +91 22 2857 5071 Fax sales.mtin@mt.com

Italy

e-mail

Mettler-Toledo S.p.A. Via Vialba 42 I-20026 Novate Milanese Phone +39 02 333 321 +39 02 356 2973 Fax e-mail customercare italia@mt.com

Japan

((

Mettler-Toledo K.K. Process Division 5F Tokyo Ryutsu Center, Annex B 6-1-1 Heiwajima, Ohta-ku JP-143-0006 Tokyo Phone +81 3 5762 07 06 +81 3 5762 09 71 Fax e-mail helpdesk.ing.jp@mt.com

Malaysia

Mettler-Toledo (M) Sdn Bhd Bangunan Electroscon Holding Lot 8 Jalan Astaka U8/84 Seksyen U8, Bukit Jelutong MY-40150 Shah Alam Selangor Malaysia Phone +60 3 78 45 57 73 +60 3 78 45 87 73 Fax e-mail ahmad rashidi@mt.com

Mexico

Mettler-Toledo S.A. de C.V. Pino No. 350, Col. Sta. MA. Insurgentes, Col Atlampa MX-06450 México D.F. Phone +52 55 55 47 57 00 +52 55 55 41 22 28 Fax e-mail mario.roca@mt.com

Poland

Mettler-Toledo (Poland) Sp.z.o.o. ul. Poleczki 21 PL-02-822 Warszawa Phone +48 22 545 06 80 +48 22 545 06 88 Fax polska@mt.com e-mail

Russia

Mettler-Toledo Vostok ZAO Sretenskij Bulvar 6/1 Office 6 RU-101000 Moscow Phone +7 095 921 92 11 +7 095 921 63 53 Fax +7 095 921 78 68 inforus@mt.com e-mail

Singapore

Mettler-Toledo (S) Pte. Ltd. Block 28 Ayer Rajah Crescent #05-01 SG-139959 Singapore Phone +65 6890 00 11 +65 6890 00 12 Fax +65 6890 00 13 e-mail ashley.kong@mt.com

Slovakia

Mettler-Toledo s.r.o. Bulharska 61 SK-82104 Bratislava Phone +421 243 42 74 96 Fax +421 243 33 71 90 predaj@mt.com e-mail

Slovenia

Mettler-Toledo d.o.o. Peske 12 SI-1236 Trzin +386 1 530 80 50 Phone +386 1 562 17 89 Fax e-mail cipot@mtslo.mt.com racman@mettler-toledo.si

South Korea

Mettler-Toledo (Korea) Ltd. Yeil Building 1 & 2 F 124-5, YangJe-Dong SeCho-Ku KR-137-130 Seoul Phone +82 2 3498 3500 +82 2 3498 3556 Fax Sales_MTKR@mt.com e-mail

Spain

Mettler-Toledo S.A.E. C/ Miguel Hernández, 69-71 ES - 08908 L'Hospitalet de Llobregat (Barcelona) Phone +34 93 223 76 00 Fax +34 93 223 76 01 bcn.centralita@mt.com e-mail

Sweden

Mettler-Toledo AB Virkesvägen 10 Box 92161 SE-12008 Stockholm Phone +46 8 702 50 00 +46 8 642 45 62 Fax e-mail sales.mts@mt.com

Switzerland

Mettler-Toledo (Schweiz) AG Im Langacher Postfach CH - 8606 Greifensee Phone +41 44 944 45 45 +41 44 944 45 10 Fax info.ch@mt.com e-mail info.ola.ch@mt.com

Thailand

Mettler-Toledo (Thailand) Ltd. 272 Soi Soonvijai 4 Rama 9 Rd., Bangkapi Huay Kwang TH-10320 Bangkok +66 2 723 03 00 Phone +66 2 719 64 79 Fax e-mail mettler@samart.co.th

USA/Canada

Mettler-Toledo Ingold, Inc. 36 Middlesex Turnpike Bedford, MA 01730, USA Phone +1 781 301 8800 Toll free +1 800 352 8763 +1 781 271 0681 Fax mtprous@mt.com e-mail ingold@mt.com

Subject to technical changes. © Mettler-Toledo GmbH, Process Analytics 01/06 Printed in Switzerland. 52 403 126

Mettler-Toledo GmbH, Process Analytics Industrie Nord CH-8902 Urdorf Tel. + 41 44 729 62 11, Fax +41 44 729 66 36

www.mtpro.com