# Electronic data recording and protection according to FDA 21 CFR Part 11

(Additional function SW 700-107 for M 700 modular process analysis system)



Latest product information: www.mtpro.com

Order number: 52 121 241





#### Warranty

Defects occurring within 1 year from delivery date shall be remedied free of charge at our plant (carriage and insurance paid by sender). Sensors, fittings and accessories: 1 year ©2003 Subject to change without notice

#### **Return of products under waranty**

Please contact your local METTLER TOLEDO representative before returning a defective device. Ship the cleaned device to the address you have been given. If the device has been in contact with process fluids, it must be decontaminated/disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.



**Disposal** (Directive 2002/96/EC of January 27, 2003) Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

#### **Registered trademarks**

The following registered trademarks are used in this instruction manual without further marking

CalCheck Calimatic Sensocheck Sensoface ServiceScope VariPower

SMARTMEDIA® is a registered trademark of Toshiba Corp., Japan

 $\ensuremath{\mathsf{InPro}}^{\circledast}$  is a registered trademark of Mettler Toledo GmbH, Switzerland

Mettler-Toledo GmbH, Process Analytics, Industrie Nord, CH-8902 Urdorf, Tel. +41 (01) 736 22 11 Fax +41 (01) 736 26 36 Subject to technical changes. Mettler-Toledo GmbH, 12/03. Printed in Germany.

### AuditTrail and signature

Electronic data recording and protection according to FDA 21 CFR Part 11 Additional function SW 700-107

#### Provided by: TAN and AuditTrail card

#### Contents

Latest product information	4
How to order an additional function	5
AuditTrail	6
Signature	7
Activating the AuditTrail	8
Signature entry	9
AuditTrail: Parameter setting	10
AuditTrail: User management	11
AuditTrail: Recording	12
AuditTrail log	13

### Latest product information

Additional functions

The M 700 is an expandable modular process analysis system. For latest product information, please refer to:

#### www.mtpro.com

#### **Additional functions**

A current version of the respective function description is available as pdf for download.

### How to order an additional function

Device-specific TAN (transaction number)

Additional functions expand the device capabilities.

The additional functions are device-specific. When ordering an additional function, you therefore have to specify the serial number and hardware version of your M 700 FRONT in addition to the respective order number. The manufacturer then supplies a TAN (transaction number) to release the additional function.

Menu	Display	Serial number FRONT module	
	Image: Constraint of the selection	Menu selection Call up diagnostics. From the measuring mode: Press <b>menu</b> key to select menu. Select diagnostics using arrow keys, confirm with <b>enter</b> .	
diag	Image: Constraint of the second se	<b>Diagnostics</b> Select Device description using arrow keys, confirm with <b>enter</b> .	
	Image: Constraint of the second se	<b>Device description</b> Please specify the <u>serial number/hardware version</u> when ordering an additional function.	

#### Serial number of FRONT module

METTLER TOLEDO

Function description

## Electronic data recording and protection according to FDA 21 CFR Part 11

Pharmaceutical applications not only require a hygienic design, but also the protection and complete recording of data (FDA regulation: 21 CFR Part 11). The M 700 with stainless steel enclosure and SmartMedia card meets these demands. The FDA regulation basically consists of two parts: "Audit Trail" and "Electronic Signature". Both parts interact directly: For every change of parameters, the device records what has been changed and who has changed it.

Functions like calibration (adjustment) often may only be performed after written approval! This approval is given by an electronic signature.

#### AuditTrail card

The AuditTrail card is a specially prepared SmartMedia card. The M 700 records the AuditTrail log on an "AuditTrail card" only. Readability from PC is not affected.

#### Note:

Only specially coded SmartMedia cards supplied from Mettler Toledo can be used as AuditTrail cards.

As delivered, the AuditTrail card includes a PC program which makes the recorded AuditTrail data readable. When log coding has been activated, it also allows decoding of the data. The program supports the data export to other programs, e.g. Microsoft Excel.

### Signature

Function description

#### Signature

A signature consists of two entries, the user name and the corresponding passcode. For each change at the device (calibration/adjustment, parameter setting, maintenance, if required also diagnostics), the M 700 records in the "AuditTrail Log" and in addition in the logbook who has performed the change. For that purpose, access authorizations are introduced.

#### Access authorization

Every access authorization consists of a signature:

- User name and
- Passcode.

The administrator and the user enter the access data in the user management menu. The administrator enters the user name (short name) and the user adds his/her personal passcode.

Each time he/she calls up a protected menu, the user must enter his/her name and passcode. The M 700 only gives access to the function if the access codes are valid. The recordings (AuditTrail log, logbook) on the AuditTrail card include the user names so that all actions can be clearly traced back.

#### Coded storage

The entries in the AuditTrail log can be coded using an algorithm and protected with checksums so that a special evaluation program can prove the authenticity of the data records on a PC.

### Activating the AuditTrail

Select menu: Parameter setting/System control/Release of options Note: The TAN for releasing an additional function is only valid for the device with the corresponding serial number

Display 7.00 pH  $\overline{7}$ ш 25.6 ℃ Menu selection cal Select: Inter Return to meas 🛍 Lingua 25.6 ℃ m 7 0.003 mS/cm Parameter setting Uiewing level (All Data) view Operator level (Operation Data) opl Admin, level (All Data) adm Return 16.9 %Air Π A 闤 Ш 25.6 ℃ Parameter setting (Administrator) System control FRONT M 700-011 BASE M 700-021 🗅 🗉 Module O2 4700 🗅 🗉 Module PH 2700 🗀 🗉 Module Cond Ind 7700 Return 0.003 mS/cm 🛕 🗎 □ 25.6 °C Release of options (Administrator) 102 Parameter sets Active 103 Measurement recorder Thactive 104 Logbook Inactive 106 Software update Inactive 107 AuditTrail Active Inactive Abort

#### Activate additional function

#### Menu selection

Call up parameter setting. From the measuring mode: Press **menu** key to select menu. Select parameter setting using arrow keys, confirm with enter.

#### **Parameter setting**

Select Administrator level using arrow keys, confirm with enter. Enter passcode and confirm (Passcode as delivered: 1989).

Select System control using arrow keys, confirm with enter. Then select Release of options using arrow keys, confirm with enter.

#### **Release of options**

Select the additional function to be released.

Set option to "active". Enter the TAN at the prompt. (Note: The TAN is only valid for the device with the corresponding serial number, see previous page.) The option is available after the TAN has been entered.

Menu

(User name, passcode)

#### Signature entry (user name, passcode) to open a menu

To protect the data against unauthorized query or modification, access data must be entered before a main menu from the menu selection is opened. Two-stage security is provided by input of user name and passcode.



The passcode is not visible. It is only displayed as 4 dots. The administrator can delete each user's passcode (with entry in AuditTrail log). The user is identified by his/her access data.

After entry of the valid signature all menus are available with the defined user rights. Only in measuring mode, the rights will be erased again. This allows accessing several menus without having to enter the access data each time.

When an invalid combination of user name and passcode is entered, access is denied with the text "Wrong passcode". The entry "Wrong passcode" with the entered (wrong) user ID is recorded in the AuditTrail log.

#### **Default setting for signature**

Signature User: ADMIN Passcode: 1989

### AuditTrail: Parameter setting

Adapt function: Parameter setting/System control/AuditTrail

#### **Parameter setting**

In the System Control menu the Administrator can adjust the AuditTrail system to the requirements of his company. For example, he/she can select the protection against unauthorized query and the coding of the data in the AuditTrail log.

	□ 7.20 pH □ 25.6 ℃				
AuditTrail (Administrator)					
🗅 User management					
AutoComplete	On Off				
Access diag menu	▼ Free				
Record AuditTrail	On Off				
Coding Log	On Off				
Return					

- "AutoComplete" is an input assistance for the user name as known from some Windows programs.
- Access to the Diagnostics menu can be "free" or require input of a "Signature" according to the access authorization granted the User Management menu.
- For the recording on SmartMedia card, a coding can be activated so that the data cannot be read or queried without corresponding programs.

#### User management

An administrator and up to 5 users can be entered in the M 700. Each user has a user name (max. 8 characters) and a passcode (4 characters). Each user can be granted access and usage rights for each main menu level. The combination of user name and passcode is the user's "Electronic Signature". Every logged modification is signed and can be clearly traced back. Different access rights for each menu can be defined for each user.

### AuditTrail: User management

Create, edit, erase user data

Menu	Display	AuditTrail: User management
unter State	Image: Constraint of the selection	Menu selection From the measuring mode: Press menu key to select menu. Select parameter setting using arrow keys, confirm with enter. Signature for first access: User: ADMIN Passcode as delivered: 1989
	Image: Second system      Image: Second system        Image: Second	Select: Parameter setting/System control/ AuditTrail/User management. Access data and rights can be created for the administrator and up to 5 users. The administrator has always admin- istrator rights for parameter setting. The user enters his/her passcode after being asked to do so by the administrator. The passcode is only displayed as a series of 4 dots.

#### Access rights for the main menu levels:

- Calibration:
- [ No access | Check | Adjustment ]
- Maintenance:
- [No access | Maintenance ]
- Parameter setting:
- Diagnostics:
- [ No access | Viewing | Operator | Administrator ] [ No access | Diagnostics ]

### METTLER TOLEDO

### Record AuditTrail

Functions

#### AutoComplete

The "AutoComplete" function facilitates the input of the user name. You only have to enter as many letters as the M 700 requires to identify the name.

This function can be switched off in the AuditTrail menu because the automatic display of the user name reduces system security, which then might not be sufficient for FDA documentation.

#### Access to diag menu

To protect all data against unauthorized query, it is also possible to protect the diagnostics level with a passcode (this will also affect the Favorites menu). As delivered, access to the diagnostics level is free.

#### Record AuditTrail

The AuditTrail log is an additional - protected if required - recording of all data on the AuditTrail card (SmartMedia card). The log is recorded in the currently selected language.

The data can be coded (to protect against unauthorized query) and secured with a checksum (to prove the authenticity of the data). The coding can be switched on or off in the AuditTrail menu of the Administrator level.

A four-digit line numbering (0000 ... 9999 P 0000 ..., cyclic) ensures uninterrupted recording. The AuditTrail log cannot be deleted by the M 700. A buffer storage is provided to store the recorded data while the AuditTrail card is being replaced, for example. An error message "Fail AuditTrail card" is generated when the buffer is full. The oldest entries in the buffer will be recorded. They are then missing in the log (the corresponding line numbers will be missing).

When an empty AuditTrail card is inserted, the buffer storage and the current entries are transferred to the empty card.

### Record AuditTrail

Parameters

The AuditTrail log cannot be queried by the M 700. It can only be evaluated using the suitable PC program. According to FDA 21 CFR Part 11 the user must ensure the transport of the data from the M 700 to evaluation and further processing or archiving. Furthermore, the user must ensure that the log data cannot be deleted from the memory card. The PC evaluation program is able to decode the data and verify the security parameters and thus to prove the authenticity of the data on the card. It allows printing the encoded data and exporting them to other PC applications.

#### Logbook

The logbook is provided in addition. In AuditTrail systems it runs simultaneously with the audit trail log in contrast to which it can be queried provided that the diagnostics menu is accessible. The administrator is permitted to delete the logbook data since all security-relevant data are stored in the AuditTrail log.

#### **Recordings in the AuditTrail log**

#### **Caution!**

#### Data loss (even total destruction) if the card is not closed before it is pulled out. (Maintenance menu, "Close card")

The following entries and messages are recorded:

- Measuring point, device description, module equipment and serial numbers
- Menu system log-ins
- Movements within menus (menu headlines)
- Modified parameters at the press of **enter**
- Calibration record at the end of a calibration or adjustment incl. user name, sensor type and serial number
- Messages such as failure and maintenance request

Audi Measu FRONT BASE M pH 270 O2 470	tTrail log rement poir M 700-011 M 700-021 00 00	nt I	BIGBLUE vessel 0006123 0006458 0007221 0006045	METTLER	TOLEDO
6653 6656 6657 6658 6659 6660 6661 6662 6663 6664 6665 6666 6667 6668 6669 6670 6671 6672 6673 6674 6675 6674 6675 6676 6677 6678 6679	07/01/04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 07.01.04 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	10:43:02Login: A 10:43:02Login: > 10:43:02Login: > 10:43:02Login: < 10:43:02Login: < 10:43:02Login: > 10:43:02Login: > 10:43:02Login: > 10:43:02Login: < 10:43:02Login: < 10:43:02Login: < 10:43:02Login: <	Admin (ADMIN) Parameter setting (A [1] Module pH 2700 Parameter setting (A Menu selection Admin (ADMIN) Calibration [1] Module pH 2700 [1] Calimatic Warn Identical Buf [1] Calimatic [1] Calimatic [1] Calimatic [1] Calimatic [1] Calimatic [1] Calimatic [2] Cal record Last calibration: 20. User: ADMIN Sensor type: El 204 Serial number: 0819 Cal mode: Calimatic Zero: 7.02 pH Slope: 58.0 mV/pH Impedance glass (29 Impedance ref (25%) 1st nominal buffer v Electrode potential: Cal temperature: 29 Response time: 19 st [1] Module pH 2700 Calibration	Administrator) (Administrator) Administrator) 11.03 10:44 5 5 5°C): 825.9 MOh C): 119.4 kOhm value: 7.00 pH 1 mV 5.0 °C	ок ок ок ок ок ок ок ок ок ок ок ок ок о
		2			

Typical AuditTrail log