
# Indholdsfortegnelse

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## 5 Operating the Instrument

The SevenExcellence™ display is equipped with a protective foil. Removing this foil increases the sensitivity of the touch screen.

SevenExcellence is delivered with default parameters for calibration and measurement as well as for methods to ensure an immediate deployment of the instrument after installation.

This quick guide leads you through the first steps of using SevenExcellence.

For a detailed guidance in using the extensive features of SevenExcellence the user is referred to the operating instructions.

### 5.1 Direct calibration

- 1 Select a suitable sensor for the calibration and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**], enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap the colored tag of the module you use for calibration, **Module settings** is displayed.
- 5 Tap the list field **Sensor name** and select the sensor you use, if necessary select the external temperature probe.
- 6 Adjust the **General settings** and **Calibration settings** if necessary.
- 7 Immerse the electrode/sensor in the sample, tap [**Calibrate**].

### 5.2 Direct measurement

- 1 Select a suitable sensor for the measurement and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**], enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap the colored tag of the module you use for measurement, **Module settings** is displayed.
- 5 Tap the list field **Sensor name** and select the sensor you use, if necessary select the external temperature probe.
- 6 Adjust **General settings** and **Measurement settings** if necessary.
- 7 Tap [**Save**].
- 8 Immerse the electrode/sensor in the sample, tap [**Read**].

### 5.3 Starting a method

- 1 Select a suitable sensor for the method and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**] enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap [**Methods**], select a method or tap [**New**] to use a **Method type** template. Use the default parameters or modify them according to your requirements.
- 5 Tap [**Save**].
- 6 Immerse the electrode/sensor in the sample, tap [**Start**].

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## Quick Check

| Content of delivery                                              |                                     | Active | Inactive |
|------------------------------------------------------------------|-------------------------------------|--------|----------|
| SevenExcellence™ instrument                                      |                                     |        |          |
| Fixing screw for the electrode arm                               |                                     |        |          |
| Country-specific power supply                                    |                                     |        |          |
| Semi-transparent module cover                                    |                                     |        |          |
| Measurement module(s) (depending on kit configuration)           | pH/Ion                              |        |          |
|                                                                  | Conductivity                        |        |          |
|                                                                  | pH/mV                               |        |          |
|                                                                  | DO/BOD                              |        |          |
| Sensor(s) (depending on kit configuration)                       | All sensors delivered               |        |          |
|                                                                  | All cables delivered                |        |          |
| uPlace™ electrode arm                                            | Base with plastic covers (inserted) |        |          |
|                                                                  | Movable electrode holder            |        |          |
|                                                                  | Extension shaft part                |        |          |
|                                                                  | Wrench                              |        |          |
|                                                                  | Two screws                          |        |          |
| Installation guide (print version, all languages)                |                                     |        |          |
| Declaration of conformity                                        |                                     |        |          |
| Calibration certificate                                          |                                     |        |          |
| Location of the instrument                                       |                                     | Active | Inactive |
| Device is not exposed to powerful vibrations                     |                                     |        |          |
| Device is not exposed to direct sunlight                         |                                     |        |          |
| Device is not exposed to atmospheric humidity greater than 80%   |                                     |        |          |
| Device is not exposed to corrosive gas atmosphere                |                                     |        |          |
| Device is not exposed to temperatures below 5 °C and above 40 °C |                                     |        |          |
| Device is not exposed to powerful electric or magnetic fields    |                                     |        |          |
| Power supply                                                     |                                     | Active | Inactive |
| Power supply is connected                                        |                                     |        |          |
| Main power is turned on                                          |                                     |        |          |
| Instrument is switched on                                        |                                     |        |          |
| Display shows the home screen                                    |                                     |        |          |
| Modules                                                          |                                     | Active | Inactive |
| Module(s) is/are attached                                        |                                     |        |          |
| Module(s) is/are displayed on the home screen                    |                                     |        |          |
| Check box(es) of the module(s) is/are activated                  |                                     |        |          |
| Semi-transparent module cover is attached                        |                                     |        |          |
| Sensors                                                          |                                     | Active | Inactive |
| Sensor(s) is/are attached                                        |                                     |        |          |
| Mini-LTW socket fits precisely (digital ISM® sensors)            |                                     |        |          |
| BNC socket fits precisely (pH, ORP and ISE electrodes)           |                                     |        |          |
| RCA socket fits precisely (temperature probes)                   |                                     |        |          |
| Mini-DIN socket fits precisely (conductivity sensors)            |                                     |        |          |
| Recognition of peripheral devices                                |                                     | Active | Inactive |
| ISM electrode recognized                                         |                                     |        |          |
| Stirrer recognized                                               |                                     |        |          |
| Printer recognized (only US-P25)                                 |                                     |        |          |
| PC recognized with <b>LabX</b> / <b>EasyDirect pH</b>            |                                     |        |          |
| USB stick recognized                                             |                                     |        |          |
| Rondolino / InMotion recognized                                  |                                     |        |          |
| uPlace™ electrode arm                                            |                                     | Active | Inactive |
| Extension shaft part used                                        |                                     |        | *        |
| Mounted as stand alone                                           |                                     |        | *        |

|                                                                                                                                                                                                                                                                                                          |       |  |   |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|--|---|
| Attached to the instrument                                                                                                                                                                                                                                                                               | Right |  | * |
|                                                                                                                                                                                                                                                                                                          | Left  |  | * |
| Fitted correctly and stable                                                                                                                                                                                                                                                                              |       |  |   |
| <p><b>i Bemærk</b></p> <ul style="list-style-type: none"> <li>If you have any questions answered in the negative, read the corresponding chapter of the installation instructions. If the problem persists, contact your vendor.</li> </ul> <p>* Irrelevant for the functionality of the instrument.</p> |       |  |   |

|                                             |           |
|---------------------------------------------|-----------|
| Instrument:                                 |           |
| Serial Number:                              |           |
| Date:                                       | Date:     |
| Company:                                    | Company:  |
| Name:                                       | Name:     |
| Surname:                                    | Surname:  |
| Signature                                   | Signature |
| <p><b>Retain on file for reference.</b></p> |           |

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► [www.mt.com/service](http://www.mt.com/service)

**www.mt.com**

For mere information

**Mettler-Toledo GmbH**

Im Langacher 44  
8606 Greifensee, Switzerland  
Tel. +41 22 567 53 22  
Fax +41 22 567 53 23  
[www.mt.com](http://www.mt.com)

Forbehold for tekniske ændringer.  
© 01/2025 METTLER TOLEDO. Alle rettigheder forbeholdes.  
30996963I

30996963