## Operating Instructions for Dissolved Oxygen Sensors LE 611/612 LE 621/622

### Unpacking

Carefully remove the protective bottle from the tip of the sensor by unscrewing the lid then remove the bottle. Remove the white shorting plug from the BNC connector and store in a safe place.

### Calibration

Connect the sensor to the instrument, rinse with deionized water and carefully dry the tip of the sensor with a tissue. Then follow your instrument instruction manual for detailed calibration instructions.

It is recommended to calibrate the sensors daily.

### Cleaning and Maintenance

This galvanic oxygen sensor has no customer serviceable parts. Cleaning is important. The sensor should be rinsed thoroughly with deionized water before, between and after measurements are made. Never allow that sample dry on the electrode. If dirt or sample deposits on the membrane, the tip of the sensor can be immersed in deionised water and ultrasonically cleaned for 30 seconds.

### Storage

For short term storage: rinse the electrode in deionised water and refit the protective bottle filled with storage solution.

For longer term storage: rinse the sensor with deionised water, refit the storage bottle filled with storage solution and the shorting plug, and store in a cool dry place.

Storage solution: 10% NaCl

### Specification

Type designation: galvanic dissolved oxygen sensor

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