
METTLER TOLEDO

DS500

Surfactant sensitive electrode

Tensidsensitive Elektrode

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General information

The DS500 surfactant sensitive electrode is suitable for the indication of titrations in aqueous systems of a cationic surfactant with a known anionic surfactant, or of an anionic surfactant with a known cationic surfactant. Use of the DS500 is described in the application brochure 22.

Specifications

Type of ion selective membrane: Polymer matrix-liquid membrane with ion carrier molecules

Resistance of the membrane: <2 MΩ

Measurement range: <10⁻⁵ mol/L to critical micelle concentration (cmc)

Measurement temperature range: 0...50 °C

Optimum pH range: 2...10, depending on the surfactant measured

Interferences: In surfactant mixtures the most sorbing surfactant is the dominant species. Lipophilic inorganic and organic ions can interfere (perchlorate, periodate and tetraphenylborate ions, cationic and anionic dyes).

Precision

In the absence of interfering substances, a precision of better than ±0.5 mV can be achieved.

Response times

When changing from lower to higher surfactant concentrations, the response time is approximately one minute; on change from higher to lower concentrations several minutes.

Materials

Electrode body of POM copolymer, membrane module of PVC with plasticizer and ion carrier substances, seals of silicon.

Electrolyte for DX200 reference electrode

3 M KCl (for anionic and cationic surfactants).

Disposal

Dispose of the metallic, plastic and glass parts separately.

Maintenance and service

Deposits on the surface of the membrane can be removed by thoroughly rinsing with deionized water. After rinsing with deionized water, conditioning in a dilute NaCl solution during approx. 10 minutes is required. For short breaks in between measurements, the electrode should be kept in deionized water. For longer breaks, e.g. over night, it is recommended to place the electrode in deionized water.

Storage

The DS500 electrode can be stored in air. It is best stored disassembled (dry) for longer periods. Before titration the electrode is conditioned in NaCl (0.01 mol/L) for at least 10 minutes. Eventually shake the electrode once like a clinical thermometer. Reference electrode DX200: allow electrolyte to run out. For storage, close filling port with the rubber cap.

Changing the membrane module

Membrane modules with a damaged membrane, excessive response times or unstable potentials (high membrane resistance) which can no longer be regenerated must be changed.

Procedure for changing a membrane module:

- Place electrode upside down and unscrew old membrane module from the electrode body.
- Refill electrode body with internal solution up to the edge ensuring absence of bubbles.
- Also fill the new membrane module with internal solution up to the sealing ring.
- Screw the membrane module on tightly.
- Give the electrode a single powerful shake (like a clinical thermometer).
- Rinse the entire electrode thoroughly with deionized water. Then condition for at least 10 minutes in NaCl (0.01 mol/L).

Accessories

Reference electrode	DX200
Cable S7-Lemo (DS500)	89601
Membrane module replacement	51107669
Electrolyte	51107899

Literature

Application brochure	51725015
Guide to ION Selective Measurement	51300075

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