
DispenSix

English	User Manual Six-Port Burette
Deutsch	Benutzerhandbuch Sechs-Port-Bürette
Español	Manual de usuario Bureta de seis puertos
Français	Guide de l'utilisateur Burette à six voies
Português	Manual do usuário Bureta de seis portas
Italiano	Manuale per l'utente Buretta a sei vie
Polski	Podręcznik użytkownika Biureta sześcioportowa
简体中文	用户手册 六通道滴定管

METTLER TOLEDO

User Manual **Six-Port Burette**

English

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1 Introduction

Thank you for choosing a METTLER TOLEDO DispensSix six-port burette.

This document provides you with the information you need to get started with your burette.

Read the Reference Manual of the main instrument and the User Manual of the burette drive for a full description of the interface between the burette, the burette drive, and the main instrument. The manuals are available online. See [Download manuals ▶ Page 7].

If you have any additional questions, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

1.1 Explanation of conventions and symbols

Refers to an external document.

Elements of instructions

Instructions always contain action steps and can contain prerequisites, intermediate results and results. If an instruction contains more than one action step, the action steps are numbered.

- Prerequisites that must be fulfilled before the individual action steps can be executed.

1 Action step 1

⇒ Intermediate result

2 Action step 2

⇒ Result

1.2 Compliance information

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ www.mt.com/ComplianceSearch

2 Safety information

The burette is designed to be used as an accessory. The burette is compatible with all main instruments that list it as an accessory in their Reference Manual.

- This User Manual contains the information you need to set up the accessory, put it into operation and maintain it.
- Keep the User Manual for future reference.
- Include the User Manual if you transfer the accessory to other parties.
- The Reference Manual of the main instrument contains additional information on the interface between the accessory and the main instrument.

Only use the accessory according to this User Manual and the Reference Manual of the main instrument. If you do not use the accessory according to both documents or if it is modified, the safety of the main instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

The User Manual and Reference Manual are available online. See [Download manuals ▶ Page 7].

2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

NOTICE

A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

Notice

2.2 Product specific safety notes

Intended use

The burette is intended to be used by trained staff. The burette, together with a burette drive, is designed to aspirate and dispense liquids. The liquids must have the following properties:

- Viscosity: <10 mPa*s at ambient temperature
- No highly concentrated crystallizing reagents
- No particles
- No degassing when underpressure is applied
- Compatible with the materials the liquids come into contact with

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the accessory and who uses the accessory and the main instrument or authorizes any person to use them, or the person who is deemed by law to be the operator of the accessory and the main instrument. The instrument owner is responsible for the safety of all users of the accessory and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the accessory and the main instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

NOTICE

Damage to the accessory or the main instrument due to the use of unsuitable parts

Using unsuitable parts with the accessory can damage the accessory or the main instrument or cause them to malfunction.


- Only use parts from METTLER TOLEDO that are intended to be used with your accessory.

See also

- [Technical data](#) ▶ Page 12
- [Burette overview](#) ▶ Page 5
- [Six-port valve](#) ▶ Page 5

3 Design and function

3.1 Burette overview

No.	Name	Function
1	Ports	<ul style="list-style-type: none"> • Openings to dispense and aspirate liquids
2	Valve head	<ul style="list-style-type: none"> • Part of the six-port valve
3	Housing	<ul style="list-style-type: none"> • Attaches the burette to the burette drive
4	Rotary encoder	<ul style="list-style-type: none"> • Protection of glass cylinder • The position determines which port is open
5	Port number	<ul style="list-style-type: none"> • Shows which port is open
6	Gear	<ul style="list-style-type: none"> • Connects the rotary encoder to the burette drive • Enables the burette drive to open and close ports
7	RFID chip	<ul style="list-style-type: none"> • Stores burette data
8	Piston	<ul style="list-style-type: none"> • Connects the burette to the piston rod on the burette drive • Enables the burette drive to fill and empty the burette
9	Glass cylinder	<ul style="list-style-type: none"> • Holds the liquid
10	Cylinder adapter	<ul style="list-style-type: none"> • Part of the six-port valve


See also

Technical data ▶ Page 12


3.2 Six-port valve

The cylinder adapter and the valve head form the six-port valve.

- The lower surface of the valve head (1) and the upper surface of the cylinder adapter (2) are pressed together and form a seal.

- From the inlet port (6), liquid passes through the central opening (5) in the valve head (1). The liquid flows down into the groove (4) in the cylinder adapter, through the opening (2) and into the glass cylinder (3).


Open and closed ports

Only one port is open at a time. The orientation of the groove on the cylinder adapter determines, which port is open. To open a port, the rotary encoder and the cylinder adapter are rotated, aligning the groove with the port.

Port 1 open



Port 3 open

4 Installation

4.1 Scope of delivery

Part	Order number
 Liquid Handler DispenseSix Six-port burette	–
 Flat seal GL45	30673280
 Burette parking station M6 GL45	–
 Stopper NS14	–
 Drying tube NS14	30673119
 Tube tip holder	–

Part		Order number
	Blind plug set M6 (6 pcs)	30869353
	Suction tube M6 87 cm (2 pcs) 87 cm, two M6 connectors	30673121
	Dispensing tube M6 100 cm 100 cm, M6 connector and tapered connector	30673120
	Tube clip set	–
	Glass bottle Brown, 1 L	–
	Cable M8/F, RJ12/M 60 cm	30634408
	User Manual	–

See also

 Accessories, spare parts and consumables ▶ Page 13

4.2 Download manuals

- 1 Go to the website www.mt.com/library.
- 2 Select the **Technical Documentation** tab.
- 3 Find the product type on the housing of the main instrument and enter it into the search field.
- 4 Start the search.
- 5 Select the manual from the result list.
- 6 Select the link.
 - ⇒ The manual is either opened or downloaded depending on the browser settings.
- 7 Check which software version is installed on your main instrument.
- 8 If the manual is not written for the installed software version, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

4.3 Unpack the burette

- 1 Remove the burette from the protective packaging.
- 2 Store the protective packaging for later transport over long distances.
- 3 Check that you have received all parts listed in the scope of delivery.
- 4 Inspect the parts visually for flaws or damage.

- 5 If parts are missing or damaged, report it to your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

See also

🔗 Burette overview ► Page 5

5 Operation

Read the Reference Manual of the main instrument for more information about how to use the burette.

See [Download manuals ► Page 7].

6 Maintenance

In this chapter, you find descriptions of the maintenance tasks you can perform yourself. Any other maintenance tasks need to be performed by a service technician who has been qualified by METTLER TOLEDO. If you experience problems with your burette, contact your authorized METTLER TOLEDO dealer or service representative.

METTLER TOLEDO recommends that a preventive maintenance and calibration certification is done at least once a year through your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

6.1 Clean the burette

NOTICE


Damage to the burette due to inappropriate cleaning methods

The burette leaks if the sealing surfaces or the sealing lips are scratched.

- 1 Make sure you only use soft cleaning materials.
- 2 Make sure the cleaning agent is compatible with the material of the part you want to clean.


Sealing surfaces

- Lower surface of the valve head (1)
- Upper surface of the cylinder adapter (2)

Sealing lips

- The five sealing lips (1) of the piston

NOTICE

Damage to the glass cylinder due to improper handling

The glass cylinder can break if it slides out of the burette and drops.

- 1 Do not tilt the burette or turn it upside down when the valve head is not installed.
- 2 Hold the glass cylinder in place with a finger on the cylinder adapter when the valve head is not installed.
- 3 Check that the valve head is installed correctly before you fill the burette or turn it upside down.

NOTICE



Incorrect data on the RFID chip due to incorrect assembly

The serial number of the glass cylinder is stored on the RFID chip. If a glass cylinder with a different serial number is installed, the data on the RFID chip is no longer correct.


- If you have more than one burette, make sure that the correct glass cylinder is installed after cleaning.

6.1.1 Disassemble the burette

- The burette is removed from the burette drive.
 - No blind plugs are connected to the valve head.
 - No tubes are connected to the valve head.
 - The burette is empty.
- 1 Hold the housing (2) firmly with one hand and grip the sides (1) of the valve head.
 - 2 **NOTICE: Damage to the valve head. Pressure on the port area in the center damages the valve head. Only exert force on the sides and top outer ring of the valve head.**
Press the valve head onto the housing, turn it counter-clockwise and remove it.
 - 3 Lift the rotary encoder (3) out of the housing (2) without tilting it.
- 4 Lift the assembled cylinder (1) out of the rotary encoder (2) and place it on the work surface.
 - 5 To remove the centering ring (3), tilt the rotary encoder (2) until the centering ring (3) slides out.

- 6 Firmly hold the glass cylinder (2) with one hand and the edge of the cylinder adapter (1) with the other hand.
- 7 Twist the cylinder adapter (1) while pulling it out along the central axis. You might need to exert some force because the parts fit together tightly.
- 8 Place the cylinder adapter right-side-up (3) on the work surface.
- 9 Hold the lower end of the piston (4) and pull the piston out along the central axis. You might need to exert some force because the parts fit together tightly.

See also

 Burette overview ▶ Page 5

6.1.2 Clean the parts

Cleaning agents

METTLER TOLEDO recommends the following cleaning agents:

- Isopropanol
- Ethanol
- Water with a mild detergent

If you have questions about the compatibility of cleaning agents, contact your authorized METTLER TOLEDO service representative or dealer.

 www.mt.com/contact

Clean by hand

- 1 Wipe the burette parts with a soft cloth moistened with the cleaning agent.
- 2 Rinse the parts with deionized water.
- 3 Air-dry the parts or dry them with a soft tissue.

Clean in the dishwasher

- 1 Place the burette parts in the dishwasher.
- 2 Check that the sealing surfaces and the sealing lips cannot be scratched during the wash cycle.
- 3 Run the dishwasher.

See also


 Six-port valve ▶ Page 5

 Technical data ▶ Page 12

6.1.3 Assemble the burette


Assemble the cylinder

- 1 Align the piston (2) along the central axis (1) of the glass cylinder.
- 2 Push the piston into the glass cylinder until the lower edge of the piston and the lower edge of the glass cylinder are aligned (3). You might need to exert some force because the parts fit together tightly.
- 3 Align the cylinder adapter (4) along the central axis of the glass cylinder and push it into the glass cylinder (5).


Assemble the rotary encoder

- 1 Insert the centering ring (1) into the rotary encoder (2).
- 2 Make sure the centering ring (6) lies flat on the bottom of the rotary encoder.
- 3 Align the narrow alignment gap (4) of the cylinder adapter with the narrow alignment tooth (5) of the rotary encoder.
- 4 Insert the assembled cylinder (3) into the rotary encoder.


Install rotary encoder and valve head

- 1 Slide the rotary encoder (1) into the housing (2).
- 2 If the top of the rotary encoder (3) juts out of the housing, rotate the rotary encoder until it slides down into place.

- 3 Align the alignment tooth (3) with the alignment gap (4) and place the valve head on the housing.
- 4 Hold the housing firmly with one hand and grip the sides (5) of the valve head.
- 5 **NOTICE: Damage to the valve head. Pressure on the port area (1) in the center damages the valve head. Only exert force on the sides (5) and top outer ring (2) of the valve head.**

Press the valve head onto the housing and turn it clockwise until it locks into place (6). You might need to exert some force because the parts fit together tightly.

See also

[Burette overview](#) ▶ Page 5

6.2 Prepare the burette for storage or shipping

- 1 Clean the burette.
- 2 To store the burette, keep it in a dry and clean place.
- 3 To ship the burette, use the original packaging.

See also

[Technical data](#) ▶ Page 12

6.3 Dispose of the burette

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

7 Technical data

7.1 Burette

Burette

Characteristic		Value
Dimensions	Width	56 mm
	Depth	59 mm
	Height	139 mm
Weight		278 g
Volume range		0.010...50 mL
Materials	Housing	PET (polyethylene terephthalate) and PBT (polybutylene terephthalate)
	Rotary encoder	PET (polyethylene terephthalate)
	Centering ring	PE-HD (high-density polyethylene)
	Valve head	PVDF (polyvinylidene difluoride) PET (polyethylene terephthalate) and PBT (polybutylene terephthalate)
	Cylinder adapter	PTFE (polytetrafluoroethylen) with graphite EPDM (ethylene propylene diene monomer (M-class) rubber)
	Glass cylinder	Borosilicate glass
	Piston	PTFE (polytetrafluoroethylen) Stainless steel (1.4301)

Site requirements

Characteristic		Value
Ambient conditions	Ambient temperature	5...40 °C
	Recommended operational temperature ¹⁾	18...28 °C
	Relative humidity	Non-condensing, max. 80 % for temperatures up to 31 °C, decreasing linearly to 50 % at 40 °C
	Altitude	Up to 5000 m above sea level
	Use	In interior spaces
	Pollution degree	2
Storage conditions	Temperature	-20...+70 °C
	Relative humidity	0...95 %

¹⁾ METTLER TOLEDO produces and tests the equipment with test tools that are certified for this temperature range. Usage outside the given range could lead to inferior performance such as burette leakage.

7.2 Burette parking station

Characteristic	Value	
Dimensions	Width	68 mm
	Depth	118 mm
	Height	47 mm
Screw thread	Tube thread type	M6
	Bottle thread type	GL45
Materials	Burette parking station	PP (polypropylene)
	Drying tube	PP (polypropylene)
	Flat seal	EPDM (ethylene propylene diene monomer)
	Tube tip holder	PP (polypropylene)
	Stopper NS14	PP (polypropylene)

7.3 Tubes and blind plug

Characteristic	Value	
Dispensing tube	Outer tube	PE-LD (low-density polyethylene)
	Inner tube	PFA (perfluoroalkoxy alkanes)
	Connector	PVDF (polyvinylidene difluoride)
	Thread type	M6
Suction tube	Inner tube	PFA (perfluoroalkoxy alkanes)
	Outer tube	PE-LD (low-density polyethylene)
	Connector	PVDF (polyvinylidene difluoride)
	Thread type	M6
Blind plug	Material	ETFE (ethylene tetrafluoroethylene)
	Thread type	M6


8 Accessories, spare parts and consumables

All accessories, spare parts and consumables are specified with their order number.

If you have any questions, contact your authorized METTLER TOLEDO service representative or dealer.


► www.mt.com/contact

Tubes


Dispensing tube InMotion M6 200 cm
200 cm, M6 connector and M8 connector

30673123


Dispensing tube M6 200 cm
200 cm, M6 connector and tapered connector

30673125


Dispensing tube M6 100 cm
100 cm, M6 connector and tapered connector

30673120


Suction tube M6 200 cm
200 cm, two M6 connectors

30673124


Suction tube M6 87 cm
87 cm, two M6 connectors

30673121

Tube syphon tip set (5 pcs)


30869304

Tube ferrules (5 pcs)


30673126

Burette


DispenSix valve head

30673133


DispenSix rotary encoder

30673128


DispenSix housing


30673127


Burette cylinder adapter 50 mL

30673131


Burette glass cylinder 50 mL

30673130


Piston 50 mL

30673132


Burette centering ring 50 mL

30673129**Miscellaneous parts**


Blind plug set M6 (6 pcs)

30869353

Drying tube NS14

30673119

Bottle adapter for suction tube M6 GL45

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© 02-2024 METTLER TOLEDO. All rights reserved.