
dDrive

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Português

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Table of Contents

1	Introduction	3
	1.1 Explanation of conventions and symbols	3
	1.2 Compliance information	3
2	Safety information	4
	2.1 Definition of signal words and warning symbols	4
	2.2 Product specific safety notes	5
3	Design and function	6
	3.1 Burette drive overview	6
	3.2 Overview of functions	6
4	Installation	7
	4.1 Scope of delivery	7
	4.2 Download manuals	7
	4.3 Unpack the burette drive	8
	4.4 Position the burette drive	8
	4.5 Install the burette	9
	4.6 Remove the burette	9
5	Operation	9
6	Maintenance	10
	6.1 Maintenance schedule	10
	6.2 Clean the burette drive	10
	6.3 Prepare the burette drive for storage or shipping	10
	6.4 Dispose of the burette drive	11
7	Technical data	11
8	Accessories, spare parts and consumables	12

1 Introduction

Thank you for choosing a METTLER TOLEDO dDrive burette drive. The burette drive can be used with the following burettes:

- DispenSix six-port burette
- 2 mL two-port burette
- 5 mL two-port burette
- 10 mL two-port burette

This document provides you with the information you need to get started with your burette drive.

Read the User Manuals for the six-port burette and the two-port burettes, and the Reference Manual for the main instrument for more information on how to install and operate the burette drive. See [Download manuals ▶ Page 7].

For third party licenses and open source attribution files, see the following link:

▶ www.mt.com/licenses

If you have any additional questions, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

1.1 Explanation of conventions and symbols

Refers to an external document.

Elements of instructions

Instructions always contain action steps and can contain prerequisites, intermediate results and results. If an instruction contains more than one action step, the action steps are numbered.

- Prerequisites that must be fulfilled before the individual action steps can be executed.
 - 1 Action step 1
 - ➔ Intermediate result
 - 2 Action step 2
 - ➔ Result

1.2 Compliance information

The instrument complies with the directives and standards listed on the Declaration of Conformity.

▶ <https://www.mt.com/doc>

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ www.mt.com/ComplianceSearch

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

▶ www.mt.com/contact

European Union

This product may contain SVHC candidate substances according to Article 33 of the EU regulation no. 1907/2006 (REACH). SVHC candidate substances are listed on the Declaration of Conformity (DoC).

▶ <https://www.mt.com/doc>

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual,

may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

2 Safety information

The burette drive is designed to be used as an accessory. The burette drive is compatible with all main instruments that list it as an accessory in their Reference Manual.

- This User Manual contains the information you need to set up the accessory, put it into operation and maintain it.
- Keep the User Manual for future reference.
- Include the User Manual if you transfer the accessory to other parties.
- The Reference Manual of the main instrument contains additional information on the interface between the accessory and the main instrument.

Only use the accessory according to this User Manual and the Reference Manual of the main instrument. If you do not use the accessory according to both documents or if it is modified, the safety of the main instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

The User Manual and Reference Manual are available online. See [Download manuals ▶ Page 7].

2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard: read the User Manual or the Reference Manual for information about the hazards and the resulting measures.

Notice

2.2 Product specific safety notes

Intended use

The burette drive is intended to be used by trained staff. Together with a burette, the burette drive is designed to aspirate and dispense liquids. The liquids must have the following properties:

- Viscosity dependent on burette:
 - Six-port burette: <10 mPa*s at ambient temperature
 - Two-port burettes with the following volumes 2 mL, 5 mL, and 10 mL: <100 mPa*s at ambient temperature
- No highly concentrated crystallizing reagents
- No particles
- No degassing when underpressure is applied
- Compatible with the materials the liquids come into contact with

The burette drive needs to be attached to the main instrument or another accessory. If not properly attached, the burette drive can fall over. As a result, the burette drive and the installed burette can be damaged.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the accessory and who uses the accessory and the main instrument or authorizes any person to use them, or the person who is deemed by law to be the operator of the accessory and the main instrument. The instrument owner is responsible for the safety of all users of the accessory and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the accessory and the main instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

NOTICE

Leakage at the valve head due to overpressure

If overpressure in the burette exceeds 4 bar, the content leaks out along the valve head.

- 1 When using a six-port burette, only use liquids with a viscosity <10 mPa*s at ambient temperature.
- 2 When using a 2 mL, 5 mL, or 10 mL two-port burette, only use liquids with a viscosity <100 mPa*s at ambient temperature.
- 3 Make sure that all ports are either connected to a tube or closed by a blind plug.

NOTICE

Damage to the accessory or the main instrument due to the use of unsuitable parts

Using unsuitable parts with the accessory can damage the accessory or the main instrument or cause them to malfunction.


- Only use parts from METTLER TOLEDO that are intended to be used with your accessory.

See also

-  Technical data ▶ Page 11
-  Burette drive overview ▶ Page 6
-  Overview of functions ▶ Page 6

3 Design and function

3.1 Burette drive overview

No.	Name	Function
1	Lock wheel	<ul style="list-style-type: none"> Locks the burette in place
2	Right side internal magnets	<ul style="list-style-type: none"> Attaches the burette drive to the main instrument or another accessory
3	Gear	<ul style="list-style-type: none"> Rotates the rotary encoder to open and close ports
4	Piston rod	<ul style="list-style-type: none"> Connects the piston of the burette to the motor in the burette drive
5	Left side internal magnets	<ul style="list-style-type: none"> Attaches the burette drive to another accessory
6	IN socket	<ul style="list-style-type: none"> Actor Bus IN socket to connect the main instrument or another burette drive
7	OUT socket	<ul style="list-style-type: none"> Actor Bus OUT socket to connect another burette drive
8	Socket cover	<ul style="list-style-type: none"> Protects the Actor Bus OUT socket from corrosion when not in use.
9	Signal LED	<ul style="list-style-type: none"> Green: burette drive connected to the main instrument and ready for use LED off: burette drive not connected to the main instrument

See also

Technical data ▶ Page 11

3.2 Overview of functions

Read the User Manuals of the six-port burette and the two-port burettes for more information on the design and function of the burettes. See [Download manuals ▶ Page 7].


Read the Reference Manual of the main instrument for more information on the programming of aspiration and dispensing. See [Download manuals ▶ Page 7].

The function of the burette drive depends on the installed burette.

- Burette drive and two-port burette (1) are designed to function as a dispensing unit.
- Burette drive and six-port burette (2) are designed to function as a liquid handler.

The burette drive is connected to a main instrument. The main instrument controls movement of the piston in the glass cylinder and the rotation of the rotary encoder.

- The motor in the burette drive moves the piston with a resolution of 120000 steps.

- The position of the rotary encoder determines, which port is open.

Aspiration

- The burette drive pulls the piston down.
- Underpressure builds in the cylinder.
- Liquid is aspirated into the cylinder through the open port.

The effect of underpressure on the liquids limits the types of liquids that can be used.

- No liquids with dissolved gases: when the pressure on the liquid is lowered, bubbles form in the liquid, which reduces the accuracy of the aspiration.
- No liquids with high viscosity:
 - The underpressure in the cylinder might not be enough to pull the liquid into the cylinder.
 - The cohesive force in the liquid might not be enough to hold it together.

Dispensing


- The burette drive pushes the piston up.
- Pressure builds in the cylinder.
- Liquid is pushed out through the open port.

See also

 Operation ▶ Page 9

4 Installation

4.1 Scope of delivery

Part		Order number
	dDrive	-
	Cable ACT M8/F M8/M 20 cm	30634406
	User Manual	-
	Declaration of conformity	-
	Test report	-

See also

 Accessories, spare parts and consumables ▶ Page 12

4.2 Download manuals

- 1 Go to the website www.mt.com/library.
- 2 Select the **Technical Documentation** tab.
- 3 Find the product type on the housing of the main instrument and enter it into the search field.
- 4 Start the search.

- 5 Select the manual from the result list.
- 6 Select the link.
 - ➔ The manual is either opened or downloaded depending on the browser settings.
- 7 Check which software version is installed on your main instrument.
- 8 If the manual is not written for the installed software version, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

4.3 Unpack the burette drive

- 1 Remove the burette drive from the protective packaging.
- 2 Store the protective packaging for later transport over long distances.
- 3 Check that you have received all parts listed in the scope of delivery.
- 4 Inspect the parts visually for flaws or damage.
- 5 If parts are missing or damaged, report it to your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

See also

🔗 [Scope of delivery](#) ▶ Page 7

4.4 Position the burette drive

The burette drive has been developed for indoor operation in a well-ventilated area.

The following site requirements apply:

- Ambient conditions within the limits specified in the technical data
- No powerful vibrations
- No direct sunlight
- No corrosive gas atmosphere
- No explosive atmosphere
- No powerful electric or magnetic fields

Procedure

- 1 Place the burette drive on a level surface.
- 2 Wait at least three hours before you connect the burette drive to the main instrument if the temperature during transport or storage was more than 20 °C lower than the temperature in the laboratory.
- 3 Install the cables in such a way that they cannot be damaged or interfere with operation.

Read the Reference Manual of the main instrument for more information about how to connect the burette drive to the main instrument.


See [[Download manuals](#) ▶ Page 7].

See also


🔗 [Technical data](#) ▶ Page 11

4.5 Install the burette

- The piston rod (2) is in its lower position.
- The lock wheel points to the open-lock symbol (1).

- 1 Slide the burette (4) into its mounting position (3).
- 2 Turn the lock wheel (1) clockwise to the closed-lock symbol (2).

See also

 Burette drive overview ▶ Page 6

4.6 Remove the burette


Read the Reference Manual of the main instrument for more information about how to move the piston rod to its lower position.

See [Download manuals ▶ Page 7].

- The piston rod of the burette drive is in its lower position.
- The burette and the connected tubes are empty.

- 1 Disconnect the tubes.
- 2 Turn the lock wheel (2) counterclockwise to the open-lock symbol (1).
- 3 Slide the burette (3) off the burette drive.

See also

 Burette drive overview ▶ Page 6

5 Operation

Read the Reference Manual of the main instrument for more information about how to control the burette drive.

See [Download manuals ▶ Page 7].

6 Maintenance

In this chapter, you find descriptions of the maintenance tasks you should perform on your burette drive. Any other maintenance tasks need to be performed by a service technician who has been qualified by METTLER TOLEDO.

If you experience problems with your burette drive, contact your authorized METTLER TOLEDO service representative or dealer.

METTLER TOLEDO recommends that a preventive maintenance and calibration certification is done at least once a year through your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

6.1 Maintenance schedule

Follow this maintenance schedule, unless otherwise required by the standard operating procedures of your company.

Every day

Task	Link
Check the burette drive and the burette visually for leakage.	–

6.2 Clean the burette drive

NOTICE

Damage to the burette drive due to incorrect cleaning methods

Incorrect cleaning agents can damage the housing or other parts of the burette drive. If liquids enter the housing, they can damage the burette drive.

- 1 Make sure the cleaning agent is compatible with the material of the part you want to clean.
- 2 Make sure that no liquid enters the interior of the burette drive.

METTLER TOLEDO recommends the following cleaning agents:

- Isopropanol
- Ethanol
- Water with a mild detergent

If you have questions about the compatibility of cleaning agents, contact your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

Procedure

- Wipe the housing with a cloth moistened with the cleaning agent.

See also

🔗 Technical data ► Page 11

6.3 Prepare the burette drive for storage or shipping

- 1 Remove the burette.
- 2 Disconnect the burette drive from the main instrument.
- 3 Clean the burette drive.
- 4 To store the burette drive, keep it in a dry and clean place.
- 5 To ship the burette drive, use the original packaging.

See also

🔗 Technical data ► Page 11

6.4 Dispose of the burette drive

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

7 Technical data

Power supply

Characteristic		Value
Input (IN)	Power consumption	24 V DC ± 10 %, 0.5 A
	Input rating	24 V DC ± 10 %, 3 A
	Socket	M8, 4 pin, male
Output (OUT) ¹⁾	Output rating	24 V DC ± 10 %, 2.5 A
	Socket	M8, 4 pin, female

¹⁾ The output has been evaluated for connections to non-hazardous safety extra-low voltage (SELV) circuits. The output must only be connected with non-hazardous safety extra-low voltage (SELV) circuits.

Instrument

Characteristic		Value
Dimensions	Width	81.2 mm
	Depth	167 mm
	Height	190.5 mm
Weight		approx. 1.75 kg
Resolution		120000 steps
Materials	Housing	PBT (polybutylene terephthalate)
	Piston rod	PEEK (polyether ether ketone)
	Socket cover	PA (polyamide)

Site requirements

Characteristic		Value
Ambient conditions	Ambient temperature	5...40 °C
	Recommended operational temperature ¹⁾	18...28 °C
	Relative humidity	Non-condensing, max. 80 % for temperatures up to 31 °C, decreasing linearly to 50 % at 40 °C
	Altitude	Up to 5000 m above sea level
	Use	In interior spaces
	Overvoltage category	II
	Pollution degree	2
Storage conditions	Temperature	-20...+70 °C
	Relative humidity	0...95 %


¹⁾ METTLER TOLEDO produces and tests the equipment with test tools that are certified for this temperature range. Usage outside the given range could lead to inferior performance such as burette leakage.

8 Accessories, spare parts and consumables

All accessories, spare parts and consumables are specified with their order number.


If you have any questions, contact your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact


Burette 2 mL
Two-port burette

30869286


Burette 5 mL
Two-port burette

30869287


Burette 10 mL
Two-port burette

30869288



Liquid Handler DispenseSix
Six-port burette

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© 02-2024 METTLER TOLEDO. All rights reserved.