METTLER **TOLEDO**

Sensor Head InSUS H60i for Single-Use O₂ Sensors

Instruction manual
Bedienungsanleitung
Instructions d'utilisation

InSUS H60i
30 802 188 A

English **3**

Deutsch **17**

Français **31**

Sensor Head InSUS H60i for Single-Use O₂ Sensors

Instruction manual

InSUS, iSense and ISM are trademarks of the METTLER TOLEDO Group.

Contents

	Page
1 Introduction.....	5
2 Important Notes.....	5
2.1 Notes on Operating Instructions	5
2.2 Intended Use	5
2.3 Safety Instructions	6
2.4 Examples of Some Typical Applications	6
3 Product Description	7
3.1 General Information.....	7
3.2 Principle	7
3.3 Scope of Delivery	8
4 Installation.....	8
4.1 Connection	8
4.1.1 Connection of the Single-Use Sensor to the InSUS H60i.....	8
4.1.2 Digital Connection of the InSUS H60i to a Transmitter	8
4.1.3 Analog Connection of the InSUS H60i to a Transmitter or Controller.....	9
4.1.4 Power Connection of InSUS H60i	9
4.1.5 Modbus Installation	10
5 Operation.....	10
5.1 Configuration	10
5.1.1 Sensor Head Detection	10
5.1.2 Sampling Rate	10
5.1.3 LED Mode	10
5.2 Calibration	11
5.2.1 Purpose of Calibration	11
5.2.2 Factory Calibration	12
5.2.3 Single Point Calibration (Slope or Process Calibration)	12
5.2.4 Process Calibration	12
5.2.5 Dual Point Calibration with M400/M800 Transmitter.....	13
5.2.6 Calibration When Connected with Analog Signal	13
6 Maintenance	13
6.1 Sensor Head Inspection	13
6.1.1 Visual Inspection	13
6.1.2 Testing the Sensor/Sensor Head with the Transmitter	14
6.1.3 Intelligent Sensor Management: ISM.....	14
7 Storage, Warranty and Disposal.....	14
8 Product Specification.....	15
8.1 Specifications	15
9 Ordering information	15
9.1 Sensor Heads	15
9.2 Sensors	15
9.3 Accessories	16
9.4 Recommended Transmitters	16

1 Introduction

Thank you for buying the **optical oxygen sensor head from METTLER TOLEDO Ingold**.

The InSUS™ H60i is operated together with optical single-use O₂ sensors of the InSUS series. InSUS O₂ sensors, for example InSUS 607 sensors, are typically provided pre-installed and pre-sterilized in single-use process devices.

The construction of METTLER TOLEDO's optical oxygen sensors employs leading-edge technology and complies with safety regulations currently in force. However, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment.

Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor head.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor head.

If you have questions that are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Important Notes

2.1 Notes on Operating Instructions

These operating instructions contain all the information needed for safe and proper use of the optical sensor head.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended Use

METTLER TOLEDO optical O₂ sensors and sensor heads are intended solely for inline measurement of the oxygen partial pressure, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose. The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- Compliance with the instructions, notes and requirements set out in this instruction manual.
- Acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- Compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (single-use bags, transmitters, etc.).
- Observance of all safety regulations governing the equipment in which the sensor is installed.

- Correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- Consultation with METTLER TOLEDO Process Analytics in the event of any uncertainties.

2.3 Safety Instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. single-use bags, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor head must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see chapter 9).
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorized modifications. The risk is borne entirely by the user.

2.4 Examples of Some Typical Applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

Measurement in liquids:

- Fermentation
- Bio-Tech

3 Product Description

3.1 General Information

The optical single-use oxygen sensor and the InSUS H60i sensor head (with integrated temperature probe) are used for measurement of oxygen in single-use process devices.

The single-use sensors are gamma and X-ray irradiation sterilizable and are typically integrated in single-use devices such as single-use bags and other similar applications. Please also refer to the documentation provided by your single-use device supplier.

3.2 Principle


The optical oxygen sensors are based on an optical detection method: fluorescence quenching. Here is a short summary of the principle. In contrast to the polarographic Clark-electrode, which detects a redox reaction of oxygen at the electrode, the optical method is based on an energy transfer between a chromophore and oxygen.

- A chromophore, embedded in the sensor is illuminated with blue light. This chromophore absorbs the energy and if no oxygen is present emits red fluorescence light with a specific lifetime. This emitted light is being detected by a detector in the sensor head.
- In the presence of oxygen, the chromophore transfers the energy to the oxygen molecule. Oxygen is then able to transfer this energy as heat to the surrounding area and no fluorescence is emitted.
- The total intensity of the fluorescence and the lifetime of the fluorescence is related to the oxygen partial pressure in the medium.
- To analyze the lifetime of the fluorescence, the excitation light is pulsed with a constant frequency, the emitted light shows the same course but with a time delay to the excitation. This delay is called Phase shift or Phase angle (Φ). The phase shift is dependent on the oxygen level and follows the Stern-Vollmer correlation.
- The sensor head detects this phase shift and calculates the oxygen concentration.
- The oxygen value is digitally transferred from the sensor head to the transmitter.

3.3 Scope of Delivery

Each sensor head is supplied fully assembled and factory tested for correct function together with:

- a protective cap and a cap clip
- this manual and a Quick Setup Guide
- a Quality Control Certificate
- a Declaration of Conformity
- a 3.1 material certificate

InSUS H60i with protective cap.

InSUS H60i connected with an InSUS O₂ single-use sensor, protective cap on clip holder.

- 1: VP8 connector
- 2: InSUS H60i
- 3: Clip
- 4: Sensor connector
- 5: Protective cap
- 6: InSUS O₂ sensor – not supplied with InSUS H60i

4 Installation

4.1 Connection

4.1.1 Connection of the Single-Use Sensor to the InSUS H60i

Remove the rubber cap from the InSUS O₂ sensor and unscrew the protective cap from the InSUS H60i. Mount the protective cap on the clip holder. Screw the InSUS H60i sensor connector to the InSUS O₂ sensor.

Note: For protection of the optical tip of the InSUS H60i from damage and dirt, reattach the protective cap to the sensor head whenever the InSUS H60i is not connected to an InSUS O₂ sensor.

4.1.2 Digital Connection of the InSUS H60i to a Transmitter

Transmitter M400 or M800

The sensor head is connected to the transmitter via a VP8 cable. The cables are available from METTLER TOLEDO in different lengths. The cable ensures a secure connection between the transmitter and the sensor head under harsh industrial conditions. To connect the cable to the sensor head, align the slit of the connector with the pin in the plug, then tightly screw the plug to fasten the two parts.

Transmitter M400

Note: For connecting the cable to the terminals of the transmitter, please refer also to the instructions given in the METTLER TOLEDO transmitter manual.

Connect the cable to the transmitter as described in the tables below.

VP8 cable for InSUS H60i

		M400	M800 1-channel	M800 2/4-channel
Color	Function	TB3	TB3	TB2
gray	24 V	18	7	9
blue	GND 24 V	17	8	10
green/yellow	Shield	13	2	12
brown	RS485B	14	3	13
pink	RS485A	15	4	14

4.1.3 Analog Connection of the InSUS H60i to a Transmitter or Controller

The InSUS H60i is equipped with two different interfaces:

1. Digital communication to a METTLER TOLEDO transmitter or a Modbus RTU system.
2. Analog output (nA or 4–20 mA) plus a simulated NTC 22 kΩ temperature signal. The mA version includes HART communication.

Attention: Never connect an InSUS H60i with mA output to an nA input. This could damage the electronics of the sensor head and/or instrument the sensor head is connected to.

VP8 cable for InSUS H60i nA or InSUS H60i mA

	InSUS H60i nA	InSUS H60i mA
Color	Function	Function
black/transp	Cathode (nA)	4–20 mA+ HART
red	Anode (nA)	4–20 mA– HART
gray	24 V	24 V
blue	GND 24 V	GND 24 V
white	NTC 22 kΩ	NTC 22 kΩ
green	GND NTC 22 kΩ	GND NTC 22 kΩ
green/yellow	Shield	Shield

Additional to this manual, the following documents and tools are available:

- Modbus Programmers Guide
- Optical O₂ HART communication manual
- mA-Output configuration tool

4.1.4 Power Connection of InSUS H60i

If the sensor head is used with a 4-wire transmitter (M400/M800), the sensor head is powered by the transmitter. In this case, no additional power supply is needed. For installation with 2-wire and BUS transmitters (M400/2(X)H; M400 FF; M400 PA) and for direct integration with a nA, mA and Modbus transmitter an additional power supply is required.

Power Supply Specification: please refer to chapter 8.1.

4.1.5 Modbus Installation

The InSUS H60i offers a Modbus RTU connection that is identical to METTLER TOLEDO's InPro 6860i optical oxygen sensor. For correct installation, please follow the instructions of the control system. For programming, please follow the information given in the "Modbus Programmers' Guide InPro 6860i" which is available at your local METTLER TOLEDO Partner.

The Modbus interface allows using the full functionality of the sensor head including all Intelligent Sensor Management (ISM™) functions.

5 Operation

5.1 Configuration

5.1.1 Sensor Head Detection

Before installing an optical sensor head, please refer to the manual for the transmitter and configure the transmitter for automatic sensor detection. Date and time must be set correctly in the transmitter, otherwise diagnostic functions such as DLI/ACT will not work properly.

If the wrong date and time are set, the calibration and setup might be corrupted.

Note: Each time the sensor head is connected with a new single-use oxygen sensor and diagnostic functions are activated, DLI/ACT settings must be reset on the transmitter. Please refer to the manual of the transmitter.

5.1.2 Sampling Rate

Optical oxygen sensors do not measure permanently. Each measurement cycle has a duration of approx. 1 second. The measurement interval can be set to any value between 1 and 60 seconds. Please choose the appropriate setting. The default setting is 3 seconds, which is sufficient for most applications.

5.1.3 LED Mode

The measurement can be switched off if the system is not needed.

When the sensor head is not measuring, the sensor head LED is off. In this state the sensor head sends a constant measurement value of -1 % air to the transmitter and the transmitter is set to the "Hold mode". To configure the "Hold mode" please refer to the transmitter manual.

Automatic switch off at high temperature

If the LED mode is set to "Auto" (default setting) the sensor head LED will be switched off as soon as a specific process temperature is reached.

Temperatures

	Maximum operating temperature	Default switch off temperature
InSUS H60i	60 °C (140 °F)	60 °C (140 °F)

This limit can be set to an individual value by the user. Using the transmitter (M400 or M800) or with iSense™. These settings are also active if the sensor head is connected to a controller with an analog output. The switch off temperature should be set at least 5° higher than the highest process temperature. For example, if the process temperature is 37 °C (99 °F), 42 °C (104 °F) should be the minimum setpoint. In this situation, as soon as the temperature exceeds 42 °C (104 °F) the sensor head will stop measuring and the LED will be switched off. For the switch on, a hysteresis of 3° is implemented, meaning that the sensor head (and LED) will be switched on as soon as the temperature drops below 39 °C (101 °F).

Manual switch off of the sensor head (M400/M800 transmitter)

The sensor head can be switched off manually via the transmitter menu (see the transmitter manual) by setting the LED mode to "off". To restart the measurement, the LED mode needs to be set manually to "on" via the transmitter menu, or via a remote signal (digital input).

Remote switch off of the sensor head (M400/M800 transmitter)

The M400 transmitter can be set to "Hold" by applying an external digital signal (see the transmitter manual). In this situation the sensor head and the sensor head LED are switched off. As soon as the "Hold Mode" is off, the optical sensor head will continue to measure using the previous settings.

5.2 Calibration

5.2.1 Purpose of Calibration

The manual of the transmitter also contains information about calibration.

Calibration must be performed every time a new single-use sensor is connected to the InSUS head.

Since the correlation between the measured phase and the oxygen value is not linear, a calibration of an optical sensor must be performed very accurately. Wrong calibrations may significantly reduce the measurement accuracy.

Each single-use oxygen sensor has its own individual phase angle at zero oxygen (phi 0) and hundred percent air saturation (phi 100).

Several methods for calibration are available for the optical oxygen sensors. The highest measurement accuracy is achieved by performing a 2-point calibration with air and a zero gas e.g. N₂ or CO₂ with a purity of at least 99.995 %.

Note: Please take into account the correct air pressure and humidity. Small deviations in air ($\pm 3\%$) are due to differences in humidity and process pressure settings. The sensor head calculates for 100 % humidity if it is set to dissolved oxygen measurement.

General remarks:

- **Calibration must be performed after irradiation sterilization** of the single-use device.
- **For calibration in gas** the single-use device in which the single-use sensor is integrated **must be filled with air (phi 100) or zero point gas (phi 0)**.
- Make sure that the settings for **oxygen saturation** of the calibration is **correct and remains constant** during calibration.
- In the event of calibration in water or process medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is very slow. Therefore it takes quite a long time until water is saturated with atmospheric oxygen.

- Make sure that all other parameters, such as temperature and pressure, are constant.
- Calibration always needs accurate pressure and temperature measurement. Only process scaling is independent of those parameters (see chapter 5.2.4).
- Make sure that the correct calibration pressure, humidity and salinity values are set in the transmitter before the calibration is started.
- Please refer also to the transmitter or biocontroller manual for detailed information.

5.2.2 Factory Calibration

The sensor head is delivered pre calibrated and ready for use.

The factory calibration data are stored in the sensor head and cannot be changed by the user.

5.2.3 Single Point Calibration (Slope or Process Calibration)

For most applications, a single point calibration should be sufficient.

By carrying out a single point calibration, the factual phase at the desired oxygen value e.g. at hundred percent oxygen (phi 100) of the sensor can be established. The corresponding calibration curve is calculated.

The calibration medium can be either air or a calibration gas with a known O₂ concentration or process medium with a known oxygen concentration.

Before starting a calibration in gas, the correct pressure and the correct humidity have to be set in the transmitter.

Note: Wrong pressure values are the most common reasons for bad measurement accuracy.

E.g. A 50 mbar difference between the ambient pressure and the value set in the transmitter result in a 5 % measurement error in air.

For calibration in gas, it is important that the temperature reading of the sensor is stable and represents the real gas temperature.

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O₂, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

5.2.4 Process Calibration

For detailed information please also refer to the transmitter manual.

Two different routines for process calibration are possible:

- **Process calibration**
- **Process scaling (M400/M800 transmitter)**

Process calibration is performed when a reliable control value is available and process pressure is known. Process pressure is only needed if the system is measuring in saturation (% air or % O₂) or gas (ppm gas) units. During this calibration the phase values of the calibration curve are adjusted.

Process scaling is performed mainly in biopharma applications when the user desires to set the system to an initial value. During this calibration the phase values of the sensor are not adjusted, only the displayed values and the 4–20 mA output are rescaled to the desired value.

Note: For process calibration, the operator can use either the process pressure or the calibration pressure, depending on how the reference value is taken.

After the sensor signal has stabilized, the complete measurement system can be calibrated to the desired variable, e.g. % air, % O₂, ppm or ppb (see instruction manual for the transmitter).

 **Note:** For this type of calibration an accurate reference value and correct pressure settings are essential.

5.2.5 Dual Point Calibration with M400/M800 Transmitter

To receive maximum accuracy of the measured values over the full measuring range, a dual point calibration is required.

By carrying out a dual point calibration, both phase angles, at zero oxygen (phi 0) and at hundred percent oxygen (phi 100), of the single-use sensor can be established.

Point 1: Slope correction (with air or other calibration media with known O₂ value)

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O₂, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Point 2: Zero point

After the sensor signal has stabilized, the sensor can be calibrated to the 0 % value of the desired measurable variable, e.g. 0 % air, 0.0 % O₂, or 0 ppm at 25 °C (77 °F) (see instruction manual for the transmitter).

 **Note:** Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.995 %.

5.2.6 Calibration When Connected with Analog Signal

As soon as a sensor head is connected with an analog communication, only rescaling of the sensor head output is available.

The procedures can be performed the same way as with amperometric sensors:

1. Slope correction
2. Offset correction

6 Maintenance

 **Note:** All maintenance work can be done without any tools.

6.1 Sensor Head Inspection

6.1.1 Visual Inspection

To check your sensor head, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before connection to a single-use O₂ sensor always examine both components optically for signs of damage. The two parts must be clean and intact. Dirty surfaces of sensor and sensor head should be wiped clean using a soft, moist tissue.

 **Attention!** Do not use any cleaning agents containing alcohol or any solvents. This could damage the products.

6.1.2 Testing the Sensor/Sensor Head with the Transmitter

If the measured values differ from the expected value, an air calibration should be performed.

Appropriate values for PhiO and Phi100 are:

InSUS O₂ sensors

PhiO	Phi100
≥48°	≥20°

The phase values of the sensor are stored in the calibration history. The actual phase value can be checked in the "Calibration – Verify" menu.

If after such procedures the above mentioned values are still not reached, the single-use sensor should not be used.

Zero oxygen measurement can be done by using CO₂ or nitrogen (N₂), or alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 5 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If after such procedures the above mentioned values are still not reached, the single-use sensor should not be used. For a detailed inspection of the InSUS H60i please contact your local METTLER TOLEDO representative.

6.1.3 Intelligent Sensor Management: ISM

Note: Each time the sensor head is connected with a new single-use oxygen sensor and diagnostic functions are activated, DLI/ACT settings must be reset on the transmitter. Please refer to the manual of the transmitter.

Dynamic Lifetime Indicator: DLI

The DLI provides information about the remaining lifetime of an optical InSUS O₂ sensor. As long as the DLI is above zero days, the system is within the specified accuracy after a calibration. If the DLI is zero, the InSUS O₂ sensor should not be used.

Contributing factors for aging of the optical InSUS O₂ sensor:

- number of measurements
- temperature during measurement
- oxygen concentration during measurement

Adaptive Calibration Timer: ACT

The ACT provides information as to when the next calibration is required to ensure measurements will remain within the specified accuracy. This calculation is based on the DLI information.

7 Storage, Warranty and Disposal

For storage, the sensor head should be clean and dry. The protection caps have to be placed on the sensor head and the cable connectors. In the event of manufacturing faults, 12 months warranty from date of production is granted.

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority for disposal or recycling advice.

8 Product Specification

8.1 Specifications

InSUS H60i	
Measurement principle	Optical
Power requirements	$U_{\min} = 19.5 \text{ VDC}$ $U_{\max} = \leq 25 \text{ VDC}$ $P_{\max} = 0.75 \text{ W (33...40 mA)}$
Analog/digital interface	
nA output (InSUS H60i nA)	60 nA at air (25 °C, 1013 mbar)
mA output (InSUS H60i mA)	4...20 mA (active)
Analog temperature output	NTC 22 kΩ
Digital interface	Modbus RTU on RS485 <ul style="list-style-type: none"> • Baudrate: 4800, 9600, 19200 (default), 38400, 56600, 115000 <ul style="list-style-type: none"> • Device address: 1 (default) to 32
Design features	
Temperature compensation	Automatic with built-in RTD
Cable connection	VP8
Certificates/Compliances	
Quality certificate	Yes
Declaration of Conformity (CE, UKCA)	Yes
Material certificate 3.1	Yes
Performance together with optical InSUS O₂ single-use sensor	
Measuring range	0...250 % air
Accuracy under defined laboratory conditions	< 2.5 % for the range 50–100 % air after 1-point calibration in 100 % air <1 % after 2-point calibration in 100 % air and 0 % oxygen
Response time	25 °C (77 °F) air → nitrogen, t98 % <30 s
Operating temperature	5...60 °C (41...140 °F)
Sensor connection	M16x1

9 Ordering information

For more detailed information ask your local distributor.

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All rights reserved