**METTLER TOLEDO**

**InPro® 6000  
Optical O<sub>2</sub> Sensors**

**Instruction manual  
Bedienungsanleitung  
Instructions d'utilisation**

**English** **Page** **3**

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# **InPro<sup>®</sup> 6000 Series Optical O<sub>2</sub> Sensors**

## **Instruction manual**

## Contents

<b>1</b>	<b>Introduction</b> .....	<b>5</b>
<b>2</b>	<b>Important notes</b> .....	<b>6</b>
2.1	Notes on operating instructions .....	6
2.2	Intended use .....	6
2.3	Safety instructions .....	7
2.4	Examples of some typical applications .....	8
2.5	Use in Ex-zones .....	8
2.6	Ex-classification ATEX .....	8
2.6.1	Introduction .....	8
2.6.2	Rated data .....	9
2.6.3	Special conditions .....	9
<b>3</b>	<b>Product description</b> .....	<b>11</b>
3.1	General information .....	11
3.2	Principle .....	11
3.3	Scope of delivery .....	12
3.4	Equipment features .....	12
<b>4</b>	<b>Installation</b> .....	<b>15</b>
4.1	Mounting the sensor .....	15
4.1.1	Retrofit kit for InPro 6870 i/InPro 6960 i/InPro 6970 i .....	15
4.2	Connection .....	16
4.2.1	Connecting the optical sensor to a cable .....	16
4.2.2	Digital connection to a transmitter .....	16
4.2.3	Analog connection to transmitter or controller .....	17
4.2.5	MODBUS Installation .....	18
4.2.6	Power Connection of sensor .....	18
<b>5</b>	<b>Operation</b> .....	<b>19</b>
5.1	Start-up .....	19
5.2	Configuration .....	19
5.2.1	Sensor detection .....	19
5.2.2	Sampling rate .....	19
5.2.3	LED mode .....	19
5.3	Calibration .....	20
5.3.1	Purpose of calibration .....	20
5.3.2	Factory calibration .....	22
5.3.3	Single point calibration (Slope or Process Calibration) .....	22
5.3.4	Slope calibration (InPro 6870 i and InPro 6960 i) .....	22
5.3.5	Process calibration .....	23
5.3.6	Dual point calibration .....	23
5.3.7	Calibration when connected with analog signal .....	24
5.3.8	Reset to factory calibration .....	24
<b>6</b>	<b>Maintenance</b> .....	<b>25</b>
6.1	Sensor inspection .....	25
6.1.1	Visual inspection .....	25
6.1.2	Testing the sensor with the transmitter .....	25
6.1.3	ISM .....	26
6.2	Replacing the OptoCap .....	27
<b>7</b>	<b>Storage</b> .....	<b>28</b>
<b>8</b>	<b>Product specification</b> .....	<b>29</b>
8.1	Certificates .....	29
8.2	Specifications .....	29
<b>9</b>	<b>Ordering information</b> .....	<b>31</b>
9.1	Sensors .....	31
9.2	Accessories .....	31
9.3	Spare parts .....	32
9.4	Recommended transmitters .....	32
9.5	Recommended housings .....	33
<b>10</b>	<b>Theory of the optical oxygen measurement</b> .....	<b>34</b>
10.1	Introduction .....	34
10.2	Principle .....	34
10.3	Principle of the design of the optical oxygen sensor .....	34
10.4	Temperature .....	35
10.5	Dependence on flow .....	35
10.6	Oxygen partial pressure – oxygen concentration .....	35

## 1 Introduction

Thank you for buying the **optical oxygen sensor from METTLER TOLEDO**.

The construction of INGOLD's optical oxygen sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment.

**Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

## 2 Important notes

### 2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the optical sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

#### Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

### 2.2 Intended use

**METTLER TOLEDO optical O<sub>2</sub> sensors are intended solely for inline measurement of the oxygen partial pressure, as described in this instruction manual.**

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose. The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor ( housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.
- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

## 2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, "Section 9.3").
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

## 2.4 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

### Measurement in liquids:

- Fermentation
- Bio-Tech
- Food & Beverage

## 2.5 Use in Ex-zones

### Attention!

For an installation in Ex-zones please read the guidelines following hereafter:

### Ex-classification ATEX:

II 1/2G Ex ia/ib IIC T6 Ga/Gb

II 1/2D Ex ia/ib IIc T 83 °C Da/Db

### Marking and number of the test certificate:

SEV 14 ATEX 0127X

IECEX SEV 14.0007X

## 2.6 Ex-classification ATEX

### 2.6.1 Introduction

According to RL 94/9/EG (ATEX 95) Appendix I, InPro 6XXX\*/\*/\*/\* oxygen sensors are devices group II, category 1/2G and according to RL 99/92/EG (ATEX 137) may be used in zones 0/1 or 0/2 and gas groups IIA, IIB and IIC that are potentially explosive due to combustible substances in the temperature class T6.

For use/installation, the requirements of EN 60079-14 must be observed.

According to RL 94/9/EG (ATEX 95) Appendix I, InPro 6XXX\*/\*/\*/\* oxygen sensors are devices group II, category 1/2D and according to RL 99/92/EG (ATEX 137) may also be used in zones 20/21 that contain combustible dusts.

The digital sensor circuit is part of a common intrinsically safe system and is for operation connected to a separately certified transmitter.

The digital sensor circuit as part of an intrinsically safe system is isolated from the not-intrinsically safe electric circuits up to a maximum rated voltage of 375 V and from grounded parts up to a maximum rated voltage of 30 V.

## 2.6.2 Rated data

### Measuring circuit:

Fail-safe ignition protection class Ex ia IIC only for connection to a certified fail-safe circuit.

### Maximum values:

$$U_i \leq 25 \text{ V}$$

$$I_i \leq 60 \text{ mA}$$

$$P_i \leq 1.5 \text{ W}$$

$$L_i = 0$$

$$C_i = 0$$

### O<sub>2</sub> sensor InPro 6860 i safe parameters

#### 1. RS485 interface

$$U_i = 15 \text{ V}, I_i = 100 \text{ mA}, P_i = 1 \text{ W}, L_i = 0 \text{ mH}, \\ C_i = 2 \text{ microF}$$

$$U_o = 4.6 \text{ V}, I_o = 91 \text{ mA}, P_o = 0.3 \text{ W}, L_o = 0 \text{ mF}, \\ C_o = 100 \text{ pF}$$

#### 2. mA HART output

$$U_o = 13.93 \text{ V}, I_o = 25 \text{ mA}, L_o = 0 \text{ mF}, \\ C_o = 100 \text{ nF}$$

#### 3. nA output

$$U_o = 7.5 \text{ V}, I_o = 1.46 \text{ mA}, L_o = 0 \text{ mF}, \\ C_o = 0 \text{ nF}$$

NTC output

$$U_o = 6.7 \text{ V}, I_o = 60 \text{ mA}, L_o = 0 \text{ mF}, \\ C_o = 1.01 \text{ microF}$$

### Note:

The above maximum values are each the total of all individual circuits of the associated intrinsically safe power supply and transmitter.

## 2.6.3 Special conditions

- The maximum permissible environment resp. medium temperature for the zone 0 (combustible gases or combustible liquids) is in accordance with the temperature class shown in the table below:

Temperature class	Max. environment resp. media temperature
T 6	60 °C

- The maximum surface temperature for the zone 20 (combustible dusts) is in accordance with the environment resp. medium temperatures shown in the table below:

Surface temperature	Max. environment resp. media temperature
T 83 °C	60 °C

- The capacitance and inductance of the connecting cable must be taken into account in the design.
- The oxygen sensors (O<sub>2</sub> sensors) can be used in/with the housings InFit 76\*-\*\*\* resp. InTrac 7\*\*-\*\*\* or in/with other suitable housings in hazardous areas.

- The metal body of the O<sub>2</sub> sensors resp. the safety weld-in-sockets resp. the independent housing are, if necessary, to be included into the periodic pressure test of the unit.
- The metal body of the O<sub>2</sub> sensors resp. the safety weld-in-sockets resp. the independent housing must be electrically connected to the potential equalizing system of the installation.

## 3 Product description

### 3.1 General information

The **optical oxygen sensors** with integrated temperature probe **are used for measurement of oxygen.**

The sensors are sterilizable and compatible with CIP (cleaning in place). The sensor InPro 6860i is also autoclavable (max. 284 °F).

### 3.2 Principle

The optical oxygen sensors are based on an optical detection method, the so called fluorescence quenching. Here is a short summary of the principle. In contrast to the polarographic Clark-electrode, which detects a redox reaction of oxygen at the electrode, the optical method is based on an energy transfer between a chromophore and oxygen.

- A chromophore, embedded in the sensor tip is illuminated with blue light. This chromophore absorbs the energy and if no oxygen is present emits red fluorescence light with a specific lifetime. This emitted light is being detected by a detector in the sensor head.
- In the presence of oxygen, the chromophore transfers the energy to the oxygen molecule. Oxygen is then able to transfer this energy as heat to the surrounding area and no fluorescence is emitted.
- The total intensity of the fluorescence and the lifetime of the fluorescence is related to the Oxygen partial pressure in the medium.
- To analyze the lifetime of the fluorescence, the excitation light is pulsed with a constant frequency, the emitted light shows the same course but with a time delay to the excitation. This delay is called Phase shift or Phase angle (Phi). The phase shift is dependent on the oxygen level and follows the Stern-Vollmer correlation.
- The sensor detects this phase shift and calculates the oxygen concentration.
- The oxygen value ist digitally transferred to the transmitter.

### **3.3 Scope of delivery**

Each sensor is supplied fully assembled and factory tested and calibrated for correct function together with:

- a quality control certificate
- inspection certificate 3.1  
(complying with EN 10204.3/1 B)


### **3.4 Equipment features**


Two different types of optical sensors are available for the fermentation application.

The InPro 6860i can be used in applications where autoclaving is performed. The InPro 6860i is fully autoclavable.

The InPro 6870i is intended for use in fermenter, where only sterilizations are performed. Here, the sensor head is not exposed to high temperatures.

The InPro 6960i/6970i is intended for use in food and beverage industry for measuring high and low ppb level of oxygen respectively.

**InPro 6860 i**

## 4 Installation

The "Plug & Measure"-concept allows the user to measure oxygen values immediately after installation.

- The sensor is recognized automatically and all relevant data are sent to the transmitter, the oxygen values are displayed.

The sensor is factory calibrated. The factory calibration data are stored in the sensor and do not need to be entered by the user.

### 4.1 Mounting the sensor


**Important!** Remove the protection cap before mounting the sensor.

#### Mounting the sensor in a housing

Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

#### Mounting the sensor directly on a pipe or a vessel

The **12 mm sensors** can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve. The sensor can be mounted in any orientation.


**Attention! Do NOT turn the sensor anti-clockwise when installing or removing it from a housing/process connection.** Use the turning nut to secure or unsecure the sensor from the housing/process connection and, if required, turn the sensor clockwise.

Turning the sensor anti-clockwise can result in loosening or removal of the sensor tip allowing ingress of fluid into the sensor body.

#### 4.1.1 Retrofit kit for InPro 6870 i/InPro 6960 i/InPro 6970 i

For protection of the sensor head against environmental influences (temperature, humidity) the use of the retrofit kit (52 403 811) is recommended. For installation, see the manual of the retrofit kit.

Retrofit kit

## 4.2 Connection

### 4.2.1 Connecting the optical sensor to a cable

The sensor is connected to the transmitter via a 5 pin data cable for InPro 6870i, 6960i and 6970i or a VP6/VP8 cable for InPro 6860i. The data cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP 67 connector housing guarantees maximum process safety.

To connect the data cable to the sensor align the slit of the connector with the pin in the plug. Then tightly screw the plug to fasten the two parts.

### 4.2.2 Digital connection to a transmitter

Transmitter M400

 **Note:** Cable assignment can be found in the METTLER TOLEDO cable instruction manual.

 **Note:** For connecting the cable to the terminals of the transmitter, please refer also to the instructions given in the METTLER TOLEDO transmitter manual.

The cables are available from METTLER TOLEDO in various lengths:

Connect the data cable to the transmitter as described in the tables below.

#### RS485 cable for InPro 6870 i/InPro 6960 i/InPro 6870 i

Color	Function	M400	M800	M800	M400
		TB4	1/2 channel TB2	4 channel TB2 or TB4	2w/FF/PA TB2
brown	24 VDC+	1	9	9	n.c.
black	24 VDC-	2	10	10	n.c.
gray	shield	6	12	12	M
yellow	shield	6	15	15	n.c.
blue	RS485-	7	13	13	N
white	RS485+	8	14	14	O

**VP8 cable InPro 6860 i/nA**

Color	Function	M400	M800	M800	M400
		TB4	1/2 channel TB2	4 channel TB2 or TB4	2w/FF/PA TB2
gray	24 VDC+	1	9	9	n.c.
blue	24 VDC-	2	10	10	n.c.
pink	RS485+	8	14	14	O
brown	RS485-	7	13	13	N
green/ yellow	shield	4	12	12	M

For M400 2-w, M400FF and M400PA, the optical oxygen sensor needs to be powered separately with an appropriate power supply.

**4.2.3 Analog connection to transmitter or controller**

The InPro 6860 i is equipped with two different interfaces:

1. Digital communication to the transmitter and any Modbus system.
2. Analog output (nA or 4–20 mA) plus a simulated NTC 22 k $\Omega$  temperature signal. The mA version includes HART communication.

**Attention:** Never connect an InPro 6860 i with mA output to a nA input. This could damage the electronics of the sensor and/or the instrument the sensor is connected to.

**VP6 cable for InPro 6860 i/nA or InPro 6860 i/mA**

Color	InPro 6860 i/nA	InPro 6860 i/mA
	Function	Function
black/transp	cathode (nA)	4–20 mA+ HART
red	anode (nA)	4–20 mA- HART
gray	24 VDC, 100 mA+	24 VDC+
blue	24 VDC, 100 mA-	24 VDC-
white	NTC 22 k $\Omega$	NTC 22 k $\Omega$
green	NTC 22 k $\Omega$ (GND)	NTC 22 k $\Omega$ (GND)
green/yellow	shield	shield

Additionally to this manual the following documents and tools are available:

- Modbus Programmers Guide
- Optical O<sub>2</sub> HART communication manual
- mA-Output configuration tool

The materials can be downloaded here:

[www.mt.com/pro-ODO-documentation](http://www.mt.com/pro-ODO-documentation)

#### 4.2.5 MODBUS Installation

The InPro 6860i offers a Modbus RTU connection. For correct installation please follow the instructions of the control system. For programming please follow the information given in the "Modbus Programmers Guide" which is available at your local METTLER TOLEDO Partner.

The Modbus interface allows to use the full functionality of the sensor including all ISM functions.

##### Wiring:

##### Protocol definitions, as implemented in InPro 6860 i

Modbus mode: RTU

Start bits: 1

Data bits: 8

Stop bits: 2

Parity: None

Baudrate: 19200 (default), 4800, 9600, 38400, 56600, 115000

Device address: 1 (default) to 32

##### Modbus RTU function codes implemented in InPro 6860 i

#3 Read Holding Registers

#4 Read Input Registers

#16 Write Multiple Registers

For detailed information please refer to the document: Modbus Programmers guide ([link](#))

#### 4.2.6 Power Connection of sensor

If the sensor is used with a 4-wire transmitter (M400/M800), the Sensor is powered over the transmitter. In this case no additional Power Supply is needed.

For installation with 2-wire and BUS-Transmitters (M400/2(X)H; M400 FF; M400 PA) and for direct integration with nA, mA and Modbus an additional power supply is required.

Power Supply Specification: 24 VDC; 500 mA

In Explosion proof applications the power supply must comply the specifications accordingly.

Mettler Toledo offers a current limited Powers Supply for use in Exi environment (see chapter 9).

## 5 Operation

### 5.1 Start-up

Each sensor is supplied ready to use. Before using remove the protecting cap.

**Note:** No polarization or calibration is necessary. "Plug and measure" conception.

### 5.2 Configuration

#### 5.2.1 Sensor detection

Before installing an optical sensor, please refer to the manual for the transmitter and configure the transmitter for automatic sensor detection. Date and time must be set correctly in the transmitter, otherwise diagnostic functions like DLI/ACT will not work properly.

In the case wrong date and time are set, the calibration and setup might be corrupted

#### 5.2.2 Sampling rate

Optical oxygen sensors do not measure permanently. Each measurement cycle has a duration of approx. 1 second. To prolong the lifetime of an OptoCap, the measurement interval can be set to any value between 1 and 60 seconds. Please choose the appropriate setting. Default setting is 10 seconds which is sufficient for most applications.

#### 5.2.3 LED mode

One contributing factor for the ageing of an OptoCap is the measurement itself. To prolong the lifetime of the OptoCap, the measurement can be switched off if the system is not needed. Especially measurement during CIP cycles or when the sensor is exposed to high oxygen levels during standby of the plant must be avoided. Measurement during CIP cycles can cause high sensor drift

When the sensor is not measuring, the sensor LED is off. In this state the sensor sends a constant measurement value of -1 % air to the transmitter and the transmitter is set to the "Hold mode". To configure the "Hold mode" please refer to the transmitter manual.

##### **Automatic switch off at high temperature**

If the LED mode is set to "Auto" (default setting) the sensor LED will be switched off as soon as a specific process temperature is reached. The default temperature set points are different for the different optical sensors.

<b>Temperatures</b>	<b>Maximum operating temperature</b>	<b>Default switch off temperature</b>
InPro 6860 i	60 °C / 140 °F	60 °C / 140 °F
InPro 6870 i	60 °C / 140 °F	60 °C / 140 °F
InPro 6960 i	40 °C / 104 °F	40 °C / 104 °F
InPro 6970 i	40 °C / 104 °F	40 °C / 104 °F

This limit can be set to an individual value by the user. Using the transmitter (M400 or M800) or with iSense. These settings are also active if the sensor is connected to the process with an analog connection (nA). The switch off temperature should be set at least 5° higher than the highest process temperature. For example, if the process temperature is 37 °C/99 °F, 42 °C/104 °F should be the minimum set-point. In this situation, as soon as the temperature exceeds 42 °C/104 °F the sensor will stop measuring and the LED will be switched off. For the switch on, a hysteresis of 3° is implemented, meaning that the sensor (and LED) will be switched on as soon as the temperature drops below 39 °C/101 °F.

#### **Manual switch off of the sensor (M400/M800 transmitter)**

The sensor can be switched off manually via the transmitter menu (see the transmitter manual) by setting the LED mode to "off". To restart the measurement, the LED mode needs to be set manually to "on" via the transmitter menu, or via a remote signal (digital input).

#### **Remote switch off of the sensor (M400/M800 transmitter)**

The M400 transmitter can be set to "Hold" by applying an external digital signal (see the transmitter manual). In this situation the sensor and the sensor LED are switched off. As soon as the "Hold Mode" is off, the optical sensor will continue to measure using the previous settings.

## **5.3 Calibration**

### **5.3.1 Purpose of calibration**

Information about the calibration, you find also in the manual of the transmitter.

Calibration should be performed at least after each change of the OptoCap, sterilizing or after autoclaving the sensor.

Since the correlation between the measured phase and the oxygen value is not linear, a calibration of an optical sensor must be performed very accurately. Wrong calibrations may significantly reduce the measurement accuracy and result in incorrect calculation of the Dynamic Lifetime Indicator (DLI) and the Adaptive Calibration Timer (ACT).

Each oxygen sensor has its own individual phase angle at zero oxygen ( $\phi_0$ ) and hundred percent air saturation ( $\phi_{100}$ ). Both values are subject to change, for example, after exchange of OptoCap or because of normal ageing of the OptoCap.

Depending on the sensor type several methods for calibration are available for the optical oxygen sensors. The highest measurement accuracy is achieved by performing a 2-point calibration with air and a zero gas e.g. N<sub>2</sub> or CO<sub>2</sub> with a purity of at least 99.9%. For the InPro 6970 i a purity of 99.99% is required.

Sensor type	Type of calibration		
	1-point (Air)	2-point	Process
InPro 6860i	•	•	•
InPro 6870i	•	•	•
InPro 6960i	•	•	•
InPro 6970i	–	•	•

In contrast to amperometric sensors, a 1-point zero calibration is not sufficient for high accuracy over the whole measurement range and is thus not possible (only for InPro 6970i).

**Note: A calibration of the optical sensor needs to be performed with a digital transmitter or iSense.**

**Note: To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100%. If not, the sensor needs a new calibration. Please take into account the correct air pressure and humidity. Small deviations in air ( $\pm 3\%$ ) are due to differences in humidity and process pressure settings. The sensor calculates for 100% humidity if it is set to dissolved oxygen measurement.**

**General remarks:**

- **For calibration in gas (air), the OptoCap must be dry**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the settings for **oxygen saturation** of the calibration is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.
- **Calibration in a fermenter should be performed after sterilization as process calibration (see chapter 5.3.5).**
- **Make sure that all other parameters, such as temperature and pressure, are constant.**
- **When sensor is removed from process, please let the sensor acclimatize to the temperature of the calibration environment before calibrating otherwise especially the zero point calibration is not correct**
- Calibration always needs accurate pressure and temperature measurement. Only process scaling is independent of those parameters (see chapter 5.3.5).
- Make sure that the correct calibration pressure, humidity and salinity values are set in the transmitter before the calibration is started.
- Please refer also to the transmitter manual for detailed informations

### 5.3.2 Factory calibration

The sensor is delivered pre calibrated and ready for use.

The factory calibration data are stored in the sensor and can not be changed by the user. During this calibration all sensor specific parameters are determined.

For continuous applications, we recommend **periodic recalibration in line with your requirements on accuracy, the type of process in operation and your own experience**. The frequency of the need for re-calibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

### 5.3.3 Single point calibration (Slope or Process Calibration)

For most applications, a single point calibration should be sufficient, as long as not the whole measuring range of the sensor is used.

By carrying out a single point calibration, the factual phase at the desired oxygen value e.g. at hundred percent oxygen (phi 100) of the sensor can be established. The corresponding calibration curve is calculated.

### 5.3.4 Slope calibration (InPro 6870 i and InPro 6960 i)

For most applications where oxygen levels between 10 % and 200 % Air saturation are measured a single point slope calibration is sufficient to reach the required accuracy.

The calibration medium can be either air or a calibration gas with known O<sub>2</sub> concentration or water with a known oxygen concentration.

Before starting the calibration in gas, the correct pressure and the correct humidity have to be set in the transmitter.

**Note: Wrong pressure values are the most common reasons for bad measurement accuracy.**

E.g. 50 mbar difference between the ambient pressure and the value set in the transmitter result in 5% measurement error at air.

For calibration in gas it is important that the temperature reading of the sensor is stable and represents the real gas temperature.

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

### 5.3.5 Process calibration

A process calibration is needed in a situation when the sensor can not be removed from the process.

For detailed information please refer also to the transmitter manual.

Two different routines for process calibration are possible:

- **Process calibration**
- **Process scaling**

**Process calibration** is performed when a reliable control value is available and process pressure is known. Process pressure is only needed if the system is measuring in saturation (% air or % O<sub>2</sub>) or gas (ppm gas) units. During this calibration the phase values of the calibration curve are adjusted.

**Process scaling** is performed mainly in biopharma applications after sterilization (autoclaving) when the user desires to set the system to an initial value. During this calibration the phase values of the sensor are not adjusted, only the displayed values and the nA outputs are rescaled to the desired value.

 **Note:** For process calibration the operator can use either the process pressure or the calibration pressure, depending on how the reference value is taken.

After the sensor signal has stabilized, the complete measurement system can be calibrated to the desired variable, e.g. % air, % O<sub>2</sub>, ppm or ppb (see instruction manual for the transmitter).

 **Note:** For this type of calibration an accurate reference value and correct pressure settings are essential.

### 5.3.6 Dual point calibration

To receive a maximum accuracy of the measured values over the full measuring range, a dual point calibration is required.

A dual point calibration is required after replacement of the OptoCap.

By carrying out a dual point calibration both phase angles at zero oxygen (phi 0) and at hundred percent oxygen (phi 100) of the sensor can be established.

Point 1: Slope correction (with air or other calibration media with known O<sub>2</sub> value)

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Point 2: Zero point

After the sensor signal has stabilized, the sensor can be calibrated to the 0% value of the desired measurable variable, e.g. 0% air, 0.0% O<sub>2</sub>, or 0 ppm at 25 °C/77 °F (see instruction manual for the transmitter).

**Note: Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.995 %.**

### 5.3.7 Calibration when connected with analog signal

As soon as a sensor is connected with an analog communication, only rescaling of the sensor output is available.

The procedures can be performed the same way as with amperometric sensors:

1. Slope correction
2. Offset correction

These procedures do not replace a real calibration of an optical sensor.

Before installing the sensor in a fermenter a one point air calibration (see chapter 5.3.2) should be performed using iSense Asset Suite or iSenseLight or a digital transmitter.

A two point calibration (see 5.3.6) should be performed after each OptoCap replacement and depending on the process stress after 5 to 10 batches.

### 5.3.8 Reset to factory calibration

In case of a wrong calibration of an optical sensor, e.g. by using wrong calibration values during slope or process calibration, a new 2-point calibration is necessary. If it is not possible to perform a good 2-point calibration a reset of the calibration data is possible followed by a good 1 point air calibration (see instruction manual for the transmitter).

## 6 Maintenance

**Note: All maintenance work can be done without any tools.**

### 6.1 Sensor inspection

#### 6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the OptoCap optically for signs of damage. The OptoCap must be intact and clean. Dirty surface should be wiped clean using a soft, moist tissue.

**Attention! Do not use any cleaning agents containing alcohol or any solvents. This could damage the sensor.**

#### 6.1.2 Testing the sensor with the transmitter

If the measured values differ from the expected value, a air calibration should be performed.

##### **Appropriate phase values after a correct calibration:**

Depending on the age of the OptoCap the phase values typically decrease over time compared with a new OptoCap (see table). Note that for the InPro 6970 i the phase at air increases.

Sensor type	New OptoCap		Limit for old OptoCap	
	Phi0	Phi100	Phi0	Phi100
InPro 6860 i	60°±2°	32°±3°	<50°	<15°
InPro 6870 i	60°±2°	32°±3°	<50°	<15°
InPro 6960 i	68°±3°	30°±2°	<55°	<15°
InPro 6970 i	82°±3°	13°±2°	<65°	>18°

The OptoCap needs to be replaced if the phase values exceed the limits.

The phase values of the sensor are stored in the calibration history. The actual phase value can be checked in the "Calibration – Verify" menu or with iSense.

If after such procedures the above mentioned values are still not reached, replace the OptoCap. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Zero oxygen measurement can be done by using CO<sub>2</sub> or nitrogen (N<sub>2</sub>), alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 5 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1%.

If after such procedures the above mentioned values are still not reached, replace the OptoCap. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Many sample media contain volatile substances which, even at very low concentrations, have a clearly perceptible smell. Similarly to oxygen, these substances are able to invade the OptoCap. Accordingly, they become noticeable when interacting with the pigment. In most cases, such substances have absolutely no influence on the measuring properties of the sensor.

### 6.1.3 ISM

#### **DLI: Dynamic Lifetime Indicator**

The DLI provides information about the remaining lifetime of the OptoCap. As long as the DLI is above zero days the system is within the specified accuracy after a calibration. If the DLI is zero the OptoCap needs to be replaced.

Contributing factors for ageing of the OptoCap:

- number of measurements
- temperature during measurement
- oxygen concentration during measurement
- number of CIP cycles
- number of SIP cycles
- number of autoclavings

The DLI is calculated in two different ways.

Continuously: With the above parameters an actual sensor stress is calculated. With each measurement the sensor load is increased. The accumulated sensor load divided by the elapsed time is the basis of the calculation of the remaining lifetime.

During calibration: The phase values are compared to the phase values of the last calibration. Using the above calculated sensor load and the elapsed measurement time since the last calibration, the remaining lifetime of the OptoCap is calculated. The calculation after a calibration gives a higher accuracy of the DLI compared to the continuous calculation. Thus the DLI value can be significantly different after a calibration.

**Note:** For a correct DLI calculation an accurate calibration is essential.

#### **ACT: Adaptive Calibration Timer**

The ACT provides information as to when the next calibration is required to ensure measurements will remain within the specified accuracy. This calculation is based on the DLI information.

#### **Calibration history**

The last four calibrations and the factory calibration data are stored in the sensor memory. These data can be read out with a transmitter or with the iSense Asset Suite software.

The calibration history gives valuable information regarding the quality of the calibration and the ageing of the OptoCap.

## 6.2 Replacing the OptoCap

To replace the OptoCap you first have to unscrew the cap-sleeve.

**Attention!** If the cap sleeve is detached, take care to the inner part of the sensor shaft. Damage and soiling of the inner parts and the optical fiber may influence the Signal or destroy the sensor. Small soiling can be removed with a lint-free cloth.

When changing the OptoCap, please observe the following instructions:

**Attention!** Make sure that this maintenance step is carried out in **clean place**.


**Attention!** For the **InPro 6860 i** and **InPro 6970 i** the cap sleeve and the OptoCap are together in one piece. To replace the OptoCap, the whole unit needs to be replaced.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the OptoCap.
3. Place the new OptoCap on top of the inner body of the sensor shaft.
4. Carefully slip the cap sleeve over the fitted OptoCap and screw it down. The cap sleeve must be clean and dry.
5. After each exchange of the OptoCap, the DLI must be resetted manually using the transmitter or iSense.
6. After each exchange of the OptoCap, the sensor has to be recalibrated with a dual-point calibration.

**Attention!** The quality of this calibration is critical for sensor performance and accuracy of the diagnostics.

### Replacement of the OptoCap

**InPro 6870 i**

**InPro 6860 i**  
**InPro 6960 i**  
**InPro 6970 i**

## **7 Storage**

For storage, the sensor should be clean and dry. The protection caps have to be placed on the sensor and the cable connectors. If the sensor shaft is stored without the sensor head, the maintenance cap has to be used to protect the fiber optics.

## 8 Product specification

### 8.1 Certificates

Each sensor is delivered with a set of **3.1 certificates** (complying with EN 10204.3/1.B).

All wetted metal parts (sensor shaft, cap sleeve and OptoCap) are identified with a engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and OptoCap) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320: 1992).

### 8.2 Specifications

	<b>InPro 6860i/6870i/ 6960i/6970i</b>	
<b>Measurement principle</b>	optical	
<b>Working conditions</b>		
Permissible pressure range during measurement	InPro 6860i/6870i: 0.2 ... 6 bar absolute 2.9 ... 87 psi absolute InPro 6960i: 0.2 ... 9 bar absolute 2.9 ... 130.5 psi absolute InPro 6970i: 0.2 ... 12 bar absolute 2.9 ... 174 psi absolute	
Mechanical pressure resistance	InPro 6860i/6870i: max. 6 bar absolute max. 87 psi absolute InPro 6960i/6970i: 12 bar absolute 174 psi absolute	
Permissible temperature range during measuring	InPro 6860i/ InPro 6870i: 0 ... 60 °C 32 ... 140 °F	InPro 6960i/ InPro 6970i: 0 ... 40 °C 32 ... 104 °F
Mechanical temperature resistance of the sensor shaft	InPro 6860i: -20...140 °C -4...266 °F steam-sterilizable autoclavable	InPro 6870i/ InPro 6960i/ InPro 6970i: -20...121 °C -4...250 °F steam-sterilizable
<b>Analog/digital Interface</b>		
nA Output: (InPro 6860i /../ nA)	60 nA at air (25 °C, 1013 mbar)	
mA Output: (InPro 6860i /../mA)	4 to 20 mA (active)	
Analog temperature Output:	NTC 22 kΩ	
Digital Interface:	Modbus RTU on RS485 - Baudrate: 4800,9600, 19200 (default), 38400, 56600, 115000 - Device address: 1 (default to 32)	
<b>Power requirements</b>	U <sub>min</sub> = 19.5 VDC U <sub>max</sub> = ≤25 VDC P <sub>max</sub> = 0.75 W → 33 ... 40 mA	

**Sensor performance**

Operating range	InPro 6860i/6870i/6960i: 8 ppb to 60 % O <sub>2</sub> saturation InPro 6970i: 2 ppb to 2000 ppb
Accuracy in aqueous media	InPro 6860i/6870i/6960i: ≤ ±[1%+8 ppb] InPro 6970i: ≤ ±[1%+2 ppb]
Response time at 25 °C (77 °F) (air > N <sub>2</sub> )	t <sub>98%</sub> < 20 s, < 70 s (InPro 6860i)

**Design features**

Temperature compensation	automatic with built-in RTD
Cable connection (digital)	data cable 5 pin, VP8 (InPro 6860i)
Wetted O-rings	EPDM FDA approved (other materials on request)
O <sub>2</sub> selective membrane material	silicone, PTFE (InPro 6860i)
Wetted sensor parts	s/steel 1.4404, [AISI 316L] PPS with material certificate 3.1 (other materials on request)

**Surface roughness of wetted s/steel parts**

(EN 1320: 1996)	Ra < 0.4 μm [16 μin]
-----------------	----------------------

**Certificates (MaxCert™)**

Quality (final inspection certificate)	yes
FDA/USP Class VI	yes
Material certificate 3.1	yes
Surface finish certificate 2.1	yes
ATEX Certificate	InPro 6860i/12/xxx/ma only


## 9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

### 9.1 Sensors

#### Designation

InPro 6860i/12/\_\_\_/\_\_\_

Analog signal (nA/mA)  
 Immersion length (a)  
 Sensor diameter 12 = 12 mm

Sensor	Order no.
InPro 6860i/12/120 mm/nA (4.7")	30 014 100
InPro 6860i/12/220 mm/nA (8.7")	30 014 101
InPro 6860i/12/320 mm/nA (12.6")	30 014 102
InPro 6860i/12/420 mm/nA (16.5")	30 014 103
InPro 6860i/12/120 mm/mA (4.7")	30 129 734
InPro 6860i/12/220 mm/mA (8.7")	30 129 735
InPro 6860i/12/320 mm/mA (12.6")	30 129 736
InPro 6860i/12/420 mm/mA (16.5")	30 129 737
InPro 6870i/12/120 mm (4.7")	52 206 380
InPro 6870i/12/220 mm (8.7")	52 206 381
InPro 6870i/12/320 mm (12.6")	52 206 382
InPro 6870i/12/420 mm (16.5")	52 206 383
InPro 6960i/12/120 mm (4.7")	52 206 500
InPro 6960i/12/220 mm (8.7")	52 206 501
InPro 6960i/12/320 mm (12.6")	52 206 502
InPro 6970i/12/120 mm (4.7")	52 206 393
InPro 6970i/12/220 mm (8.7")	52 206 394
InPro 6970i/12/320 mm (12.6")	52 206 395

### 9.2 Accessories

Accessories	Order no.
iSense Asset Suite	52 900 336
CalBox	52 300 400
iLink RS485	52 300 399
iLink RS485-VP (InPro 6860i)	30 014 134
Data cable (5 pin) for InPro 6870i/6960i/6970i Temperature range -30 ... 80 °C (-22 ... 176 °F)	
2 m (6.6 ft)	52 300 379
5 m (16.4 ft)	52 300 380
10 m (32.8 ft)	52 300 381
15 m (49.2 ft)	52 206 422
25 m (82.0 ft)	52 206 529
50 m (164 ft)	52 206 530

<b>Accessories</b>	<b>Order no.</b>
Power supply for InPro 6860i	30 014 119
Power adapter ODO T82	30 083 984
Power adapter ODO VP6	30 083 985
Ex i power supply 9143/10-244-060-20s	30 215 910
Signal converter 4/20 Hart (active/passive)	30 212 436
Cables for InPro 6860 i	
Temperature range –30 ... 80 °C ( –22 ... 176 °F)	
VP8 -ST (digital and analog connection)	
1 m (3.3 ft)	52 300 353
3 m (9.9 ft)	52 300 354
5 m (16.4 ft)	52 300 355
10 m (32.8 ft)	52 300 356
15 m (49.2 ft)	52 206 457
20 m (65.6 ft)	52 206 558
35 m (114.8 ft)	52 206 559
VP6-ST (only analog connection)	
1 m (3.3 ft)	52 300 107
3 m (9.9 ft)	52 300 108
5 m (16.4 ft)	52 300 109
10 m (32.8 ft)	52 300 110
15 m (49.2 ft)	52 300 141
20 m (65.6 ft)	52 300 144
35 m (114.8 ft)	52 300 184
VP6 cables for IP 6860 i preconfigured with connector	
Connector cable BNC, 1 m (3.3 ft)	30 032 730
Connector cable BNC, 3 m (9.9 ft)	30 032 731
Connector cable LEMO, 1 m (3.3 ft)	30 032 732
Connector cable LEMO, 3 m (9.9 ft)	30 032 733
Connector cable Lumberg, 1 m (3.3 ft)	30 032 734
Connector cable Lumberg, 3 m (9.9 ft)	30 032 735

### 9.3 Spare parts

<b>Spare parts</b>	<b>Order no.</b>
OptoCap BT01 (InPro 6870i)	52 206 225
OptoCap BT02 T (InPro 6860 i), incl. O-ring	30 018 857
OptoCap BR01 (InPro 6970 i)	52 206 403
OptoCap BW01 (InPro 6960 i)	52 206 509
O-ring set	52 206 252
Cap sleeve	52 206 232
Maintenance cap	52 206 251

### 9.4 Recommended transmitters

<b>Transmitter</b>	<b>Order no.</b>
M400, Type 2	52 121 349
M400, Type 3 (InPro 6960 i/InPro 6970 i)	52 121 350
M400 FF	30 026 616
M400 PA	30 026 617
M800 1-channel	30 026 633
M800 2-channel	52 121 813
M800 4-channel	52 121 853

## 9.5 Recommended housings

<b>Static housings</b>	<b>Order no.</b>
InFit 761	–
<b>Retractable housings</b>	
InTrac 777 e	–
InTrac 797 e	–
<b>Sensor protection</b>	
Retrofit kit optical	52 403 811

**Note:** The housings are available in different versions. Please contact your distributor to get the right ordering information.

## 10 Theory of the optical oxygen measurement

### 10.1 Introduction

Optical oxygen measurement is a non-invasive method. No electrochemical reaction occurs during measurement.

### 10.2 Principle

In contrast to the amperometric and potentiometric method the optical measurement is not based on a chemical reaction and current measurement.

A chromophore in the Sensor is illuminated with blue light. The chromophore absorbs this energy and is transferred to a higher energy level. A part of the energy is transferred as heat. After a short time the chromophore emits a red fluorescence light and returns to its ground state.

If an oxygen molecule collides with the chromophore in its excited state, the energy can be transferred to oxygen (Dynamic Quenching). In this case, no fluorescence light is emitted. Oxygen itself can transfer this energy as heat without light emission.


The emission of fluorescence light is therefore dependent on the oxygen partial pressure at the chromophore layer.

The emission of fluorescence shows a short time delay to the excitation. This time delay can be measured, if the excitation light is being modulated. In this case, the fluorescence shows the same modulation.

The phase shift between excitation and fluorescence decreases with increasing oxygen concentration.

The oxygen concentration is calculated and digitally transferred to the transmitter.

### 10.3 Principle of the design of the optical oxygen sensor

## 10.4 Temperature

The optical oxygen sensor contains a RTD and the temperature is compensated automatically.

## 10.5 Dependence on flow

The optical oxygen sensors show a strongly reduced flow dependence compared to the amperometric sensors.

The reason is that these optical sensors do not reduce oxygen. No change in the oxygen concentration occurs at the sensor.

## 10.6 Oxygen partial pressure – oxygen concentration

The sensor signal depends on the oxygen partial pressure and the O<sub>2</sub> permeability of the membrane – but not on the O<sub>2</sub> solubility in the solutions. The oxygen concentration in mg O<sub>2</sub>/L (CL) cannot therefore be determined directly with an electrode.

According to Henry's law the oxygen concentration is proportional to its partial pressure (PO<sub>2</sub>).

$$CL = pO_2 \cdot a$$

$a$  = solubility factor

If " $a$ " is constant and known, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

Medium, sat. with air	Solubility at 20 °C (68 °F) and 1013 mbar (760 mm Hg)
Water	9.2 mg O <sub>2</sub> /L
4 mol/L KCl	2 mg O <sub>2</sub> /L
50% Methanol-water	21.9 mg O <sub>2</sub> /L

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions. Thus, determination of the oxygen concentration is only possible with constant and known solubility factors " $a$ ".

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

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01/2016 © Mettler-Toledo GmbH
- France**                    **Mettler-Toledo Analyse Industrielle S.A.S.**
- Germany**                **Mettler-Toledo GmbH**
- Switzerland**           **Mettler-Toledo (Schweiz) GmbH**
Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, CH-8902 Urdorf