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Ex classification  **II 1/2G Ex ia IIB T6/T5/T4/T3 Ga/Gb in accordance with test certification SEV 10 ATEX 0137X and CML 22 UKEX 2384X**

1 Introduction

In accordance with Directive RL2014/34/EU (ATEX 114)¹⁾ Annex I, the InPro™ 7100i/12/***/**** conductivity sensor is classified in Device Group II Category 1/2 G, which according to Directive RL 99/92/EC (ATEX 137)²⁾ may be used in zones 0/1 or 0/2 or 1 or 2 as well as gas groups IIA and IIB, which have potentially explosive atmospheres due to combustible materials in the range of temperature classes T1 to T6.

It must be installed and used in compliance with the requirements of EN 60079-14.

The ISM™ measurement circuit is part of a common intrinsically safe system and is connected to and operated with a separately certified transmitter.

As part of an intrinsically safe system, the ISM measurement circuit has safe galvanic isolation at the transmitter from circuits that are not intrinsically safe up to the peak value of a nominal voltage of 375 V and from earthed parts up to the peak value of a nominal voltage of 30V.

2 Rated parameters

Measuring and supply circuit:

With type of protection intrinsic safety Ex ia IIB.
Only for connection to certified intrinsically safe circuits

Maximum values:

$U_i \leq 16 \text{ V}$, $I_i \leq 30 \text{ mA}$, $P_i \leq 50 \text{ mW}$

$L_i = 0.002 \text{ } \mu\text{H}$ (effective internal inductance)

$C_i = 0.0015 \text{ } \mu\text{F}$ (effective internal capacitance)

3 Special conditions

- 3.1 The maximum permissible process temperatures with corresponding temperature classes are shown in the following table:

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR)

Ex classification  **II 1/2G Ex ia IIB T6/T5/T4/T3 Ga/Gb in accordance with test certification SEV 10 ATEX 0137X and CML 22 UKEX 2384X**

| Temperature class | Max. permissible process temperature |
|--------------------------|---|
| T6 | 62 °C |
| T5 | 74 °C |
| T4 | 102 °C |
| T3 | 130 °C |

- 3.2 The permissible external ambient temperature range in the area of the connector head is –20 °C to +50 °C.
- 3.3 The conductivity sensor InPro 7100 i/12/***/**** must only be used in suitable process connectors of METTLER TOLEDO or other manufacturers in potentially explosive areas.
- 3.4 The capacitance and inductance of the connecting cable must be taken into account in the dimensioning.
- 3.5 In areas with potentially explosive atmospheres, InPro 7100 i/12/***/**** conductivity sensor can be used in/with the InFit™ 76*-*** or InTrac™ 7**-*** fittings or in/with other suitable fittings.
- If necessary, this fitting must be included in the periodic pressure testing of the system or facility.
- 3.6 The independent process fitting used for the installation of the conductivity sensor must be electrically connected to the potential equalization system of the system or facility, and if necessary it must be included in the periodic pressure testing of the system or facility.
- 3.7 To ensure safe working conditions in areas with potentially explosive atmospheres, the minimum conductivity of the media must always be greater than 1 nS/cm.

InPro, InFit, InTrac and ISM are trademarks of the METTLER TOLEDO Group.

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All rights reserved