# InPro 3030

## pH Combination Electrode

### Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

InPro 3030 pH combination electrode  
52 000 362 B

# InPro 3030

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
## Instruction Manual

### Contents

1	Introduction .....	15
2	Safety instructions .....	15
3	Product description .....	15
4	Installation and commissioning .....	15
5	Operation .....	16
5.1	Calibration of the electrode and pH transmitter .....	16
5.2	Sterilization of the electrode .....	16
6	Maintenance.....	16
7	Trouble-shooting (sluggish response time, insufficient sensitivity etc.) .....	17
8	Environmental protection.....	17
9	Warranty.....	17

InPro and Viscolyt are trademarks of the METTLER TOLEDO Group.

## 1 Introduction

The METTLER TOLEDO pH electrode type InPro™ 3030 is a low-maintenance, autoclavable and pressure-resistant measuring electrode with a gel-type reference electrolyte. The reference element incorporates a silver-ion trap to prevent silver sulfide precipitation on the diaphragm. Please read these instructions carefully before putting the electrode into operation, in order to ensure faultless operation.

The InPro 3030 electrode is mainly designed for the measurement of pH value in biotechnical processes at moderate pressures (< 6 bar). It can be used under industrial fermentation, pilot and laboratory conditions, and is particularly intended for those processes where intermediate sterilization or autoclaving is indispensable.

## 2 Safety instructions

The electrode is intended only for the use described in Chapter 1. The electrode consists partly of glass. If an electrode is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn.

## 3 Product description

The inscription on each electrode contains the following information:

<b>METTLER TOLEDO</b>	Manufacturer of the electrode
<b>InPro 3030/yyy</b>	Type designation/yyy = shaft length in mm
<b>combination pH</b>	Type of electrode
pH 0...14	pH measurement range
0...80 °C/140 °C	Temperature range for operation/sterilization
Order No. 52 00X XXX	Order number

Additionally, the electrode is supplied with a serial number on the connector cap for identification possibilities for each electrode produced.

## 4 Installation and commissioning

1. On unpacking, check the electrode for mechanical damage of the pH sensitive glass and shaft. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and briefly rinse the electrode with de-ionized water. After rinsing, the electrode should only be dabbed dry with a tissue. Do not rub the pH-sensitive glass, since this can lead to electrostatic charging and sluggish response times.

3. Check the space behind the pH-sensitive glass for air bubbles and remove any air bubbles by shaking the electrode slightly in a vertical direction.
4. Place the electrode into the housing as described in the instructions for the housing.
5. Connect the electrode and the pH transmitter with the corresponding cable (AK9 for InPro 3030). Observe the connections scheme that accompanies the cable or the transmitter.

## 5 Operation

### 5.1 Calibration of the electrode and pH transmitter

For the InPro 3030 electrode a 2-point pH-calibration is recommended. Before calibration, first remove the watering cap. To speed up the response time of the electrode after long transport and storage periods, and before calibration, dip the electrode into buffer pH 7, buffer pH 4 and buffer pH 7 again, each time for about one minute. The electrode is then dipped successively into two different buffer solutions with given pH reference values and the transmitter calibrated to these buffer values. Normally you start with buffer pH 7 for determining the 0-point and then buffer pH 4 or 9 for determining the slope. In the case of microprocessor-based transmitters with buffer recognition, the order of the two chosen buffer solutions is irrelevant. Whenever possible, a temperature sensor should be connected to the transmitter and «automatic temperature measurement» selected there.

Please refer to the instruction manual of the transmitter for further details.

### 5.2 Sterilization of the electrode

The InPro 3030 electrode can be sterilized either through autoclaving or through in-situ sterilization with steam or superheated fermentation medium. Sterilization temperatures above 140 °C are not allowed, and excessive sterilization time (for instance t>30 min at 140 °C) could notably shorten the life time of the electrode. Repeated sterilization and/or autoclaving of the electrode can lead to a slight coloration of the electrolyte. This color change does not effect the functioning of the electrode.

## 6 Maintenance

1. Carefully rinse the electrode tip and diaphragm with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the electrode is not in operation, store it with electrode tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt™ electrolyte.

3. If the electrode is stored mounted in its housing, the same rule as described in point 2 applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts
4. If a electrode is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with de-ionized water or alcohol, and afterwards dry carefully.

## 7 Trouble-shooting (sluggish response time, insufficient sensitivity etc.)

The glass membrane can be reactivated using a HF solution (order no. 51 319 053). The electrode tip (only pH sensitive glass part) is immersed into the HF solution for a period of about 5 minutes, is then immediately and thoroughly rinsed with de-ionized water and thereafter stored in buffer solution pH 4 (see Chapter 6) up to the diaphragm for 12 hours.

Possible protein precipitation and blockage of the diaphragm can be treated with our special electrode cleaner (order no. 51 340 068). However, use of the diaphragm cleaner (order no. 51 340 070) is normally not necessary for this electrode type, as the built-in silver-ion trap prevents the formation of silver sulfide contamination.

Also, check the connector for possible traces of moisture. If necessary, clean well using de-ionized water or alcohol, and afterwards dry carefully.

## 8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

## 9 Warranty

In the event of manufacturing faults, 12 months warranty from date of production is granted.

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本表依据 SJ/T 11364 的规定编制。本产品符合以下标志规范：  
**Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:**

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○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

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×：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

## Notes

---


For addresses of METTLER TOLEDO  
Market Organizations please go to:

**[www.mt.com/contacts](http://www.mt.com/contacts)**

**METTLER TOLEDO Group**  
Process Analytics  
Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

Subject to technical changes  
© 04/2023 METTLER TOLEDO  
All rights reserved  
Printed in Switzerland. 52 000 362 B

Management System  
certified according to  
ISO 9001 / ISO 14001

[www.mt.com/pro](http://www.mt.com/pro)

\* 5 2 0 0 0 3 6 2 B \*

**UK  
CA**

**CE**

**EAC**