Type HA405-DPA-SC pre-pressurized gel-reference pH Electrodes

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction manual	8
es	Manual de instrucciones	10
fr	Instructions d'utilisation	12
fi	Käyttöohje	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือการแนะนำ	34
zh	说明书	36

HA405-DPA-SC pH Electrodes
53 500 086 E

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Instructions for Type HA405-DPA-SC pre-pressurized gel-reference pH Electrodes

1. Introduction

The METTLER TOLEDO Type HA405-DPA-SC combination pH electrode is a low-maintenance, pre-pressurized electrode with a gel-type reference. It has a specially formulated HA (high alkalinity) pH membrane to minimize sodium error in extremely alkaline media. The reference element incorporates a silver-ion trap to prevent silver sulfide contamination of the reference junction. The gel-type reference is also available in a combination ORP and as a separate reference electrode.

2. Preparation

After unpacking, examine the electrode for any mechanical damage.

METTLER TOLEDO electrodes are factory tested and are shipped ready-to-use. To maintain factory pressurization, the Type HA405-DPA-SC electrodes have a silicone seal covering the ceramic diaphragm; this must be removed with the enclosed blade prior to use. Remove the hydration cap and remove any electrolyte crystals by rinsing the electrode with distilled water. Pat dry; do not rub as this may build up a static charge and lengthen the response time.

Eliminate any air bubbles around the internal element in the pH-sensitive membrane.

3. Operation

METTLER TOLEDO pH electrodes can be used with most pH/mV meters. Follow manufacturer's instructions for use.

A two-point calibration is recommended, using fresh buffers of known pH values which bracket the sample pH; one buffer should be close to pH 7 to set the zero point. For fast routine analysis, a standardization using one buffer at a value near the sample pH can be done.

4. Maintenance

When properly cared for, METTLER TOLEDO electrodes should provide accurate results over thousands of uses.

Storage: Electrodes should be rinsed thoroughly after use and stored in 3M KCl or pH 4 buffer solution with added KCl. DO NOT store the electrode dry or in distilled or tap water; this dehydrates the pH membrane and the electrode must be rehydrated before use by immersing for several hours in 3M KCl.

Refilling: Type HA405-DPA-SC gel-reference electrodes are not refillable.

Protein Contamination: Clean the electrode with pH Electrode Cleaner for Proteins following instructions on the bottle.

Reactivating pH membrane: Use Reactivating Solution following directions CAREFULLY.

Servicing a clogged or dry reference junction: Immerse the electrode in a beaker of 3M KCl. Heat the solution to approximately 70 °C (160 °F) and allow electrode to soak in this solution as it cools to room temperature.

5. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

6. Troubleshooting

No response (symptom)

- Poor meter connection
 - Check connectors at the meter and at the electrode
- Cracked glass
 - Replace electrode

Slow response (symptom)

- Dehydrated pH membrane
 - Rehydrate (see Maintenance)
- Clogged reference junction
 - Soak in warm KCl (see Maintenance)
- Protein contamination
 - Use pH Electrode Cleaner (see Maintenance)
- Depleted gel layer
 - Use Reactivating Solution (see Maintenance)

Noise (symptom)

- Poor meter connection
 - Check connector
- Clogged reference junction
 - (See above)

Low slope (symptom)

- Dehydrated pH membrane
 - Rehydrate (see Maintenance)
- Contaminated membrane
 - (See above)

For further assistance, call our Customer Service Center at process.hotline@mt.com.

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All rights reserved