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
# Instructions for Redox and Metal Electrodes

## Contents

1	Type summary.....	15
2	Starting up .....	15
2.1	Inspection .....	15
2.2	Preparation .....	15
2.3	Measuring.....	15
3	Maintenance.....	16
3.1	Reference electrolyte .....	16
3.2	Metal .....	16
3.3	Diaphragm.....	16
4	Notes.....	17
5	Safety remarks.....	17
6	Trouble-shooting .....	17
7	Factors that shorten life-time.....	17
8	Environmental protection.....	17
9	Literature.....	17

## Description of the electrodes

- A Screw cap
- B Connector
- C Filler opening
- D KCl reference electrolyte
- E Ag/AgCl reference system
- F Ceramic diaphragm
- G Metal electrode
- H Electrolyte reservoir

### Markings

- Order number on connector or shaft
- Temperature range, reference electrolyte and metal on shaft

### Application

- Measurement of redox potential in ion-containing liquids

Argenthal, Xerolyt and Friscolyt are trademarks of the METTLER TOLEDO Group.

## 1 Type summary

Type	Designation	Relevant section
Combination redox electrodes	Pt 4805-S7, Au 4805-S7	all
	Pt 4805-DPA-S8, Pt 4805-DXK-S8, Pt 4805-DPAS-SC-K8S	2, 3.2/3.3, 4–8
	Pt 4865-..., Pt 4805-KN-S7, Pt 4825-60-S7	all
	Pt 4805-60-88TE-S7	2, 3.1/3.2, 4–8
Metal electrodes	(requiring reference electrode, e.g. 363-S7)	
	Pt 805-S7, Au 805-S7	2, 3.2, 4–8


## 2 Starting up

### 2.1 Inspection

On unpacking the electrode, inspect the shaft for any sign of mechanical defects.


### 2.2 Preparation

- a Take off the watering cap
- b Rinse with distilled water
- c Remove rubber cap
- d Condition and check operation of the redox electrode with 250 ml of METTLER TOLEDO redox buffer 220 mV (pH 7, No. 51 340 065) at 25 °C.
  - Connect electrode to pH meter (set to mV range)
  - Immerse electrode in METTLER TOLEDO buffer
  - Reading on mV meter must agree with the figure on the flask (permitted discrepancy  $\pm 5$  mV)
  - Discrepancies of more than 5 mV must be allowed for in the measured result
  - If discrepancy exceeds 10 mV, check the reference electrode (see chap. 3 below)

### 2.3 Measuring

- e Rinse with distilled water
- f Connect electrode to a pH/mV meter or transmitter. Dip electrode in measuring medium

Industry: Observe instructions for housing

### Important

- Surface of reference electrolyte higher than surface of measured medium
- Diaphragm must be immersed

### General

- Observe instructions for pH/mV meter or transmitter

### g Interpretation of measurement:

The reading on the pH/mV meter or transmitter depends on the reference system used.

Conversion to redox voltage related to standard hydrogen electrode (SHE)

$$E_H = E + E_{ref}$$

$E$  = Instrument reading

$E_H$  = Standard voltage of reference electrode (correction)

$E_{ref}$  = Redox voltage rel. to SHE

Table 1

Reference electrolyte	3 mol/l KCl	Friscolyt™	Viscolyt™
$E_{ref}$ [mV] at 25 °C	207.0	207.0	207.0

**Note:** With changes of temperature,  $E_{ref}$  varies differently for each reference system.

## 3 Maintenance

### 3.1 Reference electrolyte




**h** Discoloration precipitates: suck electrolyte by siphon and refill.

**i** Level must be at least 2 cm below filler inlet: refill up to max. 1 cm below filler inlet.

### 3.2 Metal

**k** Clean and remove deposits with toothpaste or 2 mol/l HCl.

**l** Rinse with distilled water.

### 3.3 Diaphragm

Clean diaphragm with appropriate solution

Table 2


Nature of deposit	Cleaning agent	Order No.
Proteins	Pepsin/HCl	51 340 068
CaCO <sub>3</sub>	Dilute HCL	
Silver sulphide Ag <sub>2</sub> S	Diaphragm cleaner	51 340 70
Oil, grease	Acetone or ethanol	

### Storage

**m** In reference electrolyte, in 3 mol/l KCl or Friscolyt.

Do not store electrode in distilled water or deionized water!

#### 4 Notes

- Measurement in mV absolute, make note of temperature.
- Temperature compensation of absolute voltage readings is unnecessary
- Follow the instructions for your pH/mV meter or transmitter

#### 5 Safety remarks

- Glass: breakage, glass splinters in measured medium risk of injury
- Wear safety glasses and gloves when handling cleaning agents
- We recommend that the electrodes only be put into operation in conjunction with original products/accessories from METTLER TOLEDO. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions

#### 6 Trouble-shooting

Table 3

Symptom	Possible cause	Remedy
Reading drifts, reading unstable, response time > 1 min	Diaphragm dirty (discoloured)	Cleaning and conditioning
Voltage $U = 0$	Short circuit	Change cable, replace electrode
Voltage $U = \infty$	Broken contact	
Reading unstable, voltage $U =$	Connector socket dirty	Remove deposits mechanically

#### 7 Factors that shorten life-time

- For Xerolyt electrodes (type Pt4805-DXK): pH < 2, organic solvents, vacuum and high temperature
- Solutions containing HF
- For electrodes Pt4805-DPA/DPAS: pressure > 2.5 bar

#### 8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

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All rights reserved