Quick Setup Guide

M300/M4002G Quick Setup Guide

For InSUS pH Single-Use Sensors





Content

1.	Introduction	3
2.	Installation	3
3.	Data Entry of Pre-Calibrated pH Single-Use Sensor: Slope & Zero Point (Offset)	4
4.	Process Calibration	6
5.	pH Membrane Glass Breakage – Alarm Activation	9

1. Introduction

METTLER TOLEDO InSUS™ pH sensors are analog, gamma irradiation-sterilizable, pre-calibrated single-use pH sensors with an integrated Pt 1000 temperature probe. Please read through these instructions carefully before commissioning in order to ensure trouble-free use. Operation should be carried out only by trained personnel and staff who have read and understood the sensors' operating instructions.



Fig. 1: InSUS pH sensor (example)

Pos. Description

1 VP connector

2 Slope, Zero Point (offset) and sensor serial number

2. Installation

Connect the sensor to the pH transmitter with a VP6 cable. Observe the connection information that accompanies the cable or the transmitter and follow the wiring instructions for pH sensors with integrated **Pt 1000** temperature probe and **without solution ground.** Follow the instructions in the M300 or M400 transmitter manual to configure the measurement channel for an analog pH/ORP sensor with a Pt 1000 temperature probe.

3

3.

Data Entry of Pre-Calibrated pH Single-Use Sensor: Slope & Zero Point (Offset)

- 3.1 While the M300/M400 transmitter is in measurement mode, go to the calibration menu by selecting the calibration icon.
- 3.2 Press Calibrate Sensor



3.3 Select the measurement channel – only applicable for 2-channel transmitters.

습\CAL\Calibrate Sensor				
Chan	✓ CHAN_1			
Unit	CHAN_2			
Method	1-Point			
Options	Options			
Verify	Edit	Cal		

3.4 Press the Edit key.

\CAL \ Calibrate Sensor				
Chan	Ch1 pH/ORP Edit			
Unit	Slope	-59.100	m∿/pH	
Unit	Offset	0.0000	mV	
Metho	Temp Slope	1.0000		
Option	Temp Offset	0.0000	ĺ	
\ \	Cancel <	(1/2 > [Save	
			4	

3.5 Locate the Slope and Zero Point (Offset) values found on the sensor label (Fig. 1) or on the sensor's Quality Certificate and use these for data entry. Press the slope value and use the keys to edit the value in mV/pH.

Note: If the slope value shows as %, press the U key to change the unit to mV/pH.



3.6 Press OK and continue with the Zero Point (Offset) by pressing the offset value.

岱\CAL\Calibrate Sensor				
Chan	Ch1 pH/ORP Edit			
LL-3	Slope	-59.100	mì√/pH	
Unit	Offset	0.0000	рН	
Metho	Temp Slope	1.0000	ĺ	
Option	Temp Offset	0.0000	Ì	
	Cancel <	(1/2 > [Save	
			1	5

3.7 Edit the offset value.

Note: If the offset value shows as mV, press the U key to change the unit to pH.

Edit Offset				
7 . 1 2 pH				
1	2	3	+/-	U
4	5	6		+
7	8	9	0	ОК

3.8 Press **OK** and then **Save** to accept the data for slope and offset and to overwrite the previous calibration data.

☆ \CAL \Calibrate Sensor				
Chan	Ch1 pH/ORP Edit			
ULA	Slope	-59.100	m∿/pH	
Unit	Offset	7.1200	рН	
Metho	Temp Slope	1.0000		
Option	Temp Offset	0.0000	ĺ	
\ \	Cancel <	(1/2 > [Save	
			4	

If a warning message appears on the screen, please check the ISM /Sensor alarm settings as described in chapter 5 of this guide

Process Calibration

For highest measurement accuracy, a process calibration must be performed after the entry of the factory calibration data as described above. In a process calibration the pH value of an off-line grab sample is used to adjust the in-line measurement to that of the grab sample. This is a two-step procedure: Step 1 initiates the process calibration and stores the current pH value in the transmitter while a grab sample is taken. Step 2 is for entering the off-line value into the transmitter.

Sensor type	Minimum wetting time
InSUS 307	20 minutes
InSUS 307 XSL	120 minutes
InSUS 310	120 minutes

Table 1

Important: Prior to the process calibration, these sensors must be wetted in process liquid for a time equal to or greater than the values indicated in Table 1.

4.

4.1 While in Measurement mode, go the calibration menu, select **Calibrate Sensor** and then press **Process.**

습\CAL\Calibrate Sensor				
Chan	✓ 1-Point			
Unit	2-Point			
Method	Process			
Ortiona	Ontions			
Options				
Verify	Edit	Cal		
		ţ		

4.2 Press the **Cal** button, take a sample from the process and then press the **Enter** key to store the current measuring value.

습\CAL\Calibrate Sensor				
Chan	Ch1 pH/ORP Process			
Chair	Press "Enter" to capture the measured			
Unit	value			
Metho	7.12 рн			
\ \	Cancel 🗸			
	4			

- 4.3 To show the ongoing calibration process, "P" blinks in the measurement and menu screen if the related channel is selected in the display.
- 4.4 After determining the pH value of the sample, press the calibration icon in the measurement screen again. Enter the pH value of the sample and press **OK**.

Edit Cal Point				
7.15				
1 2 3 +/-				U
4	5	6		+
7	8	9	0	ок

4.5 Press the **Next** button to start the calculation of the calibration results.



4.6 The display shows the value for the slope and the offset resulting from the calibration. Press **SaveCal** to accept and to overwrite the previous calibration.

	•				
0	Ch1 pH/ORP Process				
	Slope	99.9 %	-59.1 mMpH		
	Offset	7.15 pH	9.9 mV		
1	Cancel	SaveCal	lack		

pH Membrane Glass Breakage – Alarm Activation

The detection of a pH membrane glass breakage can be linked to one of the transmitter's alarm relays. If set, the alarm will be activated if the glass membrane's resistance falls below $5 M\Omega$. A sensor with a broken pH membrane glass cannot be used for measurement purpose!

To set the alarm:

- 5.1 While the transmitter is in measurement mode, go to the Configuration menu by selecting the configuration icon.
- 5.2 Press "ISM / Sensor Alarm"

A CONFIG	
ISM / Sensor Alarm	•
Clean	•
Display Setup	•
Digital Inputs	•
System	•
< 214 >	1

5.3 Select the measurement channel – only applicable for 2-channel transmitters. Press **"Events"**.

	ISM / Sensor	Alarm
Options	CHAN_1	Events
Alarm Relay	#2	Normal
Delay	1	sec
Hold Relay	None	
		IJ

5.4 Activate Rg Diagnostics (glass resistance diagnostics) and press the Enter key.

Note: Do not activate the Rr Diagnostics for InSUS 307 and InSUS 310 pH sensors!

尚 10	CONFIG \ISM / Sensor Alarm	
Ontion	Events Option	
Сраон	Rg Diagnostics	
Alarm	Rr Diagnostics	
Delay		
Hold R		
		-
	L	
		ţ

5.5 Press the return arrow two times and press Yes to save the change.

CONFIG\ISM / Sensor Alarm				
Options	CHAN_1	Events		
Alarm Rel Delay Hold Relay	Save Changes f	iormal		
		5		

Note: If the Rg Diagnostics is activated, a warning indication will appear on top of the measurement screen directly after the entry of Factory Slope and Offset data. **This warning will disappear after a Process Calibration has been performed.**

		Δ	
< CHA	N_1		>
	3.91	рН	
	25.7	°C	
	200	mV	
	314	MΩ Rg	
ISM	×	1	*0

Example: Warning indication at top of the measurement screen.

Notes

The information you want is at www.mt.com/pro

The METTLER TOLEDO Process Analytics website contains a vast amount of up-to-date information on all our products and services. Content is localized for your country and tailored to suit your selections. Simple layout allows you to quickly find the information and features you are looking for.



- Learn about our most recent product developments
 Register for free webinars
- Request further information on products and services
- Obtain a quote quickly and easily
- Download our latest white papers
- Read case studies relevant to your industry
- Access buffer and electrolyte solution certificate
- and more ...

www.mt.com/pro

For more information

METTLER TOLEDO Group

Process Analytics Im Hackacker 15 CH-8902 Urdorf

Local contacts: www.mt.com/pro-MOs

Subject to technical changes © 01/2022 METTLER TOLEDO. All rights reserved UR1000en C. eVersion only MarCom Urdorf, CH