METTLER TOLEDO

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InPro X1 HLS pH Sensor

Content

1	Introduction	16
2	Safety Instructions	16
3	Product Description	16
3.1	Scope of Delivery	17
4	Installation and Commissioning	18
5	Operation: Calibration of the Sensor	19
6	Maintenance	19
7	Disposal	20
8	Warranty	20
9	Storage Conditions	20
10	Product Specifications	21
11	Drawings	23
12	Ordering Information	24
13	Accessories	24

1 Introduction

Thank you for buying the InPro™ X1 HLS pH sensor from METTLER TOLEDO Ingold. The construction of the METTLER TOLEDO InPro X1 pH sensor employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment.

Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor. If you have questions that are not answered sufficiently in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. We recommend that you only operate the sensor in conjunction with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

2 Safety Instructions

When cleaning or calibrating the sensor using solutions containing acids or alkalis, both protective spectacles and protective gloves should be worn.

3 Product Description

The METTLER TOLEDO InPro X1 HLS pH single rod measuring cell is a low maintenance, clean-in-place (CIP) resistant and pressure resistant pH sensor based on our innovative composite X-Chip™ pH-sensing technology.

The inscription on each sensor contains the following information:

METTLER TOLEDO	Manufacturer of the sensor
InPro X1 ABC-Dxxx-Eyyy	Type designation: A = Certifications B = Reference system C = Tip design D = Redox yes/no E = Cable connection xxx = Shaft materials yyy = Rod length in mm
pH 0...14	pH measurement range
X...Y°C (Z°C)	Operating temperature range (for cleaning)
Order No. XX XXX XXX	Order number

Identification and number of certificates and code for FM testing are also found on the sensor body. In addition each sensor is provided with a serial number on the red part of the plug-in head in order to enable identification. This product meets the requirements of 3A, EHEDG, EC1935/2004 and ATEX guideline 2014/34/EU. The complete declaration can be downloaded under www.mt.com/InProX1HLS.

3.1 Scope of Delivery


Each sensor is supplied fully assembled, factory tested and calibrated for correct function together with:

- Quality control certificate
- Quick setup guide
- Ex classification for pH sensors
- EC Declaration of Conformity

4 Installation and Commissioning

1. Remove the protective cap by holding the lower part and turning the upper part counter-clockwise. Rinse the sensor for a short time with deionized water. If the sensor has already been used, the sensitive part can be cleaned if necessary with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionized water. After rinsing, dab the sensor tip dry.
2. Place the sensor into the process adapter as described in its instruction manual. The InPro X1 HLS is compatible with METTLER TOLEDO InFit™, InTrac™ and InFlow™ housings. For hygienic applications, we recommend the use of an InFit761/NC type housing.

Important: InPro X1 HLS sensors must have an angle of installation of 15° or more above the horizontal. This is required due to the design of the reference system with the pre-pressurized liquid electrolyte. Additionally, in applications with little or no flow, air bubbles can become trapped at the sensor tip. Trapped bubbles lead to interference and negatively influence the pH measurement performance. This can be avoided by an installation angle 15° from the vertical. **The proper angle of installation for the InPro X1 HLS is between 15° and 75° from the vertical.**

3. Connect the sensor to the installed METTLER TOLEDO ISM™ transmitter with an AK9 cable of the desired length. Note the connections scheme that accompanies the cable or the transmitter.
4. ISM sensors from METTLER TOLEDO such as the InPro X1 HLS, support Plug and Measure and predictive diagnostics. For installation, commissioning and operation of the ISM system please observe the instructions for the transmitter and cable.

5 Operation: Calibration of the Sensor

Before a calibration, immerse the sensor for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connecting it to a METTLER TOLEDO ISM transmitter or ISMCore™ software. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. The second buffer solution (besides the pH 7.00 buffer) for the 2-point calibration should be chosen depending on the desired range of pH measurement. A clean and carefully executed calibration is decisive for the accuracy and reliability of the pH measurement. We recommend using only original METTLER TOLEDO pH buffers and solutions for calibrations and storage. For further details please consult the Accuracy Guide: Hints and Tips to Improve Accuracy In DO and pH Measurement (available on mt.com), as well as the operating instructions for the METTLER TOLEDO ISM transmitter and ISMCore software.

6 Maintenance

After each working cycle the sensor tip and diaphragm are to be carefully rinsed with a suitable solvent, such as clean water, 3M KCl, pH 4 or 7 buffer. Any drying out of residues from the process media on these parts must be avoided at all costs.

With a soft toothbrush and water any solid residues can easily be removed from the sensor. Mild washing agents can also be used. Protein contaminations on the diaphragm can be removed with METTLER TOLEDO Cleaning Solution Pepsin-HCl 6×250 ml, 51 340 069.

Once installed, the InPro X1 pH sensor is expected to be CIP cleaned. When defining the CIP procedures, care must be taken to ensure that the cleaning and disinfection chemicals used, apart from being suitable for the intended removal of contamination, must not damage the sensor materials.

Warning: Solutions containing hydrofluoric acid could damage the pH-sensitive part of the sensor and the protective titanium cap.

7 Disposal

Environmental Protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

8 Warranty

On manufacturing defects 12 months after delivery.

9 Storage Conditions

We recommend that the sensor is stored at room temperature. If the sensor is not in use it has to be stored with the sensor tip and the diaphragm immersed in a container or protective cap filled with 3 M KCl or pH 7 buffer solution, in order to avoid the X-Chip drying out.

The InPro X1 HLS should never be stored in Friscolyt electrolyte solution. Friscolyt should also never be filled in the protective cap.

If a sensor is inadvertently left in a dry condition for too long, it must be soaked in 3 M KCl or buffer solution for several hours until proper measurement function is restored.

10 Product Specifications


The InPro X1 HLS is a food-safe, in-line pH sensor that withstands CIP

It is designed to deliver consistent and continuous pH measurement in food and beverage applications. The composite X-Chip pH-sensing technology eliminates the risk of food contamination by glass fragments.

pH Range	0...14
Operating Temperature Range (for cleaning)	0...80°C (100°C)/32...176°F (212°F)
Tested Clean-In-Place (CIP) Conditions	3 % NaOH: 30 min at 90 °C (194 °F)
	H ₂ O: 30 min at 90 °C (194 °F)
	3 % HNO ₃ : 30 min at 45 °C (113 °F)
	H ₂ O: 30 min at 90 °C (194 °F)
Operating Pressure	Up to 4 bar (58 psi) overpressure
Cable Connection	K8S
Process Connection	Pg 13.5 thread
pH-Sensing Technology	Composite X-Chip
Reference System	Pre-pressurized liquid electrolyte, ceramic diaphragm, Argenthal™ cartridge with silver-ion trap
Redox Measurement	No
Lengths	120 mm, 225 mm
Shaft Material	Virgin Food Grade PEEK
O-Ring/Gasket Material	Fluoroprene XP – fluorinated elastomer
Solution Ground	Titanium Grade 2
Sterilizable/Autoclavable	No

ISM (digital)	Yes
EHEDG	Yes
3A	Yes
ATEX Compatibility	 II 1 / 2 G Ex ia HC T6/T5/T4/T3 Ga/Gb (SEV 14 ATEX 0168 X) II 1 / 2 G Ex ia HC T6/T5/T4/T3 Ga/Gb (IECEX SEV 14.0025X)
FM Compatibility	 Class I/II/III, Div1, Groups A, B, C, D, E, F, G T6; Tamb. = 0...+60 °C; (FM16US0034X, FM18CA0021X)

11 Drawings

12 Ordering Information

X-Chip pH Sensors for Food and Beverage Applications

Designation	Shaft Length (a) in mm	Order No.
InPro X1 HLS-N200-K120	120	30 389 700
InPro X1 HLS-N200-K225	225	30 389 701

13 Accessories

Designation	Order No.
Calibration on PC	
Software ISMCore Advanced	30 846 307
Cable AK9/1 m	59 902 167
Cable AK9/3 m	59 902 193
Cable AK9/5 m	59 902 213
Cable AK9/10 m	59 902 230
Cable AK9/20 m	52 300 204
Cable AK9/30 m	52 300 393
Cable AK9/50 m	52 300 394
Cable AK9/80 m	52 300 395
Cable DS AK9-RJ12 (connects ISM sensors to ISMCore)	52 300 383

Designation	Order No.
pH Buffers	
Technical Buffer pH 4.01, 6 × 250 ml	51 340 058
Technical Buffer pH 7.00, 6 × 250 ml	51 340 060
Technical Buffer pH 9.21, 6 × 250 ml	51 350 022
Cleaning/Storage Solutions	
Cleaning Solution Pepsin-HCl, 6 × 250ml	51 340 069
Electrolyte 9823 KCl 3M/L 6 × 250ml storage solution	51 340 050
Cleaning Solution Thiourea, 6 × 250ml	51 340 082
Verification Kit	
ISM Simulator pH Kit	52 300 410

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10/2023 © METTLER TOLEDO. All rights reserved.