Français

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Literature

Application brochure	51724646
Guideline to ION selective measurement	51300075

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
Subject to technical changes.
04/10 © Mettler-Toledo AG. Printed in Switzerland. 51 709 420

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English

General information

The DX223-Na⁺ sodium electrode is suitable for measurements of sodium activities and sodium concentrations in aqueous systems. Use of the DX223-Na⁺ is described in the application brochures.

Specifications

Type of ion selective membrane:
glass membrane

Resistance of the membrane:
<100 MΩ

Measurement range:
1x10⁻⁷ mol/L to 1 mol/L sodium

Measurement temperature range:
0–80 °C

Optimum pH range:
8...11

Interferences:

Silver ions (Ag⁺) must be absent. The sodium electrode is sensitive to the following ions. The ratio of this ion to sodium (X⁺/Na⁺) should thus not be larger than the value in brackets: H⁺ (<0.001), Li⁺ (<1), K⁺ (<5), NH₄⁺ (<50), Mg²⁺ (<2000).

Precision

In the absence of interfering substances, a precision of better than ±1 mV corresponding to ±4% of the measured Na⁺ concentration can be achieved.

Response times

When changing from lower to higher concentrations below 10⁻⁴ mol/L, responding takes approx. 30 seconds, above 10⁻⁴ mol/L less than 30 seconds; when changing from higher to lower concentrations several minutes.

Materials

Electrode body and membrane made of glass.

Recommended ISA solution

ISA (Ionic strength adjuster) solution for pre-conditioning both calibration and sample solutions. Prepare as follows: Weigh out 200 g NH₄Cl and dissolve in 500 mL deionized water, add 50 mL concentrated NH₃ solution and make up to 1 L. The ration of standard or sample to ISA should be 100:2.

Bridge electrolyte for reference electrode

0.1 mol/L NH₄Cl or as described in the application brochures. The ration of sample to ISA should be 100:2.

Maintenance and service

Deposits on or contamination of the surface of the membrane can be removed by rinsing the sodium electrode with deionized water or a lab cleaning agent. After every treatment, the sodium electrode must be rinsed with deionized water and then conditioned for several hours in saturated sodium solution.

Storage

The sodium electrode should be stored in 0.1 mol/L NaCl. It is an advantage to store the electrode in the protective tube supplied (protection against damage).

Reference electrode: Allow bridge electrolyte to run out. For storage, close filling port with the rubber cap. (The service life of the reference electrode is shortened if the bridge electrolyte chamber constantly remains filled above the internal ceramic frit.)

Conditioning

If required the sodium electrode might need to be conditioned in saturated NaCl solution for 48 hours.

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