**Gel electrodes** contain a high-viscous gelatine electrolyte in a sealed chamber. The electrolyte cannot be refilled; leak rates are however very low.

### 1 Pre-operation treatment

- 1.1 The internal buffer solution (3) must cover the whole inside of the glass membrane (1). Eliminate air bubbles by gently shaking in vertical direction.
- 1.2 Unscrew the glass vessel that is filled with the gel electrolyte and rinse the electrode with deion. water. The electrode is ready for use.

### 2 Adjustment and response time

The parameters zero-point (offset) and slope are measured by adjustment according to the RC1 or LabMax online help. A new electrode should have parameter values according to 2.1 - 2.2.

- 2.1 The zero point (pH value at  $E = 0$  mV) is approx. pH 7. The electrode potential, measured in a buffer solution of pH 7, must be within  $\pm 30$  mV.
- 2.2 The electrode slope (potential change per pH unit) at the pH range 4 - 9 must be within  $-55$  mV/pH and the theoretical value of  $-59.2$  mV/pH at  $25^\circ\text{C}$ .
- 2.3 Response time: After 30 s the indicated electrode potential must not change more than 2 mV within the following 30 s.
- 2.4 If testing a used electrode, the zero point should be within pH 6 and 8, the slope not less than  $-52$  mV/pH, and the potential change after 60 s should not exceed 3 mV/30 s.

### 3 Storage and lifetime

- 3.1 Store the electrode in the gel electrolyte, see 1.2 (always screw the cover onto the vessel). The frit (4) must always be covered; if necessary top up with a solution of 3 mol/L KCl.
- 3.2 The lifetime of the electrode is 6 months to 3 years, depending on its use and treatment.

### 4 Maintenance instructions

- 4.1 No reaction should take place between the sample solution and the reference electrolyte: Some cations ( $\text{Ag}^+$ ,  $\text{Hg}^{2+}$ ,  $\text{Au}^{3+}$ ,  $\text{Pb}^{2+}$ ) react with  $\text{Cl}^-$  of the reference electrolyte and cause a blockage of the frit. Solutions containing protein quickly block the frit, because protein precipitates on it.

Blocked frit by AgCl: Place the electrode overnight in concentrated ammonia, rinse and condition it for 1 hour in the gel electrolyte.

Blocked frit by protein: place the electrode at least for 1 hour in a pepsin solution (e.g. METTLER TOLEDO order No. 51 340 068).

Other contaminations: place the electrode for a few minutes in an ultrasonic bath with deion. water, ethanol or a mixed acid solution (e.g. aqua regia).

- 4.2 If the electrode does not meet the qualifications according to 2.4 anymore, try an etching of the glass membrane: Place the membrane 1 - 2 minutes in  $\text{NH}_4\text{HF}_2$  solution (e.g. METTLER TOLEDO order No. 51 340 073), rinse thoroughly with deion. water and place it for several hours in the gel electrolyte (see 1.2).
- 4.3 Just wipe the electrode after rinsing, never rub it dry.
- 4.4 Do not let the cable connector plug corrode.
- 4.5 If the glass membrane of the electrode is dried out, dissolve the KCl crystals in/on the frit with deion. water and reactivate the membrane by placing it for several hours in the gel electrolyte (see 1.2).

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## "pH-01" Electrode

Combined pH glass electrode with gel electrolyte

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