InPro 4010 pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	5
de	Betriebsanleitung	8
en	Instruction manual	11
es	Manual de instrucciones	14
fi	Käyttöohje	17
fr	Instructions d'utilisation	20
hu	Használati utasítás	23
it	Istruzioni per l'uso	26
ja	取扱説明書	29
ko	사용 매뉴얼	32
nl	Gebruikershandleiding	35
pl	Instrukcja obsługi	38
pt	Manual de instruções	41
ru	Инструкция по эксплуатации	44
sv	Bruksanvisning	47
th	คู่มือคำแนะนำ	50
zh	说明书	53

InPro 4010 pH elektroda

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Obsah

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Instruction Manual

Contents

1.	Introduction	12
2.	Safety instructions	12
3.	Product description.....	12
4.	Installation and commissioning	12
5.	Calibration and measurement.....	13
6.	Maintenance	13
7.	Removal of dirt and residue	13
8.	Environmental protection.....	13
9.	Guarantee.....	13

InPro is a trademark of the METTLER TOLEDO Group.

1. Introduction

METTLER TOLEDO pH sensors Type InPro™ 4010 are low-maintenance measuring electrodes with a polymer electrolyte and built-in temperature sensor. InPro 4010 sensors are designed for the simultaneous measurement of pH value and temperature in industrial wastewater.

Please read these instructions carefully before putting the sensor into operation, in order to ensure faultless function. InPro 4010 sensors are non-sterilizable/non-autoclavable and unsuitable for measurements at pH values <2, in media containing chlorine or in organic solutions.

2. Safety instructions

The sensor is intended only for the use described in §1 above. The sensor consists partly of glass. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. In addition, local and in-house-safety regulations must be observed.

We recommend that the electrodes only be put into operation in conjunction with original products/accessories from METTLER TOLEDO. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions.

3. Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro 4010/xxx/yyy	Type designation/ xxx = shaft length in mm/yyy = temperature sensor
combination pH	Type of electrode
pH 2 ...12	pH measurement range
0 ...60 °C (32...140 °F)	Temperature range for operation
Order No. 52 00X XXX	Order number

4. Installation and commissioning

1. On unpacking, and again prior to installation, check the sensor for mechanical damage of the glass membrane and shaft. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and rinse the sensor with deionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the glass membrane, since this can lead to electrostatic charging and sluggish response times.

3. Check the membrane chamber of the electrode for the presence of any air bubbles and remove same by lightly swinging the sensor in a vertical plane.
4. Please refer to the instruction manual of the relative housing in order to install the sensor.
5. Connect the sensor to the pH transmitter, using a cable with VP (Vario Pin) connector. Please refer to the wiring and terminal diagram accompanying the cable.

5. Calibration and measurement

InPro 4010 sensors require 2-point calibration. The sensor is then dipped in succession into two different buffer solutions with given pH reference values and the transmitter calibrated to these buffer values. Please refer to the instruction manual of the transmitter for details.

6. Maintenance

The sensor with polymer electrolyte cannot be refilled. The sensor must not be stored dry. In the event that the sensor dries out, it is necessary to recondition it before use by immersion in FRISCOLYT™ (Order No. 51 340 053) for at least 24 hours. InPro 4010 sensors should be stored with fitted watering cap containing FRISCOLYT. After each working cycle, the sensor tip and the apertures (diaphragms) must be thoroughly cleaned using deionized water.

7. Removal of dirt and residue

1. **Sluggish response time.** The glass membrane can be reactivated using a HF solution (Order No. 51 319 053). Only that section of the sensor tip below the two (diaphragms) apertures may be immersed in the HF solution, and this during a period of about 5 minutes. Immediately afterwards, rinse the sensor thoroughly and immerse in FRISCOLYT for at least 12 hours.
2. **Insufficient sensitivity (slope).** Check the apertures for contamination and rinse with warm water. Check the VP connector for traces of moisture. If necessary, clean well using deionized water or alcohol, and afterwards dry carefully.

8. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9. Guarantee

On manufacturing defects, 12 months after delivery.

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本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 04/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 000 522 D

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 5 2 0 0 0 5 2 2 0 *