METTLER TOLEDO

DM1147-SC

Platinum ring electrode with Sensor Chip

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**Plug & Play combined platinum ring
electrode**

for redox titrations without change of the pH
value.

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
1. Start-up maintenance

- 1.1 No start-up maintenance is necessary. The electrode is immediately ready for use.

2. Plug & Play performance

The Sensor Chip which stores the sensor parameter is located in the screw cap of the sensor.

The sensor is recognized after connection to the sensor input of the titrator as well as after each restart of the system by reading the sensor parameter on the chip. A popup appears on the titrator terminal when the sensor was successfully recognized. This popup allows you to directly enter the sensor setup. You find further information concerning the functionality of the Plug & Play sensor.

3. Testing

Testing a used electrode, a comparison between titration curves is suitable: perform a titration with sodium thiosulfate (0.1 mol/L) and iodine (0.1 mol/L) with the new electrode. The titration goes fast with a big steep potential change at the equivalence point. A used electrode performing the same titration should present both these parameters.

4. Storage and lifetime

- 4.1 Store the electrode in the electrolyte.
- 4.2 Never let the electrode get dry.
- 4.3 If testing a used electrode the requirements indicated under 3 should hold true.

5. Maintenance Instructions

- 5.1 The tolerable temperature range extends from 0 to 80 °C (32–176 °F).
- 5.2 In strong oxidizing solutions an oxygen adsorption takes place on the platinum surface of the electrode, in strong reducing solutions hydrogen adsorption. This can cause a change in the redox potential as well as an increased response time.
- 5.3 Do not let the cable connector plug corrode.

6. Cleaning procedures

- 6.1 Contamination of the platinum surface: depending on the contamination several cleaning procedures can be carried out.
 - a) Place the electrode for a few minutes in an ultrasonic bath with deion. water or ethanol.
 - b) Place the electrode for one minute in chromium sulfuric acid or aqua regia and rinse it thoroughly with deion. water.
 - c) Place the electrode for 30 minutes in a pH 4 buffer solution (Potassium hydrogen phthalate) saturated with quinhydrone to remove the absorbed oxygen or hydrogen.
 - d) Polish the platinum surface with Al_2O_3 powder.

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