## Instructions for welding in the threaded sockets

### 1. Safety instructions

Weld-in sockets may only be welded into place by personnel fully qualified for such work. It is the responsibility of the customer to ensure that in-house regulations or the currently valid standards/norms for welding and glueing work are fully observed.

### 2. Normal use

Weld-in sockets are intended exclusively for use in connection with specified METTLER TOLEDO housings/sensors. The function of the sockets is only ensured if the installed housing or 25 mm O<sub>2</sub> sensor is employed in combination with a captive ring nut.

### 3. Purpose

The inclined socket is supplied finish-machined with a precision-ground bore,  $\varnothing$  25-H7. It is suitable as a side entry port for reactors.

The straight socket also has a precision-ground bore,  $\varnothing$  25-H7, but it must be machined down to length, to suit the specific installation requirements.

(see illustration on the reverse)

Safety sockets are suitable for vertical arrangement through reactor headplates, or for horizontal and vertical installation into pipes.

Due to the material, the inner diameter of the threaded sockets tends to change slightly when the adaptor is welded into the counter-piece (reactor wall or top plate). To prevent the bore of the threaded socket from becoming too small or too large, we recommend making a small welding seam ( $\alpha = \text{max. } 3 \text{ mm}$ ) in one operation. You can obtain the best results with the TIG<sup>1)</sup> welding process.

For the polymers PE, PP and PVDF, the best results are to be obtained by using a hot air rapid-welding nozzle and the triangular wire welding process. For PVC sockets, Tangit adhesive is re-commended. The surfaces to be bonded, should be cleaned beforehand with Tangit pre-cleaner. After completing welding or adhesive bonding, the diameter of the bore needs to be checked using a gauge –H7<sup>2)</sup> and, if necessary, slightly reamed out with a finishing handreamer.

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11/2012 © Mettler-Toledo AG.