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### Uso conforme

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## Operating instructions

### Intended Use

METTLER-TOLEDO InLab® redox electrodes are intended for precise ORP measurement in the laboratory and in the field at atmospheric pressure.

### Technical Data

The exact specifications for the temperature range are printed on the electrode shaft and can be found at [www.mt.com/electrode-guide](http://www.mt.com/electrode-guide).

### General Safety Information

When working with chemicals, all relevant safety regulations provided by the manufacturer and the laboratory must be observed.

### Product-specific Safety Information

The electrode is only intended for the purpose specified above. The manufacturer does not assume liability for any damage that results from use outside of the intended purpose.

The electrode is made partly of glass; there is therefore a risk of injury in the event of breakage. Failure to observe the handling instructions below may result in distorted results and electrode damage. Do not use the electrode to stir the measurement medium or to touch the bottom or walls of the measuring vessel.

### Further Information

Further information is available at [www.mt.com/electrode-guide](http://www.mt.com/electrode-guide).

### Manufacturer

Mettler-Toledo GmbH  
Im Langacher 44, 8606 Greifensee, Switzerland  
Phone +41 22 567 53 22, Fax +41 22 567 53 23  
Internet [www.mt.com](http://www.mt.com)

### Commissioning

The electrode is supplied ready for use. When unpacking, please carry out a visual inspection to check that the glass is intact. Use a suitable cable to connect the electrode to the pH/mV meter (see Fig. 1). Now remove the watering cap (see Fig. 2) and rinse the electrode with distilled water. You must open SafeLock™ prior to use (see Fig. 3), check the fill level of the electrolyte, and refill if necessary. The electrode is then ready to take measurements. Platinum and gold electrodes can be checked with redox buffer solution (order no. 51350060).

Note: The mV signal depends on the reference system used. The Ag/AgCl reference system present in the electrode has an offset of 207 mV at 25 °C compared with a standard hydrogen electrode. This applies only on condition that the recommended Ag/AgCl reference electrolyte is used:  
3 mol/L KCl (order no. 51350072) /  
1 mol/L KNO<sub>3</sub> (order no. 51350078).

## Operating instructions

### Additional Information for Commissioning the InLab® Redox Flow

Turn the ground glass junction to loosen it (see Fig. 5). Allow a drop of electrolyte to flow out, then exert a small amount of pressure to reposition the ground glass junction. (Note: No electrolyte must leak out.) If the ground glass junction cannot be loosened, immerse the electrode in warm water between 60 and 80 °C for a few minutes.

### Refilling Electrolyte

Open SafeLock™ (see Fig. 3) and refill the electrode with the recommended electrolyte (see electrode shaft). The electrolyte level should be approx. 10 mm below the filling port.

### Storage

After use, rinse the electrode well with distilled water and close SafeLock™ (see Fig. 4). Replace the watering cap (see Fig. 2), filled with electrolyte or InLab storage solution (order no. 30111142). Store the half-cell dry. The electrode should be stored upright and at room temperature.

### Troubleshooting

Make sure that ...  
– all cables are connected correctly and the connections are clean.  
– the electrode is correctly filled with electrolyte.

### Disposal

In accordance with the requirements of the European Directive 2002/96/EC on waste electrical and electronic appliances (WEEE), this product must not be disposed of as household refuse. This also applies equally for countries outside of the EU in accordance with their applicable national regulations. Please dispose of this product in accordance with the local regulations in a separate container for electrical and electronic devices. If you have questions regarding this, please contact your responsible authority or the dealer from which you purchased the product. If you pass the product on to a third party (e.g., for private or commercial/industrial use), then this regulation must also be passed on accordingly. Thank you for contributing to the care of the environment.

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© Mettler-Toledo GmbH 12/2015