Content Determination of Sodium Starch Glycolate

Method for the acidimetric content determination of sodium in the sodium salt of starch glycolate in a non-aqueous medium


Sample	Sodium Starch Glycolate	Preparation and Procedures Dry and crush the resulting residue from the determination of sodium chloride in sodium starch glycolate. Weigh 0.49 to 0.51 g of the crushed residue, accurately to 0.1 mg, in a round bottomed flask. Add 80 mL acetic acid anhydrous to the flask as well as to an additional 3 empty flasks for the blank determinations, and reflux for 2 hours. Transfer the contents of the round bottomed flasks into a titration beaker, add 20 mL anhydrous acetic acid and titrate.
Substance	Sodium MG= 22.99 g/mol	
Chemicals	Acetic acid, anhydrous R	
Titrant	$c(\text{HClO}_4) = 0.1 \text{ mol/L}$	
Standard	THAM Tris-hydroxymethyl	
Instruments	DL70ES/DL77 and DL58/55+/53+	
Accessories	Titration beaker ME-101974 SP250 Peristaltic pump Temperature sensor DT120	Remarks As a result of the high expansion coefficient of the perchloric acid solution, its temperature is measured during standardisation (H16) and titration (T) and a correction factor (H5) is then calculated.
Indication	DG 113-SC	
Chemistry		
Calculation	$R(\%) = (Q[2] - x) \cdot H5 \cdot M/10 \cdot m$ Q[2] = Titrant consumption in mmol x = Mean value of the blank determination in mmol M = Molar mass Sodium m = Sample size in g H5 = Volume correction factor	
Waste disposal	Halogen-free organic solvents waste.	
Author	L. Pludra	

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