

Electroless Nickel Bath: Determination of Nickel

Method for determination of nickel content in electroless nickel bath.

Sample	Electroless nickel bath, 5 mL
Compound	Nickel, Ni M= 58.69 g/mol, z = 1
Chemicals	Deionized water, 50 mL Indicator : 0.2 g Murexide trituration with NaCl (1 : 500). Buffer pH10, 10 mL
Titrant	Ethylenediaminetetraacetic acid disodium, $C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$ $c(EDTA-Na_2) = 0.1 \text{ mol/L}$
Standard	Zinc Sulfate, $ZnSO_4$ $c(ZnSO_4) = 0.1 \text{ mol/L}$
Indication	DP5 Phototrode TM (555 nm) (Yellow to blue-violet)
Chemistry	$Ni^{2+} + \text{Murexide}^- \rightarrow \text{Ni-Murexide}^+$ $\text{Ni-Murexide}^+ + EDTA^{4-} \rightarrow \text{Ni-EDTA}^{2-} + \text{Murexide}^-$
Calculation	$R = Q \cdot C / m \cdot d, \text{ g/L}$ Q = Titrant consumption in mmol. C = M/z. M = Molar mass of sample in g. z = Equivalent no. sample, z=1 d = density of sample in g/mL. m = mass of sample in mL.
Waste disposal	Nickel solutions : If necessary, neutralize the solution before final disposal as special waste.
Author, Version	Ruby Das, MSG AnaChem, V2.0 Revised: C. De Caro, MSG AnaChem

Preparation and Procedures

CAUTION

- Use safety goggles, a lab coat and wear gloves. If possible, work in a fume hood.
- Ensure accurate cleaning of sensor is sufficient after each titration.

Sample Preparation:

- **Electroless nickel bath** : Pipette 20 mL Ginplate Ni 426-A and 20 mL Ginplate Ni 426-B in 200 mL volumetric flask and dilute it upto the mark with deionized water.

Sample titration:

- Add 0.25 g murexide trituration with NaCl (1:500) in the beaker placed on sample changer.
- Dispense 5 mL sample from dosing unit.
- Add 50 mL of deionized water from dosing unit.
- Add 10 mL of buffer pH10 from dosing unit.
- Titrate with 0.1 mol/L EDTA.
- After completion of each sample sensor, stirrer and titration tubes are rinsed by deionized water by means of membrane pump.
- Sensor is cleaned with deion. water in the conditioning beaker placed on sample changer after each sample.

Remarks

- Prior to use, adjust the output signal of the DP5 PhototrodeTM to approx. 1000 mV in deion. water before starting titration (100% transmission) by turning the small knob on the housing.
- Rinsing and conditioning of the Phototrode is crucial to achieve accurate and precise results.
- Avoid formation of bubbles during titration by low speed rate of stirrer, as they disturb photometric indication.
- This method allows a fully automated analysis procedure. This method can be easily modified for manual operation. Select "Manual stand" in the method function "Titration stand".
- Sample may be added manually using a pipette instead of using an additional dosing unit.

Literature :

- Ginplate NI 426, (<http://www.growel.com/tds/549.pdf>), a trademark of Grauer & Weil India Ltd, www.growel.com
- Mettler-Toledo Application M066 and M007

Instruments

- Titration Excellence T50/T70/T90
(Other Titrators: depending on instrument type, manual operation and method changes are necessary)
- XP205 Balance (MT-1106024)
- Rondo 20 with PowerShower™ (MT-51108003)

Accessories

- 3 x Additional dosing unit (MT-51109030)
- 1 x 20 mL DV1020 glass burette (MT-51107502)
- 2 x 10 mL DV1010 glass burette (MT-51107501)
- 1 x 5 mL DV 1005 glass burette (MT-51107500)
- 100 mL Propylene titration beakers (MT-00101974)

Results

All results

Method-ID	Nickel determination
Sample	5 mL (1/6)
R1 (Nickel)	5.44 g/L
Sample	5 mL (2/6)
R2 (Nickel)	5.44 g/L
Sample	5 mL (3/6)
R3 (Nickel)	5.45 g/L
Sample	5 mL (4/6)
R4 (Nickel)	5.43 g/L
Sample	5 mL (5/6)
R5 (Nickel)	5.43 g/L
Sample	5 mL (6/6)
R6 (Nickel)	5.44 g/L

Statistics

Method-ID	Nickel determination
R1	Nickel
Samples	6
Mean	5.44
s	0.01
srel	0.138%

Titration curve




Table of measured values


	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature °C
EQP1	0.000	NaN	757.2	NaN	NaN	0	25.0
	0.5715	0.5715	748.2	-9.0	NaN	5	25.0
	0.8570	0.2855	739.0	-9.2	NaN	13	25.0
	1.0000	0.1430	738.4	-0.6	NaN	16	25.0
	1.1000	0.1000	734.0	-4.4	NaN	20	25.0
	1.2000	0.1000	732.1	-1.9	-25.96	24	25.0
	1.3000	0.1000	729.6	-2.5	-24.64	28	25.0
	1.4000	0.1000	727.1	-2.5	-26.16	32	25.0
	1.5000	0.1000	723.3	-3.8	-24.72	37	25.0
	1.6000	0.1000	721.7	-1.6	-24.40	42	25.0
	1.7000	0.1000	720.6	-1.1	-23.28	44	25.0
	1.8000	0.1000	716.7	-3.9	-24.36	50	25.0
	1.9000	0.1000	713.5	-3.2	-28.28	55	25.0
	2.0000	0.1000	713.0	-0.5	-29.78	58	25.0
	2.1000	0.1000	708.5	-4.5	-30.43	62	25.0
	-----	-----	-----	-----	-----	-----	-----
	4.3440	0.0500	518.0	-17.4	-358.17	184	25.0
	4.3940	0.0500	499.4	-18.6	-421.44	190	25.0
	4.4440	0.0500	476.6	-22.8	-538.07	196	25.0
	4.4940	0.0500	446.3	-30.3	-658.13	202	25.0
	4.540268	NaN	420.7	NaN	-659.56	NaN	NaN
	4.5440	0.0500	418.6	-27.7	-657.52	210	25.0
	4.5940	0.0500	379.1	-39.5	-544.67	218	25.0
	4.6440	0.0500	339.4	-39.7	-393.50	230	25.0
	4.6940	0.0500	322.0	-17.4	-269.40	239	25.0
	4.7600	0.6600	320.1	-1.9	-148.20	243	25.0
	4.8600	0.1000	319.3	-0.8	NaN	245	25.0
	4.9600	0.1000	319.6	0.3	NaN	248	25.0
	5.0600	0.1000	320.4	0.8	NaN	250	25.0
	5.1600	0.1000	321.2	0.8	NaN	253	25.0
	5.2600	0.1000	321.9	0.7	NaN	255	25.0

Comments

- Titer determination of 0.1 mol/L EDTA-Na₂ is done as per the Mettler –Toledo method application M007 and mean value found is 0.99163.
- The mean value of the titer is automatically stored as part of the setup by the function TITER.
- The buffer pH 10 is prepared by dissolving 64 g NH₄Cl in 600mL 25% Ammonia solution and diluting it upto the mark with deionized water in a 1L volumetric flask..
- The shape of the titration curve is somewhat affected by the concentration of the indicator. The results, however do not differ significantly (tested range: 25-500 mg of Murexide trituration with NaCl(1 :500).
- Add the indicator before starting analysis. Keep sample free of air bubbles during titration. Air bubbles and undissolved impurities affect the photometric indication.
- Due to the steep signal change, an EQP titration with fixed increments is used. The low threshold value allows for different amounts of indicator.

Principle :

- Nickel ions forms yellow complex with murexide in alkaline solution:



- By adding EDTA , Ni forms a more stable complex with EDTA:



At the equivalence point, all Ni ions have been complexed by EDTA and murexide is free in the alkaline solution. There is a colour change from yellow to blue-violet.

Method

001 Title

Type General titration
Compatible with T50 / T70 / T90
ID Nickeldetermination
Title Ni determination
. . .

002 Sample

Number of IDs 1
ID 1 Nickel
Entry type Fixed volume
Volume 5.0 mL
Density 1.03 g/mL
Correction factor 1.0
Temperature 25.0°C
Entry Arbitrary

003 Titration stand (Rondo/TowerA)

Type Rondo/TowerA
Titration stand Rondo60/1A
Lid handling No

004 Dispense (normal) [1]

Titrant NI SAMPLE
Concentration 1
Volume 5.0 mL
Dosing rate 60.0 mL/min
Condition No

005 Dispense (normal) [2]

Titrant Water..
Concentration 1
Volume 50.0 mL
Dosing rate 60.0 mL/min
Condition No

006 Dispense (normal) [3]

Titrant BUFFER 10PH
Concentration 1
Volume 10.0 mL
Dosing rate 60.0 mL/min
Condition No

007 Stir

Speed 10%
Duration 60 s
Condition No

008 Titration (EQP) [1]

Titrant

Titrant EDTA(0.1M)
Concentration 0.1 mol/L

Sensor

Type Phototrode
Sensor DP5
Unit mV

Temperature acquisition

Temperature measurement No

Stir

Speed 10%

Predispense

Mode Volume
Volume 1.0
Wait time 0 s

Control

Control User
Titrant addition Dynamic
dE (set value) 10 mV
dV (min) 0.05 mL
dV (max) 0.1 mL
Mode Equilibrium controlled
dE 1.0 mV
dt 2 s
t (min) 2 s
t (max) 12 s

Evaluation and recognition

Procedure Standard
Threshold 200.0 mV/mL
Tendency None
Ranges 0
Add. EQP criteria No

Termination

At Vmax 10.0 mL
At potential No
At slope No
After number of recognized EQPs Yes

Number of EQPs 1
Combined termination criteria No

Accompanying stating

Accompanying stating No

Condition

Condition No

009 Calculation R1

Result Ni content
Result unit g/L
Formula $R1 = Q \cdot C / m \cdot d$
Constant $C = M / z$
M $M[\text{Nickel}]$
z $z[\text{Nickel}]$
Decimal places 2
Result limits No
Record statistics Yes
Extra statistical func. No
Send to buffer No
Condition No

010 Rinse

Auxillary reagent WATER
Rinse cycles 1
Vol.per cycle 10 mL
Position Current position
Drain No
Condition No

011 Condition

Type Fix
Interval 1
Position Conditioning beaker
Time 20 s
Speed 10 %
Condition No

012 End of sample