

## Sodium Determination with Ion Selective Electrode

Method for the sodium determination by direct measurement with a sodium ion selective electrode.

Sample	50 mL Sample solution with 0.1 to 0.001 mol/L sodium	<b>Preparation and Procedures</b>  <b>Milk:</b> Dilute sample 1:10. To 50 mL of this add 1 mL ISA solution and measure. <b>Apple juice:</b> To 50 mL undiluted sample add 1 mL ISA sol. and 1.5 mL NH <sub>4</sub> OH 5% . The pH of the sample is adjusted to 7.5 with excess NH <sub>4</sub> OH (For correct measurement with Na-ISE, the pH must be over 7). <b>Wine, grape juice:</b> To 50 mL undiluted sample add 1 mL ISA sol. and 1.5 mL NH <sub>4</sub> OH 5% . The pH of the sample is adjusted to 8 with excess NH <sub>4</sub> OH (the pH must be over 7). <b>Mineral water:</b> To 50 mL sample add 1 mL ISA solution. <b>Drinking water:</b> To 50 mL sample add 1 mL ISA solution and measure. Due to the low sodium content, the standards 10-3 and 10-4 mol/L were used for calibration for this sample. <b>Bouillon:</b> Place a bouillon cube in 300 mL deion.water and warm till it has dissolved. Let cool and use deion.water to adjust volume to 500 mL. Dilute this solution 1:20 resp. 1:50 with deion.water. To 50 mL add 1 mL ISA sol. and measure. <b>Herb patè:</b> Mix entire contents of can (125g) with 700 mL deion. water for about 5 minutes. Decant fat. Fill to 1000 mL with deion. water. To 50 mL of this add 1 mL ISA solution and measure. <b>Infusion solution:</b> Dilute sample 1:100 and to 50 mL of this add 1 mL ISA solution and measure.																												
Substance	Sodium M = 22.99																													
Chemicals	1 mL ISA solution (4 m NH <sub>4</sub> Cl and 1 m NH <sub>4</sub> OH)																													
Titrant																														
Standard																														
Instruments	DL67, DL70ES, DL77 Printer (HP Desk Jet 500)																													
Accessories	Titration beaker ME101974 DT120 (Temp. sensor Pt100)																													
Indication	DX223 sodium ion selective electrode DX200 reference electrode																													
Chemistry		<b>Remarks</b>  <b>Reproducibility and Recovery</b> Aqueous Na solutions made from NaCl were used (titrimetric standard). The Na-ISE was recalibrated daily with standards 0.1 and 0.001 g Na <sup>+</sup> /L. It was stored in a concentrated sodium solution. To 50 mL sample add 1 mL ISA sol. then measure. The solutions were produced daily. <table><tr><th>Conc.</th><th>Recovery</th><th colspan="2">RSD of sev. series(n=6)</th></tr><tr><th>mol/L</th><th>%</th><th>Method A</th><th>Method B</th></tr><tr><td>0.1</td><td>98-102%</td><td>0.2-0.5%</td><td>0.2-0.6%</td></tr><tr><td>0.01</td><td>95-100%</td><td>0.2-0.6%</td><td>0.3-0.7%</td></tr><tr><td>0.001</td><td>97-101%</td><td>0.2-0.8%</td><td>0.2-0.9%</td></tr><tr><td>0.0005</td><td>101-106%</td><td>0.2-0.4%</td><td>0.2-0.5%</td></tr><tr><td>0.0001</td><td>131-129%</td><td>0.2-0.9%</td><td>0.2-0.9%</td></tr></table> <b>Result</b> The differences between A and B are not decisive. Below 0.001 mol/L longer stir times (15-20 min) were necessary to get a good reproducibility. The concentration 0.001 mg/L leaves the linear range (see cal. curve). Thus, the useful concentration range was limited from 0.1 mol/L to 0.001 mol/L. <b>General remarks:</b> All measurements were performed using the simple direct measurement (method B). All measurements should be performed using the same stirring conditions. This means the same speed, stirrer type, distance from stirrer to electrode etc.	Conc.	Recovery	RSD of sev. series(n=6)		mol/L	%	Method A	Method B	0.1	98-102%	0.2-0.5%	0.2-0.6%	0.01	95-100%	0.2-0.6%	0.3-0.7%	0.001	97-101%	0.2-0.8%	0.2-0.9%	0.0005	101-106%	0.2-0.4%	0.2-0.5%	0.0001	131-129%	0.2-0.9%	0.2-0.9%
Conc.	Recovery		RSD of sev. series(n=6)																											
mol/L	%		Method A	Method B																										
0.1	98-102%		0.2-0.5%	0.2-0.6%																										
0.01	95-100%	0.2-0.6%	0.3-0.7%																											
0.001	97-101%	0.2-0.8%	0.2-0.9%																											
0.0005	101-106%	0.2-0.4%	0.2-0.5%																											
0.0001	131-129%	0.2-0.9%	0.2-0.9%																											
Calculation																														
Waste disposal																														
Author	A. Aichert/Anke Stock																													

## Results

### Results

Sample	Dilution	Direct measurement with sodium electrode			Content (Producer spec.)
		n	Mittelwert	RSD	
Milk, full cream, pasteurized	1:10	6	400.8 mg/L	0.20%	430-490 mg/L*
Milk, full cream, raw	1:10	6	433.8 mg/L	0.20%	430-490 mg/L*
Apple Juice	none	4	8.2 mg/L	0.64%	[30 mg/L*
Grape Juice, red	none	5	23.8 mg/L	0.66%	[30 mg/L*
Wine, red Italian	none	5	12.2 mg/L	0.63%	10-80 mg/L*
Wine, White, Swiss	none	5	8.4 mg/L	2.0%	10-80 mg/L*
Mineral Water	none	6	8.17 mg/L	0.22%	8.2 mg/L
Drinking Water, Schwerzenbach	none	5	3.81 mg/L	0.35%	--
Bouillon sodium reduced	1:20	5	468.5 mg/L	0.35%	[1.2 g/L
Bouillon normal	1:50	5	4110 mg/L	0.41%	-4 g/L
Herb patè sodium reduced	1:8	6	175.9 mg/kg	0.52%	[0.4 g/kg
Infusion Solution	1:100	5	142.5 mmol/L	0.72%	141.5 mmol/L

\* Obtained from: Schweizerisches Lebensmittelbuch

**Method A** stirs for 5 minutes then takes 6 measurements with an interval of at least 10 seconds. The sodium concentration is calculated from the mean of these 6 measured values.

**Method B** measures once after a stir time of 5 minutes. The sodium concentration is calculated from this. If the sample is diluted for the measurement, the dilution factor can be entered as a correction factor (f) for each sample. This will be used in the calculation.

### Table of measured values

### Titration curve

## Method

Method A  
(Direct measurement with mean value)

Title  
Method ID . . . . . Na0A  
Title . . . . . Na+ -content  
Date/time . . . . . 15-Mar-1994 19:48

Sample  
Number samples . . . . . 6  
Titration stand . . . . . Stand 1  
Entry type . . . . . Fixed volume U  
Volume [mL] . . . . . 50.0  
ID1 . . . . .  
Molar mass M . . . . . 0.0  
Equivalent number z . . . . . 1  
Temperature sensor . . . . . TEMP A

Stir  
Speed [%] . . . . . 50  
Time [s] . . . . . 600

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.2  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 10.0  
t (max) [s] . . . . . 90.0

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.2  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 10.0  
t (max) [s] . . . . . 90.0

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.2  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 10.0  
t (max) [s] . . . . . 90.0

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.2  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 10.0  
t (max) [s] . . . . . 90.0

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.2  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 10.0  
t (max) [s] . . . . . 90.0

Calculation  
Result name . . . . .  
Formula . . . . .  $R1=E[1]$   
Constant . . . . .  
Result unit . . . . . p(Na+)  
Decimal places . . . . . 4

Calculation  
Result name . . . . . Na+ -single  
Formula . . . . .  $R2=pw(-E[1])*1000$   
Constant . . . . .  
Result unit . . . . . mmol/L  
Decimal places . . . . . 4

Calculation  
Result name . . . . .  
Formula . . . . .  $R3=(C3+E[5]+E[6])/6$   
Constant . . . . .  $C3=E[1]+E[2]+E[3]+E[4]$   
Result unit . . . . .  
Decimal places . . . . . 3

Calculation  
Result name . . . . . Na+ -x of 6  
Formula . . . . .  $R4=pw(-R3)*1000$   
Constant . . . . .  
Result unit . . . . . mmol/L  
Decimal places . . . . . 4

Statistics  
Ri (i=index) . . . . . R2  
Standard deviation % . . . . . Yes  
Rel. standard deviation srel . . . . . Yes  
Outlier test . . . . . Yes

Statistics  
Ri (i=index) . . . . . R4  
Standard deviation % . . . . . Yes  
Rel. standard deviation srel . . . . . Yes  
Outlier test . . . . . Yes

Record  
Output unit . . . . . Printer  
All results . . . . . Yes

Method B  
(Simple direct measurement)

Title  
Method ID . . . . . Na0B  
Title . . . . . Na+ -content  
Date/time . . . . . 22-Mar-1994 13:07

Sample  
Number samples . . . . . 6  
Titration stand . . . . . ST20  
Entry type . . . . . Fixed volume U  
Volume [mL] . . . . . 50.0  
ID1 . . . . .  
Molar mass M . . . . . 22.99  
Equivalent number z . . . . . 1  
Temperature sensor . . . . . TEMP A

Stir  
Speed [%] . . . . . 50  
Time [s] . . . . . 300

Measure  
Sensor . . . . . Na+ -sensor  
Unit of meas. . . . . As installed  
dE [mV] . . . . . 0.5  
dt [s] . . . . . 5.0  
t(min) mode . . . . . Fix  
t (min) [s] . . . . . 60.0  
t (max) [s] . . . . . 120.0

Calculation  
Result name . . . . .  
Formula . . . . .  $R=E$   
Constant . . . . .  
Result unit . . . . . p(Na+)  
Decimal places . . . . . 3

Calculation  
Result name . . . . . Na+ -content  
Formula . . . . .  $R2=pw(-E)*f*1000$   
Constant . . . . .  
Result unit . . . . . mg/L  
Decimal places . . . . . 2

Statistics  
Ri (i=index) . . . . . R2  
Standard deviation % . . . . . Yes  
Rel. standard deviation srel . . . . . Yes  
Outlier test . . . . . Yes

Statistics  
Ri (i=index) . . . . . R3  
Standard deviation % . . . . . Yes  
Rel. standard deviation srel . . . . . Yes  
Outlier test . . . . . Yes

Record  
Output unit . . . . . Printer  
All results . . . . . Yes