Content Determination of Magnesium in Magnesium Stearate

Method for the complexometric determination of magnesium in Magnesium Stearate

Sample	Magnesium stearate, $C_{36}H_{70}O_4Mg$	Preparation and Procedures Weigh a sample of 0.74 to 0.76 g magnesium stearate, accurately to 0.1 mg, in a 100 mL titration beaker and add 50 mL of a mixture of equal volumes 1-butanol and anhydrous ethanol. Pipette into the sample 5 mL ammonium solution 26%, 3 mL ammonium chloride buffer solution pH10 and 30 mL 0.1M Na_2 -EDTA. Place the titration beaker in an ultrasonic bath at 45 - 50 °C, and leave for approximately five minutes until the sample is completely dissolved. Prepare three blank determination samples in the same way.
Substance	Magnesium $M = 24.31 \text{ g/mol}$	
Chemicals	1-Butanol R Ammonium chloride buffer solution pH 10 R 0.1M Na_2 -EDTA	
Titrant	$c(ZnSO_4) = 0.1 \text{ mol/L}$	
Standard	EDTA-Dihydrate (M539)	
Instruments	DL77/DL70ES and DL58/55+/53+	
Accessories	Titration beaker ME-101974, SP250 Peristaltic pump, Rondo60 Sample changer	
Indication	DM 140-SC (amalgamated)	Remarks The previously determined loss on drying content is saved as the auxiliary value H20.
Chemistry		
Calculation	$R(\%) = (x-Q) \cdot M \cdot F/10 \cdot m$ Q = Titrant consumption in mmol x = Mean value of the blank determination in mmol F = $100/(100-H20)$ Formula for the loss on drying of Magnesium Stearate M = Molar mass Magnesium m = Sample size in g	
Waste disposal	Halogen-free organic solvents waste, zinc containing.	
Author	L. Pludra	

Results

METTLER DL70ES Titrator V3.1 Dr. Kade Pharm. Fabrik GmbH
Onkel Tuka Quality Control Laboratory


Method: G061 Cont. Magnesium stearate Date: 29-10-1999 10:02
User: Kuntze/Karasch Measured: 29-10-1999 13:19

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