

Water Determination in Ketones

The water content in acetone (ketone) is determined by volumetric Karl Fischer titration.

Sample	Acetone 1.5-2 g
Compound	Water, H ₂ O M = 18.01 g/mol
Chemicals	60 mL solvent HYDRANAL [®] KetoSolver or apura [®] CombiSolvent 5 Keto
Titrant	HYDRANAL [®] Composite 5 K (5 mg H ₂ O/mL) or apura [®] CombiTitant 5 Keto
Standard	Water Standard 10.0 mg/g (HYDRANAL [®] Water Standard 10.0)
Indication	DM143-SC Electrode
Chemistry	$\text{CH}_3\text{OH} + \text{SO}_2 + 3 \text{RN} + \text{I}_2 + \text{H}_2\text{O} \rightarrow (\text{RNH}) \cdot \text{SO}_4\text{CH}_3 + 2 (\text{RNH})\text{I}$
Calculation	<p>Water content in ppm:</p> $\frac{(\text{VEQ} \cdot \text{CONC} - (\text{TIME} \cdot \text{DRIFT} / 1000)) \cdot C}{m}$ <p>VEQ: Titrant consumption to EP TIME: Total time since sample request C: 1000 (constant for ppm calc.)</p>
Waste disposal	Organic solvents
Author, Version	Market Support Group Anachem Maria-José Schmid, August 2009

Preparation and Procedures

- Drawn out approximately 2-3 mL of sample into a 10 mL syringe with needle to rinse it before sample analysis.
- Discard the sample into disposal container and fill the syringe again with the sample.
- The sample is then injected into the KF titration vessel in aliquots of approx. 2-3 mL.
- The weight is determined by back-weighing technique.
- If the amount of water is very small, the titration can also be started with the parameter "Cautious" (see Control Parameters).

Remarks

- First, the concentration determination is performed by pressing "Start conc." (Water Standard 10.0 mg/g, approx. 0.8 g).
- The concentration determination is calculated according to the formula:

$$\text{CONC} = \text{CONT} \cdot m / (\text{VEQ} - (\text{DRIFT} / \text{CONC}(\text{alt})) \cdot \text{TIME} / 1000)$$

CONT = Conc. of liquid water standard (mg/g)
CONC(alt) = CONC(alt) refers to the Setup value current at the time of calculation.
- Subsequently, the sample analysis is started by pressing "Start sample" on the display.

Instruments

- KF Compact Line Volumeters V20/V30 (V2.0)
- Titration Excellence T70/T90 (V3.1.1)
- XP205 Balance

Accessories

- LabX Titration Software
- 10 mL syringe with needle

Results**HYDRANAL® Composite 5 K / KetoSolver**


No.	Comment / ID	Start time	Sample size and results			
1/6	--	20.08.2009 09:53:25	R1 =	0.209	%	Content
2/6	--	20.08.2009 11:06:25	R1 =	0.216	%	Content
3/6	--	20.08.2009 11:48:21	R1 =	0.218	%	Content
4/6	--	20.08.2009 12:02:55	R1 =	0.205	%	Content
5/6	--	20.08.2009 12:09:48	R1 =	0.209	%	Content
6/6	--	20.08.2009 12:19:20	R1 =	0.214	%	Content

Statistics

Rx	Name	n	Mean value	Unit	s	srel[%]
R1	Content	6	0.212	%	0.005	2.340

apura® CombiTitrant 5 Keto/CombiSolvent 5 Keto**Statistics**

Rx	Name	n	Mean value	Unit	s	srel[%]
R1	Content	6	0.209	%	0.005	2.182

Titration curve

Sample 1/6 20.08.2009 09:53:25


Table of measured values

Volume mL	Increment mL	H ₂ O mg	Online drift µg/min	Signal mV	Change mV	Time s
0.00000	NaN	0.0000	0.0	671.1	NaN	0
0.00000	0.00000	0.0000	0.0	671.1	0.0	0
0.00000	0.00000	0.0000	0.0	671.7	0.6	1
0.00100	0.00100	0.0049	153.2	672.7	1.0	2
0.00550	0.00450	0.0267	552.6	672.2	-0.5	3
0.01400	0.00850	0.0680	1045.9	659.5	-12.7	4
0.02675	0.01275	0.1299	1590.5	617.2	-42.3	5
0.02700	0.00025	0.1311	1333.2	612.0	-5.2	6
0.02750	0.00050	0.1335	1161.1	630.3	18.3	7
0.03075	0.00325	0.1493	1134.0	634.6	4.3	8
0.03800	0.00725	0.1845	1243.9	621.2	-13.4	9
0.04950	0.01150	0.2404	1456.7	579.7	-41.5	10
...
0.97900	0.00025	4.7538	50.5	147.5	-0.4	232
0.97900	0.00000	4.7538	43.3	145.6	-1.9	234
0.97900	0.00000	4.7538	28.8	149.2	3.6	234
0.97900	0.00000	4.7538	21.6	147.0	-2.2	236
0.97925	0.00025	4.7550	28.8	147.5	0.5	236
0.97950	0.00025	4.7562	36.1	146.9	-0.6	238
0.97950	0.00000	4.7562	36.1	147.8	0.9	238
0.97975	0.00025	4.7575	36.1	146.7	-1.1	240
0.97975	0.00000	4.7575	36.1	147.8	1.1	241
0.97975	0.00000	4.7575	28.8	147.9	0.1	242
0.98000	0.00025	4.7587	28.8	147.7	-0.2	243
0.98050	0.00050	4.7611	43.3	146.1	-1.6	244
0.98050	0.00000	4.7611	43.3	149.1	3.0	245
0.98050	0.00000	4.7611	43.3	146.3	-2.8	246
0.98075	0.00025	4.7623	43.3	149.4	3.1	246

Sample 1/6

Comments

- Ketones such as acetone react with methanol forming ketals. In this reaction water is released:



Therefore, methanol-free solvent and titrant are used to avoid this side-reaction. For this, special reagents (the so-called K-reagents) are commercially available.

- However, side-reactions can not be completely suppressed. The methanol-free reagents also react slowly with acetone. This leads to an increased drift value at the end of each sample titration. Thus, it is necessary to wait for a constant drift value (e.g. approx. 5 minutes) after each sample determination for the side reaction to subside.
- The drift value is also increasing with increasing number of titrated samples. Thus, the solvent must be replaced after 2-3 samples, depending on the sample size. This can be done by defining in the setup the solvent exchange after 2-3 titrated samples (Setup > Global > Solvent control > Monitoring no. of samples).
- To achieve precise results it is necessary to weigh the sample accurately, and in particular, to first clean the syringe with a few mL of sample which is then discarded.
- To ensure a more efficient cleaning, the syringe is gently shaken to allow for absorption of the moisture on the inner wall of the syringe. Subsequently, the syringe is completely filled.

Method (V2.0)

001 Title		
Type	Karl Fischer titration	
Compatible with	T70/T90/V20/V30	
ID	M303	
Title	Acetone Dry	
Author	METTLER TOLEDO	
Date/Time	18.08.2009 11:49:10	
Modified at	18.08.2009 11:49:10	
Modified by	METTLER TOLEDO	
Protect	No	
SOP	None	
002 Sample		
Number of IDs	1	
ID 1	--	
Entry type	Weight	
Lower limit	3.0 g	
Upper limit	8.0 g	
Density	1.0 g/mL	
Correction factor	1.0	
Weight per piece	1.0 g	
Temperature	25.0°C	
Autostart	No	
Entry	After addition	
Concentration		
Titrant	KF 1-comp 5K	
Nominal concentration	5 mg/mL	
Standard	Water-Standard 10.0	
Entry type	Weight	
Lower limit	0.6 g	
Upper limit	1.8 g	
Temperature	25.0°C	
Mix time	10 s	
Autostart	Yes	
Entry	After addition	
Concentration limits	Yes	
Lower limit	4.5 mg/mL	
Upper limit	5.8 mg/mL	
003 Titration stand (KF stand)		
Type	KF stand	
Titration stand	KF stand	
Source for drift	Determination	
Max. start drift	25 µg/min	
004 Mix time		
Duration	5 s	
005 Titration (KF Vol) [1]		
Titrant		
Titrant	KF 1-comp 5K	
Nominal concentration	5 mg/mL	
Reagent type	1-comp	
Sensor		
Type	Polarized	
Sensor	DM143-SC	
Unit	mV	
Indication	Voltametric	
Ipol	10.0 µA	
Stir		
Speed	35 %	
Predispense		
Mode	None	
Wait time	0 s	
Control		
End point	150 mV	
Control band	300.0 mV	
Dosing rate(max)	8 mL/min	
Dosing rate(min)	4 µL/min	
Start	Normal	
Termination		
Type	Delay time	
Delay time	10	
At VMax	10.0 mL	
Min. time	0 s	
Max. time	infinity	
006 Calculation R1		
Result tape	Predefined	
Result	Content	
Result unit	%	
Formula	$R1 = (VEQ * CONC - TIME * DRIFT / 1000) * C / m$	
Constant C=	0.1	
Decimal places	3	
Result limits	No	
Record statistics	Yes	
Extra statistical functions	No	
007 Record		
Summary	Per sample	
Results	No	
Raw results	No	
Table of measured values	No	
Sample data	No	
Resource data	No	
E-V curve	No	
E-t curve	No	
V-t curve	No	
H ₂ O-t	No	
Drift-t	No	
H ₂ O-t & Drift-t curves	No	
V-t & drift-t	No	
Method	No	
Series data	No	
008 End of sample		
Open series	Yes	