

Hydroxyl Number of Raw Polyols according to ASTM E1899-02

Polyols are basic raw materials used to produce different products like e.g. polyurethanes. To assure product constancy and quality of the final product, the hydroxyl number of the raw polyol has to be known. It is determined according to ASTM E1899-02.

Sample	Raw polyol, 0.5-1 g Pure octanol, C ₈ H ₁₆ OH, M = 130.23 g/mol
Compound	Hydroxyl end group, -OH
Chemicals	Tetrahydrofuran, THF p-Toluenesulfonyl isocyanate, TSI Water Acetonitrile, CH ₃ CN, 2-propanol
Titrant	Tetrabutylammonium hydroxide (TBAH) in 90% 2-propanol C(TBAH) = 0.1 mol/L
Standard	Potassium hydrogenphthalate, KHP n-octanol
Indication	DGi116-Solvent Electrolyte: Lithium nitrate saturated in acetonitrile
Chemistry	$\text{ROH} + \text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-NCO} \rightarrow \text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-NH-COOR}$ $\text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-NCO} + \text{H}_2\text{O} \rightarrow \text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-NH}_2 + \text{CO}_2$ $\text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-NH-COOR} + \text{OH}^- \rightarrow \text{CH}_3\text{-C}_6\text{H}_4\text{-SO}_2\text{-N}^-\text{-COOR} + \text{H}_2\text{O}$
Calculation	Hydroxyl number (mg KOH/g): $R1 = Q2 \cdot C / m$; C=56.1 Hydroxyl equivalent weight (g/equiv.) $R2 = (56.1 \cdot 1000) / R1$ Pure subst.: Molecular weight (g/mol) $R3 = (56.1 \cdot 1000) \cdot C / R1$ C = 1 (nr. OH-group/molecule, e.g. 1 for n-octanol)
Waste disposal	Organic solvent waste
Author, Version	Chris Hynes / Ross C. Koile MT-NA, June 2009

Preparation and Procedures

- All reagents must be of high purity (e.g. HPLC grade) since the amount of water contained interferes with the analysis.
- Isocyanate solution (TSI):
Add 25 g of p-toluenesulfonyl isocyanate to 450 mL of anhydrous acetonitrile.
- Titrant TBAH:
100 mL 1M TBAH in methanol are added in a 1 L volumetric flask with 500 mL 2-propanol. Mix and fill up to the mark with 2-propanol.
- Special beaker 1 on Rondo Sample Changer is filled with acetonitrile to eliminate water traces.
- Method is best if run with two peristaltic pumps (THF, water) and a dosing unit (TSI) to automate the sample preparation (Step 1):
 - 10 mL THF are added with the 1st pump
 - 10 mL TSI solution are dispensed (dosing unit)
 - stir the solution
 - 2 mL water are added by the 2nd pump
 - Stir the sample solution.
- Subsequently, 40 mL THF are added and the titration is started (Step 2).

Remarks

- A blank value is dependent on the amount of alcohol in the THF. In fact, alcohol affects the chemical reaction. If a high quality, anhydrous THF supply is used, there should rarely show a detectable blank.
- The following points have to be performed to achieve a relative standard deviation below 1%:
 1. The titration tip needs to be cut above the sample line of step 1 (30 mL) but below it after final addition of THF (Step 2).
 2. All other sample tips should be cut high enough to never be below the sample line.
- Sample size is critical to this test method:
 - Use 0.5 g for 1000 g/equiv
 - Use 1 g for 2000 g/equiv

Great care should be taken not to get any sample on the sides of the beaker. This will not allow for reaction of polyol with isocyanate resulting in low hydroxyl numbers.

Instruments	<ul style="list-style-type: none"> - Titration Excellence T70 / T90 - 2 SP250 peristaltic pumps - Rondo 20 Sample Changer
Accessories	<ul style="list-style-type: none"> - 1 x Dosing Unit - 2 DV1020 glass burettes - Balance XS204 - LabX titration software

Results

OH-Number of pure n-octanol

Results

No.	Comment / ID	Start time	Sample size and results			
1/3	OCTANOL	2/25/2009 1:11:40 PM	0.1310	g		
			R1 = 431.633	mg KOH/g	Hydroxyl number	
			R2 = 129.971	g/equiv	Equivalent weight	
			R3 = 259.943	g/mol	Molecular weight	
2/3	OCTANOL	2/25/2009 1:40:04 PM	0.1321	g		
			R1 = 429.911	mg KOH/g	Hydroxyl number	
			R2 = 130.492	g/equiv	Equivalent weight	
			R3 = 260.984	g/mol	Molecular weight	
3/3	OCTANOL	2/25/2009 2:08:06 PM	0.1310	g		
			R1 = 426.361	mg KOH/g	Hydroxyl number	
			R2 = 131.579	g/equiv	Equivalent weight	
			R3 = 263.157	g/mol	Molecular weight	

Statistics

Rx	Name	n	Mean value	unit	s	srel[%]
R1	Hydroxyl number	3	429.302	mg KOH/g	2.688	0.626
R2	Equivalent weight	3	130.681	g/equiv	0.820	0.628
R3	Molecular weight	3	261.361	g/mol	1.640	0.627

Additional results

Statistics

Rx	Name	n	Mean value	unit	s	srel[%]
R1	Hydroxyl number	3	437.653	mg KOH/g	2.969	0.678
R2	Equivalent weight	3	128.188	g/equiv	0.873	0.681
R3	Molecular weight	3	256.376	g/mol	1.746	0.681


Remarks:

- The titration of OH end groups of polyols can be tested using pure n-octanol, $C_8H_{16}OH$.
- Note that n-Octanol has only **one** OH group at the end of its aliphatic chain. For this reason, the equivalent weight corresponds to the molecular weight, i.e. $R2 = R3 = 130.681$ g/mol.
- However, method m453A has been developed to be performed with any polyol, and therefore both end groups have to be taken into account into calculation R3.
- For this reason, calculation R3 is exactly two times R2 in the case of octanol. Obviously, calculation R3 is not valid for n-octanol and does not have to be considered.

Titration curve

E - V curve EQP titration [1]

Sample 1/3



Sample 1/3, series 2/25/2009, 1:11:38 PM

Table of measured values

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature °C
EQP1	0.000	NaN	217.8	NaN	NaN	0	25.0
	0.040	0.040	186.4	-31.4	NaN	6	25.0
	0.080	0.040	166.5	-19.9	NaN	12	25.0
	0.120	0.040	152.7	-13.8	NaN	18	25.0
	0.160	0.040	142.0	-10.7	NaN	24	25.0
	0.200	0.040	133.8	-8.2	-174.32	30	25.0
	0.240	0.040	127.0	-6.8	-150.74	36	25.0

	1.641	0.040	10.2	-6.2	-155.22	250	25.0
	1.681	0.040	4.2	-6.0	-164.36	256	25.0
	1.711412	NaN	-1.4	NaN	-166.84	NaN	NaN
	1.721	0.040	-3.2	-7.4	-166.71	262	25.0
	1.761	0.040	-9.4	-6.2	-161.43	268	25.0
	1.801	0.040	-15.8	-6.4	-150.08	274	25.0
	1.841	0.040	-21.6	-5.8	-135.19	280	25.0
	1.881	0.040	-26.6	-5.0	-118.36	286	25.0
EQP2
	11.375	0.040	-197.8	-6.3	-196.06	737	25.0
	11.415	0.040	-206.0	-8.2	-234.50	744	25.0
	11.455	0.040	-216.5	-10.5	-255.80	751	25.0
	11.468579	NaN	-220.5	NaN	-256.63	NaN	NaN
	11.495	0.040	-228.2	-11.7	-253.28	760	25.0
	11.535	0.040	-238.1	-9.9	-227.90	767	25.0
	11.575	0.040	-246.1	-8.0	-187.95	775	25.0
	11.615	0.040	-252.0	-5.9	-146.99	781	25.0
	11.655	0.040	-257.1	-5.1	-116.57	787	25.0
	11.695	0.040	-261.4	-4.3	NaN	793	25.0
	11.735	0.040	-265.2	-3.8	NaN	799	25.0
	11.775	0.040	-268.6	-3.1	NaN	805	25.0
	11.815	0.040	-271.7	-3.1	NaN	812	25.0
	11.855	0.040	-274.5	-2.8	NaN	818	25.0

Sample 1/3, series 2/25/2009, 1:11:38 PM

Chemistry and Calculations

- 1) The acid carbamate is titrated with TBAH yielding two inflection points. The first equivalence point is the total amount of acid in the system, whereas the *second equivalence point* represents the amount of acid carbamate formed during the sample preparation portion of the method.
- 2) The Hydroxyl Number is given as the milligrams of potassium hydroxide equivalent to the hydroxyl content of 1 g sample:

$$OH - Number = \frac{(Q2 - B[OH - number]) * 56.1}{m}$$

where

m	=	Sample mass (g)
Q2	=	Titration consumption between first and second equivalence point (mmolL)
B[OH-number]	=	Titration consumption of the blank value (mmolL)
56.1	=	Molar mass of KOH

- 3) For a **pure compound**, the OH-number is inversely proportional to the hydroxyl equivalent weight

$$Hydroxyl\ equivalent\ weight\ (g/equivalent) = 56100 / Hydroxyl\ Number$$

and the molecular weight

$$Molecular\ weight\ (g/mol) = (56100 \times Number\ of\ OH-group\ per\ molecule) / Hydroxyl\ Number$$

- 4) The average molecular weight (g/mol) of the polyol can be calculated as follows:

$$MW = \frac{m \cdot 1000 \cdot 2}{Q2 - B[OH_number]}$$

where

m	=	Sample mass (g)
Q2	=	Titration consumption between first and second equivalence point (mmolL)
B[OH-number]	=	Titration consumption of the blank value (mmolL)
1000	=	Conversion factor that converts grams to milligrams
2	=	Two end groups per polyol chain

If we combine equations 1) with equation 3) we find that the product of the Molecular Weight and the Hydroxyl Number is always 112200:

$$MW * OH - Number = 112200$$

This relationship can be used to calculate the average molecular weight of the polyols.

Tips and hints for the titration

- This test works best with an anhydrous, stabilized with anti-oxidant 2,6-Di-tert-butylphenol (BHT) and alcohol free tetrahydrofuran.
- Acetonitrile used for the preparation of the isocyanate solution (TSI) must also be anhydrous and alcohol free.
- Special care must be taken in order to minimize the amount of water that is introduced during the sample preparation portion of the test method. Thus, appropriate conditioning is required.
- A molecular sieve trap must be placed above the isocyanate solution since TSI is very sensitive to moisture. The titrant TBAH must be protected against the intake of CO₂ from air.
- The electrode will be removed from the 50:50 water:2-propanol conditioning beaker and placed in a beaker with anhydrous, alcohol free THF or acetonitrile to remove all water from the system. The electrode will then move to the sample beaker for the titration.
- Adding too much water for the neutralization of remaining isocyanate will change the matrix of the titration. This is not recommended for maximum repeatability.

Results

Raw Polyol No. 1

Results

No.	Comment / ID	Start time	Sample size and results			
1/8	Polyol No. 1	10/22/2008 12:12:08 PM	0.2239	g		
			R1 = 539.4	mg KOH/g	Hydroxyl number	
			R2 = 203	g/equiv	Equivalent weight	
2/8	Polyol No. 1	10/22/2008 12:41:34 PM	0.2151	g		
			R1 = 547.3	mg KOH/g	Hydroxyl number	
			R2 = 205	g/equiv	Equivalent weight	
3/8	Polyol No. 1	10/22/2008 1:11:16 PM	0.2183	g		
			R1 = 558.2	mg KOH/g	Hydroxyl number	
			R2 = 201	g/equiv	Equivalent weight	
4/8	Polyol No. 1	10/22/2008 1:40:49 PM	0.1821	g		
			R1 = 558.2	mg KOH/g	Hydroxyl number	
			R2 = 201	g/equiv	Equivalent weight	
5/8	Polyol No. 1	10/22/2008 2:09:54 PM	0.2283	g		
			R1 = 547.3	mg KOH/g	Hydroxyl number	
			R2 = 205	g/equiv	Equivalent weight	
6/8	Polyol No. 1	10/22/2008 2:38:33 PM	0.2013	g		
			R1 = 563.8	mg KOH/g	Hydroxyl number	
			R2 = 199	g/equiv	Equivalent weight	
7/8	Polyol No. 1	10/22/2008 3:07:38 PM	0.1935	g		
			R1 = 561.0	mg KOH/g	Hydroxyl number	
			R2 = 200	g/equiv	Equivalent weight	
8/8	Polyol No. 1	10/22/2008 3:37:18 PM	0.2494	g		
			R1 = 550.0	mg KOH/g	Hydroxyl number	
			R2 = 204	g/equiv	Equivalent weight	

Statistics


Rx	Name	n	Mean value	unit	s	srel[%]
R1	Hydroxyl number	8	553.2	mg KOH/g	8.4	1.518
R2	Equivalent weight	8	203	g/equiv	3	1.524

Additional results

Statistics

Rx	Name	n	Mean value	unit	s	srel[%]
Polyol No. 2						
R1	Hydroxyl number	18	271.9	mg KOH/g	4.3	1.585
R2	Equivalent weight	18	206	g/equiv	3	1.568
3 decimal digits:						
Polyol No. 3						
R1	Hydroxyl number	4	56.898	mg KOH/g	0.460	0.808
R2	Equivalent weight	4	986.024	g/equiv	7.976	0.809
Polyol No. 4						
R1	Hydroxyl number	3	57.821	mg KOH/g	0.268	0.463
R2	Equivalent weight	3	970.257	g/equiv	4.490	0.463
Polyol No. 5						
R1	Hydroxyl number	3	461.388	mg KOH/g	4.155	0.901
R2	Equivalent weight	3	121.596	g/equiv	1.101	0.902
Polyol No. 6						
R1	Hydroxyl number	3	59.947	mg KOH/g	0.285	0.476
R2	Equivalent weight	3	935.834	g/equiv	4.454	0.476
Polyol No. 7						
R1	Hydroxyl number	3	116.002	mg KOH/g	0.884	0.762
R2	Equivalent weight	3	483.630	g/equiv	3.688	0.763

Titration curve



Polyol No. 1, sample series 10/22/2008 12:12:07, sample 8/8

Table of measured values

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s	Temperature °C
EQP1	0.000	NaN	22.7	NaN	NaN	0	25.0
	0.114	0.114	10.1	-12.6	NaN	5	25.0
	0.171	0.057	4.0	-6.1	NaN	12	25.0
	0.200	0.029	1.2	-2.8	NaN	18	25.0
	0.234	0.034	-2.1	-3.3	NaN	44	25.0
	0.264	0.030	-4.3	-2.2	-85.24	51	25.0
	0.321	0.057	-9.0	-4.7	-77.44	57	25.0

	1.012	0.031	-65.2	-3.1	-102.72	173	25.0
	1.041	0.029	-68.9	-3.1	-103.03	179	25.0
	1.048733	NaN	-69.7	NaN	-103.29	NaN	NaN
	1.068	0.027	-71.6	-2.7	-103.51	185	25.0
	1.099	0.031	-74.8	-3.2	-102.88	191	25.0
	1.127	0.028	-77.8	-3.0	-101.59	197	25.0
	1.155	0.028	-80.4	-2.6	-99.19	203	25.0
	1.191	0.036	-84.1	-3.7	-96.44	209	25.0
EQP2
	12.653	0.024	-306.9	-2.8	-125.22	662	25.0
	12.679	0.026	-310.1	-3.4	-129.25	669	25.0
	12.700	0.021	-313.0	-2.9	-132.64	675	25.0
	12.717160	NaN	-315.1	NaN	-133.07	NaN	NaN
	12.721	0.021	-315.6	-2.6	-132.55	681	25.0
	12.747	0.026	-319.1	-3.5	-131.58	687	25.0
	12.768	0.021	-322.0	-2.9	-129.55	693	25.0
	12.790	0.022	-324.7	-2.7	-127.59	699	25.0
	12.816	0.026	-327.8	-3.1	-124.80	705	25.0
	12.842	0.026	-331.2	-3.4	NaN	713	25.0
	12.863	0.021	-333.7	-2.5	NaN	719	25.0
	12.890	0.027	-336.7	-3.0	NaN	725	25.0
	12.919	0.029	-339.8	-3.1	NaN	732	25.0
	12.948	0.029	-342.6	-2.8	NaN	738	25.0

Polyol No. 1, sample series 10/22/2008 12:12:07, sample 8/8

Method

001 Title		Control	
Type	General titration	Control	User
Compatible with	T70 / T90	Titration addition	Dynamic
ID	m453A	dE (set value)	6.0 mV
Title	Hydroxyl number	dV (min)	0.02 mL
Author	Administrator	dV (max)	0.5 mL
Date/Time	11.02.2009 13:20:04	Mode	Equilibrium controlled
Modified at	11.02.2009 14:25:58	dE	0.5 mV
Modified by	admin	dt	3 s
Protect	No	t (min)	6 s
SOP	None	t (max)	30 s
002 Sample		Evaluation and recognition	
Number of IDs	1	Procedure	Standard
ID 1	Polyol No. 1	Threshold	50
Entry type	Weight	Tendency	Negative
Lower limit	0.2 g	Ranges	0
Upper limit	0.3 g	Add. EQP criteria	No
Density	1.0 g/mL	Termination	
Correction factor	1.0	At Vmax	25 mL
Temperature	25.0°C	At potential	No
Entry	Arbitrary	At slope	No
003 Titration stand (Rondo/Tower A)		After number of recognized EQPs	Yes
Type	Rondo/Tower A	Number of EQPs	2
Titration stand	Rondo60/1A	Combined termination criteria	No
Lid handling	No	Accompanying stating	
004 Conditioning		Accompanying stating	No
Type	Fix	Condition	No
Interval	1	013 Calculation R1	
Position	Special beaker 1	Result	Hydroxyl number
Time	10 s	Result unit	mg KOH/g
Speed	30%	Formula	$R1 = (Q2 - B[OH \text{ number}]) * C/m$
Condition	No	Constant	M/z
005 Pump		M	M[KOH]
Auxiliary reagent	THF	z	z[KOH]
Volume	10.0 mL	Decimal places	1
Condition	No	Result limits	No
006 Dispense (normal) [1]		Record statistics	Yes
Titration	Isocyanate	Extra statistical func.	No
Concentration	1.0 mol/L	Send to buffer	No
Volume	10.0 mL	Condition	No
Dosing rate	60.0 mL/min	014 Calculation R2	
Condition	No	Result	Equivalent weight
007 Stir		Result unit	g/equiv.
Speed	30%	Formula	$R2 = (56.1 * 1000) / R1$
Duration	600 s	Constant	0
Condition	No	M	M[None]
008 Pump		z	z[None]
Auxiliary reagent	H ₂ O	Decimal places	1
Volume	2.0 mL	Result limits	No
Condition	No	Record statistics	Yes
009 Stir		Extra statistical func.	No
Speed	30%	Send to buffer	No
Duration	120 s	Condition	No
Condition	No	015 Conditioning	
010 Pump		Type	Fix
Auxiliary reagent	THF	Interval	1
Volume	40.0 mL	Position	Conditioning beaker
Condition	No	Time	30 s
011 Stir		Speed	30%
Speed	30%	Condition	No
Duration	30 s	016 End of sample	
Condition	No		
012 Titration (EQP) [1]			
Titration			
Titration	TBAH		
Concentration	0.1 mol/L		
Sensor			
Type	pH		
Sensor	DG1116-Solvent		
Unit	mV		
Temperature acquisition			
Temperature acquisition	No		
Stir			
Speed	30%		
Predispense			
Mode	None		
Waiting time	0 s		
		Conditioning beaker:	50:50 water:2-propanol
		Special rinse beaker 1:	anhydrous, alcohol free THF or acetonitrile to remove all water from the electrode

Blank value determination

001 Title

Type General titration
Compatible with T70 / T90
ID m453B
Title Blank Hydroxyl number
Author Administrator
Date/Time 11.02.2009 13:15:34
Modified at 11.02.2009 14:18:03
Modified by admin
Protect No
SOP None

002 Sample

Number of IDs 1
ID 1 Blank value OH-No.
Entry type Fixed volume
Volume 60 mL
Density 1.0 g/mL
Correction factor 1.0
Temperature 25.0°C
Entry Arbitrary

003 Titration stand (Rondo/Tower A)

Type Rondo/Tower A
Titration stand Rondo60/1A
Lid handling No

004 Conditioning

Type Fix
Interval 1
Position Special beaker 1
Time 10 s
Speed 30%
Condition No

005 Pump

Auxiliary reagent THF
Volume 10.0 mL
Condition No

006 Dispense (normal) [1]

Titrant Isocyanate
Concentration 1.0 mol/L
Volume 10.0 mL
Dosing rate 60.0 mL/min
Condition No

007 Stir

Speed 30%
Duration 600 s
Condition No

008 Pump

Auxiliary reagent H₂O
Volume 2.0 mL
Condition No

009 Stir

Speed 30%
Duration 120 s
Condition No

010 Pump

Auxiliary reagent THF
Volume 40.0 mL
Condition No

011 Stir

Speed 30%
Duration 30 s
Condition No

012 Titration (EQP) [1]

Titration
Titration
Titration TBAH
Concentration 0.1 mol/L
Sensor
Type pH
Sensor DG116-Solvent
Unit mV
Temperature acquisition
Temperature acquisition No
Stir
Speed 30%
Predispense
Mode None
Waiting time 0 s

Control

Control
Titration addition User
dE (set value) Dynamic
dV (min) 2.0 mV
dV (max) 0.005 mL
Mode 0.1 mL
dE Equilibrium controlled
dt 0.5 mV
t (min) 3 s
t (max) 6 s
t (max) 30 s

Evaluation and recognition

Procedure Standard
Threshold 50
Tendency Negative
Ranges 0
Add. EQP criteria No

Termination

At Vmax 5 mL
At potential No
At slope No
After number of recognized EQPs Yes
Number of EQPs 1
Combined termination criteria No
Accompanying stating Accompanying stating No
Condition Condition No

013 Calculation R1

Result Blank value
Result unit mmol
Formula R1=Q
Constant 1
M M[none]
z z[none]
Decimal places 1
Result limits No
Record statistics Yes
Extra statistical func. No
Send to buffer No
Condition No

014 Conditioning

Type Fix
Interval 1
Position Conditioning beaker
Time 30 s
Speed 30%
Condition No

015 End of sample

016 Blank

Name OH number
Value B = Mean[R1]
Unit mmol
Limits No
Condition No

Conditioning beaker: 50:50
water:2-propanol

Special rinse beaker 1: anhydrous, alcohol
free THF
or
acetonitrile to
remove all water from
the electrode