

## Content Determination of Sodium Metabisulfite

Via a iodometric back-titration with sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3$ ) the content of sodium metabisulfite ( $\text{Na}_2\text{S}_2\text{O}_5$ ) in a solution is determined.

<b>Sample</b>	Metabisulfite, aqueous solution. (Used here: 12 mL of ~0.01 mol/L solution)
<b>Compound</b>	Sodium metabisulfite, $\text{Na}_2\text{S}_2\text{O}_5$ $M = 190.10 \text{ g/mol}$ , $z = 4$
<b>Chemicals</b>	10 mL $\text{KIO}_3$ ; 0.01 mol/L 10 mL $\text{KI}$ ; 0.06 mol/L 10 mL $\text{H}_2\text{SO}_4$ ; 0.35 mol/L
<b>Titrant</b>	Sodium thiosulfate, $\text{Na}_2\text{S}_2\text{O}_3$ $c(\text{Na}_2\text{S}_2\text{O}_3) = 0.01 \text{ mol/L}$
<b>Standard</b>	Potassium iodate, $\text{KIO}_3$
<b>Indication</b>	DM140-SC
<b>Chemistry</b>	See comments section for complete chemical reactions.
<b>Calculation</b>	$R1 = \text{VEQ}$ Consumption (mL)  $R2 = (\text{QENDDi}[2]^*6 - \text{Q}) * \text{C} / \text{m}$ $\text{C} = \text{M}/z$ Content $\text{Na}_2\text{S}_2\text{O}_5$ (g/L)
<b>Waste disposal</b>	Neutralize acidic waste before disposal
<b>Author, Version</b>	Melanie Nijman, MSG Anachem, June 2009

## Preparation and Procedures

An exact amount of metabisulfite sample solution is pipetted into a titration beaker.

Added to this are (either manually or by additional dosing units):


- 10 mL 0.01 mol/L  $\text{KIO}_3$
- 10 mL 0.06 mol/L  $\text{KI}$
- 10 mL 0.35 mol/L  $\text{H}_2\text{SO}_4$

The amounts are chosen in such a way that an excess amount of iodine is generated and that consecutively a back-titration with sodium thiosulfate can be performed, to find out the iodine consumption of the sample (metabisulfite).

If larger amounts of sample are used or if the expected concentration is higher than ~ 0.01 M the amounts of  $\text{KIO}_3$ ,  $\text{KI}$  and  $\text{H}_2\text{SO}_4$  have to be adapted such that still an excess of iodine is present for a correct back-titration with sodium thiosulfate.

## Remarks

When dissolved in water, sodium metabisulfite leads immediately to the formation of sulfur dioxide:



Literature:

[http://en.wikipedia.org/wiki/Sodium\\_metabisulfite](http://en.wikipedia.org/wiki/Sodium_metabisulfite)  
<http://de.wikipedia.org/wiki/Natriumdisulfid>

<b>Instruments</b>	T50/T70/T90 with optionally 3 extra dosing units for dispensing of auxiliary reagents
--------------------	---

<b>Accessories</b>	DV1010 burettes (3x) Titration beakers ME-101974 LabX pro titration software
--------------------	--

## Results

### Samples


1/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
2/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
3/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
4/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
5/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
6/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL
7/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	12.0 mL

Results	Comment / ID	Rx	Result	Unit	Name
1/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.408	mL	Consumption
		R2 =	1.7660	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
2/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.462	mL	Consumption
		R2 =	1.7639	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
3/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.412	mL	Consumption
		R2 =	1.7659	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
4/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.517	mL	Consumption
		R2 =	1.7617	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
5/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.529	mL	Consumption
		R2 =	1.7612	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
6/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.494	mL	Consumption
		R2 =	1.7626	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
7/7	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	R1 =	15.489	mL	Consumption
		R2 =	1.7628	g/L	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>

### Statistics

Rx	Name	n	Mean value	Unit	s	srel [%]
R1	Consumption	7	15.473	mL	0.048	0.31
R2	Content Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>	7	1.7634	g/L	0.0019	0.109

## Titration curve



## Table of measured values

	Volume mL	Increment mL	Signal mV	Change mV	1st deriv. mV/mL	Time s
EQP1	0	NaN	432.2	NaN	NaN	0
	7.143	7.143	420.8	-11.4	NaN	16
	10.714	3.571	411.6	-9.2	NaN	52
	12.500	1.786	403.6	-8.0	NaN	58
	12.900	0.400	401.8	-1.8	NaN	72
	13.300	0.400	399.3	-2.5	-9.87	76
	13.700	0.400	395.7	-3.6	-11.85	79
	14.100	0.400	391.9	-3.8	-14.41	82
	14.500	0.400	386.4	-5.5	-20.67	86
	14.900	0.400	377.8	-8.6	-43.57	90
	15.091	0.191	371.0	-6.8	-80.60	93
	15.196	0.105	367.0	-4.0	-111.55	96
	15.340	0.144	355.8	-11.2	-134.60	100
	15.390	0.050	348.4	-7.4	-136.58	104
	15.411644	NaN	341.1	NaN	-137.66	NaN
	15.440	0.050	331.6	-16.8	-131.42	110
	15.490	0.050	304.4	-27.2	-112.35	127
	15.540	0.050	300.9	-3.5	-123.37	132
	15.664	0.124	296.3	-4.6	-99.09	139
	15.968	0.304	289.6	-6.7	NaN	147
16.368	0.400	283.3	-6.3	NaN	156	
16.768	0.400	279.4	-3.9	NaN	161	
17.168	0.400	275.2	-4.2	NaN	167	
17.568	0.400	272.2	-3.0	NaN	172	

## Comments


When sodium metabisulfite is dissolved in water, the following reaction occurs:




The formed bisulfite,  $\text{NaHSO}_3$ , can be measured quantitatively via a iodometric back titration. In this titration, first an excess  $\text{I}_2$  is generated with the chemicals added to the  $\text{Na}_2\text{S}_2\text{O}_5$  sample,  $\text{KIO}_3$ ,  $\text{KI}$  and  $\text{H}_2\text{SO}_4$ :



Part of this excess  $\text{I}_2$  is used by the  $\text{Na}_2\text{S}_2\text{O}_5$  sample that is present in the titration beaker as  $\text{NaHSO}_3$  after its reaction with water:



The rest of iodine can be determined quantitatively by back-titration with sodium thiosulfate,  $\text{Na}_2\text{S}_2\text{O}_3$ :



## Method

### 001 Title

Type	General titration
Compatible with	T50 / T70 / T90
ID	Metabisulfite
Title	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> content with EQP redox
Author	admin
Date/Time	23.06.2009 08:30:27
Modified at	23.06.2009 15:23:07
Modified by	admin
Protect	No
SOP	None

### 002 Sample

Number of IDs	1
ID 1	Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub>
Entry type	Volume
Lower limit [mL]	5.0
Upper limit [mL]	15.0
Density [g/mL]	1.0
Correction factor	1.0
Temperature [°C]	25.0°C
Entry	Before

### 003 Titration stand (Manual stand)

Type	Manual stand
Titration stand	Manual stand 1

### 004 Dispense (normal) [1]

Titrant	KI
Concentration [mol/L]	0.06
Volume [mL]	10.0
Dosing rate [mL/min]	60.0
Condition	No

### 005 Dispense (normal) [2]

Titrant	KIO <sub>3</sub>
Concentration [mol/L]	0.01
Volume [mL]	10.0
Dosing rate [mL/min]	60.0
Condition	No

### 006 Dispense (normal) [3]

Titrant	H <sub>2</sub> SO <sub>4</sub>
Concentration [mol/L]	0.35
Volume [mL]	10.0
Dosing rate [mL/min]	60.0
Condition	No

### 007 Stir

Speed [%]	30
Duration [s]	10
Condition	No

### 008 Titration (EQP) [1]

#### Titrant

Titrant	Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>
Concentration [mol/L]	0.01

#### Sensor

Type	mV
Sensor	DM140-SC
Unit	mV

#### Temperature acquisition

Temperature acquisition	No
-------------------------	----

#### Stir

Speed [%]	30
-----------	----

#### Predispense

Mode	Volume
Volume [mL]	150/m
Waiting time [s]	10

#### Control

Control	User
Titrant addition	Dynamic
dE (set value) [mV]	6
dV (min) [mL]	0.05
dV (max) [mL]	0.4
Mode	Equilibrium controlled
dE [mV]	1.0
dt [s]	2
t (min) [s]	3
t (max) [s]	30

#### Evaluation and recognition

Procedure	Standard
Threshold [mV/mL]	50
Tendency	None
Ranges	0
Add. EQP criteria	No

### Termination

At Vmax [mL]	40
At potential	No
At slope	No
After number of recognized EQPs	Yes
Number of EQPs	1
Combined termination criteria	No
Accompanying stating condition	No

### 009 Calculation R1

Result	Consumption
Result unit	mL
Formula	R1=VEQ
Constant C=	1
M	M[None]
z	z[None]
Decimal places	3
Result limits	No
Record statistics	Yes
Extra statistical func.	No
Send to buffer	No
Condition	No

### 010 Calculation R1

Result	Consumption
Result unit	mL
Formula	R2=(QENDDi [2] *6-Q) *C/m
Constant C=	M/z
M	M[Metabisulfite]
z	z[Metabisulfite]
Decimal places	4
Result limits	No
Record statistics	Yes
Extra statistical func.	No
Send to buffer	No
Condition	No

### 011 End of sample