
**Operating Instruction
for NH_4^+ Ion-Selective Electrode**

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METTLER TOLEDO

A graphic element consisting of a series of parallel diagonal lines forming a triangular shape, positioned behind the company name.

Operating Instructions for Type 15 223 3000 Ammonium Electrode

Contents

Pages 3–1

1. General information

- 1.1 Technical data
- 1.2 Scope of delivery
- 1.3 Apparatus and accessories required
- 1.4 Chemical solutions required

2. Use of Electrode

- 2.1 Preparation for use

2.2 Calibration and measurement

- 2.3 Particular requirements for determining low ammonium concentrations
- 2.4 Particular requirements for determining ammonium in solutions with high ionic strengths
- 2.5 Standard addition method for ammonium-ion determination

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
6. Literatur

7. Bestellinformationen

Figure 1: Structural arrangement of an ammonium-selective electrode assembly.

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Structural arrangement of
type 15 223 3000 **Ammonium
Electrode**

- 1 S7 screw cap
- 2 plastic electrode body (POM copolymer)
- 3 Ag/AgCl internal lead-out
- 4 electrical screening
- 5 air-filled tube
- 6, 8 silicone sealing washers
- 7 plastic membrane module body (PVC)
- 10 ion-selective polymer membrane

Structural arrangement of type
373-90-WTE-ISE-S7 **Reference
Electrode**

- 1 S7 screw cap
- 11 filling port for intermediate electrolyte
- 12 internal reference electrolyte (maintenance-free)
- 13 Ag/AgCl reference element
- 14 internal (ceramic) diaphragm
- 15 ground-joint diaphragm

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1. General information

The Type 15 223 3000 Ammonium Electrode is intended for measuring the ammonium-ion content of aqueous sys-

tems. Its preferred applications are to be found in the field of food-products, soil, fertilizer and biological-fluid analysis.

1.1 Technical data

Type of ion-selective membrane:	PVC polymer-matrix
Membrane impedance:	<10 MΩ
Primary measured value:	potential in mV, proportional to $\log a\text{NH}_4^+$
Measuring range:	<10 ⁻⁶ mol/l to 1 mol/l NH_4^+
Operating temperature range:	0 to 50 °C
Sample medium:	aqueous
Sample volume:	preferable greater than 5 ml, although smaller samples may be determined in special measuring vessels.
Optimum pH range:	4 to 7
Interfering ions ¹⁾ :	cationic surface-tension agents (detergents) must be absent; the molar proportions of the following ions to the ammonium ion must lie below the following levels (ratios): $\text{K}^+/\text{NH}_4^+ < 0.05$; $\text{Na}^+, \text{Li}^+ < 5$; $\text{H}^+ < 0.3$; $\text{Rb}^+ < 0.2$; $\text{Cs}^+ < 2$; $\text{Mg}^{++} < 1000$; Ca^{++} , $\text{Ba}^{++} < 350$; $\text{Sr}^{++} < 100$.
Precision:	in the absence of interfering ions a precision of better than ± 0.5 mV corresponding to $\pm 2\%$ of the measured NH_4^+ concentration is obtainable.
Response time:	when passing from lower to higher concentrations, a few seconds.
Electrode dimensions:	length from lower edge of screw cap: 120 mm; diameter: 12 mm
Materials:	electrode body in chemically resistant plastic (POM copolymer); membrane module body in PVC; silicone sealing washers.
Working life:	the working life of the membrane module depends on the operating conditions, but is generally longer than 4 months.
Connections:	AS7 cable with appliance coupler as specified by the manufacturer of the pH/mV- or ion-meter (see leaflet).

1.2 Scope of Delivery

The electrode is supplied with fitted membrane module and filled with internal filling solution. It is thus ready for use without any special preparation.

1.3 Apparatus and accessories required

- Type 373-90-WTE-ISE-S7 Reference electrode with Ag/AgCl reference system and intermediate electrolyte chamber (other types of reference electrodes for special applications may be supplied to order).
- Ion-meter or pH/mV-meter.
- Magnetic stirrer with stirring bars.
- Electrode holder
- Measuring vessels
- Various pipettes, preferably microliter pipettes.

The following accessories may be useful, particularly when no ion-meter is available:

- Log/linear graph paper and/or pocket calculator with scientific and engineering functions.

1.4 Chemical solutions required

- Deionised water for preparing solutions when determining very low ammonium concentrations (i.e. <10⁻⁴ mol/l), preferably double-distilled water.
- Intermediate electrolyte for the reference electrode and ISA (Ionic Strength Adjuster) solution for preconditioning both calibration and sample solutions; $\text{ISA} = 0.9 \text{ mol/l Al}_2(\text{SO}_4)_3$. Prepare as follows:
Weigh out 600 g pure $\text{Al}_2(\text{SO}_4)_3 \cdot 18 \text{ H}_2\text{O}$ and dissolve with the aid of heating in about 800 ml deionised water. Allow to cool to room temperature and make up to 1000 ml in a calibrated flask.
- Standard ammonium solution (0.1 mol/l NH_4^+ or 1 g/l NH_4^+):
Dry AR grade ammonium chloride for 2 hours at 120 °C. Weigh out the following quantities:
5.349 g NH_4Cl for a 0.1 mol/l NH_4^+ solution, or
2.965 g NH_4Cl for a 1.0 g/l NH_4^+ solution, or
3.819 g NH_4Cl for a 1.0 g/l N solution
Dissolve in deionised water in a 1000 ml calibrated flask and make up to the mark.
The above solutions have a storage life of several months.

2. Use of ammonium electrode

2.1 Preparation for use

First withdraw the electrode from the protective tube, taking care not to touch the ion-selective membrane. Connect the electrode with the cable to the socket of the ion-meter or pH/mV-meter provided for ion-selective or pH electrodes, i. e. the high-impedance input.

Connect the reference electrode to the socket provided for it (4 mm banana plug or US 2 mm pin).

2.2 Calibration and measurement

Please note!

- Calibration and sample solutions should always be at the same temperature (the measuring vessel should be thermally insulated from possible heating by the magnetic stirrer – generally a sheet of cardboard or plastic foam material will suffice).
- All measurements should be carried out under the same stirring conditions, i. e. identical stirring speed, identical stirring bars, identical distance between electrodes and stirring bars, etc.
- Damage (scratches, etc.) and deposits on the surface of the ion-selective membrane will impair the response characteristics (for remedies, see § 3.).
- For total ammonium determinations the ionic strength should be the same in all calibration and sample solutions. The pH value should lie between 4 and 7. This can generally be ensured by adding ISA solution (see above). For particular requirements, see §§ 2.3 to 2.5).

Preparation of calibration solutions

- 1) Prepare solutions with different ammonium concentrations by diluting the standard stock solution with deionised water as appropriate.
- 2) To 100 parts of each of these solutions add 2 parts of ISA solution (see § 1.4)²⁾.

In each case, prepare two, and for high-

accuracy measurements at least five different calibration solutions. Ammonium solutions with concentrations below 10⁻⁴ mol/l or 2 mg/l should always be freshly prepared before each calibration. Other solutions have a storage life of about 1 month.

Preparation of sample solutions

In each case, add 2 parts of ISA solution to 100 parts of the sample solution with

an unknown ammonium content²⁾.

General instructions for determining total ammonium content by direct measurement

a) Preparatory operations

- 1) Prepare calibration solutions with different NH_4^+ concentrations. Ideally, each of these solutions should be placed in a different measuring vessel. If, however, only one such vessel is available, it should be thoroughly rinsed when changing solutions.
- 2) Rinse NH_4^+ electrode and wipe off electrode body with a paper tissue; the ion-selective membrane should only be carefully 'dabbed' off.
- 3) Reference electrode

- Fill with intermediate electrolyte up to the level of the filling port⁴⁾.
- Slightly loosen the ground-joint diaphragm (by careful turning) so that a little electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with intermediate electrolyte as necessary.
- Rinse reference electrode and wipe off with paper tissue
- 4) Fix electrode in holder and connect to measuring instrument.

b) Calibration

- Calibrate with NH_4^+ solutions of different concentrations one after another (preferably proceeding from dilute towards more concentrated solutions).
- Thoroughly rinse and wipe off electrodes (see above) before immersion in a new solution.
- 1) For calibration, immerse both ammonium and reference electrodes in the solution, and switch on the magnetic stirrer. Allow the measured value to settle down and note the final constant value.

- 2) Plot measured mV values on a calibration diagram (see Fig. 2) against the logarithm of the NH_4^+ concentrations of the corresponding solution (use log/linear graph paper):
 NH_4^+ concentration → logarithmic axis
mV value → linear (mm) axis³⁾
(The kind and frequency of recalibration operations depends on measuring conditions, and particularly on the degree of accuracy desired).

c) Measurement

For measurement, proceed as for calibration and note the final constant values. Use the calibration diagram (see above) to read off the ammonium con-

centration of the sample solutions from the measured mV values, or evaluate results with a pocket calculator³⁾.

d) After measurement

- 1) Rinse ammonium electrode. Store dry preferably in the protective tube.
- 2) Drain electrolyte out of the reference

electrode⁴⁾. Replace cap over filling port for storage.


2) The addition of ISA will produce an ionic strength of about 0.1 mol/l

3) Particularly when more than two calibration points are being used, it may be advantageous to employ a scientific/engineering pocket calculator (with least squares function) for evaluation of the calibration diagram. (See for example: S. Ebel, E. Glaser, A. Seuring, Z. Anal. Chem., 1978, 291, 108-112).

4) The life of the reference electrode will be shortened if the intermediate electrolyte chamber is always filled over the level of the internal diaphragm (leaching out of KCl).

¹⁾ At the ratios mentioned, the interference is <1%

Figure 2:
Typical calibration curve for an ammonium-selective electrode assembly in ISA solution (100 parts of deionised water + 2 parts ISA) at 20 °C.

2.3 Particular requirements for determining low ammonium concentrations

The detection limit of the ammonium electrode is determined by the low solubility in water of the ammonium-selective component in the membrane. This manifests itself in a flattening of the slope of the curve in the region of the detection limit (see figure 2).

For measurements in the curved region of the calibration graph (i. e. below about $5 \cdot 10^{-5}$ mol/l NH_4^+ or 2 mg/l NH_4^+), a calibration curve with at least 5 points

should be plotted, or an enrichment method with several standard additions, or an addition/dilution method (see literature) should be used.

For plotting the calibration curve it is advisable to make stepwise additions of a dilute standard solution to an initially ammonium-free, dilute ISA solution (standard addition method). For this purpose, the use of a strip chart recorder is recommended.

2.4 Particular requirements for determining ammonium in solutions with high ionic strengths

Sample solutions with ionic strengths markedly greater than 0.01 mol/l require very considerable additions of the conditioning solution in order to obtain comparable NH_4^+ activity coefficients in

both calibration and sample solutions. With very high salt concentrations an enrichment method (see below) should be used.

2.5 Standard addition method for ammonium-ion determination

The unknown ammonium concentration of the sample solution may be determined by adding a known quantity of ammonium to this solution.

In principle, the following procedure should be applied (taking account of the general hints given in § 2.2):

- 1) Measure the electrode potential E_P in the sample solution (preconditioned if necessary).

- 2) Add a definite quantity of ammonium (preferably in the form of a standard ammonium solution) (enrichment). This addition should produce a potential change of app. 10 to 30 mV.
- 3) Measure the electrode potential E_{PS} in the enriched solution.

The unknown ammonium concentration in the original sample solution can then be calculated from the equation:

$$C_P = \frac{C_S \cdot \left(\frac{V_S}{V_P + V_S} \right)}{10^{\frac{E_{PS} - E_P}{S}} - \left(\frac{V_P + V_S}{V_P} \right)}$$

If the dilution occasioned by the standard addition does not need to be taken

into account (i. e. if $V_S < 1/100 V_P$), equation (1) above may be simplified to:

$$C_P = \frac{C_S \cdot \frac{V_S}{V_P}}{10^{\frac{E_{PS} - E_P}{S}} - 1}$$

The terms of equations (1) and (2) above have the following significance:

C_P = NH_4^+ concentration in sample solution

C_S = NH_4^+ concentration in standard solution

V_P = volume of sample solution (ml)

V_S = volume of standard solution added (ml)

E_P = electrode potential in sample solution (mV)

E_{PS} = electrode potential in enriched solution (mV)

The slope S of the ammonium-selective electrode assembly depends primarily on the measuring temperature. Table 1

below gives the theoretical slopes for various temperatures:

Table 1:
NERNST'S factors S_{th} for different temperatures expressed in mV per 10-fold change in concentration:

T (°C)	S_{th}	T (°C)	S_{th}	T (°C)	S_{th}
0	54.20	25	59.16	50	64.12
5	55.19	30	60.15	55	65.11
10	56.18	35	61.14	60	66.10
15	57.17	40	62.13	65	67.09
20	58.17	45	63.13	70	68.09

The literature (see § 6) describes many variations of the standard addition method, in which particular attention is

paid to the exact slope of the electrode-assembly over the ranges of NH_4^+ concentration being measured.

3. Care and maintenance of the ammonium electrode

Deposits on the active surface of the membrane module may be removed by immersion the NH_4^+ electrode in deionised water for a few minutes. After rins-

ing, recondition the electrode by immersing for a few minutes in 0.01 mol/l standard NH_4^+ solution.

3.1 Storage

The ammonium electrode should be stored dry in the air, preferably in the pro-

tective tube supplied with the electrode.

3.2 Replacing membrane module

Modules with a damaged sensitive membrane, characterized by sluggish response and/or unstable electrode potential, and which are beyond regeneration by the means described in § 3 above must be replaced.

To fit a new membrane module, proceed as follows⁵⁾:

- 1) Shake the electrode vigorously with the head outwards (as with a clinical thermometer) so as to remove air bubbles from the upper part.
- 2) Place electrode with the head downwards and unscrew module from the electrode body.

- 3) Refill electrode body completely with internal filling solution ensuring that there are no air bubbles.

- 4) Fill the new membrane module with internal filling solution up to the level of the sealing washer.

- 5) Screw the membrane module right home.

- 6) Rinse the entire electrode thoroughly and then immerse it for 5 to 10 minutes in a dilute ammonium solution (about 0.01 mol/l NH_4^+).

After rinsing, the electrode is ready for use again.

3.3 Replacing electrode body

Faulty electrode bodies often produce unstable indications and then require replacement.

To fit a new electrode body, proceed as follows⁵⁾:

- 1) Unscrew and remove membrane module from the old electrode body.
- 2) Fill a certain quantity of internal filling solution into the new electrode body.
- 3) Insert the air-filled tube in the new

electrode body and then fill the latter completely with internal filling solution (without air bubbles), and carefully wipe off any superfluous solution.

- 4) Wipe dry the threaded portion of the membrane module. Refill the module with internal filling solution if necessary (see § 3.2, Item 4) above).

Proceed further as per § 3.2, Item 5) above.

4. Troubleshooting procedures

Fault:

Possible cause:	Remedy:
Value outside measuring range:	
Fault in instrument	Check with input short-circuited.
Faulty contact in electrode plug	Insert and withdraw repeatedly.
Reference electrode empty	Refill, and check ground-joint diaphragm.
Ground-joint diaphragm dry	Loosen the outer taper so as to allow a few drops of electrolyte to flow out.
Unstable indication:	
Fault in instrument	Check with input short-circuited.
Inadequate preconditioning of sample solution	Add ISA.
Faulty earthing	Check earthing of instrument and magnetic stirrer.
Static charges	Treat all plastic parts of the instrument with a suitable antistatic or wetting agent.
Ground-joint diaphragm blocked	Clean ground-joint diaphragm.
Large air bubbles in electrode body (after replacing membrane module)	Shake electrode vigorously with the membrane module outwards (as with a clinical thermometer) refit completely in the case of inverted mounting (see § 3).
Faulty electrode body	Replace electrode body (see § 3.3).

⁵⁾ The proper functioning of the electrode, even in an inverted position, is ensured if the internal filling solution is filled without air bubbles. Pressure fluctuations are absorbed by the air-filled tube inside the electrode body.

Indication drifts:

Temperature not constant

Membrane contaminated or damaged

Slope incorrect:

Calibration solution(s) contaminated

Sample preconditioning inadequate

Membrane module not screwed up tight enough

Membrane module damaged

Incorrect result despite correct calibration values:

Defective calibration curve

Incorrect sign

Incorrect calibration solution

Incorrect ISA addition

Incorrect units

Use a thermostatic jacket, or wait until all solutions are at room temperature.

Clean or replace membrane module (see § 3.2 above).

Prepare new calibration solution(s) (see § 2.2.).

Add ISA.

Tighten membrane module.

Replace membrane module (see § 3.2).

Plot mV on the linear axis and check plotting direction along the logarithmic axis.

Note sign when reading off mV.

Repeat dilution, and prepare new calibration solution if necessary.

Add ISA to calibration and sample solutions in the same proportions.

 $10^{-3} \text{ mol/l} = 18.04 \text{ mg/l as } \text{NH}_4^+$ $10^{-3} \text{ mol/l} = 14.01 \text{ mg/l as N.}$ **5. Conversion tables****The following conversion factors apply:**

Given value:	Multiply by:	To obtain:
mol/l NH_4^+	$18.04 \cdot 10^3$	mg/l NH_4^+
mol/l NH_4^+	$14.01 \cdot 10^3$	mg/l N
Given value:	Divide by:	To obtain:
mg/l NH_4^+	$18.04 \cdot 10^3$	mol/l NH_4^+
mg/l N	$14.01 \cdot 10^3$	mol/l NH_4^+

Relationships between concentration values:

mol/l NH_4^+	mg/l NH_4^+	mg/l NH_4^+
10^{-1}	1804	1401
10^{-2}	180.4	140.1
10^{-3}	18.04	14.01
10^{-4}	1.804	1.401
10^{-5}	0.1804	0.1401
10^{-6}	0.01804	0.01401

6. Literature

K. Cammann, Working with Ion-Selective Electrodes. Springer-Verlag, 1979.
D. Midgley, K. Torrance, (Potentiometric Water Analysis) John Wiley and Sons, 1978.

J. Vesely, D. Weiss, K. Stulik, (Analysis with Ion-Selective Electrodes). Ellis Horwood Ltd./John Wiley and Sons, 1978.
Further literature for special applications can be supplied on request.

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J. Vesely, D. Weiss, K. Stulik, Analysis with Ion-Selective Electrodes. Ellis Horwood Ltd./John Wiley and Sons, 1978.
Mettler-Toledo AG, Process