
**Operating Instructions
for F- ion-selective Electrode**

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METTLER TOLEDO

Operating Instructions for Type 15 215 3000
Fluoride Electrode

Contents:	
1. General information	Pages 3–8
1.1 Technical data	
1.2 Scope of Delivery	
1.3 Apparatus and accessories required	
1.4 Chemical solutions required	
2. Use of electrode	
2.1 Preparation for use	
2.2 Calibration and measurement	
2.3 Particular requirements for determining low fluoride concentrations	
2.4 Particular requirements for measuring fluoride-ion activities	
2.5 Particular requirements for determining fluoride in strongly alkaline or strongly acid solutions, as well as in solutions with high ionic strength	
3. Care and maintenance of the fluoride electrode	
3.1 Storage	
3.2 Replacing membrane module	
3.3 Replacing electrode body	
4. Troubleshooting procedures	
5. Conversion tables	
6. Literature references	
7. Ordering information	Page 20

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
Structural arrangement of Type 15 215 3000 Fluoride Electrode	Structural arrangement of Type 373-90-WTE-ISE-S7 Reference Electrode
1 S7 screw cap	1 S7 screw cap
2 plastic electrode body (POM copolymer)	11 filling port for intermediate electrolyte
3 Ag/AgCl internal lead-out	12 internal reference electrolyte (maintenance-free)
4 electrical screening	13 Ag/AgCl reference element
5 air-filled tube	14 internal (ceramic) diaphragm
6, 8, 9 silicone sealing washers	15 ground-joint diaphragm
7 plastic module body (POM copolymer)	
10 ion-selective crystal membrane (LaF ₃)	

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1. General information

The Type 15 215 3000 Fluoride Electrode is intended for measuring fluoride ion concentrations and activities in aqueous, partially aqueous and non-aqueous solutions. Its preferred applications are to be

found in the fields of drinking-water, effluent-gas, bone, cement, fish-protein, glass, phosphate-mineral and tooth-paste analysis, as well as for the determination of aluminium and phosphate.

1.1 Technical data

Type of ion-selective membrane:	LaF ₃ single crystal
Membrane impedance:	<1 MΩ
Primary measured value:	potential in mV, proportional to log a _{F⁻} (F ⁻ is also indicated as SiF ₆ ²⁻).
Measuring range:	<10 ⁻⁶ mol/l F ⁻ up to fluoride saturation.
Operating temperature range:	0 to 80 °C
Sample media:	aqueous, partially aqueous or non-aqueous.
Sample volume:	preferably greater than 5 ml, although smaller samples may be determined in special measuring vessels.
Interfering ions:	OH ⁻ (indicated along with fluoride when OH ⁻ > 0.001 F ⁻); high concentrations of lanthanum complexing anions (such as citrate, phosphate, bicarbonate) will impair the response characteristics.
Optimum pH range:	5 to 8
Precision:	in the absence of interfering ions a precision better than ±0.25 mV, corresponding to ±1% of the measured fluoride concentration is obtainable.
Response time:	when passing from lower to higher concentrations below 10 ⁻⁴ mol/l F ⁻ , about 1 minute; above 10 ⁻⁴ mol/l F ⁻ less than 1 minute; when passing from higher to lower concentrations, several minutes.
Electrode dimensions:	length from lower edge of screw cap: 120 mm; diameter: 12 mm
Materials:	electrode body and membrane module body in chemically resistant plastic (POM copolymer); silicone sealing washers.
Working life:	the working life depends on the operating conditions, but is generally longer than 1 year.
Connections:	AS7 cable with appliance coupler as specified by the manufacturer of the pH/mV or ion-meter manufacturer (see leaflet).

1.2 Scope of Delivery

The electrode is supplied with fitted membrane module and filled with internal filling solution. It is thus ready for use without any special preparation.

1.3 Apparatus and accessories required

- Type 373-90-WTE-ISE-S7 Reference Electrode with Ag/AgCl reference system and intermediate electrolyte chamber (other types of reference electrodes for special applications may be supplied to order).
 - Ion meter or pH/mV meter.
 - Magnetic stirrer with stirring bars.
 - Electrode holder.
 - Measuring vessels.
 - Various pipettes, preferably microlitre pipettes.
- The following accessories may be useful, particularly when no ion-meter is available:
- Log/linear graph paper and/or a pocket calculator with scientific and engineering functions.

1.4 Chemical solutions required

- Deionized water for preparing solutions; when determining very low fluoride concentrations (i.e. <10⁻⁴ mol/l), preferably double-distilled water.
- Intermediate electrolyte (3 mol/l KCl) for the reference electrode: (INGOLD Ordering number 20 9823 250)
- INGOLD TISAB solution for conditioning of solutions: (Ordering No. 20 9870 500)
- Standard fluoride solution 0.1 mol/l F⁻ or 1 g/l F⁻:

Prepare as follows:

Dry AR grade sodium fluoride at 120 °C for 2 hours. Weigh out 4.199 g NaF (for a 0.1 mol/l F⁻ solution), or 2.210 g NaF (for a 1.0 g/l F⁻ solution), dissolve in deionized water in a 1000 ml calibrated flask and make up to the mark.

These standard solutions may be stored for several months in plastic bottles.

2. Use of electrode

2.1 Preparation for use

First withdraw the electrode from the protective tube, taking care not to touch the ion-selective membrane. Connect the electrode with the cable to the socket of the ion-meter or mV/pH meter provided for ion-selective or pH

electrodes, i.e. the high-impedance input.

Connect the reference electrode to the socket provided for it (4 mm banana plug or US 2 mm pin).

2.2 Calibration and measurement

Please note!

- Calibration and sample solutions should always be at the **same temperature** (the measuring vessel should be thermally insulated from possible heating by the magnetic stirrer – generally a sheet of cardboard or plastic foam material will suffice).
- All measurements should be carried out under the **same stirring conditions**, i.e. identical stirring speed, identical stirring bars, identical distance between electrodes and stirring bars, etc.
- **Damage (scratches, etc.)** and deposits on the surface of the sensor crystal will impair the response characteristics (for remedies, see § 3).
- The F⁻ electrode should not be left immersed in solutions containing TISAB for long periods, since the crystal membrane will be attacked by complexing agents such as EDTA (the same applies to citrates, etc.).
- For total fluoride determinations the **ionic strength should be the same** in all calibration and sample solutions, and all **residual bound fluoride** in the sample should be liberated. The **pH value should lie between 5 and 8**. This can generally be ensured by adding TISAB (see above). (For particular requirements, see §§ 2.3 to 2.6).

Preparation of calibration solutions

- 1) Prepare solutions with different fluoride concentrations by diluting the standard stock solution with deionized water as appropriate.
- 2) To 10 parts of each of these solutions add 1 part of INGOLD TISAB¹⁾ (see § 7).

In each case, prepare two, and for high-

Preparation of sample solutions

In each case, add 1 part of INGOLD TISAB to 10 parts of the sample solution

accuracy measurements at least five different calibration solutions.

Solutions with fluoride concentrations below 10⁻⁴ mol/l or 1.9 mg/l should always be freshly prepared before each calibration. Stronger solutions can be stored for about a month in plastic bottles.

with an unknown fluoride content¹⁾.

General instructions for determining total fluoride content

a) Preparatory operations

- 1) Prepare calibration solutions with different fluoride concentrations. Ideally, each of these solutions should be placed in a different measuring vessel. If, however, only one such vessel is available, it should be thoroughly rinsed when changing solutions.
- 2) Rinse F⁻ electrode and wipe off electrode body with paper tissue; the ion-selective membrane should only be carefully 'dabbed' off.
- 3) Reference electrode:
 - Remove cap from electrolyte filling port.
 - Slightly loosen the ground-joint diaphragm (by careful turning) so that a little electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with 3 mol/l KCl electrolyte as necessary.
 - Rinse reference electrode and wipe off with paper tissue.
- 4) Fix electrodes in holder and connect to measuring instrument.

b) Calibration

- Calibrate with F⁻ solutions of different concentrations one after another (preferably proceeding from dilute towards more concentrated solutions).
- Thoroughly rinse and wipe off electrodes (see above) before immersion in a new solution.
- 1) For calibration, immerse both fluoride and reference electrodes in the solution, and switch on the magnetic stirrer. Allow the measured value to settle down and note the final constant value.
- 2) Plot measured mV values on a calibration diagram (see Fig. 2) against the logarithm of the fluoride concentration of the corresponding solution (use log/linear graph paper):
F⁻ concentration → logarithmic axis
mV value → linear (mm) axis²⁾
(The kind and frequency of recalibration operations depends on measuring conditions, and particularly on the degree of accuracy desired).

c) Measurement

For measurement, proceed as for calibration and note the final constant values. Use the calibration diagram (see above) to read off the fluoride concentra-

tions of the sample solutions from the measured mV values, or evaluate results with a pocket calculator²⁾.

¹⁾ The addition of TISAB will adjust the pH value to between 5 and 8 and set the ionic strength to about 0.5 mol/l.


²⁾ Particularly when more than two calibration points are being used, it may be advantageous to employ a scientific/engineering pocket calculator (with least squares function) for evaluation of the calibration diagram. (See, for example, S. Ebel, E. Glaser, A. Seuring, *Z. Anal. Chem.*, 1978, 291, 108–112).

c) After measurement

- 1) Rinse fluoride electrode. Store dry or immersed in a dilute solution (not in TISAB!), preferably in the protective tube.
- 2) Rinse reference electrode. Replace cap over filling port for storage.

Figure 2:

Typical calibration curve for a fluoride-selective electrode assembly in INGOLD TISAB solution (10 parts of deionized water + 1 part TISAB) at 20 °C.

2.3 Particular requirements for determining low fluoride concentrations

The theoretical detection limit of the fluoride electrode is determined by the very low solubility of the lanthanum fluoride crystal. The practical detection limit is generally somewhat higher, and is primarily due to fluoride traces in the reagents used, as well as to the increasingly sluggish response at greater dilutions. If INGOLD TISAB is used for adjusting the pH value and ionic strength of the calibration and sample solutions, linear calibration curves down to 10⁻⁶ mol/l or 0.02 mg/l F⁻ can be obtained, the detection limit being at slightly above 10⁻⁷ mol/l or 0.002 mg/l F⁻. To avoid erroneous results due to impurities, it is advisable when carrying out fluoride determinations in the range below 10⁻⁵ mol/l or 0.2 mg/l F⁻ to apply a

standard addition method (see literature) involving several standard additions, and to use a strip chart recorder. When doing this, it is best to start with a fluoride-free (blank solution) and to renew this as required until the measured electrode potential no longer increases with a further change of solution. This solution should then be enriched by successively adding definite quantities of the (liquid or solid) sample substance, and of standard fluoride solution, reading off the electrode potential after each addition. The measured value from the (blank solution) will then enable the practical detection limit and hence the accuracy of the fluoride determination to be checked.

2.4 Particular requirements for measuring fluoride-ion activities

Fluoride-ion activities can be measured directly in the sample solution, without the need for preconditioning by adding TISAB. Under these conditions the only ions acting on the fluoride-selective membrane are free fluoride ions which are neither combined with hydrogen ions (as HF or HF₂⁻) nor with fluoride

complexing ions such as aluminium or iron. If the electrode assembly is first calibrated in dilute standard solutions of sodium fluoride with reference to the activity coefficients given in the table below, the measuring results will give the free fluoride-ion activities:

NaF concentration (mol/l)	Fluoride ion activity (mol/l)
10 ⁻¹	0.77 · 10 ⁻¹
10 ⁻²	0.91 · 10 ⁻²
10 ⁻³	0.97 · 10 ⁻³
10 ⁻⁴	1.00 · 10 ⁻⁴


Thus, the procedure for determining fluoride-ion activity is distinguished from that for total fluoride determination by the following features:

- 1) The electrode assembly is calibrated by means of diluted standard NaF solutions (without TISAB preconditioning), with reference to the activity coefficients.

- 2) Direct measurement in the sample solution (also without TISAB preconditioning). However, the pH value of the sample solution should be less than 9, so as to prevent hydroxyl ion interference.

Figure 3:

Fluoride-ion activity as a function of total ionic strength. The data in the region below 10⁻¹ mol/l have been calculated according to the Debye-Hückel theory using the values given by Kielland in the J. Amer. Chem. Soc. 1937, 59, 1675. Data in the range from 10⁻¹ mol/l to 1 mol/l are based on actual electrode measurements.

2.5 Particular requirements for determining fluoride in strongly alkaline or strongly acid solutions, as well as in solutions with high ionic strength

In order to ensure error-free determinations of total fluoride content, the pH value of both calibration and sample solutions must lie between 5 and 8, and the ionic strength of all solutions concerned should be approximately equal. However, in strongly alkaline or strongly acid media, the TISAB additions mentioned in § 2.2 above may sometimes be inadequate to ensure that such conditions obtain in the sample solution.

In such cases it is permissible to increase the addition of TISAB to the calibration and sample solutions. However, it must be borne in mind that the ionic strengths of the solutions as well as the concentrations of complexing agents (Titriplex IV, CDTA) will be modified in accordance with the following table:

Volumetric parts of solution	+	Volumetric parts of INGOLD TISAB	Ionic strength produced by TISAB (mol/l)	Resultant CDTA concentration (mol/l) ¹⁾
10	+	1	env. 0.45	0.0045
5	+	1	env. 0.83	0.0083
2	+	1	env. 1.65	0.0165
1	+	1	env. 2.50	0.0250

¹⁾ Hereby, metallic ions will be complexed up to the same total concentration.

Note: The detection limit of the fluoride electrode may be shifted up to a higher level of F⁻ concentration with increasing addition of TISAB.

3. Care and maintenance of the fluoride electrode

Deposits on the crystal may be removed by careful wiping with a dampened paper tissue. When doing this, the greatest care should be taken to avoid mechanical damage (scratches, etc.) to the surface of the sensitive crystal. Finally, conditioning for a few minutes in a dilute

fluoride solution is recommended. Modules whose crystal surface has become dulled or scratched beyond regeneration by this method should be sent to our works or your local INGOLD Service Agent for repolishing.

3.1 Storage

The fluoride electrode may be stored dry in the air or immersed in a dilute fluoride solution. For that purpose it is best to make use of the protective tube supplied with the electrode on delivery.

Electrodes should never be stored immersed in solutions containing TISAB! (See § 2.2 above).

3.2 Replacing membrane module

Membrane modules with a damaged crystal, characterised by sluggish response and/or unstable electrode potential, which are beyond regeneration by the means described in § 3 above, must be replaced.

To fit a new membrane module, proceed as follows³⁾:

- 1) Shake the electrode vigorously with the head outwards (as with a clinical thermometer) so as to remove air bubbles from the upper part.
- 2) Place electrode with the head downwards and unscrew the old module from the electrode body.

- 3) Refill electrode body completely with internal filling solution, ensuring that there are no air bubbles.
- 4) Fill the new membrane module with internal filling solution up to the level of the sealing washer.
- 5) Screw the membrane module right home.

- 6) Rinse the entire electrode thoroughly and then immerse it for 5 to 10 minutes in a dilute fluoride solution (about 0.01 mol/l F⁻).

After rinsing, the electrode is ready for use again.

3.3 Replacing electrode body

Faulty electrode bodies often produce unstable indications and then require replacement.

To fit a new electrode body, proceed as follows³⁾:

- 1) Unscrew and remove membrane module from the old electrode body.

- 2) Fill a certain quantity of internal filling solution into the new electrode body.
- 3) Insert the air-filled tube in the new electrode body and then fill the latter completely with internal filling solution (without air bubbles), and carefully wipe off any superfluous solution.

4) Wipe dry the threaded portion of the membrane module. Refill the module with internal filling solution if necessary (see § 3.2 Item 4) above).

Proceed further as per § 3.2, Item (5) onwards.

4. Troubleshooting procedures

Fault:

Possible cause:	Remedy:
Value outside measuring range:	
Fault in instrument	Check with input short-circuited.
Faulty contact in electrode plug	Insert and withdraw repeatedly.
Reference electrode empty	Refill, and check ground-joint diaphragm.
Ground-joint diaphragm dry	Loosen the outer taper so as to allow a few drops of electrolyte to flow out.
Unstable indication:	
Fault in instrument	Check with input short-circuited.
Inadequate preconditioning of sample solution	Add TISAB.
Faulty earthing	Check earthing of instrument and magnetic stirrer.
Static charges	Treat all plastic parts of the instrument with a suitable antistatic or wetting agent.
Ground-joint diaphragm blocked	Clean ground-joint diaphragm.
Large air bubbles in electrode body (after replacing membrane module)	Shake electrode vigorously with the membrane module outwards (as with a clinical thermometer); refit completely in the case of inverted mounting (see § 3.2 above).
Faulty electrode body	Replace electrode body (see § 3.3 above).
Temperature not constant	Use a thermostatic jacket, or wait until all solutions are at room temperature.
Sample preconditioning inadequate	Add TISAB.
Adsorption on surfaces of glass vessels	Use plastic vessels.
Slope incorrect:	
Calibration solution(s) contaminated or incorrectly diluted	Prepare new calibration solution(s) (see § 2.2 above).
Sample preconditioning inadequate	Add TISAB.
Membrane module not screwed up tight enough	Tighten membrane module.
Membrane module damaged	Replace membrane module (see § 3 above).
Incorrect result despite correct calibration values:	
Defective calibration curve	Plot mV on the linear axis and check plotting direction along the logarithmic axis.
Incorrect sign	Note sign when reading off mV.
Incorrect calibration solution	Repeat dilution, and prepare new calibration solution if necessary.
Incorrect TISAB addition	Add TISAB to calibration and sample solutions in the same proportion.
Incorrect units	$10^{-3} \text{ mol/l} = 19 \text{ mg/l F}^-$ $= 41.99 \text{ mg/l NaF}$

5. Conversion tables

The following conversion factors apply:

Given value:	Multiply by:	To obtain:
mol/l F^-	$19.00 \cdot 10^3$	mg/l F^-
mol/l F^-	$41.99 \cdot 10^3$	mg/l NaF
Given value:	Divide by:	To obtain:
mg/l F^-	$19.00 \cdot 10^3$	mol/l F^-
mg/l NaF	$41.99 \cdot 10^3$	mol/l F^-

Relationships between concentration values:

mol/l F^-	mg/l F^-	mg/l NaF
10^{-1}	1900	4199
10^{-2}	190	419.9
10^{-3}	19	41.99
10^{-4}	1.9	4.199
10^{-5}	0.19	0.4199
10^{-6}	0.019	0.04199

6. Literature

K. Cammann, Working with Ion-Selective Electrodes. Springer-Verlag, 1979.
D. Midgley, K. Torrance, Potentiometric Water Analysis. John Wiley and Sons, 1978.

J. Vesely, D. Weiss, K. Stulik, Analysis with Ion-Selective Electrodes. Ellis Horwood Ltd./John Wiley and Sons, 1978.
Further literature for special applications can be supplied on request.

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7. Ordering information

The following components of the fluoride electrode can be supplied separately.

- Membrane kit (replacement F⁻ membrane module + 12 ml internal filling solution): Ordering No. **15 215 2000**
 - 500 ml INGOLD TISAB solution: Ordering No. **20 9870 500**
 - Type 373-90-WTE-ISE-S7 Reference Electrode: Ordering No. **10 373 3145**
 - 250 ml intermediate electrolyte solution (3 mol/l KCl): Ordering No. **20 9823 250**
 - Replacement electrode body: Ordering No. **15 200 1000**
- AS7 connecting cables for both fluoride and reference electrodes can be supplied on request to suit the measuring instrument in use.

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20 606 1091

Mettler-Toledo AG, Process, Industrie Nord, CH-8902 Urdorf
Tel. (01) 736 22 11, Fax (01) 736 26 36
Subject to technical changes. 12/95.
© Mettler-Toledo AG. Printed in Switzerland

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- CH** Mettler-Toledo (Schweiz) AG
- D** Mettler-Toledo Prozeßanalytik GmbH
- UK** Mettler-Toledo Ltd., Analytical
- USA** Mettler-Toledo Process Analytical Inc.