**Operating Instruction
for Cl^- Ion-Selective Electrode**

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METTLER TOLEDO

A graphic element consisting of a series of parallel diagonal lines forming a triangular shape, positioned behind the company name.

Operating Instructions for Type 15 213 3000 Chloride Electrode

Contents

Pages 3-8

1. General information
 - 1.1 Technical data
 - 1.2 Scope of delivery
 - 1.3 Apparatus and accessories required
 - 1.4 Chemical solutions required
2. Use of electrode
 - 2.1 Preparation for use
 - 2.2 Calibration and measurement
 - 2.3 Particular requirements for determining low chloride concentrations
 - 2.4 Particular requirements for determining chloride in ground waters, tap water and effluents
 - 2.5 Particular requirements for determining chloride in the presence of sulphide
 - 2.6 Particular requirements for determining chloride in solutions with high ionic strengths
 - 2.7 Standard addition method for chloride determination
3. Care and maintenance of the chloride electrode
 - 3.1 Storage
 - 3.2 Replacing membrane module
 - 3.3 Replacing electrode body
4. Troubleshooting procedures
5. Conversion tables
6. Literature references
7. Ordering information

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
Seite 20

Figure 1: Structural arrangement of a chloride-selective electrode assembly.

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Structural arrangement of Type 15 213 3000 chloride electrode

- | | |
|---------|---|
| 1 | S7 screw cap |
| 2 | plastic electrode body (POM copolymer) |
| 3 | Ag/AgCl internal lead-out |
| 4 | electrical screening |
| 5 | air-filled tube |
| 6, 8, 9 | silicone sealing washers |
| 7 | plastic membrane module body (POM copolymer) |
| 10 | ion-selective membrane (compacted AgCl/Ag ₂ S) |

Structural arrangement of Type 373-90-WTE-ISE-S7 reference electrode

- | | |
|----|---|
| 1 | S7 screw cap |
| 11 | filling port for intermediate electrolyte |
| 12 | internal reference electrolyte (maintenance-free) |
| 13 | Ag/AgCl reference element |
| 14 | internal ceramic diaphragm |
| 15 | ground-joint diaphragm |

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1. General information

The Type 15 213 3000 Chloride Electrode is intended for measuring the chloride ion concentrations in aqueous, partially aqueous and non-aqueous solutions. Its preferred applications are to be found in

the fields of ground-water, drinking water and effluent analysis, as well as that of highly purified water, food products, soil samples and biological fluids (such as urine and perspiration, etc.).

1.1 Technical data

Type of ion-selective membrane:	compacted AgCl/Ag ₂ S
Membrane impedance:	<100 k Ω
Primary measured value:	potential in mV, proportional to $\log a_{Cl^-}$
Measuring range:	<5 · 10 ⁻⁶ mol/l to >1 mol/l Cl ⁻
Operating temperature range:	0 to 80 °C
Sample media:	aqueous, partially aqueous or non-aqueous.
Sample volume:	preferably greater than 5 ml, although smaller samples may be determined in special measuring vessels.
Optimum pH range:	2 to 11
Interfering ions ¹⁾ :	S ²⁻ must be absent; the molar proportions of the following ions to the chloride ion must lie below the following levels (ratios): CN ⁻ /Cl ⁻ <2 · 10 ⁻⁷ , I ⁻ <5 · 10 ⁻⁷ , Br ⁻ <3 · 10 ⁻³ , S ₂ O ₃ ²⁻ <0.01, NH ₃ <0.12, OH ⁻ <80.
Precision:	in the absence of interfering substances a precision of better than ±0.5 mV corresponding to ±2% of the measured Cl ⁻ concentration is obtainable.
Response time:	when passing from lower to higher concentrations, below 10 ⁻³ mol/l Cl ⁻ , several minutes, and above 10 ⁻³ mol/l Cl ⁻ , less than 1 minute; when passing from higher to lower concentrations, several minutes.
Electrode dimensions:	Length from lower edge of screw cap: 120 mm; diameter: 12 mm
Materials:	electrode body and membrane module body in chemically resistant plastic (POM copolymer); silicone sealing washers.
Working life:	the working life depends on the operating conditions, but is generally longer than 1 year.
Connections:	AS7 cable with appliance coupler as specified by the manufacturer of the pH/mV- or ion-meter (see leaflet).

1.2 Scope of Delivery

The electrode is supplied with fitted membrane module and filled with internal filling solution. It is thus ready for use without any special preparation.

1.3 Apparatus and accessories required

- Type 373-90-WTE-ISE-S7 Reference Electrode with Ag/AgCl reference system and intermediate electrolyte chamber (other types of reference electrodes for special applications may be supplied to order).
- Ion-meter or pH/mV-meter.
- Magnetic stirrer with stirring bars.
- Electrode holder.
- Measuring vessels.
- Various pipettes, preferably microliter pipettes.

The following accessories may be useful, particularly when no ion-meter is available:

- Log/linear graph paper and/or pocket calculator with scientific and engineering functions.

1.4 Chemical solutions required

- Deionized water for preparing solutions; when determining very low chloride concentrations (i.e. <10⁻⁴ mol/l), preferably double-distilled water.
- Intermediate electrolyte (1 mol/l KNO₃) for the reference electrode (Ordering No. 20 9813 250).
- ISA (Ionic Strength Adjuster) solution (5 mol/l NaNO₃) for conditioning sample solutions:
Prepare as follows:
Weigh out 42.5 g AR grade sodium nitrate, dissolve in deionized water and make up to 100 ml.
- Standard chloride solution (0.1 mol/l Cl⁻ or 1 g/l Cl⁻):
Dry AR grade sodium chloride for 2 hours at 120 °C.
Weigh out the following quantities:
5.844 g NaCl for a 0.1 mol/l Cl⁻ solution, or:
1.648 g NaCl for a 1.0 g/l Cl⁻ solution.
Dissolve in deionized water in a 1000 ml calibrated flask and make up to the mark.

The above solutions have a storage life of several months.

2. Use of electrode

2.1 Preparation for use

First withdraw the electrode from the protective tube, taking care not to touch the ion-selective membrane. Connect the electrode with the cable to the socket of the ion-meter or pH/mV-meter provided for ion-selective or pH electrodes, i.e. the high-impedance input.

Connect the reference electrode to the socket provided for it (4 mm banana plug or US 2 mm pin).

2.2 Calibration and measurement

Please note!

- Calibration and sample solutions should always be at the same temperature (the measuring vessel should be thermally insulated from possible heating by the magnetic stirrer – generally a sheet of cardboard or plastic foam material will suffice).
- All measurements should be carried out under the same stirring conditions, i.e. identical stirring speed, identical stirring bars, identical distance between electrodes and stirring bars, etc.
- All measurements should be carried out under the same conditions of illumination, preferably in light protected vessels.
- Damage (scratches, etc.) and deposits on the surface of the

ion-selective membrane will impair the response characteristics (for remedies, see § 3 below).

- The Cl⁻ electrode should not be left immersed in concentrated chloride solutions for long periods, since the membrane surface will be attacked (leaching out of AgCl by complexing of silver).

- For total chloride determinations the ionic strength should be the same in all calibration and sample solutions. The pH value should lie between 2 and 11. This can generally be ensured by adding ISA solution (see above). (For particular requirements, see §§ 2.3 to 2.5).

Preparation of calibration solutions

- 1) Prepare solutions with different chloride concentrations by diluting the standard stock solution with deionized water as appropriate.
- 2) To 100 parts of each of these solutions add two parts of ISA solution (see § 1.4 above)²⁾.

In each case, prepare two, and for high-

accuracy measurements at least five different calibration solutions.

Solutions with chloride concentrations below 10⁻⁴ mol/l or 3.5 mg/l should always be freshly prepared before each calibration. Other solutions can be stored for about a month.

Preparation of sample solutions

In each case, add 2 parts of ISA solution to 100 parts of the sample solution with

an unknown chloride content²⁾.

General instructions for determining total chloride content by direct measurement

a) Preparatory operations

- 1) Prepare calibration solutions with different chloride concentrations. Ideally, each of these solutions should be placed in a different measuring vessel. If, however, only one such vessel is available, it should be thoroughly rinsed when changing solutions.
- 2) Rinse Cl⁻ electrode and wipe off electrode body with a paper tissue; the ion-selective membrane should only be carefully 'dabbed' off.
- 3) Reference electrode
- Fill with intermediate electrolyte up

to the lower edge of the filling port⁴⁾.

- Slightly loosen the ground-joint diaphragm (by careful turning) so that a little electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with 1 mol/l KNO₃ intermediate electrolyte as necessary.

- Rinse reference electrode and wipe off with paper tissue.

- 4) Fix electrode in holder and connect to measuring instrument.

b) Calibration

- Calibrate with Cl⁻ solutions of different concentrations one after another (preferably proceeding from dilute towards more concentrated solutions).

- Thoroughly rinse and wipe off electrodes (see above) before immersion in a new solution.

- 1) For calibration, immerse both chloride and reference electrodes in the solution, and switch on the magnetic stirrer. Allow the measured value to settle down and note the final constant value.

- 2) Plot measured mV values on a calibration diagram (see Fig. 2) against the logarithm of the chloride concentration of the corresponding solution (use log/linear graph paper):
Cl⁻ concentration → logarithmic axis
mV value → linear (mm) axis³⁾

(The kind and frequency of recalibration operations depends on measuring conditions, and particularly on the degree of accuracy desired).

c) Measurement

For measurement, proceed as for calibration and note the final constant values. Use the calibration diagram (see above) to read off the chloride concen-

trations of the sample solutions from the measured mV values, or evaluate results with a pocket calculator³⁾.

²⁾ The addition of ISA will produce an ionic strength of about 0.1 mol/l.

³⁾ Particularly when more than two calibration points are being used, it may be advantageous to employ a scientific/engineering pocket calculator (with a least-squares function) for evaluation of the calibration diagram. (See, for example: S. Ebel, E. Glaser, A. Seuring, *Z. Anal. Chem.*, 1978, 291, 108–112).


¹⁾ At the ratios mentioned, the interference is <1%.

d) After measurement

- 1) Rinse chloride electrode. Store dry and in the dark (not immersed in chloride solution), preferably in the protective tube.
- 2) Drain electrolyte out of the reference electrode⁴⁾. Replace cap over filling port for storage.

Figure 2:

Typical calibration curve for a chloride-selective electrode assembly in ISA solution (100 parts of deionized water + 2 parts ISA) at 20°C.

2.3 Particular requirements for determining low chloride concentrations

The detection limit of the chloride electrode is determined by the low solubility of the silver chloride contained in the membrane. This manifests itself in a flattening off of the slope of the curve in the region of the detection limit (see Figure 2 above).

For measurements in the curved region of the calibration graph (i. e. below about $2 \cdot 10^{-4}$ mol/l Cl^- or 10 mg/l Cl^-), a cali-

bration curve with at least 5 points should be plotted, or an enrichment method with several standard additions, or again an addition/dilution method (see literature) should be used.

To plot the calibration curve it is best to add several increments of a standard solution to a dilute, chloride-free solution of ISA (standard addition).

Lowering the detection limit:

The solubility of AgCl and hence the detection limit may be lowered by adding, for example, AR grade isopropanol in a proportion of 1:1. At the same time, the straight portion of the calibration curve is extended. The same effects may be ob-

tained by cooling all calibration and sample solutions to below room temperature. A combination of both methods enables the range of measurement to be extended down to Cl^- concentrations of below 10^{-6} mol/l.

2.4 Particular requirements for determining chloride in ground-waters, tap water and effluents

Interference is particularly likely to be caused with such waters, which contain alkaline-earth ions, owing to deposits of hydroxides and carbonates on the surface of the compacted chloride selective membrane. Such interference may be

prevented by using a slightly acid citrate buffer as sample conditioning solution instead of the usual ISA solution.

This solution should be added to calibration and sample solutions in a proportion of 1:1.

X^- conditioning solution:

Weigh out 106 g AR grade (or low-chloride) citric acid, place in a 2 l beaker and dissolve in about 800 ml deionized water, stirring as required, and then carefully add 40 g AR grade (or low-chloride) sodium hydroxide (wear protective

goggles!). After cooling to room temperature, pour into a 1000 ml calibrated flask and make up to the mark with deionized water. The pH value should be about 5.

2.5 Particular requirements for determining chloride in the presence of sulphide

To prevent sulphide from interfering with the chloride determination, it must first be oxidized before immersing the electrodes. This can be achieved by adding the special $\text{X}^-(\text{S}^{2-})$ conditioning solution described below in a proportion of 1:1 to

all calibration and sample solutions instead of the usual ISA solution. Allow the solutions so treated to stand for 10 minutes before measurement so as to ensure complete oxidation of the sulphide.

$\text{X}^-(\text{S}^{2-})$ conditioning solution:

Weigh out 62 g boric acid, place in a 2 l beaker, add about 800 ml deionized water and stir (the boric acid will not completely dissolve). Now add 10% AR grade sodium hydroxide solution, bit by bit, with stirring, until the boric acid is com-

pletely dissolved and a pH value of 8.5 is obtained. Finally add 100 ml of 30% (= 100 volume) hydrogen peroxide, pour into a 1000 ml calibrated flask and make up to the mark with deionized water.

2.6 Particular requirements for determining chloride in solutions with high ionic strengths

Sample solutions with ionic strengths markedly greater than 0.01 mol/l require correspondingly greater additions of the conditioning solution in order to obtain comparable Cl^- activity coefficients in

both calibration and sample solutions. With very high salt concentrations an enrichment method (see below) should be used.

2.7 Standard addition method for chloride determination

The unknown chloride concentration of the sample solution may be determined by adding a known quantity of chloride to this solution.

In principle, the following procedure should be applied (taking account of the general hints given in § 2.2 above):

- 1) Measure the electrode potential E_P in the sample solution (preconditioned if necessary).
- 2) Add a definite quantity of chloride, preferably in the form of a standard chloride solution (enrichment). This addition should produce a potential change of 10 to 30 mV.
- 3) Measure the electrode potential E_{PS} in the enriched solution.

The unknown chloride concentration in the original sample solution can then be calculated from the following equation:

$$C_P = \frac{C_S \cdot \left(\frac{V_S}{V_P + V_S} \right)}{10^{\frac{E_{PS} - E_P}{S}} - \left(\frac{V_P + V_S}{V_P} \right)}$$

If the dilution occasioned by the stan-

dard addition needs not be taken into account (i. e. if $V_S < 1/100 V_P$), equation (1) above may be simplified to:

$$C_P = \frac{C_S \cdot \frac{V_S}{V_P}}{10^{\frac{E_{PS} - E_P}{S}} - 1}$$

The terms of equations (1) and (2) above have the following significance:

C_P = Cl^- concentration in sample solution

C_S = Cl^- concentration in standard solution

V_P = volume of sample solution (ml)

V_S = volume of standard solution added (ml)

E_P = electrode potential in sample solution (mV)

E_{PS} = electrode potential in enriched solution (mV)

The slope S of the chloride-selective electrode assembly depends primarily on the measuring temperature. Table 1 below gives the theoretical slopes for various temperatures.

Table 1:

NERNST'S factors S_{th} for different temperatures expressed in mV per 10-fold change in concentration:

T (°C)	S_{th}	T (°C)	S_{th}	T (°C)	S_{th}
0	54.20	25	59.16	50	64.12
5	55.19	30	60.15	55	65.11
10	56.18	35	61.14	60	66.10
15	57.17	40	62.13	65	67.09
20	58.17	45	63.13	70	68.09

The literature (see § 6 below) describes many variations of the standard addition method, in which particular attention is paid to the exact slope of the electrode-assembly over the ranges of Cl^- concentration being measured.

3. Care and maintenance of the chloride electrode

If necessary, the active surface of the membrane module may be rubbed wet with the finest grade of wet-and-dry abrasive paper and then polished with a suitable abrasive paste, in order to remove deposits and damage due to chemical corrosion. After being cleaned in this way, there may be a small shift (a

few mV) in the electrode zero point, which, however, will regain its initial position within a few days. For this reason, during the period immediately after cleaning the Cl^- electrode should be recalibrated at frequent intervals and kept stored immersed in deionized water.

3.1 Storage

The Cl^- electrode may either be stored dry in the air or immersed in deionized water. For that purpose it is best to make use of the protective tube supplied with

the electrode on delivery. Do not store the electrode in solutions containing chloride (see § 2.2 above).

3.2 Replacing membrane module

Modules in which the compacted AgCl/Ag₂S membrane has suffered damage, characterized by sluggish response and/or unstable electrode potential, which are beyond regeneration by the means described above, must be replaced. Scrapped modules should be disposed off in the special dustbin reserved for spent batteries, mercury cells, etc., to avoid pollution of the environment.

To fit a new membrane module, proceed as follows⁵⁾:

- 1) Shake the electrode vigorously with the head outwards (as with a clinical thermometer) so as to remove air bubbles from the upper part.
- 2) Place electrode with the head down-

wards and unscrew module from the electrode body.

- 3) Refill electrode body completely with internal filling solution ensuring that there are no air bubbles.
- 4) Fill the new membrane module with internal filling solution up to the level of the sealing washer.
- 5) Screw the membrane module right home.
- 6) Rinse the entire electrode thoroughly and then immerse it for 5 to 10 minutes in a dilute chloride solution (about 0.01 mol/l Cl^-).

After rinsing, the electrode is ready for use again.

⁴⁾ The life of the reference electrode will be shortened if the intermediate electrolyte chamber is always filled over the level of the internal diaphragm (leaching out of KCl).

3.3 Replacing electrode body

Faulty electrode bodies often produce unstable indications and then require replacement.

To fit a new electrode body, proceed as follows⁵⁾:

- 1) Unscrew and remove membrane module from the old electrode body.
- 2) Fill a certain quantity of internal filling solution into the new electrode body.
- 3) Insert the air-filled tube in the new

electrode body and then fill the latter completely with internal filling solution (without air bubbles), and carefully wipe off any superfluous solution.

- 4) Wipe dry the threaded portion of the membrane module. Refill the module with internal filling solution if necessary (see § 3.2, Item 4) above).

Proceed further as per § 3.2, Item 5) above.

4. Troubleshooting procedures

Fault:

Possible cause:	Remedy:
Value outside measuring range:	
Fault in instrument	Check with input short-circuited.
Faulty contact in electrode plug	Insert and withdraw repeatedly.
Reference electrode empty	Refill, and check ground-joint diaphragm.
Ground-joint diaphragm dry	Loosen the outer taper so as to allow a few drops of electrolyte to flow out.
Unstable indication:	
Fault in instrument	Check with input short-circuited.
Inadequate preconditioning of sample solution	Add ISA.
Faulty earthing	Check earthing of instrument and magnetic stirrer.
Static charges	Treat all plastic parts of the instrument with a suitable antistatic or wetting agent.
Ground-joint diaphragm blocked	Clean ground-joint diaphragm.
Large air bubbles in electrode body (after replacing membrane module)	Shake electrode vigorously with the membrane module outwards (as with a clinical thermometer); refit completely in the case of inverted mounting (see § 3.2 above).
Faulty electrode body	Replace electrode body (see § 3.3 above).
Indication drifts:	
Temperature not constant	Use a thermostatic jacket, or wait until all solutions are at room temperature.
Sample preconditioning inadequate	Add ISA.
Membrane contaminated or damaged	Clean or replace membrane module (see § 3.2 above).
Slope incorrect:	
Calibration solution(s) contaminated or incorrectly diluted	Prepare new calibration solution(s) (see § 2.2 above).
Sample preconditioning inadequate	Add ISA.
Membrane module not screwed up tight enough	Tighten membrane module.
Membrane module damaged	Replace membrane module (see § 3.2 above).
Incorrect result despite correct calibration values:	
Defective calibration curve	Plot mV on the linear axis and check plotting direction along the logarithmic axis.
Incorrect sign	Note sign when reading off mV.
Incorrect calibration solution	Repeat dilution, and prepare new calibration solution if necessary.
Incorrect ISA addition	Add ISA to calibration and sample solutions in the same proportions.
Incorrect units	10 ⁻³ mol/l = 35.45 mg/l as Cl ⁻ 10 ⁻³ mol/l = 58.44 mg/l as NaCl.

5. Conversion tables

The following conversion factors apply:

Given value:	Multiply by:	To obtain:
mol/l Cl ⁻	35.45 · 10 ³	mg/l Cl ⁻
mol/l Cl ⁻	58.44 · 10 ³	mg/l NaCl
Given value:	Divide by:	To obtain:
mg/l Cl ⁻	35.45 · 10 ³	mol/l Cl ⁻
mg/l NaCl	58.44 · 10 ³	mol/l Cl ⁻

Relationships between concentration values:

mol/l Cl ⁻	mg/l Cl ⁻	mg/l NaCl
10 ⁻¹	3545	5844
10 ⁻²	354.5	584.4
10 ⁻³	35.45	58.44
10 ⁻⁴	3.545	5.844
10 ⁻⁵	0.3545	0.5844
10 ⁻⁶	0.03545	0.05844

6. Literature

K. Cammann, Working with Ion-Selective Electrodes. Springer-Verlag, 1979.
D. Midgley, K. Torrance, Potentiometric Water Analysis. John Wiley and Sons, 1978.

J. Vesely, D. Weiss, K. Stulik, Analysis with Ion-Selective Electrodes. Ellis Horwood Ltd./John Wiley and Sons, 1978.
Further literature for special applications can be supplied on request.

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Mettler-Toledo AG, Process