
General information

The DX287 fluoroborate electrode is suitable for measurements of fluoroborate concentrations and fluoroborate activities in aqueous and nonaqueous systems. General information on the use of ISEs can be found in the guide to ion selective measurement.

Specifications

Type of ion selective membrane: Polymer matrix-liquid membrane
Resistance of the membrane: -
Measurement range: >10⁻⁵ mol/L to approx. 1 mol/L fluoroborate
Measurement temperature range: 0 - 40 °C
Optimum pH range: 2.5 - 11
Interferences: Anionic surfactants (detergents) must be absent. The fluoroborate electrode is sensitive to the following ions. The ratio of these ions to fluoroborate (X/BF₄⁻) should thus not be larger than the value shown in brackets: ClO₄⁻ (<0.005), I⁻ (<0.05), ClO₃⁻ (<0.5), CN⁻ (<5), Br⁻ (<10), NO₂⁻ (<10), NO₃⁻ (<50), HCO₃⁻ (<333.3), CO₃²⁻ (<5555.6), Cl⁻ (<500), H₂PO₄⁻ (<9091), PO₄³⁻ (<500000), OAc⁻ (<2000), F⁻ (<5882.4), SO₄²⁻ (<31250).

Precision

In the absence of interfering substances, a precision of better than ±0.5 mV corresponding to ±2% of the measured BF₄⁻ concentration can be achieved.

Response times

When changing from lower to higher concentrations below 10⁻⁴ mol/L BF₄⁻, response time approx. half a minute, above 10⁻⁴ mol/L BF₄⁻ less than half a minute; when changing from higher to lower concentrations several minutes.

Materials

Electrode body of POM copolymer, membrane module of PVC, seals of silicon.

Bridge electrolyte for DX200 reference electrode
Inner electrolyte: 3 M KCl, saturated with AgCl.
Bridge electrolyte: 0.1 M NH₄NO₃.

Maintenance and service

Deposits on the surface of the membrane can be removed by placing the DX287 fluoroborate electrode in deionized water for a few minutes. After rinsing with deionized water, conditioning in dilute fluoroborate solution (approx. 0.01mol/L) during 5 to 10 minutes is required before the next calibration.

Storage

The fluoroborate electrode can be stored in air. Storage in a 10⁻³ molar solution of the appropriate ion maintains the electrode ready for use.
The electrode is best stored in the supplied tube during longer periods.
DX200 reference electrode: allow bridge electrolyte to run out. For storage, close filling port with the rubber cap. The service life of the reference electrode is shortened if the bridge electrolyte chamber constantly remains filled above the internal ceramic frit.

Changing the membrane module

Membrane modules with a damaged membrane (excessive response times or unstable potentials) which can no longer be regenerated must be changed.

Procedure for changing a membrane module:

- Place electrode upside down and unscrew old membrane module from the electrode body.
- Fill electrode body with internal solution up to the edge ensuring absence of bubbles.
- Also fill the new membrane module with internal solution up to the sealing ring.
- Screw the membrane module on tightly.
- Give the electrode a single powerful shake towards the electrode head (like a clinical thermometer).
- Rinse the entire electrode thoroughly with deionized water. Then condition during 10 minutes in dilute fluoroborate solution (approximately 0.01 mol/L).

Accessories

Reference electrode	DX200
Cable for ISE	ME-51089954
Membrane module replacement	ME-51107690

Literature

Guide to Ion Selective Measurement	ME-51300075
------------------------------------	-------------

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