**Instruction manual
for Ammonia Electrode (NH₃)**

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NH3 MSG SA 001

METTLER TOLEDO

A large, stylized graphic element at the bottom of the page, consisting of a series of parallel diagonal lines forming a triangular shape pointing upwards, with the company name 'METTLER TOLEDO' overlaid on it.

English **Page** **3**

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Ammonia Electrode (NH₃)

Instruction manual

Contents

	Page
1. General information	5
1.1 Technical data	
1.2 Scope of delivery	
1.3 Apparatus and accessories required	
1.4 Chemical solutions required	
2. Use of electrode	7
2.1 Preparation for use	
2.2 Calibration and measurement	
2.3 Particular requirements for determining low NH ₃ concentrations	
3. Care and maintenance of the NH ₃ electrode	10
3.1 Storage	
3.2 Replacing membrane module	
4. Trouble-shooting procedures	10
5. Conversion tables	12
6. Literature	12
7. Ordering information	12
8. Appendix	13

1. General information

The Type 15 230 3000 Ammonia Electrode is intended for measuring the ammonia content (NH_3) of aqueous systems. Its preferred applications are to be found in the fields of water and effluent, soil and microbiological analysis, as well as that of Kjeldahl digestions (e.g. proteins).

1.1 Technical data

Primary measured value:	potential in mV, proportional to $\log(\text{NH}_3)$.
Measuring range:	10^{-6} mol/l to $5 \cdot 10^{-2}$ mol/l NH_3
Operating temperature range:	0 to 50 °C
Sample medium:	aqueous only (surface-active agents will shorten the working life of the membrane module)
Maximum permissible ionic strength of the sample medium:	0.2 mol/l
Sample volume:	preferably greater than 20 ml (at concentration $\leq 10^{-5}$ mol/l, still larger volumes will be required)
Interfering substances:	volatile bases, such as amines
Electrode poisons:	surface-active agents (e.g. detergents, etc.)
Precision:	in the absence of interfering substances and over the linear portion of the characteristic curve it is possible to obtain a precision of better than $\pm 2\%$ of the measured NH_3 concentration, corresponding to ± 0.5 mV
Response time:	when passing from 10^{-5} mol/l NH_3 to higher concentrations, about 1 minute, and when passing from high concentrations down to 10^{-5} mol/l, several minutes (deviation from final constant value = 1 mV)
Electrode dimensions:	length from lower edge of screw cap: 120 mm diameter: 12 mm
Materials:	electrode body and membrane module body in chemically resistant plastic (PPS)
Working life:	at least 1 year for the pH electrode assembly. The life of the membrane module depends on the operating conditions, but is generally several months

Connections: AS7 cable with appliance coupler as specified by the manufacturer of the pH/mV- or ion-meter (see leaflet)

1.2 Scope of delivery

The electrode, or more precisely the electrode body, is supplied fitted with a hydrating cap, and thus is not immediately ready for use. Three membrane modules (of which two are spares), and a 50 ml bottle of electrolyte are supplied in a separate package.

1.3 Apparatus and accessories required

- Ion-meter or pH/mV-meter
- Magnetic stirrer with stirring bars
- Covered measuring vessels (titrating vessels) (pressure compensation)
- Suitable electrode/measuring vessel adapter
- Various pipettes, preferably microlitre pipettes
- Calibrated flasks

The following accessories may be useful, particularly when no ion-meter is available:

- Log/linear graph paper and/or a pocket calculator with scientific and engineering functions.

1.4 Chemical solutions required

- Deionized water for preparing solutions; when determining very low NH_3 concentrations (i.e. $\leq 10^{-5}$ mol/l), preferably double-distilled water (N.B.: water from an ion-exchanger may contain interfering amines)
- Electrolyte solution as supplied with the membrane kit (Ordering No. 15 230 2000)
- NH_3 sample conditioning solution

Prepare as follows:

Dissolve 40.0 g sodium hydroxide (NaOH) in about 80 ml deionized water (caution: wear safety goggles) and make up to 100 ml (resulting concentration: 10 mol/l).

- NH_3 electrode storage solution

Prepare as follows:

Dilute the sample conditioning solution described above 1:100 with deionized water.

- 0.1 mol/l standard NH_3 solution

Prepare as follows:

Weigh out 5.35 g AR grade ammonium chloride (NH_4Cl), dissolve in deionized water and make up to 1000 ml.

or:

- 1000 mg/l standard NH_3 solution

Prepare as follows:

Weigh out 3.14 g AR grade ammonium chloride (NH_4Cl), dissolve in deionized water and make up to 1000 ml.

2. Use of electrode

2.1 Preparation for use

Before preparing the NH_3 electrode for use, make up a 0.1 mol/l NaOH electrode storage solution as described in § 1.4.

Remove the hydrating cap and rinse the pH electrode assembly with deionized water.

Fill app. 15 drops or 0.8 to 1.0 ml respectively of electrolyte solution into the membrane module and remove any air bubbles adhering to the sides by tapping gently with the fingers. Now screw the membrane module to the electrode body and immerse the NH_3 electrode in the storage solution.

Neither the gas-permeable membrane of the NH_3 electrode nor the glass membrane of the pH electrode should be touched with the fingers.

Connect the electrode with the cable to the input socket of the ion-meter or pH/mV-meter.

Before calibration or measurement, immerse the NH_3 electrode for at least 30 minutes, before its very first use for several hours, in the storage solution to precondition it.

2.2 Calibration and measurement

Please note:

- Calibration and sample solutions should always be **at the same temperature** (the measuring vessel should be thermally insulated from possible heating by the magnetic stirrer – generally a sheet of cardboard or plastic foam material will suffice).
- All measurements should be carried out under the **same stirring conditions**, i.e. identical stirring speed, identical stirring bars, identical distance between electrodes and stirring bars, etc.
- When not in use, the NH_3 electrode should always be kept immersed in the storage solution (see also § 3).
- In all calibration and sample solutions the **ionic strength should be the same**, and should not exceed a maximum value of 0.2 mol/l. After addition of the sample preconditioning solution, the pH value should be higher than 11.0.
- The electrode should not be immersed to a depth greater than the length of the membrane module, i.e. not more than 6 cm.

General instruction for determining ammonia content

a) Calibration

- 1) Prepare a series of dilutions of the stock standard ammonia solution (0.1 mol/l or 1000 mg/l), preferably choosing concentrations which bracket that expected from the sample solution. When calibrating for the first time, at least three calibration points should be determined.
- 2) Run equal volumes of the calibration solutions, beginning with the lowest concentrations, into measuring vessels and stir.
- 3) Immerse the NH_3 electrode in a calibration solution and then add to this solution $1/100$ th of its volume of sample conditioning solution. Remove any gas bubbles adhering to the membrane by tapping with the fingers against the electrode body.
- 4) Once the mV indication has settled down, read off the final constant value.
- 5) Remove the electrode, rinse it briefly with deionized water and dab off with a paper tissue.
- 6) Repeat steps 3) to 5) above for all the other calibration solutions.
- 7) Plot the measured mV values against the logarithms of the NH_3 concentrations of the corresponding calibration solutions on log/linear graph paper to obtain a calibration curve:

NH_3 concentration	→ logarithmic axis
mV value	→ linear axis

b) Measurement

For measurement, proceed as for calibration. Read off the NH_3 content of the sample solution from the measured mV value with the aid of the calibration curve (see Figure 2), or calculate it with the pocket calculator.

The frequency of recalibration operations depends on the degree of accuracy required and on the measuring conditions.

2.3 Particular requirements for determining low NH_3 concentrations


The measurement of low NH_3 concentrations (i.e. less than 10^{-5} mol/l) is affected by the following factors:

- flatter NH_3 electrode slope
- sluggish response
- interfering substances or contamination

Use the largest possible volumes of both calibration and sample solutions. If the electrode has been previously used for determining high NH_3 concentrations, it is recommended to drain out and refill with internal electrolyte. It is advisable to use a standard addition method with several NH_3 addition steps, employing strip chart recorder. Start with a NH_3 -free "blank" solution (such as the storage solution), and renew this repeatedly until the measured electrode potential does not increase with a further renewal. Then add a definite quantity of the sample solution, followed by standard NH_3 solution additions, noting the electrode potential each time.

The measured value in the "blank" solution then enables the practical limit of detection to be determined, and hence the accuracy of the subsequent NH_3 determinations to be checked (see Figure 2).

Figure 2: Typical calibration curve for NH_3 electrode

3. Care and maintenance of the NH_3 electrode

3.1 Storage

The NH_3 electrode should not be kept immersed in the storage solution for more than a few days. For longer periods, the membrane module is removed and the pH electrode assembly and the membrane module thoroughly rinsed with deionized water. Fill the hydrating cap with deionized water and fit it to the pH electrode assembly, this, together with the membrane module, then being placed for storage in the original delivery package. Properly cleaned membrane modules may be kept stored dry for an indefinite period.

3.2 Replacing membrane module

Place app. 15 drops or 0.8 to 1.0 ml respectively of electrolyte solution in the new module and remove any air bubbles adhering to the sides by gently tapping with the fingers. Now screw the new module to the electrode (the latter previously rinsed with deionized water) and immerse the complete NH_3 electrode in the storage solution.

4. Trouble-shooting procedures

Fault:

Possible cause:	Remedy:
Value outside measuring range:	
Fault in instrument	Check with input short-circuited
Fault contact in electrode plug	Plug in and out several times
Insufficient electrolyte	Top up electrolyte as necessary
Fault in pH electrode assembly	Remove membrane module, measure potential in pH 7.0 buffer solution containing 0.1 mol/l NaCl (0.58 g NaCl to 100 ml buffer solution), and then in pH 4.0 buffer solution containing an equivalent amount of NaCl (measure small volumes). After settling down, the difference in potential should be 160 to 180 mV.
Unstable indication:	
Fault in instrument	Check with input short-circuited
No sample preconditioning	Add sample conditioning solution
Faulty earthing	Check earthing of instrument and magnetic stirrer
Static charges	Treat all plastic parts of the instrument with a suitable antistatic or wetting agent
Air bubble between plastic membrane and glass pH membrane	Loosen membrane module, tap it with the fingers and then screw it up again
Faulty electrode body	Replace electrode body

Possible cause:	Remedy:
Indication drifts:	
Temperature not constant	Use a thermostatic jacket, or wait until all solutions are at room temperature.
No sample preconditioning	Add sample conditioning solution
Membrane module not screwed up tight enough	Tighten membrane module.
Insufficient electrolyte	Top up or renew electrolyte.
Damage to gas-permeable membrane due to surface-active agents (detergents), ageing or mechanical action	Check for damaged areas, and replace membrane module if necessary.
Escape of gases	Check closure of measuring vessel.
Ionic strength of sample too high (>0.2 mol/l)	Dilute sample solution.
pH electrode assembly stored dry too long	Immerse pH electrode assembly in water.
Inadequate preconditioning	Immerse electrode for at least 30 minutes in storage solution before use.
Sample solution unstable or contaminated with interfering substances	Modify standard measuring procedure as appropriate and if this does not work, contact your METTLER TOLEDO agent for advice.
Measurement of low NH_3 concentrations after previous measurement of high NH_3 concentrations ("memory effect")	Immerse in storage solution before the next measurement.
Slope incorrect:	
Calibration solution(s) incorrectly prepared	Prepare new calibration solution(s).
No sample preconditioning	Add sample conditioning solution
Membrane module not screwed up tight enough	Tighten membrane module.
Membrane damaged	Replace membrane module.
Electrode not properly cleaned when passing from high to low concentration	Rinse electrode thoroughly, and immerse in storage solution if necessary.
pH electrode assembly damaged	Check in pH 7.0 and pH 4.0 buffer solutions as described under "Value outside measuring range" above.

5. Conversion tables

The following conversion factors apply:

Given value:	Multiply by:	To obtain:
mol/l NH_3	$17.0 \cdot 10^3$	mg/l NH_3
mol/l NH_3	$14.0 \cdot 10^3$	mg/l N
Given value:	Divide by:	To obtain:
mg/l NH_3	$17.0 \cdot 10^3$	mol/l NH_3

Relationship between concentration values:

mol/l NH_3	mg/l NH_3	mg/l N
10^{-1}	1700	1400
10^{-2}	170	140
10^{-3}	17	14
10^{-4}	1.7	1.4
10^{-5}	0.17	0.14

6. Literature

K. Cammann, H. Galster, "Working with Ion-Selective Electrodes". Springer-Verlag, 1996.

D. Migley, K. Torrance, "Potentiometric Water Analysis". John Wiley and Sons, 1978.

J. Vesely, D. Weiss, K. Stulik, "Analysis with Ion-Selective Electrodes".

Ellis Horwood Ltd./John Wiley and Sons, 1978.

Further literature for special applications can be supplied on request.

7. Ordering information


The following components of the Ammonia Electrode can be supplied separately:

- | | |
|--|---------------------------------|
| – Membrane kit (3 replacement membrane modules + 50 ml electrolyte solution) | Ordering No. 15 230 2000 |
| – 50 ml electrolyte solution | Ordering No. 15 230 1000 |
| – AS7 connecting cable consisting of: | on request |
| – AS7 cable socket | |
| – ST-Li-Koax 2.8 coax cable and plug to suit the instrument in use | |

8. Appendix

Structural arrangement of a potentiometric ammonia electrode.

AS7 Connecting cable

Structural arrangement of type 15 230 3000 NH₃ Electrode:

- 1) S7 screw cap
 - 2) plastic electrode body (PPS)
 - 3) Ag/AgCl reference electrode
 - 4) glass pH membrane
 - 5) hydrating cap (only fitted on delivery or for long-term storage of electrode)
 - 6) plastic module body (PPS)
 - 7) gas-permeable membrane (micro-porous)
- } pH electrode assembly

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- DE** **Mettler-Toledo GmbH, Prozeßanalytik,**
- USA** **Mettler-Toledo Ingold, Inc.,**
Mettler-Toledo GmbH, Process, Industrie Nord, 8902 Urdorf, Switzerland
Subject to technical changes. Mettler-Toledo GmbH, 11/04. 20 606 1096