

7. Ordering information

The following components of the NO_x electrode can be supplied separately:

- Membrane kit (3 replacement membrane modules + 50 ml electrolyte solution) Ordering No. **15 231 2000**
 - 50 ml electrolyte solution Ordering No. **15 231 1000**
 - Replacement electrode body Ordering No. **15 200 2000**
- AS7 connecting cable consisting of:
- AS7 cable socket
 - ST-Li-Koax 2.8 coax cable
 - and plug to suit the instrument in use on request

7. Informations pour commande de matériel

Les parties suivantes de l'électrode à NO_x peuvent être commandées séparément:

- (Membrane-Kit) (3 modules à membrane de remplacement et 50 ml d'électrolyte) No. de cde. **15 231 2000**
 - 50 ml d'électrolyte No. de cde. **15 231 1000**
 - Corps de remplacement de l'électrode No. de cde. **15 200 2000**
- Cable de raccordement AS7 comportant:
- une prise AS7 pour le câble
 - un câble coax ST-Li-Koax 2.8
 - une fiche correspondant à votre appareil de mesure sur demande

7. Bestellinformationen

Folgende Teile der NO_x-Elektrode können separat bestellt werden:

- Membran-Kit (3 Membranmodule und 50 ml Elektrolyt) Bestell-Nr. **15 231 2000**
 - 50 ml Elektrolyt Bestell-Nr. **15 231 1000**
 - Elektrodenkörper Bestell-Nr. **15 200 2000**
- AS7-Verbindungskabel bestehend aus:
- Kabelbuchse AS7
 - Koax-Kabel ST-Li-Koax 2.8
 - und Stecker passend zum Messgerät auf Anfrage

**Operating Instruction
for NO_x Electrode**

**Manuel d'Instruction
pour l'Electrode à NO_x**

**Anleitung
zur NO_x Elektrode**

Mettler-Toledo AG, Process
Im Hackacker 15
8902 Urdorf/Switzerland
Phone: ++1 736 22 11
Fax: ++1 736 26 36

20 606 1009
231-O2-CH

© Mettler-Toledo AG, 8606 Greifensee, Switzerland

METTLER TOLEO

Operating Instructions for Type 15 231 3000 NO_x Electrode

Contents

Pages 3-7

- General information
 - 1.1 Technical data
 - 1.2 Scope of delivery
 - 1.3 Apparatus and accessories required
 - 1.4 Chemical solutions required
- Use of electrode
 - 2.1 Preparation for use
 - 2.2 Calibration and measurement**
 - 2.3 Particular requirements for determining low NO_x concentrations in the presence of CO₂
- Care and maintenance of the NO_x electrode
 - 3.1 Storage
 - 3.2 Replacing membrane module
- Troubleshooting procedures
- Conversion tables
- Literature
- Ordering information

Page 20

Instructions pour l'Electrode à NO_x Type 15 231 3000

Sommaire

Pages 3 et 8-12

- Informations générales
 - 1.1 Spécifications
 - 1.2 Etendue de la livraison
 - 1.3 Appareillage et accessoires
 - 1.4 Solutions nécessaires
- Emploi de l'électrode
 - 2.1 Mise en service
 - 2.2 Etalonnage et mesurage**
 - 2.3 Particularités lors de la détermination de faibles teneurs en NO_x en présence de CO₂
- Maintenance de l'électrode
 - 3.1 Conservation
 - 3.2 Remplacement du module à membrane
- Dépannage
- Tables de conversion
- Bibliographie
- Informations pour commande de matériel

Page 20

Einleitung zur NO_x-Elektrode Typ 15 231 3000

Inhalt

Seiten 3 und 13-16

- Allgemeine Informationen
 - 1.1 Spezifikationen
 - 1.2 Lieferumfang
 - 1.3 Geräte und Hilfsmittel
 - 1.4 Erforderliche Lösungen
- Gebrauch der Elektrode
 - 2.1 Inbetriebnahme
 - 2.2 Eichen und Messen**
 - 2.3 Besonderheiten bei der Bestimmung niedriger NO_x-Gehalte in Gegenwart von CO₂
- Wartung und Unterhalt der Elektrode
 - 3.1 Lagerung
 - 3.2 Ersetzen des Membranmoduls
- Fehlersuche
- Umrechnungstabellen
- Literatur
- Bestellinformationen

Seite 20

Figure 1: Structural arrangement of a potentiometric NO_x electrode.


Figure 1: Construction de l'électrode potentiométrique à NO_x.

Abbildung 1: Aufbau einer potentiometrischen NO_x-Elektrode.

AS7 Connecting cable
AS7 Câble de raccordement
AS7 Verbindungskabel

Electrode body
Corps de l'électrode
Elektrodenkörper

Membrane module
Module à membrane
Membranmodul



Structural arrangement of type 15 231 3000 NO_x Electrode

- 1 S7 screw cap
- 2 plastic electrode body (POM copolymer)
- 3 Ag/AgCl reference electrode
- 4 glass pH membrane
- 5 hydrating cap (only fitted on delivery or for long-term storage of electrode)
- 6 plastic module body (POM copolymer)
- 7 gas-permeable membrane (micro-porous)

Construction de l'électrode à NO_x type 15 231 3000

- 1 Tête enfichable S7
- 2 Corps en matière synthétique (copolymère POM)
- 3 Electrode de référence à Ag/AgCl
- 4 Membrane de verre sensible au pH
- 5 Capuchon de mouillage (seulement pour l'envoi et pour la conservation à long terme de l'électrode)
- 6 Corps modulaire en matière synthétique (copolymère POM)
- 7 Membrane perméable aux gaz (microporeuse)

Aufbau der NO_x-Elektrode Typ 15 231 3000

- 1 S7-Steckkopf
- 2 Kunststoffkörper (POM Copolymer)
- 3 Ag/AgCl-Referenzelektrode
- 4 pH-Glasmembran
- 5 Wässerungskappe (nur beim Versand und längerem Lagern der Elektrode)
- 6 Kunststoff-Modulkörper (POM Copolymer)
- 7 Gaspermeable Membran (mikroporös)

1. General information

English

The Type 15 2313000 NO_x Electrode is intended for measuring the NO_x content of aqueous systems. Its preferred applications are to be found in the fields of water and waste water, soil, food analysis etc.

1.1 Technical data

Primary measured value:	potential in mV, proportional to log [NO _x]
Measuring range:	≤10 ⁻⁶ mol/l to 5 · 10 ⁻³ mol/l NO _x
Operating temperature range:	0 to 50 °C
Sample medium:	aqueous only (surface-active agents will shorten the working life of the membrane module).
Maximum permissible ionic strength of the sample medium:	0.8 mol/l
Sample volume:	preferably greater than 20 ml (at concentrations ≤10 ⁻⁵ mol/l, even larger volumes will be required).
Interfering substances:	SO ₂ , CO ₂ , H ₂ S, volatile carbonic acids.
Electrode poisons:	Surface-active agents (e. g. detergents), H ₂ S, etc.
Precision:	in the absence of interfering substances and over the linear portion of the characteristic curve it is possible to obtain a precision of better than ±2% of the measured NO _x concentration, corresponding to ±0.5 mV.
Response time:	when passing from 10 ⁻⁵ mol/l NO _x to higher concentrations, about 1 minute, and when passing from high concentrations down to 10 ⁻⁵ mol/l, several minutes (deviation from final constant value = 1 mV).
Electrode dimensions:	length from lower edge of screw cap: 120 mm diameter: 12 mm
Materials:	electrode body and membrane module body in chemically resistant plastic (POM copolymer).
Working life:	at least 1 year for the pH electrode assembly. The life of the membrane module depends on the operating conditions, but is generally a few months.
Connections:	AS7 cable with appliance coupler as specified by the manufacturer of the pH/mV- or ion-meter (see leaflet).

1.2 Scope of delivery

The electrode, or more precisely the electrode body, is supplied fitted with a hydrating cap, and thus is not immediately ready for use. Three membrane

modules (of which two are spares), and a 50 ml bottle of electrolyte are supplied in a separate package.

1.3 Apparatus and accessories required

- Ion-meter or pH/mV-meter
- Magnetic stirrer with stirring bars
- Closed measuring vessels (titrating vessels) (pressure equalization with the atmosphere)
- Suitable electrode/vessel adapter
- Various pipettes, preferably microlitre pipettes
- Calibrated flasks
- The following accessories may be useful, particularly when no ion-meter is available:
- Log/linear graph paper and/or a pocket calculator with scientific and engineering functions.

1.4 Chemical solutions required

- Deionized water for preparing solutions
- Electrolyte solution as supplied with the membrane kit (Ordering No. 15 2312000)
- NO_x electrode storage solution
Prepare as follows:
Dissolve 190 g sodium sulphate (Na₂SO₄) in about 800 ml deionized water and make up to 1000 ml. Take 100 ml of this stock solution and make up to 1000 ml with deionized water to obtain the storage solution.
- NO_x sample conditioning solution
Prepare as follows:
Take the remaining 900 ml of the so-

dium sulphate solution prepared as described above and carefully add 48 ml concentrated sulphuric acid (96 to 97% H₂SO₄) (caution: wear safety goggles).

- 0.1 mol/l standard NO_x solution
Prepare as follows:
Weigh out 6.9 g AR grade sodium nitrite (NaNO₂), dissolve in deionized water and make up to 1000 ml.
- 1000 mg/l standard NO_x solution
Prepare as follows:
Weigh out 4.93 g AR grade sodium nitrite (NaNO₂), dissolve in deionized water and make up to 1000 ml.

2. Use of electrode

2.1 Preparation for use

Before preparing the NO_x electrode for use, make up an electrode storage solution as described in § 1.4.
Remove the hydrating cap and rinse the pH electrode assembly with deionized water.
Fill approx. 15 drops or 0.8 to 1.0 ml of electrolyte solution into the membrane module and remove any air bubbles adhering to the sides by tapping gently with the fingers. Now screw the membrane module to the electrode body and immerse the NO_x electrode in the storage solution.


Neither the gas-permeable membrane of the NO_x electrode nor the glass membrane of the pH electrode should be touched with the fingers.
Connect the electrode with the cable to the input socket of the ion- or pH/mV-meter.
Prior to carrying out the first measurement with a new NO_x electrode immerse it for several hours in the storage solution to precondition it. Afterwards, before calibration or measurement, immerse the NO_x electrode at least for 30 minutes in the storage solution.

2.2 Calibration and measurement

Please note!

- Calibration and sample solutions should always be at the same temperature (the measuring vessel should be thermally insulated from possible heating by the magnetic stirrer – generally a sheet of cardboard or plastic foam material will suffice).
- All measurements should be carried out under the same stirring conditions, i. e. identical stirring speed, identical stirring bars, identical distance between electrodes and stirring bars, etc.
- When not in use, the NO_x electrode should always be kept immersed in the storage solution (see also § 3).
- In all calibration and sample solutions the ionic strength should be the same and should not exceed a maximum value of 0.8 mol/l. After addition of the sample preconditioning solution, the pH value should be less than 1.8 (see also § 2.5).
- The electrode should not be immersed to a depth greater than the length of the membrane module, i. e. not more than 6 cm.
- The calibration and measuring vessels should be well closed but allowing a pressure equalization with the atmosphere. They shall have a small volume ratio of the gas phase over the liquid phase.

Figure 2:
Typical calibration curve for an NO_x electrode



General instructions for determining NO_x content

a) Calibration

- 1) Prepare a series of dilutions of the stock standard NO_x solution (0.1 mol/l or 1000 mg/l), preferably choosing concentrations which bracket that expected from the sample solution. When calibrating for the first time, at least three calibration points should be determined.
- 2) Run equal volumes of the calibration solutions, beginning with the lowest concentrations, into measuring vessels and stir.
- 3) Immerse the NO_x electrode in a calibration solution and then add to this solution 1/10th of its volume of sample conditioning solution and immediately stopper the measuring vessel. Gas bubbles adhering to the membrane

are removed by tapping the electrode body with the fingers.

- 4) Once the mV indication has settled down, read off the final constant value.
- 5) Remove the electrode, rinse it briefly with deionized water and dab off with a paper tissue.
- 6) Repeat steps 3) to 5) above for all the other calibration solutions.
- 7) Plot the measured mV values against the logarithms of the NO_x concentrations of the corresponding calibration solutions on log/linear graph paper to obtain a calibration curve:

NO_x concentration → logarithmic axis
mV value → linear axis

b) Measurement

For measurement, proceed as for calibration. Read off the NO_x content of the sample solution from the measured mV value with the aid of the calibration curve (see Figure 2), or calculate it with the pocket calculator.

The frequency of recalibration operations depends on the degree of accuracy required and on the measuring conditions.

2.3 Particular requirements for determining low NO_x concentrations in the presence of CO_2

The measurement of low NO_x concentrations (i. e. less than 10^{-5} mol/l) is affected by the following factors:

- flatter NO_x electrode slope
- sluggish response
- interfering substances, particularly CO_2

Use the largest possible volumes of both calibration and sample solutions. If the electrode has been previously used for determining high NO_x concentrations, it should be left immersed in the storage solution until stable values are obtained. If necessary, the electrolyte should be renewed.

It is advisable to use an addition or enrichment method with several NO_x addition steps when measured in the bent region of the calibration curve, employing a strip chart recorder. Start with an NO_x -free (blank solution) (such as the storage solution), and renew this repeatedly until


the measured electrode potential does not decrease with a further renewal. Then add a definite quantity of the sample solution, followed by standard NO_x solution additions, noting the electrode potential each time.

The measured value in the (blank solution) then enables the practical limit of detection to be determined, and hence the precision of the subsequent NO_x determinations to be checked (see figure 2).

The interfering substance most frequently encountered is CO_2 . This is liable to be present in considerable concentrations in samples such as effluent water, tap water, meat, etc. If the display indicates an overshoot when sample conditioning solution is added, followed by a slow fall-off, this is generally due to CO_2 interference. The response behaviour can be seen from figure 3.

Figure 3

Typical response behaviour of the NO_x electrode in the presence of CO_2 or HCO_3^- . Addition of sample conditioning solution to a 10^{-5} mol/l nitrite solution. For higher nitrite concentrations the CO_2 interference is smaller.



3. Storage and maintenance of the NO_x electrode

3.1 Storage

The NO_x electrode should not be kept immersed in the storage solution for more than a few days. For longer periods, remove the membrane module. Thoroughly rinse the pH electrode assembly and the membrane module with deionized water. Fill the hydrating cap with deionized water and fit it to the pH

electrode assembly, this, together with the membrane module, then being placed for storage in the original delivery package. Properly cleaned membrane modules may be kept stored dry for an indefinite period.

3.2 Replacing membrane module

Place app. 15 drops or 0.8 to 1.0 ml of electrolyte solution in the new module and remove any air bubbles adhering to the sides by gently tapping with the fingers. Now screw the new module to the

electrode (the latter previously rinsed with deionized water) and immerse the complete NO_x electrode in the storage solution.

4. Troubleshooting procedures

Fault:

Possible cause:	Remedy:
Value outside measuring range:	
Fault in instrument	Check with input short-circuited.
Faulty contact in electrode plug	Insert and withdraw repeatedly.
Insufficient electrolyte	Top up electrolyte as necessary.
Fault in pH electrode assembly	Remove membrane module, measure potential in pH 7.0 buffer solution containing 0.1 mol/l NaCl (0.58 g NaCl to 100 ml buffer solution), and then in pH 4.0 buffer solution containing an equivalent amount of NaCl (use small volumes). After settling down, the difference in potential should be 160 to 180 mV.
Unstable indication:	
Fault in instrument	Check with input short-circuited.
No sample preconditioning	Add sample conditioning solution.
Faulty earthing	Check earthing of instrument and magnetic stirrer.
Static charges	Treat all plastic parts of the instrument with a suitable antistatic or wetting agent.
Air bubble between plastic NO_x membrane and glass pH membrane	Loosen membrane module, tap it with the fingers and then screw it up again.
Faulty electrode body	Replace electrode body.
Indication drifts:	
Temperature not constant	Use a thermostatic jacket, or wait until all solutions are at room temperature.
No sample preconditioning	Add sample conditioning solution.
Membrane module not screwed up tight enough	Tighten membrane module.
Insufficient electrolyte	Top up or renew electrolyte.
Damage to gas-permeable membrane due to surface-active agents (detergents), ageing or mechanical action	Check for damaged areas, and replace membrane module if necessary.
Escape of gases	Check closure of measuring vessel.
Ionic strength of sample too high (>0.8 mol/l)	Dilute sample solution.
pH electrode assembly stored dry too long	Immerse pH electrode assembly in water.
Inadequate preconditioning	Immerse electrode for at least 30 minutes in storage solution before use.
Sample solution unstable or contaminated with interfering substances	Modify standard measuring procedure as appropriate and if this does not work, contact your INGOLD Agent for advice.
Measurement of low NO_x concentrations after previous measurement of high NO_x concentrations (memory effect)	Immerse in storage solution before the next measurement.
Slope incorrect:	
Calibration solution(s) incorrectly prepared	Prepare new calibration solution(s).
No sample preconditioning	Add sample conditioning solution.
Membrane module not screwed up tight enough	Tighten membrane module.
Membrane damaged	Replace membrane module.
Electrode not properly cleaned when passing from high to low concentration	Rinse electrode thoroughly, and immerse in storage solution if necessary.
pH electrode assembly damaged	Check in pH 7 and pH 4 buffer solutions as described under (Value outside measuring range) above.

5. Conversion tables

The following conversion factors apply:

Given value:	Multiply by:	To obtain:
mol/l NO_2^-	$46.0 \cdot 10^3$	mg/l NO_2^-
mol/l NO_2^-	$14.0 \cdot 10^3$	mg/l N- NO_2^-
Given value:	Divide by:	To obtain:
mg/l NO_2^-	$46.0 \cdot 10^3$	mol/l NO_2^-

Relationship between concentration values:

mol/l NO_2^-	mg/l NO_2^-	mg/l N- NO_2^-
10^{-1}	4600	1400
10^{-2}	460	140
10^{-3}	46	14
10^{-4}	4.6	1.4
10^{-5}	0.46	0.14

6. Literature

K. Cammann, (Working with Ion-Selective Electrodes), Springer-Verlag, 1979.
D. Midgley, K. Torrance, (Potentiometric Water Analysis), John Wiley and Sons, 1978.

J. Vesely, D. Weiss, K. Stulik, (Analysis with Ion-Selective Electrodes), Ellis Horwood Ltd./John Wiley and Sons, 1978.
Further literature for special applications can be supplied on request.