
English
Français
Español
Polski

User Manual **SevenExcellence™**
Guide de l'utilisateur **SevenExcellence™**
Manual de usuario **SevenExcellence™**
Podręcznik użytkownika **SevenExcellence™**

METTLER TOLEDO

User Manual **SevenExcellence™**

English

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Polski

Table of Contents

1	Introduction	3
2	Safety information	4
2.1	Definitions of signal words and warning symbols	4
2.2	Product specific safety notes	4
3	Design and Function	6
3.1	Description of the instrument	6
3.2	Overview	7
3.3	Module connections	9
3.4	User interface	12
3.4.1	Homescreen	12
4	Putting into Operation	16
4.1	Scope of delivery	16
4.2	Connecting and removing modules	16
4.3	Mounting uPlace™ electrode arm	18
4.4	Installing power supply	19
4.5	Switching the instrument on and off	20
4.6	Connecting sensors	20
4.7	Connecting autosampler	20
4.8	Connecting stirrers	20
4.9	Connecting PnP devices	21
4.10	Connecting computers	21
5	Operating the Instrument	22
5.1	Direct calibration	22
5.2	Direct measurement	22
5.3	Starting a method	22
6	Maintenance and Care	23
6.1	Cleaning the Instrument	23
6.2	Maintenance of electrodes	23
6.3	Transporting the instrument	24
6.4	Disposal	24
7	Accessories	25
8	Technical Data	28

1 Introduction

The METTLER TOLEDO SevenExcellence™ is a modern, professional instrument suitable for use in a wide variety of application areas and segments. It can be used, for example, in quality control as well as in research and development and meets the most demanding requirements.

This flexible meter perfectly combines easy-to-understand operation with a high level of measurement accuracy and outstanding reliability. Thanks to the plug & play capability, the detection of external devices and ISM® sensors occurs automatically, in this way supporting the comfortable handling of SevenExcellence™.

Straightforward user guidance on the large color touch screen makes the operation particularly intuitive. User-definable shortcuts allow methods to be activated directly from the home screen with a single tap.

About this document

This document provides you with the information you need to get started with your METTLER TOLEDO instrument. The instructions in this document refer to running firmware version **4.3** or higher.

For a full description of the instrument and its functions, refer to the Reference Manual, supplied as PDF file online.

► www.mt.com/library

► <http://mt.com/pHLabsupport>

If you have any additional questions, contact your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

Conventions and symbols

Note for useful information about the product.

Elements of instructions

- Prerequisites
- 1 Steps
- 2 ...
 - ⇒ Intermediate results
- ⇒ Results

2 Safety information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

User Manual and Reference Manual are available online.

► www.mt.com/library

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

WARNING A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

Electrical shock

2.2 Product specific safety notes

Intended use

This instrument is designed to be used by trained staff. The SevenExcellence™ is intended for measuring pH, mV, ion, conductivity, DO/BOD.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

METTLER TOLEDO assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. METTLER TOLEDO assumes that the instrument owner provides the necessary protective gear.

Safety notes

WARNING

Danger of death or serious injury due to electric shock!

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.

NOTICE

Risk of damage to the instrument due to the use of unsuitable parts!

Using unsuitable parts with the instrument can damage the instrument or cause it to malfunction.

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

FCC Rules

This device complies with Part 15 of the FCC Rules and Radio Interference Requirements of the Canadian Department of Communications. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

3 Design and Function

3.1 Description of the instrument

The SevenExcellence™ bench meters are based on precise electrochemical measuring technique. Modular expansion units, automatically recognized ISM® sensors and the plug & play function for external devices make SevenExcellence a comfortable instrument for reliable measurement.

For the optimal use of your SevenExcellence™ bench meters, the following devices can be connected:

- **Modules**


- METTLER TOLEDO pH/mV module
- METTLER TOLEDO pH/Ion module
- METTLER TOLEDO conductivity module
- METTLER TOLEDO DO/BOD module

- **Sensors**

- Select from the wide range of METTLER TOLEDO electrodes, the electrode value box supports you to find the right tool. www.electrodes.net.
- ISM® sensors are detected automatically.
- Temperature probes can be connected separately.
- Rondolino sample changer allows fully automatic analysis of up to 9 samples. The optional PowerShower™ perfectly rinses contaminated electrodes.
- InMotion sample changer allows fully automatic analysis of up to 303 samples in different beaker sizes (only if operated by **LabX**).
- Computer for operation with **LabX** or **EasyDirect pH** software.
- uMix™ magnetic stirrer to optimize workflows and to enhance the reproducibility of results.
- Barcode reader for scanning sample data.
- Logstraight™ fingerprint reader for user identification.
- Compact printer (USB-P25), to print out results.
- USB stick to store and transfer data.
- LAN for printing with network printer.
- Additional compact printers, including RS-P25, RS-P26, RS-P28.


3.2 Overview

Front view

1	On/Off button	2	USB A interface
3	USB A interface	4	Touch screen
5	Info Opens a screen with information about the instrument and the attached modules.	6	Home Brings you back to the homescreen.


Rear view

1	LAN	Network connection
2	USB	USB B interface for PC connection (LabX or EasyDirect pH software)
3	Data	RS232 interface
4	Stir	Socket for METTLER TOLEDO stirrer (Mini-DIN)
5	DC	Power supply socket

Pin Assignments


Below the PIN assignments for the RS-232 interface are shown. METTLER TOLEDO printers such as RS-P25 can be connected to this interface.

Pin 1	NC	Pin 6	NC
Pin 2	TxD (out)	Pin 7	NC
Pin 3	RxD (in)	Pin 8	NC
Pin 4	NC	Pin 9	NC
Pin 5	RSGND		

3.3 Module connections


Instrument connection

In the lower area you find the socket for connecting the module to the instrument. A guide rail on the backside of the housing supports attaching the module.

Connections on the pH/mV module

The type of module is identified by a colored tag. Yellow for pH/mV.

Digi

Mini-LTW socket for digital sensors

ISFET

Mini-DIN socket for ISFET sensors

ATC int


RCA (Cinch) socket for internal temperature input (NTC30k)

pH

BNC socket for mV/pH sensors

Connections on the conductivity module


The type of module is identified by a colored tag. Purple for conductivity.

- Digi** Mini-LTW socket for digital sensors
- ATC ext** RCA (Cinch) socket for external temperature input (NTC30k or PT1000)
- Cond** Mini-DIN socket for conductivity sensors

Connections on the pH/Ion module


The type of module is identified by a colored tag. Orange for pH/Ion.

- Digi** Mini-LTW socket for digital sensors
- ATC ext** RCA (Cinch) socket for external temperature input (NTC30k or PT1000)
- ref** 2 mm banana socket for reference electrode
- ATC int** RCA (Cinch) socket for internal temperature input (NTC30k)
- pH** BNC socket for mV/pH sensors

Connections on the DO/BOD module

The type of module is identified by a colored tag. Blue for DO/BOD.

Digi

Mini-LTW socket for digital sensors

ATC ext

RCA (Cinch) socket for external temperature input (NTC22k)

ATC int

RCA (Cinch) socket for internal temperature input (NTC22k)

DO

BNC socket for the DO sensorss


3.4 User interface

NOTICE

Danger of damaging the touch screen with pointed or sharp objects!
Pressing on the touch screen with pointed or sharp objects may be damaged it.
– Operate the touch screen by applying gentle pressure with the pad of your finger.

3.4.1 Homescreen

Name	Explanation
1 Shortcuts	User-specific shortcuts for frequently used methods. Shortcuts are saved in the user profile and can be defined, changed and deleted by the user.
2 Modules	Attached modules and connected electrodes are shown.
3 Status bar	The status bar contains the current menu item, user name as well as date and time.
4 News	Indicates important changes of the system like expiration of resources etc. or disconnected PnP devices.
5 Tasks	Shows the current working status of the instrument. The list of tasks can be accessed.
6 Read	Initiate a direct measurement from the Home Screen, without the need to first set up a measurement method.
7 User data	Access an overview of the currently logged in user, the user group, and the time since the current user has been logged in.
8 Log out	Allows you to log out of the instrument, change user levels, or lock the instrument with a personalized password.

Name	Explanation
9 Menus	Methods Create and handle methods for every measurement type.
	Series Create and manage series of individual samples
	Results Display all measurement results, print out or export them. Find detail information about every single result.
	Setup Define all system settings in this menu, e.g. hardware settings, user management or user preferences. These settings are usually made during installation of the instrument.
	Sensors Modify and manage sensors.

Footer buttons

Depending on the selected submenu there are specific buttons in the footer.

AddToHome	Create a shortcut of a quick analysis.
Back	Go back one step in the menu structure.
Cancel	Cancel the current input without saving.
Delete method	Delete the selected method. Shortcuts and series which reference to this method will also be deleted.
Expired resources	Shows the expired resources connected to this instrument.
Insert	Insert a method function in an existing method.
Log out	Log out the current user.
New	Creates a new method, sample series, sensor etc.
OK	Confirm the entered settings.
Preview	When entering a text with an alphanumeric keypad, you can preview your input.
Print	Print out the currently shown screen in tabular form. A printer must be connected to the instrument.
Samples	Open the list with your defined samples.
Read	Start a quick analysis directly from the homescreen.
Start	Start a selected method.
User data	Access an overview of the currently logged in user

Input field types

In the single menus you have different kinds of parameter-fields to enter information, values or names or to choose an option in a drop-down list. Depending on the input field (value or name) you will get a numeric or alphabetic keypad for entering your values.

Text input field

Any text comprised of letters (up to 30 characters), numbers and symbols can be entered into these fields.

Text input field (extended)

Any text comprised of letters (up to 500 characters), numbers and symbols can be entered into these fields.

Number input field

Numbers and formulas can be entered into these fields.

Drop-down list

A drop-down list opens from which you can select an entry.

Shortlist

A shortlist opens from which you can select an entry.

Formula field

A formula must be entered in these fields.


Info field

The displayed content is only for information (read only).

Keypads

Alphabetic keypad

- Tap (1) to preview the entry.
- Tap (2) to enter capital letters.
- Tap (3) to enter lowercase letters.
- Tap (4) to switch to a numeric keypad and (2) to turn back to alphanumeric.
- Tap (5) to delete all entered letters or numbers.
- Tap (6) to delete the last entered letter or number.

Numeric keypad

- Tap (1) to delete all entered numbers.
- Tap (2) to delete the last entered number.

If the user language is set to Chinese, an additional button **CN** appears in the alphabetic keypad at the end of the input field. Tap it to enter Chinese character in Pinyin and select the character from the list below the input field. Be aware that Chinese characters cannot be printed on compact printers.

4 Putting into Operation

4.1 Scope of delivery

Unpack the instrument and check the scope of delivery. Keep the calibration certificate in a safe place.

SevenExcellence™ is delivered with:

- uPlace™ electrode arm
- Measurement modules and blank modules (according to your order)
- Sensors (kit version or according to your order)
- Universal AC adapter
- CD-ROM with Operating Instructions (all available languages)
- User Manual
- Declaration of conformity
- Calibration certificate

4.2 Connecting and removing modules

Note

- Before connecting or removing modules switch off the instrument.


METTLER TOLEDO offers different types of modules. The type of module is identified by a colored tag.

The instrument can be equipped with up to three modules. The modules can be combined in any composition and order according to your requirements.


With less than three modules installed, protect the open sockets against moisture or contamination by installing a blank module.

Connecting modules

- 1 Insert the expansion unit into the slot. The guide rails simplify the connecting.


- 2 Attach the cover.


Removing modules

- 1 To remove the cover, press the clips on both sides and lift it up.


- 2 To remove a module, pull it out of the slot.


4.3 Mounting uPlace™ electrode arm

The electrode arm can be used as stand alone or it can be attached to the instrument on the left or right side, according to your preferences. The height of the electrode arm can be varied by using the extension shaft part. Use the wrench to attach the extension part.



Assembly of the electrode arm


- 1 Use the wrench to attach the base to the electrode arm by tightening the screws. The electrode arm can now be used in the stand alone mode.

- 2 Then insert the foot of the meter to the arm base and shift the meter in the direction of the arrow to make the foot fit.

- 3 Use the lock screw to attach the meter to the base of the arm.

4.4 Installing power supply

WARNING

Danger of death or serious injury due to electric shock!

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.

NOTICE


Danger of damage to the AC adapter due to overheating!

If the AC adapter is covered or in a container, it is not sufficiently cooled and overheats.

- 1 Do not cover the AC adapter.
- 2 Do not put the AC adapter in a container.

The instrument is operated using an AC adapter. The AC adapter is suitable for all supply line voltages ranging from 100...240 V AC $\pm 10\%$ and 50-60 Hz.

- 1 Insert the correct connector plug into the AC adapter until it is completely inserted.
- 2 Connect the cable of the AC adapter with the DC socket of the instrument.
- 3 Install the cables in such a way that they cannot be damaged or interfere with operation.

- 4 Insert the plug of the AC adapter in a power outlet that is easily accessible.
- ⇒ To remove the connector plug, push the release button and withdraw the connector plug.

4.5 Switching the instrument on and off

Switching the instrument on

- 1 Press the **On/Off** button.
 - ⇒ The power button lights up green.
- 2 The instrument is initialized and a system test is performed. This process will take some time.
- ⇒ If only one user is defined, the instrument starts directly with the home screen without showing the login screen.

Switching the instrument off

- 1 Tap **Home > Log out > Shut down**.
 - or -
 - Press the **On/Off** button.
- 2 The instrument stops running tasks and shuts down. This process will take some time.
- ⇒ When the light goes out, the instrument is in the standby mode.

Note

- In the standby mode, the control circuit for the **On/Off** switch is energized. The rest of the instrument is no longer energized.

4.6 Connecting sensors

Connect the sensor and make sure that the plugs are properly inserted. For sensors with a built-in temperature probe, connect the second cable to the appropriate ATC socket (**ATC int**). When using separate temperature sensors, connect the cable to the appropriate ATC socket (**ATC ext**). Twist the RCA (Cinch) plug to ease the attachment of the sensor.

ISM® sensor

ISM® sensors are automatically recognized upon connection.

When connecting an ISM® sensor to the instrument, verify that no analysis is running. In case an analysis is running, the sensor identification will be displayed after the endpoint of the analysis is reached.

4.7 Connecting autosampler

InMotion autosampler

With the InMotion autosampler the automatic measurement of up to 303 samples is possible. InMotion is only supported when the instrument is connected to the **LabX** software. The InMotion autosampler is connected to the instrument with an USB cable and can only be used in **LabX** connected mode. Several different configurations are available for the InMotion autosampler depending on your requirements.

Rondolino sample changer

With the Rondolino sample changer the automatic measurement of up to 9 samples is possible. To connect the Rondolino sample changer to the instrument, a "Rondolino USB-TTL box" is required. Use the USB cable to connect the USB A interface of the instrument to the USB-TTL box. Use the 9-pin cable to connect the Rondolino sample changer with the TTL box. Connect the compact stirrer to the stirrer adapter cable and then to the instrument. The compact stirrer is an overhead stirrer.

4.8 Connecting stirrers

Connect uMix™, the METTLER TOLEDO external magnetic stirrer, or the METTLER TOLEDO Compact Stirrer (overhead stirrer, adapter cable needed) to the instrument. The stirrers are powered by the instrument and will be automatically switched on/off during an analysis according to the settings. Before connecting a stirrer, switch off the instrument. Use the stirrer socket (**Stir**) on the rear side, note the designation. Note the corresponding chapter of the operating instructions concerning the settings of peripheral devices.

4.9 Connecting PnP devices

USB sticks, barcode reader, US-P25 printer and fingerprint reader can be connected to the instrument via the USB connections on the left side of the instrument for immediate use. The instrument has a plug & play detection to support the comfortable handling of peripheral devices.

4.10 Connecting computers

The instrument can be connected to a PC with **LabX** or **EasyDirect pH** software.

Connected to the PC software **LabX** via USB or Ethernet, SevenExcellence profits from additional functionalities such as data and resource management, report creation and ERP/LIMS integration. The instrument can be controlled via touchscreen and PC software. In addition, **LabX** offers regulation support (e.g. 21 CFR part 11) and allows combined analysis with other METTLER TOLEDO laboratory instruments.

Find out more and visit us on:

 www.mt.com/LabX

The PC software **EasyDirect pH** allows for easy and efficient transfer of results to a PC. SevenExcellence is connected via USB, one instrument per PC.

Furthermore, the previous PC software **LabX direct pH** can be used to transfer results to Microsoft Excel or store them in a text file.

For simplified data transmission to a PC or network via RS232, the RS232 Compact Printer can be selected in the printer settings.

5 Operating the Instrument

The SevenExcellence™ display is equipped with a protective foil. Removing this foil increases the sensitivity of the touch screen.

SevenExcellence is delivered with default parameters for calibration and measurement as well as for methods to ensure an immediate deployment of the instrument after installation.

This quick guide leads you through the first steps of using SevenExcellence.

For a detailed guidance in using the extensive features of SevenExcellence the user is referred to the operating instructions.

5.1 Direct calibration

- 1 Select a suitable sensor for the calibration and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**], enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap the colored tag of the module you use for calibration, **Module settings** is displayed.
- 5 Tap the list field **Sensor name** and select the sensor you use, if necessary select the external temperature probe.
- 6 Adjust the **General settings** and **Calibration settings** if necessary.
- 7 Immerse the electrode/sensor in the sample, tap [**Calibrate**].

5.2 Direct measurement

- 1 Select a suitable sensor for the measurement and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**], enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap the colored tag of the module you use for measurement, **Module settings** is displayed.
- 5 Tap the list field **Sensor name** and select the sensor you use, if necessary select the external temperature probe.
- 6 Adjust **General settings** and **Measurement settings** if necessary.
- 7 Tap [**Save**].
- 8 Immerse the electrode/sensor in the sample, tap [**Read**].

5.3 Starting a method

- 1 Select a suitable sensor for the method and connect it to the module (ISM sensors are detected automatically, tap [**Home**] and go on with step 4).
- 2 Tap [**Sensors**], tap [**New**] enter **Parameters**. If an external temperature probe is used, repeat this step and enter the parameters of the probe.
- 3 Tap [**Save**], tap [**Home**].
- 4 Tap [**Methods**], select a method or tap [**New**] to use a **Method type** template. Use the default parameters or modify them according to your requirements.
- 5 Tap [**Save**].
- 6 Immerse the electrode/sensor in the sample, tap [**Start**].

6 Maintenance and Care

Do not open the housing of the instrument; it does not contain any parts that can be maintained, repaired or replaced by the user. If you experience problems with your instrument, contact your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

6.1 Cleaning the Instrument

NOTICE

Danger of damage to the instrument due to inappropriate cleaning agents!

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is sensitive to some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). If liquids enter the housing they can damage the instrument.

- 1 Use only water and a mild detergent to clean the housing.
- 2 Wipe off any spills immediately.
- 3 The instrument is IP54 splash water proof: Do not immerse the instrument in liquid.

- The instrument is turned off and disconnected from the electrical outlet.
- Clean the housing of the instrument using a cloth dampened with water and a mild detergent.

6.2 Maintenance of electrodes

The instrument monitors the condition of the attached pH electrodes.

Slope: 95-105%
and offset: $\pm$ (0-20) mV
Electrode is in good condition

Slope: 90-94%
or offset: $\pm$ (20-35) mV
Electrode needs cleaning

Slope: 85-89%
or offset: $\pm$ (>35) mV
Electrode is defective or too old

When cleaning, always follow the instructions in the manual of the electrodes used. Make sure the pH electrode is always kept filled with the appropriate filling solution. For maximum accuracy, any filling solution that may have "crept" and encrusted the outside of the electrode should be removed with deionized water. Always store the electrode according to the manufacturer's instructions and do not allow it to dry out.

If the electrode slope falls rapidly, or if the response becomes sluggish, the following procedures may help. Try one of the following, depending on your sample.

Problem	Action
Fat or oil build-up	Either rinse the membrane with soap solution or acetone/ethanol or shortly soak the tip of the electrode in hot water. When rinsed with organic solvent, place the membrane overnight in 0.1 mol/L HCl.
pH electrode membrane has dried out	Soak the tip of the electrode overnight in 0.1 mol/L HCl. If this procedure has no effect, soak the tip of the electrode for a few minutes in reactivation solution for pH electrodes.
Protein build-up in the diaphragm of a pH electrode	Remove deposits by soaking the electrode in an HCl/pepsin solution.
Silver sulfide contamination of pH electrode	Remove deposits by soaking the electrode in a thiourea solution.

Run a new calibration after treatment.

Note

- Cleaning and filling solutions should be handled with the same care as that given to toxic or corrosive substances.
- The condition of the pH electrode can also be checked, using the provided METTLER TOLEDO method Sensor Test.

6.3 Transporting the instrument

Note the following instructions when transporting the instrument to a new location:

- Transport the instrument with care to avoid damage! The instrument may be damaged if not transported correctly.
- Unplug the instrument and remove all connected cables.
- Remove the electrode arm.
- To avoid damage to the instrument when transporting it over long distances, please use the original packaging.
- If the original packaging is no longer available, choose packaging that will ensure safe handling.

6.4 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

7 Accessories

Accessories	Order No.
SevenExcellence™ pH/mV module	30034472
SevenExcellence™ conductivity module	30034473
SevenExcellence™ pH/Ion module	30034471
SevenExcellence™ DO/BOD module	30034474
SevenExcellence™ blank module	30034475
SevenExcellence™ protective film (2 pcs.)	30041155
uPlace™ electrode arm, complete with base and extension	30019823
SevenExcellence™ semi-transparent cover	30041154
InMotion sample changer bundle: Flex base and 100 mL rack (18 samples)	30094120
Rondolino sample changer, complete for SevenExcellence™	51302888
Rondolino USB-TTL box	30046261
Magnetic stirrer uMix™	30040002
Compact stirrer, incl. 2 propeller stirrers	51109150
Cable for compact stirrer	30098212
Disposable beakers for InMotion and Rondolino (PP, 100 mL, 1400 pcs.)	00101974
PowerShower™	51108219
LogStraight™ fingerprint reader	51192107
Barcode reader	21901297
USB cable for barcode reader	21901309
Printer USB-P25	11124301
EasyDirect pH PC software	-
USB A-B cable 1.8 m for EasyDirect pH software (not included with EasyDirect pH software)	51191926
Adapter for auxiliary instruments (USB-RS232 adapter)	51105856

Buffers and standards	Order No.
pH 4.01 buffer sachets, 30 x 20mL	51302069
pH 4.01 buffer solution, 6 x 250mL	51350018
pH 7.00 buffer sachets, 30 x 20mL	51302047
pH 7.00 buffer solution, 6 x 250mL	51350020
pH 9.21 buffer sachets, 30 x 20mL	51302070
pH 9.21 buffer solution, 6 x 250mL	51350022
pH 10.01 buffer sachets, 30 x 20mL	51302079
pH 10.00 buffer solution, 6 x 250mL	51350024
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)	51302068
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)	51302080
10 µS/cm conductivity standard solution, 250 mL	51300169
84 µS/cm conductivity standard solution, 250 mL	51302153
500 µS/cm conductivity standard solution, 250 mL	51300170
1413 µS/cm conductivity standard solution, 30 x 20 mL	51302049
1413 µS/cm conductivity standard solution, 6 x 250 mL	51350096
12.88 mS/cm conductivity standard solution, 30 x 20 mL	51302050
12.88 mS/cm conductivity standard solution, 6 x 250 mL	51350098
Zero oxygen tablets (24 pcs.)	51300140
pH sensors with fixed cable	Order No.
InLab® Expert Pro-ISM, rugged 3-in-1 pH sensor, PEEK shaft, ATC	30014096
ISM® sensors with MultiPin™ head	Order No.
InLab® Routine Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, refillable	51344055
InLab®Micro Pro-ISM, 3-in-1 pH sensor, glass shaft, 5 mm shaft diameter, ATC, refillable	51344163
InLab® Power Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, pressurized SteadyForce™ reference system	51344211
InLab® Pure Pro-ISM, 3-in-1 pH sensor, glass shaft, immovable glass sleeve, ATC, refillable	51344172
InLab® Science Pro-ISM, 3-in-1 pH sensor, glass shaft, movable glass sleeve, ATC, refillable	51344072
InLab® Solids Pro-ISM, 3-in-1 pH sensor, glass shaft, open junction, sharp membrane, ATC	51344155
Electrode cable 1.2 m, BNC/RCA-MultiPin™	30281896
Electrode cable 3 m, BNC/RCA-MultiPin™	30281897
ATC probe, temperature sensor	12997876
Digital ISM® pH sensors with diagnostic functions	Order No.
InLab® Smart Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, pressurized SteadyForce™ reference system	30027775
Cable InLab® Digital, 1.2 m	30027776

Conductivity sensors with fixed cable	Order No.
InLab® 731-ISM, 4-electrode graphite conductivity sensor, ATC	30014092
InLab® 741-ISM, 2-electrode steel conductivity sensor, ATC	30014094
InLab® Trace, high-accuracy sensor for low conductivities, ATC	30014097
Flow-Cell for InLab® Trace	30014098
InLab® Trace Kit (Sensor and Flow-Cell)	30014099
Polarographic Dissolved Oxygen sensors with fixed cable	Order No.
InLab® 605-ISM-2 m	51344611
InLab® 605-ISM-5 m	51344612
InLab® 605-ISM-10 m	51344613
Optical Dissolved Oxygen sensors with fixed cable	Order No.
InLab® OptiOx, 1.8 m	51344621
InLab® OptiOx, 5 m	51344622
InLab® OptiOx, 10 m	51344623
OptiOx parts	Order No.
OptiOx replacement cap	51344630
OptiOx calibration tube	51344631
OptiOx protective guard	51344632
OptiOx BOD adapter	51344633
Solutions	Order No.
HCl/Pepsin solution (removes protein contamination), 1 x 250 mL	51350100
Thiourea solution (remove silver sulfide contamination), 1 x 250 mL	51350102
Reactivation solution for pH electrodes, 1 x 25 mL	51350104
InLab Storage solution, 1 x 250 mL	30111142

8 Technical Data

Screen		Color TFT
Interfaces		RS232, USB A, USB B, Ethernet
Ambient conditions	Ambient temperature	5 °C to 40 °C
	Relative humidity	5% to 80% (non-condensing)
	Overtoltage category	Class II
	Pollution degree	2
	Range of application	For indoor use only
	Maximum operating altitude	Up to 2000 m
Dimensions	Width	235 mm
	Depth	188 mm
	Height	75 mm
Weight	Basic device	1120 g
	1 module	111 - 130 g
Power rating instrument	Input voltage	12 V $\overline{\text{=}}$
	Power consumption	10 W
Power rating AC adapter	Line voltage	100 - 240 V $\sim \pm 10\%$
	Input frequency	50/60 Hz
	Input current	0.3 A
	Output voltage	12 V $\overline{\text{=}}$
	Output current	0.84 A
Materials	Housing	ABS/PC
	Expansion units	ABS/PC
	Electrode arm	ABS/PC
	Touch panel	Tempered glass
pH/mV module		
pH mode	Parameter	Measurement range
	pH	-2.000...20.000
	mV	± 2000.0
ISFET mode	Temperature °C	-30.0...130.0
	pH	0.000...14.000
	Temperature °C	-30.0...130.0
Conductivity module		
	Parameter	Measurement range
	Conductivity	0.001...999999 $\mu\text{S/cm}$
	TDS	0.001...1000 ppt, g/L 0.001...999999 mg/L, ppm
	Salinity	0.0...80 psu, ppt
	Resistivity	0.01...100 M Ω -cm 0.01...999999 Ω -cm
	Temperature °C	-30.0...130.0

pH/Ion module	Parameter	Measurement range
pH mode	pH	-2.000...20.000
	mV	± 2000.0
	Temperature °C	-30.0...130.0
Ion mode	Ion	0...999999 mg/L, ppm
		0...100 mol/L, %
		0...100000 mmol/L
		-2.000...20.000 pX
	Temperature °C	-30.0...130.0
DO/BOD module	Parameter	Measurement range
Optical sensor (digital)	Dissolved oxygen	0.000...50 mg/L (ppm)
	DO Saturation	0.0...500%
	Temperature °C	0.0 °C...50.0 °C
	Pressure	500 to 1100 mbar
Polarographic sensor (analog)	Dissolved oxygen	0.000...99 mg/L (ppm)
	DO Saturation	0.0...600%
	Temperature °C	0.0 °C...60.0 °C
	Pressure	500 to 1100 mbar

Quick Check

Content of delivery		Active	Inactive
SevenExcellence™ instrument			
Fixing screw for the electrode arm			
Country-specific power supply			
Semi-transparent module cover			
Measurement module(s) (depending on kit configuration)	pH/Ion		
	Conductivity		
	pH/mV		
	DO/BOD		
Sensor(s) (depending on kit configuration)	All sensors delivered		
	All cables delivered		
uPlace™ electrode arm	Base with plastic covers (inserted)		
	Movable electrode holder		
	Extension shaft part		
	Wrench		
	Two screws		
CD-ROM with operating instructions and installation instructions (all languages)			
Installation guide (print version, all languages)			
Declaration of conformity			
Calibration certificate			
Location of the instrument		Active	Inactive
Device is not exposed to powerful vibrations			
Device is not exposed to direct sunlight			
Device is not exposed to atmospheric humidity greater than 80%			
Device is not exposed to corrosive gas atmosphere			
Device is not exposed to temperatures below 5 °C and above 40 °C			
Device is not exposed to powerful electric or magnetic fields			
Power supply		Active	Inactive
Power supply is connected			
Main power is turned on			
Instrument is switched on			
Display shows the home screen			
Modules		Active	Inactive
Module(s) is/are attached			
Module(s) is/are displayed on the home screen			
Check box(es) of the module(s) is/are activated			
Semi-transparent module cover is attached			

Sensors		Active	Inactive
Sensor(s) is/are attached			
Mini-LTW socket fits precisely (digital ISM® sensors)			
BNC socket fits precisely (pH, ORP and ISE electrodes)			
RCA socket fits precisely (temperature probes)			
Mini-DIN socket fits precisely (conductivity sensors)			
Recognition of peripheral devices		Active	Inactive
ISM electrode recognized			
Stirrer recognized			
Printer recognized (only US-P25)			
PC recognized with LabX / EasyDirect pH			
USB stick recognized			
Rondolino / InMotion recognized			
uPlace™ electrode arm		Active	Inactive
Extension shaft part used			*
Mounted as stand alone			*
Attached to the instrument	Right		*
	Left		*
Fitted correctly and stable			
Note			
If you have any questions answered in the negative, read the corresponding chapter of the installation instructions. If the problem persists, contact your vendor.			
* Irrelevant for the functionality of the instrument.			

Instrument:	
Serial Number:	
Date:	Date:
Company:	Company:
Name:	Name:
Surname:	Surname:
Signature	Signature
Retain on file for reference.	

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© Mettler-Toledo GmbH 04/2018