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Acetic acid content in vinegar

The acetic acid content is determined in vinegar by titration with sodium hydroxide.

Sample	- Vinegar, 0.7 – 1g
Preparation and procedures	- Approx. 1g vinegar in a titration beaker. - Add 50mL of deionized water. - Titrate it with 0.1 mol/L sodium hydroxide.
Compound	- Acetic acid, CH ₃ COOH - M = 60.05 g/mol, z = 1
Chemicals	- Deionized water - Sodium hydroxide
Titrant	- Sodium hydroxide, NaOH - c(NaOH) = 0.1 mol/L
Standard	- Potassium hydrogen phthalate,
Instruments, Accessories	- T50/T70/T90 Titration Excellence, G20 Compact Titrator, - Rondo or Rondolino Sample Changer, - DV1010 10 mL glass burette, - 100 mL PP titration beakers, XS205 Balance, LabX [®] pro titration software
Indication	- DG115-SC (Combined pH glass electrode)
Results	- Acid content: 4.738 ± 0.007% srel = 0.152% n = 6
Waste	- Neutralization before final disposal as aqueous solution.
Comments	- Use safety goggles, wear gloves and a lab coat. - The method parameter has been optimized for the sample used in this application. It may be necessary to slightly adapt the method to your specific sample. - Adjustment of the DG115-SC pH glass electrode. The electrode is adjusted using pH buffers 4.01, 7.00 and 9.21 before starting the analysis. - The titer of sodium hydroxide is determined using potassium hydrogen phthalate as a primary standard; 80mg is dissolved into 50mL deionized water and titrated. - The rinsing time was defined to 2s. In this way, the electrode is cleaned with deionized water before starting the subsequent sample.
Literature	- Mettler-Toledo Titration Application M400-2009
Author	Thomas Hitz, MSG, July 2006. Ralph Egli, MSG, May 2009. Compiled by: Vineesh Pallath, IMSG AnaChem, June 2012 – v1.0 Revised: V. Gärtner – MSG AnaChem, August 2012

Titration method

001 Title
 Type General titration
 Compatible with T50/T70/T90
 ID M400
 Title Acetic acid content in vinegar

002 Sample
 Number of IDs 1
 ID 1 Vinegar
 Entry type Weight
 Lower Limit 0.7g
 Upper Limit 1.0g

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Stir
 Speed 35%
Predispense
 Mode Volume
 Volume 1 mL
 Wait time 5 s
Control
 Control User
 Titrant addition Dynamic
 dE (set value) 12 mV
 dV (min) 0.002 mL
 dV (max) 0.5 mL
 Meas. Value acquisition Equilibrium controlled
 dE 0.5 mV
 dt 1 s
 t (min) 3 s
 t (max) 30 s
Evaluation and recognition
 Procedure Standard
 Threshold 5 pH/mL
 Tendency Positive
 Ranges 0
 Add EQP criteria No


Termination
 At Vmax 10.0 mL
 At potential No
 At slope No
 After number of recognized EQPs Yes
 Number of EQPs 1
 Combined termination criteria No
Accompanying stating
 Accompanying stating No
Condition
 Condition No

006 Calculation R1
 Result Content
 Result unit %
 Formula $R1=Q \cdot C/m$
 Constant C $M/(10 \cdot z)$
 M $M[\text{Acetic acid}]$
 z $z[\text{Acetic acid}]$
 Decimal places 3
 Result limits No
 Record statistics Yes
 Extra statistical func. No
 Send to buffer No
 Condition No

007 Record
 Results Per series
 Raw results Per series
 Table of meas. Values Last titration function
 Sample data Per series
 Resource data No
 E - V Last titration function
 dE/dV - V Last titration function
 log dE/dV - V No
 dE/dV - V No
 Beta - V No
 E - t No
 V - t No
 dV/dt - t No
 T - t No
 E - V & dE/dV - V No
 V - t & dV/dt - t No
 Method No
 Series data No
 Condition No

008 End of sample

Titration curve


Vitamin C Content in Orange Juice: Voltametric Indication	
Content determination of vitamin C (ascorbic acid) in orange juice by voltametric titration with DPI. The titration is monitored by a polarized DM143 electrode.	
Sample	- Orange Juice, 5 g
Preparation and procedures	<ul style="list-style-type: none"> - The juice is well homogenised and is weighed directly into the titration beaker and diluted with 50mL deionised water. - After addition of 50mL deionised water the sample can be purged with nitrogen gas to avoid oxidation of vitamin C - The pH value is adjusted to pH 3 with 2% oxalic acid by titrating the sample using an EP titration function.
Compound	<ul style="list-style-type: none"> - L(+) ascorbic acid, $C_6H_8O_6$ - $M = 176.13 \text{ g/mol}$, $z = 2$
Chemicals	<ul style="list-style-type: none"> - 50mL Deionized water, - 2% Oxalic acid for adjustment to pH 3
Titrant	<ul style="list-style-type: none"> - 2,6-dichlorophenolindophenol, DPI - $c(1/2 \text{ DPI}) = 0.01 \text{ mol/L}$
Standard	- Ascorbic acid, 0.01 mol/L
Instruments, Accessories	<ul style="list-style-type: none"> - T50/T70/T90 Titration Excellence, G20 Compact Titrator, - Rondo Sample Changer, Cover Up, T-Box with connecting cable - DV1010 10 mL glass burette, - 100 mL PP titration beakers, XS205 Balance, LabX[®] pro titration software
Indication	- DM143-SC (Double pin platinum electrode)
Results	- Vitamin C: 34.816 ± 0.334 mg/100g srel = 0.960% n = 5
Waste	- Neutralization before final disposal as aqueous solution.
Comments	<ul style="list-style-type: none"> - Use safety goggles, wear gloves and a lab coat. - The method parameters have been optimized for the sample used in this application. It may be necessary to slightly adapt the method to your specific sample. - Both DPI and ascorbic acid standard solution are not stable. To get accurate results, it is recommended to prepare them fresh before use. - The DPI titrant has to be stored in brown glass bottles to protect it from light. - Ascorbic acid is sensitive to light, temperature and oxygen. Titration beakers with light protection as well as use of purge gas are recommended. Keep in the dark. - The sample size depends on the amount of vitamin C: 1 -15 mg vitamin C corresponds to 0.5 – 8.5mL 0.01 mol/L ½ DPI - This method describes the voltametric determination of L(+) ascorbic acid with a polarized sensor DM143-SC. An alternative the current of 10 -12 µA is applied to the double pin platinum electrode and the potential is measured.
Literature	- Mettler-Toledo Titration Application M411-2006
Author	<p>Thomas Hitz, MSG, July 2006.</p> <p>Compiled by: Vineesh Pallath, MSG AnaChem, June 2012 – v1.0</p> <p>Revised: V. Gärtner – MSG AnaChem, August 2012</p>

Titration method

001 Title			
Type	General titration		
Compatible with	T50/T70/T90		
ID	M411		
Title	Vitamin C content		
.			
002 Sample			
Number of IDs	1		
ID 1	Orange Juice		
Entry type	Weight		
Lower Limit	0.0g		
Upper Limit	10.0g		
.			
003 Titration stand (Manual stand)			
Type	Manual stand		
Titration stand	Manual stand1		
004 Stir			
Speed	35%		
Duration	10 s		
Condition	No		
005 Titration (EQP) [1]			
 Titrant			
Titrant	¼ DPI		
Concentration	0.01 mol/L		
 Sensor			
Type	Polarized		
Sensor	DM143-SC		
Unit	mV		
Indication	Voltametric		
Ipol	12 µA		
 Stir			
Speed	35%		
 Predispense			
Mode	Volume		
Volume	1 mL		
Wait time	2 s		
 Control			
Control	User		
Titrant addition	Incremental		
dV	0.05 mL		
Meas. Value acquisition	Equilibrium controlled		
dE	2.0 mV		
dt	1 s		
t (min)	2 s		
t (max)	30 s		
 Evaluation and recognition			
Procedure	standard		
Threshold	500 mV/mL		
Tendency	Negative		
Ranges	0		
Add EQP criteria	No		
Termination			
At Vmax	10.0 mL		
At potential	No		
At slope	No		
After number of recognized EQPs	Yes		
Number of EQPs	1		
Combined termination criteria	No		
Accompanying stating			
Accompanying stating	No		
 Condition			
Condition	No		
006 Calculation R1			
Result	Content		
Result unit	mg/100g		
Formula	$R1=Q*C/m$		
Constant C	$100*M/z$		
M	$M[\text{Ascorbic acid}]$		
z	$z[\text{Ascorbic acid}]$		
Decimal places	3		
Result limits	No		
Record statistics	Yes		
Extra statistical func.	No		
Send to buffer	No		
Condition	No		
007 Calculation R2			
Result	Content		
Result unit	%		
Formula	$R2=Q*C/m$		
Constant	$C=M/(10*z)$		
M	$M[\text{Ascorbic acid}]$		
Z	$z[\text{Ascorbic acid}]$		
Decimal places	4		
Result limits	No		
Record statistics	Yes		
Extra statistical function	No		
Send to buffer	No		
Condition	No		
008 Record			
Results	Per series		
Raw results	Per series		
Table of meas. Values	Last titration function		
Sample data	Per series		
Resource data	No		
E - V	Last titration function		
dE/dV - V	Last titration function		
log dE/dV - V	No		
dE/dV - V	No		
Beta - V	No		
E - t	No		
V - t	No		
dV/dt - t	No		
T - t	No		
E - V & dE/dV - V	No		
V - t & dV/dt - t	No		
Method	No		
Series data	No		
Condition	No		
009 End of sample			

Titration curve

Acid Content in Wines and fruit juices	
General method for the determination of acidity of wines and juices. The method could easily be adapted for other acid determinations.	
Sample	- Wines and Juices, 10mL
Preparation and procedures	- Place 10mL of well homogenized juice or wine into a titration beaker. - Add 40mL of deionized water. - Titrate with 0.1 mol/L NaOH.
Compound	- Malic acid, M = 134.09 g/mol, z = 2 - Citric acid, M = 192.13 g/mol, z = 3 - Tartaric acid, M = 150.09 g/mol, z = 2
Chemicals	- 40mL Deionized water,
Titrant	- Sodium hydroxide, NaOH - c(NaOH) = 0.1 mol/L
Standard	- Potassium hydrogen phthalate
Instruments, Accessories	- T50/T70/T90 Titration Excellence, G20 Compact Titrator, - Rondo or Rondolino Sample Changer, - DV1010 10 mL glass burette, - 100 mL PP titration beakers, XS205 Balance, LabX [®] pro titration software
Indication	- DG111-SC (Combined pH glass electrode)
Results	- Acid content: 6.380 ± 0.01076 g/L srel = 0.169% n = 6
Waste	- Neutralization before final disposal as aqueous solution.
Comments	- Use safety goggles, wear gloves and a lab coat. - The method parameters have been optimized for the sample used in this application. It may be necessary to slightly adapt the method to your specific sample. - The sample was degassed in an ultrasonic bath for 5 minutes, in order to get rid of CO ₂ . - Temperature was measured before the titration and the endpoint values were automatically corrected accordingly.
Literature	- Mettler-Toledo Titration Application M561
Author	A.Aichert, N.Spiru, MSG, Compiled by: Vineesh Pallath, IMSG AnaChem, June 2012 – v1.0 Revised: V. Gärtner – MSG AnaChem, August 2012

Titration method

001 Title	Type	General titration		006 Calculation R1	Result	Volume
	Compatible with	T50/T70/T90			Result unit	mL
	ID	M561			Formula	R1=VEQ
	Title	Acid content			Constant C	1
				M	M[None]
002 Sample	Number of IDs	1			z	z[None]
	ID 1	Apple Juice			Decimal places	3
	Entry type	Fixed volume			Result limits	No
	Volume	10 mL			Record statistics	Yes
				Extra statistical func.	No
003 Titration stand (Manual stand)	Type	Manual stand			Send to buffer	No
	Titration stand	Manual stand1			Condition	No
004 Stir	Speed	50%		007 Calculation R2	Result	Acid content
	Duration	60 s			Result unit	g/L
	Condition	No			Formula	R2=Q*C/m
005 Titration (EP) [1]					Constant	C=M/z
	Titrant				M	M[Malic acid]
	Titrant	NaOH			Z	z[Malic acid]
	Concentration	0.1 mol/L			Decimal places	3
	Sensor				Result limits	No
	Type	pH			Record statistics	Yes
	Sensor	DG111-SC			Extra statistical function	No
	Unit	pH			Send to buffer	No
	Stir				Condition	No
	Speed	35%		008 Record	Results	Per series
	Predispense				Raw results	Per series
	Mode	Volume			Table of meas.Values	Last titration function
	Volume	3 mL			Sample data	Per series
	Wait time	0 s			Resource data	No
	Control				E - V	Last titration function
	Mode	Absolute			dE/dV - V	Last titration function
	Tendency	Positive			log dE/dV - V	No
	End point value	8.1 pH			dE/dV - V	No
	Control band	1.0 pH			Beta - V	No
	Dosing rate [max]	1.0 mL/min			E - t	No
	Dosing rate [min]	10 µL/min			V - t	No
	Termination				dV/dt - t	No
	At EP	Yes			T - t	No
	Termination delay	0 s			E- V & dE/dV - V	No
	At V max	20.0 mL			V - t & dV/dt - t	No
	Max.time	infinite			Method	No
	Accompanying stating				Series data	No
	Accompanying stating	No			Condition	No
	Condition					
	Condition	No		009 End of sample		

Titration curve

Not Available

Free Sulphur Dioxide (SO ₂) content in Wine	
The content of free sulphur dioxide (SO ₂) in wine is determined by redox titration with iodine as a titrant. The titration is monitored using a Pt double pin electrode DM143-SC at a fixed polarization current (voltametric indication)	
Sample	- 50 mL White wine
Preparation and procedures	<ul style="list-style-type: none"> - 50mL wine is transferred into a titration beaker with a pipette. - 5mL of 10% potassium iodide solution is added. This step has been performed using an additional dosing unit. It can be performed manually using a pipette. - 5mL 20% sulfuric acid is added to the titration beaker. This step has been performed using an additional dosing unit. If this is not possible, this step can be done manually with a pipette.
Compound	<ul style="list-style-type: none"> - Sulphur dioxide, SO₂, - M = 64.06 g/mol, z = 2
Chemicals	<ul style="list-style-type: none"> - 5mL 10% Potassium iodide, KI - 5mL 20% Sulfuric acid, H₂SO₄
Titrant	<ul style="list-style-type: none"> - Iodine, I₂ - c(1/2 I₂) = 0.02 mol/L
Standard	- Ascorbic acid, C ₆ H ₈ O ₆ ,
Instruments, Accessories	<ul style="list-style-type: none"> - T50/T70/T90 Titration Excellence, G20 Compact Titrator, - 2 Additional dosing units. - DV1010 10 mL glass burette, - 100 mL PP titration beakers, XS205 Balance, LabX[®] pro titration software
Indication	- DM143-SC (Double pin platinum electrode)
Results	- Sulphur dioxide content: 21.185 ± 0.067 mg/L srel = 0.315% n = 6
Waste	- Neutralize the acid waste before final disposal.
Comments	<ul style="list-style-type: none"> - Use safety goggles, wear gloves and a lab coat. - The method parameters have been optimized for the sample used in this application. It may be necessary to slightly adapt the method to your specific sample. - The titer determination is performed using pure ascorbic acid as a primary standard. - In order to avoid loss of SO₂, the sample must be taken from a freshly opened bottle. After opening the bottle, free SO₂ can evaporate with CO₂ or be oxidized while in contact with air resulting in too low values. - Working with the sample changer leads to SO₂ losses from the sample beakers already prepared on the rack. Therefore, it is recommended to work without a sample changer to achieve a higher accuracy. - Iodine is also reduced by other wine components. These competing reactions can partly be delayed by the addition of 5 mL 10% potassium iodide solution. - The reaction only takes place in acidic medium. Therefore, 5 mL 20% H₂SO₄ must be added immediately before titration.
Literature	- Mettler-Toledo Titration Applications M419-2006, M563 and M564
Author	<p>Claudia Schreiner, MSG, July 2006. Compiled by: Vineesh Pallath, IMSG AnaChem, June 2012 – v1.0 Revised: V. Gärtner – MSG AnaChem, August 2012</p>

Titration method

001 Title
 Type General titration
 Compatible with T50/T70/T90
 ID M419
 Title Free SO₂ content in wine

002 Sample
 Number of IDs 1
 ID 1 White wine
 Entry type Fixed volume
 Volume 50mL
 Density 1.0 g/mL

003 Titration stand (Manual stand)
 Type Manual stand
 Titration stand Manual stand1

004 Dispense (normal) [1]
 Titrant 10%KI
 Concentration 1 mol/L
 Volume 5.0 mL
 Dosing rate 60.0 mL/min

005 Dispense (normal) [2]
 Titrant 20% H₂SO₄
 Concentration 1 mol/L
 Volume 5.0 mL
 Dosing rate 60.0 mL/min

006 Stir
 Speed 50%
 Duration 10 s
 Condition No

007 Titration (EP) [1]
Titrant
 Titrant $\frac{1}{2}$ I₂
 Concentration 0.01 mol/L
Sensor
 Type Polarized
 Sensor DM143-SC
 Unit mV
 Indication Voltametric
 Ipol 10 μ A
 Temperature acquisition
 Temperature acquisition No
Stir
 Speed 50%
Predispense
 Mode None
 Wait time 10 s
Control
 Control User
 End point type Absolute
 Tendency Negative
 Endpoint value 100 mV
 Control band 30 mV
 Dosing rate (max) 1.25 mL/min
 Dosing rate (min) 100 μ L/min

Termination
 At EP Yes
 Termination delay 5 s
 At Vmax 20.0 mL
 Max. time 600 s

008 Calculation R1
 Result Content
 Result unit mg/L
 Formula $R1=Q \cdot C / m$
 Constant C $M \cdot 1000 / z$
 M $M[\text{Sulfur dioxide}]$
 z $z[\text{Sulfur dioxide}]$
 Decimal places 3
 Result limits No
 Record statistics Yes
 Extra statistical func. No
 Send to buffer No
 Condition No

009 Record
 Results Per series
 Raw results Per series
 Table of meas. Values Last titration function
 Sample data Per series
 Resource data No
 E - V Last titration function
 dE/dV - V Last titration function
 log dE/dV - V No
 dE/dV - V No
 Beta - V No
 E - t No
 V - t No
 dV/dt - t No
 T - t No
 E- V & dE/dV - V No
 V - t & dV/dt - t No
 Method No
 Series data No
 Condition No


010 End of sample

G20 Compact Titrator Method

The G20 method includes slightly modifications

008 Calculation R1
 Result type Predefined
 Calculation type Direct titration
 Formula $R1=Q \cdot C / m$
 Result Content
 Result unit mg/L
 Constant C $M \cdot 1000 / z$
 M $M[\text{Sulfur dioxide}]$
 z $z[\text{Sulfur dioxide}]$
 Decimal places 3
 Result limits No
 Record statistics Yes

Titration curve


Citrate content determination using a Cu-Ion selective electrode	
Determination of citrate content in mineral water drinks by Potentiometric titration with copper sulphate and a copper ion selective electrode (Cu-ISE) as indicating electrode.	
Sample	- Mineral water drinks, containing citrate or citric acid.
Preparation and procedures	<ul style="list-style-type: none"> - Degas the mineral water for 5 mins first then pass it through a cation exchange resin (Fill a glass tube with approx. 20 g strongly acidic cation exchange resin i.e. Dowex 50 WX8 and rinse it with deionised water, eg. 10 mL. Mineral water drink is run through the column containing the exchange resin, and collected into a beaker. The column is then rinsed with 10 mL deionised water). - Transfer 25mL of sample into a titration beaker. - Methanol, borate buffer and NaOH solutions are added by additional dosing units.
Compound	<ul style="list-style-type: none"> - Citric acid, $C_3H_5O(COOH)_3$ - $M = 192.12 \text{ g/mol}$ $z = 1$ (1:1 soluble complex in solution)
Chemicals	<ul style="list-style-type: none"> - Deionised water, Methanol. - 0.025 mol/L borax ($Na_2B_4O_7 \cdot 10H_2O$), pH adjusted to 9.0 with 0.1 mol/L HCl. - Strongly acidic cation exchange resin (Dowex 50 WX8) - 1.0 mol/L Sodium Hydroxide (NaOH), 0.1 mol/L Hydrochloric acid (HCl)
Titrant	<ul style="list-style-type: none"> - Copper sulphate, $CuSO_4$ - $c(CuSO_4) = 0.05 \text{ mol/L}$
Standard	- Sodium thiosulphate
Instruments, Accessories	<ul style="list-style-type: none"> - T50/T70/T90 Excellence Titrator. - 2 X DV1020 20 mL glass burette, 2 X DV1010 10 mL glass burette, 3 additional dosing units, Rondo 20 sample changer with powershower™ SP250 Peristaltic pump. - 100 mL PP titration beakers (ME-101974), LabX® pro titration software.
Indication	- DG115-SC pH glass electrode, perfectION™ combination Cu-ISE
Results	- Citrate content : 3.068 ± 0.013 g/L srel = 0.437 % n =6
Waste	- Comply with Federal, State and local regulations. Sweep up crystals, powder or insoluble copper carbonate and dispose as special waste.
Comments	<ul style="list-style-type: none"> - Use safety goggles, wear gloves and a lab coat. - It may be necessary to adapt the method to your specific sample. - Maintain the sample solution at pH 9.4 with 1 mol/L sodium hydroxide. Borate buffer is used to avoid interferences. - Copper(II) ions Cu^{2+} forms 1:1 stable complex with citrate anions. The free copper(II) concentration is monitored with the combined Cu-ISE. - The electrode is conditioned in deionised water after every titration. If necessary gently polish the Cu-ISE solid state membrane.
Literature	- Mettler-Toledo Titration Application M472-2010
Author	<p>Vineesh Pallath/Claudia Schreiner/Cosimo DeCaro.</p> <p>Compiled by: Sohail Ansari, IMSG AnaChem, April 2012 - v1.0</p> <p>Revised: V. Gärtner – MSG AnaChem, August 2012</p>

Titration method

001 Title	Type	General titration	Control	Control	User
	Compatible with	T50/T70/T90		Titrant addition	Dynamic
	ID	VP0076dispNAOH		dE(set value) [mV]	6.0
	Title	Citrate Content		dV(min) [mL]	0.02
			dV(max) [mL]	0.2
002 Sample	Number of IDs	1		Mode	Equilibrium controlled
	ID 1	7UP with extraction		dE[mV]	0.5
	Entry type	Fixed Volume		dt[s]	2.0
	Volume [mL]	25		t(min)[s]	10
			t(max)[s]	30
003 Titration stand (Rondo 60)	Type	Rondo/Tower A	Evaluation and recognition	Procedure	Standard
	Titration stand	Rondo60/1A		Threshold[mV/mL]	35
	Lid handling	No		Tendency	Positive
004 Dispense (normal) [1]	Titrant	Methanol		Ranges	0
	Concentration	100		Add. EQP criteria	No
	Volume [mL]	25		
005 Dispense (normal) [2]	Titrant	Borax Buffer 9	Termination	At Vmax[mL]	10
	Concentration	100		At potential	No
	Volume [mL]	25		At slope	No
006 Dispense (normal) [3]	Titrant	NaOH		After number of recognized EQPs	No
	Concentration	1		Combined termination criteria	No
	Volume [mL]	2.5	Accompanying stating	Accompanying stating	No
007 Stir	Speed	35%	Condition	No
	Duration	10 s		010 Calculation R1	
008 Measure (normal) [1]	Sensor			Result	Content
	Type	pH		Result unit	g/L
	Sensor	DG1115-SC		Formula	R1=EQE*C/m
	Unit	pH		Constant C=	M/z
	Temperature acquisition			M	M[Citrate]
	Temperature Measurement	No		z	z[Citrate]
	Stir			Decimal places	3
	Speed	35%		
	Acquisition of measured values	Equilibrium controlled	011 Calculation R2	Result	pH value
	Acquisition			Result unit	
	dE [mV]	0.5		Formula	R2=E[1]
	dt [s]	1		Constant C=	1
	t(min) [s]	2		M	M[None]
	t(max) [s]	30		z	z[None]
	Mean value	Yes		Decimal places	2
	No of measured values	4		
	dt [s]	1	012 Rinse	Auxiliary reagent	Water
	Condition			Rinse cycles	1
	Condition	No		Vol. per cycle[mL]	10
009 Titration (EQP) [1]	Titrant			Position	Current position
	Titrant	CuSO ₄		Drain	Yes
	Concentration[mol/L]	0.5		Drain pump	Drain pump 1
	Sensor			Condition	No
	Type	ISE	013 Conditioning	Type	Fix
	Sensor	Cu ISE		Interval	1
	Unit	mV		Position	Conditioning beaker
	Ion Charge	2		Time[s]	60
	Temperature acquisition			Speed[%]	30
	Temperature Measurement	No		Condition	No
	Stir		014 End of sample		
	Speed	30%	015 Park	Titration stand	Rondo60/1A
	Predispense			Position	Conditioning beaker
	Mode	Volume		Condition	No
	Volume [mL]	6			
	Wait time [s]	10			

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- ▶ www.mt.com/gwp 用於稱量
- ▶ www.mt.com/gtp 用於滴定
- ▶ www.mt.com/gpp 用於移液
- ▶ www.mt.com/gdrp 用於測定密度和折射率

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© 09/2012 Mettler-Toledo AG
Mettler-Toledo AG