# Quick Setup Guide Transmitter M400

---

## Multilingual Quick Setup Guide

en	Quick Setup Guide	3
de	Quick-Setup-Leitfaden	27
fr	Guide de paramétrage rapide	51
it	Guida alla configurazione rapida	75
es	Guía de configuración rápida	99
pt	Guia de Configuração Rápida	123
ru	Руководство по быстрой настройке	147
ja	クイック セットアップ ガイド	171

**METTLER TOLEDO**

© 02/2023 METTLER TOLEDO.

Subject to technical changes. Printed in Switzerland. 52 121 395 I

# Quick Setup Guide

## Transmitter M400

---

### Content

1 Operation	4
2 Menu Structure	5
3 Terminal Block (TB) definitions	6
4 Wiring example for pH Transmitter	9
5 General Setup (applies for all parameters)	10
6 pH Calibration incl. ISFET	11
7 O <sub>2</sub> Calibration	13
8 CO <sub>2</sub> Calibration (InPro 5000)	14
9 CO <sub>2</sub> Hi (High) Calibration (InPro 5500i)	17
10 M400 Type 1 Cond Ind only: Conductivity Calibration for inductive Sensors	19
11 Conductivity Quick Setup	20
12 Conductivity Calibration for 2-e and 4-e Sensors	21
13 TDL GPro™ 500	22
14 Environmental protection	26

# 1 Operation

## Entry of data values, selection of data entry options

Use the ▲ key to increase or the ▼ key to decrease a digit. Use the same keys to navigate within a selection of values or options of a data entry field.

**Note:** Some screens require configuring multiple values via the same data field (ex: configuring multiple setpoints). Be sure to use the ► or ◀ key to return to the primary field and the ▲ or ▼ key to toggle between all configuration options before entering to the next display screen.

## Navigation with ↑ in Display

If a ↑ is displayed on the bottom right hand corner of the display, you can use the ► or ◀ key to navigate to it. If you click [ENTER] you will navigate backwards through the menu (go back one screen). This can be a very useful option to move back up the menu tree without having to exit into the measuring mode and re-enter the menu.

## Exit menu

**Note:** Exit the menu at any time by pressing the ◀ and ► key simultaneously (ESCAPE). The transmitter returns to the Measurement mode.


## “Save changes” dialog

Three options are possible for the “Save changes” dialog:

- “Yes & Exit”: Save changes and exit to measuring mode.
- “Yes & ↑”: Save changes and go back one screen.
- “No & Exit”: Don’t save changes and exit to measuring mode.

The “Yes & ↑” option is very useful if you want to continue configuring without having to re-enter the menu.

## 2 Menu Structure


\* Only available in combination with ISM sensors.

\*\* Do not use after configuration

### 3 Terminal Block (TB) definitions

Power connections are labeled –N for Neutral and +L for Line, for 100 to 240 VAC or 20–30 VDC.

TB2 for ½DIN	
1	AO1+
2	AO1–/AO2–
3	AO2+
4	AO3+
5	AO3–/AO4–
6	AO4+
7	DI1+
8	DI1–/DI2–
9	DI2+

#### Note:

This is a 4-wire-product with an active 4–20 mA analog output.

Please do not supply to Pin1–Pin6 of TB2.

TB2 for ½DIN			
1	NO1	8	NC5
2	COM1	9	COM6
3	NC1	10	NO6
4	NO2	11	NO3
5	COM2	12	COM3
6	NC2	13	NO4
7	COM5	14	COM4

**NO** = normally open (contact is open if un-actuated).

**NC** = normally closed (contact is closed if un-actuated).

#### TB4 – ISM (digital) Sensors

		TDL GPro™ 500	Optical Oxygen, InPro™ 5500 i (CO <sub>2</sub> )	InPro 6860 i (VP8 cable)	pH, amp. Oxygen, Cond 4-e, InPro 5000 i
Term.	Function	Color	Color	Color	Color
1	24 VDC	–	brown	grey	–
2	GND (24 VDC)	–	black	blue	–
3	1-Wire	–	–	–	transparent (cable core)
4	GND (5 VDC)	–	green/yellow	green/yellow	red (shield)
5	–	–	–	–	–
6	GND (5 VDC)	brown	–	–	–
7	RS485-B	yellow	blue	brown	–
8	RS485-A	green	white	pink	–
9	5 VDC	–	–	–	–

### TB3 – Conventional (analog) Sensors

Cond 2-E or 4-E			amp. Oxygen	
Term.	Function	Color*	Function	Color
1	Cnd inner1	white	–	–
2	Cnd outer1	white/blue	Anode	red
3	Cnd inner2	blue	–	–
4	Cnd outer2/Shield	black	Shield/GND	green/yellow
5	–	–	Cathode	transparent
6	RTD ref/GND	bare shield	GND/NTC	green
7	RTD sense	red	–	–
8	RTD	green	NTC	white
9	+5V	–	+5V	–

Terminal 4 and 6 are internally connected, either terminal can be used to connect a wire.

\* Transparent not connected.

### TB3 – Conventional (analog) Sensors

pH			Redox (ORP)	
Term.	Function	Color*	Function	Color
1	Glass	transparent	Platinum	transparent***
2	–	–	–	–
3	Reference**	red	Reference	Shield
4	Solution GND/Shield**	green/yellow and blue	–	–
5	–	–	–	–
6	RTD ref/GND	white	–	–
7	RTD sense	–	–	–
8	RTD	green	–	–
9	+5V	–	–	–

Terminal 4 and 6 are internally connected, either terminal can be used to connect a wire.

\* Grey wire not used.  
 \*\* For pH without solution ground, install jumper 3 to 4.  
 \*\*\* Remove black coating.

### TB3 – Conventional (analog) Sensors

Cond Ind			
Term.	Color. InPro 7250 ST/PFA	Color. InPro 7250 HT	Function
1	Coax inner/transparent	Coax inner/transparent	receive hi
2	red	yellow	receive lo
3	green/yellow	green/yellow	shield/GND
4	brown, white	violet	send lo
5	blue	black	send hi
6	white	white	RTD ref/GND
7	grey	grey	RTD sense
8	green	green	RTD
9	–	–	not used

Terminal 4 and 6 are internally connected, either terminal can be used to connect a wire.

### TB3 – Conventional (analog) Sensors

Dissolved carbon dioxide		
Term.	Color*	Function
1	Coax inner/transparent	Glass
2	–	not used
3**	Coax shield/red	Reference
4**	green/yellow	GND/Shield
5	–	not used
6	white	RTD ref/GND
7	–	not used
8	green	RTD
9	–	+5V

Terminal 4 and 6 are internally connected, either terminal can be used to connect a wire.

- \* Grey wire not used.
- \*\* Jumper 3 to 4 has to be installed.

### TB3 – Conventional (analog) Sensors

ISFET		
Term.	Color	Function
1	Coax inner/pink	FET
2	–	not used
3*	yellow	Reference
4*	green/yellow	GND/Shield
5	–	not used
6	white	RTD ref/GND
7	–	not used
8	grey	RTD
9	brown	+5V


Terminal 4 and 6 are internally connected, either terminal can be used to connect a wire.

- \* Jumper 3 to 4 has to be installed.

## 4 Wiring example for pH Transmitter (using TB3)

pH measurement with monitoring of glass electrode.

- 9 +5V
- 8 RTD
- 7 RTD sense
- 6 RTD ret/GND
- 5 not used
- 4 Solution GND/Shield
- 3 Reference
- 2 not used
- 1 Glass

**Note:** Wire colors only valid for connection with VP cable, grey not connected.

## 5 General Setup (applies for all parameters) (PATH: Menu/Quick Setup)


While in Measurement mode press the [MENU] key to bring up the Menu selection. Select QUICK SETUP and press the [ENTER] key.

### Display Convention:

1st line on display → a      3rd line on display → c  
2nd line on display → b      4th line on display → d

Select the units of measurement for a and b. Only lines a and b can be configured in QUICK SETUP. Go to the Configuration Menu to configure lines c and d.

### Channel Selection

### Please select the type of Sensor:

Analog: For conventional analog sensors (will be displayed on channel "A").


ISM: For ISM sensors (will be displayed on channel "B").

### Please select the parameter requirement:

The choice of parameter depends on the level of transmitter. If an ISM sensor is selected, the setting "Auto" means, all possible ISM sensors will be recognized and accepted. If a special parameter is chosen, only this parameter will be recognized and accepted on the transmitter.

**By using analogue 2-e or 4-e conductivity sensors,  
see "Conductivity Quick Setup" below for intermediate steps.**

### Analog outputs


By selecting YES the linear 4–20 mA analog output Aout1 will be set up for measurement a when [ENTER] is pressed. Selecting NO means that no analog output is set up.

Aout1 min, Aout1 max are the minimum and maximum measurement values for the 4 and 20 mA values respectively.


**Note for multi channel devices:** The user can configure output A out 3 to measurement c by going back to the previous menu and selecting c.

## Setpoints


After configuring the Analog Output a setpoint can be configured for that output. If NO is selected and [ENTER] is pressed then the quick setup is done and the menus are exited without setting up any setpoint.

Selecting YES and pressing [ENTER] means a setpoint can be configured for channel a.

Following setpoint Types can be selected:

- High (High value has to be set).
- Low (Low value has to be set).
- Between (High and Low value has to be set).
- Outside (High and Low value has to be set).

After setting the setpoint value(s) a Relay (none, 1, 2, 3, 4) can be configured for that setpoint. The Relay delay is set to 10 seconds and the Hysteresis is set to 5%. If a relay is configured as CLEAN relay, it will not be selectable in this menu.


**For multi-channel devices:** Also setpoints 5 and 6 can be configured and the user can configure a setpoint to measurement c by going back through previous menus and selecting c.

## 6 pH Calibration incl. ISFET (PATH: Cal)

For Display setup, analog outputs and setpoints see the section “General Setup”.

### pH (resp. ISFET) Two point Calibration


**Remark:** When using ISFET sensors, you have to adjust the zero point first.


While in Measurement mode press the ► key. If the display prompts you to enter the calibration security code, press the ▲ or ▼ key to set the calibration security code, then press the [ENTER] key to confirm the calibration security code.

Press the ▲ or ▼ key to select the pH (resp. ISFET) calibration sub function.


A flashing “H” shows the ongoing calibration process.


Select 2 POINT CALIBRATION by pressing the [ENTER] key.

Place the electrode in the first buffer solution and then press the [ENTER] key.


As soon as the stabilisation criteria have been fulfilled (or [ENTER] was pressed in manual mode) the display changes and prompts you to place the electrode in the second buffer solution.


As soon as the stabilisation criteria have been fulfilled (or [ENTER] was pressed in manual mode) the display changes to show the slope calibration factor S and the offset calibration factor Z. Select ADJUST to save the calibration values and the successful Calibration is confirmed on the display. If you select CALIBRATE, the values will not be taken. If an ISM Sensor is used, the Calibration data will be stored in the calibration history.

## Process calibration


Select PROCESS CALIBRATION by pressing the ▲ key once followed by the [ENTER] key. To show the ongoing Calibration Process an "H" is displayed in the top left hand corner.

The "H" changes to "A" (or "B", depends on type of sensor) if Process Calibration is selected to show the user the ongoing calibration on Channel "A" or "B".

Take a sample and press the [ENTER] key again to store the current measuring Value.

After determining the pH Value of the Sample press the ► key again to proceed with the calibration. If the display prompts you to enter the calibration security code, press the ▲ or ▼ key to set the calibration security code, then press the [ENTER] key to confirm the calibration security code.


Enter the pH value of the sample then press the [ENTER] key to start calibration.

After the calibration the slope calibration factor S and the offset calibration factor Z are displayed. Select ADJUST to save the new calibration values and the successful Calibration is confirmed on the display. If you select CALIBRATE, the values will not be taken. If an ISM Sensor is used, the Calibration data will be stored in the calibration history. The "A" in the top left hand corner disappears.

## 7 O<sub>2</sub> Calibration (PATH: Cal)

For Display setup, analog outputs and setpoints see the section "General Setup".


### O<sub>2</sub> Calibration

**For multi-channel devices:** Using the ▲ and ▼ key on the "Channel A" field lets the user change the channel to be calibrated.


Enter Calibration mode by pressing the ► key.

A flashing "H" (alternating with "A" or "B" to show the channel being calibrated for multi-channel devices) in the top left hand corner shows the ongoing calibration process.

A DO sensor calibration is always a one point calibration either in Air (Slope) or a zero (Offset) calibration.

It is possible to select a SLOPE or ZERO CALIBRATION. A one point slope calibration is done in air and a one point offset calibration is done at 0 ppb DO. Press the [ENTER] key after selecting SLOPE or OFFSET.

Enter the value for Point 1 including a decimal point. DO is the value being measured by the transmitter and sensor in the units set by the user. Press [ENTER] when this value is stable to perform the calibration.

After the calibration the slope calibration factor S and the offset calibration factor Z are displayed.

H 21.7 %sat  
A 25.0 °C  
Calibration Successful

Select ADJUST to save the calibration values and the successful calibration is confirmed on the display. If you select CALIBRATE, the values will not be taken. If an ISM Sensor is used, the Calibration data will be stored in the calibration history.

## 8 CO<sub>2</sub> Calibration (InPro 5000) (PATH: Cal)

For Display setup, analog outputs and setpoints see the section "General Setup".

### CO<sub>2</sub> One point Calibration

A 180.4 hPa  
A 26.1 °C  
Calibrate Sensor  
Channel A CO2 ↑

While in measurement mode press the key [CAL]. Press the ▲ or ▼ key to select the CO<sub>2</sub> CALIBRATION and press the [ENTER] key to call the sub function.

H 180.4 hPa  
A 26.1 °C  
CO2 Calibration  
Type = 1 Point ↑

Select 1 POINT and press [ENTER].

H 137.5 hPa  
A 26.1 °C  
Press ENTER when  
Sensor is in Buffer 1 ↑

Place the sensor in the buffer solution and press the [ENTER] key to start the calibration.

A 154.5 hPa  
A 26.1 °C  
A Point1 = 7.00 pH  
A CO2 = 7.07 pH ↑

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

H 154.5 hPa  
A 26.1 °C  
pH S=100.0% Z=7.04pH  
Save Adjust ↑

As soon as the drift conditions have stabilized the display changes to show the slope calibration factor S and the offset calibration factor Z.

After a successful calibration, the calibration values are taken over (Adjust) or were aborted (Calibrate or Abort).

## CO<sub>2</sub> Two point Calibration

Start as in chapter CO<sub>2</sub> ONE POINT CALIBRATION described.

```
A 154.5 hPa
A 26.1 °C
CO2 Calibration
Type = 2 Point ↑
```

Select 2 POINT CALIBRATION.

```
H 137.5 hPa
A 26.1 °C
Press ENTER when
Sensor is in Buffer 1 ↑
```

Place the sensor in the first buffer solution and press the [ENTER] key to start the calibration.

```
154.5 hPa
A 26.1 °C
A Point1 = 7.00 pH
A CO2 = 7.07 pH ↑
```

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

```
122.4 hPa
A 26.1 °C
Press ENTER when
Sensor is in Buffer 2 ↑
```

As soon as the drift conditions have stabilized, the display changes and prompts you to place the electrode in the second buffer.

```
2.8 hPa
A 26.1 °C
A Point2 = 9.21 pH ...
A CO2 = 8.80 pH ↑
```

Place the sensor in the second buffer solution and press the [ENTER] key to go on with the calibration.

The display shows the second buffer the transmitter has recognized (Point 2) and the measured value.

```
2.8 hPa
A 26.1 °C
pH S=74.21% Z=6.948pH
Save Adjust ↑
```

As soon as the drift conditions have stabilized the display changes to show the slope calibration factor S and the offset calibration factor Z.

After a successful calibration, the calibration values are taken over (Adjust) or were aborted (Calibrate or Abort).

## CO<sub>2</sub> Process Calibration

Start as in chapter CO<sub>2</sub> ONE POINT CALIBRATION described.

```
A 17.3 hPa
A 27.3 °C
CO2 Calibration
Type = Process ↑
```

Select PROCESS CALIBRATION.

```
A 17.3 hPa
A 27.3 °C
A Point1 = 00000 hPa
A CO2 = 17.3 hPa ↑
```

Take a sample and press the [ENTER] key again to store the current measuring value. To show the ongoing calibration process, A or B (depending on the channel) is blinking in the display. After determining the CO<sub>2</sub> value of the sample, press the ► key again to proceed with the calibration.

```
A 17.3 hPa
A 27.3 °C
A Point1 = 16.90 hPa
A CO2 = 17.3 hPa ↑
```

Enter the CO<sub>2</sub> value of the sample then press the [ENTER] key to start calibration.

```
A 17.3 hPa
A 27.3 °C
pH S=100.0% Z=7.009pH
Save Adjust ↑
```

The display shows the slope calibration factor S and the offset calibration factor Z.

After a successful calibration, the calibration values are taken over (Adjust) or were aborted (Calibrate or Abort).

## 9 CO<sub>2</sub> Hi (High) Calibration (InPro 5500 i) (PATH: Cal)

For Display setup, analog outputs and setpoints see the section “General Setup”.

### One point calibration

```
B 189.0 hPa
B 25.0 °C
Calibrate Sensor
Channel B CO2 Hi ↑
```

While in measurement mode press the key [CAL]. Press the ▼ or ▲ key to select the CO<sub>2</sub> Hi and press the [ENTER] key to call the sub function.

```
B 189.0 hPa
H 25.0 °C
CO2 Hi Calibration
Type = 1 Point ↑
```

Select “1 point” and press [ENTER].

```
B 189.0 hPa
H 25.0 °C
Press ENTER when
Sensor is in Gas (CO2) ↑
```

Place the sensor in the calibration gas (e.g. 100% CO<sub>2</sub>) resp. solution. Press [ENTER].

```
B 189.0 hPa
H 25.0 °C
B Point1=1013. hPa .
B CO2=189.0 hPa ↑
```

Enter the value for Point 1 including a decimal point and units. The value in the second text line is the value being measured by the transmitter and sensor in the units selected by the user.

Depending on the Drift control (see User Manual chapter 8.2.3.10 “CO<sub>2</sub> Hi Parameters”), one of the two following modes is active.

```
B 189.0 hPa
H 25.0 °C
CO2 S= 9.28 mV EL= 253 mV
Save Adjust ↑
```

If Drift is auto, as soon as the signal of stabilization criteria have been fulfilled, transmitter will display the calibration result automatic.

If Drift is Manual, press [ENTER] to proceed as soon as the value is stable, transmitter will display the calibration result after manual press [ENTER]

If “Adjust” or “Calibrate” are chosen, the message “Calibration successful” is displayed. In any case you will get the message “Re-install sensor” and “Press ENTER” on the display. After pressing “ENTER” the M400 returns to the measuring mode.

## Process calibration

Start as in chapter One point Calibration described.

```
B 189.0 hPa
H 25.0 °C
CO2 Hi Calibration
Type Process ↑
```

Select Process Calibration

```
B 189.0 hPa
B 25.0 °C
Press ENTER to Capture
B CO2=189.0 hPa ↑
```

Take a sample and press the [ENTER] key again to store the current measuring value. To show the ongoing calibration process, A or B (depending on the channel) is blinking in the display.

```
B 189.0 hPa
H 25.0 °C
B Point1=1013. hPa
B CO2=189.0 hPa ↑
```

After determining the CO<sub>2</sub> value of the sample press the [CAL] key again to proceed with the calibration.

Enter the CO<sub>2</sub> value of the sample then press the [ENTER] key to start the calculation of the calibration results.

```
B 189.0 hPa
H 25.0 °C
CO2 S= 9.28 mV EL= 253 mV
Save Adjust ↑
```

After the calibration the slope calibration factor S and the offset calibration factor Z are displayed.

In case of a successful calibration, the calibration values are stored in the cal history\* and taken over (Adjust), stored in the cal history\* and not taken over (Calibrate) or discarded (Abort).

If "Adjust" or "Calibrate" are chosen, the message "Calibration successful" is displayed. The M400 returns to the measuring mode.

## 10 M400 Type 1 Cond Ind only: Conductivity Calibration for inductive Sensors (PATH: Cal)

For Display setup, analog outputs and setpoints see the section "General Setup".

### Cond Ind Zero point Calibration

A 1.25  $\mu\text{S}/\text{cm}$   
A 25.00  $^{\circ}\text{C}$   
Calibrate Sensor  
Channel A Conductivity ▲

While in measurement mode press the key [CAL]. Press the ▲ or ▼ key to select the CONDUCTIVITY CALIBRATION and press the [ENTER] key to call the sub function.

A 1.25  $\mu\text{S}/\text{cm}$   
A 25.00  $^{\circ}\text{C}$   
Cal Compensation  
Standard ▲

Select the kind of calibration (standard, linear or n. water).

H 40.5  $\text{mS}/\text{cm}$   
A 23.9  $^{\circ}\text{C}$   
Conductivity Calibration  
Type = Zero Point ↑

Select ZERO POINT and press [ENTER]. Remove the sensor and dry it.

A 1.035  $\text{mS}/\text{cm}$   
A 21.9  $^{\circ}\text{C}$   
A Point 1 = 0.000  $\text{mS}/\text{cm}$   
A C = 1.035  $\text{mS}/\text{cm}$  ↑

Press the [ENTER] key again to go on with the calibration.

As soon as the conditions have stabilized (or [ENTER] pressed in the manual mode) the display changes to show the multiplier and adder.

After a successful calibration, the calibration values are taken over and stored in the cal history (Adjust), only stored in the cal history (Calibrate) or the calibration can be aborted (Abort).

## Cond Ind One point Slope Calibration

Start as in chapter COND IND ZERO POINT CALIBRATION described. Select the kind of calibration.

Select 1 POINT SLOPE and press [ENTER].

	<b>217.4</b>	$\mu\text{S}/\text{cm}$
A	<b>25.0</b>	$^{\circ}\text{C}$
Conductivity Calibration		
Type = 1 Point Slope ↑		

	<b>217.4</b>	$\mu\text{S}/\text{cm}$
A	<b>25.0</b>	$^{\circ}\text{C}$
A Point 1 =	215.0	$\mu\text{S}/\text{cm}$
A	C =	217.4 $\mu\text{S}/\text{cm}$ ↑

Enter the value of calibration point 1 and then press the [ENTER] key to start calibration.

The cell multiplier and the adder are displayed. After a successful calibration, the calibration values are taken over and can be stored in the cal history (Adjust), only stored in the cal history (Calibrate) or the calibration can be aborted (Abort).

## 11 Conductivity Quick Setup (PATH: Menu/Quick Setup)

For Display setup, analog outputs and setpoints see the section "General Setup".

### Sensor Type Selection

	<b>1.25</b>	$\mu\text{S}/\text{cm}$
	<b>25.0</b>	$^{\circ}\text{C}$
Sensor Type = Cond (2) ▲		

Select the type of sensor to be used with the M400 transmitter. Choices are "Cond(2)", used for all 2-Electrode type sensors and "Cond (4)" for all 4-electrode sensors. Press the [ENTER] key.

### Cell Constant

	<b>1.25</b>	$\mu\text{S}/\text{cm}$
	<b>25.0</b>	$^{\circ}\text{C}$
p	M=0.1003	A=0.0000
s	M=1.0000	A=0.0000

Enter the appropriate cell constant(s): from the sensor label or certificate (M) for 2-electrode sensors, leaving (A) at 0.000; or (M) and (A) values for 4-electrode sensors. Press the [ENTER] key.

### Measurement units

	<b>1.25</b>	$\mu\text{S}/\text{cm}$
	<b>25.0</b>	$^{\circ}\text{C}$
a	S/cm	
Analog Output? Yes ▲		

Select the measurement (conductivity or temperature) and units for measurement.

If using analog output, select YES. Refer back to section "General Setup" to continue setup.

## 12 Conductivity Calibration for 2-e and 4-e Sensors (PATH: Cal)

For Display setup, analog outputs and setpoints see the section "General Setup".

### Cond One point Calibration

A 1.25  $\mu\text{S}/\text{cm}$   
A 25.0  $^{\circ}\text{C}$   
Calibrate Sensor  
Channel A Conductivity ▲

While in Measurement mode press the [CAL] key. If the display prompts you to enter the calibration security code, press the ▲ or ▼ key to set the calibration security code, then press the [ENTER] key to confirm the calibration security code.

Press [ENTER] to select the CONDUCTIVITY CALIBRATION SUB FUNCTION.

A flashing "H" (alternating with "A" or "B" to show the channel being calibrated for multi-channel devices) in the top left hand corner shows the ongoing calibration process.

1.25  $\mu\text{S}/\text{cm}$   
25.0  $^{\circ}\text{C}$   
Cal Compensation  
Standard ▲

Choose the compensation mode by using the ▲ or ▼ key.

Press [ENTER] to confirm the compensation mode.

1.25  $\mu\text{S}/\text{cm}$   
25.0  $^{\circ}\text{C}$   
Conductivity Calibration  
Type = 1 point ▲

Select 1 POINT instead of 2 POINT CALIBRATION by pressing the ▼ key once followed by the [ENTER] key.

**CAUTION:** Rinse sensors with a high-purity water solution before every calibration to prevent contamination of the reference solutions.

1.25  $\mu\text{S}/\text{cm}$   
25.0  $^{\circ}\text{C}$   
A Point1 = 0.055  $\mu\text{S}/\text{cm}$   
A C = 0.057  $\mu\text{S}/\text{cm}$  ▲

Place the electrode in the reference solution.

Enter the value of Point 1 and press the [ENTER] key.

<b>1.25</b>	$\mu\text{S}/\text{cm}$
<b>25.0</b>	$^{\circ}\text{C}$
C M=0.1000 A=0.0000	
Save Calibration? Yes	▲

After the calibration of the cell multiplier or slope calibration factor "M" i.e. cell constant and the Adder or offset calibration factor "A" are displayed.

Select YES to save the calibration values and the successful Calibration is confirmed on the display.

The user gets the message RE-INSTALL SENSOR and "Press [ENTER]" on the display.

After pressing [ENTER] the M400 returns to the measuring mode.

## 13 TDL GPro™ 500

**For Display setup, analog outputs and setpoints see the section "General Setup"**

**Note:** TDL is calibrated at the factory and does not require calibration at installation & Startup.

### TDL Installation (PATH: Quick Setup/TDL/Installation)

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
MENU		
Quick Setup		↑

While in measurement mode press the key [MENU]. Press the ▲ or ▼ key to select the TDL and then the Installation menu item.

B	<b>0.0</b>	%V G1
B	<b>0.0</b>	%V G2
Transmission: 000.00%		
NSL:000%		↑

In this mode, the current live value of the % transmission and the Noise Signal Level (NSL) are displayed during 5 minutes until it automatically returns to the Measurement mode. Best performance of the TDL is obtained when the % transmission value is maximized (50% or more) and the NSL is minimized (50 or less). Use these values to rotate the blue sensor head attached with a loose clamp connection on to the probe until the maximum transmission is found. In this position, and secure the blue sensor head into position and tighten the clamp.

### TDL Commissioning (PATH: Quick Setup/TDL/Commissioning)

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
O2-TDL		
Commissioning		↑

While in measurement mode press the key [MENU]. Press the ▲ or ▼ key to select the TDL and then the Commissioning menu item.

First, select the type of pressure compensation selected:

- External: current external pressure value coming from a pressure transducer of 4...20 mA analog output

– Fixed: pressure compensation uses a fixed value to be set manually.

**Note:** if this pressure compensation mode is selected, a considerable gas concentration measurement error resulting from a non- realistic pressure value can take place.

B	<b>20.9</b>	%V O <sub>2</sub>
B	<b>20.9</b>	%Trm
Pressure-External		↑

If External compensation is selected, then the minimum (4 mA) and maximum (20 mA) analog output signals from the pressure transducer must be mapped to the corresponding Analog input of the TDL. Key in the minimum and maximum values of the pressure in the following units:

- hPa
- mmHg
- mbar
- psi
- kPa

B	<b>20.9</b>	%V O <sub>2</sub>
B	<b>20.9</b>	%Trm
Ain	4mA = 940.0	mbar
Ain	20mA = 2000.	mbar
		↑

In general, METTLER TOLEDO recommends the use of absolute pressure transducers for more accurate signal compensation over a broad pressure range. If, however, small pressure variations around atmospheric pressure are to be expected, relative pressure sensors will produce better results; but the variations of the underlying barometric pressure will be ignored.

For relative pressure sensors, the minimum and maximum values must be mapped so that the TDL can interpret the analog pressure signal as "absolute", i.e. a fixed barometric pressure of 1013 mbar (for example) has to be added to the mapped values.

B	<b>20.9</b>	%V O <sub>2</sub>
B	<b>20.9</b>	%Trm
Pressure=fixed		↑

If Fixed compensation is selected, the fixed pressure value with which the measurement signal will be calculated has to be keyed in manually. For the fixed pressure, the following units can be used:

- hPa
- mmHg
- mbar
- psi
- kPa

B	<b>20.9</b>	%V O <sub>2</sub>
B	<b>20.9</b>	%Trm
Pressure=1013. mbar		↑

If External compensation is selected, then the minimum (4 mA) and maximum (20 mA) analog output signals from the temperature transducer must be mapped to the corresponding Analog input of the TDL. Key in the minimum and maximum values of the temperature in °C.

B	<b>20.9</b>	%V O <sub>2</sub>
B	<b>20.9</b>	%Trm
Temperature=External		↑

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
Ain 4mA = 0.000 °C		
Ain 20mA = 250.0 °C		
		↑

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
Temperature=Fixed		
		↑

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
Temperature=320.0 °C		
		↑

B	<b>20.9</b>	%V O2
B	<b>20.9</b>	%Trm
Pathlength=00200 mm		
		↑

If Fixed compensation is selected, the fixed temperature value with which the measurement signal will be calculated has to be keyed in manually. For the fixed temperature, only °C can be used.

Last, select the initial optical path length corresponding to the probe length installed:

- 290 mm probe: 200 mm
- 390 mm probe: 400 mm
- 590 mm probe: 800 mm


This initial value is valid when instrument purging on the instrument and on the process side is running. Depending on the process conditions and after the optimum of the process purging flow has been found (see next chapter), this value may have to be slightly adapted.

### Setting the correct process side purging

The flow rate of the purging will affect the effective path length and consequently the measurement value.

Therefore the following procedure should be used. Start with a very high flow rate and gradually decrease it. The measurement value will then start at a low value and increase with decreasing purge flow. At some point it will level out and stay constant for a while and then again start increasing. Choose a purge flow in the middle of the constant region.

### Optimizing the purge flow

On the x-axis there is purge flow and on the y-axis there is the instrument concentration reading.

- 1 Concentration reading with high purge flow. The path length is now shorter than the effective path length since the purge tubes is completely filled with purging gas and some of the purging gas is flowing into the measurement path.
- 2 Concentration reading with optimized purge flow. The path length is now equal to the effective path length since the purge tubes are completely filled with purge gas. See the illustration below.
- 3 Concentration reading with no purge flow. The path length is now equal to the nominal path length since the probe is completely filled with process gas.
- 4 The optimized purge flow.

**WARNING:** Always start purging at maximum flow before starting the process.

**WARNING:** Purging must always be switched on in order to avoid dust deposition onto the optical surfaces.

## 14 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

GPro, InPro and ISM are trademarks of the METTLER TOLEDO Group.  
All other trademarks are the property of their respective holders.

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All rights reserved