
InPro® 6900 Series O₂ Sensors

**Instruction manual
Bedienungsanleitung
Instructions d'utilisation**

English **Page** **3**

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InPro[®] 6900 Series O₂ Sensors

Instruction manual

Contents

1	Introduction.....	5
2	Important notes	6
2.1	Notes on operating instructions	6
2.2	Intended use	6
2.3	Safety instructions	7
2.4	Examples of some typical applications	8
2.5	Use in Ex-zones.....	8
2.6	Ex-classification ATEX	9
2.6.1	Introduction.....	9
2.6.2	Rated data	9
2.6.3	Special conditions.....	10
2.7	Ex-classification FM approved	11
3	Product description.....	12
3.1	General information.....	12
3.2	Principle	12
3.3	Scope of delivery	12
3.4	Equipment features	13
4	Installation.....	14
4.1	Mounting the sensor.....	14
4.2	Connection.....	14
4.2.1	Connecting the InPro 6900 to a VP cable	14
4.2.2	Connecting the VP cable to the transmitter	15
5	Operation	16
5.1	Start-up and polarizing	16
5.2	Calibration	17
5.2.1	Purpose of calibration	17
5.2.2	What you have to know for calibration	17
5.2.3	Single point calibration	19
5.2.4	Dual point calibration	19
6	Maintenance	20
6.1	Inspection of the sensor.....	20
6.1.1	Visual inspection	20
6.1.2	Testing the METTLER TOLEDO O ₂ Sensor-Master InPro 6900	21
6.1.3	Testing the sensor via a transmitter	22
6.2	Changing the electrolyte, the membrane body or the interior body.....	23
7	Storage	26
8	Product specification	26
8.1	Certificates	26
8.2	Specifications	27
9	Ordering information.....	28
9.1	Sensors	28
9.2	Accessories	28
9.3	Spare parts	29
9.4	Recommended transmitters	29
9.5	Recommended housings	29
10	Theory of the polarographic sensor	30
10.1	Introduction.....	30
10.2	Principle of the design of an oxygen electrode	31
10.3	Parameters determining current	32
10.4	Polarization voltage	33
10.5	Temperature	33
10.6	Dependence on flow.....	34
10.7	Oxygen partial pressure – oxygen concentration	35

1 Introduction

Thank you for buying the **InPro 6900 sensor from METTLER TOLEDO**.

The construction of the InPro 6900 sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment. **Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro 6900.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Important notes

2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 6900 sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended use

Mettler-Toledo InPro 6900 sensors are intended solely for inline measurement of the oxygen partial pressure in liquids and gases, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, «Section 9.3»).

- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

2.4 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

Measurement in liquids:

- Brewing
- Beverage filtration
- Filling stations

2.5 Use in Ex-zones

Attention!

For an installation in Ex-zones please read the guidelines following hereafter:

CE 1258

Ex-classification ATEX:

 Ex ia IIC T6/T5/T4/T3 Ga/Gb

 Ex ia IIIC T69°C/T81°C/T109°C/T161°C Da/Db

Number of the test certificate:

SEV 14 ATEX 0169 X

IECEX SEV 14.0026X

Ex-classification FM approved:

IS/I, II, III / 1 / ABCDEFG / T6 Ta = 60 °C
- 53 800 002; Entity

2.6 Ex-classification ATEX

2.6.1 Introduction

According to Directive 94/9/EC (ATEX 95) Appendix I, the O₂ oxygen sensors type InPro 6XXX is a devices of equipment group II, category 1/2G which, according to Directive 99/92/EC (ATEX 137) can be used in zones 0/1 or 1/2 or 1 or 2 as well as gas groups IIA, IIB and IIC, which are potentially explosive due to combustible substances in the temperature T3 to T6.

The requirements specified in EN 60079-14 must be observed during use/installation.

According to Directive 94/9/EC (ATEX 95) Appendix I, the O₂ oxygen sensors type InPro 6XXX is a devices of equipment group II, category 1/2D which, according to Directive 99/92/EC (ATEX 137) can be used in zones 20/21 or 21/22 or 21 or 22, which are potentially explosive due to combustible dust.

The requirements specified in EN 60079-14 must be observed during use/installation.

For the analog version of the O₂ oxygen electrode, the O₂ measurement circuit, temperature measurement circuit and data chip circuit are part of the common intrinsically safe system and are jointly connected to and operated by a separately certified transmitter.

The digital version of the O₂ oxygen sensor is connected to and operated by two-wire cable to the certified transmitter.

The intrinsically safe circuits are galvanically isolated from the non-intrinsically safe circuits up to a nominal voltage peak value of 375 V and from the earthed parts up to a nominal voltage peak value of 30 V.

2.6.2 Rated data

Analog O₂ Oxygen sensor

With type of protection: intrinsic safety to Ex ia IIC

O₂ measuring circuit, temperature measuring circuit and data chip circuit

Only for connection to certified intrinsically safe circuits. Maximum values:

$U_i \leq 16 \text{ V}$, $I_i \leq 190 \text{ mA}$, $P_i \leq 200 \text{ mW}$

$L_i = 0$ (effective internal inductance)

$C_i = 900 \text{ pF}$ (effective internal capacitance)

The values above apply, each as the sum of all the individual circuits of the associated intrinsically safe supply and evaluation unit (transmitter).

Digital O₂ Oxygen sensor

With type of protection: intrinsic safety to Ex ia IIC

Two-wire current circuit

Only for connection to certified intrinsically safe circuits. Maximum values:

$$U_i \leq 16 \text{ V}, I_i \leq 30 \text{ mA}, P_i \leq 50 \text{ mW}$$

$$L_i = \text{negligible}$$

$$C_i = \text{negligible}$$

2.6.3 Special conditions

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1G applications, zone 0, is shown in the following table:

Temperature class	Max. ambient or media temperature
T 6	68 °C
T 5	80 °C
T 4	108 °C
T 3	160 °C

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1D applications, zone 20, is shown in the following table:


Temperature class	Max. ambient or media temperature
T 69 °C	68 °C
T 81 °C	80 °C
T 109 °C	108 °C
T 161 °C	160 °C

- The capacitance and inductance of the connecting cable has to be considered.
- The O₂ Oxygen sensor type InPro 6XXX can be used in /with the fittings InFit 76*-*** or InTrac 7**-***, or in /with other suitable fittings in potentially explosive areas.
- The metal body of the O₂ Oxygen sensors, or the fittings InFit 76*-*** or InTrac 7**-***, or other appropriate fitting is optionally included in the routine pressure test of the system.
- The independent fitting used for installation of O₂ Oxygen sensor must be conductively connected to the equipotential bonding system.

2.7

Ex-classification FM approved

3 Product description

3.1 General information

The **sensor InPro 6900** with integrated temperature probe is used for measurement of oxygen at low concentrations.

The sensor is **sterilizable** and **compatible with CIP** (cleaning in place).

InPro 6900 sensors with ISM functionality offer Plug and Measure as well as enhanced diagnostics features.

3.2 Principle

Here is a short summary of the principle of polarographic measurement on which this sensor is based (Clark 1961).

- a) The Clark polarographic sensor basically consist of a working electrode (cathode), a counter/reference electrode (anode), and an oxygen permeable membrane which separates the electrodes from the sample medium.
- b) The transmitter supplies a constant polarization voltage to the cathode, needed to reduce oxygen.
- c) The oxygen molecules which migrate through the permeable membrane are reduced at the cathode. At the same time, oxidation takes place at the anode and oxidized anode metal (silver) is liberated as silver ions into the electrolyte. The electrolyte closes the electric circuit between anode and cathode (ion conductivity).
- d) The current produced by the reactions described above is measured by the transmitter and is proportional to the partial pressure of oxygen (pO₂) in the sample medium.

Additionally, the InPro 6900 features a guard ring cathode which guarantees rapid response by the sensor and good signal stability at low DO concentrations.

Please refer to «Section 10 - Theory of the polarographic sensor» for further information.


3.3 Scope of delivery

Each sensor is supplied fully assembled and factory-tested for correct function together with:

- an electrolyte bottle (52 201 005)
- a quality control certificate
- inspection certificates 3.1 B (complying with EN 10204.3/1 B)

3.4 Equipment features

12 mm Sensor

METTLER TOLEDO DO sensors are supplied with fitted membrane body filled with electrolyte and covered with the watering cap and have been checked for proper function.

4 Installation

4.1 Mounting the sensor


Important! Remove the watering cap before mounting the sensor.

Mounting the sensor in a housing

Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

Mounting the sensor directly on a pipe or a vessel

The **12 mm sensors** can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.

4.2 Connection

4.2.1 Connecting the InPro 6900 to a VP cable


Sensors with ISM functionality require the use of a special VP-8 cable as well as an ISM-compatible O₂ transmitter. The sensor is connected to the transmitter via a VP cable. The VP cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP68 connector housing guarantees maximum process safety.

If you have chosen to use a sensor with ISM functionality, the VP connector will incorporate an EEPROM chip. It is absolutely necessary to protect this electronic component against any electrical charge.


Do not touch the sensor at the VP connector plug.

To connect the VP cable to the sensor align the slit of the VP connector with the pin in the plug. Then tightly screw the plug to fasten the two parts.

4.2.2 Connecting the VP cable to the transmitter

Note: Cable assignment can be found in the METTLER TOLEDO VP cable instruction manual.

Note: For connecting the cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual.

5 Operation

Attention! Before using the sensors for the first time, the electrolyte should be replaced (see «Chapter 6.2»). Due to possible adverse conditions during transport and storage (e.g. airfreight; pressure and temperature variations), the quality of the electrolyte may become impaired. Poor electrolyte quality can lead to erroneous measurement values.

5.1 Start-up and polarizing

Attention! The watering cap must be removed before mounting the sensor in the process.

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O₂ Sensor-Master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O₂ transmitter or to a sensor master. After 6 hours, the sensor is fully polarized and ready for operation. During this time, it is recommended to keep the watering cap filled with conditioning solution on the sensor (see «Section 7 - Storage»), especially if the sensor is being polarized for longer than 6 hours.

A shorter polarization period is sufficient if the sensor has been disconnected for only a few minutes. The following table serves to establish the correct polarization time in relation to the depolarization time.

Depolarization time ¹⁾ t_{depol} [Min.]	Minimum required polarization time ²⁾ [Min.]
$t_{\text{depol}} > 30$	360
$30 > t_{\text{depol}} > 15$	$6 * t_{\text{depol}}$
$15 > t_{\text{depol}} > 5$	$4 * t_{\text{depol}}$
$t_{\text{depol}} < 5$	$2 * t_{\text{depol}}$

1) Depolarization time:

Time span in which the polarization voltage is cut off from the sensor. This is the case during:

- change of electrolyte
- change of membrane body
- the time the cable is disconnected or no transmitter or sensor master is connected to the cable

2) Polarization time:

Time span during which the sensor is under a polarization voltage.

Attention! Setting of the polarization voltage on the transmitter for correct measurements:

- Standard applications (e.g. measurement in biotechnology): **-675 mV**
- Measurement of permanently low oxygen concentrations <500 ppb in the presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries): **-500 mV**

Important! To ensure the supply of the correct polarization voltage the transmitter must be set accordingly.

5.2 Calibration

5.2.1 Purpose of calibration

Each oxygen sensor has its own individual slope and own individual zero point. Both values are subject to change, for example, through electrolyte consumption or after exchange of electrolyte or membrane body. To ensure high measurement accuracy of the sensor, a calibration must be carried out after each change of electrolyte or membrane. Prior to calibration, the sensor has to be polarized for at least 6 hours.

Please remove the watering cap from the sensor, rinse the sensor with water and dry it for at least 10 minutes before starting calibration.

A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.


To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100 %. If not, then the sensor needs a new calibration.

5.2.2 What you have to know for calibration

Attention! If dissolved oxygen is being measured at low concentrations (<500 ppb) in the presence of volatile acidic components (e.g. carbon dioxide (CO₂) measurements in breweries), the calibration procedure described below should be followed in order to achieve the best measurement performance.

No special calibration procedure is required if dissolved oxygen is being measured in standard applications (e.g. measurement in biotechnology).

First step: Adjust the transmitter polarization voltage to –675 mV. Detailed information is given in the transmitter instruction manual.

Second step: After adjusting the polarization voltage to –675 mV, equilibration time is required. Wait for about 5 minutes before carrying out step 3.

Third step: Perform calibration according to the transmitter' instruction manual.

Fourth step: Adjust the transmitter polarization voltage back to –500 mV.

Steps 1 to 4 can be accomplished automatically if you use the M 700 transmitter. Thus the routine can be minimized to a few keystrokes. In addition the M 700 transmitter must be equipped with the software function SW-700-011 «high CO₂ compensation». Please ask your local Mettler-Toledo distributor.

General remarks:

- **For calibration in air, the sensor membrane must be dry**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the **oxygen saturation index** of the calibration medium is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is only very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.
- For correct calibration, a minimum flow rate of the calibration medium is necessary.
- **Calibration in a fermenter should be performed after sterilization** (as sterilization may alter the sensor slope), but prior to inoculation. If it is not possible to perform the calibration after sterilization, the use of an existing membrane body that has been pre-sterilized while mounted in the sensor is recommended. A slope alteration of some per cent can occur with new membrane bodies, particularly after a first sterilization, as the tension of the membrane may be altered by the sterilization process.
- Make sure that all other parameters, such as temperature and pressure, are constant.

For continuous applications, we recommend **periodic recalibration in line with your requirements on accuracy, the type of process in operation and your own experience**. The frequency of the need for recalibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

5.2.3 Single point calibration

By carrying out a single point calibration, the factual slope of the sensor can be established. The calibration medium can be water with known oxygen saturation index (e.g. air-saturated water) or air with known water-vapor saturation (e.g. water-vapor saturated air).

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O₂, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Single point calibration should be sufficient for almost all process applications.

5.2.4 Dual point calibration

By carrying out a dual point calibration both slope and zero point of the sensor can be established.

Important! In case of a **dual point calibration, always start by the zero point calibration** before calibrating the slope.

Due to the very low zero current of METTLER TOLEDO sensors, dual point calibration is normally not necessary for standard applications.

As a rule, the zero point should be adjusted to zero manually, or it is automatically performed by the transmitter (see instruction manual for the transmitter).

A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.

Attention! Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.995%.

After the sensor signal has stabilized (after 20...30 minutes), the sensor can be calibrated through the relevant transmitter to the 100 % value of the desired measurable variable, e.g. 0 % air, 0.0 % O₂, or 0.0 ppm (see instruction manual for the transmitter).

6 Maintenance

6.1 Inspection of the sensor

6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the membrane foil optically for signs of damage. The foil must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.
Note: An undulated membrane has no influence on the sensor performance, assuming the membrane is intact.
- The membrane body must be replaced if the sensor has too long a response time, the reading is unstable or subject to drift, and if the sensor cannot be calibrated or the membrane shows sign of mechanical damage.
- Check the cathode area for discoloration, contamination or cracks in the glass. If necessary rinse with demineralized water and clean with a clean soft brush or soft paper tissue.

Attention! Do not use any cleaning agents containing alcohol. This could damage the sensor or lead to fault current.


Attention! The glass body is fragile and sensitive to vibration.

6.1.2 Testing the METTLER TOLEDO O₂ Sensor-Master InPro 6900

We recommend to use the METTLER TOLEDO O₂ Sensor-Master InPro 6900 to check the quality of your sensor as follows:

- Connect the sensor to the O₂ Sensor-Master.

As soon as the sensor is connected to the O₂ Sensor-Master, the polarization function is automatically activated. Please note: if the sensor was disconnected from the transmitter for longer than 5 minutes, the sensor must be polarized first (polarizing time see «Section 5.1») to get representative test results.

- **Battery Check:**

Push the switch to the left. If the battery is ok and the O₂ Sensor-Master is operational the green LED lights up. Otherwise, please consult the instruction manual of the O₂ Sensor-Master.

- **Sensor Check:**

For this test the O₂ sensor must be **fully polarized** and the **membrane of the sensor must be dry and clean**.

Expose the sensor (connected to O₂ Sensor-Master) to the air. By pushing the switch to the right to the position «2 – Sensor check», the O₂ Sensor-Master checks whether the electrode current for air measurement delivered by the sensor is within the admissible range (250 to 500 nA for InPro 6900).

If the green LED lights up the current for air measurement is within the admissible range.

If the LED does not light up, you should check the battery of the O₂ Sensor-Master (see instruction manual «Accessories»). If the battery is working, than there is probably a problem with your sensor. You should change the electrolyte and/or the membrane body of your sensor. If after a membrane change the LED still does not light up, this means that there is maybe something wrong with the interior body of the sensor. You should than change it (see «Section 6.2»).

Important! The Sensor Check function only verifies the correctness of the electrode current for air measurement. In order to be absolutely sure of the functionality of the sensor, the residual signal in an oxygen free medium should also be controlled (see «Section 6.1.3»).

6.1.3 Testing the sensor via a transmitter

A periodic zero current measurement (**no zero point calibration**) is recommended for verification of proper sensor function.

Note: At the time you carry out the zero current measurement, the sensor must be polarized.

Zero current measurement can be done by using zeroing gel (order no. 34 100 1032) or nitrogen (N₂) or carbon dioxide (CO₂) calibration gases with a purity of at least 99.995 %, alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 10% of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If the measured values are too high, this suggests a depleted electrolyte or a defective membrane. In the first instance replace the electrolyte, and in the second case exchange both the membrane body and the electrolyte accordingly.

If after such procedures the above mentioned values are still not reached, replace the interior body. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Many sample media contain volatile substances which, even at very low concentrations, have a clearly perceptible smell. Similarly to oxygen, these substances are able to invade the electrolyte through the gas-permeable membrane. Accordingly, they become noticeable when changing the electrolyte. In most cases, such substances have absolutely no influence on the measuring properties of the sensor. Slight discoloration of the electrolyte also has no effect on the measuring properties.

6.2 Changing the electrolyte, the membrane body or the interior body

Note: The InPro 6900 uses a special electrolyte which contains an oxygen scavenger. This electrolyte ensures fast response time and together with the guard ring increases the signal stability of the sensor. The electrolyte must be exchanged on a regular basis or when the sensor has been exposed to the air for more than 24 hours without a watering cap or if the sensor has been stored for several months.

Note: The content of the electrolyte bottle has to be used within 3 months after opening.

If the membrane exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body has to be replaced.

Warning! The O₂ electrolyte has an alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.**

If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.

When changing the electrolyte, the membrane body or the interior body, please observe the following instructions (see also the following illustration):

Attention! Make sure that this maintenance step is carried out in **clean place**.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.

Note: steps 4 to 7 may only be carried out when changing the interior body.

4. Unscrew the retainer nut of the interior body with an adjustable wrench or with a 3/8" wrench.
5. Remove the interior body by pulling it out of the sensor shaft. If necessary use a plier.


Warning! Do not twist the interior body. Otherwise the connection pins can be damaged.

6. Insert the new interior body in the sensor shaft. Turn the interior body in the shaft until the slit of the interior body is aligned with the pin placed in the shaft.
7. Press the body in the shaft and screw the new retainer nut in place.
8. Examine the O-rings visually for mechanical defects, and replace if necessary.
9. Half-fill the new membrane body with O₂ electrolyte.

Note: The electrolyte bottle is equipped with a special pouring system. To ensure proper functioning, hold the bottle vertically, upside-down.

Note: make sure that all air bubbles are removed from the membrane body. Air bubbles can be removed by carefully tapping on the membrane body.


10. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue.

Important! No electrolyte, sample media or contamination may be present between the membrane body and the cap sleeve. Please check carefully.

11. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
12. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated.

Replacement of the interior body

7 Storage

For storage periods longer than 24 hours, we recommend you to use a watering cap filled with the same cleaning and conditioning solution as that used for our portable DO measurement system (order no. 52 200 255). This solution contains an oxygen scavenger which ensures that the electrolyte of the sensor does not come into contact with oxygen when not in use.

Prepare the cleaning and conditioning solution as follows: Dissolve one tablet in 40 ml of deionized water and wait 5 minutes for the tablet to be completely dissolved. Fill the watering cap with this solution and place it over the tip of the sensor. This solution has some cleaning properties which will keep the membrane free of microorganisms. In case you do not have any cleaning and conditioning set, you may also use checking gel or deaerated water in the watering cap. Before mounting the sensor in the process, always remove the watering cap and rinse the tip of the sensor with water.

Attention! The electrolyte should be replaced before use if the storage period exceeds 3 months.

Attention! If the storage period exceeds 6 months, the sensor should be **stored dry**, i.e. without any electrolyte in the membrane body. A sensor being stored dry (without electrolyte in the membrane body) **may on no account be connected to the O₂ Sensor-Master InPro 6900.**

8 Product specification

8.1 Certificates

Each sensor is delivered with a set of **3.1 B certificates** (complying with EN 10204.3/1.B).

All wetted metal parts (sensor shaft, cap sleeve and membrane body) are identified with a engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and membrane body) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

8.2 Specifications


InPro 6900	
Measurement principle	Amperometric/Polarographic
Working conditions	
Pressure resistance measurement	0.2...6 bar absolute [2.9...87.0 psi absolute] with standard membrane body 0.2...9 bar absolute [2.9...130 psi absolute] with reinforced membrane body
Mechanical pressure resistance	max. 12 bar absolute [174.0 psi absolute]
Measuring temperature range	0...80°C [32...176°F]
Temperature range	-5...121°C [23...250°F] (sterilizable)
Construction	
Temperature compensation	Automatic with built-in RTD
Cable connection	VarioPin (IP 68) straight or angled
O-ring material	Silicone FDA and USP VI approved
Membrane material	PTFE/Silicone/PTFE (reinforced with steel mesh)
Wetted metal parts	Stainless steel Special material on request
Surface roughness of wetted metal parts (ISO 1320:1992)	N5 (RA < 0.4 µm [16 µin])
Quick disconnect interior body	Low ppb
Cathode	Pt
Anode	Ag
Guard ring	Yes
Dimensions	
Sensor diameter	12 mm [0.47"]
Immersion length (a) for 12mm sensor	70, 120, 220 mm [2.8, 4.7, 8.66"]
Performances	
Detection limit	1 ppb (water), 3 ppb (water/CO ₂)
Accuracy	± [1 % + ± 1 ppb] of reading in liquids ± [1 % + ± 3 ppb] of reading in CO ₂ containing liquids
Response time at 25°C/77°F (air → N ₂)	98 % of final value < 90 s
Sensor signal in ambient air (25°C/77°F)	250...500 nA
Residual signal in oxygen-free medium	< 0.03% of the signal in ambient air
Maximum flow error	≥ 5%
Certification	
EHEDG, 3A	Yes
3.1 B (EN 10204.3/1.B)	Yes
ATEX certificate	Yes
FM Approval	Yes
FDA / USP VI	Yes
Quality control	Yes
Compatibility	
with METTLER TOLEDO transmitters	see «Section 9.4»
with METTLER TOLEDO housings	see «Section 9.5»

9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

9.1 Sensors

Intelligent Sensor Management (ISM) O₂ sensor configuration

Example configuration:
InPro 6900/12/120/ISM
→ Sensor **with** ISM functionality
InPro 6910/25/220
→ Sensor **without** ISM functionality

Sensor without ISM functionality	Order No.
InPro 6900 / 12 / 070 (2.8")	52 200 944
InPro 6910 / 12 / 070 (2.8")	52 200 947
InPro 6900 / 12 / 120 (4.7")	52 200 945
InPro 6910 / 12 / 120 (4.7")	52 200 948
InPro 6900 / 12 / 220 (8.7")	52 200 946
InPro 6910 / 12 / 220 (8.7")	52 200 949
InPro 6900 / 12 / 060 (2.4") V80	52 200 951

9.2 Accessories

Accessory	Order No.
O ₂ Sensor-master InPro 6900	52 200 893
O ₂ Sensor-simulator	52 200 891
Cap sleeve with protective cage P-Type	52 200 038
VP cable VP6-ST/3 m	52 300 108
VP cable VP6-ST/5 m	52 300 109
VP cable VP6-HT/3 m	52 300 112
VP cable VP6-HT/5 m	52 300 113
VP cable VP8-ST/3 m	52 300 354
VP cable VP8-ST/5 m	52 300 355
VP cable VP8-HT/3 m	52 300 361
VP cable VP8-HT/5 m	52 300 362
Zeroing gel	34 100 1032

For other cable lengths or types, please contact your local Mettler-Toledo representative.

9.3 Spare parts

Spare part	Order No.
Membrane body, Single T-6900	52 201 049
Membrane kit T-6900 (4 membrane bodies, 1 O-ring set, 2 x 5 ml electrolyte bottle, watering cap)	52 201 003
Membrane body T-6900 (20 pcs.)	52 201 050
Membrane body, single T-6900 R	52 201 108
Membrane kit T-6900 R (4 membrane bodies R, 1 O-ring set, 2 x 5 ml electrolyte bottle, watering cap)	52 201 109
Membrane body T-6900 R (20 pcs.)	52 201 181
Electrolyte bottle InPro 6900 (5 ml)	52 201 005
Interior body InPro 6900 (with quick disconnect)	52 200 943
Upgrade kit 6800/6900	52 201 004

9.4 Recommended transmitters

Transmitter	Order No.
O ₂ Transmitter 4100 PA	Ask your local distributor or refer to the technical data sheet.
O ₂ Transmitter 4100 ppb	
O ₂ Transmitter 4220 X	
O ₂ Transmitter M 700	

9.5 Recommended housings

Housings	12 mm Ø	25 mm Ø	
Static housings			Please ask your local Mettler-Toledo representative.
INGOLD "safety weld-in socket"	–	✓	
InFit 761 CIP	✓	–	
Retractable housings			
InFit 777 e	✓	–	
InFit 797 e	✓	–	
Immersion housing			
InDip 550	✓	–	

Note: The housings are available in different versions. Please contact your distributor to get the right ordering information.

10 Theory of the polarographic sensor

10.1 Introduction

Two types of electrodes are employed in analytical work: **potentiometric** and **amperometric** electrodes.


- Potentiometric electrodes develop a voltage generated by the activity of a particular ion. Examples of such electrodes are glass electrodes (like pH electrodes) and most ion-selective electrodes. Their individual potentials cannot be determined. The measurable quantity is the difference of potential between the measuring electrode and an inert reference electrode. The potential of the reference electrode must be constant.

All potentiometric electrodes are subject to Nernst's law and for this reason electrodes and measuring instruments are in most cases interchangeable. An important requirement of potentiometric measurements is the virtually currentless determination of the electrode voltage. During measurement no chemical reaction occurs and the solution remains in equilibrium.


- In the case of **amperometric electrodes**, such as the **oxygen electrode**, activity measurement is based on a current measurement.

The oxygen electrode consists of a cathode and an anode conductively connected by an electrolyte. A suitable polarization voltage between the anode and the cathode selectively reduces the oxygen at the cathode.

Cathode reaction:

Anode reaction:

These chemical reactions result in an electric current which is proportional to the oxygen partial pressure (PO₂). The oxygen electrode consumes oxygen which is continuously extracted from the solution. The viscosity and flow rate of the solution are therefore important parameters.

The electrode current of an oxygen electrode is determined not only by the oxygen partial pressure but by many other electrode parameters. The electrode currents of different electrode types may differ by several powers of ten. For this reason oxygen electrodes and amplifiers cannot be freely interchanged.


10.2 Principle of the design of an oxygen electrode

There are two main types of oxygen electrodes:

- Electrodes **without** membrane
- Electrodes **with** gas-permeable membrane (Clark Principle)

The membrane electrode according to Clark is today most widely used. As compared to the electrode without membrane it possesses the following advantages:


- Oxygen measurement in gases and solutions
- No mutual contamination of electrode and solution
- No or little dependency on flow

In the case of the Clark electrode, geometrical configuration is very important. In particular, the thickness of the electrolyte film between the cathode and the membrane must be within narrow tolerances so as to ensure good linearity and a low zero current (current in a nitrogen atmosphere). The figure above exemplifies the principle of the design of a Clark-type oxygen electrode.

Principle of the design of the InPro 6900

The InPro 6900 design is based on a Clark electrode. It features an additional cathode ring connected to a separate polarization voltage. This additional cathode, also designated as “guard ring”, consumes the oxygen diffusing from the bulk electrolyte toward the measuring cathode. This together with our electrolyte ensures exceptional signal stability and a very fast response time of the sensor.

10.3 Parameters determining current

The quantity of oxygen diffused in and the magnitude of the electrode current are influenced by the following parameters:

- Oxygen partial pressure of the solution
- Membrane material and thickness
- Size of cathode
- Polarization voltage
- Temperature
- Flow conditions in the solution


Fick's law gives the mathematical relationship between these parameters:

$$I = k \cdot D \cdot a \cdot A \cdot \frac{pO_2}{X}$$

- I = electrode current
 k = constant
 D = coefficient of diffusion of O₂ in the membrane
 a = solubility of O₂ in the membrane material
 A = cathode surface
 pO₂ = oxygen partial pressure in the solution
 X = thickness of gas-permeable membrane

10.4 Polarization voltage

The voltage between the anode and the cathode is so selected that the oxygen is fully (> A, see polarogram) reduced while other gases are unaffected (< D). The ideal voltage for the Pt/Ag/AgCl system is between -500 and -750 mV.

The polarization voltage should remain as constant as possible. Besides a constant source of voltage the following prerequisites must be met:

The electrical resistance of electrolyte film must not exceed a particular value so as to avoid a voltage drop. The anode must have a large surface so as to prevent polarization of the anode by the electrode current.

10.5 Temperature

The temperature dependence of the current passing through an oxygen electrode when referred to a constant O₂ partial pressure is determined mainly by the properties of the gas-permeable membrane.

10.6 Dependence on flow

With most oxygen electrodes the electrode current is smaller in stagnant solutions than in agitated ones. The consumption of oxygen by the electrode results in extraction of oxygen from the solution in the close vicinity of the cathode outside the membrane. The oxygen is replaced by diffusion. If the electrode current is strong, the solution cannot fully restore the oxygen by diffusion. This results in an electrode current weaker than would correspond to conditions in the solution. In agitated solutions the oxygen is conveyed to the surface of the membrane not only by diffusion but additionally by the flow (convection). In that case no oxygen impoverishment occurs at the membrane surface.

A high degree of flow dependence occurs mainly with large cathodes, thin and highly permeable membranes, i.e. where electrode currents are large.

The problem of flow dependence is often solved by prescribing a minimum flow rate.

In METTLER TOLEDO InPro 6900 electrodes, the thin PTFE membrane determining the electrode current (i.e. the actual measuring signal) is separated from the sample solution by a relatively thick silicone membrane. This latter is highly permeable to oxygen molecules and thus acts as an oxygen reservoir. The diffusion of oxygen out of the sample solution into the silicone membrane is spread over a wide area. Since this results in less oxygen being extracted from the sample solution per unit area, the double PTFE/silicone membrane forms an effective buffer against disturbances due to hydrodynamic flow.

This membrane coupled with the guard ring and the special electrolyte ensures excellent signal stability even when the hydrodynamic flow stops (as on a beer filling line for example)

10.7 Oxygen partial pressure – oxygen concentration

The electrode current depends on the oxygen partial pressure and the O₂ permeability of the membrane – but not on the O₂ solubility in the solutions. The oxygen concentration in mg O₂/l (CL) cannot therefore be determined directly with an electrode.

According to Henry's law the oxygen concentration is proportional to its partial pressure (PO₂).

$$CL = pO_2 \cdot a$$

a = solubility factor

If « a » is constant, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

Medium, sat. with air	Solubility at 20 °C (68 °F) and 760 mm Hg
Water	9.2 mg O ₂ /l
4 mol/l KCl	2 mg O ₂ /l
50 % Methanol-water	21.9 mg O ₂ /l

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions.

Thus, determination of the oxygen concentration is only possible with constant and known solubility factors « a ».

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

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- CH Mettler-Toledo (Schweiz) GmbH,** Im Langacher, CH-8606 Greifensee, Switzerland
- D Mettler-Toledo GmbH,** Prozeßanalytik, Ockerweg 3, D-35396 Gießen, Germany
- F Mettler-Toledo Analyse Industrielle S.A.S.,** 30, Boulevard de Douaumont, F-75017 Paris