

# Operation Manual Multi-parameter Transmitter M800





# **Operation Manual Multi-parameter Transmitter M800**



## Content

<b>1</b>	<b>Introduction</b>	<b>9</b>
<b>2</b>	<b>Safety instructions</b>	<b>10</b>
2.1	Definition of equipment and documentation symbols and designations	10
2.2	Correct disposal of the unit	11
<b>3</b>	<b>Unit overview</b>	<b>12</b>
3.1	Overview	12
3.2	Display	13
3.2.1	Start Screen	13
3.2.2	Activation Menu Screen	13
3.3	Graphic Trend Measurement	14
3.3.1	Activation Trend Display Screen	14
3.3.2	Settings for Trend Display Screen	15
3.3.3	Deactivation Trend Display Screen	15
3.4	Control/Navigation	16
3.4.1	Menu Structure	16
3.4.2	Navigation	17
3.4.2.1	Enter the Main Menus	17
3.4.2.2	Navigating the Menu Tree	18
3.4.2.3	Exit a Menu	18
3.4.2.4	Confirm Data and Values	18
3.4.2.5	Return to the Menu Screen	18
3.4.3	Entry of Data	18
3.4.4	Selection Menus	18
3.4.5	"Save changes" Dialog	18
3.4.6	Security Passwords	19
3.4.7	Display	19
<b>4</b>	<b>Installation instruction</b>	<b>20</b>
4.1	Unpacking and inspection of equipment	20
4.1.1	Panel cutout dimensional information	20
4.1.2	Installation procedure	21
4.1.3	Assembly	21
4.1.4	Dimension drawings	22
4.1.5	Pipe mounting	22
4.2	Connection of power supply	23
4.3	Terminal Definition	24
4.3.1	M800 2-channel	24
4.3.2	M800 4-channel	25
4.3.3	TB2 and TB4 – Terminal Assignment for Optical Oxygen Sensor, CO2 hi, and UniCond2e	26
4.3.4	TB2 and TB4 – Terminal Assignment for pH, Amp. Oxygen, Cond 4e, CO2 and O3 Sensors	26
4.3.5	TB3 – Terminal Assignment for Flow Sensors	27
4.4	Connection of Flow Sensor	27
4.4.1	Flow Sensor Input Wiring Kit	27
4.4.2	Kit Contents	28
4.4.3	Flow sensor wiring for Compatible Sensors	28
4.4.4	Wiring for "HIGH" type flow sensors	28
4.4.5	Wiring for "LOW" type flow sensors	31
4.4.6	Wiring for "TYPE 2" flow sensors	31
<b>5</b>	<b>Placing transmitter in, or out, of service</b>	<b>32</b>
5.1	Placing transmitter in service	32
5.2	Placing transmitter out of service	32
<b>6</b>	<b>Guided Setup</b>	<b>33</b>
<b>7</b>	<b>Calibration</b>	<b>34</b>
7.1	Sensor Calibration	34
7.2	Calibration of UniCond2e Sensors	34
7.2.1	Conductivity Calibration of UniCond2e Sensors	34
7.2.1.1	One-Point Calibration	36
7.2.1.2	Two-Point Calibration	37
7.2.1.3	Process Calibration	38
7.2.2	Temperature Calibration of UniCond2e Sensors	39
7.2.2.1	One-Point Calibration	39
7.2.2.2	Two-Point Calibration	40

7.3	Calibration of Cond4e Sensors	42
7.3.1	One-Point Calibration	42
7.3.2	Two-Point Calibration	43
7.3.3	Process Calibration	43
7.4	pH Calibration	44
7.4.1	One-Point Calibration	44
7.4.2	Two-Point Calibration	45
7.4.3	Process Calibration	45
7.5	ORP Calibration of pH Sensors	46
7.6	Calibration of Amperometric Oxygen Sensors	46
7.6.1	One-Point Calibration	47
7.6.2	Process Calibration	47
7.7	Calibration of Optical Oxygen Sensors	48
7.7.1	One-Point Calibration	48
7.7.2	Two-Point Calibration	49
7.7.3	Process Calibration	50
7.8	Calibration of Dissolved Carbon Dioxide Sensors	50
7.8.1	One-Point Calibration	50
7.8.2	Two-Point Calibration	51
7.8.3	Process Calibration	52
7.9	Calibration of Thermal Conductivity CO <sub>2</sub> (CO <sub>2</sub> high) Sensors	52
7.9.1	One-Point Calibration	53
7.9.2	Process Calibration	53
7.10	Calibration of O <sub>3</sub> Sensors	54
7.10.1	One-Point Calibration	54
7.10.2	Process Calibration	55
7.11	Calibration of Flow Sensors	56
7.11.1	One-Point Calibration	56
7.11.2	Two-Point Calibration	57
7.12	Sensor Verification	58
7.13	Edit Calibration Constants for Flow Sensors	59
7.14	UniCond2e Electronics Calibration	59
7.15	Flow Meter Calibration	60
7.16	Flow Meter Verification	61
7.17	Analog Output Calibration	61
7.18	Analog Input Calibration	62
7.19	Maintenance	62
<b>8</b>	<b>Configuration</b>	<b>63</b>
8.1	Measurement	63
8.1.1	Channel Setup	63
8.1.2	Derived Measurements	64
8.1.2.1	% Rejection measurement	65
8.1.2.2	Calculated pH (Power Plant Applications only)	65
8.1.2.3	Calculated CO <sub>2</sub> (Power plant applications only)	66
8.1.3	Display Mode	66
8.1.4	Parameter related Settings	67
8.1.4.1	Conductivity Settings	67
8.1.4.2	pH Settings	68
8.1.4.3	Settings for Oxygen Measurement Based on Amperometric Sensors	69
8.1.4.4	Settings for Oxygen Measurement Based on Optical Sensors	70
8.1.4.5	Dissolved Carbon Dioxide Settings	71
8.1.4.6	Settings for Thermal Conductivity Dissolved CO <sub>2</sub> Measurement (CO <sub>2</sub> hi)	72
8.1.4.7	Settings for TOC Measurement	72
8.1.4.8	Settings for Flow Measurement	73
8.1.4.9	Deionization Capacity (DI-Cap™)	73
8.1.5	Concentration Curve Table	74
8.2	Analog Outputs	75
8.3	Set Points	76
8.4	ISM Setup	77
8.4.1	Sensor Monitor	77
8.4.2	CIP Cycle Limit	78
8.4.3	SIP Cycle Limit	79
8.4.4	AutoClave Cycle Limit	80
8.4.5	DLI Stress Adjustment	81
8.4.6	SAN Cycle Parameters	81
8.4.7	Reset Counters for UniCond2e Sensors	82
8.4.8	Set Calibration Interval for UniCond2e Sensors	82

8.5	General Alarm	83
8.6	ISM / Sensor Alarm	83
8.7	Clean	84
8.8	Display Setup	84
8.9	Digital Inputs	85
8.10	System	85
8.11	PID Controller	86
8.12	Service	90
8.12.1	Set Analog Outputs	90
8.12.2	Read Analog Outputs	90
8.12.3	Read Analog Inputs	90
8.12.4	Set Relay	90
8.12.5	Read Relay	90
8.12.6	Read Digital Inputs	90
8.12.7	Memory	91
8.12.8	Display	91
8.12.9	Calibrate TouchPad	91
8.12.10	Channel Diagnostic	91
8.13	Technical Service	91
8.14	User Management	92
8.15	Reset	92
8.15.1	System Reset	92
8.15.2	Reset Sensor Calibration for Optical DO Sensors	93
8.15.3	Reset Sensor Calibration for UniCond2e Sensors	93
8.15.4	Reset Total Flow	93
8.15.5	Reset for CO2 hi Measurement	94
8.16	RS485 Printer Output Configuration	94
8.17	USB Measurement Interface	95
<b>9</b>	<b>ISM</b>	<b>96</b>
9.1	iMonitor	96
9.2	Messages	97
9.3	ISM Diagnostics	97
9.3.1	pH/ORP, Oxygen, O3 and Cond4e Sensors	98
9.3.2	UniCond2e Sensors	98
9.4	Calibration Data	99
9.4.1	Calibration Data for All ISM Sensors excluding UniCond2e	99
9.4.2	Calibration Data for UniCond2e Sensors	100
9.5	Sensor Info	100
9.6	HW / SW Version	101
9.7	Log Book	101
<b>10</b>	<b>Wizards</b>	<b>102</b>
10.1	Set Wizard	102
10.2	Access to Wizards	102
<b>11</b>	<b>Maintenance</b>	<b>103</b>
11.1	Front panel cleaning	103
<b>12</b>	<b>Troubleshooting</b>	<b>104</b>
12.1	Warning- and Alarm Indication	104
12.1.1	Warning Indication	104
12.1.2	Alarm Indication	105
<b>13</b>	<b>Accessories and Spare Parts</b>	<b>106</b>
<b>14</b>	<b>Specifications</b>	<b>107</b>
14.1	General specifications	107
14.2	Electrical specifications	109
14.3	Mechanical specifications	109
14.4	Environmental specifications	110
<b>15</b>	<b>Warranty</b>	<b>111</b>
<b>16</b>	<b>Buffer tables</b>	<b>112</b>
16.1	Standard pH buffers	112
16.1.1	Mettler-9	112
16.1.2	Mettler-10	113
16.1.3	NIST Technical Buffers	113
16.1.4	NIST standard buffers (DIN and JIS 19266: 2000-01)	114
16.1.5	Hach buffers	114

---

16.1.6	Ciba (94) buffers	115
16.1.7	Merck Titrisole, Riedel-de-Haën Fixanale	115
16.1.8	WTW buffers	116
16.1.9	JIS Z 8802 buffers	116
16.2	Dual membrane pH electrode buffers	117
16.2.1	Mettler-pH/pNa buffers (Na+ 3.9M)	117



# 1 Introduction

Statement of Intended Use – The M800 multiparameter transmitter is a multi-channel online process instrument for measuring various properties of fluids and gases. These include Conductivity, Dissolved Oxygen, O<sub>2</sub> gas, dissolved Ozone, dissolved carbon dioxide, pH/ORP and Flow. The M800 is available in four different versions. The version indicates the amount of measurement parameters which can be covered and the kind of parameter. The version are indicated through there part numbers on the label of the transmitter.

The M800 transmitter is compatible with (digital) ISM and flow sensors.

## M800 parameter fit guide

Version Part no.	Water 2-ch 58 000 802	Water 4-ch 58 000 804	Process 2-ch 52 121 813	Process 4-ch 52 121 853
pH/ORP	•	•	•	•
pH/pNa	–	–	•	•
UniCond2e	•	•	•	•
Cond4e	•	•	•	•
Amp. DO ppm/ppb/trace	•/–/–*	•/–/–*	•/•/•**	•/•/•**
Amp. O <sub>2</sub> gas ppm/ppb/trace	•/–/–*	•/–/–*	•/•/•**	•/•/•**
Optical DO	–	–	•***	•***
Dissolved Carbon Dioxed	–	–	•	•
CO <sub>2</sub> hi (thermal conductivity)	–	–	•***	•***
TOC	•	•	–	–
Dissolved O <sub>3</sub>	•	•	–	–
Flow	•	•	–	–

\* THORNTON sensors

\*\* INGOLD sensors

\*\*\* One (two) optical DO or thermal conductivity CO<sub>2</sub> sensor(s) can be used together with 2-channel (4-channel) transmitter M800.

A colored touch screen conveys measuring data and setup information. The menu structure allows the operator to modify all operational parameters by using the touch screen. A menu-lock-out feature, with password protection, is available to prevent the unauthorized use of the meter. The M800 Multiparameter transmitter can be configured to use up to eight analog and/or up to eight relay outputs for process control.

The M800 Multiparameter transmitter is equipped with a USB communication interface. This interface provides up- and download capabilities of the transmitter configuration via a Personal Computer (PC).

This description corresponds to the firmware release, version 1.2 for the transmitter M800 Water 2-channel, M800 Process 2-channel, M800 Water 4-channel and M800 Process 4-channel. Changes are taking place constantly, without prior notification.

## 2 Safety instructions

This manual includes safety information with the following designations and formats.

### 2.1 Definition of equipment and documentation symbols and designations



**WARNING:** POTENTIAL FOR PERSONAL INJURY.



**CAUTION:** possible instrument damage or malfunction.



**NOTE:** Important operating information.



On the transmitter or in this manual text indicates: Caution and/or other possible hazard including risk of electric shock (refer to accompanying documents)

The following is a list of general safety instructions and warnings. Failure to adhere to these instructions can result in damage to the equipment and/or personal injury to the operator.

- The M800 Transmitter should be installed and operated only by personnel familiar with the transmitter and who are qualified for such work.
- The M800 Transmitter must only be operated under the specified operating conditions (see chapter 14 “Specifications”).
- Repair of the M800 Transmitter must be performed by authorized, trained personnel only.
- With the exception of routine maintenance, cleaning procedures or fuse replacement, as described in this manual, the M800 Transmitter must not be tampered with or altered in any manner.
- Mettler-Toledo accepts no responsibility for damage caused by unauthorized modifications to the transmitter.
- Follow all warnings, cautions, and instructions indicated on and supplied with this product.
- Install equipment as specified in this instruction manual. Follow appropriate local and national codes.
- Protective covers must be in place at all times during normal operation.
- If this equipment is used in a manner not specified by the manufacturer, the protection provided by it against hazards may be impaired.



**WARNINGS:**

- Installation of cable connections and servicing of this product require access to shock hazard voltage levels.
- Main power and relay contacts wired to separate power source must be disconnected before servicing.
- Switch or circuit breaker shall be in close proximity to the equipment and within easy reach of the OPERATOR; it shall be marked as the disconnecting device for the equipment.
- Main power must employ a switch or circuit breaker as the disconnecting device for the equipment.
- Electrical installation must be in accordance with the National Electrical Code and/or any other applicable national or local codes.

**NOTE: RELAY CONTROL ACTION**

the M800 transmitter relays will always de-energize on loss of power, equivalent to normal state, regardless of relay state setting for powered operation. Configure any control system using these relays with fail-safe logic accordingly.

**NOTE: PROCESS UPSETS**

Because process and safety conditions may depend on consistent operation of this transmitter, provide appropriate means to maintain operation during sensor cleaning, replacement or sensor or instrument calibration.



**NOTE:** This is a 4-wire-product with an active 4–20 mA analog output.

Please do not supply to terminal 3 to 10 of TB1 and terminal 1 to 8 of TB3.


## 2.2 Correct disposal of the unit

When the transmitter is finally removed from service, observe all local environmental regulations for proper disposal.

### 3 Unit overview

The M800 models are available in 1/2DIN case size. The M800 models provide an integral IP66 housing for wall- or pipe mount.


#### 3.1 Overview



- 1: Hard Polycarbonate case
- 2: VGA Screen
- 3: Power Supply Terminals
- 4: Relay Output Terminals
- 5: Analog Output/Digital Input Terminals
- 6: Sensor Input Terminals


## 3.2 Display

### 3.2.1 Start Screen

After starting the M800, the following Start Screen (logout screen) is shown automatically. To return from the Menu Screen to the Start Screen press . The M800 will return automatically after 240 seconds from the Menu Screen or any configuration screen to the Start Screen if the user has not pressed the touch screen.



### 3.2.2 Activation Menu Screen

While the M800 shows the Start Screen (logout screen) touch the display to activate the Menu Screen. To return to the Menu Screen from other menus press .




### 3.3 Graphic Trend Measurement

Any single measurement may be displayed as a trend measurement over time. Measurement values will be indicated by a value on the Y axis and time elapsed on the X axis of the graph displayed. An actual measurement for the selected value will also display numerically above the graphic trend display. The measurement value is refreshed once per second.

Graphic trending will only display the data within maximum/minimum range. Out of range values or invalid values will not be displayed. The Y axis will display the maximum value unit with its range; X axis unit uses "mins" for minutes for measurements less than one hour and "hrs" for one day. 4 scales for X/Y axis. The maximum value on Y-axis is one decimal place.

#### 3.3.1 Activation Trend Display Screen

While the M800 displays the Menu Screen, touch any measurement value line of the display screen twice (1-chan, 2-chan, 4-meas, 8-meas) to activate the trend display for that measurement.



If a sensor is disconnect/connect pop-up window come up, after close it will go back to the Menu Screen.

Red/yellow bar on top line will display for any message occurring during trending. 'H', 'P', 'AB' will display when this channel is in hold or process.

### 3.3.2 Settings for Trend Display Screen

For setting configurations, touch any area of the graphic trend display to go to the pop-up window of this measurement parameter. Settings are at the default values. However, these settings may be changed when options are available, as needed.



**Time:** option button. For graphic display time (X axis)  
 1-h (default value)  
 1-day

**NOTE:** 1 h means: 1 meas storage/15 seconds, totally 240 measurements for 1h. 1 day means: 1 meas storage/6 minutes, totally 240 measurements for 1 day;

**Range:** option button  
 Default(default value)  
 Individual

When "Default" modes are set for the maximum or minimum value, this indicates the full measurement range for this unit. A Max or Min button is not displayed. If setting is selectable, the user can set maximum and minimum settings manually.

**Max:** Edit button.  
 Maximum value of this unit on Y axis. xxxxxx, floating decimal point.

**Min:** Edit button.  
 Minimum value of this unit on Y axis. xxxxxx, floating decimal point.  
 Max Value > Min Value



**NOTE:** Settings for y- and x-axis and the corresponding measurement values are stored the transmitters memory. A power down returns to default settings.

### 3.3.3 Deactivation Trend Display Screen

Press  in activated graphic trend screen to return to Menu Screen.




**NOTE:** If a sensor is disconnect/connect pop-up window come up, after close it will go back to the Menu Screen.

### 3.4 Control/Navigation

#### 3.4.1 Menu Structure

Below is the structure of the M800 menu tree:





## 3.4.2 Navigation

### 3.4.2.1 Enter the Main Menus

Activate the Menu Screen and press one of icons to enter the different main menus:



ISM menu



Wizard menu



Calibration menu



Configuration menu

### 3.4.2.2 Navigating the Menu Tree

To browse through the menus, press the open arrows > and/or <. To access a menu touch the corresponding arrow ► in the same line.

### 3.4.2.3 Exit a Menu

Press ⏪ to exit the menu. Press 🏠 to return to the Menu Screen (see chapter 3.2.2 "Activation Menu Screen").

### 3.4.2.4 Confirm Data and Values

Use the ↵ key to confirm values. Press ESC and the values will not be taken over.

### 3.4.2.5 Return to the Menu Screen

Press 🏠 to return to the Menu Screen (see chapter 3.2.2 "Activation Menu Screen"). To return from the Menu Screen to the Start Screen press 🔒.

## 3.4.3 Entry of Data

The M800 displays a keypad for modifying values. Press the ↵ button and the transmitter will take over the value. Press the ESC button to exit the keypad without changing data.



**NOTE:** For some value the unit can be modified. In this case the keypad shows a button with a u. To select another unit for the entered value on the keypad press the u button. To return again press the 0–9 button.



**NOTE:** For some entries letters and/or numbers can be used. In that case the keypad shows a button 'A,a,0'. Press this button to change between capital letters, small letters and numbers on the keypad.

## 3.4.4 Selection Menus

Some menus require a selection of a parameter / data. In this case the transmitter displays a pop up window. Press the according field to select the value. The pop window will be closed and the selection will be taken over.

## 3.4.5 "Save changes" Dialog

If the M800 brings up the "Save changes" dialog there are the following options. No will discard the entered values, Yes will save changes made and Cancel will bring you back to continue configuring.

### 3.4.6 Security Passwords

The M800 transmitter allows a security lock-out of various menus. If the security lock-out feature of the transmitter has been enabled, a security password must be entered to allow access to the menu. See chapter 8.14 "User Management".

### 3.4.7 Display



**NOTE:** In the event of an alarm or other error condition the M800 Transmitter will display a flashing bar graph on the display. This bar graph will remain until the condition that caused it has been cleared (see chapter 12.1 "Warning- and Alarm Indication").



**NOTE:** During calibrations, clean, Digital In with Analog Output/Relay/USB in Hold state, a flashing "H" (Hold) will appear in the upper right corner of the display for the corresponding channel. This symbol will remain for 20 sec., after end of calibration. This symbol will remain for 20 seconds until after the calibration or clean is completed. This symbol will also disappear when Digital In is deactivated.

## 4 Installation instruction

### 4.1 Unpacking and inspection of equipment


Inspect the shipping container. If it is damaged, contact the shipper immediately for instructions. Do not discard the box.

If there is no apparent damage, unpack the container. Be sure all items shown on the packing list are present.

If items are missing, notify Mettler-Toledo immediately

#### 4.1.1 Panel cutout dimensional information

Below are cut-out dimensions required by the transmitter when mounted within a flat panel or on a flat enclosure door. This surface must be flat and smooth. Textured or rough surfaces are not recommended and may limit the effectiveness of the gasket seal provided.



Optional hardware accessories are available that allow for panel- or pipe-mount.

Refer to chapter 13 "Accessories and Spare Parts" for ordering information.

## 4.1.2 Installation procedure


### General:

- Orient the transmitter so that the cable grips face downward.
- Wiring routed through the cable grips shall be suitable for use in wet locations.
- In order provide IP66 enclosure ratings, all cable glands must be in place. Each cable gland must be filled using a cable, or suitable Cable Gland Hole Seal.

### For Pipe Mount:


- Use only manufacturer-supplied components for pipe-mounting the M800 transmitter and install per the supplied instructions. See chapter 13 "Accessories and Spare Parts" for ordering information.

## 4.1.3 Assembly




- 1: 3 M20x1.5 cable glands  
2: 2 M25x1.5 cable glands  
3: 4 screws

### 4.1.4 Dimension drawings



### 4.1.5 Pipe mounting



## 4.2 Connection of power supply

All connections to the transmitter are made on the inside of all models.



Be sure power to all wires is turned off before proceeding with the installation

A three-terminal connector on TB6 of all M800 models is provided for power connection. All M800 models are designed to operate from a 20–30 VDC or a 100 to 240 VAC power source. Refer to specifications for power requirements and ratings and size power wiring accordingly (16–24 AWG, wire cross-section between 0.2 mm<sup>2</sup> and 1.5 mm<sup>2</sup>).

The terminal block for power connections is labeled TB6 on the rear panel of the transmitter. One terminal is labeled **N** (–) for the Neutral wire and the other **L** (+) for the Line (or Load) wire. For DC power, use the polarity shown in parentheses.

The terminals are suitable for single wires and flexible leads with a wire cross-section from 0.2 mm<sup>2</sup> up to 1.5 mm<sup>2</sup>, (16–24 AWG).


### 4.3 Terminal Definition

#### 4.3.1 M800 2-channel

Power connections:

**N** for Neutral and **L** for Line, for 100 to 240 VAC or 20–30 VDC.

Terminal number	TB1	TB2 (ISM Ch1,2)	TB3	TB4	TB5	TB6
1	DI1+	DI2+	Aout5+	<b>Not installed</b>	AI1+	L
2	DI1–	DI2–	Aout5–		AI1–	N
3	Aout1+	1-Wire_Ch1	Aout6+		DI4+	Ground
4	Aout1–	GND5V_Ch1	Aout6–		DI4–	Relay1_NC
5	Aout2+	RS485B_Ch1	Aout7+		DI5+	Relay1_COM
6	Aout2–	RS485A_Ch1	Aout7–		DI5–	Relay2_NO
7	Aout3+	GND5V_Ch1	Aout8+		DI6+	Relay2_COM
8	Aout3–	5V_Ch1	Aout8–		DI6–	Relay3_NO
9	Aout4+	24V_Ch2	Ain_Ch5		Relay5_NO	Relay3_COM
10	Aout4–	GND24V_Ch2	AJ_Ch5		Relay5_COM	Relay4_NO
11	n. a.	1-Wire_Ch2	5V_Ch5		Relay6_NO	Relay4_COM
12	n. a.	GND5V_Ch2	GND5V_Ch5		Relay6_COM	n. a.
13	n. a.	RS485B_Ch2	Bin_Ch6		Relay7_NO	n. a.
14	n. a.	RS485A_Ch2	BJ_Ch6		Relay7_COM	n. a.
15	n. a.	GND5V_Ch2	5V_Ch6		Relay8_NC	n. a.
16	n. a.	5V_Ch2	GND5V_Ch6		Relay8_COM	n. a.



NO: normally open (contact open if un-actuated).  
 NC: normally closed (contact closed if un-actuated).  
 Ain:  
 Bin:  
 n.a. not available

AO: Analog Output  
 DI: Digital Input  
 AJ:  
 BJ:



**NOTE:** This is a 4-wire-product with an active 4–20 mA analog output.  
 Please do not supply to terminal no. 3 to 10 of TB1 and 1 to 8 of TB3.




### 4.3.2 M800 4-channel

Power connections:

**N** for Neutral and **L** for Line, for 100 to 240 VAC or 20–30 VDC.

Terminal number	TB1	TB2 (ISM Ch1,2)	TB3	TB4 (ISM Ch3,4)	TB5	TB6
1	DI1+	DI2+	Aout5+	DI3+	AI1+	L
2	DI1–	DI2–	Aout5–	DI3–	AI1–	N
3	Aout1+	1-Wire_Ch1	Aout6+	1-Wire_Ch3	DI4+	Ground
4	Aout1–	GND5V_Ch1	Aout6–	GND5V_Ch3	DI4–	Relay1_NC
5	Aout2+	RS485B_Ch1	Aout7+	RS485B_Ch3	DI5+	Relay1_COM
6	Aout2–	RS485A_Ch1	Aout7–	RS485A_Ch3	DI5–	Relay2_NO
7	Aout3+	GND5V_Ch1	Aout8+	GND5V_Ch3	DI6+	Relay2_COM
8	Aout3–	5V_Ch1	Aout8–	5V_Ch3	DI6–	Relay3_NO
9	Aout4+	24V_Ch2	Ain_Ch5	24V_Ch4	Relay5_NO	Relay3_COM
10	Aout4–	GND24V_Ch2	AJ_Ch5	GND24V_Ch4	Relay5_COM	Relay4_NO
11	n. a.	1-Wire_Ch2	5V_Ch5	1-Wire_Ch4	Relay6_NO	Relay4_COM
12	n. a.	GND5V_Ch2	GND5V_Ch5	GND5V_Ch4	Relay6_COM	n. a.
13	n. a.	RS485B_Ch2	Bin_Ch6	RS485B_Ch4	Relay7_NO	n. a.
14	n. a.	RS485A_Ch2	BJ_Ch6	RS485A_Ch4	Relay7_COM	n. a.
15	n. a.	GND5V_Ch2	5V_Ch6	GND5V_Ch4	Relay8_NC	n. a.
16	n. a.	5V_Ch2	GND5V_Ch6	5V_Ch4	Relay8_COM	n. a.



NO: normally open (contact open if un-actuated).

NC: normally closed (contact closed if un-actuated).

Ain:

Bin:

n.a. not available

AO: Analog Output

DI: Digital Input

AJ:

BJ:



**NOTE:** This is a 4-wire-product with an active 4–20 mA analog output. Please do not supply to terminal no. 3 to 10 of TB1 and 1 to 8 of TB3.

### 4.3.3 TB2 and TB4 – Terminal Assignment for Optical Oxygen Sensor, CO2 hi, and UniCond2e

	<b>TB2 (ISM Ch1, 2)</b>	<b>TB4 (ISM Ch3,4)</b>	<b>Optical Oxygen*, CO2 hi*,</b>	<b>UniCond2e**</b>
<b>Terminal no.</b>	<b>Function</b>	<b>Function</b>	<b>Sensor wire color</b>	<b>Sensor wire color</b>
1	DI2+	DI6+	–	–
2	DI2–	DI6–	–	–
3	1-Wire_Ch1	1-Wire_Ch3	–	–
4	GND5V_Ch1	GND5V_Ch3	–	–
5	RS485B_Ch1	RS485B_Ch3	–	black
6	RS485A_Ch1	RS485A_Ch3	–	red
7	GND5V_Ch1	GND5V_Ch3	–	white
8	5V_Ch1	5V_Ch3	–	blue
9	24V_Ch2	24V_Ch4	brown	–
10	GND24V_Ch2	GND24V_Ch4	black	–
11	1-Wire_Ch2	1-Wire_Ch4	–	–
12	GND5V_Ch2	GND5V_Ch4	grey	–
13	RS485B_Ch2	RS485B_Ch4	blue	black
14	RS485A_Ch2	RS485A_Ch4	white	red
15	GND5V_Ch2	GND5V_Ch4	yellow	white
16	5V_Ch2	5V_Ch4	–	blue

\* Always one O2 optical or thermal conductivity CO2 sensor can be connected to plug TB2 and TB4.

\*\* Transparent wire not connected

### 4.3.4 TB2 and TB4 – Terminal Assignment for pH, Amp. Oxygen, Cond 4e, CO2 and O3 Sensors

	<b>TB2 (ISM Ch1, 2)</b>	<b>TB4 (ISM Ch3,4)</b>	<b>pH, amp. Oxygen, Cond 4e, CO2, and O3</b>
<b>Terminal no.</b>	<b>Function</b>	<b>Function</b>	<b>Sensor wire color</b>
1	DI2+	DI6+	–
2	DI2–	DI6–	–
3	1-Wire_Ch1	1-Wire_Ch3	transparent (cable core)
4	GND5V_Ch1	GND5V_Ch3	red
5	RS485B_Ch1	RS485B_Ch3	–
6	RS485A_Ch1	RS485A_Ch3	–
7	GND5V_Ch1	GND5V_Ch3	–
8	5V_Ch1	5V_Ch3	–
9	24V	24V	–
10	GND24V	GND24V	–
11	1-Wire_Ch2	1-Wire_Ch4	transparent (cable core)
12	GND5V_Ch2	GND5V_Ch4	red
13	RS485B_Ch2	RS485B_Ch4	–
14	RS485A_Ch2	RS485A_Ch4	–
15	GND5V_Ch2	GND5V_Ch4	–
16	5V_Ch2	5V_Ch4	–

### 4.3.5 TB3 – Terminal Assignment for Flow Sensors

	<b>TB3</b>	<b>Flow hi, Flow lo, Flow Type2</b>
<b>Terminal no.</b>	<b>Transmitter</b>	<b>Function</b>
1	Aout5+	–
2	Aout5-	–
3	Aout6+	–
4	Aout6-	–
5	Aout7+	–
6	Aout7-	–
7	Aout8+	–
8	Aout8-	–
9	Ain_Ch5	Flow Pulse Input
10	AJ_Ch5	+ 10 VDC
11	5V_Ch5	+ 5 VDC
12	GND5V_Ch5	Ground
13	Ain_Ch6	Flow Pulse Input
14	AJ_Ch6	+ 10 VDC
15	5V_Ch6	+ 5 VDC
16	GND5V_Ch6	Ground

## 4.4 Connection of Flow Sensor

The M800 transmitter is designed to operate with various types of sensors. These sensors require different wiring configurations. Listed below are instructions for wiring the various types of sensors offered by Mettler-Toledo THORNTON for use with this transmitter. Please consult the factory for assistance if attempting to wire sensors not offered by Mettler-Toledo Thornton as some sensors may not be compatible.

### 4.4.1 Flow Sensor Input Wiring Kit

This kit contains components that may be needed at input terminals to condition sensor signals. Refer to the following sections or to the instruction manual for wiring details.

## 4.4.2 Kit Contents

This kit contains the following items:

- 4x Wire nuts
- 4x 10K ohm resistors for use with Burkert 8020 and 8030 type sensors, and GF Signet 2500-series sensors.
- 4x 1K ohm resistors for use with Data Industrial 200-series and Fluidyne insertion type sensors.
- 4x 0.33uF, 50 V capacitors for use with Burkert 8020 and 8030 type sensors, Data Industrial 200-series and 4000-series sensors, GF Signet 2500-series sensors, Sanitary Turbine-Type sensors, Fluidyne insertion type sensors and Racine Federated (Formerly Asahi/America) vortex-style sensors.

## 4.4.3 Flow sensor wiring for Compatible Sensors


The following sections provide wiring information to properly connect various compatible flow sensors to the M800 transmitter. When using the Configuration menu of the transmitter to setup the flow sensor, the first prompt asks to select the TYPE of flow sensor being connected.

There are three choices as follows:

- High: All flow sensors described in Section 4.4.4  
 Low: P515 Signet flow sensors only, described in section 4.4.5  
 Type 2: Asahi flow sensors, described in Section 4.4.6

## 4.4.4 Wiring for "HIGH" type flow sensors


The following wiring information is used when connecting (Burkert 8020 and 8030 type) inline Hall effect 5VDC, flow sensors. **Thornton models 33901 thru 33935.**



Extension cable not provided. Use 2-conductor twisted pair with shield, 22 AWG (Belden 8451 or equivalent), 1,000 ft (305 m) maximum length.

The following wiring information is used when connecting Badger (formerly Data Industrial 200-Series) forward-swept paddlewheel type flow sensors.


**Thornton models 33142 thru 33145 and 33159 thru 33162 and 33273.**



Extension cable provided with sensor. Use 2-conductor twisted pair with shield 20 AWG (Belden 9320 or equivalent) to extend length to 2000 ft (610 m) max.


The following wiring information is used when connecting Badger (formerly Data Industrial 4000-Series) forward-swept paddlewheel type flow sensors.

**Thornton models 33174 thru 33177 and 33171 and 33172.**




20 ft (6.1 m) extension cable provided with sensor. Use 3-conductor with shield, 20 AWG (Belden 9364 or equivalent) to extend length to 2000 ft (610 m) maximum.

The following wiring information is used when connecting (GF Signet 2500-Series) Hall Effect paddlewheel type flow sensors. **Thornton models 33282, 33285, 33287, 33298 thru 33305.**




25 ft (7.6 m) extension cable provided with sensor. Use 2-conductor with shield, 22 AWG (Belden 8451 or equivalent) to extend length to 1000 ft (305 m) maximum.

The following wiring information is used when connecting Sanitary Turbine type flow sensors. **Thornton models 33336 thru 33377 (Hoffer) and 33441 thru 33450 (Sponsler).**




20 ft (6.1 m) extension cable provided with sensor. Use 3-conductor with shield, 20 AWG (Belden 9364 or equivalent) to extend length to 3000 ft (915 m) maximum.



20 ft (6.1 m) extension cable provided with sensor. Use 3-conductor with shield, 20 AWG (Belden 9364 or equivalent) to extend length to 3000 ft (915 m) maximum.


The following wiring information is used when connecting Spirax Sarco/Emco flow (formerly Fluidyne) insertion type flow sensors. **Thornton models 33358 thru 33375.**



Extension cable not provided. Use 2-conductor twisted pair with shield, 20 AWG (Belden 9320 or equivalent), 2000 ft (610 m) maximum length.

### 4.4.5 Wiring for "LOW" type flow sensors


The following wiring information is used when connecting (GF Signet 515) type flow sensors. **Thornton models 33189, 33193, 33195, 33196, and 33229.**



Extension cable not provided. Use 2-conductor twisted pair with shield, 22 AWG (Belden 8451 or equivalent), 200 ft (61 m) maximum length.

### 4.4.6 Wiring for "TYPE 2" flow sensors

The following wiring information is used when connecting Racine Federated (formerly Asahi/America) vortex flow sensors. **Thornton models 33308 to 33335.**



Extension cable not provided. Use 3-conductor with shield, 20 AWG (Belden 9364 or equivalent), 1000 ft (305 m) maximum length.

## **5 Placing transmitter in, or out, of service**

### **5.1 Placing transmitter in service**



After connecting the transmitter to power supply circuit, it will be active as soon as the circuit is powered.

### **5.2 Placing transmitter out of service**

First disconnect the unit from the main power source, then disconnect all remaining electrical connections. Remove the unit from the panel. Use the installation instruction in this manual as reference for dis-assembling mounting hardware.

All transmitter settings stored in memory are non volatile.




## 6 Guided Setup

PATH:  \ CONFIG \ Guided Setup



**NOTE:** Please do not use Guided Setup menu after configuration of the transmitter, because some of the settings i.e. analog output configuration will may be set to default values again.

See the following explanation to get more details about the different settings for the guided setup.



Select the desired **Channel** for the guided setup and in the same line the parameter.

If Auto is selected, M800 transmitter automatically recognizes the type of sensor. The channel can also be fixed to a certain measurement parameter (parameter = pH/ORP, UniCond2e, Cond4e, O2 hi, O2 lo, O2 Trace, O2 opt, O3 and Flow hi, Flow lo, Flow Type 2), depending on the type of transmitter. For detailed information refer to chapter 8.1.1 "Channel Setup"

Press the corresponding button to measurement **M1** to configure the measurement. For detailed information about the configuration options refer to chapter 8.1.1 "Channel Setup".



**NOTE:** If the guided setup for a flow sensor has been selected, the calibration factor of the sensor from the sensor label or certificate can be entered. Press therefore the Cal Factor button. For sensor types High and Low slope and offset can be entered. In case of sensor Type 2 the slope followed by a table of K and F values can be entered.

Assign the corresponding output signal **Aout'X'** to the measurement through pressing Yes. For detailed information about the configuration of the analog output signal refer to chapter 8.2 "Analog Outputs"

Enter the **Min Value**, that corresponds with start point of the analog output range.

Enter the **Max Value**, that corresponds with end point of the analog output signal.

Additional settings can be done by navigating to the next page of the menu.



Assigning the corresponding **Set Point'X'** to the measurement through pressing Yes. For detailed information about the configuration of the set point refer to chapter 8.3 "Set Points"

Select the **Type** for the setpoint.



The type of the setpoint can be High, Low, Between, Outside or Off. An "Outside" setpoint will cause an alarm condition whenever the measurement goes above its high limit or below its low limit. A "Between" setpoint will cause an alarm condition to occur whenever the measurement is between its high and low limits.



**NOTE:** If the type of set point is not Off additional settings can be done. See the following description.


According to the selected type of set point, value(s) according to the limit(s) can be entered.

Select the desired relay that will be activated if the defined conditions are reached through the parameter **SP Relay**. If the chosen relay is used for another task, the transmitter shows a message on the screen that there is a Relay Conflict.

To escape the menu of the settings for Guided Setup press . To return to the Menu Screen (see chapter 3.1 "Display") press . The M800 will bring up the Save Changes dialog.

# 7 Calibration

PATH: \ Cal



**NOTE:** During calibration, the outputs for the corresponding channel will default to be held at their current values until 20 seconds after the calibration menu is exited. A flashing H appears in the upper right corner of the display while outputs are held. Refer to chapter 8.2 "Analog Outputs" and chapter 8.3 "Set Points" to change the hold output status.

## 7.1 Sensor Calibration

PATH: \ Cal \ Calibrate Sensor



Select the desired channel (**Chan**) for calibration.

**NOTE:** During sensor calibration, the outputs will default to be held at their current values until 20 seconds after the calibration menu is exited. A flashing H appears in the upper right corner of the display while outputs are held. Refer to chapter 8.2 "Analog Outputs" and chapter 8.3 "Set Points" to change the hold output status.

See the following explanation to get more details about the calibration options and procedure.

## 7.2 Calibration of UniCond2e Sensors

### 7.2.1 Conductivity Calibration of UniCond2e Sensors

The M800 provides the ability to perform a one-point, two-point or process conductivity or resistivity calibration for 2e- sensors.




**NOTE:** When performing calibration on a conductivity sensor, results will vary depending on the method, calibration apparatus and/or quality of reference standards used to perform the calibration.



**NOTE:** For measuring tasks the temperature compensation for the application as defined through the parameter settings for conductivity will be considered and not the temperature compensation selected through the calibration procedure (see also chapter 8.1.4.1 "Conductivity Settings"; PATH:  $\mathcal{H}$ \CONFIG\Meas\Parameter Setting).

Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration"; PATH:  $\mathcal{H}$ \Cal\Calibrate Sensor) and choose the desired channel for calibration.



The following menus can be called up:

**Unit:** Choose between the units for conductivity (S/cm) and resistivity ( $\Omega$ -cm).

**Method:** Select the desired calibration procedure. Available are 1-point, 2-point or process calibration.

**Options:** The desired compensation mode for the calibration process can be selected.

Choices are "None", "Standard", "Light 84", "Std 75 °C", "Linear 25°C", "Linear 20°C", "Glycol.5", "Glycol1", "Cation", "Alcohol" and "Ammonia".

None does not make any compensation of the measured conductivity value. The uncompensated value will be displayed and proceeded.

Standard compensation includes compensation for non-linear high purity effects as well as conventional neutral salt impurities and conforms to ASTM standards D1125 and D5391.

Light 84 compensation matches the high purity water research results of Dr. T.S. Light published in 1984. Use only if your institution has standardized on that work.

Std 75 °C compensation is the Standard compensation algorithm referenced to 75 °C. This compensation may be preferred when measuring Ultrapure Water at an elevated temperature (Resistivity of ultrapure water compensated to 75 °C is 2.4818 Mohm-cm.)

Linear 25 °C compensation adjusts the reading by a coefficient or factor expressed as %/°C (deviation from 25 °C). Use only if the solution has a well-characterized linear temperature coefficient. The factory default setting is 2.0%/°C.

Linear 20 °C compensation adjusts the reading by a coefficient or factor expressed as %/°C (deviation from 20 °C). Use only if the solution has a well-characterized linear temperature coefficient. The factory default setting is 2.0%/°C.

Glycol.5 compensation matches the temperature characteristics of 50% ethylene glycol in water. Compensated measurements using this solution may go above 18 Mohm-cm.

Glycol1 compensation matches the temperature characteristics of 100% ethylene glycol. Compensated measurements may go well above 18 Mohm-cm.

Cation compensation is used in power industry applications measuring the sample after a cation exchanger. It takes into account the effects of temperature on the dissociation of pure water in the presence of acids.

Alcohol compensation provides for the temperature characteristics of a 75% solution of isopropyl alcohol in pure water. Compensated measurements using this solution may go above 18 Mohm-cm.

Ammonia compensation is used in power industry applications for specific conductivity measured on samples using ammonia and/or ETA (ethanolamine) water treatment. It takes into account the effects of temperature on the dissociation of pure water in the presence of these bases.




**NOTE:** If compensation mode "Linear 25 °C" or "Linear 20 °C" has been chosen, the coefficient for the adjustment of the reading can be modified. In this case an additional input field will be displayed.


The changes are valid until the calibration mode has been escaped. After the values defined in the configuration menu are valid again.

### 7.2.1.1 One-Point Calibration


Select calibration procedure 1-Point (see chapter 7.2.1 "Conductivity Calibration of UniCond2e Sensors"). With 2e-sensors a one point calibration is always performed as a slope calibration.



Press the button Cal for starting calibration.



Place the electrode in the reference solution and press Next button.




The second value displayed on the screen is the value being measured by the transmitter and sensor in units selected by the user.

Press the input field for **Point1** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button and the transmitter will take over the value.




**NOTE:** To select another unit for the entered value on the keypad press the u button. To return again press the 0-9 button.



The screen shows the entered value for the reference solution (1st line) and the measured value of the M800 (2nd line).


Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (press button SaveCal) or discarded (press button Cancel).

Use the Back button to go one step back in the calibration procedure.




If "SaveCal" is chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor". After pressing the Done button the M800 returns to the calibration menu for the sensor.

### 7.2.1.2 Two-Point Calibration

Select calibration procedure 2-Point. With 2e-sensors a two point calibration is always performed as an offset and slope calibration.

Press the button Cal for starting calibration.




Place the electrode in the first reference solution and press Next button.

**CAUTION:** Rinse sensors with a high-purity water solution between calibration points to prevent contamination of the reference solutions.



The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.


Press the input field for **Point1** to enter the calibration point. The M800 displays a keypad for modifying the value. Press the ← button and the transmitter will take over the value.




**NOTE:** To select another unit for the entered value on the keypad press the u button. To return again press the 0-9 button.

The screen shows the entered value for the first reference solution (1st line) and the measured value of the M800 (2nd line).

Press the Next button to go on with the calibration.



Place the electrode in the second reference solution and press Next button.






The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.

Press the input field for **Point2** to enter the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button and the transmitter will take over the value.




**NOTE:** To select another unit for the entered value on the keypad press the u button. To return again press the 0–9 button.



The screen shows the entered value for the second reference solution (1st line) and the measured value of the M800 (2nd line).

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (press button SaveCal) or discarded (press button Cancel).

Use the Back button to go one step back in the calibration procedure.




If "SaveCal" is chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor". After pressing the Done button the M800 returns to the calibration menu for the sensor.

### 7.2.1.3 Process Calibration

Select calibration procedure Process (see chapter 7.2.1 "Conductivity Calibration of UniCond2e Sensors"). With 2e-sensors a process calibration is always performed as a slope calibration.



Press the button Cal for starting calibration.




Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu screen if the related channel is selected in the display.



After determining the conductivity value of the sample, press the calibration icon in the Menu Screen (see chapter .... Enter Calibration Mode) again.




Press the input field for **Point1** and enter the conductivity value of the sample. Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (press button SaveCal) or discarded (press button Cancel).

Through the Back button it is possible to go one step back in the calibration procedure.



If "SaveCal" is chosen, the message "Calibration Saved Successfully!" is displayed. After pressing the Done button the M800 returns to the Menu Screen.

## 7.2.2 Temperature Calibration of UniCond2e Sensors

The M800 provides the ability to perform a one-point or two-point calibration for the temperature sensor of the UniCond2e.

Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration"; PATH: \ Cal \ Calibrate Sensor) and choose the desired channel for calibration.




The following menus can be called up:

**Unit:** Choose between the units °C and °F.

**Method:** Select the desired calibration procedure. Available are 1-point and 2-point calibration.

### 7.2.2.1 One-Point Calibration

Select calibration procedure 1-Point. With UniCond2e-sensors a one point temperature calibration can be performed as a slope or offset calibration.



Press the right input field for the parameter **Method**. Choose Slope or Offset calibration through pressing the corresponding field.



Press the button Cal for starting calibration.




Place the electrode in the reference solution and press Next button.




The second value displayed on the screen is the value being measured by the transmitter and sensor.

Press the input field for **Point1** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button and the transmitter will take over the value.



The screen shows the entered value for the reference solution (1st line) and the measured value of the M800 (2nd line).

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (press button SaveCal) or discarded (press button Cancel).

Use the Back button to go one step back in the calibration procedure.




If "SaveCal" is chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor". After pressing the Done button the M800 returns to the calibration menu for the sensor.

## 7.2.2.2 Two-Point Calibration

Select calibration procedure 2-Point (see chapter 7.2.2 "Temperature Calibration of UniCond2e Sensors"). With 2e-sensors a two point calibration is always performed as an offset and slope calibration.




Press the button Cal for starting calibration.




Place the electrode in the first reference solution and press Next button.






The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.

Press the input field for **Point1** to enter the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button and the transmitter will take over the value.




The screen shows the entered value for the first reference solution (1st line) and the measured value of the M800 (2nd line).

Press the Next button to go on with the calibration.



Place the electrode in the second reference solution and press Next button.



The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.

Press the input field for **Point2** to enter the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button and the transmitter will take over the value.



The screen shows the entered value for the second reference solution (1st line) and the measured value of the M800 (2nd line).

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (press button SaveCal) or discarded (press button Cancel).

Use the Back button to go one step back in the calibration procedure.



If "SaveCal" is chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor". After pressing the Done button the M800 returns to the calibration menu for the sensor.

## 7.3 Calibration of Cond4e Sensors

PATH:  \ Cal \ Calibrate Sensor


The M800 provides the ability to perform a one-point, two-point or process conductivity or resistivity calibration for 4e sensors.



**NOTE:** When performing calibration on a conductivity sensor, results will vary depending on the method, calibration apparatus and/or quality of reference standards used to perform the calibration.



**NOTE:** For measuring tasks the temperature compensation for the application as defined through the parameter settings for conductivity will be considered and not the temperature compensation selected through the calibration procedure (see also chapter 8.1.4.1 "Conductivity Settings").



The following menus can be called up:

**Unit:** Between the units for conductivity and resistivity can be chosen.

**Method:** Select the desired calibration procedure, 1-point, 2-point or process calibration.

**Options:** Select the desired temperature compensation mode for the calibration process.




**NOTE:** If compensation mode "Linear 25 °C" or "Linear 20 °C" has been chosen, the coefficient for the adjustment of the reading can be modified.

The changes are valid until the calibration mode has been exited. After the values defined in the configuration menu are valid again.

### 7.3.1 One-Point Calibration

With 4e-sensors a one point calibration is always performed as a slope calibration.



Press the button Cal for starting calibration.

Place the electrode in the reference solution and press Next button.

Enter the value for the calibration point (**Point1**).

Press the Next button to start the calculation of the calibration results.


The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

### 7.3.2 Two-Point Calibration

With 4e-sensors a two point calibration is always performed as an offset and slope calibration.



Press the button Cal for starting calibration.

Place the electrode in the first reference solution and press Next button.

**CAUTION:** Rinse sensors with a high-purity water solution between calibration points to prevent contamination of the reference solutions.

Enter the value for the first calibration point (**Point1**).

Press the Next button to go on with the calibration.

Place the electrode in the second reference solution and press Next button.

Enter the value for the second calibration point (**Point2**).

Press the Next button to start the calculation of the calibration results.


The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

### 7.3.3 Process Calibration

With 4e-sensors a process calibration is always performed as a slope calibration.



Press the button Cal for starting calibration.

Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu screen if the related channel is selected in the display.

After determining the conductivity value of the sample, press the calibration icon in the Menu Screen again.

Enter the conductivity value of the sample. Press the Next button to start the calculation of the calibration results.

The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.


## 7.4 pH Calibration

PATH:  \ Cal \ Calibrate Sensor

For pH sensors, the M800 transmitter features one-point, two-point or process calibration with 9 preset buffer sets or manual buffer entry. Buffer values refer to 25 °C. To calibrate the instrument with automatic buffer recognition, you need a standard pH buffer solution that matches one of these values. Please select the correct buffer table before using automatic calibration (see chapter 16 "Buffer tables"). The stability of the sensor signal during calibration can be checked by the user or automatically by the transmitter (see chapter 8.1.4.2 "pH Settings").



**NOTE:** For dual membrane pH electrodes (pH/pNa) only buffer Na+ 3.9M (see section 16.2.1 "Mettler-pH/pNa buffers") is available.



The following menus can be called up:

**Unit:** Select pH.

**Method:** Select the desired calibration procedure, 1-point, 2-point or process calibration.

**Options:** The buffer used for the calibration and the required stability of the sensor signal during the calibration can be selected (see also chapter. 8.1.4.2 "pH Settings"). The changes are valid until the calibration mode has been escaped. After the values defined in the configuration menu are valid again.

### 7.4.1 One-Point Calibration

With pH sensors a one point calibration is always performed as an offset calibration.



Press the button Cal for starting calibration.

Place the electrode in the buffer solution and press the Next button.

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.


The transmitter shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will get the message "Please re-install sensor".

## 7.4.2 Two-Point Calibration

With pH sensors a two point calibration is always performed as calibration of slope and offset.



Press the Cal button to start calibration.

Place the electrode in buffer solution 1 and press Next button.

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter prompts you to place the electrode in the second buffer solution.

Press Next button to proceed with the calibration.

The display shows the buffer the transmitter has recognized (Point 2) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

## 7.4.3 Process Calibration

With pH sensors a process calibration is always performed as an offset calibration.



Press the Cal button to start calibration.

Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu Screen if the related channel is selected in the display.

After determining the pH value of the sample, press the calibration icon in the Menu Screen again.

Enter the pH value of the sample. Press the Next button to start the calculation of the calibration results.

The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (press button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.5 ORP Calibration of pH Sensors

PATH:  \ Cal \ Calibrate Sensor

For pH sensors with solution ground based on ISM technology the M800 transmitter gives the option to make in addition to the pH calibration an ORP calibration.



**NOTE:** In case of choosing ORP calibration the parameters defined for pH (see chapter 8.1.4.2 "pH Settings") will not be considered. For pH sensors, the M800 transmitter features one-point calibration for ORP.



The following menu can be called up:

**Unit:** Select ORP through pressing the corresponding field.

**Method:** 1-Point calibration is displayed.

Press the button Cal for starting calibration.

Enter the value for calibration point 1 (**Point1**).

Press the Next button to start the calculation of the calibration results.

The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".


## 7.6 Calibration of Amperometric Oxygen Sensors

PATH:  \ Cal \ Calibrate Sensor

The M800 provides the ability to perform a one-point or process calibration for amperometric oxygen sensors.



**NOTE:** Before air calibration, for highest accuracy, enter the barometric pressure and relative humidity, as described in chapter 8.1.4.3 "Settings for Oxygen Measurement Based on Amperometric Sensors".



The following menu can be called up:


**Unit:** Between several units for DO and O2gas can be chosen.

**Method:** Select the desired calibration procedure, 1-point or process calibration.

**Options:** In case the method 1-point has been chosen the calibration pressure, relative humidity and - for slope calibration - the stability mode for the sensor signal during the calibration can be selected. For the method Process the values for the process pressure, calibration pressure and the parameter ProcCalPress can be modified. See also chapter 8.1.4.3 "Settings for Oxygen Measurement Based on Amperometric Sensors". The changes are valid until the calibration mode has been escaped. After the values defined in the configuration menu are valid again.

### 7.6.1 One-Point Calibration

A one-point calibration of oxygen sensors is always either a one point slope (i.e. with air) or a zero (offset) calibration. A one point slope calibration is done in air and a one point offset calibration is done at 0 ppb oxygen. A one-point zero dissolved oxygen calibration is available but not normally recommended since zero oxygen is very hard to achieve. A zero-point calibration is only recommended if high accuracy at low oxygen level (below 5% air) is needed.



Choose Slope or Offset calibration through pressing the corresponding field.

Press the button Cal for starting calibration.



**NOTE:** If the polarization voltage for the measuring mode and calibration mode is different, the transmitter will wait 120 seconds before starting the calibration. In this case the transmitter will also go after the calibration for 120 seconds to the HOLD Mode, before returning to the measuring mode again.

Place the sensor in air or the calibration gas and press Next button

Enter the value for the calibration point (**Point1**).

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.



**NOTE:** For an offset calibration the Auto mode is not available. If Auto mode has been chosen and afterwards slope calibration has been changed to offset calibration, the transmitter will perform the calibration in Manual mode.

The transmitter shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

### 7.6.2 Process Calibration

A process calibration of oxygen sensors is always either a slope or an offset calibration.



Choose Slope or Offset calibration through pressing the corresponding field.

Press the Cal button to start calibration.

Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu screen if the related channel is selected in the display.

After determining the oxygen value of the sample, press the calibration icon in the Menu Screen again.

Enter the oxygen value of the sample. Press the Next button to start the calculation of the calibration results.

The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.7 Calibration of Optical Oxygen Sensors

PATH:  \ Cal \ Calibrate Sensor

Oxygen calibration for optical sensors can be performed as a two-point, process or, depending on the sensor model connected to the transmitter, also as a one-point calibration.



**NOTE:** Before air calibration, for highest accuracy, enter the barometric pressure and relative humidity, as described in chapter 8.1.4.4 "Settings for Oxygen Measurement Based on Optical Sensors".




The following menus can be called up:

- Unit:** Between several units can be chosen. The units are displayed during the calibration.
- Method:** Select the desired calibration procedure, 1-point, 2-point or process calibration.
- Options:** In case the method 1-point has been chosen the calibration pressure, relative humidity and the stability mode for the sensor signal during the calibration can be selected. For the method Process the values for the process pressure, calibration pressure, the parameter ProcCalPress and the mode of the process calibration can be modified. See also chapter 8.1.4.4 "Settings for Oxygen Measurement Based on Optical Sensors". The changes are valid until the calibration mode has been escaped. After the values defined in the configuration menu are valid again.

### 7.7.1 One-Point Calibration

Typically a one point calibration is done in air. Nevertheless other calibration gases and solutions are possible.

The calibration of an optical sensor is always a calibration of the phase of the fluorescence signal towards the internal reference. During a one point calibration the phase in this point is measured and extrapolated over the measuring range.



Press the button Cal for starting calibration.

Place the sensor in air or the calibration gas and press Next button

Enter the value for the calibration point (**Point1**).

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.



The transmitter shows the value for the phase of the sensor at 100% air (P100) and at 0% air (P0) as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

## 7.7.2 Two-Point Calibration

The calibration of an optical sensor is always a calibration of the phase of the fluorescence signal towards the internal reference. A two-point calibration is a combination of first a calibration in air (100%) where a new phase P100 is measured and then a calibration in nitrogen (0%) where a new phase P0 is measured. This calibration routine gives the most accurate calibration curve over the whole measuring range.



Press the Cal button to start calibration.

Place the sensor in air or the calibration gas and press Next button

Enter the value for the first calibration point (**Point1**).

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter prompts you to change the gas.

Press Next button to proceed with the calibration.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter shows the value for the phase of the sensor at 100% air (P100) and at 0% air (P0) as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

## 7.7.3 Process Calibration



Press the Cal button to start calibration.

Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the start and Menu Screen if the related channel is selected in the display.

After determining the oxygen value of the sample, press the calibration icon in the Menu Screen.

Enter oxygen value of the sample. Press the Next button to start the calculation of the calibration results.

The display shows now the values for the phase of the sensor at 100% air (P100) and at 0% (P0) air.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).



**NOTE:** If for process calibration Scaling has been chosen (see chapter 8.1.4.4 "Settings for Oxygen Measurement Based on Optical Sensors") the calibration values are not stored in the calibration history.

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.8 Calibration of Dissolved Carbon Dioxide Sensors

For dissolved carbon dioxide (CO<sub>2</sub>) sensors, the M800 transmitter features one-point, two-point or process calibration. For the one-point or two-point calibration the solution with pH = 7.00 and/or pH = 9.21 of the Mettler – 9 standard buffer can be used (see also chapter 8.1.4.5 "Dissolved carbon dioxide settings") or the buffer value can be entered manually.



The following menus can be called up:

**Unit:** Between several units for partial pressure, and dissolved carbon dioxide can be selected.

**Method:** Select the desired calibration procedure, 1-point or process calibration.

**Options:** The buffer used for the calibration and the required stability of the sensor signal during the calibration can be selected (see also chapter 8.1.4.5 "Dissolved carbon dioxide settings"). The changes are valid until the calibration mode has been escaped. After the values defined in the configuration menu are valid again.

### 7.8.1 One-Point Calibration


With CO<sub>2</sub> sensors a one point calibration is always performed as an offset calibration.

Press the button Cal for starting calibration.

Place the electrode in the buffer solution and press the Next button.

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.





**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will get the message "Please re-install sensor".

## 7.8.2 Two-Point Calibration

With CO<sub>2</sub> sensors a two point calibration is always performed as calibration of slope and offset.



Press the Cal button to start calibration.

Place the electrode in buffer solution 1 and press Next button.

The display shows the buffer the transmitter has recognized (Point 1) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter prompts you to place the electrode in the second buffer solution.

Press Next button to proceed with the calibration.

The display shows the buffer the transmitter has recognized (Point 2) and the measured value.

The M800 checks the stability of the measuring signal and proceeds as soon as the signal is sufficiently stable.



**NOTE:** If **option** Stability is set to Manual press 'Next' after the measuring signal is stable enough to go on with the calibration.

The transmitter shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor".

### 7.8.3 Process Calibration

With CO<sub>2</sub> sensors a process calibration is always performed as an offset calibration.



Press the Cal button to start calibration.

Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu Screen if the related channel is selected in the display.

After determining the corresponding value of the sample, press the calibration icon in the Menu Screen again.

Enter the value of the sample. Press the Next button to start the calculation of the calibration results.

The display shows the value for the slope and the offset as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (press button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.9 Calibration of Thermal Conductivity CO<sub>2</sub> (CO<sub>2</sub> high) Sensors

PATH: H \ Cal \ Calibrate Sensor

The M800 provides the ability to perform an one-point calibration using a reference gas (CO<sub>2</sub>) with a known carbon dioxide partial pressure value. It offers also to perform a process calibration based on a analyzed process sample.



**NOTE:** The sensor is designed to measure CO<sub>2</sub> partial pressure or concentration values accurately in liquid phase only! In gas phase the sensor will only show correct CO<sub>2</sub> gas partial pressure values in the 1-point calibration menu.

The following menus can be called up:



**Unit:** Between the units of CO<sub>2</sub> pressure or concentration can be chosen.

**Method/Options:** Select the desired calibration procedure, 1-point or process calibration and stability option (manual/auto).

In case the method 1-point has been chosen only the calibration pressure and the option stability mode for the sensor signal during the calibration can be selected (Sensor expects to be in a calibration gas).

For the method Process only concentration values can be chosen as pressure or concentration values (Sensor expects to be in liquids).




**NOTE:** With reference Gas (CO<sub>2</sub>) --> Use 1-point calibration. With liquids --> Use process. When changing MembraCap always first perform a 1-point gas calibration. The changes are valid until the calibration mode has been exited. After the values defined in the Configuration menu are valid again.

## 7.9.1 One-Point Calibration



With the thermal conductivity sensor a one point calibration is always performed as a slope calibration. Press the button Cal for starting calibration.

Expose the TC-Sensor to a reference gas of a known CO2 concentration and press Next button. Enter the value for the calibration point (Point1) in mbar or hPa.



Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the baseline as the result of the calibration. The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.9.2 Process Calibration



With the thermal conductivity sensor a process calibration is always performed as a slope calibration.

Choose process calibration and desired unit in the calibration menu  
Press the button Cal for starting calibration.



Take a sample and press the  $\leftarrow$  button to store the current measuring value. To show the ongoing calibration process, P is blinking in the Start and Menu screen if the related channel is selected on the display.

After determining the CO2 value of the sample, press the calibration icon in the Menu Screen again. Enter the CO2 value of the sample.

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the baseline as the result of the calibration.

The calibration values are stored in the calibration history and taken over (button Adjust), stored in the cal history and not taken over (button Calibrate) or discarded (button Cancel).

If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed.

## 7.10 Calibration of O3 Sensors

The M800 provides the ability to perform a 1-Point or process calibration for O3 sensors. Dissolved Ozone must be performed quickly because O3 decays rapidly into oxygen, especially at warm temperatures.

Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration"; PATH:  $\mathcal{H}$ \CAL\Calibrate Sensor) and choose the desired channel for calibration.



The following menus can be called up:

**Unit:** Several units for dissolved O3 can be chosen.

**Method:** Select the desired calibration procedure, 1-Point or process calibration.

### 7.10.1 One-Point Calibration

Select the 1-Point calibration method. A 1-Point calibration of O3 sensors is always a zero (offset) calibration



Press the button Cal for starting calibration.



Place the sensor in the calibration gas, such as air, and press the Next button.



The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.

Press the input field for **Point1** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button to accept value.



The screen shows the entered value for the reference solution (1st line) and the measured value of the M800 (2nd line).

When the measuring signal is stable, press Next to continue with the calibration



The display shows the value for the slope and the offset as result of the calibration.

Press Adjust to perform the calibration and store the calibration values. Press Calibrate to calculate the calibration factors and view the results without saving them. Press Cancel to terminate the calibration.


Use the Back button to go one step back in the calibration procedure



If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. In either case you will see the message "Please re-install sensor". After pressing the Done button the M800 returns to the calibration menu.

### 7.10.2 Process Calibration

Select the Process calibration method. A Process calibration of O3 sensors can be performed as a slope or offset calibration.



Select the desired calibration **Method**.



Press Cal to start the calibration.



Take a sample and press the  $\leftarrow$  button to store the current measuring value. "P" will blink in the measurement screen indicating a Process calibration is active.



After determining the O3 value of the sample, press the calibration icon to complete the Process calibration.



Press the input field for **Point1** and enter the O3 value of the sample. Press the  $\leftarrow$  button to accept the value.

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as the result of the calibration.

Press Adjust to perform the calibration and store the calibration values. Press Calibrate to calculate the calibration factors and view the results without saving them. Press Cancel to terminate the calibration.

Use the Back button to go one step back in the calibration procedure.



If "Adjust" or "Calibrate" are chosen, the message "Calibration Saved Successfully!" is displayed. After pressing the Done button the M800 returns to the Menu Screen

## 7.11 Calibration of Flow Sensors

The M800 transmitter provides the ability to perform a 1-Point or 2-Point Sensor calibration for flow, Edit of saved calibration constants, and Verify of the flow signal. The most common method of calibration for flow sensors is to enter the calibration constants appropriate for the sensor using the Edit function. Some users may choose to perform an in-line calibration using a 1-point or 2-point Sensor flow calibration. This requires an external reference system. When performing an in-line calibration on a flow sensor, results will vary depending on the methods and calibration apparatus used to perform the calibration.

Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration"; PATH: \Cal\Calibrate Sensor) and choose the desired channel for calibration.



**NOTE:** The channel for Flow Type 2 can not be selected. The M800 transmitter allows during the Guided Setup (see chapter 6 "Guided Setup") to enter a table of K and F factors

Select the channel (four channel models only) and the desired calibration option. Choices are GPM, liters/minute meters<sup>3</sup>/hour, ft/sec, or meters/sec (for a 1-Point or 2-Point flow calibration), Edit and Verify. Press [ENTER].



The following menus can be called up:

- Unit:** Several units for Flow can be chosen.
- Method:** Select the desired calibration procedure, 1-point or 2-point calibration.

### 7.11.1 One-Point Calibration

Select calibration method 1-Point . A 1-Point calibration of a Flow sensor is always a slope calibration.


Press Cal to start the calibration.



Set the desired flow rate and press Next.









The second value displayed on the screen is the value being measured by the transmitter and sensor in the units selected by the user.

Press the input field for **Point1** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button to accept the value.



The screen shows the entered value for the reference system (1st line) and the measured value of the M800 (2nd line).

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as result of the calibration.

Selecting Cancel will discard the entered values and the M800 will return to the calibration menu.

Use the Back button to go one step back in the calibration procedure.

Press SaveCal to save the calibration factors.



If "SaveCal" is chosen, "Calibration Saved Successfully" and "Please re-install sensor" is displayed. After pressing the Done button the M800 returns to the calibration menu.

### 7.11.2 Two-Point Calibration


Select calibration method 2-Point . A 2-point calibration of a Flow sensor calculates a new slope and offset.



Press Cal to start the calibration




Set the desired flow rate for the first point and press Next.



The second value displayed on the screen is the value being measured by the transmitter and sensor in units selected by the user.

Press the input field for **Point1** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button to accept the value.




The screen shows the entered value for the reference system (1st line) and the measured value of the M800 (2nd line).

Press Next to continue the calibration.



Set the desired flow rate for the second point and press Next.



The second value displayed on the screen is the value being measured by the transmitter and sensor in units selected by the user.

Press the input field for **Point2** to enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press the  $\leftarrow$  button to accept the value.



The screen shows the entered value for the reference system (1st line) and the measured value of the M800 (2nd line).

Press the Next button to start the calculation of the calibration results.



The display shows the value for the slope and the offset as result of the calibration.

Selecting Cancel will discard the entered values and the M800 will return to the calibration menu.

Use the Back button to go one step back in the calibration procedure.

Press SaveCal to save the calibration factors.



If "SaveCal" is chosen, "Calibration Saved Successfully" and "Please re-install sensor" is displayed. After pressing the Done button the M800 returns to the calibration menu.

## 7.12 Sensor Verification

Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration";  
 PATH:  $\leftarrow$  \ Cal \ Calibrate Sensor) and choose the desired channel for verification




Press the Verify button to start verification.

The measured signal of the primary and the secondary measurement in basic (mostly electrical) units are shown. The meter calibration factors are used when calculating these values.

Press the  $\leftarrow$  button and the transmitter returns to the calibration menu.


## 7.13 Edit Calibration Constants for Flow Sensors


This function is the most commonly used calibration method for flow sensors


Enter the menu Calibrate Sensor (see chapter 7.1 "Sensor Calibration";  
PATH:  \ Cal \ Calibrate Sensor) and choose the desired channel.



Press the Edit button.



Press the input field for **Slope** to modify the slope value. The M800 displays a keypad for modifying the value. Press the  button and to accept the value.

Press the input field for **Offset** to modify the offset value. The M800 displays a keypad for modifying the value. Press the  button to accept the value.

Selecting Cancel will discard the entered values and the M800 will return to the calibration menu.

Press Save to save the calibration factors.



If "Save" is chosen "Calibration Saved Successfully" and "Please re-install sensor" is displayed. After pressing the Done button the M800 returns to the calibration menu.

Exit Menu Calibrate Sensor.

Press . To return to the Menu Screen press .

## 7.14 UniCond2e Electronics Calibration

The M800 provides the ability to calibrate or verify the electronic circuits of UniCond2e conductivity sensors. UniCond2e sensors have 3 resistance range circuits that require individual calibration. These measuring circuits are calibrated using the Thornton ISM Conductivity Sensor Calibration Module part number 58 082 305 and supplied Y-connector. Before calibration, remove the sensor from the process, rinse with deionized water and allow to completely dry. Power the transmitter and sensor at least 10 minutes prior to calibration to assure stable operating temperature of the circuitry.



Press the Cal button.

Enter menu Calibrate Electronics.

Press button Chan\_x and select the desired channel for calibration.

Choose **Verify** or **Cal**.

Reference Thornton ISM Conductivity Sensor Calibration Module (part number 58 082 305) for detailed calibration and verification instructions.

## 7.15 Flow Meter Calibration

Although it is not normally necessary to perform meter re-calibration unless extreme conditions cause an out of spec operation shown by Calibration Verification, periodic verification/re-calibration may be necessary to meet Q.A. requirements. The frequency calibration requires a 2-point calibration. It is recommended that point one be at the low end of the frequency range and point two at the high end.

Press the Cal button.



Enter menu Calibrate Meter.

Press button Chan\_x and select the desired channel for calibration.



Press Cal to start the calibration process.

Connect a frequency generator to the xIN and xGND terminals and press Next.

Adjust the frequency generator to the first calibration point and press Next.

Press the input field for **Point1** and enter the value for the calibration point. The M800 displays a keypad for modifying the value. Press  $\leftarrow$  to accept the value.



**NOTE:** To select another unit for the entered value on the keypad press the U button, to return to the keypad press 0–9.

Press Next

Adjust the frequency generator to the second calibration point and press Next.

Press the input field for **Point2** and enter the value for the second calibration point. The M800 displays a keypad for modifying the value. Press  $\leftarrow$  to accept the value.

The display shows the slope and offset values as the result of the calibration.

**SaveCal:** save new calibration values and store in the calibration history

**Cancel:** cancel the calibration procedure and return to Calibrate Sensor Menu.

**Back:** go one step back in the calibration procedure.

Press Done to return to Calibrate menu.

## 7.16 Flow Meter Verification

Enter Calibration Mode as described in chapter 7 "Enter Calibration Mode"

Enter menu Calibrate Meter.

Press button Chan\_x and select the desired channel for calibration.

Press Verify to start the verification process.

Connect a frequency generator to xIn and xGND and press Next.

The measured frequency is displayed.

Press ← to return to the calibration menu.

## 7.17 Analog Output Calibration

PATH:  \ CAL \ Calibrate Analog Outputs



Each analog output can be calibrated at 4 and 20 mA. Select the desired output signal for calibration by pressing the #1 button for output signal 1, #2 for output signal 2, etc.

Connect an accurate milliamp meter to the analog output terminals and then adjust the five digit number in the display until the milliamp meter reads 4.00 mA and repeat for 20.00 mA.

As the five digit number is increased the output current increases and as the number is decreased the output current decreases. Thus coarse changes in the output current can be made by changing the thousands or hundreds digits and fine changes can be made by changing the tens or ones digits.

After adjusting both values press Next button to start the calculation of the calibration results.


The display shows the calibration slope and zero point as the result of the output signal calibration.

Selecting Cancel will discard the entered values. Pressing SaveCal will making the entered values the current ones.

If "SaveCal" is chosen, "Calibration Saved Successfully" is displayed.

## 7.18 Analog Input Calibration

PATH:  \ CAL \ Calibrate Analog Outputs



Each analog input can be calibrated at 4 and 20 mA. Select the input signal for calibration by pressing the #1 button.

Connect an 4 mA signal to the analog input terminals. Press next button.

Enter the right value for the input signal (**Point1**).

Press the Next button go on with the calibration.

Connect an 20 mA signal to the analog input terminals. Press next button.

Enter the right value for the input signal (**Point2**).

Press the Next button go on with the calibration.

The display shows the calibration slope and zero point as the result of the input signal calibration.


Selecting Cancel will discard the entered values. Pressing SaveCal will making the entered values the current ones.

If "SaveCal" is chosen, "Calibration Saved Successfully" is displayed.

## 7.19 Maintenance

PATH:  \ CAL \ Calibrate Analog Outputs

The different channels of the transmitter M800 can be switched manually into Hold state. Furthermore a cleaning cycle can be started / stopped manually.




Select the channel, which should be set to Hold manually.

Press Start button for **Manual Hold** to activate the Hold state for the selected channel. To deactivate the Hold state again, press the Stop button, which is now displayed instead of the Start button.

Press the Start button for **Manual Clean** to switch the cleaning relay to the state for starting a cleaning cycle. To switch back the relay press the Stop button, which is now displayed instead of the Start button.

## 8 Configuration




### 8.1 Measurement

PATH: \ CONFIG \ Meas

#### 8.1.1 Channel Setup

PATH: \ CONFIG \ Meas \ Channel Setup



Select the **Channel** for the setup through pressing the button #1 for channel 1, #2 for channel 2 etc.

Press the right input field in the line of the setting for **Channel**. A parameter for the corresponding channel is chosen through pressing the according field.

If Auto is selected, M800 transmitter automatically recognizes the type of sensor. The channel can also be fixed to a certain measurement parameter (parameter = pH/ORP, pH/pNa, UniCond2e, Cond4e, O2 hi, O2 lo, O2 Trace, O2 opt, CO2, CO2 hi, O3 and Flow hi, Flow low, Flow Type2), depending on the type of transmitter.

Measurement parameter		Type
pH/ORP	= pH and ORP	Water, Process
pH/pNa	= pH and ORP (with pH/pNa electrode)	Process
UniCond2e	= 2 electrode conductivity	Water, Process
Cond4e	= 4 electrode conductivity	Water, Process
O2 hi	= Dissolved oxygen or oxygen in gas (ppm)	Process
O2 lo	= Dissolved oxygen or oxygen in gas (ppb)	Process
O2 Trace	= Dissolved oxygen or oxygen in gas	Process
O2 Opt	= Dissolved oxygen optical	Process
O2 lo Thornton	= Dissolved oxygen	Water
CO2	= Dissolved carbon dioxide	Process
CO2 hi	= Dissolved carbon dioxide for beverages	Process
TOC	= Total organic carbon	Water
O3	= Dissolved O3	Water
Flow hi, low, Type2	= Flow	Water

Enter the name with a maximum length of 6 characters for the channel through pressing the input field in the line **Descriptor**. The name of the channel will always be displayed, if the channel has to be selected. The name will also be displayed on the in the Start Screen and Menu Screen if the Display Mode (see chapter 8.1.3 "Display Mode") has been set to 1-Channel or 2-Channel.

Choose one of the measurements **M1 to M6** (e.g. for measuring value M1 the left button, for measuring M2 the right button in the corresponding line).

Select in the input field for **Measurement** the desired parameter to show.



**NOTE:** Beside the parameters pH, O2, T, etc. also the ISM values DLI, TTM and ACT can be linked to the measurements.

Choose **Range factor** of the measuring value. Not all parameters allow a modification of the range.

The menu **Resolution** allows the setting of the resolution for the measurement. The accuracy of the measurement is not effected by this setting. Possible setting are 1, 0.1, 0.01, 0.001.

Selected the menu **Filter**. The averaging method (noise filter) for the measurement can be selected. The options are None (default), Low, Medium, High and Special.

None = no averaging or filtering  
 Low = equivalent to a 3 point moving average  
 Medium = equivalent to a 6 point moving average  
 High = equivalent to a 10 point moving average  
 Special = averaging depending on signal change (normally High averaging but Low averaging for large changes in input signal)

## 8.1.2 Derived Measurements

The M800 enables the setup of derived measurements (total, difference, ratio) based on two measuring values like pH, conductivity, etc.. To get the derived measurements, first set up the two primary measurements, which will be used to calculate the derived measurement. Define the primary measurements as if they were stand-alone readings. Then chose the corresponding unit for the derived measurement for the first channel. The transmitter M800 will displays an additional menu **Other Channel** to select the second channel with the corresponding measurement.

There are three additional derived measurements available for configuration with two conductivity sensors: %Rej (% Rejection), pH Cal (Calculated pH) and CO<sub>2</sub> Cal (Calculated CO<sub>2</sub>).



### 8.1.2.1 % Rejection measurement

For reverse osmosis (RO) applications, percent rejection is measured with conductivity to determine the ratio of impurities removed from product or permeate water to the total impurities in the incoming feed water. The formula for obtaining Percent Rejection is:

$$[1 - (\text{Product/Feed})] \times 100 = \% \text{ Rejection}$$

Where Product and Feed are the conductivity values measured by the respective sensors.

Figure a shows a diagram of an RO installation with sensors installed for Percent Rejection.




Figure a: % Rejection



**NOTE:** The product monitoring sensor must be on the channel that will measure percent rejection. If the product conductivity sensor is installed in channel 1, then percent rejection must be measured in channel 1.

### 8.1.2.2 Calculated pH (Power Plant Applications only)

Calculated pH may be obtained very accurately from specific and cation conductivity values on power plant samples when the pH is between 7.5 and 10.5 due to ammonia or amines and when the specific conductivity is significantly greater than the cation conductivity. This calculation is not suitable where significant levels of phosphates are present. The M800 uses this algorithm when pH Cal is selected as a measurement.

The calculated pH must be configured on the same channel as specific conductivity. For example, set up measurement M1 on CHAN\_1 to be specific conductivity, measurement M1 on CHAN\_2 to be cation conductivity, measurement M2 on CHAN\_1 to be calculated pH and measurement M3 on CHAN\_1 to be temperature. Set the temperature compensation mode to "Ammonia" for measurement M1 on CHAN\_1 and to "Cation" for measurement M1 on CHAN\_2.



**NOTE:** If operation goes outside the recommended conditions, a glass electrode pH measurement is needed to obtain an accurate value. On the other hand, when sample conditions are within the ranges noted above, the calculated pH provides an accurate standard for one-point trim calibration of the electrode pH measurement.

### 8.1.2.3 Calculated CO<sub>2</sub> (Power plant applications only)

Carbon dioxide may be calculated from cation conductivity and degassed cation conductivity measurements on power plant samples using tables from ASTM Standard D4519. The M800 has these tables stored in memory, which it uses when units of CO<sub>2</sub> CAL are selected.

The calculated CO<sub>2</sub> measurement must be configured to the same channel as cation conductivity. For example, set up measurement M1 on CHAN\_1 to be cation conductivity, measurement M1 on CHAN\_2 to be degassed cation conductivity, measurement M2 on CHAN\_1 to be calculated CO<sub>2</sub> and measurement M2 on CHAN\_2 to be temperature. Set the temperature compensation mode to "Cation" for both conductivity measurements.

### 8.1.3 Display Mode

PATH:  \ CONFIG \ Meas \ Display Mode



Press the input field in the line of the setting for **Disp Mode** and choose the measuring values, which are displayed on the Start Screen and Menu Screen.

Choose between the display of the measuring values for one channel, the measuring values for two channels, 4 measuring values (4-meas) or 8 measuring values (8-meas).

**NOTE:** If 1-Channel or 2-Channel has been chosen the measuring values, that will be displayed are defined in the menu Channel Setup (see chapter 8.1.1 "Channel Setup"). If 1-Channel has been chosen, M1 to M4 of every channel will be displayed. In case of 2-Channel M1 and M2 of every channel will be displayed.

Additional settings can be done if 4-meas or 8-meas has been selected.

Select the **Page** of the Start Screen or Menu Screen the measuring value will be displayed.

Choose the **Line** of the according page the measuring value will be displayed.

Select the **Channel** which should be displayed in the according line of the page through pressing the corresponding field.


Choose the measured value of the selected channel which should be displayed through the parameter **Measure**.

**NOTE:** Beside the measurement values pH, O<sub>2</sub>, T, etc. also the ISM values DLI, TTM and ACT can be displayed.

## 8.1.4 Parameter related Settings

PATH:  \ CONFIG \ Meas \ Parameter Setting

Measuring and calibration parameters can be set for the parameters pH, conductivity, oxygen, and flow




Access the menu **Channel** and select the channel.

Depending on the selected channel and assigned sensor the measuring and calibration parameters are displayed.

See the following explanation to get more details about the different parameter settings.

### 8.1.4.1 Conductivity Settings



Select measurement (M1-M6). For more information regarding measurements see chapter 8.1.1 "Settings for Channel Setup".

If the selected measurement can be temperature compensated, the compensation method may be selected.



**NOTE:** During calibration, the compensation method must also be selected. (see chapter 7.2 "Calibration of UniCond2e Sensors" and chapter 7.3 "Calibration of Cond4e Sensors").

Press **Compens.** to select the desired temperature compensation method. Choices are "None", "Standard", "Light 84", "Std 75 °C", "Linear 25°C", "Linear 20°C", "Glycol.5", "Glycol1", "Cation", "Alcohol" and "Ammonia".

None does not make any compensation of the measured conductivity value. The uncompensated value will be displayed and proceeded.

Standard compensation includes compensation for non-linear high purity effects as well as conventional neutral salt impurities and conforms to ASTM standards D1125 and D5391.

Light 84 compensation matches the high purity water research results of Dr. T.S. Light published in 1984. Use only if your institution has standardized on that work.

Std 75 °C compensation is the Standard compensation algorithm referenced to 75 °C. This compensation may be preferred when measuring Ultrapure Water at an elevated temperature (Resistivity of ultrapure water compensated to 75 °C is 2.4818 Mohm-cm.)

Linear 25 °C compensation adjusts the reading by a coefficient or factor expressed as %/°C (deviation from 25 °C). Use only if the solution has a well-characterized linear temperature coefficient. The factory default setting is 2.0%/°C.

Linear 20 °C compensation adjusts the reading by a coefficient or factor expressed as %/°C (deviation from 20 °C). Use only if the solution has a well-characterized linear temperature coefficient. The factory default setting is 2.0%/°C.

Glycol.5 compensation matches the temperature characteristics of 50% ethylene glycol in water. Compensated measurements using this solution may go above 18 Mohm-cm.

Glycol1 compensation matches the temperature characteristics of 100% ethylene glycol. Compensated measurements may go well above 18 Mohm-cm.

Cation compensation is used in power industry applications measuring the sample after a cation exchanger. It takes into account the effects of temperature on the dissociation of pure water in the presence of acids.

Alcohol compensation provides for the temperature characteristics of a 75% solution of isopropyl alcohol in pure water. Compensated measurements using this solution may go above 18 Mohm-cm.


Ammonia compensation is used in power industry applications for specific conductivity measured on samples using ammonia and/or ETA (ethanolamine) water treatment. It takes into account the effects of temperature on the dissociation of pure water in the presence of these bases.



**NOTE:** If compensation mode "Linear 25 °C" or "Linear 20 °C" has been chosen, the coefficient for the adjustment of the reading can be modified. In this case an additional input field will be displayed.

Press the input field for **Coef.** and adjust the coefficient or factor for the compensation.

### 8.1.4.2 pH Settings



If a pH sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameters Buffer Tab, Stability, IP, STC and calibration temperature as well as the displayed units for slope and/or zero point can be set or adjusted. The same parameters will be displayed if during the channel setup not Auto but pH/ORP has been set.

Select the buffer through the parameter **Buffer Tab**.

For automatic buffer recognition during calibration, select the buffer solution set that will be used: Mettler-9, Mettler-10, NIST Tech, NIST Std = JIS Std, HACH, CIBA, MERCK, WTW, JIS Z 8802 or None. See chapter 16 "Buffer tables" for buffer values. If the auto buffer feature will not be used or if the available buffers are different from those above, select None.



**NOTE:** For dual membrane pH electrodes (pH/pNa) buffer Na+ 3.9M (see section 16.2.1 "Mettler-pH/pNa buffers") is available.

Select the required **Stability** of the measuring signal during the calibration procedure. Choose manual if the user will decide when a signal is stable enough to complete the calibration. Select Low, Medium or Strict if an automatic stability control of the sensor signal during calibration through the transmitter should be done.

If the parameter stability is set to medium (default) the signal deviation has to be less than 0.8 mV over a 20 second interval to be recognized by the transmitter as stable. The calibration is done using the last reading. If the criteria is not met within 300 seconds then the calibration times out and the message "Calibration Not Done" is displayed.

Adjust the parameter **IP pH**.

**IP** is the isothermal point value (Default = 7.000 for most applications). For specific compensation requirements or non standard inner buffer value, this value can be changed.


Adjust the value of the parameter **STC pH/°C**.

STC is the solution temperature coefficient in units of pH/°C referenced to the defined temperature. (Default = 0.000 pH/°C for most applications). For pure waters, a setting of -0.016 pH/°C should be used. For low conductivity power plant samples near 9 pH, a setting of -0.033 pH/°C should be used.

If the value for STC is  $\neq$  0.000 pH/°C an additional input field for the reference temperature will be displayed.

The value for **pH Ref Temperature** indicates to which temperature the solution temperature compensation is referenced. The displayed value and the output signal is referenced to this temperature. Most common reference temperature is 25°C.

### 8.1.4.3 Settings for Oxygen Measurement Based on Amperometric Sensors



If an amperometric oxygen sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameters CalPressure, ProcPressure, ProcCalPress, Stability, Salinity, RelHumidity, UpolMeas and UpolCal can be set or adjusted. The same parameters will be displayed if during the channel setup not Auto but O2 hi, O2 lo or O2 trace has been set.

Enter the value for the calibration pressure through the parameter **CalPressure**.

**NOTE:** For a modification of the unit for the calibration pressure press 'u' on the displayed keypad.

Press the button Option for the parameter **ProcPressure** and select the how to get applying process pressure through choosing the **Type**.


The applied process pressure can be entered by choosing Edit or measured over the analog input of the M800 by choosing Ain\_1.

If Edit has been chosen an input field for entering the value manually is displayed on the screen. In case that Ain\_1 has been selected two input fields are displayed to enter the start value (4 mA) and the end value (20 mA) of the range for the 4 to 20 mA input signal.

For the algorithm of the process calibration the applied pressure has to be defined. Select the pressure through the parameter **ProcCalPress**. For the process calibration the value of the process pressure (ProcPress) or the calibration pressure (CalPress) can be used.

Select the required **Stability** of the measuring signal during the calibration procedure. Choose Manual if the user will decide when a signal is stable enough to complete the calibration. Select Auto and an automatic stability control of the sensor signal during calibration through the transmitter will be done.

Additional settings can be done by navigating to the next page of the menu.



The **Salinity** of the measured solution can be modified.

In addition the relative humidity (button **Rel.Humidity**) of the calibration gas can also be entered. The allowed values for relative humidity are in the range 0% to 100%. When no humidity measurement is available, use 50% (default value).

The polarization voltage of amperometric oxygen sensors in the measuring mode can be modified through the parameter **UpolMeas**. For entered values 0 mV to -550 mV the connected sensor will be set to a polarization voltage of -500mV. If the entered value is less then -550 mV, the connected sensor will set to a polarization voltage of -674 mV.

The polarization voltage of amperometric oxygen sensors for calibration can be modified through the parameter **UpolCal**. For entered values 0 mV to -550 mV the connected sensor will be set to a polarization voltage of -500mV. If the entered value is less then -550mV, the connected sensor will set to a polarization voltage of -674mV.

**NOTE:** During a process calibration, the polarization voltage UpolMeas, defined for the measuring mode, will be used.



**NOTE:** If a one point calibration is executed, the transmitter sends the polarization voltage, valid for the calibration, to the sensor. If the polarization voltage for the measuring mode and calibration mode is different, the transmitter will wait 120 seconds before starting the calibration. In this case the transmitter will also go after the calibration for 120 seconds to the HOLD Mode, before returning to the measuring mode again.

### 8.1.4.4 Settings for Oxygen Measurement Based on Optical Sensors



If an optical oxygen sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel setup") Auto has been chosen the parameters CalPressure, ProcPressure, ProcCalPress, Stability, Salinity, RelHumidity, Sample Rate, LED Mode and Toff can be set or adjusted. The same parameters will be displayed if during the channel setup not Auto but Optical O2 has been set.

Enter the value for the calibration pressure through the parameter **CalPressure**.

Press the button Option for the parameter **ProcPress** and select the how to get applying process pressure through pressing the according button in the line **Type**.


The applied process pressure can be entered by choosing Edit or measured over the analog input of the M800 by choosing AIN\_1.

If Edit has been chosen an input field for entering the value manually is displayed on the screen. In case that AIN\_1 has been selected two input fields are displayed to enter the start value (4mA) and the end value (20 mA) of the range for the 4 to 20 mA input signal.

For the algorithm of the process calibration the applied pressure has to be defined. Select the pressure through the parameter **ProcCal**. For the process calibration the value of the process pressure (ProcPress) and the value of the calibration pressure (CalPress) can be used. Select between Scaling and Calibration for the process calibration. If Scaling has been chosen, the calibration curve of the sensor will be untouched, but the output signal of the sensor will be scaled. In case of calibration value <1%, the offset of the sensor output signal will be modified during scaling, for value >1% the slope of the sensor output will be adjusted. For further information about scaling refer to the sensor manual.

Selecting the required **Stability** of the measuring signal during the calibration procedure. Choose Manual if the user will decide when a signal is stable enough to complete the calibration. Select Auto and an automatic stability control of the sensor signal during calibration through the transmitter will be done.

Additional settings can be done by navigating to the next page of the menu.



The **Salinity** of the measured solution can be modified.

In addition the relative humidity (button **Rel.Humidity**) of the calibration gas can also be entered. The allowed values for relative humidity are in the range 0% to 100%. When no humidity measurement is available, use 50% (default value).

Adjust the required **Sample Rate** of the optical sensor during measurement. The time interval from one measuring cycle of the sensor to the next can be adjusted i.e. adapted to the application. A higher value will increase the life time of the OptoCap of the sensor.

Select the **LED Mode** of the sensor. There are the following options.

Off: LED is permanently switched off.

On: LED is permanently switched on.

Auto: The LED is switched on as long as the measured media temperature is smaller then Toff (see next value) or switched off through a digital input signal (see chapter 8.9 "Digital Inputs").



**NOTE:** If the LED is switched off, no oxygen measurement is performed.

Enter the limit for the measuring temperature to switch off the LED of the sensor automatically by the M800 through the parameter **Toff**.

If the media temperature is higher then Toff, the LED will be switched off. The LED will be switched on as soon as the media temperature falls below Toff – 3K. This function give the option to increase the lifetime of the OptoCap by switching off the LED during SIP or CIP cycles.



**NOTE:** This function is only active if the LED Mode is set to "Auto".

### 8.1.4.5 Dissolved Carbon Dioxide Settings



If an dissolved carbon dioxide sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel setup") Auto or CO<sub>2</sub> has been chosen, the buffer used for calibration and the parameters stability, salinity, HCO<sub>3</sub>, TotPres can be set resp. adjusted.

Select the buffer through the parameter **Buffer Tab**. For automatic buffer recognition during calibration, select buffer solution Mettler-9 if it will be used. If the auto buffer feature will not be used or if the available buffer are different from Mettler-9 select None.

Select the required **Stability** of the measuring signal during the calibration procedure. Choose manual if the user will decide when a signal is stable enough to complete the calibration. Select Low, Medium or Strict if an automatic stability control of the sensor signal during calibration through the transmitter should be done.

If the unit for the measured dissolved carbon dioxide is %sat, the pressure during the calibration resp. measurement has to be considered. This will be done by setting the parameter **TotPres**. If another unit then %sat has been selected, the result will not be influenced by this parameter.

The **Salinity** describes the total amount of solved salts in the CO<sub>2</sub> electrolyte of the sensor connected to the transmitter. It is a sensor specific parameter. The default value (28.00 g/L) is valid for the InPro 5000i. Do not change this parameter if the InPro 5000i will be used.


Additional settings can be done by navigating to the next page of the menu.



The parameter **HCO<sub>3</sub>** describes the concentration of hydrogen carbonate in the CO<sub>2</sub> electrolyte of the sensor connected to the transmitter. It is also a sensor specific parameter. The default value 0.050 Mol/L is valid for the InPro 5000i. Do not change this parameter if the InPro 5000i will be used.

### 8.1.4.6 Settings for Thermal Conductivity Dissolved CO<sub>2</sub> Measurement (CO<sub>2</sub> hi)

If during the channel setup (see chapter 8.1.1 "channel setup") the parameter CO<sub>2</sub> Hi has been chosen, the parameters stability (manual/auto) and CO<sub>2</sub> solubility (CO<sub>2</sub>-solubility and Temperature Factor), be set resp. adjusted.



Select the required **Stability** of the measuring signal during the calibration procedure. Choose manual if the user will decide when a signal is stable enough to complete the calibration. Select Auto if an automatic stability control of the sensor signal during calibration through the transmitter should be done.

The sensor offers a choice of CO<sub>2</sub> **Solubility's** for measurement in beer, water and cola. The cola setting is to be used with carbonated soft drinks. For other beverages the user has the possibility to enter individual values for CO<sub>2</sub> solubility and temperature factors.

Default values for measurement in beer (valid for temperatures – 5 ... 50 °C):

CO<sub>2</sub> solubility (A): 1.420 g/L  
Temp. factor (B): 2485

Values for pure water:


CO<sub>2</sub> solubility (A): 1.471 g/L  
Temp. factor (B): 2491

Values for cola:

CO<sub>2</sub> solubility (A): 1.345 g/L  
Temp. factor (B): 2370



**NOTE:** The sensor is delivered factory calibrated and is set up to measure in beer as the default.



For beverages where the user knows the exact CO<sub>2</sub> solubility and the temperature factor the values can be changed **individually**.

If the user desires to evaluate the solubility (**CO<sub>2</sub>-solub.**) and temperature factors (**Temp.-Factor**) they can be evaluated with the following formulas.

$$HCO_2 = A \cdot 3 \exp(B \cdot 3 (1/T - 1/298.15))$$

$$cCO_2 = HCO_2 \cdot 3 pCO_2$$

HCO<sub>2</sub> : Calculated CO<sub>2</sub> Solubility (Henry constant) at measured process temp.

A: Solubility of CO<sub>2</sub> (g / L at 25 °C)

B: Temperature factor (valid for – 5 ... 50 °C)


cCO<sub>2</sub>: Calculated CO<sub>2</sub> concentration in g/l or V/V

### 8.1.4.7 Settings for TOC Measurement

For information on how to configure parameter related settings associated with TOC measurement, refer to the 5000TOCi operating manual provided with the 5000TOCi Total Organic Carbon sensor.



### 8.1.4.8 Settings for Flow Measurement



For the dedicated flow channels, Pipe ID and Flow Cut can be adjusted



Press the input field for the value of the parameter **Pipe ID**. The M800 displays a keypad for modifying the value. Flow velocity measurements in ft/sec require the inside diameter of the pipe (in which the sensor is installed) for calculations. Enter the precise inside diameter in inches. . Press  $\leftarrow$  to accept the value.



To select the low flow cutoff, select **Flow Cut**. The M800 displays a keypad for modifying the value. If the measured value is less than the entered value for Flow Cut, the M800 transmitter sets the measured value for flow to 0. Press  $\leftarrow$  to accept the value.

### 8.1.4.9 Deionization Capacity (DI-Cap™)

**Total equivalents, ppm-gallons or grains:** The M800 can monitor the flow rate and mineral concentration entering a deionization vessel and infer the extent of resin capacity consumption. By multiplying the total dissolved solids (TDS) based on conductivity times flow rate and integrating the result over time, total mineral content that has entered the deionization vessel can be monitored.

The M800 can provide this in real time by setting measurement units to equivalents, ppm-gallons or grains. From this and knowledge of the total ion exchange capacity of the vessel, the “% of run” and/or anticipated time for next regeneration can be determined. This measurement requires installation of a flow sensor and a conductivity sensor.

To set up a deionization capacity measurement:

1. Set up the measurement for the flow sensor.
2. Set up the measurement for the conductivity sensor.
3. Set up an additional measurement for deionization capacity for the flow channel with units of equiv, ppm-G or Grains.
4. Following the path <Config/Measurement/Parameter Setting> for the flow measurement channel, select the channel of the conductivity measurement being used for the calculation and set the TDS factor appropriately for the water being measured as determined below.

Deionization capacity can be reset to zero, similar to total flow, using the path <Config/Reset/Configure/Chan X/DiCap/Reset>.

**Total Dissolved Solids (TDS)** may be inferred and displayed based on conductivity/resistivity data. TDS is the concentration of sodium chloride (or other conductive substance) corresponding to the measured conductivity. Salinity is the same as TDS, specifically for sodium chloride. Both are given in units of parts per billion (ppb), parts per million (ppm) or parts per thousand (ppk, as abbreviated on M800).

The default setting of 1.0 for the TDS Factor provides conversion based on the conductivity of sodium chloride at 0.462 ppm per  $\mu\text{S}/\text{cm}$ , with non-linear corrections at very low and very high


conductivities. The TDS factor may be changed to provide accurate conversion for other compositions. It is a multiplier on the sodium chloride conversion. Values for other materials are given in the table below (normalized to NaCl). These values adjust the TDS value for the actual conductivity of the materials in the table. Different values are needed for measurements involving ion exchange calculations - see below.

Material	TDS factor
NaCl	1.0000
KCl	1.0786
CaCl2	0.8839
CaCO3	0.8407
NaOH	0.3480

**Total Dissolved Solids for ion exchange calculations** are based on the conductivity and weight of the materials present expressed as their ion exchange equivalent as calcium carbonate. For a defined composition of neutral minerals with conductivity the same as sodium chloride, a TDS factor of 0.856 will give readout as ppm NaCl expressed as CaCO3. For conditions of strong base exchange, a TDS factor of 0.435 will give readout as ppm NaOH expressed as CaCO3.

### 8.1.5 Concentration Curve Table


To specify a concentration curve for customers-specific solutions, up to 5 concentration values can be edited in a matrix together with up to 5 temperatures. To do so the desired values are edited under the concentration curve table menu. Beside the temperature values, the conductivity and concentration values for the corresponding temperature are edited. The concentration curve can be selected resp. used in combination with conductivity sensors.



Enter the name with a maximum length of 6 characters for the concentration curve through pressing the input field in the line **Descriptor**.

Enter the amount of desired temperature points (**TempPoint**) and concentration points (**ConcPoint**).

The different values can be entered by navigating to the next page of the menu.



Enter the values for temperature (**T1...T5**), concentration (**Conc1...Conc5**) and the corresponding conductivity through pressing the according input field. The unit for the value of the conductivity can be adjusted as well in the according input field.



**NOTE:** The values for the temperature have to increase from T1 to T2 to T3 etc.. The values for the concentration have to increase from Conc1 to Conc2 to Conc3 etc..



**NOTE:** The conductivity values at the different temperatures have to increase or decrease from Conc1 to Conc2 to Conc3 etc.. Maxima and/or minima are not permitted. If the conductivity values at T1 are increasing with the different concentrations, they have to increase also at the other temperatures. If the conductivity values at T1 are decreasing with the different concentrations, they have to decrease also at the other temperatures.

## 8.2 Analog Outputs

PATH:  \ CONFIG \ Analog Outputs

See the following explanation to get more details about the different settings for the analog outputs.



Press the input field in the line of the setting for **Aout** and select the desired output signal for configuration by pressing button #1 for output signal 1, #2 for output signal 2 etc. Press the related button for the assignment of the channel (**Chan**). Select the channel, which has to be linked to the output signal.

Press the button for the assignment of the measuring parameter – based on the selected channel – that has to be linked to the output signal.



**NOTE:** Beside the measurement values pH, O<sub>2</sub>, T, etc. also the ISM values DLI, TTM and ACT can be linked to the output signal.

Select the **Range** for the output signal.

To adjust the value for the analog output signal if an alarm occurs, press the input field in the line for the setting of **Alarm**. Off means, that an alarm has now influence on the output signal.




**NOTE:** Not only the alarms occurred on the assigned channel will be considered, but every alarm coming up on the transmitter.

The value for the output signal if the transmitter goes into hold mode can be defined. It can be chosen between the last value (i.e. the value before the transmitter switched to the hold mode) or an fixed value.

Press the input field in the line for the setting of the **Hold Mode** and select the value. If a fixed value is chosen, the transmitter shows an additional input field.

Additional settings can be done by navigating to the next page of the menu.



The **Aout Type** can be Normal, Bi-Linear, Auto-Range or Logarithmic. The range can be 4–20 mA or 0–20 mA. Normal provides linear scaling between the minimum and maximum scaling limits and is the default setting. Bi-Linear will also prompt for a scaling value for the mid-point of the signal and allows two different linear segments between the minimum and maximum scaling limits.

Press the button for the **Min Value**, that corresponds with start point of the analog output range.

Press the button for the **Max Value**, that corresponds with end point of the analog output signal.

Depending on the chosen Aout type additional values can be entered.

**Bi-Linear** will also prompt for a scaling value for the Mid Value of the signal and allows two different linear segments between the defined Min and Max Values.

**Auto-Range** scaling provides two ranges of output. It is designed to work with a PLC to provide a wide measurement range at the high end of the scale, and a narrower range with high resolution at the low end. Two separate settings are used, one for the maximum limit of the high range and one for the maximum limit of the low range, for the single 0/4-20 mA signal.

Max1 is the maximum limit of the low range on auto-range. The maximum value for the high range on auto-range is set with the Max Value. Both ranges have the same minimum value that

is set through Min Value. If the input value is higher than value of Max1, the transmitter switches automatically to the second range. To indicate the currently valid range a relay can be assigned. The relay will be switched if the transmitter changes from on range to the other.

If **Logarithmic** Range was selected, it will prompt for the Max Value and also for the number of decades.

## 8.3 Set Points

PATH:  \ CONFIG \ Set Points

See the following explanation to get more details about the different settings for the set points.



Press the input field in the line of the setting for **Set Point** and select the desired set point for configuration through pressing the button #1 for set point 1, #2 for set point 2 etc..

Press the related button for the assignment of the channel (**Chan**). Select the channel, which has to be linked to the set point.

Press the button for the assignment of the measuring parameter – based on the selected channel – that has to be linked to the set point.

Mx in the display indicates the measurement assigned to the set point. (see chapter 8.1.1 "Channel Setup").



**NOTE:** Beside the parameters pH, O<sub>2</sub>, T, mS/cm, %EP WFI etc. also the ISM values DLI, TTM and ACT can be linked to the set point.

The **Type** of the setpoint can be High, Low, Between, Outside or Off. An "Outside" setpoint will cause an alarm condition whenever the measurement goes above its high limit or below its low limit. A "Between" setpoint will cause an alarm condition to occur whenever the measurement is between its high and low limits.



**NOTE:** If the type of set point is not Off additional settings can be done. See the following description.

According to the selected type of setpoint, value(s) regarding the limit(s) can be entered.

Additional settings can be done by navigating to the next page of the menu.



Once configured a relay could be activated if a sensor **Out of Range** condition is detected on the assigned input channel.

To select the desired relay that will be activated if the defined conditions are reached press the input field in the line for the setting of **SP Relay**. If the chosen relay is used for another task, the transmitter shows the message on the screen that there is a Relay Conflict.

The operation mode of the relay can be defined.

Relay contacts are in normal mode until the associated setpoint is exceeded, then the relay is activated and the contact states change. Select Inverted to reverse the normal operating state of the relay (i.e. normally open contacts are in a closed state, and normally closed contacts are in an open state, until the setpoint is exceeded).

Enter the **Delay** time in seconds. A time delay requires the setpoint to be exceeded continuously for the specified length of time before activating the relay. If the condition disappears before the delay period is over, the relay will not be activated.


Enter the value for the **Hysteresis**. A hysteresis value requires the measurement to return within the setpoint value by a specified percentage before the relay is deactivated.

For a high setpoint, the measurement must decrease more than the indicated percentage below the setpoint value before the relay is deactivated. With a low setpoint, the measurement must rise at least this percentage above the setpoint value before the relay is deactivated. For example, with a high setpoint of 100, when this value is exceeded, the measurement must fall below 90 before the relay is deactivated.

Enter the relay **Hold Mode** of "Off", "Last Value" or "On". This is the state of the relay during hold status.

## 8.4 ISM Setup

PATH:  \ CONFIG \ ISM Setup




See the following explanation to get more details about the different parameter settings for the ISM Setup.

### 8.4.1 Sensor Monitor

If a pH/ORP, O2 hi, O2 lo, O2 trace, O3 or O2 optical sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameter Sensor Monitor can be set or adjusted. The menu Sensor Monitor will also be displayed if during the channel setup not Auto but one of the mentioned sensors has been set.

Press the button Sensor Monitor.



Enter the value for the initial Time To Maintenance interval (**TTM Initial**) in days. The initial value for TTM can be modified according to the application experience.

For pH/ORP sensor the timer estimates when the next cleaning cycle should be performed to keep the best possible measurement performance. The timer is influenced by significant changes on the DLI parameters.

For amperometric oxygen and ozone sensors, the time to maintenance indicates a maintenance cycle for the membrane and electrolyte.

Press the input field for **TTM Reset**. Select Yes if Time To Maintenance (TTM) for the sensor should be reset to the initial value.

Time To Maintenance needs to be reset after the following operations.

pH sensors: manual maintenance cycle on the sensor.  
 Oxygen or ozone sensor: manual maintenance cycle on the sensor or exchanging of the membrane of the sensor



**NOTE:** The menu TTM Initial and TTM Reset is for O2 optical sensors not available.



**NOTE:** By connecting a sensor, the actual value for TTM of the sensor is read out from the sensor.

Enter the **ACT Initial** value in days. The new value will be loaded down to the sensor after saving the changes.

The Adaptive Calibration Timer (ACT) estimates when the next calibration should be performed to keep the best possible measurement performance. The timer is influenced by significant changes on the DLI parameters. The ACT will be reset to his initial value after a successful calibration. The initial value for the ACT can be modified according to the application experience and loaded down to the sensor.



**NOTE:** By connecting a sensor, the actual value for the ACT of the sensor is read out from the sensor.

Press the input field for **DLI Reset**. Select Yes if Dynamic Lifetime Indicator (DLI) for the sensor should be reset to the initial value. The reset will be done after saving the changes.

The DLI allows an estimation, when the pH electrode, the inner body of an amperometric oxygen or ozone sensor or the OptoCap of an optical oxygen sensor is at the end of his lifetime, based on the actual stress he is exposed to. The sensor permanently takes the averaged stress of the past days into consideration and is able to increase/decrease the lifetime accordingly.

The following parameters affect the lifetime indicator:

**Dynamic parameters:**

- Temperature
- pH or oxygen value
- Glass impedance (only pH)
- Reference impedance (only pH)

**Static parameters:**

- Calibration history
- Zero and Slope
- Phase 0 and phase 100 (only optical DO)
- Illumination time (only optical DO)
- Sampling rate (only optical DO)
- CIP/SIP/Autoclaving cycles

The sensor keeps the information stored in the built in electronics and can be retrieved via a transmitter or the iSense asset management suite.

For amperometric oxygen sensors, the DLI is related to the inner-body of the sensor. After exchanging the inner-body perform DLI Reset.

For optical DO sensors the lifetime indicator is related to the OptoCap. After exchanging the OptoCap perform DLI Reset.



**NOTE:** By connecting a sensor, the actual values for the DLI of the sensor are read out from the sensor.




**NOTE:** The menu DLI Reset for pH sensors not available. If the actual value for the DLI of a pH sensor is 0 the sensor has to be replaced.

## 8.4.2 CIP Cycle Limit

If a pH/ORP, oxygen or conductivity sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameter CIP Cycle Limit can be set or adjusted. The menu CIP Cycle Limit will also be displayed if during the channel setup not Auto but one of the mentioned sensors has been set.

Press the button CIP Cycle Limit.



Press the button in the input field for the parameter **Max Cycles** and enter the value for the maximum CIP cycles. The new value will be written to the sensor after saving the changes.

The CIP cycles are counted by the transmitter. If the limit (value for Max Cycles) is reached, an alarm can be indicated and set to a certain output relays.

If the Max Cycles setting is on 0, the counter functionality is turned off.



**NOTE:** In case of an optical oxygen sensor, the value for Max Cycles will also be written to the sensor. The transmitter M800 uploads the value Max Cycles from an optical oxygen sensor after the connection.

Press the button in the input field for the parameter **Temp** and enter the temperature, which has to be exceeded, that the a CIP cycle will be counted.

CIP Cycles will be automatically recognized by the transmitter. Since CIP cycles will vary in intensity (duration and temperature) for each application the algorithm of the counter recognizes an increase of the measurement temperature above the level defined through the value for Temp. If the temperature does not decrease below the defined temperature level - 10°C within the next 5 minutes after the first temperature was reached, the counter in question will be incremented by one and also locked for the next two hours. In the case the CIP would last longer than two hours the counter would be incremented by one once more.

Press the input field for **Reset**. Select Yes if CIP counter for the sensor should be reset to 0. The reset will be done after saving the changes.

If an oxygen sensor is connected, the reset should be performed after the following operations.  
 optical sensor: exchanging of the OptoCap  
 amperometric sensor: exchanging of the inner-body of the sensor.




**NOTE:** For pH/ORP sensor the menu Reset is not available. A pH/ORP sensor should be replaced if the number for Max Cycles has been exceeded.

### 8.4.3 SIP Cycle Limit

If a pH/ORP, oxygen or conductivity sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameter SIP Cycle Limit can be set or adjusted. The menu SIP Cycle Limit will also be displayed if during the channel setup not Auto but one of the mentioned sensors has been set.

Press the button SIP Cycle Limit.



Press the button in the input field for the parameter **Max Cycles** and enter the value for the maximum SIP cycles. The new value will be written to the sensor after saving the changes.

The SIP cycles are counted by the transmitter. If the limit (value for Max Cycles) is reached, an alarm can be indicated and set to a certain output relays.

If the Max Cycles setting is on 0, the counter functionality is turned off.



**NOTE:** In case of an optical oxygen sensor, the value for Max Cycles will also be written to the sensor. The transmitter M800 uploads the value Max Cycles from an optical oxygen sensor after the connection.

Press the button in the input field for the parameter **Temp** and enter the temperature, which has to be exceeded, that the a SIP cycle will be counted.

SIP Cycles will be automatically recognized by the transmitter. Since SIP cycles will vary in intensity (duration and temperature) for each application the algorithm of the counter recognizes an increase of the measurement temperature above the level defined through the value for Temp. If the temperature does not decrease below the defined temperature level - 10°C within the next 5 minutes after the first temperature was reached, the counter in question will be incremented by one and also locked for the next two hours. In the case the SIP would last longer than two hours the counter would be incremented by one once more.

Press the input field for **Reset**. Select Yes if SIP counter for the sensor should be reset to 0. The reset will be done after saving the changes.

If an oxygen sensor is connected, the reset should be performed after the following operations.  
 Optical sensor: exchanging of the OptoCap  
 Amperometric sensor: exchanging of the inner-body of the sensor.



**NOTE:** For pH/ORP sensor the menu Reset is not available. A pH/ORP sensor should be replaced if the number for Max Cycles has been exceeded.

## 8.4.4 AutoClave Cycle Limit

If a pH/ORP, amperometric oxygen or, depending on the model, optical oxygen sensor is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameter AutoClave Cycle Limit can be set or adjusted. The menu AutoClave Cycle Limit will also be displayed if during the channel setup not Auto but one of the mentioned sensors has been set.

Press the button AutoClave Cycle Limit.



Press the button in the input field for the parameter **Max Cycles** and enter the value for the maximum AutoClave cycles. The new value will be written to the sensor after saving the changes.

If the Max Cycles setting is on 0, the counter functionality is turned off.

Since during the autoclaving cycle the sensor is not connected to the transmitter, you will be asked after every sensor connection, whether the sensor was autoclaved or not. According to your selection, the counter will be incremented or not. If the limit (value for Max Cycles) is reached, an alarm can be indicated and set to a certain output relay.



**NOTE:** In case of an optical oxygen sensor, the value for AutoClave Max will be written to the sensor. The transmitter M800 uploads the value Max Cycles from an optical oxygen sensor after plugging in.

Press the input field for **Reset**. Select Yes if the AutoClave counter for the sensor should be reset to 0. The reset will be done after saving the changes.

If an oxygen sensor is connected, the reset should be performed after the following operations.  
 Optical sensor: exchanging of the OptoCap  
 Amperometric sensor: exchanging of the inner-body of the sensor.



**NOTE:** For pH/ORP sensor the menu Reset is not available. A pH/ORP sensor should be replaced if the number for Max Cycles has been exceeded.



## 8.4.5 DLI Stress Adjustment

If a pH/ORP is connected to the selected channel while during the channel setup (see chapter 8.1.1 "Channel Setup") Auto has been chosen the parameter DLI Stress Adjustment can be adjusted. With this setting the user can adjust the sensor sensitivity to the stress of his specific application for the DLI calculation.



Browse to page 2 of "ISM Setup".

Press the button **DLI Stress Adjustment**.

Select between low / medium / high for the **Type** of DLI Stress Adjustment.

LOW: DLI extended (-30% sensitivity )

MEDIUM: standard DLI (default)

HIGH: DLI reduced (+30% sensitivity)

Press **↵** to accept the setting.

## 8.4.6 SAN Cycle Parameters

If an ozone sensor is connected , values for the following SAN Cycle Parameters can be set, Max Cycles (the maximum number of sanitization cycles), Conc. Max (the maximum allowed O3 concentration), Conc. Min (the minimum allowed O3 concentration), Cycle Time (length of cycle), and Reset.

Press the button SAN Cycle Parameters.



Press the input field next to Max Cycles and enter the value for the maximum SAN cycles. Press **↵** to accept the value. The new value will be written to the sensor after saving the changes.

The SAN cycles are counted by the transmitter. If the limit (value for Max Cycles) is reached, an alarm can be configured. If the Max Cycles setting = 0, the counter functionality is turned off.

Press the input field next to Conc. Max and enter the ozone concentration above which a sanitization cycle is to be detected. Press **↵** to accept the value.

Press the input field next to Conc. Min. Enter the value for the ozone concentration below which a sanitization cycle is no longer detected. Press **↵** to accept the value

Press the input field next to Cycle Time. Enter the value for the time, the ozone concentration has to be higher than the Conc. Min value after the Conc. Max value has been exceeded to count a sanitization cycle. Press  $\leftarrow$  to accept the value.

Press the input field next to Reset. Select Yes to reset the sanitization counter to zero. This is typically performed after sensor replacement. The reset will be done after saving the changes

Press  $\leftarrow$  to exit the menu SAN Cycle Parameters.

## 8.4.7 Reset Counters for UniCond2e Sensors

For UniCond2e sensors, the following counters can be reset: High Temp and High Conductivity.

Press the button Reset Counters.




Select Yes for the desired counter to be reset and press enter. The reset will be done after saving the changes.

Press  $\leftarrow$  to exit the menu Reset Counters.

## 8.4.8 Set Calibration Interval for UniCond2e Sensors

For UniCond2e sensor the Cal Interval (calibration interval) can be set.

Press the button Cal Interval.



Press the input field next to **Cal Interval** and enter the value for the calibration interval. Based on this value the Time To Calibration (TTCal) will be calculated by the transmitter. Press  $\leftarrow$  to accept the value. The new value will be written to the sensor after saving the changes.

Press  $\leftarrow$  to exit the menu Cal Interval.

## 8.5 General Alarm

PATH:  \ CONFIG \ General Alarm

See the following explanation to get more details about the different settings for General Alarm.



Press the button **Event** in the line of the settings for **Option** and select the events, that should be considered for an alarm.

To activate a relay if the defined conditions are reached press the input field in the line for the settings of **Relay**. Only relay 1 can be assigned to general alarm. For general alarms the operation mode of the assigned relay is always inverted.

Enter the **Delay** time in seconds. A time delay requires the setpoint to be exceeded continuously for the specified length of time before activating the relay. If the condition disappears before the delay period is over, the relay will not be activated.

## 8.6 ISM / Sensor Alarm

PATH:  \ CONFIG \ ISM / Sensor Alarm

See the following explanation to get more details about the different settings for ISM / Sensor Alarm.



Select the channel by pressing the related button in the line of the settings for **Option**.

Depending on the selected channel or assigned sensor the **Events** that will be considered for generating an alarm can be selected. Some alarms will be considered in any case and not have to be selected or deactivated.

To select the desired relay that will be activated if an event has taken place press the input field in the line for the settings for **Relay**.

The operation mode of the relay can be defined.


Relay contacts are in normal mode until one of the selected events has taken place. Then the relay is activated and the contact states change. Select **Inverted** to reverse the normal operating state of the relay (i.e. normally open contacts are in an open state, and normally closed contacts are in a closed state if an event has taken place).

Enter the **Delay** time in seconds. A time delay requires the event to be occurred continuously for the specified length of time before activating the relay. If the condition disappears before the delay period is over, the relay will not be activated.

## 8.7 Clean

PATH:  \ CONFIG \ Clean

See the following explanation to get more details about the different settings for Clean



Enter the cleaning **Interval** time in hours. The cleaning interval can be set from 0.000 to 99999 hours. Setting it to 0 turns the clean cycle off.

Enter the **Clean Time** in seconds. The clean time can be 0 to 9999 seconds and must be smaller than the cleaning interval.


**Assign** the channel(s) for cleaning cycles. The assigned channels will be in HOLD state during the cleaning cycle.

Choose a **Relay**. Relay contacts are in normal mode until the cleaning cycle starts, then the relay is activated and the contact states change. Select Inverted to reverse the normal operating state of the relay (i.e. normally open contacts are in a open state, and normally closed contacts are in a closed state when the cleaning cycle starts).

## 8.8 Display Setup

PATH:  \ CONFIG \ Display Setup

See the following explanation to get more details about the different settings for Display Setup



Enter the name for the transmitter M800 (**Instrument Tag**). The instrument tag will also be displayed on the line at the top of the Start Screen and Menu Screen.

Choose with parameter **BackLight** to switched off or dim the transmitter screen after a defined time period without interaction. The transmitter screen will automatically come back after pressing the display.

Enter the **Light Time** in minutes. The light time is the time period without interaction before the transmitter screen will be dimmed or switched off.

**NOTE:** In case of an unacknowledged warning or alarm the transmitter screen will not be dimmed or switched off even if the light time has been elapsed


The parameter **Max** allows the setting of the backlight during operation. With the parameter **Dim** the backlight of the transmitter screen during the dimmed state can be adjusted. Press the + or - buttons in the corresponding line to adjust the parameters.



## 8.9 Digital Inputs

PATH:  \ CONFIG \ Digital Inputs

See the following explanation to get more details about the different settings for the digital inputs



Press the related button for the assignment of the **Channel** (Chan\_). Select the channel, which has to be linked to the digital input signal.

Press the input field in the line of the setting for **Mode** and select the impact of an active digital input signal. Choose 'Hold' to lead the assigned channel in hold state. For flow sensor the digital input signal can be used to reset the totalized flow value (Reset T-flow) for the channel. If an optical DO sensors is connected, the digital input signal can be used for LED controlling.

Press the related button for the assignment of the **Digital Inputs** (#1 for DI1, #2 for DI2 etc.) and select the digital input signal, which has to be linked to the channel.


An additional setting can be done, if an digital input signal has been selected.

Press the input field in the line for the setting of the **State** and select if the digital input is active at high or low level of the voltage input signal.

## 8.10 System

PATH:  \ CONFIG \ System

See the following explanation to get more details about the different settings for the System.



Select the desired **Language**. The following languages are available: English, French, German, Italian, Spanish, Portuguese, Russian, Chinese, Korean or Japanese.

Enter **Date&Time**.

The automatic change-over from summertime to wintertime and vice-versa frees the users from having to correct the time twice a year.

The winter to summer time-change are carried out automatically using the 12-month clock integrated in the transmitter. The date for the time-change can be set with the parameter **Summer**.

Provided it is a Sunday, the time-change would take place on the day that equates with the value, otherwise on the following Sunday. The winter/summer time-change takes place at 02:00 h.

The summer to winter time-change are carried out automatically using the 12-month clock integrated in the transmitter. The date for the time-change can be set through the parameter **Winter**.

Provided it is a Sunday, the time-change would take place on the day that equates with the value, otherwise on the following Sunday. The winter/summer time-change takes place at 03:00 h.

The number of hours, the clock will be shifted through the winter to summer and summer to winter time-change can be chosen. Press the related button for the setting of the **Shift Hour**.

## 8.11 PID Controller

PATH:  \ CONFIG \ PID Controller

PID control is proportional, integral and derivative control action that can provide smooth regulation of a process. Before configuring the transmitter, the following process characteristics must be identified.

Identify the **control direction** of the process

– **Conductivity:**

Dilution – direct acting where increasing measurement produces increasing control output such as controlling the feed of low conductivity diluting water to rinse tanks, cooling towers or boilers

Concentrating – reverse acting where increasing measurement produces decreasing control output, such as controlling chemical feed to attain a desired concentration

– **Dissolved Oxygen:**

Deaeration – direct acting where increasing DO concentration produces increasing control output such as controlling the feed of a reducing agent to remove oxygen from boiler feedwater

Aeration – reverse acting where increasing DO concentration produces decreasing control output, such as controlling an aerator blower speed to maintain a desired DO concentration in fermentation or wastewater treatment

– **pH/ORP:**

Acid feed only – direct acting where increasing pH produces increasing control output, also for ORP reducing reagent feed

Base feed only – reverse acting where increasing pH produces decreasing control output, also for ORP oxidizing reagent feed

Both acid and base feed – direct and reverse acting

Identify the **control output type** based on the control device to be used:

Pulse frequency – used with pulse input metering pump


Pulse length – used with solenoid valve

Analog – used with current input device such as electric drive unit, analog input metering pump or current-to-pneumatic (I/P) converter for pneumatic control valve


Default control settings provide linear control, which is appropriate for conductivity, dissolved oxygen. Therefore, when configuring PID for these parameters (or simple pH control) ignore settings of deadband and corner points in the tuning parameter section below. The non-linear control settings are used for more difficult pH/ORP control situations.

If desired, identify the non-linearity of the pH/ORP process. Improved control can be obtained if the non-linearity is accommodated with an opposing non-linearity in the controller. A titration curve (graph of pH or ORP vs. reagent volume) made on a process sample provides the best information. There is often a very high process gain or sensitivity near the setpoint and decreasing gain further away from the setpoint. To counteract this, the instrument allows for adjustable non-linear control with settings of a deadband around the setpoint, corner points further out and proportional limits at the ends of control as shown in the figure below.

Determine the appropriate settings for each of these control parameters based on the shape of the pH process titration curve.



See the following explanation to get more details about the different settings for PID Controller.



The M800 provides to 2 PID controllers. Press the input field in the line of the setting for **PID** and select the desired PID controller for configuration through pressing the button #1 for PID controller 1 and #2 for PID controller 2.

Press the related button for the assignment of the channel (**Chan**). Select the channel, which has to be linked to the PID Controller. To deactivate the PID controller press None.

Press the button for the assignment of the measuring parameter – based on the selected channel – that has be linked to the PID controller. Chose the measuring parameter through pressing the according field. Mx in the display indicates the measurement assigned to the PID Controller. (see chapter 8.1.1 “Channel Setup”).

The M800 offers the display of control output (%PID) of the PID controller in the Start Screen and Menu Screen. Press the related button for **Display For** and select the line, the control output should be displayed through pressing the corresponding field.



**NOTE:** The control output of the PID controller will be displayed instead of the measurement, that has been defined to be shown in the corresponding line (see chapter 8.1.1 “Channel Setup”).


Select with the parameter **PID Hold** the state of the control output for the PID controller if the transmitter M800 is in hold mode. Off means that the control output will be 0%PID if the transmitter is in hold mode. If Last Value has been chosen, the value for the control output signal before the transmitter went into hold mode will be used.

The parameter **PID A/M** allows selection of auto or manual operation for the PID controller. If auto has been chosen, the transmitter calculates the output signal based on the measured value and the settings of the parameters for the PID controller. In case of manual the transmitter shows in the Menu Screen at the line where the output signal is displayed two additional arrow buttons. Press the arrows buttons to increase or decrease the PID output signal.



**NOTE:** If Manual has been chosen the values for the time constants, gain, corner points, proportional limits, setpoint and deadband do not have any influence on the output signal.

Additional settings can be done by navigating to the next page of the menu.



The **PID Mode** assigns a relay or analog output for PID control action. Based on the control device being used, select one of the three options Relay PL, Relay PF and Aout through pressing the corresponding field

- Relay PL: If using a solenoid valve, select Relays PL (Pulse Length).
- Relay PF: If using a pulse input metering pump, select Relays PF (Pulse Frequency)
- Aout: For using an analog control select Aout.

Link the output signal **Out1,2** of the PID controller to the desired output of the transmitter. Press the related button for Out 1 and Out 2 and select the corresponding number for the output through pressing the according field. #1 means relay 1 or Aout 1, #2 means relay 2 our Aout 2 etc.



**NOTE:** Take care if reed type relays are linked to the controlling function. The reed type relays could be used for pulse frequency control devices and light duty applications. The current is limited to 0.5 amps and 10 watts (see also chapter 14.2 "Electrical Specifications"). Do not connect to this relays higher current devices.

If the PID Mode is set to Relay PL, the Puls Length for the output signal of the transmitter can be adjusted. Press the button for **Pulse Length** and the M800 displays a keypad for modifying the value. Enter the new value in the unit seconds according to the table below and press **←**.



**NOTE:** A longer pulse length will reduce wear on the solenoid valve. The % "on" time in the cycle is proportional to the control output.

	1 <sup>st</sup> Relay Position (Out 1)	2 <sup>nd</sup> Relay Position (Out 2)	Pulse Length (PL)
Conductivity	Controlling concentrating reagent feed	Controlling dilution water	Short (PL) provides more uniform feed. Suggested start point = 30 sec
pH/ORP	Feeding base	Feeding acid	Reagent addition cycle: short PL provides more uniform addition of reagent. Suggested start point = 10 sec
Dissolved Oxygen	Reverse control action	Direct acting control action	Feed cycle time: short PL provides more uniform feed. Suggested start point = 30 sec

If the PID Mode is set to Relay PF, the Puls Frequency for the output signal of the transmitter can be adjusted. Press the button for **Pulse Freq** and enter the new value in the unit pulse / minute according to the table below.



**NOTE:** Set the pulse frequency to the maximum frequency allowed for the particular pump being used, typically 60 to 100 pulses/minute. Control action will produce this frequency at 100% output.



**CAUTION:** Setting the pulse frequency too high may cause the pump to overheat.



	1 <sup>st</sup> Relay Position = #3	2 <sup>nd</sup> Relay Position = #4	Pulse Frequency (PF)
Conductivity	Controlling concentrating chemical feed	Controlling dilution water	Max allowed for the pump used (typically 60–100 pulses/minute)
pH/ORP	Feeding base	Feeding acid	Max allowed for the pump used (typically 60–100 pulses/minute)
Dissolved Oxygen	Reverse control action	Direct acting control action	Max allowed for the pump used (typically 60–100 pulses/minute)

If the PID Mode is set to **Aout**, the type for the analogue output signal of the transmitter can be selected. Press the corresponding button and choose between 4 to 20 mA and 0 to 20 mA for the output signal through pressing the according field.

For the assignment of the analogue output signal consider the table below.

	1 <sup>st</sup> Analogout Position = Out 1	2 <sup>nd</sup> Analogout Position = Out 2
Conductivity	Controlling concentrating chemical feed	Controlling dilution water
pH/ORP	Feeding base	Feeding acid
Dissolved Oxygen	Reverse control action	Direct acting control action


Press the input field for the parameter **Gain** to enter the gain of the PID controller as a unitless value. Gain represents the maximum value of the output signal of the PID controller in per cent (value 1 corresponds to 100%).

Press the corresponding input field in the line of **min** to adjust the Parameter integral or reset time **Tr** (left button) and/or rate of derivate time **Td** (right button).



**NOTE:** Gain, integral and derivate time are usually adjusted later by trial and error on process response. It is recommended to start with the value  $Td = 0$ .

Further settings can be done by navigating to the next page of the menu.



The display shows PID controller curve with input buttons for the corner points, setpoint and proportional limit for 100%.

Press the button **CP** to enter the menu for adjusting the corner points.

Page 1 shows the Corner Limit Low settings. Press the corresponding button to modify the value for the process parameter and the related output signal in %.

Browse to page 2 and the Corner Limit High settings are displayed. Press the corresponding button to modify the value for the process parameter and the related output signal in %.

Press the button **SP** to enter the menu for adjusting the setpoint and the dead band.

Press the button **Lim** to enter the menu for adjusting the proportional limit high and the proportional limit low, the range over which control action is required.

## 8.12 Service

PATH:  \ CONFIG \ Service

This menu is a valuable tool for troubleshooting and provides diagnostic functionality for the following items: Calibrate TouchPad, Set Analog Outputs, Read Analog Outputs, Read Analog Inputs, Set Relays, Read Relays, Read Digital Inputs, Memory, Display and optical DO sensors.



Select through the parameter **System** the desired item for diagnostic by pressing the according field.

Select through **Chan** the channel for diagnostic information of the optical DO sensor. This menu is only displayed if an optical DO sensor is connected.

The provided diagnostic functionality can now be called up through pressing the button **Diagnostic**.

### 8.12.1 Set Analog Outputs

The menu enables the user to set all analog outputs to any mA value within the 0–22 mA range. Use the + and – button to adjust the mA output signal. The transmitter will adjust the output signals according to the measurement and configuration of the analog output signals.

### 8.12.2 Read Analog Outputs

The menu shows the mA value of the analog outputs.

### 8.12.3 Read Analog Inputs

The menu shows the mA value of the analog input signals.

### 8.12.4 Set Relay

The menu allows the user to open or close each relay manually. If the menu will be exited, the transmitter will switch the relay according configuration.

### 8.12.5 Read Relay

The menu shows state of every relay. On indicates the relay is closed, Off indicates that the relay is open.

### 8.12.6 Read Digital Inputs

The menu shows the state of the digital input signals.

## 8.12.7 Memory

If Memory is selected the transmitter will perform a memory test of all connected transmitter boards and ISM sensors.

## 8.12.8 Display

The transmitter shows every 5 seconds red, green, blue, grey and dark grey display and returns afterwards to the menu Service. If within the 5 seconds for every color the screen will be pressed the transmitter goes to the next step.

## 8.12.9 Calibrate TouchPad

Press during the 4 calibrations steps always the center of the shown circle in the 4 corner of the display. The transmitter shows the calibration result.

## 8.12.10 Channel Diagnostic

If errors has been occurred on the optical DO sensor, the corresponding messages are displayed.

## 8.13 Technical Service

PATH:  \ CONFIG \ Technical Service

With this menu the calibration factors for the analog input and output signals can be shown.




Select through the parameter **Options** the signal(s), the calibration factors should be displayed for.

## 8.14 User Management

PATH:  \ CONFIG \ User Management

This menu allows for the configuration of different user and administrator passwords, as well as setting up a list of allowed menus for the different users. The administrator has rights to access all menus. All default passwords for new transmitters are "00000".



Press the input field in the line of **Protection** and select the desired kind of protection. Available are the following options

**Off:** no protection

**Active:** Activation of the Menu Screen (see chapter 3.2.2 "Activation Menu Screen") has to be confirmed

**Password:** Activation of the Menu Screen is only possible with a password

Press the according button for **Option** to select the profile for the administrator (Admin) or one of the users.

**NOTE:** The administrator has always the rights to access all menus. For different users the access rights can be defined.

Press the input button for **UserID** to enter the name for the user or administrator. The name for the user or administrator will be displayed if the protection via password is selected for activation of the Menu Screen.

For changing the password of the selected user or administrator press the input field for **Password**. Enter the old password in the field Old PW, the new one in the field New PW and confirm it in the field confirm PW. The default password is "00000000" for the administrator and all users.

If the profile for an user has been selected an additional input field to define the access rights will be displayed.

To assign access rights the according button for the menu has to be pressed. In case of an assignment of the access rights,  is displayed in the related button.

## 8.15 Reset

PATH:  \ CONFIG \ Reset

Depending on the transmitter version and configuration different options for a reset are available.

See the following explanation to get more details about the different option to reset data and / or configurations.

### 8.15.1 System Reset

This menu allows the reset of the transmitter M800 to the factory default settings (setpoints off, analog outputs off, passwords, etc.). Furthermore the calibration factors for analog in- and outputs, meter etc. can be set to the last factory values.

Press the input field for **Options** and select System.

Press the input field for **Items** (Configure button) and select the different parts of the configuration that will be reset.

If an item has been selected the Action menu is displayed. Press the Reset button.

## 8.15.2 Reset Sensor Calibration for Optical DO Sensors

If an optical oxygen sensor is connected to the transmitter, a menu is available that allows the reset of the calibration data of the sensor to the factory settings.

Press the input field for **Options** and select the channel the optical DO sensor is connected to.

Press the input field for **Items** (Configure button). Select SensorCal to Factory through pressing the according button.

If SensorCal to Factory has been selected the Action menu is displayed. Press the Reset button.



**NOTE:** Through a reset of the calibration data the Adaptive Calibration Timer (see chapter 9.1 "iMonitor") will set to 0.



**NOTE:** To ensure best measuring results, a new calibration of the sensor is recommended after a reset of the calibration data to factory settings. Depending on the application and sensor, the calibration should be performed as a one point calibration or two point calibration (see chapter 7.7 "Calibration of Optical Oxygen Sensors").

## 8.15.3 Reset Sensor Calibration for UniCond2e Sensors

For UniCond2e sensors, the SensorCal (sensor calibration) and ElecCal (sensor electronics calibration) can be restored to factory settings.

Press the input field for **Options** and select the channel the UniCond2e sensor is connected to.

Press the input field for **Item** (Configure button). Select SensorCal to Factory and/or ElecCal to Factory by checking the adjacent box. Press **↵** enter to accept the value.

If an item has been selected the Action menu is displayed. Press the Reset button.

The M800 will bring up the confirmation dialog. Select Yes and the reset will be executed. Press No to go back to menu Reset without performing the reset.

## 8.15.4 Reset Total Flow

For models of the M800 transmitter that accept flow signals the totalized flow value for each channel can be reset.

Press the input field for **Options** and select the channel for which the totalized flow value should be reset.

Press the input field for **Items** (Configure button). Select Total Flow by checking the adjacent box. Press **↵** to accept the value.

If Total Flow has been selected the Action menu is displayed. Press the Reset button.

The M800 will bring up the confirmation dialog. Select Yes and the reset will be executed. Press No to go back to menu Reset without performing the reset.

## 8.15.5 Reset for CO2 hi Measurement

If a thermal conductivity dissolved CO2 sensor is connected to the transmitter, a menu is available that allows the reset of the measurement circuit from the sensor.

Under the circumstances the sensor detects an error, the sensor will run into the sensor protection mode. The electronic measurement circuit will be shut down for sensor protection and needs to be restarted after the remedy of the failure for accurate CO2 measurement.

Press the input field for **Options** and select the channel with the CO2 sensor, which should be reset.

Press the input field for **Items** (Configure button). Select CO2 Measurement by checking the adjacent box. Press  $\leftarrow$  to accept the value.

If CO2 Measurement has been selected the Action menu is displayed. Press the Reset button.

The M800 will bring up the confirmation dialog. Select Yes and the reset will be executed. Press No to go back to menu Reset without performing the reset.

## 8.16 RS485 Printer Output Configuration

PATH:  \ CONFIG \ Printer

The Printer menu option allows configuring the M800 RS485 output to send data to a suitable printer. The printer output may be configured to print up to 6 configure measurements on separate lines, for each available sensor input, including pulsed input channels. At each print cycle, the output will include a header line with data and time based on the M800 internal clock, and one line for each configured measurement including channel, measurement descriptor, measurement value and unit of measure.

The output will appear as follows:

11/May/2012 15:36

Ch Label Measurement

1 CHAN\_1 302 ppbC

2 CHAN\_2 0.54 uS/cm


3 CHAN\_3 7.15 pH

To configure the printer output, access the printer menu and configure the following options:

**Lines to Print** will configure the number of measurements that will be printed for each print cycle. Enter the total number of measurements to be configured for output.

**Output Time** defines the time in minutes between each print cycle. Output time may be set from 1 to 1000 minutes.





Once the output time and print lines has been established, press the Configure button to format the printer output. The number at the left of the window shows the order in which the lines will appear on the printer output. From the first dropdown, select the channel with which the desired sensor is connected. This dropdown will list the labels associated with each channel as configured under Channel Setup. Usint the second dropdown, select the unit associated with the measurement to be displayed. Note that if more the 4 lines of output has been selected, use the < and > icons to navigate through the pages to be configured.

## 8.17 USB Measurement Interface

The user may access measurement values via the USB. The user provides a command and the M800 responds using the following format.

Command: [0x02][0x02]"Dx"(x is the channel index: 1-6)

Response: "XXXXXXXXuuuuuu XXXXXXXXuuuuuu XXXXXXXXuuuuuu XXXXXXXXuuuuuu"<cr>



**NOTE:** The first instance of 0x02 is the ID for M800, which must be 0x02 only. The second instance of 0x02 is the length, which must be 0x02 only. The response provides M1~M4 only. XXXXXXXX is measurement floating value in ASCII.

uuuuuu is the unit in ASCII, if current unit is less than 6 characters, the format is right aligned, e.g. if unit is pH, response " pH".

<cr> means carriage return (0x0D, 0x0A)

If the sending command is not correct an error message is generated.

Error response format: "ERROR #xx"

xx is the error code


01: Invalid opcode ---- if it is not D.

02: Parameter error ---- if x is not 1-6

07: Length error ---- if length is not 2.

## 9 ISM


PATH:  \ ISM



### 9.1 iMonitor

PATH:  \ ISM \ iMonitor

The iMonitor gives an overview of the current state of the complete loop at a glance.



The iMonitor of the first channel is displayed on the screen. To browse through the iMonitor for the different channels press > at line at the bottom of the display.

The values DLI, TTM and ACT as well as TTCal in combination with UniCond2e sensors are shown as bar graph. If the values falls below 20% of the initial value the bar graph changes from green to yellow color. If the value falls below 10% the color changes to red.

For Cond4e sensors the days in operation of the sensor are displayed.

Furthermore SIP-, CIP-, AutoClave-, SAN-cycles as well as the values for Rg and Rref can be displayed and assigned to a colored button if the values are provided by the sensor.

The color for the related button of SIP-, CIP-, Autoclave- and SAN-cycles will change from green to yellow if less then 20% of the defined maximum quantity for the cycle remain and to red if less then 10% remain. For configuration of the maximum quantity see chapter 8.4 "ISM Setup"

The buttons for Rg and Rref change to yellow if the conditions for a warning messages are fulfilled and to red if the conditions for a alarm message are fulfilled. The buttons remain grey if the corresponding ISM alarm is not configured (see chapter 8.6 "ISM / Sensor Alarm").

Depending on the measured parameter (connected sensor) the following data are available in the menu iMonitor:

pH: DLI, TTM, ACT, CIP, AutoClave, SIP\*, Rg\*\*, Rref\*\*  
 Amperometric O2: DLI, TTM, ACT, CIP, AutoClave, SIP\*, Electrolyte\*\*\*  
 Optical O2: DLI, ACT, CIP, AutoClave, SIP\*  
 O3: DLI, TTM, ACT, SAN  
 Cond: Days in operation, TTCal\*\*\*\*, CIP, SIP


- \* if AutoClave has not been activated (see chapter 8.6 "ISM / Sensor Alarm" )
- \*\* if the alarm for Rg and/or Rref has been activated (see chapter 8.6 "ISM / Sensor Alarm")
- \*\*\* if the alarm for Electrolyte Level Error has been activated (see chapter 8.6 "ISM / Sensor Alarm")
- \*\*\*\* if UniCond2e sensor is connected



## 9.2 Messages

PATH:  \ ISM \ Messages

The messages for occurred warnings and alarms are listed in this menu. Up to 100 entries will be listed.



5 messages per page are listed. If more than 5 messages are available additional pages can be accessed.

Not acknowledged alarms or warning will be listed at the beginning. Then the acknowledged but still existing alarm or warning are listed. At the end of the list the already solved warning and alarms are described. Between this groups the messages are listed chronological.

The state of the warning or alarm is indicated through the following signs:

Red button blinking	Alarm exists and has not been acknowledged
Red button not blinking	Alarm exists and has been acknowledged
Yellow button blinking	Warning exists and has not been acknowledged
Yellow button not blinking	Warning exists and has been acknowledged
Grey button not blinking	Warning or alarm has been solved

A unacknowledged warning or alarm will be acknowledged by pressing the **Info** button in the corresponding line.

For every message the corresponding **Info** button can be pressed. Message information, date and time the warning or alarm has been occurred and the status of the alarm or message are displayed.

If warning or alarm has already been solved the pull up window for the message shows an additional button to clear the message i.e. to delete it from the message list.


## 9.3 ISM Diagnostics

PATH:  \ ISM \ ISM Diagnostics

The transmitter M800 provides for all ISM sensors a diagnostic menu. Access the menu Channel and select the channel by pressing the related input field.

Depending on the selected channel and assigned sensor different diagnostic menus are displayed. See the following explanation to get more details about the different diagnostic menus.

### 9.3.1 pH/ORP, Oxygen, O3 and Cond4e Sensors



If an pH/ORP, oxygen, O3 or Cond4e sensor is connected to the selected channel, the diagnostic menus cycles, sensor monitor and max. temperature are available.

Press the **Cycle** button and the information for CIP, SIP and Autoclave cycles of the connected sensor are displayed. The displayed information shows the amount of cycles the sensor has been exposed and the max. limitation for the corresponding cycle as defined in the menu 'ISM Setup' (see chapter 8.4 "ISM Setup").



**NOTE:** For Cond4e and optical DO sensors, which are not autoclavable the menu AutoClave Cycles is not displayed.



**NOTE:** For O3 sensors the SAN cycles are displayed.

Press the **Sensor Monitor** button and the information for DLI, TTM and ACT of the connected sensor are displayed. The values DLI, TTM and ACT are shown as bar graph. If the values falls below 20% of the initial value the bar graph changes from green to yellow color. If the value falls below 10% the color changes to red.



**NOTE:** For Optical DO sensors TTM does not exist.




**NOTE:** For Cond4e sensors the operating hours are displayed.

Press the **Max. Temperature** button and the information about the maximum temperature, that the connected sensor has ever seen, together with a time stamp of this maximum is displayed. This value is stored on the sensor and cannot be changed. During autoclaving the max. temperature is not recorded.



**NOTE:** For Optical DO sensors the max. temperature of the board and of the spot are displayed.

### 9.3.2 UniCond2e Sensors



For UniCond2e sensors, the following diagnostic items can be viewed: Excursion Counters including High Temp and High Conductivity, Highest Measured including Highest Temp and Highest Cond, Cycles including CIP cycles and SIP Cycles.

## 9.4 Calibration Data

PATH:  \ ISM \ Calibration Data

The transmitter M800 provides for all ISM sensors a calibration history. Depending on the selected channel and assigned sensor different data of the calibration history are available.

See the following explanation to get more details about the different data of the calibration history.

### 9.4.1 Calibration Data for All ISM Sensors excluding UniCond2e



If an ISM sensor - excluding UniCond2e - is connected to the selected channel between the calibration data set of

**Actual** (Actual adjustment): This is the actual calibration dataset which is used for the measurement. This dataset moves to Cal1 position after the next adjustment.

**Factory** (Factory calibration): This is the original dataset, determined in the factory. This dataset remains stored in the sensor for reference and cannot be overwritten.

**1.Adjust** (First adjustment): This is the first adjustment after the factory calibration. This dataset remains stored in the sensor for reference and cannot be overwritten

**Cal1** (last calibration/adjustment): This is the last executed calibration/adjustment data set. This dataset moves to Cal2 and then to Cal3 when a new calibration/adjustment is performed. Afterwards, the dataset is not available anymore. Cal2 and Cal3 acting in the same way as Cal1.

**Cal2** and **Cal3** can be chosen. For the selection of the calibration data set press the corresponding field.



**NOTE:** The amperometric oxygen sensor of THORNTON and the O3 sensor do not provide the data set Cal1, Cal2, Cal3 and 1.Adjust.

Press the **Cal Data** button and the corresponding calibration data set is displayed. Furthermore the time stamp for the calibration and the User ID is listed.



**NOTE:** This function requires the correct setting of date and time during calibration and / or adjustment tasks (see chapter 8.10 "System").

## 9.4.2 Calibration Data for UniCond2e Sensors



For UniCond2e sensors the following three sets of calibration data may be selected:

**Actual (Actual calibration):** This is the actual calibration dataset which is used for the measurement.

**Factory (Factory calibration):** This is the original dataset, determined in the factory. This dataset remains stored in the sensor for reference and cannot be overwritten.

**Cal1 (last calibration/adjustment):** This is the last executed calibration/adjustment data set.

Press the Cal Data button and the corresponding calibration data set is displayed.

If the data set of the actual calibration has been chosen, on page 1, the date and time of the calibration, User ID, conductivity calibration constants, and reference conductivity values to calibrate are displayed. On page 2 the As-found conductivity values and the deviation from the reference are shown. On page 3 and 4 the same information for temperature is displayed. On page 5 the calibration cycles applied to the sensor and the next calibration date for conductivity (C) and temperature (T) are displayed.

If the data set of the factory calibration has been chosen, on page 1, the date and time of the calibration, the conductivity calibration constants, and reference conductivity values used to calibrate are displayed. On page 2, the same values for temperature are shown.

Press  $\leftarrow$  to exit the menu Cal Data.



**NOTE:** This function requires the correct setting of date and time during calibration and / or adjustment tasks (see chapter 8.10 "System").

## 9.5 Sensor Info

PATH:  $\leftarrow$  \ ISM \ Sensor Info

The model, hardware and software version, last calibration date as well as the product and serial number of the ISM sensors, that are connected to the M800 transmitter can be displayed on the screen.

Enter Sensor Info.



The data of the first channel, a sensor is connected, are displayed on the screen. Press the input field in the line of Chan. To get the data of the desired sensor select the corresponding channel through pressing the according field.

The data Model, Cal Date (date of last adjustment), S/N (serial number), P/N (product number), SW Ver (software version) and HW Ver (hardware version) of the select sensor are displayed.

**NOTE:** If a UniCond2e sensor is connected the following data is also displayed, Temp Sens. (temperature sensor) Electrode (electrode material), Body/Ins Mat: (body and/or insulator material), Inner: (inner electrode material), Outer (outer electrode material) Fitting: (fitting material), Class VI (FDA Class VI material).

To exit the menu Sensor Info press  $\leftarrow$ . To return to the Menu Screen press  $\leftarrow$ .

## 9.6 HW / SW Version

PATH:  \ ISM \ HW/SW Version

The hardware and software version as well as the product number and serial number of the M800 transmitter itself or the different boards, that are plugged can be displayed on the screen.




The data of the transmitter are displayed on the screen. Press the input field in the line of **M800**. To select the data of the desired board or the transmitter itself press the corresponding field.

The data S/N (serial number), P/N (product number), SW Ver (software version) and HW Ver (hardware version) of the select board or transmitter are displayed.

## 9.7 Log Book

PATH:  \ ISM \ Sensor Info

The transmitter M800 provides a log book with 250 entries. The log book is managed as an ring buffer, i.e. entry 251 causes the erasing of entry no. 1 etc..



The entries show time stamp and action.

## 10 Wizards

PATH:  \ WIZARD

The transmitter M800 allows to set up to 4 wizards / favorites to ensure a quick access for often used functions.

### 10.1 Set Wizard

PATH:  \ WIZARD \ Set Wizard



The main menus are displayed. Choose the menu, that contains the function, which should be defined as a wizard (favorite), e.g. ISM through pressing the corresponding arrow ► in the same line.

Choose the function, that should be set as an wizard by pressing the according button. A function, which is set as a wizard shows ★ icon.

**NOTE:** To delete the link to wizards, press the according button for the function. The wizard ★ icon is not shown any more.

### 10.2 Access to Wizards

Access the menu Set Wizards. The wizards defined are listed on this page. Press the corresponding arrow ► for the function in the same line.

## **11 Maintenance**

### **11.1 Front panel cleaning**

Clean the front panel with a damp soft cloth (water only, no solvents). Gently wipe the surface and dry with a soft cloth.

## 12 Troubleshooting

If the equipment is used in a manner not specified by Mettler-Toledo, the protection provided by the equipment may be impaired.

Review the table below for possible causes of common problems:

Problem	Possible Cause
Display is blank.	<ul style="list-style-type: none"> <li>– No power to M800.</li> <li>– Blown fuse.</li> <li>– Hardware failure.</li> </ul>
Incorrect measurement readings.	<ul style="list-style-type: none"> <li>– Sensor improperly installed.</li> <li>– Incorrect units multiplier entered.</li> <li>– Temperature compensation incorrectly set or disabled.</li> <li>– Sensor or transmitter needs calibration.</li> <li>– Sensor or patch cord defective or exceeds recommended maximum length.</li> <li>– Hardware failure.</li> </ul>
Measurement readings not stable.	<ul style="list-style-type: none"> <li>– Sensors or cables installed too close to equipment that generates high level of electrical noise.</li> <li>– Recommended cable length exceeded.</li> <li>– Averaging set too low.</li> <li>– Sensor or patch cord defective.</li> </ul>
Displayed red or yellow bar graph is flashing.	<ul style="list-style-type: none"> <li>– Setpoint is in alarm condition (setpoint exceeded).</li> <li>– Alarm has been selected (see chapter 8.6 "ISM / Sensor Alarm") and occurred.</li> </ul>
Cannot change menu settings.	<ul style="list-style-type: none"> <li>– User locked out for security reasons.</li> </ul>

### 12.1 Warning- and Alarm Indication

#### 12.1.1 Warning Indication



It will be indicated through a yellow bar graph on the display if there are conditions, that have generate a warning. If the corresponding channel is shown on the current Menu Screen or Start Screen (see chapter 3.2 "Display") the yellow bar graph is displayed in the line with the name of the channel. A warning message will be recorded and can be selected through the menu Messages (PATH: \ISM\Messages; see also chapter 9.2 "Messages").



If a channel, that is not shown on the current Menu Screen or Start Screen is concerned by a warning, a yellow bar graph is displayed on the head line of the display. A warning message will be recorded and can be selected through the menu Messages (PATH: \ISM\Messages; see also chapter 9.2 "Messages").



**NOTE:** If the warning has not been acknowledged the bar graph is blinking. If the warning has already been acknowledged, the bar graph will be displayed continuously. See also chapter 9.2 "Messages". In case of an unacknowledged warning or alarm the transmitter screen will not be dimmed or switched off even if the light time has been elapsed (see chapter 8.8 "Display Setup")



**NOTE:** If at the same time for a channel an alarm and a warning should be indicated, the indication of the alarm will have higher priority. The alarm will be indicated (see chapter 12.1 "Warning and Alarm Indication") on the Menu Screen or Start Screen, while the warning will not be shown.





Pressing the yellow bar graph on the Menu Screen will lead to the Messages. Refer to chapter 9.2 "Messages" for the description of the functionality for this menu.

**NOTE:** The detection of some warnings can be activated/deactivated through (de)activating the corresponding alarm. Refer therefore to chapter 8.6 "ISM / Sensor Alarm"

### 12.1.2 Alarm Indication



It will be indicated through a red bar graph on the display if there are conditions, that have generate an alarm. If the corresponding channel is shown on the current Menu Screen or Start Screen (see chapter 3.2 "Display") the red bar graph is displayed in the line with the name of the channel. An alarm message will be recorded and can be selected through the menu Messages (PATH:  $\text{ISM}\backslash\text{Messages}$ ; see also chapter 9.2 "Messages").




If a channel, that is not shown on the current Menu Screen or Start Screen is concerned by an alarm, a red bar graph is displayed on the head line of the display. An alarm message will be recorded and can be selected through the menu Messages (PATH:  $\text{ISM}\backslash\text{Messages}$ ; see also chapter 9.2 "Messages").



**NOTE:** If the alarm has not been acknowledged the bar graph is blinking. If the alarm has already been acknowledged, the bar graph will be displayed continuously. See also chapter 9.2 "Messages". In case of an unacknowledged warning or alarm the transmitter screen will not be dimmed or switched off even if the light time has been elapsed (see chapter 8.8 "Display Setup")



**NOTE:** If at the same time for a channel an alarm and a warning should be indicated, the indication of the alarm will have higher priority. The alarm will be indicated (see chapter 12.1 "Warning and Alarm Indication") on the Menu Screen or Start Screen, while the warning will not be shown.



Pressing the red bar graph on the Menu Screen will lead to the Messages. Refer to chapter 9.2 "Messages" for the description of the functionality for this menu.



**NOTE:** The detection of some alarms can be activated/deactivated. Refer therefore to chapter 8.6 "ISM / Sensor Alarm"



**NOTE:** Alarms which are caused by a violation of the limitation of a setpoint or the range (PATH:  $\text{CONFIG}\backslash\text{Set Points}$ ; see also chapter 8.3 "Set Points") will also be indicated on the display and recorded through the menu Messages (PATH:  $\text{ISM}\backslash\text{Messages}$ ; see also chapter 9.2 "Messages").

## 13 Accessories and Spare Parts

Please contact your local Mettler-Toledo sales office or representative for details for additional accessories and spare parts.

<b>Description</b>	<b>Order no.</b>
Pipe Mount Kit for 1/2DIN models	52 500 212
Panel Mount Kit for 1/2DIN models	52 500 213
Protective Hood for 1/2DIN models	52 500 214

# 14 Specifications

## 14.1 General specifications

<b>Conductivity/resistive Specifications for Cond4e Sensors</b>	
Conductivity range	0.01 to 650 mS/cm (1.54 $\Omega$ x cm to 0.1 M $\Omega$ x cm)
Chemical concentration curves	NaCl: 0–26% @ 0 °C to 0–28% @ +100 °C NaOH: 0–12% @ 0 °C to 0–16% @ +40 °C to 0–6% @ +100 °C HCl: 0–18% @ –20 °C to 0–18% @ 0 °C to 0–5% @ +50 °C HNO <sub>3</sub> : 0–30% @ –20 °C to 0–30% @ 0 °C to 0–8% @ +50 °C H <sub>2</sub> SO <sub>4</sub> : 0–26% @ –12 °C to 0–26% @ +5 °C to 0–9% @ +100 °C H <sub>3</sub> PO <sub>4</sub> : 0–35% @ +5 °C to +80 °C User defined concentration curve table (5x5 matrix)
TDS ranges	NaCl, CaCO <sub>3</sub>
Sensor maximum distance	ISM: 80 m (260 ft)
Cond/Res accuracy	± 1 digit
Cond/Res repeatability	± 1 digit
Cond/Res resolution	auto/0.001/0.01/0.1/1 (can be selected)
Temperature range	–40 to +200.0 °C (–40 to 392 °F)
Temperature resolution	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	± 1 digit
Temperature repeatability	± 1 digit

<b>pH Specifications</b>	
pH range	–1.00 to 15.00 pH
Sensor maximum distance	80 m (260 ft)
pH resolution	auto/0.01/0.1/1 (can be selected)
pH accuracy	± 1 digit
mV range	–1500 to 1500 mV
mV resolution	auto/0.01/0.1/1 mV
mV accuracy	± 1 digit
Temperature range	–30 to 130 °C (–22 to 266 °F)
Temperature resolution	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	± 1 digit
Temperature repeatability	± 1 digit

<b>Available Buffer Sets:</b>
<b>Standard buffers</b>
MT-9 buffers, MT-10 buffers, NIST Technical Buffers, NIST Standard Buffers (DIN 19266:2000–01), JIS Z 8802 buffers, Hach buffers, CIBA (94) buffers, Merck Titrisols-Reidel Fixanals, WTW buffers
<b>Dual membrane electrodes pH buffers (pH/pNa)</b>
Mettler-pH/pNa buffers (Na+ 3.9M)

<b>Specifications for Amperometric Oxygen Sensors</b>	
Display range for current	0 to 9999 (depending on sensor) nA
Sensor maximum distance	80 m (260 ft)
DO concentration range	0.1 ppb (µg/l) to 50.00ppm (mg/l)
DO saturation range	0 to 500% air
O2 gas concentration range	0 to 9999 ppm O2 gas
O2 gas saturation range	0 to 100% O2 gas
DO accuracy	± 1 digit
O2 gas accuracy	± 1 digit
Resolution	auto/0.001/0.01/0.1/1, (can be selected)
Temperature measuring range	-30 to 150 °C (-22 to 302 °F)
Temperature resolution	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	± 1 digit
Temperature repeatability	± 1 digit

<b>Specifications for Optical Oxygen Sensors</b>	
Sensor maximum distance	50 m (164 ft)
DO concentration range	0.1 ppb (µg/l) to 50.00ppm (mg/l)
DO saturation range	0 to 500%
DO accuracy	± 1 digit
Resolution	auto/0.001/0.01/0.1/1, (can be selected)
Temperature measuring range	-30 to 150 °C (-22 to 302 °F)
Temperature resolution	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	± 1 digit
Temperature repeatability	± 1 digit

<b>Dissolved Carbon Dioxide Specifications</b>	
CO <sub>2</sub> measuring ranges	0 ... 5000 mg/l 0 ... 200%sat 0 ... 1500 mmHg 0 ... 2000 mbar 0 ... 2000 hPa
Sensor maximum distance	80 m (260 ft)
CO <sub>2</sub> accuracy	± 1 digit
CO <sub>2</sub> resolution	auto/0.001/0.01/0.1/1, (can be selected)
mV range	-1500 to 1500 mV
mV resolution	auto/0.01/0.1/1 mV
mV accuracy	± 1 digit
Total pressure range (TotPres)	0 ... 4000 mbar
Temperature measuring range	-30 to 150 °C (-22 to 302 °F)
Temperature resolution	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	-40 to + 200.0 °C (-40 to 392 °F)
Temperature repeatability	auto/0.001/0.01/0.1/1 K (°F), (can be selected)
Temperature accuracy	± 1 digit
Temperature repeatability	± 1 digit
<b>Available Buffer Set:</b>	
MT-9 buffers with solution pH = 7.00 and pH = 9.21 @ 25 °C	

## 14.2 Electrical specifications

Power requirements	100 to 240 V AC or 20 to 30 V DC, 10W, AWG 16–24, 0.2 mm <sup>2</sup> to 1.5 mm <sup>2</sup>
Frequency	50 to 60 Hz
Analog output signals	Eight 0/4 to 20 mA outputs, galvanically isolated from input and from earth/ground
Measurement error through analog outputs	< ±0.05 mA over 0 to 22 mA range
Analog output configuration	Linear, Bi-Linear, Logarithmic, Autoranging
Load	max. 500 Ω
Connection terminals	Spring cage terminals appropriate for AWG 16–24, 0.2 mm <sup>2</sup> to 1.5 mm <sup>2</sup> wires
Digital communication	USB port, Type B connector
PID process controller	2xPID; pulse length, pulse frequency or analog control
Cycle time	Ca. 1 second
Digital inputs	6 (5 for 2-channel version) with switching limits 0.00 VDC to 1.00 VDC for low level and 2.30 VDC to 30.00 VDC for high level
Analog inputs	One 4 to 20 mA input, galvanically isolated from other signals
Measurement error through analog input	< ±0.05 mA over 0 to 22 mA range
Mains power fuse	2.0 A slow blow type FC, not replaceable
Relays	4-SPST mechanical rated at 250 VAC, 3 Amps Relay1 NC, Relay2 to 4 NO 4-SPST Type Reed 250 VAC or DC, 0.5 A (Relay5 to 8)
Alarm Relay delay	0–999 s
User interface	Color touch screen 5.7" Resolution 320 x 240 pixel 256 colors
Max. cable length	80 m (260 ft) for pH, amp. oxygen, Cond4e, ozone 15 m (50 ft) for optical DO, UniCond2e



**NOTE:** This is a 4-wire-product with an active 4–20 mA analog output.  
Please do not supply to terminal no. 3 to 10 of TB1 and terminal no. 1 to 8 of TB3.

## 14.3 Mechanical specifications

Dimensions (housing – H x W x D)*	150 x 158 x 170 mm (5.36" x 6.22" x 6.69")
Front bezel – H x W	150 x 158 mm (5.36" x 6.22")
Max. D – panel mounted	125 mm (4.92")
Weight	1.6 kg (3.5 lb)
Material	Polycarbonate / PC
Ingress rating	IP 66 (when back cover is attached)

\* H = Height, W = Width, D = Depth

## 14.4 Environmental specifications

Storage temperature	-40 to 70 °C (-40 to 158 °F)
Ambient temperature operating range	-20 to 50 °C (-4 to 122 °F)
Relative humidity	0 to 95% non-condensing
Emissions	According to EN55011 Class A
Hazardous areas	cFMus Class I Division 2*,
Ratings / Approvals	CE Compliant
Altitude, maximum	2,700 m

\* in progress

## 15 Warranty

METTLER TOLEDO warrants this product to be free from significant deviations in material and workmanship for a period of one year from the date of purchase. If repair is necessary and not the result of abuse or misuse within the warranty period, please return by freight pre-paid and amendment will be made without any charge. METTLER TOLEDO's Customer Service Dept. will determine if the product problem is due to deviations or customer abuse. Out-of-warranty products will be repaired on an exchange basis at cost.

The above warranty is the only warranty made by METTLER TOLEDO and is lieu of all other warranties, expressed or implied, including, without limitation, implied warranties of merchantability and fitness for a particular purpose. METTLER TOLEDO shall not be liable for any loss, claim, expense or damage caused by, contributed to or arising out of the acts or omissions of the Buyer or Third Parties, whether negligent or otherwise. In no event shall METTLER TOLEDO's liability for any cause of action whatsoever exceed the cost of the item giving rise to the claim, whether based in contract, warranty, indemnity, or tort (including negligence).

## 16 Buffer tables

M800 transmitters have the ability to do automatic pH buffer recognition. The following tables show different buffers that are automatically recognized.

### 16.1 Standard pH buffers

#### 16.1.1 Mettler-9

Temp (°C)	pH of buffer solutions			
0	2.03	4.01	7.12	9.52
5	2.02	4.01	7.09	9.45
10	2.01	4.00	7.06	9.38
15	2.00	4.00	7.04	9.32
20	2.00	4.00	7.02	9.26
25	2.00	4.01	7.00	9.21
30	1.99	4.01	6.99	9.16
35	1.99	4.02	6.98	9.11
40	1.98	4.03	6.97	9.06
45	1.98	4.04	6.97	9.03
50	1.98	4.06	6.97	8.99
55	1.98	4.08	6.98	8.96
60	1.98	4.10	6.98	8.93
65	1.98	4.13	6.99	8.90
70	1.99	4.16	7.00	8.88
75	1.99	4.19	7.02	8.85
80	2.00	4.22	7.04	8.83
85	2.00	4.26	7.06	8.81
90	2.00	4.30	7.09	8.79
95	2.00	4.35	7.12	8.77



## 16.1.2 Mettler-10

Temp (°C)	pH of buffer solutions				
0	2.03	4.01	7.12	10.65	
5	2.02	4.01	7.09	10.52	
10	2.01	4.00	7.06	10.39	
15	2.00	4.00	7.04	10.26	
20	2.00	4.00	7.02	10.13	
25	2.00	4.01	7.00	10.00	
30	1.99	4.01	6.99	9.87	
35	1.99	4.02	6.98	9.74	
40	1.98	4.03	6.97	9.61	
45	1.98	4.04	6.97	9.48	
50	1.98	4.06	6.97	9.35	
55	1.98	4.08	6.98		
60	1.98	4.10	6.98		
65	1.99	4.13	6.99		
70	1.98	4.16	7.00		
75	1.99	4.19	7.02		
80	2.00	4.22	7.04		
85	2.00	4.26	7.06		
90	2.00	4.30	7.09		
95	2.00	4.35	7.12		

## 16.1.3 NIST Technical Buffers

Temp (°C)	pH of buffer solutions				
0	1.67	4.00	7.115	10.32	13.42
5	1.67	4.00	7.085	10.25	13.21
10	1.67	4.00	7.06	10.18	13.01
15	1.67	4.00	7.04	10.12	12.80
20	1.675	4.00	7.015	10.07	12.64
25	1.68	4.005	7.00	10.01	12.46
30	1.68	4.015	6.985	9.97	12.30
35	1.69	4.025	6.98	9.93	12.13
40	1.69	4.03	6.975	9.89	11.99
45	1.70	4.045	6.975	9.86	11.84
50	1.705	4.06	6.97	9.83	11.71
55	1.715	4.075	6.97		11.57
60	1.72	4.085	6.97		11.45
65	1.73	4.10	6.98		
70	1.74	4.13	6.99		
75	1.75	4.14	7.01		
80	1.765	4.16	7.03		
85	1.78	4.18	7.05		
90	1.79	4.21	7.08		
95	1.805	4.23	7.11		

### 16.1.4 NIST standard buffers (DIN and JIS 19266: 2000–01)

Temp (°C)	pH of buffer solutions			
0				
5	1.668	4.004	6.950	9.392
10	1.670	4.001	6.922	9.331
15	1.672	4.001	6.900	9.277
20	1.676	4.003	6.880	9.228
25	1.680	4.008	6.865	9.184
30	1.685	4.015	6.853	9.144
35	1.694	4.028	6.841	9.095
40	1.697	4.036	6.837	9.076
45	1.704	4.049	6.834	9.046
50	1.712	4.064	6.833	9.018
55	1.715	4.075	6.834	8.985
60	1.723	4.091	6.836	8.962
70	1.743	4.126	6.845	8.921
80	1.766	4.164	6.859	8.885
90	1.792	4.205	6.877	8.850
95	1.806	4.227	6.886	8.833



**NOTE:** The pH(S) values of the individual charges of the secondary reference materials are documented in a certificate of an accredited laboratory. This certificate is supplied with the respective buffer materials. Only these pH(S) values shall be used as standard values for the secondary reference buffer materials. Correspondingly, this standard does not include a table with standard pH values for practical use. The table above only provides examples of pH(PS) values for orientation.

### 16.1.5 Hach buffers

Buffer values up to 60 °C as specified by Bergmann & Beving Process AB.

Temp (°C)	pH of buffer solutions		
0	4.00	7.14	10.30
5	4.00	60	10.23
10	4.00	7.04	10.11
15	4.00	7.04	10.11
20	4.00	7.02	10.05
25	4.01	7.00	10.00
30	4.01	6.99	9.96
35	4.02	6.98	9.92
40	4.03	6.98	9.88
45	4.05	6.98	9.85
50	4.06	6.98	9.82
55	4.07	6.98	9.79
60	4.09	6.99	9.76

### 16.1.6 Ciba (94) buffers

Temp (°C)	pH of buffer solutions				
0	2.04	4.00	7.10	10.30	
5	2.09	4.02	7.08	10.21	
10	2.07	4.00	7.05	10.14	
15	2.08	4.00	7.02	10.06	
20	2.09	4.01	6.98	9.99	
25	2.08	4.02	6.98	9.95	
30	2.06	4.00	6.96	9.89	
35	2.06	4.01	6.95	9.85	
40	2.07	4.02	6.94	9.81	
45	2.06	4.03	6.93	9.77	
50	2.06	4.04	6.93	9.73	
55	2.05	4.05	6.91	9.68	
60	2.08	4.10	6.93	9.66	
65	2.07*	4.10*	6.92*	9.61*	
70	2.07	4.11	6.92	9.57	
75	2.04*	4.13*	6.92*	9.54*	
80	2.02	4.15	6.93	9.52	
85	2.03*	4.17*	6.95*	9.47*	
90	2.04	4.20	6.97	9.43	
95	2.05*	4.22*	6.99*	9.38*	

\* Extrapolated

### 16.1.7 Merck Titrisole, Riedel-de-Haën Fixanale

Temp (°C)	pH of buffer solutions				
0	2.01	4.05	7.13	9.24	12.58
5	2.01	4.05	7.07	9.16	12.41
10	2.01	4.02	7.05	9.11	12.26
15	2.00	4.01	7.02	9.05	12.10
20	2.00	4.00	7.00	9.00	12.00
25	2.00	4.01	6.98	8.95	11.88
30	2.00	4.01	6.98	8.91	11.72
35	2.00	4.01	6.96	8.88	11.67
40	2.00	4.01	6.95	8.85	11.54
45	2.00	4.01	6.95	8.82	11.44
50	2.00	4.00	6.95	8.79	11.33
55	2.00	4.00	6.95	8.76	11.19
60	2.00	4.00	6.96	8.73	11.04
65	2.00	4.00	6.96	8.72	10.97
70	2.01	4.00	6.96	8.70	10.90
75	2.01	4.00	6.96	8.68	10.80
80	2.01	4.00	6.97	8.66	10.70
85	2.01	4.00	6.98	8.65	10.59
90	2.01	4.00	7.00	8.64	10.48
95	2.01	4.00	7.02	8.64	10.37

## 16.1.8 WTW buffers

Temp (°C)	pH of buffer solutions			
0	2.03	4.01	7.12	10.65
5	2.02	4.01	7.09	10.52
10	2.01	4.00	7.06	10.39
15	2.00	4.00	7.04	10.26
20	2.00	4.00	7.02	10.13
25	2.00	4.01	7.00	10.00
30	1.99	4.01	6.99	9.87
35	1.99	4.02	6.98	9.74
40	1.98	4.03	6.97	9.61
45	1.98	4.04	6.97	9.48
50	1.98	4.06	6.97	9.35
55	1.98	4.08	6.98	
60	1.98	4.10	6.98	
65	1.99	4.13	6.99	
70		4.16	7.00	
75		4.19	7.02	
80		4.22	7.04	
85		4.26	7.06	
90		4.30	7.09	
95		4.35	7.12	

## 16.1.9 JIS Z 8802 buffers

Temp (°C)	pH of buffer solutions			
0	1.666	4.003	6.984	9.464
5	1.668	3.999	6.951	9.395
10	1.670	3.998	6.923	9.332
15	1.672	3.999	6.900	9.276
20	1.675	4.002	6.881	9.225
25	1.679	4.008	6.865	9.180
30	1.683	4.015	6.853	9.139
35	1.688	4.024	6.844	9.102
38	1.691	4.030	6.840	9.081
40	1.694	4.035	6.838	9.068
45	1.700	4.047	6.834	9.038
50	1.707	4.060	6.833	9.011
55	1.715	4.075	6.834	8.985
60	1.723	4.091	6.836	8.962
70	1.743	4.126	6.845	8.921
80	1.766	4.164	6.859	8.885
90	1.792	4.205	6.877	8.850
95	1.806	4.227	6.886	8.833

## 16.2 Dual membrane pH electrode buffers

### 16.2.1 Mettler-pH/pNa buffers (Na+ 3.9M)

Temp (°C)	pH of buffer solutions			
0	1.98	3.99	7.01	9.51
5	1.98	3.99	7.00	9.43
10	1.99	3.99	7.00	9.36
15	1.99	3.99	6.99	9.30
20	1.99	4.00	7.00	9.25
25	2.00	4.01	7.00	9.21
30	2.00	4.02	7.01	9.18
35	2.01	4.04	7.01	9.15
40	2.01	4.05	7.02	9.12
45	2.02	4.07	7.03	9.11
50	2.02	4.09	7.04	9.10

## Sales and Service:

### Australia

Mettler-Toledo Ltd.  
220 Turner Street  
Port Melbourne  
AUS-3207 Melbourne/VIC  
Phone +61 1300 659 761  
Fax +61 3 9645 3935  
e-mail info.mtaus@mt.com

### Austria

Mettler-Toledo Ges.m.b.H.  
Südrandstraße 17  
A-1230 Wien  
Phone +43 1 604 19 80  
Fax +43 1 604 28 80  
e-mail infoprocess.mtat@mt.com

### Brazil

Mettler-Toledo Ind. e Com. Ltda.  
Alameda Araguaia, 451  
Alphaville  
BR-06455-000 Barueri/SP  
Phone +55 11 4166 7444  
Fax +55 11 4166 7401  
e-mail mettler@mettler.com.br  
service@mettler.com.br

### China

Mettler-Toledo Instruments  
(Shanghai) Co. Ltd.  
589 Gui Ping Road  
Cao He Jing  
CN-200233 Shanghai  
Phone +86 21 64 85 04 35  
Fax +86 21 64 85 33 51  
e-mail mtcs@public.sta.net.cn

### Croatia

Mettler-Toledo d.o.o.  
Mandlova 3  
HR-10000 Zagreb  
Phone +385 1 292 06 33  
Fax +385 1 295 81 40  
e-mail mt.zagreb@mt.com

### Czech Republic

Mettler-Toledo s.r.o.  
Trebohosticka 2283/2  
CZ-100 00 Praha 10  
Phone +420 2 72 123 150  
Fax +420 2 72 123 170  
e-mail sales.mtcz@mt.com

### Denmark

Mettler-Toledo A/S  
Naverland 8  
DK-2600 Glostrup  
Phone +45 43 27 08 00  
Fax +45 43 27 08 28  
e-mail info.mtdk@mt.com

### France

Mettler-Toledo  
Analyse Industrielle S.A.S.  
30, Boulevard de Douaumont  
F-75017 Paris  
Phone +33 1 47 37 06 00  
Fax +33 1 47 37 46 26  
e-mail mtpro-f@mt.com

### Germany

Mettler-Toledo GmbH  
Prozeßanalytik  
Ockerweg 3  
D-35396 Gießen  
Phone +49 641 507 333  
Fax +49 641 507 397  
e-mail prozess@mt.com

### Great Britain

Mettler-Toledo LTD  
64 Boston Road, Beaumont Leys  
GB-Leicester LE4 1AW  
Phone +44 116 235 7070  
Fax +44 116 236 5500  
e-mail enquire.mtuk@mt.com

### Hungary

Mettler-Toledo Kereskedelmi KFT  
Teve u. 41  
HU-1139 Budapest  
Phone +36 1 288 40 40  
Fax +36 1 288 40 50  
e-mail mthu@axelero.hu

### India

Mettler-Toledo India Private Limited  
Amar Hill, Saki Vihar Road  
Powai  
IN-400 072 Mumbai  
Phone +91 22 2857 0808  
Fax +91 22 2857 5071  
e-mail sales.mfin@mt.com

### Italy

Mettler-Toledo S.p.A.  
Via Vialba 42  
I-20026 Novate Milanese  
Phone +39 02 333 321  
Fax +39 02 356 2973  
e-mail customercare.italia@mt.com

### Japan

Mettler-Toledo K.K.  
Process Division  
6F Ikenohata Nisshoku Bldg.  
2-9-7, Ikenohata  
Taito-ku  
JP-110-0008 Tokyo  
Phone +81 3 5815 5606  
Fax +81 3 5815 5626  
e-mail helpdesk.ing.jp@mt.com

### Malaysia

Mettler-Toledo (M) Sdn Bhd  
Bangunan Electroscon Holding, U 1-01  
Lot 8 Jalan Astaka U8/84  
Seksyen U8, Bukit Jelutong  
MY-40150 Shah Alam Selangor  
Phone +60 3 78 44 58 88  
Fax +60 3 78 45 87 73  
e-mail  
MT-MY.CustomerSupport@mt.com

### Mexico

Mettler-Toledo S.A. de C.V.  
Ejercito Nacional #340  
Col. Chapultepec Morales  
Del. Miguel Hidalgo  
MX-11570 México D.F.  
Phone +52 55 1946 0900  
e-mail ventas.lab@mt.com

### Poland

Mettler-Toledo (Poland) Sp.z.o.o.  
ul. Poleczki 21  
PL-02-822 Warszawa  
Phone +48 22 545 06 80  
Fax +48 22 545 06 88  
e-mail polska@mt.com

### Russia

Mettler-Toledo Vostok ZAO  
Sretenskij Bulvar 6/1  
Office 6  
RU-101000 Moscow  
Phone +7 495 621 56 66  
Fax +7 495 621 63 53  
e-mail inforus@mt.com

### Singapore

Mettler-Toledo (S) Pte. Ltd.  
Block 28  
Ayer Rajah Crescent #05-01  
SG-139959 Singapore  
Phone +65 6890 00 11  
Fax +65 6890 00 12  
+65 6890 00 13  
e-mail precision@mt.com

### Slovakia

Mettler-Toledo s.r.o.  
Hattalova 12/A  
SK-831 03 Bratislava  
Phone +421 2 4444 12 20-2  
Fax +421 2 4444 12 23  
e-mail predaj@mt.com

### Slovenia

Mettler-Toledo d.o.o.  
Pot heroja Trtnika 26  
SI-1261 Ljubljana-Dobrunje  
Phone +386 1 530 80 50  
Fax +386 1 562 17 89  
e-mail keith.racman@mt.com

### South Korea

Mettler-Toledo (Korea) Ltd.  
Yeil Building 1 & 2 F  
124-5, YangJe-Dong  
SeCho-Ku  
KR-137-130 Seoul  
Phone +82 2 3498 3500  
Fax +82 2 3498 3555  
e-mail Sales\_MTKR@mt.com

### Spain

Mettler-Toledo S.A.E.  
C/Miguel Hernández, 69-71  
ES-08908 L'Hospitalet de Llobregat  
(Barcelona)  
Phone +34 902 32 00 23  
Fax +34 902 32 00 24  
e-mail mtemkt@mt.com

### Sweden

Mettler-Toledo AB  
Virkesvägen 10  
Box 92161  
SE-12008 Stockholm  
Phone +46 8 702 50 00  
Fax +46 8 642 45 62  
e-mail sales.mts@mt.com

### Switzerland

Mettler-Toledo (Schweiz) GmbH  
Im Langacher  
Postfach  
CH-8606 Greifensee  
Phone +41 44 944 45 45  
Fax +41 44 944 45 10  
e-mail salesola.ch@mt.com

### Thailand

Mettler-Toledo (Thailand) Ltd.  
272 Soi Soonvijai 4  
Rama 9 Rd., Bangkapi  
Huay Kwang  
TH-10320 Bangkok  
Phone +66 2 723 03 00  
Fax +66 2 719 64 79  
e-mail  
MT-TH.CustomerSupport@mt.com

### USA/Canada

Mettler-Toledo Ingold, Inc.  
36 Middlesex Turnpike  
Bedford, MA 01730, USA  
Phone +1 781 301 8800  
Freephone +1 800 352 8763  
Fax +1 781 271 0681  
e-mail mtpro@mt.com



Management System  
certified according to  
ISO 9001 / ISO 14001



Subject to technical changes.  
© Mettler-Toledo AG, Process Analytics  
06/12 Printed in Switzerland. 52 121 825

Mettler-Toledo AG, Process Analytics  
Im Hackacker 15, CH-8902 Urdorf, Switzerland  
Tel. +41 44 729 62 11, Fax +41 44 729 66 36