NewClassic Balances ML Models



Direct Adjustment Access

Execute predefined adjusting (calibration) procedure Factory Setting: Internal adjustment (MLxxxE models: external adjustment)

External Adjustment Procedure



Operation Keys

Tare

METTLER TOLEDO									
	lo		BCU.			_			
	+ F	On	→0/T←	ĉ	(<u>)</u> ≈ ●	ل •			
12	3	4	5	6	7	89			

Key F	Key Functions						
No.	Key	Press briefly (less than 1.5 s)	Press and hold (longer than 1.5 s)				
1		 To change display resolution (1/10d function) while application is running 	no function				
2	Δ̈́Δ	 To navigate back (scroll up) within menu topics or menu selections Decrease (numerical) parameters within menu and in applications 	 To select the weighing application Decrease (numerical) parameters quickly within menu or in applications 				
3	∔ ⊧	 To navigate forward (scroll down) within menu topics or menu selections Increase (numerical) parameters within menu and in applications 	 To select assigned application and entering the parameter settings of application Default application assignment: Piece counting Increase (numerical) parameters quickly within menu or in applications 				
4	ON/OFF	Switch on	Switch off				
5	→0/T←	• Zero/Tare	Switch off				
6	С	 Cancel and to leave menu without saving (one step back in the menu) 	no function				
7	<u>دی</u>	 With entries: scroll down To navigate through menu topics or menu selections To toggle between unit 1, recall value (if selected), unit 2 (if different from unit 1) and the application unit (if any) 	 Execute predefined adjusting (calibration) procedure 				
8	Ĵ Ĺ	To enter or leave menu selection (from / to menu topic) To enter application parameter and switch to next parameter To store parameter	Enter and leave menu (Parameter settings)				
9	.⊟.	 Printout display value Printout menu settings Transmit data 	no function				

Status Icons Application Icons ⋳**⋴**⊓ ४४०००००००००००००००००००००००००० ₽∿ MNet 🖸 FACT .W.W.W GNCtis%bahtin msgPCStbidæt V Ŵ ·M ÷ Unit Field Application Icons Status Icons Indicates stored values (memory) 0 Menu locked Μ ö Menu settings activated To mark Net weight values Net ₫∆ Application "Weighing" Adjustments (calibration) ٢ started Application "Piece counting" FACT FACT activated Application "Percent weighing" % Diagnostic mode and Applica-tion "Routine test" activated Y ⊳⊲ Application "Check weighing" Service reminder, Service mode activated ᢣ <u>......</u> Application "Statistics" Application "Formulation / Net-Total" Acoustic feedback for pressed (((•))) 四 keys activated Weighing range 1 (Dual Range models only) Σ Application "Totaling" W1 Application "Dynamic weighing" Weighing range 2 (Dual Range models only) $\overline{\mathbf{W}}$ **W2** Charge of battery; full, 2/3, 1/3, discharged (Battery operated models only) Application "Multiplication factor" F×T Application "Division factor" F÷∎ Weight Value Field Ð Application "Density" Indicates negative values Unit Field Marking of weighing units: mg milligram g gram kg kilogram ct carat ib pound oz ounce ozt troy ounce GN grain dwt pennyweight mom momme msg mesghal tih Hong Kong teels tis Singapore teels tis Singapore teels tif Taiwan teels tola tola Ο Indicates unstable values GNctls%bohtlh msgPCStbidat kgmgm mg g kg ct Ib oz ozt GN dwt mon msg tih tis tit tola baht * Indicates calculated values Indicates uncertified digits (approved models only) Weighing-in aid Marking of nominal or target weight ⊳ Marking of tolerance limit T+ tola baht ► Marking of tolerance limit T-

Display Panel



Menu Map







Legend 144V

- Factory setting Only those weight units allowed by the appropriate national legislation are selectable. 1
- Not available with MLxxxE models θ
- Not available with approved models Not available with approved models Ō
- ۲

(except OIML accuracy class I approved models) Settings are automatically set for the 2nd display

Software version V2.20

Mettler-Toledo AG

www.mt.com/newclassic

For more information



Laboratory Weighin CH-8606 Greifensee, Switzerland