


Pharmacy Balances

PHS Models



METTLER TOLEDO

Table of Contents

1	Introduction	7
	1.1 Conventions and Symbols Used in These Operating Instructions	7
2	Safety Precautions	8
3	Overview	9
	3.1 S Platform	9
	3.2 L Platform	10
	3.3 Operation Keys	11
	3.4 Display Panel	12
4	Setting up the Balance	14
	4.1 Unpacking and Delivery Inspection	14
	4.2 Installing the Components	15
	4.3 Selecting the Location and Leveling the Balance	16
	4.3.1 Selecting the Location	16
	4.3.2 Leveling the Balance	17
	4.4 Power Supply	18
	4.5 Transporting the Balance	18
	4.6 Weighing Below the Balance	18
	4.7 General Requirements	19
	4.7.1 Switching on the Balance	19
	4.7.2 Adjusting the Balance	19
	4.8 Adjustment (Calibration)	20
	4.8.1 Fully Automatic Adjustment FACT	20
	4.8.2 Manual Adjustment with Internal Weight	20
	4.8.3 Manual Adjustment with External Weight	20
	4.8.4 Customer Fine Adjustment	21
5	Weighing Made Simple	23
	5.1 Switching the Balance On and Off	23
	5.2 Performing a Simple Weighing	24
	5.3 Zero Setting / Taring	24
	5.4 METTLER TOLEDO DeltaRange Balances	25
	5.5 Switching Weight Units	25
	5.6 Recall / Recall Weight Value	25
	5.7 Weighing with the Weighing-in Aid	25
	5.8 Print / Transmit Data	26
6	The Menu	27
	6.1 What is in the Menu ?	27
	6.2 Menu Operation	28
	6.3 Description of Menu Topics	29
	6.3.1 Main Menu	29
	6.3.2 Basic Menu	30
	6.3.3 Advanced Menu	31
	6.3.4 Interface Menu	35

7	Application "Formulation" (Net Total Formulation)	43
8	Application "Piece Counting"	46
9	Application "Check Weighing"	49
10	Application "Statistics"	52
11	Application "Routine Test"	54
12	Application "Diagnostics"	57
	12.1 Repeatability Test	57
	12.2 Display Test	58
	12.3 Key Test	59
	12.4 Motor Test	60
	12.5 Balance History	60
	12.6 Calibration History	61
	12.7 Balance Information	62
	12.8 Service Provider Information	63
13	Communication with Peripheral Devices	64
	13.1 Function PC-Direct	64
	13.2 USB Device Interface	65
14	Firmware (Software) Updates	67
	14.1 Operating Principle	67
	14.2 Update Procedure	67
15	Error and Status Messages	69
	15.1 Error Messages	69
	15.2 Status Messages	70
16	Cleaning and Service	71
	16.1 Cleaning the Glass Draft Shield (0.1 mg and 1 mg Models)	71
17	Interface Specification	74
	17.1 RS232C Interface	74
	17.2 USB Device Interface	74
	17.3 MT-SICS Interface Commands and Functions	75
18	Technical Data	76
	18.1 General Data	76
	18.2 Model-Specific Data	77
	18.2.1 Balances with Readability of 0.1 mg, S Platform with Draft Shield	77
	18.2.2 Balances with Readability of 1 mg, S Platform with Draft Shield	78
	18.2.3 Balances with Readability of 0.01 g, S Platform	80
	18.2.4 Balances with Readability of 0.1 g, L Platform	81
	18.3 Dimensions	82
	18.3.1 Balances with Readability of 0.1 mg, S Platform with Draft Shield	82
	18.3.2 Balances with Readability of 1 mg, S Platform with Draft Shield	83
	18.3.3 Balances with Readability of 0.01 g, S Platform	84

	18.3.4	Balances with Readability of 0.1 g, L Platform	85
19		Accessories and Spare Parts	86
	19.1	Accessories	86
	19.2	Spare Parts	90
		Index	91


1 Introduction

Thank you for choosing a METTLER TOLEDO balance.

The precision balances of the Pharmacy line combine a large number of weighing possibilities with easy operation. These balances are designed specially for use in pharmacies.

These operating instructions apply to all balance models PHS in the Pharmacy line. However, the different models have different characteristics regarding equipment and performance. Special notes in the text indicate where this makes a difference to operation.

1.1 Conventions and Symbols Used in These Operating Instructions

Key designations are indicated by double angular brackets (e.g. «»).



This symbol indicates press key briefly (less than 1.5 s).



This symbol indicates press and hold key down (longer than 1.5 s).




This symbol indicates a flashing display.



This symbol indicates an automatic sequence.



These symbols indicate safety notes and hazard warnings which, if ignored, can cause personal danger to the user, damage to the balance or other equipment, or malfunctioning of the balance.



This symbol indicates additional information and notes. These make working with your balance easier, as well as ensuring that you use it correctly and economically.

2 Safety Precautions

Always operate and use your balance only in accordance with the instructions contained in this manual. The instructions for setting up your new balance must be strictly observed.

If the balance is not used according to these Operating Instructions, protection of the balance may be impaired and METTLER TOLEDO assumes no liability.



It is not permitted to use the balance in explosive atmosphere of gases, steam, fog, dust and flammable dust (hazardous environments).



Use only the original Universal AC adapter delivered with your balance.

The L platform has a built-in power supply unit. Hazard of electric shock if the power cable is damaged! Check the power cable for damage regularly. Unplug the power cord immediately if the power cable is damaged.



Do not use sharply pointed objects to operate the keyboard of your balance! Although your balance is very ruggedly constructed, it is nevertheless a precision instrument. Treat it with corresponding care.

Do not open the balance: It does not contain any parts which can be maintained, repaired, or replaced by the user. If you ever have problems with your balance, contact your METTLER TOLEDO dealer.

Use only balance accessories and peripheral devices from METTLER TOLEDO; they are optimally adapted to your balance.



Disposal






In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.


3 Overview

3.1 S Platform




Name and Function of Components	
1	Display
2	Operation keys
3	Glass draft shield
4	Handle for operation of the draft shield door
5	Weighing pan
6	Draft shield element
7	Level indicator
8	Kensington slot for anti-theft purposes
9	Safety feet (with 10 mg, 0.1 g S series models)
10	Leveling foot
11	USB device interface
12	RS232C serial interface
13	Socket for AC Adapter
14	Model sticker (with approved models only)
15	Product label

3.2 L Platform











Name and Function of Components	
1	Display
2	Operation keys
3	Weighing pan
4	Level indicator
5	Leveling foot
6	USB device interface
7	RS232C serial interface
8	Power cord with country-specific plug
9	Security slot for anti-theft purposes
10	Model sticker (with approved models only)
11	Product label

3.3 Operation Keys




Key Functions

No.	Key	Press briefly (less than 1.5 s) 	Press and hold (longer than 1.5 s) 
1		<ul style="list-style-type: none"> To change display resolution (1/10d function) while application is running Note: not available with approved models in selected countries. 	no function
2		<ul style="list-style-type: none"> Enter or leave menu (Parameter settings) Save parameters 	no function
3		<ul style="list-style-type: none"> Execute predefined adjusting (calibration) procedure 	no function
4		<ul style="list-style-type: none"> Printout display value Printout active user menu settings Transfer data 	no function
5		<ul style="list-style-type: none"> To navigate back (scroll up) within menu topics or menu selections Decrease (numerical) parameters within menu and in applications 	<ul style="list-style-type: none"> To select the weighing application Decrease (numerical) parameters quickly within menu and in applications
6		<ul style="list-style-type: none"> To navigate forward (scroll down) within menu topics or menu selections Increase (numerical) parameters within menu and in applications 	<ul style="list-style-type: none"> To select assigned F1 application and entering the parameter settings of application. Default F1 application assignment: Formulation Increase (numerical) parameters quickly within menu and in applications

No.	Key	Press briefly (less than 1.5 s)	Press and hold (longer than 1.5 s)
7	F2	<ul style="list-style-type: none"> With entries: scroll down To navigate through menu topics or menu selections To toggle between unit 1, recall value (if selected), unit 2 (if different from unit 1) and the application unit (if any) 	<ul style="list-style-type: none"> To select assigned F2 application and entering the parameter settings of application. Default F2 application assignment: Percent weighing
8	F3	<ul style="list-style-type: none"> To enter or leave menu selection (from / to menu topic) To enter application parameter or switch to next parameter To confirm parameter 	<ul style="list-style-type: none"> To select assigned F3 application and entering the parameter settings of application. Default F3 application assignment: Check weighing
9	ON/OFF →0/T←	<ul style="list-style-type: none"> Switch on Zero/Tare 	<ul style="list-style-type: none"> Switch off
10	C	<ul style="list-style-type: none"> Cancel and to leave menu without saving (one step back in the menu). 	no function

3.4 Display Panel











Application Icons			
	Menu locked		Application "Check weighing"
	Menu setting activated		Application "Statistics"
	Application "Weighing"		Application "Formulation / Net-Total"
	Application "Piece counting"		

Note

While an application is running, the corresponding application icon appears at the top of the display.

Status Icons			
M	Indicates stored value (Memory)		Service reminder
Net	Indicates Net weight values		Acoustic feedback for pressed keys activated
	Adjustments (calibration) started	W1	Weighing range 1 (Dual Range models only)
FACT	FACT activated	W2	Weighing range 2 (Dual Range models only)
	Applications "Diagnostics" and "Routine Test"		Charge of battery: full, 2/3, 1/3, discharged (Battery operated models only)

Weight Value Field and Weighing-in aid			
	Indicates negative values		Brackets to indicate uncertified digits (approved models only)
	Indicates unstable values		Marking of nominal or target weight
	Indicates calculated values		Marking of tolerance limit T+
			Marking of tolerance limit T-

Unit Field		
	g	gram
	kg	kilogram
	mg	milligram

4 Setting up the Balance



The balance must be disconnected from the power supply when carrying out all setup and mounting work.


4.1 Unpacking and Delivery Inspection

- 1 Open the packaging and carefully remove all components.
- 2 Check the delivered items.

The standard scope of delivery contains the following items:

Components		S platform				L platform
		0.1 mg	1 mg	0.01 g	0.1 g	0.1g
Draft shield	236 mm	✓	–	–	–	–
	168 mm	–	✓	–	–	–
Weighing pan	Ø 90 mm	✓	–	–	–	–
	127 x 127 mm	–	✓	–	–	–
	170 x 200 mm	–	–	✓	–	–
	190 x 226 mm	–	–	–	✓	–
	245 x 351 mm	–	–	–	–	✓
Draft shield element		✓	–	✓	–	–
Pan support		✓	✓	✓	✓	–
Bottom plate		✓	✓	–	–	–
Protective cover		✓	✓	✓	✓	✓
Universal AC adapter (country specific)		✓	✓	✓	✓	–
Mounted country specific power cable		–	–	–	–	✓
Operating instructions printed or on CD-ROM depending on the country		✓	✓	✓	✓	✓
Quick Guide		✓	✓	✓	✓	✓
EC declaration of conformity		✓	✓	✓	✓	✓

4.2 Installing the Components




Balances with readability of 0.1 mg, S platform with draft shield (236 mm)

Place the following components on the balance in the specified order:

Note: Push the side glass back as far as will go and grasp the draft shield (1) with both hands on the bars at the top.

- 1 Turn draft shield lock (2) to position "U" (unlock).
- 2 Place draft shield on the balance.
- 3 Turn draft shield lock to "L" (lock) and place bottom plate (3).
- 4 Place draft shield element (4) and weighing pan (6) with pan support (5).

Note: Cleaning the draft shield see section "Maintenance and cleaning".




Balances with readability of 1 mg, S platform with draft shield (168 mm)

Place the following components on the balance in the specified order:

Note: Push the side glass back as far as will go and grasp the draft shield (1) with both hands on the bars at the top.


- 1 Turn draft shield lock (2) to position "U" (unlock).
- 2 Place draft shield on the balance.
- 3 Turn draft shield lock to "L" (lock) and place bottom plate (3).
- 4 Place weighing pan (5) with pan support (4).

Note: Cleaning the draft shield see section "Maintenance and cleaning".




Balances with readability of 10 mg, S platform

- Place the following components on the balance in the specified order:
 - Pan support (1)
 - Weighing pan (2)
 - Draft shield element (3)



Balances with readability of 0.1 g, S platform

- Place the following components on the balance in the specified order:
 - Pan support (1)
 - Weighing pan (2)



Balances with readability to 1 g, L platform


- Place the weighing pan (1) on the balance.

4.3 Selecting the Location and Leveling the Balance

Your balance is a precision instrument and will thank you for an optimum location with high accuracy and dependability.

4.3.1 Selecting the Location

Select a stable, vibration-free position that is as horizontal as possible. The surface must be able to safely carry the weight of a fully loaded balance.




Observe ambient conditions (see Technical Data).

Avoid the following:


- Direct sunlight
- Powerful drafts (e.g. from fans or air conditioners)
- Excessive temperature fluctuations

4.3.2 Leveling the Balance



The balances have a level indicator and two (S Platform) or four (L Platform) adjustable leveling feet to compensate for slight irregularities in the surface of the weighing bench. The balance is exactly horizontal when the air bubble is in the middle of the level glass.


Note: The balance should be leveled and adjusted each time it is moved to a new location.



Balances with S platform and readability of 0.1 mg and 1 mg

- Adjust the two leveling feet appropriately until the air bubble comes to rest exactly in the middle of the glass:

Air bubble at	"12 o'clock"	turn both feet clockwise
Air bubble at	"3 o'clock"	turn left foot clockwise, right foot counterclockwise
Air bubble at	"6 o'clock"	turn both feet counterclockwise
Air bubble at	"9 o'clock"	turn left foot counterclockwise, right foot clockwise



Balances with S platform and readability of 10 mg and 0.1 g

- 1 Remove the clamps (A) for the safety feet by turning them outwards.
Note: Turn the clamps (A) outwards as far as they will go (~ 90°), so that the safety feet can move freely.
- 2 Now level the balance by turning both leveling screws (B) until the air bubble is in the inner circle of the level indicator (see procedure above).
- 3 Secure the safety feet by turning the clamps (A) inwards as far as they will go.

Balances with L platform

- Align the balance horizontally by turning the leveling screws of the balance housing until the air bubble is in the inner circle of the level indicator.

4.4 Power Supply

Your balance is supplied with an country-specific AC adapter or with a country-specific power cable. The power supply is suitable for all line voltages in the range: 100 - 240 VAC, 50/60 Hz (for exact specifications, see section "technical data").




First, check the local line voltage is in the range 100 - 240 VAC, 50/60 Hz and whether the power plug fits your local power supply connection. **If this is not the case, on no account connect the balance or the AC adapter to the power supply**, but contact the responsible METTLER TOLEDO dealer.




Important:

- Before operating, check all cables for damage.
- Guide the cables so that they cannot become damaged or interfere with the weighing process!
- Take care that the AC adapter cannot come into contact with liquids!
- The power plug must be always accessible.



Allow your balance to warm up for 30 minutes (0.1 mg models 60 minutes) to enable it to adapt itself to the ambient conditions.



Connect the AC adapter to the connection socket on the back of your balance (see figure) and to the power line.

4.5 Transporting the Balance

Switch off the balance and remove the power cable and any interface cable from the balance. Refer to the notes in Section "Selecting the location" regarding the choice of an optimal location.

Transporting Over Short Distances



For balances with a draft shield: Observe the following instructions to transport your balance over a short distance to a new location: **Never lift the balance by the glass draft shield. The draft shield is not sufficiently fastened to the balance.**

Transporting Over Long Distances


If you would like to transport or send your balance over long distances, **use the complete original packaging.**

4.6 Weighing Below the Balance


The balances are equipped with a hanger for carrying out weighings below the work surface (weighing below the balance).

**Attention:**

- **Do not place the balance on the pan support location bolt** (0.1 mg and 1 mg models).
- Models with a glass draft shield: Carefully lift the draft shield from the weighing platform and put it aside.

**Note:**

- For below-the-balance weighing with the L Platform models, you will need hook 11132565 from the accessories range.
- Weighing below the balance is not possible with "MS-KL" models.



- 1 Switch off the balance and remove the power cable and any interface cable from the balance.
- 2 Remove the draft shield element (10 mg models).
- 3 Remove the weighing pan and pan support.
- 4 Remove the bottom plate and unlock the draft shield (models with draft shield). Carefully lift the draft shield from the weighing platform and put it aside.
- 5 Turn the balance carefully on its side.
- 6 Remove and retain the 2 screws (A) and the cover plate (B). The hanger is now accessible.
- 7 Then turn the balance to its normal position and simply reinstall all components in the reverse order.

4.7 General Requirements

4.7.1 Switching on the Balance

Before working with the balance, it must be warmed up in order to obtain accurate weighing results. To reach operating temperature, the balance must be connected to the power supply for at least

- 30 minutes on balances with a readability of 0.001 g (0.01 ct) to 0.1 g.
- 60 minutes on balances with a readability of 0.1 mg (0.001 ct) and better.

See also

- Adjusting the Balance (page 19)

4.7.2 Adjusting the Balance

To obtain accurate weighing results, the balance must be adjusted to match the gravitational acceleration at its location and depending on the ambient conditions. After reaching the operation temperature, adjusting is necessary

- before the balance is used for the first time.
- when the balance (readability of 0.1 mg and better) was disconnected from the power or in case of power failure.

- after a change of the location.
- at regular intervals during weighing service.

See also

- Switching on the Balance (page 19)

4.8 Adjustment (Calibration)

Attention

Before adjusting the balance, it must be warmed up.


4.8.1 Fully Automatic Adjustment FACT

The **factory setting** is fully automatic adjustment **FACT** (Fully Automatic Calibration Technology) with the internal weight (see also section "The Menu").


The balance adjusts itself automatically:

- after the warm-up phase on connection to the power supply.
- when a change in the ambient conditions, e.g. the temperature, could lead to a noticeable deviation in the measurement.
- on a predefined time. (see menu topic "FACT")
- time interval. (with OIML accuracy class II approved models)

4.8.2 Manual Adjustment with Internal Weight




Requirement: To carry out this operation, in the menu topic "**CAL**" (Adjustment) of advanced menu "**ADJ.INT**" must be selected.

- 1 Unload weighing pan
- 2 Press «» to execute "Internal Adjustment".


The balance adjusts itself automatically. The adjusting is finished when the message "**ADJ DONE**" appears briefly on the display. The balance returns to the last active application and is ready for operation.

Sample adjustment printout using internal weight:



4.8.3 Manual Adjustment with External Weight

Note: Because of certification legislation, the approved models cannot be adjusted with an external weight * .
 * except OIML accuracy class I approved models.



Requirement: To carry out this operation, in the menu topic "**CAL**" (Adjustment) of advanced menu "**ADJ.EXT**" must be selected.

Note

We recommend to disable FACT.

- 1 Have required adjustment weight ready.
- 2 Unload weighing pan.
- 3 Press «» briefly to execute "External Adjustment". The required (predefined) adjustment weight value flashes on the display.
- 4 Place adjustment weight in center of pan. The balance adjusts itself automatically.
- 5 When "0.00 g" flashes, remove adjustment weight.

The adjusting is finished when the message "**ADJ DONE**" appears briefly on the display. The balance returns to the last active application and is ready for operation.

Sample adjustment printout using external weight:

```

- External Adjustment --
21.Jan 2010      12:56

METTLER TOLEDO

Balance Type     PH4002S
SNR              1234567890

Temperature      22.5 °C
Nominal          2000.00 g
Actual          1999.99 g
Diff             5 ppm

Adjustment done

Signature

.....
-----

```

4.8.4 Customer Fine Adjustment

Attention

This function should be executed only by trained personnel.


The function customer fine adjustment "**ADJ.CUST.F**" allows you to adjust the value of the internal adjustment weight with your own adjustment weight. The adjustable range of the adjustment weight is possible only in a very small range. Customer fine adjustment impacts the function of internal adjustment. The customer fine adjustment can be deactivated at any time.

Note

- This feature is available on models with internal weight only.
- Because of certification legislation, approved models cannot be adjusted with customer fine adjustment (depending on selected countries' certification legislation).
- Use certificated weights.

- Balance and test weight have to be on operating temperature.
- Observe the correct environmental conditions.

Execute customer fine adjustment



- The balance is under measuring condition.
- 1 Have required adjustment weight ready.
 - 2 Unload weighing pan
 - 3 Select in the menu "**ADVANCED**": **ADJ.CUST.F**
 - 4 Confirm "**ADJ.CUST.F**" with «←».
 - 5 To carry out this operation select "**EXECUTE**"
 - 6 Start Adjustment with «←»
 - ⇒ "**SET REF.**" appears briefly.
 - ⇒ The last saved value flashes on the display.
 - 7 Select the target adjustment weight.
 - For coarse setting, press «↶» to change the value.
 - For fine setting, press «+» to increase the value or press "-" to decrease the value.
 - 8 Press and hold «←» to confirm and execute "**ADJ.CUST.F**".
 - ⇒ The required adjustment weight value flashes in the display. This could take some time.
 - 9 Place required adjustment weight in center of pan.
 - 10 Remove adjustment weight when zero is flashing.
 - 11 Wait until "**ADJ DONE**" briefly appears.
 - ⇒ The adjusting is finished when the message "**ADJ DONE**" appears briefly on the display. The balance returns to the last active application and is ready for operation
 - ⇒ If the error message "**WRONG ADJUSTMENT WEIGHT**" appears, the weight is not within the allowed value range and could not be accepted. "**ADJ.CUST.F**" could not be executed.

Note

Storing the adjustment is not required.

Deactivate customer fine adjustment

- 1 Select in the menu "**ADVANCE.**": "**ADJ.CUST.F**".
 - 2 Confirm "**ADJ.CUST.F**" with «←».
 - 3 To carry out this operation select "**RESET**"
 - 4 Start **RESET** by pressing «←»
 - ⇒ "**NO?**" appears.
 - 5 Select "**YES?**" and confirm with «←».
- ⇒ The adjusting is finished when the message "**ADJ DONE**" appears briefly on the display. The balance returns to the last active application and is ready for operation with initial adjustment.


5 Weighing Made Simple



This section shows you how to perform simple weighings and how you can accelerate the weighing process.

5.1 Switching the Balance On and Off

Switching on




Connecting to the mains

- 1 Remove any load from weighing pan.
- 2 Connect balance via AC adapter to the mains.

The balance performs a display test (all segments in the display light up briefly), "WELCOME", Software version, **Maximum load** and **Readability** appears briefly. (Startup "FULL" mode only)


After the warm-up time, the balance is ready for weighing or for operation with the last active application, **see** General Requirements.



Mains operated (standby mode)

- Press «On».

The balance is ready for weighing or for operation with the last active application. Approved balances will execute an initial zero.




Battery operated

- 1 Remove any load from weighing pan.
- 2 Press and hold «On»

The balance performs a display test (all segments in the display light up briefly), "WELCOME", Software version, **Maximum load** and **Readability** appears briefly. (Startup "FULL" mode only)

After the warm-up time, the balance is ready for weighing or for operation with the last active application, **see** General Requirements.

Switching off




- Press and hold the «Off» key until "SHUTOFF" appears on the display. Release the key.
- ⇒ Mains operated balances switch into standby mode.
- ⇒ Battery operated balances switch off completely.

Note:

- After switching on from standby mode, your balance needs no warm-up time and is immediately ready for weighing.
- Standby mode is not possible with approved balances (only available in selected countries).
- If your balance has been switched off after a preselected time, the display is dimly lit and shows date, time, maximum load and readability.


- If your balance has been switched off manually, the display is off.
- To completely switch off mains operated balances, they must be disconnected from the power supply.

5.2 Performing a Simple Weighing



- 1 Press «→0/T←» to zero the balance.
Note: If your balance is not in the weighing mode, press and hold the « $\Delta\Delta$ » key down until "WEIGHING" appears in the display. Release the key. Your balance is in the weighing mode and set to zero.
- 2 Place weighing sample on the weighing pan.
- 3 Wait until the instability detector "O" disappears and the stability beep sounds.
- 4 Read the result.


5.3 Zero Setting / Taring



Zero setting

- 1 Unload the balance.
- 2 Press «→0/T←» to set the balance to zero. All weight values are measured in relation to this zero point (see menu topic "ZERO RNG").

Note: Use the «→0/T←» zeroing key before you start with a weighing.



Taring

If you are working with a weighing container, first set the balance to zero.


- 1 Place empty container on the balance. The weight is displayed.
- 2 Press «→0/T←» to tare the balance.

"0.00 g" and "Net" appears in the display. "Net" indicates that all weight values displayed are net values.

Note:

- If the container is removed from the balance, the tare weight will be shown as a negative value.
- The tare weight remains stored until the «→0/T←» key is pressed again or the balance is switched off.
- With METTLER TOLEDO DeltaRange balances, the fine range with its 10 times smaller display increments (depending on the model) is available again after every taring operation.

5.4 METTLER TOLEDO DeltaRange Balances




METTLER TOLEDO DeltaRange balances have a movable fine range with 10 times smaller display increments over their entire weighing range. In this fine range an additional decimal place always appears in the display.

The balance operates in the fine range

- after switching on.
- after every taring operation.

If the fine range is exceeded, the balance display automatically switches to coarser display increments.

5.5 Switching Weight Units




The $\leftarrow \rightleftarrows \rightarrow$ key can be used at any time to toggle between weight unit "UNIT 1", "RECALL" value (if selected) and weight unit "UNIT 2" (if different from weight unit 1) and the application unit (if any).

5.6 Recall / Recall Weight Value

Recall stores stable weights with an absolute display value bigger than 10d.

Requirement: The function "RECALL" must be activated in the menu.




- 1 Load weighing sample. The display shows weight value and stores stable value.
- 2 Remove weighing sample. When the weight is removed the Display shows zero.
- 3 Press $\leftarrow \rightleftarrows \rightarrow$. The display shows last stored stable weight value for 5 seconds together with asterisk (*) and Memory (M) symbols. After 5 seconds the display goes back to zero. This can be repeated unlimited times.

Delete last weight value

As soon a new stable weight value is displayed, the old recall value becomes replaced by the new weight value. When pressing $\leftarrow \rightarrow 0/T \leftarrow$, the recall value is set to 0.


Note: If the power is switched off, the recall value is lost. The recall value can not be printed.


5.7 Weighing with the Weighing-in Aid



The weighing-in aid is a dynamic graphic indicator which shows the used amount of the total weighing range. You can thus recognize at a glance when the load on the balance approaches the maximum load.

5.8 Print / Transmit Data



Pressing the «» key transmits the weighing results over the interface e.g. to a printer or a PC.

6 The Menu

6.1 What is in the Menu ?



The Menu allows you to match your balance to your specific weighing needs. In the menu you can change the settings of your balance and activate functions. The main menu has 4 different menus and these contains 45 different **topics**, each of which allows you various **selection** possibilities. For Menu "**PROTECT**" see chapter "Description of menu topics" section "Main menu".

Note: See Quick Guide for the graphical overview of the menu (Menu Map) with all setting possibilities.

Menu "BASIC"

Topic	Description
DATE	Setting the current date.
TIME	Setting the current time.
UNIT 1	Specification of the 1 st weight unit in which the balance should show the result.
UNIT 2	Specification of the 2 nd weight unit in which the balance should show the result.
KEY BEEP	Setting the key beep level.
STAB.BEEP	Setting the stability beep level.
RESET	Call up of the factory settings.

Menu "ADVANCED"

Topic	Description
WEIG.MODE	Adapting the balance to the weighing mode.
ENVIRON.	Matching the balance to the ambient conditions.
CAL	Settings for the type of adjustment (calibration).
ADJ.CUST.F	Executing customer fine adjustment.
FACT	Settings for fully automatic balance adjustment based on a selected time.
FACT PRT.	Switching the automatic FACT printout on or off.
DATE.FORM	Setting the date format.
TIME.FORM	Preselection of the time format.
RECALL	Switching the application "Recall" for storing stable weights on or off.
SHUTOFF	Setting the time after which the balance should be switched off automatically.
BCKLIGHT	Setting the time after which the display backlight should be switched off automatically.
DISPLAY	Adjusting the brightness and contrast of the display.
AUTOZERO	Switching the automatic zero correction (Autozero) on or off.
ZERO RNG	Setting the zero limit of the zero/tare key.
LANGUAGE	Setting the preferred language.
ASSIGN:F1	Selection of assigned F1 key application and entering their parameter settings.
ASSIGN:F2	Selection of assigned F2 key application and entering their parameter settings.
ASSIGN:F3	Selection of assigned F3 key application and entering their parameter settings.
DIAGNOSE	Starting a diagnostic application.
SERV.ICON	Switching the service icon (service reminder) on or off.
SRV.D.RST	Reset service date and hours (service reminder).


Menu "INT.FACE"

Topic	Description
RS232	Matching the serial interface RS232C to a peripheral unit.
HEADER	Setting the header for printout of individual values.
SINGLE	Setting the information for printout of individual values.
SIGN.L	Setting the footer for printout of individual values.




Topic	Description
LINE.FEED	Setting line feed for printout of individual values.
ZERO PRT.	Setting the auto print function for printing zero.
COM.SET	Setting the data communication format of the serial interface RS232C.
BAUDRATE	Setting the transfer speed of the serial interface RS232C.
BIT/PAR.	Setting the character format (Bit/Parity) of the serial interface RS232C.
STOPBIT	Setting the character format (stop bit) of the serial interface RS232C.
HD.SHAKE	Setting the transfer protocol (Handshake) of the serial interface RS232C.
RS E.O.L.	Setting the end of line format of the serial interface RS232C.
RS CHAR	Setting the char set of the serial interface RS232C.
USB	Matching the USB interface to a peripheral unit.
USB COM.S.	Setting the data communication format of the USB interface.
USB E.O.L.	Setting the end of line format of the USB interface.
USB CHAR	Setting the char set of the USB interface.
INTERVAL	Selection of the time interval for the simulated print key press.

6.2 Menu Operation


In this Section you will learn how to work with the menu.





Select Menu


- 1 Press «» to activate main menu. The first menu "BASIC" is displayed (except menu protection is active).
- 2 Press «» repeatedly to change menu (Scrolling down/up «+» / «-» keys).
- 3 Press «» to confirm the selection.

Note: The menu selection "BASIC", "ADVANCED" or "INT.FACE" can not be saved. The selection "PROTECT" must be saved.






Select Menu Topic

- Press «»». The next menu topic appears in the display. Each time the «» or the «+» key is pressed, the balance switches to the next menu topic; the «-» key to the previous menu topic.




Change Settings in a Selected Menu Topic

The ">>" flashing symbol in the display indicates selectable options available.

- 1 Press «»». The display shows the current setting in the selected menu topic. Each time «» or «+» is pressed, the balance switches to the next selection; press «-» to the previous selection. After the last selection, the first is shown again.
- 2 Press «»» to confirm the setting. For store the setting see section **Saving Settings and Closing the Menu**.


Change Settings in a Submenu Selection

The same procedure as for menu topics.



Input Principle of Numerical Values

- 1 Press «←» for input of numerical values.
- 2 Press «→» to select a digit or a value (depending on the application). The selected digit or the selected value is blinking.
- 3 For changing digits or values, press «+» to scroll up or «-» to scroll down.
- 4 Press «←» to confirm the input.



Saving Settings and Closing the Menu

- 1 Press «⚙» briefly to leave menu topic.
- 2 Press «←» to execute "SAVE:YES". Changes are saved.
- 3 Press «←» to execute "SAVE:NO". Changes are not saved. To toggle between "SAVE:YES" and "SAVE:NO" press «↔».



Cancel

- For leaving menu topic or menu selection without saving press «C» (one step back in the menu).

Note: If no entry is made within 30 seconds, the balance reverts to last active application mode. Changes are not saved. If changes are made, the balance asks "SAVE:NO".

6.3 Description of Menu Topics

In this Section you will find information regarding the individual menu topics and the available selections.

6.3.1 Main Menu

Selecting the menu.

"BASIC"	The small "BASIC" menu for simple weighing is displayed.
"ADVANCED"	The extended "ADVANCED" menu for further weighing settings is displayed.
"INT.FACE"	The menu "INT.FACE" for all interface parameter settings for peripheral devices e.g. printer is displayed.
"PROTECT"	Menu protection. Protection of balance configurations against unmeant manipulation.
"OFF"	Menu protection is off. (Factory setting)
"ON"	Menu protection is on. The menu BASIC , ADVANCED and INT.FACE are not displayed. This is indicated with "🔒" in the display.

Note:

- The menu selection "BASIC", "ADVANCED" or "INT.FACE" can not be saved.
- To activate "PROTECT" "ON" or "OFF", this selection must be saved.

6.3.2 Basic Menu

"DATE" – Date

Setting the current date according to date format.

Note: A reset of the balance will not change this setting.

"TIME" – Time

Setting the current time according to time format

" +1H "	Set the current time forwards by 1 hour (to adjust summer or winter time). (Factory setting)
" -1H "	Set the current time backwards by 1 hour (to adjust summer or winter time).
" SET TIME "	Enter the current time.

Note: A reset of the balance will not change this setting.

"UNIT 1" – Weight Unit 1

The balance can operate with the following units .

- With approved balances, this menu topic has a fixed setting and cannot be changed.
- Conversion table for weight units see chapter Appendix.

Units:	
g ¹⁾	Gram
kg	Kilogram
mg	Milligram
¹⁾ factory setting	

"UNIT 2" – Weight Unit 2

If it is required to show the weighing results in weighing mode in an additional unit, the desired second weight unit can be selected in this menu topic. Units see "**UNIT 1**". Select "**NO**", if you do not want to use "**UNIT 2**".

"KEY BEEP" – Key Beep

This menu topic allows you to select the volume of the key beep. The according key beep is emitted during the setting.

" MED "	Medium level (Factory setting)
" HIGH "	High level
" OFF "	Beep switched off
" LOW "	Low level

"STAB.BEEP" – Stability Beep

If the unstable symbol disappears, the stability beep becomes active. This menu topic allows you to preselect the volume of the stability beep.

" LOW "	Low level (Factory setting)
" MED "	Medium level
" HIGH "	High level
" OFF "	Beep switched off

"RESET" – Reset Balance Settings

This menu topic allows you to call-up the factory settings.

To toggle between "**YES?**" and "**NO?**" press «» (or «**+**» or «**-**»).

Note: A reset of the balance will not change the "**DATE**", "**TIME**" and "**ZERO RNG**" settings.

6.3.3 Advanced Menu

"WEIG.MODE" – weighing mode settings

This setting can be used to adapt the balance to the weighing mode.

"UNIVERS."	For all standard weighing applications. (Factory setting)
"DOSING"	For dosing liquid or powdery products. With this setting, the balance responds very quickly to the smallest changes in weight.

"ENVIRON." – Environment Settings

This setting can be used to match your balance to the ambient conditions.

"STANDARD"	Setting for an average working environment subject to moderate variations in the ambient conditions. (Factory setting)
"UNSTABLE"	Setting for a working environment where the conditions are continuously changing.
"STABLE"	Setting for a working environment which is practically free from drafts and vibrations.

"CAL" – Adjustment (calibration)

In this menu topic you can preselect the function of the « $\overline{\overline{\square}}$ » key. Your balance can be adjusted with internal or external weights by pressing the « $\overline{\overline{\square}}$ » key. If you have attached a printer to your balance, the data of the adjustment (calibration) are printed out.

"ADJ.OFF"	The adjustment is switched off . The « $\overline{\overline{\square}}$ » key has no function.
"ADJ.INT"	Internal adjustment: adjustment is performed at a keystroke with the built-in weight (depending on the model, see technical data).
"ADJ.EXT"	External adjustment: adjustment is performed at a keystroke with a selectable external weight. Note: This function is not available for approved balances * (depend on selected countries' certification legislation). * except OIML accuracy class I approved models.
"200.00 g"	Defining the external adjustment weight: define the weight of the external adjustment weight (in grams). Factory setting: depends on the model.

"ADJ.CUST.F" – Customer fine adjustment

At this menu topic you can fine-adjust the internal weights. Further information refer to chapter Customer Fine Adjustment.

"EXECUTE"	Start customer fine adjustment "ADJ.CUST.F".
"RESET"	Deactivate customer fine adjustment after confirming with YES?
NO?	No deactivation.
YES?	Confirm to deactivation.

"FACT" – Fully Automatic Adjustment

Fully automatic internal adjustment (calibration) **FACT** (Fully Automatic Calibration Technology) provides fully automatic balance adjustment based on temperature criteria and on preselected time. (depending on the model, see technical data)

"TIME"	Execute FACT (with selected time).
"12:00"	Specify the time for a fully automatic adjustment to take place every day. Factory setting: 12:00 (according to time format)
"OFF"	FACT is switched off .

"FACT.PRT." – Protocol Trigger for Fact

This setting specifies whether an adjustment report should be printed automatically.

Note: This menu topic does not affect the printing of adjustments with an internal or external adjustment weight.

"OFF"	Protocol switched off: if the balance adjusts automatically (FACT), a protocol is not printed out.
"ON"	Protocol switched on: a record is printed out after every automatic adjustment of the balance (FACT). Note: The protocol is printed out without a line for signatures.

"DATE.FORM" – Date Format

This menu topic allows you to preselect the date format.

The following date formats are available:

	Display examples	Printing examples
"DD.MM.Y"	01.02.2009	01.02.2009
"MM/DD/Y"	02/01/09	02/01/2009
"Y-MM-DD"	09-02-01	2009-02-01
"D.MMM Y"	1.FEB.09	1.FEB 2009
"MMM D Y"	FEB.1.09	FEB 1 2009

Factory setting: "DD.MM.Y"

"TIME.FORM" – Time Format

This menu topic allows you to preselect the time format.

The following date formats are available:

	Display examples
"24:MM"	15:04
"12:MM"	3:04 PM
"24.MM"	15.04
"12.MM"	3.04 PM

Factory setting: "24:MM"

"RECALL" – Recall

This menu topic allows you to switch the "RECALL" function on or off. When it is switched on recall stores the last stable weight if the absolute display value was bigger than 10d.

"OFF"	"RECALL" switched off (Factory setting)
"ON"	"RECALL" switched on

Note: The recall value is displayed with an asterisk and cannot be printed.

"STARTUP" – Startup Mode

You can set your balance such that it either immediately starts from the standby mode when you load a weight or it must be switched on with the «ON/OFF» key after which it then performs a display test.

Note: This topic is not visible with approved balances (only available in selected countries).

"QUICK"	"Quickstart": The balance can be started directly from the standby mode and is immediately ready for weighing. You can load the weight in the standby mode and the balance immediately shows the current weighing result. This is the Factory setting Note: Standby mode is available on line powered balances only.
---------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

"FULL"

Start with display test: You must switch on the balance with the «ON/OFF» key. After it has been switched on, it performs a display test for approx. 2 sec. in which all display elements lights up, it shows "WELCOME", software version, maximum load and readability. The balance is ready for weighing.

"SHUTOFF" – Automatic Shutoff

If the automatic shutoff function is activated, the balance automatically switches itself off after a preselected time of inactivity (i.e. with no key being pressed or changes of weight occurring etc.) and is switched to the standby mode.

"A.OFF 10" min	Automatic shutoff after 10 minutes of inactivity. (Factory setting)
"A.OFF –"	Automatic shutoff not activated.
"A.OFF 2" min	Automatic shutoff after 2 minutes of inactivity.
"A.OFF 5" min	Automatic shutoff after 5 minutes of inactivity.

"BCKLIGHT" – Backlight

Under this menu topic, the display backlight can be switched off automatically. If the automatic switch-off is activated, the backlight will turn off automatically after the selected period of inactivity has elapsed. The backlight is reactivated when a key is pressed or the weight is changed.

"B.L. ON"	Backlight is always on . (Factory setting)
"B.L. 30" s	Automatic switch-off after 30 seconds inactivity.
"B.L. 1" min	Automatic switch-off after 1 minute inactivity.
"B.L. 2" min	Automatic switch-off after 2 minutes inactivity.
"B.L. 5" min	Automatic switch-off after 5 minutes inactivity.

"DISPLAY" – Display Settings

This menu topic allows you to adjust brightness and contrast of the display.

"BRIGHTN"	To set the brightness in 1% steps.
"50%"	Factory setting: 50%
"CONTRAST"	To set the contrast in 1% steps.
"75%"	Factory setting: 75%

"AUTOZERO" – Automatic Zero Setting

This menu topic allows you to switch the automatic zero setting on or off.

"ON"	"AUTOZERO" switched on (Factory setting). The automatic zero setting continuously corrects possible variations in the zero point that might be caused through small amounts of contamination on the weighing pan.
"OFF"	"AUTOZERO" switched off. The zero point is not automatically corrected. This setting is advantageous for special applications (e.g. evaporation measurements).

Note: With approved balances, this setting is not available (only available in selected countries).

"ZERO RNG" – Zero Range

This menu topic allows you to set a zero limit for the «→0/T←» key. Up to and including this limit the «→0/T←» key will execute a zero. Above this limit the «→0/T←» key will execute a tare.

"21g" To set the upper limit of the zero setting range as weight in the definition unit of the balance.
(Factory setting: 0.5% of weighing range)
Note: With approved balances, this setting is not available and fixed to 3e (only available in selected countries).

Note: A reset of the balance will not change this setting.

"LANGUAGE" – Language

Factory setting: Generally, the language of the destination country (if available) or English is set.

The following languages are available:

"ENGLISH"	English	"POLSKI"	Polish
"DEUTSCH"	German	"CESKY"	Czech
"FRANCAIS"	French	"MAGYAR"	Hungarian
"ESPANOL"	Spanish	"NEDERL."	Dutch
"ITALIANO"	Italian	"BR.PORTUG."	Brazil Portuguese
"RUSSIAN" РУССКИЙ	Russian		

"ASSIGN:F1" – Assign Application Key F1

At this menu topic you can assign an application to the «**F1**» key. The following applications are available (depending on the model):

"FORMULA"	Formulation / Net-Total (Factory setting)
"COUNTING"	Piece counting
"CHECK"	Checkweighing
"STAT"	Statistics

"ASSIGN:F2" – Assign Application Key F2

At this menu topic you can assign an application to the «**F2**» key. The following applications are available (depending on the model):

"COUNTING"	Piece counting (Factory setting)
"CHECK"	Checkweighing
"STAT"	Statistics
"FORMULA"	Formulation / Net-Total

"ASSIGN:F3" – Assign Application Key F3

At this menu topic you can assign an application to the «**F3**» key. The following applications are available (depending on the model):

"CHECK"	Checkweighing (Factory setting)
"STAT"	Statistics
"FORMULA"	Formulation / Net-Total
"R. TEST"	Routine test
"COUNTING"	Piece counting

"DIAGNOSE" – Diagnostics Application

At this menu topic you can start a diagnostic application. For more information see chapter application "Diagnostics".

The following diagnostics are available:

"REPEAT.T"	Repeatability test (models with internal weights only)
"DISPLAY"	Display test

"KEYPAD T"	Key test
"CAL.MOT. T"	Motor test (models with internal weights only)
"BAL.HIST"	Balance history
"CAL.HIST"	Calibration history
"BAL.INFO"	Balance information
"PROVIDER"	Service provider information

"SERV.ICON" – Service Reminder

This menu topic allows you to switch the service reminder "☞" on or off.

"ON"	Service reminder "☞" switched on (factory setting). You will be informed after a preset time (e.g. one Year or 8000 operating hours) to call service for recalibration. This will be indicated by the flashing service icon: "☞". (Factory setting)
"OFF"	Service reminder "☞" switched off .

"SRV.D.RST" – Service Date Reset

This menu topic allows you to reset service date and hours.

Note: This menu topic is only available if "SERV.ICON" setting "ON" was selected.





To toggle between "YES?" and "NO?" press «☞» (or «+» or «-»)

6.3.4 Interface Menu

"RS232" – RS232C Interface ¹⁾

At this menu topic you can select the peripheral device connected to the RS232C interface and specify how the data is transmitted.

"PRINTER"	Connection to a printer . (Factory setting) Note: <ul style="list-style-type: none"> • Only one printer possible. • See recommended printer settings found in section "Appendix", as well as the printer-specific user's manual.
"PRT.STAB"	If the «☞» key is pressed, the next stable weight value will be printed. (Factory setting)
"PRT.AUTO"	Every stable weight value will be printed, without pressing the «☞» key.
"PRT.ALL"	If the «☞» key is pressed, the weight value will be printed regardless of stability.
"PC-DIR."	Connection to a PC : the balance can send data (as a Keyboard) to the PC used for PC applications e.g. Excel. Note: The balance sends the weight value without the unit to the PC.
"PRT.STAB"	If the «☞» key is pressed, the next stable weight value will be sent followed by an enter. (Factory setting)
"PRT.AUTO"	Every stable weight value will be sent followed by an enter, without pressing the «☞» key.
"PRT.ALL"	If the «☞» key is pressed, the weight value will be sent followed by an enter regardless of stability.
"HOST"	Connection to a PC , Barcode Reader etc.: the balance can send data to the PC and receive commands or data from the PC).
"SEND.OFF"	Send mode switched off. (Factory setting)


"SEND.STB"	If the «  » key is pressed, the next stable weight value will be sent.
"SEND.CONT"	All weight value updates will be sent regardless of stability, without pressing the «  » key.
"SEND.AUTO"	Every stable weight value will be sent, without pressing the «  » key.
"SEND.ALL"	If the «  » key is pressed, the weight value will be sent regardless of stability.
"2.DISPLAY"	Connection of an optional auxiliary display unit Note: The transmission parameters cannot be selected. Settings are automatically set.



Attention:

- If you select 2nd Display "2.DISPLAY", first make sure that no other device is connected at COM1 as a 2nd display. Other devices could be damaged because of the voltage on connector Pin 9. Necessary for powering the 2nd display (see chapter "Interface Specification")

"HEADER" – Options for the Printout Header of individual values


This menu topic allows you to specify the information that is to be printed at the top of the printout for every individual weighing results (after pressing «»).

Note: This menu topic is only available if "PRINTER" setting was selected.

"NO"	The header is not be printed (Factory setting)
"DAT / TIM"	Date and time are printed
"D / T / BAL"	Date, time and balance information (Balance type, SNR, Balance ID) are printed.

Note: Balance ID only if set.


"SINGLE" – Options for Printing out the Result of individual values

This menu topic allows you to specify the information that is to be printed for every individual weighing result (after pressing «»).

Note: This menu topic is only available if "PRINTER" setting was selected.

"NET"	The value of the Net weight from the current weighing is printed (Factory setting)
"G / T / N"	The values of the Gross weight, the Tare weight and the Net weight are printed


"SIGN.L" – Options for the Printout Footer for Signature Line of individual values

This menu topic allows you to set a footer for signature at the bottom of the printout for every individual weighing result (after pressing «»).

Note: This menu topic is only available if "PRINTER" setting was selected.

"OFF"	The signature footer is not be printed. (Factory setting)
"ON"	The signature footer is printed

"LINE.FEED" – Options for Complete the Printout of individual values

This menu topic allows you to specify the number of blank lines to complete the printout (line feed) for every individual weighing result (after pressing «»).

Note: This menu topic is only available if "PRINTER" setting was selected.

"0"	Possible numbers of blank lines: 0 to 99 (Factory setting = 0)
-----	-------------------------------------------------------------------------

"ZERO PRT." – Options for "PRT.AUTO" 1)

This menu topic allows you to specify the auto print function "PRT.AUTO" for printing zero "YES" or "NO".

"OFF"	Zero is not be printed (Zero +/- 3d) (Factory setting)
"ON"	Zero is always printed

Note: This menu topic is only available if "PRT.AUTO" fuction of the "PRINTER" or "PC-DIR." was selected.

"COM.SET" – Options for the Data Communication Format (RS232C)("HOST") 1)

This menu topic allows you to set the data format depending on which peripheral device is connected.

Note: This menu topic is only available if "HOST" setting was selected.

"MT-SICS"	The MT-SICS data transfer formats is used. (Factory setting) For more information see section "MT-SICS Interface Commands and Functions".
-----------	-------------------------------------------------------------------------------------------------------------------------------------------------------

"MT-PM" The following PM balance commands are supported:

S	Send value
SI	Send immediate value
SIR	Send immediate value and repeat
SR	Send value and repeat
SNR	Send next value and repeat
T	Tare
TI	Tare immediatelly
B	Base *)
MI	Modify ambient vibration
MZ	Modify Auto Zero
M	Modified settings reset
ID	Identify
CA	Calibrate
D	Display (only symbol N and G available)

*) Limitation:

- Negative values are limited up to the current tare value.
- B command is additive.
- The sum of the B values plus the previous tare value, before a "TA", "T" or "Z" is sent, must be less than the total weighing range.

"SART" The following Sartorius commands are supported:

K	Ambient conditions: very stable
L	Ambient conditions: stable
M	Ambient conditions: unstable
N	Ambient conditions: very unstable
O	Block keys
P	Print key (print, auto print; activate or block)
Q	Acoustic signal
R	Unblock keys
S	Restart/self-test
T	Tare key
W	Calibration/adjustment (depending on the menu setting)

*)

Z	Internal calibration/adjustment **)
f0_	Function key (F)
f1_	Function key (CAL)
s3_	C key
x0_	Perform internal calibration **)
x1_	Print balance/scale model
x2_	Print weighing cell serial number
x3_	Print software version

*) may be inaccessible on verified balances/scales

**) only on models with built-in motorized calibration weight

Functionality mapping

"HOST" settings:	Sartorius printer settings:
"SEND.OFF"	not applicable
"SEND.STB"	manually print with stability
"SEND.ALL"	manually print without stability
"SEND.CONT"	automatically print without stability
"SEND.AUTO"	similar applicable to automatically print when load is changed

"BAUDRATE" – Baud rate RS232C ¹⁾

This menu topic allows you to match the data transmission to different serial RS232C receivers. The baud rate (data transfer rate) determines the speed of transmission via the serial interface. For problem-free data transmission the sending and receiving devices must be set at the same value.

The following settings are available:

600 bd, 1200 bd, 2400 bd, 4800 bd, 9600 bd, 19200 and 38400 bd. (default: **9600 bd**)

Note:

- Not visible for 2nd display.
- Each device has separate settings.

"BIT/PAR." – Bit/Parity RS232C ¹⁾

At this menu topic you can set the character format for the attached RS232C serial peripheral device.

"8/NO"	8 data bits/no parity (Factory setting)
"7/NO"	7 data bits/no parity
"7/MARK"	7 data bits/mark parity
"7/SPACE"	7 data bits/space parity
"7/EVEN"	7 data bits/even parity
"7/ODD"	7 data bits/odd parity

Note:

- Not visible for 2nd display.
- Each device has separate settings.

"STOPBIT" – Stop Bits RS232C ¹⁾

At this menu topic you can set the stop bits of the transmitted data to different RS232C serial receivers.

"1 BIT"	1 Stop bit (Factory setting)
"2 BITS"	2 Stop bits

Note:

- Not visible for 2nd display.
- Each device has separate settings.

"HD.SHAKE" – Handshake RS232C ¹⁾

This menu topic allows you to match the data transmission to different RS232C serial receivers.

"XON/XOFF"	Software handshake (XON/XOFF) (Factory setting)
"RTS/CTS"	Hardware handshake (RTS/CTS)
"OFF"	No handshake

Note:

- Not visible for 2nd display.
- Each device has separate settings.

"RS.TX.E.O.L." – End of Line RS232C ¹⁾

At this menu topic you can set the "End of Line" character of the transmitted data to different RS232C serial receivers.

"(CR)(LF)"	<CR><LF> Carriage Return followed by Line feed (ASCII-Codes 013+010) (Factory setting)
"(CR)"	<CR> Carriage Return (ASCII-Code 013)
"(LF)"	<LF> Line feed (ASCII-Code 010)
"(TAB)"	<TAB> Horizontal tab (ASCII-Code 011), only settable if PC-DIR. is selected.

Note:

- Not visible for 2nd display.
- Each device has separate settings.

"RS CHAR" – Char Set RS232C ¹⁾

At this menu topic you can set the "Character Set" of the transmitted data to different RS232C serial receivers.

"IBM/DOS"	Char Set IBM/DOS (Factory setting)
"ANSI/WIN"	Char Set ANSI/WINDOWS

Note:



- Not visible for 2nd display.
- Each device has separate settings.







"USB" – USB Interface

At this menu topic you can select the mode of the "USB Device" interface and specify how the data is transmitted.

Note:

- DISCONNECT THE USB CONNECTION FROM THE BALANCE PRIOR TO CHANGE THE SETTINGS.
- This port is not usable for printers or displays.

"PC-DIR."	Connection to a PC : the balance can send data (as a Keyboard) to the PC used for PC applications e.g. Excel. Note: The balance sends the weight value without the unit to the PC.
"SEND.OFF"	Send mode switched off (Factory setting)
"SEND.STB"	If the «  » key is pressed, the next stable weight value will be sent.
"SEND.CONT"	All weight value updates will be sent regardless of stability, without pressing the «  » key.

"SEND.AUTO"	Every stable weight value will be sent, without pressing the «  » key.
"SEND.ALL"	If the «  » key is pressed, the weight value will be sent regardless of stability.
"HOST"	Connection to a PC , Barcode Reader etc.: the balance can send data to the PC and receive commands or data from the PC).
"SEND.OFF"	Send mode switched off. (Factory setting)
"SEND.STB"	If the «  » key is pressed, the next stable weight value will be sent.
"SEND.CONT"	All weight value updates will be sent regardless of stability, without pressing the «  » key.
"SEND.AUTO"	Every stable weight value will be sent, without pressing the «  » key.
"SEND.ALL"	If the «  » key is pressed, the weight value will be sent regardless of stability.

"USB COM.S." – Options for the Data Communication Format (USB)

This menu topic allows you to set the data format depending on which peripheral device is connected.

"MT-SICS"	The MT-SICS data transfer formats is used. (Factory setting) For more information see section "MT-SICS Interface Commands and Functions".
"MT-PM"	The following PM balance commands are supported: S Send value SI Send immediate value SIR Send immediate value and repeat SR Send value and repeat SNR Send next value and repeat T Tare TI Tare immediately B Base *) MI Modify ambient vibration MZ Modify Auto Zero M Modified settings reset ID Identify CA Calibrate D Display (only symbol N and G available) *) Limitation: <ul style="list-style-type: none"> • Negative values are limited up to the current tare value. • B command is additive. • The sum of the B values plus the previous tare value, before a "TA", "T" or "Z" is sent, must be less than the total weighing range.
"SART"	The following Sartorius commands are supported: K Ambient conditions: very stable L Ambient conditions: stable M Ambient conditions: unstable N Ambient conditions: very unstable

O	Block keys
P	Print key (print, auto print; activate or block)
Q	Acoustic signal
R	Unblock keys
S	Restart/self-test
T	Tare key
W	Calibration/adjustment (depending on the menu setting *)
Z	Internal calibration/adjustment **)
f0_	Function key (F)
f1_	Function key (CAL)
s3_	C key
x0_	Perform internal calibration **)
x1_	Print balance/scale model
x2_	Print weighing cell serial number
x3_	Print software version

*) may be inaccessible on verified balances/scales

**) only on models with built-in motorized calibration weight

Functionality mapping

"HOST" settings:	Sartorius printer settings:
"SEND.OFF"	not applicable
"SEND.STB"	manually print with stability
"SEND.ALL"	manually print without stability
"SEND.CONT"	automatically print without stability
"SEND.AUTO"	similar applicable to automatically print when load is changed

"USB E.O.L." – End of Line USB

At this menu topic you can set the "End of Line" character of the transmitted data to USB device.


"(CR)(LF)"	<CR><LF> Carriage Return followed by Line feed (ASCII-Codes 013+010) (Factory setting)
"(CR)"	<CR> Carriage Return (ASCII-Code 013)
"(LF)"	<LF> Line feed (ASCII-Code 010)
"(TAB)"	<TAB> Horizontal tab (ASCII-Code 011), only settable if PC-DIR. is selected.

"USB CHAR" – Char Set USB

At this menu topic you can set the "Character Set" of the transmitted data to USB device.

"ANSI/WIN"	Char Set ANSI/WINDOWS (Factory setting)
"IBM/DOS"	Char Set IBM/DOS

"INTERVAL" – Print Key Simulation

At this menu topic you can activate a simulation of the «» key. "INTERVAL" simulates a print key press every x seconds.

Range:	0 to 65535 seconds
0 sec:	disables the print key simulation

Factory setting: 0 sec

Note: The executed action is according to the configuration of the print key. (see interface setting)

1) Note for 2nd RS232C Interface

- If an optional 2nd interface is installed, the menu topic is displayed for each interface, e.g
"BAUDRATE.1" for standard interface
"BAUDRATE.2" for optional 2nd interface
 - Only one printer can be set if two RS232 interfaces are existing.
-


7 Application "Formulation" (Net Total Formulation)




The "Formulation" (Net Total) application allows you to

- weigh in (add and store) up to 999 individual component weights and displays the total. If a printer is connected, the component weights are printed individually and as a total.
- tare/pre-tare and store up to 999 container weights and displays the total. If a printer is connected, the tare weights are printed out individually and as a total.

Requirement: The function "FORMULA" must be assigned to an «Fx» key (see advanced menu topic "ASSIGN:Fx"). Connect a printer or a PC if present.



- 1 Activate function formulation "FORMULA" by pressing and holding the appropriate assigned «Fx» key.
- 2 Press «←» to continue formulation weighing. For a new formulation press «↶» (or «+» or «-») to select "Yes" and press «←» to clear the memory.
Note: If the memory is already cleared (sample and tare/pre-tare counter is zero) the memory clear question will be not displayed.




Tare container (if used):

- 1 Press «→0/T←» to zero or tare the balance if needed.
- 2 Place the empty container on the weighing pan.
- 3 Press «→0/T←». The container is tared and the tare count "- T 1 -" is displayed and the tare weight is printed.

Note:

- If you pre-tare via MT-SICS (e.g. bar code reader) "- PT1 -" is displayed.
- Zero range setting (menu topic "ZERO RNG") has no effect. The zero-limit is less than or equal 10d.



Weighing the first component weight:

- 1 Load the first component weight.
- 2 Press «←». The display briefly shows the component count "- 1 -", the current weight is stored as sample and the component weight is printed. The display is set back to zero.


Weighing further component weights:

The same procedure as for the first component weight with the same or new container).

- 1...999 sample values are possible.
- max 999 tare values are possible.
- max 999 pre-tare values are possible.

Results:





- If the numbers of sample are greater than or equal to 2, press «≡», the results are displayed and printed.











Displayed results:





- 1 Press «←» to show the next statistical value.
- 2 Press «C» to cancel displaying results and to continue weighing next component.

0.5 seconds

number of samples  N  *  8 

sum of all tare values (T and PT)  T.TOTAL  *  452.76 g 

sum of all component gross weight values  G.TOTAL  *  546.79 g 

sum of all component net weight values  N.TOTAL  *  94.03 g 

Printout:

```

----- Formulation -----
21.Jan 2010           12:56


METTLER TOLEDO

Balance Type      PH4002S
SNR               1234567890
-----
1 T                10.33 g
1 N                8.85 g
2 N                9.23 g
2 T                10.84 g
3 N                7.43 g
.
.
n                  8
T Total           452.76 g
G Total           546.79 g

N Total           94.03 g
-----
    
```


Function "FILL UP"

This function allows you to add an additional component weight to the total weight of all components to reach a desired target weight (Fill up).



Starting the fill up function.

- Activate function "FILL UP" by pressing «+».
- Deactivate function "FILL UP" by pressing «-».




Filling up with an additional component weight:

- ▶ The last total of the component weights is displayed.
- 1 Add component weight until the desired target weight is reached.
- 2 Press «←» to confirm.
- ⇒ The display briefly shows the next component count marked with "F", the current weight is stored as sample and the component weight is printed. The display is set back to zero.


Filling up further additional component weights:

The same procedure, beginning with starting up the "FILL UP" function.

Terminate the application

Press and hold «» to terminate the application and to return to the weighing application.

8 Application "Piece Counting"




The **"Piece Counting"** application allows you to determine the number of pieces put on the weighing pan.

Requirement: The function **"COUNTING"** must be assigned to an «Fx» key (see advanced menu topic **"ASSIGN:Fx"**, factory setting: F1).




- Activate function **"COUNTING"** by pressing and holding the appropriate assigned «Fx» key (factory setting: F1).



Piece Counting first requires the setting of a reference weight, there are 4 possibilities:


- A** Setting the reference **by multiple pieces with fix reference values.**
- B** Setting the reference **by multiple pieces with variable reference values.**
- C** Setting the reference **for 1 piece in weighing mode.**
- D** Setting the reference **for 1 piece in manual mode.**



Setting possibility

A Setting the reference by multiple pieces with fix reference values


- 1 Select a number of reference pieces by scrolling with «». Possible numbers* are 5, 10, 20 and 50.
* with approved balances in selected countries: min 10
- 2 Press «» to tare. If using: place empty container on the weighing pan first or tare again.
- 3 Add the selected number of reference pieces to container.
- 4 Press «» to confirm.



Setting possibility

B Setting the reference by multiple pieces with variable reference values

- 1 Select **"VAR.REF"** by scrolling with «». Press «» to confirm.
- 2 Select a number of reference pieces by scrolling up («+» key) or down («-» key). Speed up by press and hold. Possible numbers* are 1 to 999.
* with approved balances in selected countries: min 10
- 3 Press «» to tare. If using: place empty container on the weighing pan first or tare again.
- 4 Add the selected number of reference pieces to container.
- 5 Press «» to confirm.




Setting possibility

C Setting the reference for one piece in weighing mode

- 1 Select "PCS.WGT" by scrolling with «G».
- 2 Press «→0/T←» to tare. If using: place empty container on the weighing pan first or tare again.
- 3 Add one reference piece to container. The weight of one piece is displayed.
- 4 Press «←» to confirm.

Note: With approved balances, this setting is not available in selected countries.




Setting possibility

D Setting the reference for one piece in manual mode

- 1 Select "PCS.WGT" by scrolling with «G».
- 2 Press «←» to confirm.
- 3 Enter the final reference one piece weight by scrolling up («+» key) or down («-» key). Speed up by press and hold.
- 4 Press «←» to confirm.

Note: With approved balances, this setting is not available in selected countries.




Switching between manual mode and weighing mode

- Press «G» to switch between manual and weighing mode.

Note: By switching from weighing mode to manual mode the weight value will be transferred and can be changed manually.

Note: If without any key press within 60 seconds, the balance returns to the previous active application. Press «C» to cancel and returns to the previous active application.

On completion of the setting procedure, your balance is ready for piece counting.




Switching between piece count and weight display.

You can use the «G» key at any time to switch the display between piece display, weighing unit "UNIT 1", "RECALL" value (if activated) and weighing unit "UNIT 2" (if different from "UNIT 1").

Note:

- The "RECALL" value is displayed with an asterisk (*) and icon "M" and can not be printed.
- Take into account minimum values: min. reference weight = 10d (10 digits), min. piece weight* = 1d (1 digit)!
* with approved balances in selected countries: min 3e
- The current reference weight remains stored until the reference setting is changed.

Terminate the application


Press and hold «» to terminate the application and to return to the weighing application.

9 Application "Check Weighing"




The **"Check weighing"** application allows you to check the deviation of a sample weight within a tolerance limit to a reference target weight.

Requirement: The function **"CHECK"** must be assigned to an «Fx» key (see advanced menu topic **"ASSIGN:Fx"**, factory setting: F3).



- Activate function **"CHECK"** by pressing and holding the appropriate assigned «Fx» key (factory setting: F3).



Step 1: Check Weighing first requires the setting of a reference weight that should corresponds to the nominal weight, there are 2 possibilities:


- 1A** Setting the reference **in manual mode** (enter nominal weight).
- 1B** Setting the reference **in weighing mode** (weigh nominal weight).

Step 2: Check weighing needs the upper and lower limits, there are 2 possibilities::

- 2A** Setting the **upper and lower limits in percentage**.
- 2B** Setting the **upper and lower limits by weight**.

Step 3: Setting tolerance beep


- 3** Activate or deactivate **tolerance beep**.



Step 1, setting possibility:

1A **Setting the reference by manual mode** (enter nominal weight)


- 1 Press «←» to activate manual mode.
- 2 Select the reference target weight by scrolling up («+» key) or down («-» key). Speed up by press and hold.
- 3 Press «←» to confirm the nominal weight.



Step 1, setting possibility:

1B **Setting the reference by weighing mode** (weigh nominal weight)


- 1 Press «→0/T←» to tare the balance and to activate the weighing mode. If using: place empty container on the weighing pan first or tare again.
- 2 Load the nominal weight.
- 3 Press «←» to confirm the nominal weight.



Switching between manual mode and weighing mode

- Press «↻» to switch between manual mode and weighing mode.

Note: By switching from weighing mode to manual mode the weight value will be transferred and can be changed manually.




Step 2, setting possibility:

2A Setting the upper and lower limits (in percentage):

- 1 Press «←» to start setting.
- 2 Press «←» to confirm the default limit of +/- 2.5 % or enter the limit value by scrolling up («+» key) or down («-» key). Press «←» to confirm the limits.


Note: Press «↶» to switch between "UNIT 1" and Unit "%".



Step 2, setting possibility:


2B Setting the upper and lower limits by weight:

- 1 Press «←» to start setting.
- 2 Press «↶» to switch to **UNIT 1**.
- 3 Press «←» to confirm the default limit or enter the limit value by scrolling up («+» key) or down («-» key). Press «←» to confirm the limits.



Switching between percentage and weight unit 1

- Press «↶» to switch between setting in percentage and setting by weight .



Step 3:

3 Setting tolerance beep:

The tolerance beep indicates whether the weighing sample lies within the tolerance by beeping three times.


Note: The beep level corresponds to the setting in menu topic "STAB.BEEP" (Basic menu). If "STAB.BEEP" is set to "OFF", the tolerance beep level is medium.

- To activate tolerance beep press «←». To deactivate tolerance beep press «↶» to select "NO" and press «←».

Note:

- If without any key press within 60 seconds, the balance returns to the previous active application. Press «C» to cancel.
- The nominal weight must be at least 10 digit.

On completion of the setting procedure, your balance is ready for checkweighing.




Weighing-in-Aid

The Weighing-in-Aid helps you quickly determine the position of the sample weight regarding the tolerance.

- 1 Lower limit
- 2 Target weight
- 3 Upper limit

Terminate the application


Press and hold «» to terminate the application and to return to the weighing application.

10 Application "Statistics"



The **"Statistics"** application allows you to generate statistics of a series of weighing values. 1 to 999 values are possible.


Requirement: The function **"STAT"** must be assigned to an **«Fx»** key (see advanced menu topic **"ASSIGN:Fx"**). Connect a printer or a PC if present.



- 1 Activate function **"STAT"** by pressing and holding the appropriate assigned **«Fx»** key.
- 2 To continue the last statistics press **«←»**. For a new statistical evaluation press **«↻»** to select "Yes" and press **«←»** to clear the memory.

Note:

If the memory is already cleared (at the first start of this application or sample counter is 0) the memory clear question will be not displayed.




Weighing the first sample weight:

- 1 Press **«→0/T←»** to zero/tare the balance if needed.
- 2 Load the first sample weight.
- 3 Press **«←»**. The display shows the sample count **"- 1 -"** and the current weight is stored as sample and the weight is printed out.
Note: When the sample counter is displayed you may press **«C»** to undo (drop) this sample.
- 4 Unload the first sample weight.

Weighing further sample weights:

The same procedure as for the first sample weight.

- 1...999 samples are possible.
- The next value will be accepted if the sample weight is in the range 70% –130% of the current average value. **"OUT OF RANGE"** will be displayed if the sample is not accepted.



Results:

- If the numbers of sample are greater than or equal to 2, press **«≡»**, the results are displayed and printed.

Displayed results:


- 1 Press **«←»** to show the next statistical value.
- 2 Press **«C»** to cancel displaying results and to continue weighing next sample.

		0.5 seconds
number of samples		
average		
standard deviation		
relative standard deviation		
lowest value (minimum)		
highest value (maximum)		

Displayed results:

- 1 Press «**←**» to show the next statistical value.
- 2 Press «**C**» to cancel displaying results and to continue weighing next sample.

different between the minimum and the maximum sum of all values



Printout:

```
----- Statistics -----  
21.Jan 2010      12:56  
  
METTLER TOLEDO  
  
Balance Type      PH4002S  
SNR               1234567890  
-----  
1                46.36 g  
2                55.81 g  
3                47.49 g  
4                53.28 g  
5                49.71 g  
n                5  
x                50.530 g  
s dev           3.961 g  
s rel           7.84 g  
Min.            46.36 g  
Max.            55.81 g  
Diff            9.45 g  
Sum             252.65 g  
-----
```

11 Application "Routine Test"



The **"Routine Test"** application allows you to determine the sensitivity of the balance. More about periodic sensitivity tests (routine tests) see: **GWP®** (Good Weighing Practice) on www.mt.com/gwp.


GWP gives clear recommendation for routine testing:

- how should I test my balance?
- how often?
- where can I reduce efforts?

More about test weights see www.mt.com/weights.

Requirement:


- The function **"R. TEST"** must be assigned to **«F3»** key (see advanced menu topic **"ASSIGN:F3"**).
- It is recommended to connect a printer or a PC to the balance for showing the results.



- 1 Activate function **"R. TEST"** by pressing and holding the assigned **«F3»** key.
- 2 Select "No" (no tare weight used). If a tare weight is used during the test select "Yes" (use a tare weight). To toggle between "Yes" and "No" use **«G»** (or **«+»** or **«-»**).
- 3 Press **«←»** to confirm the selection.

Note:


- It is recommended to test the sensitivity without tare load. (factory setting "No").
- If using tare: Make sure that tare weight plus test weight is not exceeding max. load.



Setting the reference test weight value

The default value of the test weight: Next smaller OIML weight than the maximum load of your balance according to the GWP® recommendation.

- 1 For changing the value, press **«+»** to scroll up or **«-»** to scroll down. Progressing speed by press and hold.
- 2 Press **«←»** to confirm the value.




Setting the Control Limit

The default value of the control limit:

Test weight x weighing process tolerance / 2

Example: 5000 g x 0.1% / 2 = 2.50 g.

- 1 For changing the value, press **«+»** to scroll up or **«-»** to scroll down. Progressing speed by press and hold.
- 2 Press **«←»** to confirm the value.




Setting the Warning Limit

The default value of the warning limit:
 Warning limit = control limit / safety factor
 Example: 2.5 g / 2 = 1.25 g.

- 1 For changing the value, press «+» to scroll up or «-» to scroll down. Progressing speed by press and hold.
- 2 Press «←» to confirm the value.

Note: The default values of control limit and the warning limit are evaluated according the GWP recommendation. These are based under the assumption that the weighing process tolerance is 0.1% and the safety factor is 2.



On completion of the setting procedure, your balance is ready for the routine test procedure.

Note: The test weight must be acclimatized to the ambient temperature of the balance.

- 1 Press «←» to start the test.
- 2 Follow the instructions on the display. If the test weight value is flashing: Load the test weight (displayed value).

The printout starts after the weighing pan is unloaded.

Exit the current test procedure:

- Press and hold « $\Delta\Delta$ », «F1», «F2» for executing a new application.

Printout:

```

----- Routine Test -----
21.Jan 2010           12:56

METTLER TOLEDO

Balance Type      PH6002S
SNR               1234567890

Sensitivity:
Test weight      5000.00 g
Value            5000.11 g
Warning L.       1.25 g
Control L.       2.50 g
Warning L.       OK
Control L.       OK

Signature

.....
-----
  
```

What if Warning Limit or Control Limit are "FAILED"?

The "SOP for Periodic Sensitivity Tests (Routine Tests)" provides information about measures when routine tests fail. Find a download version of these SOPs on www.mt.com/gwp, link "GWP® The Program / Routine Operation".

Content of SOP:

- Preparation
- Test procedure

- Evaluation
- Deviation
 - If Warning Limit "**FAILED**"
 - If Control Limit "**FAILED**"

12 Application "Diagnostics"



The "**Diagnostics**" application allows you to carry out predefined diagnostics tests and to view or print predefined sets of balance information. This diagnostics tool helps you find errors faster and more efficiently.

Requirement: A printer or a PC is connected to the balance for showing the results.

- 1 Activate "**ADVANCED**" menu. (See section menu operation)
- 2 Activate function "**DIAGNOSE**" by pressing «←↓».
- 3 Use «←↶» to select appropriate tests.

12.1 Repeatability Test

The repeatability test allows you to repeat tests with internal weight for a given number of times.

Note: On models with internal weights only.

- 1 Press «←↓» to activate repeatability test "**REPEAT.T**". "**R. TST. 10**" appears on the Display.
- 2 Enter the number of times (blinking) by pressing «+» or «-». Possible values are 5, 10 (default), 20, 50, 100 times.
- 3 Press «←↓» to start the test. The message "**RUNNING REPEAT TEST**" is displayed till the tests are completed.
- 4 Press «←☰» to print the test information..
- 5 Press «←↓» to scroll forward through the displayed list.
- 6 Press «C» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample Information Displayed:

Displayed for 0.5 s	Display
"S DEV"	* 0.004 g
"MAX. TEMP"	21.2 °C
"MIN. TEMP"	21.0 °C
"MEAN. TEMP"	21.1 °C
"TOT.TIME"	00:01:26

Sample Printout:

```
-- Repeatability Test --
21.Jan 2010          11:34

METTLER TOLEDO

Balance Type      PH6002S
SNR               1234567890
SW               V1.00
Temperature       21.3 °C
No. of tests      10
-----
1. Time           00:00:00
1. Temp.          21.3 °C
2. Time           00:00:04
2. Temp.          21.3 °C
.
.
.
-----
s Dev.            0.004 g
Max Temp.         21.2 °C
Min Temp.         21.0 °C
Mean Temp.        21.1 °C
Total Time        00:00:44
-----
```



Examples:

Repeatability test is a tool to do functional check with the balance. It may be performed:

- **To check function of balance**
 - during installation to store print out with installation documents.
 - after preventative maintenance to store print out with installation maintenance report.
 - when remarkable decrease of weighing performance occurs, so that you can email/fax print out to service support provider for diagnose purposes.
- **To develop the optimal environment settings** (see menu topic "**ENVIRON.**"). Measure the time you need to perform repeatability test with each "**STABLE**", "**STANDARD**" and "**UNSTABLE**" setting. The setting with the fastest total time suits best for the existing environmental conditions.

12.2 Display Test

The display test allows you to test the display of the balance.

- 1 Press «» to start "**DISPLAY**".
All possible segments and icons on the display will illuminate.
- 2 Press «» to print the test information.
- 3 Press «**C**» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample Printout:

```
----- Display Test -----
21.Jan 2010      11:34

METTLER TOLEDO






Balance Type      PH204S
SNR               1234567890
SW                V1.00
Display Test      DONE
-----
```

12.3 Key Test

The key test allows you to test the keys of the balance.

- 1 Press «←↓» to start "KEYPAD T".
- 2 The message "KEY TEST - PRESS KEY TO BE TESTED" is displayed scrolling during the duration of the key test. Press every Key briefly. Each press of a key beeps and echoes with "OK" on the display.
- 3 Second press «C» key to print the test information. The test procedure will be cancelled and the balance will return to the topic "DIAGNOSE". If a key has not been tested before printing, then the test results will be indicated with a "----" line.

Sample Information Displayed:

Key	Display
«  »	1/10 D OK
«  »	MENU OK
«  »	CAL OK
«  »	PRINT OK
«-»	MINUS OK
«+»	PLUS OK
«  »	TOGGLE OK
«←↓»	ENTER OK
«C»	C OK
«→0/T←»	O/T OK

Sample Printout:

```
----- Key Test -----
21.Jan 2010          11:34


METTLER TOLEDO

Balance Type        PH204S
SNR                 1234567890
SW                  V1.00
1/10 d Key         OK
Menu Key           OK
Cal Key            OK
Print Key          OK
Minus Key          OK
Plus Key           OK
Toggle Key         OK
Enter Key          OK
Zero/Tare Key      OK
Cancel Key         OK
-----
```

12.4 Motor Test

The motor test allows you to test the calibration motor of the balance.

Note: On models with internal weight only.

- 1 Press «←|» to start "**CAL.MOT. T**".
"**RUNNING**" is displayed during the Motor Test. A motor test is deemed successful when all the motor positions have been successfully tested. At the end of the test, the test information will be printed.
- 2 Press «» for printout.
- 3 Press «C» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample Printout:


```
----- Motor Test -----
21.Jan 2010          11:34

METTLER TOLEDO

Balance Type        PH204S
SNR                 1234567890
SW                  V1.00
Motor Test          OK
-----
```

12.5 Balance History

The balance history function allows you to view and print the history of the balance.

- 1 Press «←|» to start "**BAL.HIST**".
- 2 Press «» for printout.
- 3 Press «←|» to scroll forward through the displayed list of balance history information.
- 4 Press «C» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample Information Displayed:

Information	Display
Operation Time (year:day:hour)	00:018:04
Total load kg	115.7191 kg
Number of weighings	1255
Number of key pressed	4931
Number of motor movements	1012
Backlight time (year:day:hour)	00:018:04
Next service due date	01:01:2010

Sample Printout:

```

--- Statistical Info ---
21.Jan 2010      11:34


METTLER TOLEDO

Balance Type     PH4002S
SNR              1234567890
SW               V1.00
-----
Operating time
                18d 4h
Total weight loaded
                115.7191 kg
Number of weighings
                1255
Number of key presses
                4931
Motor movements
                1012
Backlight operating time
                18d 4h
Next service due date
                01.01.2011
-----

```

12.6 Calibration History

The "Calibration History" function allows you to view and print information of the last 30 (thirty) balance adjustment. Adjustments made by a service technician and normal user are counted together.

- 1 Press «←|» to start "**CAL.HIST**".
- 2 Press «» for printout.
- 3 Press «←|» key to scroll forward through the displayed list of Adjustments history information.
- 4 Press «C» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample Information Displayed:

Note	Display	
S = External adjusted service	05:03:09S	01
	-3 PPM	
F = FACT	05:03:09F	02
	2 PPM	

Note	Display	
	.	.
	.	.
	.	.
I = Internal adjusted	04:03:09I	28
	-1 PPM	
E = External adjusted user	03:03:09E	29
	4 PPM	
F = FACT	02:03:09F	30
	1 PPM	

Sample Printout:

```

----- Calibration -----
05.Mar 2009      11:34

METTLER TOLEDO

Balance Type      PH204S
SNR              1234567890
SW               1.50
-----

01 05.Mar 2009   11:34
External ADJ SERVICE
                23.5°C
Diff            -3ppm
-----

02 05.Mar 2009   09:00
FACT
                22.4°C
Diff            2ppm
-----

.
.
.
28 03.Mar 2009   10:59
Internal ADJ
                22.6°C
Diff            -1ppm
-----


29 02.Mar 2009   16:34
External ADJ USER
                24.6°C
Diff            4ppm
-----

30 02.Mar 2009   18:36
FACT
                22.4°C
Diff            1ppm
-----

```

12.7 Balance Information

The balance information function allows you to view and print information about your balance.

- 1 Press «←» to start "BAL.INFO".
- 2 Press «» for printout.
- 3 Press «←» to scroll forward through the displayed list of Balance information.

4 Press «**C**» to cancel the test procedure. The balance will return to the topic "**DIAGNOSE**".

Sample information displayed:

Information	Display
Balance type	TYPE MS6002S
Max. load	MAX 6200 g
Software platform	PLATFORM RAINBOW
Serial number	SNR 1234567890
Type definition number	TDNR 9.6.3.411
Software version	SOFTWARE V1.00
Cell ID	CELL ID 1172400044
Cell type	CELL TYPE MMAI6000G2
Tolerance revision number	TOLERANCE NO2
Language	LANGUAGE ENGLISH

Sample Printout:

```

-- Balance Information -
05.Mar 2010      11:34

METTLER TOLEDO

Balance Type      PH6002S
SNR               1234567890
SW               V1.00
Max              6200 g
Platform         Rainbow
TDNR             9.6.3.411.2-03
Cell ID          1172400044
Cell Type        MMAI6000G2
Tolerance Rev.no.  2
Language         English
-----

```

12.8 Service Provider Information

The service provider Information function allows you to print information about your service provider.

- 1 Press «**←**» to start "**PROVIDER**". The service provider information will be displayed.
- 2 Press «**☰**». The service provider information will be printed and the balance will return to the topic "**DIAGNOSE**".

Sample Printout:

```

--- Service Provider ---
21.Jan 2009      11:34

METTLER TOLEDO
Im Langacher
CH-8606 Greifensee
Switzerland
(+41) 044 944 22 11
-----

```

13 Communication with Peripheral Devices

13.1 Function PC-Direct

The numerical value displayed at the balance can be transferred to the cursor position in Windows Applications (e.g. Excel, Word) as by typing with the keyboard.

Note: The units will not be transferred.

Requirements

- PC with one of the Microsoft Windows® operating system 32bit/64bit: XP (SP3), Vista (SP2), Win 7 (SP1) or Win 8.
- Serial interface RS232 or USB.
- Administrator rights for installing software (for USB not required).
- Windows Application (e.g. Excel).
- Balance to PC connection with cable RS232 or USB.

Settings at the balance:

Attention

- DISCONNECT THE USB CONNECTION FROM THE BALANCE PRIOR TO CHANGE THE SETTINGS.
- USB does not work with keyboards where the "Shift" key must be pressed for entering numbers.

Balance Interface Settings (see Interface Menu):

- Topic "**RS232**" or "**USB**": set "**PC-DIR.**" and select the most appropriate option for the desired weighing result.
- Topic "**RS.TX.E.O.L.**"/"**RS E.O.L.**" or "**USB E.O.L.**"/"**USB E.O.L.**":
 - set <TAB> to write into the same row (e.g. in Excel).
 - set <CR><LF> to write into the same column (e.g. in Excel).
- Save changes.

Settings at the PC:

Installing SerialPortToKeyboard

Operation of PC-Direct via serial port RS232 requires the installation of **SerialPortToKeyboard** on your host computer.

Using CD-ROM

- 1 Insert the product CD in the CD/DVD drive of the host computer.
- 2 Double click the folder **SerialPortToKeyboard**.

Using internet

- 1 Go to the site <http://www.mettler-toledo-support.com>.
- 2 Log in to the METTLER TOLEDO Balance Support Site (registration with the serial number of a METTLER TOLEDO instrument required).
- 3 Click **Customer Support**
- 4 Click appropriate product folder and save the program file **SerialPortToKeyboard.exe** on your specified storage location.

Installing procedure


- 1 Right-click on **SerialPortToKeyboard.exe** and select **Run as Administrator** from the menu.
- 2 Follow the installer's instructions.

Settings for SerialPortToKeyboard

- 1 Select the serial port (COM) to be used for connection with the balance.
- 2 Set the baud rate to 9600.
- 3 Activate "Connect"

Note

- The window can be minimized.
- Closing of the window terminates the session.



Checking Operation

- 1 Start **SerialPortToKeyboard** (RS232)
- 2 Start Excel (or another application) at the PC.
- 3 Activate a cell in Excel.

According to your selected "**PC-DIR.**" option, the displayed values will appear e.g. in the column one after the other one in the different rows.

13.2 USB Device Interface

To perform the functionality "**HOST**" with a PC equipped only with a USB Interface, you have to assign an appropriate USB Driver on the PC first.

Requirements

- Balance with USB Device Interface.
- PC with one of the Microsoft Windows® operating system 32bit/64bit: XP (SP3), Vista (SP2), Win 7 (SP1) or Win 8.
- Administrator rights for installing software.
- PC to balance USB connection cable.

Installing USB Driver on the PC:

Using CD-ROM

- 1 Insert the product CD in the CD/DVD drive of the host computer.
- 2 Double click the folder **USB Driver**.
- 3 Click **USBDriverInstaller.exe**.

Using internet

- 1 Connect to the Internet
- 2 Go to the site <http://www.mettler-toledo-support.com>.
- 3 Log in to the METTLER TOLEDO Balance Support Site (registration with the serial number of a METTLER TOLEDO instrument required).
- 4 Click **Customer Support**.
- 5 Click appropriate product folder.
- 6 Click **USB Driver**.


7 Click **USBDriverInstaller.exe**.

Installing procedure

- 1 Click **Save** to download to your specified location.
- 2 Right-click on the downloaded install program: **USBDriverInstaller.exe** and select **Run as Administrator** from the menu.
- 3 If a safety warning appears, allow Windows to install.



4 Click **Next** and follow the installer's instructions.




Installing Instrument

- 1 Switch the balance **off**.
- 2 Connect the balance to the preferred USB Port on the PC.
- 3 Switch the balance **on**.
- 4 Follow the instructions of the Wizard and install the software automatically (recommended)

Note: The wizard appears again for each USB port, either on your PC or if another balance is connected.

Warning: Do not click **Cancel** as for the connected USB port, it might not be possible anymore to perform the installation process.



14 Firmware (Software) Updates

METTLER TOLEDO is continuously improving its balance firmware (software) for the benefit of customers, so that the customer can benefit quickly and easily from further developments, METTLER TOLEDO makes the latest firmware versions available on the Internet. The firmware made available on the Internet has been developed and tested by Mettler-Toledo AG using processes that meet the guidelines of ISO 9001. Mettler-Toledo AG does not, however, accept liability for consequences that might arise from using the firmware.

14.1 Operating Principle

You will find all the relevant information and updates for your balance on the METTLER TOLEDO website at the following address:

www.mettler-toledo-support.com

A program known as the "**e-Loader II**" is loaded onto your computer together with the firmware update. You can use this program to download the firmware to the balance. The "e-Loader II" can also save the settings in your balance before the new firmware is downloaded to it. You can reload the saved settings into the balance manually or automatically after the software is downloaded.

If the selected update includes an application that is not described in these instructions (or that has been updated in the meantime) you can download the corresponding instructions in Adobe Acrobat® PDF format.

Note

New applications might not be visible unless the type data are updated by a service technician.

Requirements

The minimum requirements for obtaining applications from the Internet and downloading them into your balance are as follows:

- PC with one of the following Microsoft Windows® operating system:
 - Microsoft® Windows® XP Home or Professional with Service Pack 3 (32 bit)
 - Microsoft® Windows Vista® Home Premium, Business, Ultimate, or Enterprise with Service Pack 2 (32 bit and 64 bit)
 - Microsoft® Windows 7 with Service Pack 1 Home Premium, Professional, Ultimate, or Enterprise (32 bit and 64 bit)
- Administrator rights for installing software.
- PC to balance connection cable (e.g. No. 11101051 see chapter accessories)

14.2 Update Procedure

Installing the "e-Loader II" software from the Internet onto the PC.

- 1 Connect to the Internet.
- 2 Go to the site <http://www.mettler-toledo-support.com>.
- 3 Log in to the **METTLER TOLEDO Balance Support Site** (registration with the serial number of a METTLER TOLEDO instrument required).
- 4 Click **Customer Support**.
- 5 Click appropriate product folder.
- 6 Click the firmware version (e-Loader II) you need and save it on your specified storage location.
- 7 Right-click on the **firmware SNxxx.exe** and select **Run as Administrator** from the menu.
- 8 Follow the installer's instructions.



Loading the new firmware into the balance.

- 1 Right-click on **METTLER TOLEDO e-Loader II** and select Run as Administrator from the menu.
- 2 Follow the instructions, which will take you step-by-step through the installation.

15 Error and Status Messages


15.1 Error Messages

Error messages in the display draw your attention to incorrect operation or that the balance could not execute a procedure properly.

Error Message	Cause	Rectification
NO STABILITY	No stability.	Ensure more stable ambient conditions. If not possible, check settings for environment.
WRONG ADJUSTMENT WEIGHT	Wrong adjustment weight on pan or none at all.	Place required adjustment weight in center of pan.
REFERENCE TOO SMALL	Reference for piece counting too small.	Increase reference weight.
EEPROM ERROR - PLEASE CONTACT CUSTOMER SERVICE	EEPROM (memory) error.	Please contact METTLER TOLEDO customer service.
WRONG CELL DATA - PLEASE CONTACT CUSTOMER SERVICE	Wrong cell data.	Please contact METTLER TOLEDO customer service.
NO STANDARD ADJUSTMENT - PLEASE CONTACT CUSTOMER SERVICE	No standard calibration.	Please contact METTLER TOLEDO customer service.
PROGRAM MEMORY DEFECT - PLEASE CONTACT CUSTOMER SERVICE	Program memory defect.	Please contact METTLER TOLEDO customer service.
TEMP SENSOR DEFECT - PLEASE CONTACT CUSTOMER SERVICE	Temperature sensor defect.	Please contact METTLER TOLEDO customer service.
WRONG LOAD CELL BRAND - PLEASE CONTACT CUSTOMER SERVICE	Wrong load cell brand.	Please contact METTLER TOLEDO customer service.
WRONG TYPE DATA SET - PLEASE CONTACT CUSTOMER SERVICE	Wrong type data set.	Please contact METTLER TOLEDO customer service.
BATTERY BACKUP LOST - CHECK DATE TIME SETTINGS	Backup battery is empty. This battery ensures that the date and time are not lost when the balance is disconnected from power.	Connect the balance to the power supply for charging the battery (e.g. during the night) or contact METTLER TOLEDO customer service.
	Overload - The weight on the pan exceeds the weighing capacity of the balance.	Reduce the weight on the weighing pan.
	Underload	Check that the weighing pan is positioned correctly.
INITIAL ZERO RANGE EXCEEDED	Wrong weighing pan or pan is not empty.	Mount correct weighing pan or unload weighing pan.
BELOW INITIAL ZERO RANGE	Wrong weighing pan or pan is missing.	Mount correct weighing pan.
MEM FULL	Memory full.	Clear the memory and start a new evaluation.
FACTOR OUT OF RANGE	Factor is outside the allowed range.	Select a new factor.
STEP OUT OF RANGE	Step is outside the allowed range.	Select a new step.
OUT OF RANGE	Sample weight is outside the allowed range.	Unload the pan and load a new sample weight.

15.2 Status Messages

Status messages are displayed by means of small icons. The status icons indicate the following:

Status Icon	Signification
	Service Reminder Your balance is due for servicing. Contact your dealer's customer service department as soon as possible to have a technician service your balance. (See menu topic " SERV.ICON ")

16 Cleaning and Service

Every now and then, clean the weighing pan, draft shield element, bottom plate, draft shield (depending on the model) and housing of your balance. Your balance is made from high-quality, durable materials and can therefore be cleaned using a damp cloth or with a standard cleaning agent.

To thoroughly clean the draft shield glass panels, remove the draft shield from the balance. When reinstalling the draft shield, ensure that it is in the correct position.


Please observe the following notes:



- The balance must be disconnected from the power supply
- Ensure that no liquid comes into contact with the balance or the AC adapter.
- Never open the balance or AC adapter – they contain no components, which can be cleaned, repaired or replaced by the user.



- On no account use cleaning agents which contain solvents or abrasive ingredients, as this can result in damage to the operation panel overlay.




Please contact your METTLER TOLEDO dealer for details of the available service options. Regular servicing by an authorized service engineer ensures constant accuracy for years to come and prolongs the service life of your balance.

16.1 Cleaning the Glass Draft Shield (0.1 mg and 1 mg Models)


1 Remove the following parts:

- 1 Remove weighing pan, draft shield element (0.1 mg models) and pan support.
- 2 Remove the bottom plate.
- 3 Unlock the draft shield, lift it off the balance and place it on a clean surface.




2

- 1 Push the **glass doors (A)** back.
- 2 Turn the two **lock covers (B)** on the front as far as they will go.




3

- 1 Tilt the **front glass (C)** forward.
- 2 Remove the front glass.




4

- Pull the **top glass door (D)** out from the front.




5

- Lift the **side glass doors (E)** at (F) and pull them out from the front.




6

- Push the lock button (G) to release the **rear glass**.



7

- Remove the rear glass (H).



9

8

- 1 Turn the **draft shield lock** to the "Service" position.
- 2 Remove the draft shield lock.

After cleaning reinstall all components in the reverse order. For balance mounting see chapter "Setting up the Balance – Installing the Components".

17 Interface Specification

17.1 RS232C Interface

Each balance is equipped with an RS232C Interface as standard for the attachment of a peripheral device (e.g. printer or computer).

Schematic	Item	Specification
	Interface type	Voltage interface according to EIA RS-232C/DIN66020 CCITT V24/V.28)
	Max. cable length	15 m
	Signal level	Outputs: +5 V ... +15 V (RL = 3–7 kΩ) –5 V ... –15 V (RL = 3–7 kΩ) Inputs: +3 V ... +25 V –3 V ... –25 V
	Connector	Sub-D, 9-pole, female
	Operating mode	Full duplex
	Transmission mode	Bit-serial, asynchronous
	Transmission code	ASCII
	Baud rates	600, 1200, 2400, 4800, 9600, 19200, 38400 (software selectable)
	Bits/parity	7-bit/none, 7-bit/even, 7-bit/odd, 8-bit/none (software selectable)
	Stop bits	1 stop bit
	Handshake	None, XON/XOFF, RTS/CTS (software selectable)
	End-of-line	<CR><LF>, <CR>, <LF> (software selectable)
	Power supply for 2nd display	+ 12 V, max 40 mA (software selectable, 2nd display mode only)

17.2 USB Device Interface

Each balance is equipped with an "USB Device" Interface as standard for the attachment of a peripheral device (e.g. computer).

Note: This interface is not suitable to communicate with a Printer.

Schematic	Item	Specification										
<table border="1" style="margin-top: 10px;"> <tr> <td>1</td> <td>VBUS (+5 VDC)</td> </tr> <tr> <td>2</td> <td>D- (Data -)</td> </tr> <tr> <td>3</td> <td>D+ (Data +)</td> </tr> <tr> <td>4</td> <td>GND (Ground)</td> </tr> <tr> <td>Shield</td> <td>Shield</td> </tr> </table>	1	VBUS (+5 VDC)	2	D- (Data -)	3	D+ (Data +)	4	GND (Ground)	Shield	Shield	Standard	In conformity with USB Specification Revision 1.1
	1	VBUS (+5 VDC)										
	2	D- (Data -)										
	3	D+ (Data +)										
	4	GND (Ground)										
Shield	Shield											
Speed	Full speed 12 Mbps (requires shielded cable)											
Function	CDC (Communication Device Class) serial port emulation											
Power usage	Suspended device: Max 10 mA											
Connector	Type B											

17.3 MT-SICS Interface Commands and Functions

Many of the instruments and balances used have to be capable of integration in a complex computer or data acquisition system.

To enable you to integrate balances in your system in a simple manner and utilize their capabilities to the full, most balance functions are also available as appropriate commands via the data interface.

All new METTLER TOLEDO balances launched on the market support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS). The commands available depending on the functionality of the balance.

For further information please refer to the Reference Manual MT-SICS downloadable from the Internet under

▶ www.mt.com/sics-newclassic

18 Technical Data

18.1 General Data

Power Supply

- S Platform:

AC/DC Adapter


Primary: 100V–240V, $\pm 10\%$, 50/60Hz, 0.3 A

Secondary: 12VDC, 0.84A (with electronic overload protection)

Power supply to the balance: 11–20VDC, 10W



Use only with a tested AC Adapter with SELV output current.

Ensure correct polarity 

- L Platform:

Power supply 100V–240V, $\pm 10\%$, 50/60Hz, 0.3 A

Power cable 2-core with country-specific plug

Protection and Standards

- Overvoltage categorie:
- Degree of pollution:
- Degree of protection:
- Standards for safety and EMC:
- Range of application:

II, III

2

Protected against dust and water

See Declaration of Conformity

For use only in enclosed interior rooms

Environmental conditions

- Height above mean sea level:
- Ambient temperature range::
- Relative air humidity::
- Warm-up time

up to 4000 m

10 to 30 °C (S platform)

5 to 40 °C (L platform)

10% to 80 % up to 31 °C, linearly decreasing to 50 % at 40 °C, non-condensing

After connecting the balance to the power supply at least

- 30 minutes on balances with a readability of 0.001 g to 0.1 g.
- 60 minutes on balances with a readability of 0.1 mg and better.

Materials

- Housing:
- Weighing pan:
- Draft shield element:
- Draft shield:
- In-use-cover:

Die-cast aluminum, lacquered

Stainless steel X2CrNiMo 17-12-3 (1.4404)

245 x 351 mm: Stainless steel X5CrNiMo 18-10 (1.4301)

with 0.1 mg models: Stainless steel X2CrNiMo 17-12-2 (1.4404)

with 10 mg models: Plastic (PBT)

Plastic (PBT), glass

Plastic (PET)

18.2 Model-Specific Data

18.2.1 Balances with Readability of 0.1 mg, S Platform with Draft Shield

Technical Data

		PH204S
Limit values		
Maximum capacity		220 g
Readability		0.1 mg
Repeatability (sd)		0.1 mg
Linearity deviation		0.2 mg
Sensitivity temperature drift (10...30 °C)		1.5 ppm/°C
Typical values		
Repeatability (at nominal load)		0.08 mg
Linearity deviation		0.06 mg
Minimum sample weight (U=1%, k=2)		0.016 g
Minimum sample weight OIML		0.01 g
Settling time		2 s
Adjustment		Int. Cal, FACT
Interfaces		1 RS232, 1 USB
Balance dimensions W x D x H)		204x347x348
Weighing pan dimensions		Ø 90 mm
Usable heigh of draft shield		236 mm
Weight of balance		6.5 kg
Weights for routine testing		
OIML CarePac		#11123001
	Weights	200 g F2, 10 g F1
ASTM CarePac		#11123101
	Weights	200 g 1, 10 g 1

sd = Standard deviation

18.2.2 Balances with Readability of 1 mg, S Platform with Draft Shield

Technical Data

	PH303S	PH403S	PH603S
Limit values			
Maximum capacity	320 g	420 g	620 g
Readability	0.001 g	0.001 g	0.001 g
Repeatability (sd)	0.001 g	0.001 g	0.001 g
Linearity deviation	0.002 g	0.002 g	0.002 g
Sensitivity temperature drift (10...30 °C)	3 ppm/°C	3 ppm/°C	3 ppm/°C
Typical values			
Repeatability (at nominal load)	0.7 mg	0.7 mg	0.7 mg
Linearity deviation	0.6 mg	0.6 mg	0.6 mg
Minimum sample weight (U=1%, k=2)	0.14 g	0.14 g	0.14 g
Minimum sample weight OIML	0.02 g	0.02 g	0.02 g
Settling time	1.5 s	1.5 s	1.5 s
Adjustment	Int. Cal, FACT	Int. Cal, FACT	Int. Cal, FACT
Interfaces	1 RS232, USB	1 RS232, USB	1 RS232, USB
Balance dimensions W x D x H)	204x347x283 mm	204x347x283 mm	204x347x283 mm
Weighing pan dimensions	127x127 mm	127x127 mm	127x127 mm
Usable height of draft shield	168 mm	168 mm	168 mm
Weight of balance	6.2 kg	6.2 kg	6.2 kg
Weights for routine testing			
OIML CarePac	#11123001	#11123001	#11123007
Weights	200 g F2, 10 g F1	200 g F2, 10 g F1	500 g F2, 20 g F1
ASTM CarePac	#11123101	#11123101	#11123107
Weights	200 g 1, 10 g 1	200 g 1, 10 g 1	500 g 1, 20 g 1

sd = Standard deviation

Model	PH1003S
Limit values	
Maximum capacity	1020 g
Readability	0.001 g
Repeatability (sd)	0.001 g
Linearity deviation	0.002 g
Sensitivity temperature drift (10...30 °C)	3 ppm/°C
Typical values	
Repeatability (at nominal load)	0.7 mg
Linearity deviation	0.6 mg
Minimum sample weight (U=1%, k=2)	0.14 g
Minimum sample weight OIML	0.02 g
Settling time	1.5 s
Adjustment	Int. Cal, FACT
Interfaces	1 RS232, USB
Balance dimensions W x D x H)	204x347x283 mm
Weighing pan dimensions	127x127 mm
Usable heigh of draft shield	168 mm
Weight of balance	6.9 kg
Weights for routine testing	
OIML CarePac	#11123008
Weights	1000 g F2, 50 g F2
ASTM CarePac	#11123108
Weights	1000 g 1, 50 g 1

sd = Standard deviation

18.2.3 Balances with Readability of 0.01 g, S Platform

Technical Data

	PH3002S	PH3002SDR	PH4002S
Limit values			
Maximum capacity	3200 g	3200 g	4200 g
Maximum capacity, fine range	–	620 g	–
Readability	0.01 g	0.1 g	0.01 g
Readability, fine range	–	0.01 g	–
Repeatability (sd)	0.01 g	0.06 g	0.01 g
Repeatability (sd), fine range	–	0.01 g	–
Linearity deviation	0.02 g	0.2 g	0.02 g
Linearity deviation, fine range	–	0.02 g	–
Sensitivity temperature drift (10...30 °C)	3 ppm/°C	3 ppm/°C	3 ppm/°C
Typical values			
Repeatability (at nominal load)	0.007 g	0.007 g	0.007 g
Linearity deviation	0.006 g	0.006 g	0.006 g
Minimum sample weight (U=1%, k=2)	1.4 g	1.4 g	1.4 g
Minimum sample weight OIML	0.5 g	0.5 g	0.5 g
Settling time	1 s	1 s	1 s
Adjustment	Int. Cal, FACT	Int. Cal, FACT	Int. Cal, FACT
Interfaces	1 RS232, 1 USB	1 RS232, 1 USB	1 RS232, 1 USB
Balance dimensions (W x D x H)	194x347x99 mm	194x347x99 mm	194x347x99 mm
Weighing pan dimensions	170x200 mm	170x200 mm	170x200 mm
Weight of balance	5.5 kg	5.5 kg	5.5 kg
Weights for routine testing			
OIML CarePac	#11123009	#11123009	#11123010
Weights	2000 g F2, 100 g F2	2000 g F2, 100 g F2	2000 g F2, 200 g F2
ASTM CarePac	#11123109	#11123109	#11123110
Weights	2000 g 1, 100 g 1	2000 g 1, 100 g 1	2000 g 4, 200 g 4

sd = Standard deviation

18.2.4 Balances with Readability of 0.1 g, L Platform

Technical Data


	PH16001S	PH32001S
Limit values		
Maximum capacity	16200 g	32200 g
Readability	0.1 g	0.1 g
Repeatability (sd)	0.1 g	0.1 g
Linearity deviation	0.2 g	0.3 g
Sensitivity temperature drift (10...30 °C)	5 ppm/°C	5 ppm/°C
Typical values		
Repeatability (at nominal load)	0.07 g	0.07 g
Linearity deviation	0.06 g	0.06 g
Minimum sample weight (U=1%, k=2)	14 g	14 g
Minimum sample weight OIML	5 g	5 g
Settling time	1 s	1 s
Adjustment	Int. Cal, FACT	Int. Cal, FACT
Interfaces	1 RS232, 1 USB	1 RS232, 1 USB
Balance dimensions (W x D x H)	363x346x118 mm	363x346x118 mm
Weighing pan dimensions	351x245 mm	351x245 mm
Weight of balance	10.7 kg	10.7 kg
Weights for routine testing		
OIML Weights	10000 g F2, 500 g F2	20000 g F2, 1000 g F2
ASTM Weights	10000 g 4, 500 g 4	20000 g 4, 1000 g 4

sd = Standard deviation

18.3 Dimensions


18.3.1 Balances with Readability of 0.1 mg, S Platform with Draft Shield

Models:
PH204S




18.3.2 Balances with Readability of 1 mg, S Platform with Draft Shield

Models:
 PH303S
 PH403S
 PH603S
 PH1003S




18.3.3 Balances with Readability of 0.01 g, S Platform

Models:
PH4002S
PH6002S










18.3.4 Balances with Readability of 0.1 g, L Platform

Models:
PH16001S
PH32001S



19 Accessories and Spare Parts

19.1 Accessories

	Description	Part No.
Draft shields		
	Draft shield with sliding doors "mg" (usable height 168 mm)	12122405
	Draft shield with sliding doors "0.1 mg" (usable height 236 mm)	12122404
	Draft shield MS-DS-21 for models with readability of 0.1 g up to 0.01 g.	12121014
Printers		
	RS-P25 printer with RS232C connection to instrument	11124300
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	RS-P26 printer with RS232C connection to instrument (with date and time)	11124303
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	RS-P28 printer with RS232C connection to instrument (with date, time and applications)	11124304
	Paper roll, set of 5 pcs	00072456
	Paper roll, self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975
	P-56RUE Thermal Printer with RS232C, USB and Ethernet connections, simple printouts, Date and Time, Label printing (limited).	30094673
	Paper roll, white, set of 10 pcs	30094723
	Paper roll, white, self-adhesive, set of 10 pcs	30094724
	Paper roll, white, self-adhesive labels, set of 6 pcs	30094725



P-58RUE Thermal Printer with RS232C, USB and Ethernet connections, simple printouts, Date and Time, Label printing, Balance applications: Statistics, Formulation, Totaling,

30094674

Paper roll, white, set of 10 pcs

30094723


Paper roll, white, self-adhesive, set of 10 pcs

30094724

Paper roll, white, self-adhesive labels, set of 6 pcs


30094725

Cables for RS232C interface




RS9 – RS9 (m/f): connection cable for PC, length = 1 m

11101051



RS9 – RS25 (m/f): connection cable for PC, length = 2 m


11101052



RS232 - USB converter cable – Cable with converter to connect a balance (RS232) to a USB port

64088427

Cables for USB interface



USB (A –B) connection cable for connection to PC, length = 1 m

12130716

Cable replacement (wireless)




Bluetooth RS232 Serial Adapter ADP-BT-S for wireless connection between **printer** and Excellence balance* or between **balance** and PC*. Fits printers P-56 / P-58 and the following balance models (SW V2.20 or higher required): MS, MS-S/L, ML, PHS, JP, JS.

30086494

* Bluetooth interface required

- 1 Bluetooth RS232 Serial Adapter (slave)
- 1 MT-DB9 male to female connector
- 1 MT-DB9 male to male connector




Bluetooth RS232 Serial Adapter set ADP-BT-P for wireless connection between printer and balance. Fits printers P-56 / P-58 and the following balance models (SW V2.20 or higher required): MS, MS-S/L, ML, PHS, JP, JS.

30086495

- 2 Bluetooth RS232 Serial Adapter paired (slave/master)
- 1 MT-DB9 male to female connector
- 1 MT-DB9 male to male connector


Auxiliary displays



RS232 auxiliary display AD-RS-M7


12122381

Power supplies



AC/DC universal adapter (EU, USA, AU, UK) 100–240 VAC, 50/60 Hz, 0.3 A, 12 VDC 0.84 A


11120270



PowerPac-M-12V, for mains independent operation of balances, 12 VDC/1 A

12122363


Pan protections



Protective foils, 166x196 mm, set of 20 pcs, pan protection for weighing pans from 170x200 mm to 190x226 mm


30113800

Protective covers




Protective cover for S platform with draft shield

12121850



Protective cover for S platform without draft shield


12121851



Protective cover for L platform up to "1 g"

12121852


Anti-theft devices



Steel cable

11600361


Software



LabX direct balance (simple data transfer)

11120340


Weighing below the balance



Hook for Platform L

11132565


Transport cases



Transport case for S platform balances

11124245


Adjustment weights




OIML / ASTM Weights (with calibration certificate) see www.mt.com/weights

19.2 Spare Parts

Draft Shield

Drawing	Pos	Description	Part No.	
	5	Draft shield lock	12122013	
	6	Bottom plate	12122019	
	Draft Shield "168 mm"			
	1	Top glass with handle	12121884	
	2	Rear glass low	12122015	
	3	Side glass door left low with handle	12121881	
	4	Side glass door right low with handle	12121883	
	7	Front glass low	12122014	
	Draft Shield "236 mm"			
	1	Top glass with handle	12121884	
	2	Rear glass high	12122012	
	3	Side glass door left high with handle	12121880	
	4	Side glass door right high with handle	12121882	
	7	Front glass high	12122011	

Weighing Pans / Draft Shield Elements / Support

Drawing	Pos	Description	Part No.	
	For S platform			
	1	0.1 mg	Weighing pan Ø 90 mm	12122010
	2	0.1 mg	Pan support Ø 90 mm	11124249
	3	0.1 mg	Draft shield element	12122008
	4	1 mg	Weighing pan 127x127 mm	12122009
	5	1 mg	Pan support 127x127 mm for models up to 999 g	12122017
	5	1 mg	Pan support 127x127 mm for models from 1000 g	12122016
	6	10 mg	Draft shield element 170x200 mm	12122018
	7	10 mg	Weighing pan 170x200 mm	11124247
	8	10 mg	Pan support 170x200 mm	12121064
	7	0.1 g	Weighing pan 190x226 mm	11124248
	8	0.1 g	Pan support 190x226 mm	12121066
	9	from 10 mg	Pan support caps	11131029
	10	Leveling foot		11106323
	For L platform			
	11	Weighing pan 245x351 mm		12122020
	12	to 1 g	Pan support caps	00239104
12	from 2 g	Pan support caps	12122006	
13	Leveling foot		00230236	

Index

A

Accessories	86
Adjusting	19
Adjustment	21, 31, 31, 31
Advanced Menu	27
Ambient conditions	16
Application "Check Weighing"	49
Application "Diagnostics"	34, 57
Application "Formulation"	43
Application "Piece Counting"	46
Application "Routine Test"	54
Application "Statistics"	52
Application "Weighing"	23
Application icons	12
Assign Application	34, 34, 34
Auto print	37
Automatic adjustment	20
Automatic shutoff	33
Automatic zero setting	33
Autozero	33
Average (Statistics)	52

B

Backlight	33
Balance history	60
Balance information	62
Basic menu	27, 30
Baudrate	38
Beep	30, 30
Bit/Parity	38

C

Calibration	31, 31
Calibration history	61
Cancel	29
Change settings	28, 28
Char Set	39, 41
Check Weighing	49
Cleaning	71
Closing the menu	29
Control Limit	55
Conventions and symbols	7
Customer fine adjustment	21, 31

D

Data communication format	37, 40
Date	30
Date format	32
Delivery inspection	14
DeltaRange balances	25
Diagnose	34
Diagnostics	57
Diagnostics application	34
Dimensions	82
Display	33, 36
Display panel	12

Display test	58
Disposal	8
Dosing	31
Draft Shield	15, 90

E

End of Line	39, 41
Environment	31
Error messages	69
External weight	21

F

FACT	20, 31, 32
Fill up	44
Firmware update	67
Formulation	43
Fully automatic adjustment	20, 31
Function "Fill up"	44
Function PC-Direct	64

G

Good Weighing Practice	54
GWP	54, 55

H

Handshake	39
Header	36
Host	35, 40

I

Icons	12
Input principle	28
Installing the components	15
Interface	
MT-SICS	75
Interface menu	27, 35
Interface RS232C	35, 74
Interface USB device	39, 74
Internal weight	20
Interval	41
Introduction	7

K

Key assign	34, 34, 34
Key beep	30
Key functions	11
Key test	59

L

L platform overview	10
Language	34
Leveling the balance	17
Line feed	36
Location	16

M

Main Menu	29
Manual adjustment with external weight	21

Manual adjustment with internal weight	20	SOP	55
Menu	27, 29	Spare parts	90
Menu Advanced	27	Stability beep	30
Menu Basic	27, 30	Standard Deviation (Statistics)	52
Menu Interface	27, 35	Startup	32
Menu operation	28	Statistics	52
Menu protection	29	Status icons	12
Menu topic	28, 28, 29	Status messages	70
Motor test	60	Stop bit	38
MT-SICS	75	Submenu	28
<hr/>		Switching	
N		On	19
Net	24	Switching the balance on and off	23
Numerical values	28	Switching weight units	25
<hr/>		Symbols and conventions	7
O		<hr/>	
Operating temperature	19	T	
Operation keys	11	Taring	24
Overview	9	Technical data dimensions	82
<hr/>		Technical data general	76
P		Technical data model-specific	77
PC-DIR	35	Time	30
PC-Direct	64	Time format	32
Performing a simple weighing	24	Topic	28, 28, 29
Piece Counting	46	Transmit data	26
Power supply	18	Transporting the balance	18
Print	26	<hr/>	
Printer	35	U	
Protect	29	Unit	30, 30
Protocol trigger	32	Unpacking	14
<hr/>		USB device interface	39, 65, 74
Q		USB-Driver installing	65
Quickstart	32	<hr/>	
<hr/>		W	
R		Warm-up time	19, 76
Recall	25, 32	Warning Limit	55
Repeatability test	57	Weighing below the balance	18
Reset	30	Weighing made simple	23
Routine Test	54	Weighing mode	31
RS232C interface	35, 74	Weighing-in aid	25
<hr/>		Weight unit	25, 30, 30
S		<hr/>	
S platform overview	9	Z	
Safety precautions	8	Zero print	37
Saving settings	29	Zero range	33
Select menu	28	Zero setting	24
Select menu topic	28	Zeroing	33
Selecting the location	16		
Service	35, 35		
Service date reset	35		
Service icon	35		
Service provider information	63		
Service reminder	35		
Setting up the balance	14		
Shutoff	23, 33		
Signature line	36		
Single	36		
Software update	67		

GWP® – Good Weighing Practice™

The global weighing guideline GWP® reduces risks associated with your weighing processes and helps to

- choose the appropriate balance
- reduce costs by optimizing testing procedures
- comply with the most common regulatory requirements

► www.mt.com/GWP

www.mt.com/pharmacy

For more information

Mettler-Toledo AG, Laboratory Weighing

CH-8606 Greifensee, Switzerland

Tel. +41 (0)44 944 22 11

Fax +41 (0)44 944 30 60

www.mt.com

Subject to technical changes.

© Mettler-Toledo AG 03/2014

11781377B en

