Instruction manual
O₂ Transmitter 4100 e/2(X)H



Warranty

Defects occurring within 1 year from delivery date shall be remedied free of charge at our plant (carriage and insurance paid by sender).

Subject to change without notice.

Return of products under waranty

Please contact METTLER TOLEDO's Customer Service Dept. before returning a defective device. Ship the <u>cleaned</u> device to the address you have been given. If the device has been in contact with process fluids, it must be decontaminated/disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.



Disposal (Directive 2002/96/EC of January 27, 2003) Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".







Mettler-Toledo GmbH, Process Analytics, Industrie Nord, CH-8902 Urdorf, Tel. +41 (01) 736 22 11 Fax +41 (01) 736 26 36 Subject to technical changes. Mettler-Toledo GmbH, 04/04. Printed in Germany.

METTLER TOLEDO

Contents

Sarety Information	
Intended use	
Trademarks	7
Certificates	8
EC Declaration of Conformity	8
Support of FDA 21 CFR Part 11	9
EC-Type-Examination Certificate	10
Overview of O ₂ Transmitter 4100 e/2(X)H	.13
Assembly	.14
Package contents	
Mounting plan	
Pipe mounting, panel mounting	
Installation and connection	
Information on installation	
Division 2 wiring	
Terminal assignment	
Wiring examples	
User interface, display	.26
Operation: Keypad	.28
Safety functions	
Sensocheck, Sensoface sensor monitoring	29
GainCheck device self test	29
Automatic device self-test	
Hold mode	30
Outputs	.31
(Current output / loop current, HART communication, Alarm)	.31
Passcodes (factory setting)	.33
Configuration	.34
Menu structure of configuration	
Overview of configuration steps	
Individual settings (for copy)	

Contents

	Current output Temperature compensation Correction Calibration mode Alarm settings	.50 .52 .54
Pa	asscodes according to FDA 21 CFR Part 11	.58
	Calibration Calibration to percent saturation (SAT) Calibration to concentration (Conc) Calibration to volume concentration Product calibration Zero calibration Temp probe adjustment	.62 .64 .66 .68 .70 .72
	easurement	
Di	agnostics functions Display of output currents Display of calibration data (Cal Info) Sensor monitor Display of last error message Specify output current	.73 .73 .73 .73
	eaning	
	perating states	
Er	ror messages (error codes)	
Se	ensoface	.80
-	Ppendix Product line and accessories Specifications FM Control Drawing Explosion protection CSA Control Drawing	.83 .84 .90 .92 .94
In	dex	.97

Safety information

Be sure to read and observe the following instructions!

The device has been designed using state of the art technology and it complies with the applicable safety regulations. When operating the device, certain conditions may nevertheless lead to danger for the operator or damage to the device.

Caution!

Commissioning may only be carried out by trained experts. Whenever it is likely that protection has been impaired, the device shall be made inoperative and secured against unintended operation.

The protection is likely to be impaired if, for example:

- the device shows visible damage
- the device fails to perform the intended measurements
- after prolonged storage at temperatures above 70 °C
- after severe transport stresses

Before recommissioning the device, a professional routine test in accordance with EN 61010-1 must be performed. This test should be carried out by the manufacturer.

Caution!

Before commissioning it must be proved that the device may be connected with other equipment.

Safety precautions for installation

- The stipulations of EN 60079-10 / EN 60079-14 must be observed during commissioning.
- The O₂ Transmitter 4100 e/2H is approved for measurements in FM Class I Div 2
- The O₂ Transmitter 4100 e/2XH is approved for operation in the following locations: ATEX, FM Zone 1 with measurement in Zone 0, and FM Class I Div 1.

Connection to supply units

- O₂ Transmitter 4100 e/2H: Before connecting this device to a supply unit, make sure that its output voltage cannot exceed 30 V DC.
 - Do not use alternating current or mains power supply!
- O₂ Transmitter 4100 e/2XH: This device may only be connected to an explosion-proof power supply unit (for input ratings refer to annex of EC-Type-Examination Certificate).
 - Before commissioning it must be made sure that the connections to other equipment such as power supply unit and cables are intrinsically safe.

Terminals:

suitable for single wires / flexible leads up to 2.5 mm² (AWG 14)

Note for cleaning in a hazardous location

In hazardous locations the Transmitter may only be cleaned with a damp cloth to prevent electrostatic discharge.

Intended use

The $\rm O_2$ Transmitter 4100 e/2(X)H is used for measuring dissolved or gaseous oxygen and temperature in biotechnology, chemical and pharmaceutical industry, as well as in the field of industry, environment, food processing, and sewage treatment.

The rugged molded enclosure can be wall or pipe mounted or fixed into a control panel.

The protective hood provides additional protection against direct weather exposure and mechanical damage.

The Transmitter has been designed for application with amperometric METTLER TOLEDO sensors of the InPro 6800 and InPro 6900 series.

- The O₂ Transmitter 4100 e/2H is approved for measurements in FM Class I Div 2.
- The O₂ Transmitter 4100 e/2XH is approved for operation in the following locations: ATEX, FM Zone 1 with measurement in Zone 0, and FM Class I Div 1.

Trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

GainCheck

Sensocheck

Sensoface

InPro® is a registered trademark of Mettler-Toledo.

HART® is a registered trademark of the HART Communications Foundation (HCF).

EC Declaration of Conformity

Mettler-Toledo GmbH

Process Analytics

Adresse Im Hackacker 15 (Industrie Nord), CH-8902 Urdorf, Schweiz Briefadresse Postfach, CH-8902 Urdorf Telefon 01-736 22 11 Telefox 01-736 26 36

Bank Credit Suisse First Boston, Zürich (Acc. 0835-370501-21-90

Declaration of conformity Konformitätserklärung Déclaration de conformité



We/Wir/Nous Mettler-Toledo GmbH, Process Analytics

> Im Hackacker 15 8902 Urdorf Switzerland

declare under our sole responsibility that the product. erklären in alleiniger Verantwortung, dass dieses Produkt, déclarons sous notre seule responsabilité que le produit.

Description Beschreibung/Description 0,-4100/2(X)H

> to which this declaration relates is in conformity with the following standard(s) or other normative document(s).

auf welches sich diese Erklärung bezieht, mit der/den folgenden Norm(en) oder Richtlinie(n) übereinstimmt.

auguel se réfère cette déclaration est conforme à la (aux) norme(s) ou au(x) document(s) normatif(s).

Explosion Protection / Protection contre les explosions 94/9/EG

Explosionsschutzrichtlinie

EMC Directive/EMV-Richtlinie Directive concernant to CEM 89/336/FWG SR 734.5, VEMV

Low-voltage directve/Niederspannungs-Richtlinie/

Directive basse tension 73/23/EWG SR 734.26. NEV

Norm/Standard/Standard EN 50 014:

EN 50 020: EN 61326 / VDE 0843 Teil 20: 2002-03 EN 61010 Teil / VDE 0411 Teil 1: 2002-08

Place and Date of issue Ausstellungsort / - Datum Lieu et date d'émission

Urdorf, March 15, 2004

Mettler-Toledo GmbH. Process Analytics

General Manager PO Urdorf

Artikel Nr.: 52960316 KE Dateiname: 52980316KE-4100e-2(X)H-Inte Sitz der Gesellschaft - Mettler-Toledo GmbH, Im Langacher, CH-8606 Greifensee



Support of FDA 21 CFR Part 11

In their directive "Title 21 Code of Federal Regulations, 21 CFR Part 11, Electronic Records; Electronic Signatures" the US American health agency FDA (Food and Drug Administration) regulates the production and processing of electronic documents for pharmaceutical development and production. This results in requirements for measuring devices used for corresponding applications. The following features ensure that the measuring devices of the O₂ Transmitter 4100 e Series meet the demands of FDA 21 CFR Part 11:

Electronic Signature

Access to the device functions is regulated and limited by individually adjustable codes – "Passcodes" (for Passcode Editor see Page 58, overview of factory settings on back of manual). This prevents unauthorized modification of device settings or manipulation of the measurement results. Appropriate use of these passcodes makes them suitable as electronic signature.

Audit Trail

Every (manual) change of device settings can be automatically documented. For that purpose, each change is marked by a "Configuration Change Flag", which can be interrogated and documented via HART communication. Then the changed device settings/parameters can also be retrieved and documented via HART cmmunication.

EC-Type-Examination Certificate



Translation

(1) EC-TYPE EXAMINATION CERTIFICATE

- Equipment and protective systems intended for use in potentially explosive atmospheres - Directive 94/9/EC
- (3) EC-Type Examination Certificate Number



TÜV 04 ATEX 2431

(4) Equipment: Transducer type O₂ 4100/2XH (5) Manufacturer: Mettler Toledo GmbH

Process Analytics

(6) Address: CH-8902 Urdorf, Im Hackacker 15

- (7) This equipment or protective system and any acceptable variation thereto are specified in the schedule to this certificate and the documents therein referred to.
- (8) The TÜV NORD CERT GmbH & Co. KG, TÜV CERT-Certification Body, notified body number N° 0032 in accordance with Article 9 of the Council Directive of the EC of March 23, 1994 (94/9/EC), certifies that this equipment or protective system has been found to comply with the Essential Health and Safety Requirements relating to the design and construction of equipment and protective systems intended for use in potentially explosive atmospheres given in Annex II to the Directive.

The examination and test results are recorded in the confidential report N° 04 YEX 551230.

(9) Compliance with the Essential Health and Safety Requirements has been assured by compliance with:

EN 50014:1997 EN 50020:2002 EN 50284:1999

- (10) If the sign "X" is placed after the certificate number, it indicates that the equipment or protective system is subject to special conditions for safe use specified in the schedule to this certificate
- (11) This EC-type examination certificate relates only to the design, examination and tests of the specified equipment in accordance to the Directive 94/9/EC. Further requirements of the Directive apply to the manufacturing process and supply of this equipment. These are not covered by this certificate.
- (12) The marking of the equipment or protective system must include the following:



TÜV NORD CERT GmbH & Co. KG TÜV CERT-Certification Body Am TÜV 1 D-30519 Hannover

Tel.: 0511 986-1470 Fax: 0511 986-2555

Head of the Certification Body



Hanover, 2004-02-09

This certificate may only be reproduced without any change, schedule included. Excerpts or changes shall be allowed by the TÜV NORD CERT GmbH & Co. KG

page 1/3



SCHEDULE (13)

(14) EC-TYPE EXAMINATION CERTIFICATE N° TÜV 04 ATEX 2431

(15) Description of equipment

The transducer type O₂ 4100/2XH is preferable intended for the detection and processing of electrochemical parameters of fluids. For this purpose the transducer is equipped with one input for the oxygen measurement and one for temperature measurement.

The maximum permissible ambient temperature is 55°C.

Electrical data

Loop measuring circuit in type of protection Intrinsic Safety EEx ib IIC (KL 10, 11)

only for the connection of certified intrinsically safe circuit with

the following maximum values:

U = 30 V I = 100 mA P.= 0.8 W

effective internal capacitance C = 20 nF effective internal inductance L₁ = 0,2 mH

Oxygen measuring circuit in type of protection Intrinsic Safety EEx ia IIC Maximum values:

(KL 1, 2, 4, 5, 6) U = 10 V

I_o= 10 mA P_o= 13 mW

Characteristic line: linear effective internal capacitance $C_1 = 15 \text{ nF}$

The internal inductance is negligibly small. max. permissible outer capacitance C o = 1,5 μF max. permissible outer inductance Lo= 1,0 mH

Temperature measuring circuit in type of protection Intrinsic Safety EEx ia IIC (KL 7, 8) Maximum values:

U .= 5 V $I_0 = 1 \text{ mA}$ P = 2 mW

Characteristic line: linear

effective internal capacitance C = 120 nF

The internal inductance is negligibly small. max. permissible outer capacitance C a = 1,38 μF max, permissible outer inductance L = 1.0 mH

Oxygene/temperature measuring circuit

(Circuits interconnected) (KL 1, 2, 4, 5, 6, 7, 8)

in type of protection Intrinsic Safety EEx ia IIC

Maximum values: U 0= 10 V

1 = 11 mA P = 14 mW



Schedule EC-Type Examination Certificate Nº TÜV 04 ATEX 2431

characteristic line: linear effective internal capacitance C_1 = 135 nF The internal inductance is negligibly small. max. permissible outer capacitance C_0 = 1,38 μ F max. permissible outer inductance L_0 = 1,0 mH

or

for the connection of the oxygen sensors type InPro $6XXX/^*/^*/^*/^*$ according to SNCH 01 ATEX 3277 X

PA (KL 9) for the connection to the equipotential bonding system

The loop measuring circuit is safely galvanically separated from all other measuring circuits up to a voltage of 60 V. The oxygen measuring circuit and the temperature measuring circuit are galvanically connected.


- (16) Test documents are listed in the test report No.: 04 YEX 551230.
- (17) Special conditions for safe use

none

(18) Essential Health and Safety Requirements

no additional ones

Overview of the O₂ Transmitter 4100 e/2(X)H




Assembly

Package contents

Check the shipment for transport damage and completeness. The package should contain:

- Front unit
- Lower case
- Bag containing small parts
- Instruction manual
- Specific test report



- 1 Jumper (2 piece)
- Washer (1 piece), for conduit mounting: place washer between enclosure and nut
- 3 Cable ties (3 pieces)
- 4 Hinge pin (1 piece), insertable from either side
- 5 Enclosure screws (4 pieces)

- 6 Sealing inserts (1 piece)
- 7 Rubber reducer (1 piece)
- 8 Cable glands (3 pieces)
- 9 Filler plugs (3 pieces)
- 10 Hexagon nuts (5 pieces)
- 11 Sealing plugs (2 pieces): for sealing in case of wall mounting

Fig.: Assembling the enclosure

Mounting plan





Fig.: Mounting plan

Pipe mounting, panel mounting



- Protective hood (if required)
- 2 Hose clamps with worm gear drive to DIN 3017 (2 pieces)
- 3 Pipe-mount plate (1 piece)
- 4 For vertical or horizontal posts or pipes
- 5 Self-tapping screws (4 pieces) (All dimensions in mm.)

Fig.: Pipe-mount kit




Fig.: Protective hood for wall and pipe mounting

METTLER TOLEDO




Fig.: Panel-mount kit

- 1 Screws (4 pieces)
- 2 Gasket (1 piece)
- 3 Panel
- 4 Span pieces (4 pieces)
- 5 Threaded sleeves (4 pieces)

Panel cutout 138 x 138 mm (DIN 43700)

Installation and connection

Information on installation

- Installation may only be carried out by trained experts in accordance with this instruction manual and as per applicable local and national codes.
- Be sure to observe the technical specifications and input ratings.
- Be sure not to notch the conductor when stripping the insulation
- When commissioning, a complete configuration must be carried out by the system administrator.

Connection to supply units

- O₂ Transmitter 4100 e/2H: Before connecting this device to a supply unit, make sure that its output voltage cannot exceed 30 V DC. Do not use alternating current or mains power supply!
- O₂ Transmitter 4100 e/2XH: This device may only be connected to an explosion-proof power supply unit (for input ratings refer to annex of EC-Type-Examination Certificate).
 Terminals:
 - suitable for single wires / flexible leads up to 2.5 mm² (AWG 14)

Warning!

Additional safety precautions have to be taken for applications in hazardous locations to CSA! (See Pg 93 et seq.)

Division 2 wiring



The connections to the Transmitter must be installed in accordance with the National Electric Code (ANSI-NFPA 70) Division 2 hazardous (classified) location non-incendive wiring techniques.

FM Control Drawing: Refer to page 90.

Terminal assignments




Fig.: Terminal assignments O₂ Transmitter 4100 e/2H





Fig.: Terminal assignments O₂ Transmitter 4100 e/2XH




- 5 6 (7)
- **1** Recommended stripping lengths for multi-core cables
- 2 Pulling out the terminals using a screwdriver (also see 6)
- 3 Cable laying in the device
- 4 Connecting lines for current supply
- **5** Cover for sensor and temperature probe terminals
- **6** Area for placing the screwdriver to pull out the terminals
- 7 Connection of current supply and HART handheld terminal

Fig.: Information on installation, rear side of device

Wiring example 1


Sensors with connection via T82



Connection	Terminal	InPro6800 sensor
cathode	1	white (coax core)
	2	not connected
guard anode shield RTD RTD	4 5 6 7 8	not connected red (coax shield) yellow/green (external shield) black black

Wiring example 2


Sensors with connection via VP cable




		InPro6800/6900 sensor
Connection	Terminal	VP cable
cathode	1 2	transparent (coax core) not connected
guard anode shield RTD RTD	4 5 6 7 8	gray (InPro6900) red (coax shield) yellow/green (external shield) green white

Wiring example 2

Connecting sensor and VP cable




VP cable connector assignment



User interface and display


User interface



- 1 Display
- 2 Mode indicators (no keys), from left to right:
 - Measuring mode
 - Calibration mode
 - Alarm
 - Wash contact (Model O₂ 4100 e only)
 - Configuration mode

- 3 Keypad
- 4 Coding
- 5 Rating plate6 Model designation
- 7 Alarm LED

Display



- 1 Passcode entry
- 2 Display of measured variable*
- 3 Temperature
- 4 Current output
- 5 Limit values
- 6 Alarm
- 7 Sensocheck
- 8 Calibration
- 9 Interval/response time
- 10 Wash contact
- 11 Measurement symbols
- 12 Proceed with enter
- 13 Bar for identifying the device status, above mode indicators from left to right:
 - Measuring mode
 - Calibration mode
 - Alarm
 - Wash contact*
 - Configuration mode

- 14 Lower display
- 15 Manual temp indicator
- 16 Hold mode active
- 17 Waiting time running
- 18 Sensor data
- 19 Main display 20 Sensoface
- * Not in use

Operation: Keypad

cal	Start, end calibration
conf	Start, end configuration
•	Select digit position (selected position flashes)
A	Edit digit
enter	 Calibration: Continue in program sequence Configuration: Confirm entries, next configuration step Measuring mode: Display output current

cal	Cal Info, display of zero current and slope
conf	Error Info: Display of last error message
+ 4	Start GainCheck device self-test

Safety functions

Sensocheck, Sensoface sensor monitoring

Sensocheck continuously monitors the sensor and lines. Sensocheck can be switched off (Configuration, Pg 57). With sensor type B Sensocheck must be switched off



Sensoface provides information on the sensor condition



The slope, response time and Sensocheck during calibration are evaluated. The three Sensoface indicators provide the user with information about wear and required maintenance of the sensor.

GainCheck device self-test

A display test is carried out, the software version is displayed and the memory and measured value transfer are checked.

Start GainCheck device self-test: + +



Automatic device self-test

The automatic device self-test checks the memory and measured-value transfer. It runs automatically in the background at fixed intervals

Safety functions

Hold mode

Display: 🖺

The Hold mode is a safety state during configuration and calibration. Output current is frozen (LAST) or set to a fixed value (FIX).

If the calibration or configuration mode is exited, the Transmitter remains in the Hold mode for safety reasons. This prevents undesirable reactions of the connected peripherals due to incorrect configuration or calibration. The measured value and "HOLD" are displayed alternately. The Transmitter only returns to measuring mode after **enter** is pressed and a waiting time of 20 s has passed.

Configuration mode is also exited automatically 20 minutes after the last keystroke (timeout). The Transmitter returns to measuring mode.

Timeout is not active during calibration.

Behavior of output signal:

LAST: The loop current is frozen at its last value.

The process should not change decisively during configuration/calibration. Changes are not noticed with this setting!

FIX: The loop current is set to a value that is noticeably different from the process value in order to signal the control system that the Transmitter is being worked at.

Configuration: Pg 49.

Current outputs

Loop current

The loop current is controlled by the process variable selected in the configuration.

The current start and end can be set to represent any desired value.

To check connected peripherals (e.g. limit switches, controllers), the loop current can be manually specified (see Pg. 74).

HART communication

The O_2 Transmitter 4100 e/2(X)H can be remote-controlled via HART communication. It can be configured using a handheld terminal or from the control room. Measured values, messages and device identification can be downloaded at any time. This allows easy integration also in fully automatic process cycles.

A list of the HART commands can be found in the "O₂ Transmitter 4100 e/2(X)H Transmitter-Specific Command Specification": www.mtpro.com/transmitters.

Alarm

The alarm delay is configurable.

Error messages can also be signaled by a 22 mA loop current (see Configuration, Pg 49).

The alarm LED on the front panel can be configured as follows:

HOLD off: Alarm: LED flashing

HOLD on: Alarm: LED on. HOLD: LED flashing.

Passcodes (factory setting)

The passcodes allow fast access to the functions

Calibration

Key+p	asscode	Description	Page
cal	0000	Cal Info	73
cal	1001	Zero calibration	70
cal	1100	Slope calibration: for saturation	62
		for concentration	64
		Volume concentration (gas)	66
cal	1105	Product calibration Adjusting the zero / slope (product)	68
cal	1015	Temp probe adjustment	72

Configuration

Key+passcode		Description	
conf	0000	Error Info Display last error and erase	73
conf	1200	Configuration	34
conf	2222	Sensor monitor Display sensor current / temperature	73
conf	5555	Current source Specify output current	74

Passcode editor

Key+passcode		Description	Page
conf	1989	Administrator passcode Changing the passcodes	58

Configuration

In the Configuration mode you set the device parameters.

Activate	conf	Activate with conf
		Enter passcode "1200"*. Edit parameter with ▶ and ♠, confirm/continue with enter. (End with conf and then enter.)
Hold	HOLD A of Hold	During configuration the Transmitter remains in the Hold mode for reasons of safety. The loop current is frozen (at its last value or at a preset fixed value, depending on the configuration), Sensoface is off, mode indicator "Configuration" is on. Red LED flashes when "HOLD ON" has been set.
Input errors	Err _	The configuration parameters are checked during the input. In the case of an incorrect input "Err" is displayed for approx. 3 sec. The incorrect parameters cannot be stored. Input must be repeated.
End	conf p enter	End with conf . The measured value and Hold are displayed alternately, "enter" flashes. End Hold mode with enter . The display shows the measured value. The output current remains frozen for another 20 s (HOLD icon on, "hourglass" flashes).

^{*} Factory setting, for passcode editing, see Pg 58

Menu structure of configuration

The configuration steps are assigned to different menu groups:

- Current output (code: o1.)
- Temperature compensation (code: tc.)
- Correction (code: Co.)
- Calibration mode (code: CA.)
- Alarm settings (code: AL.)

With the arrow keys you can jump between the individual menu groups. Each menu group contains menu items for setting the parameters.




Example:

- "o1." is displayed with all menu items of the
- "Current output" menu group.

Pressing **enter** opens a menu item. The values are edited using the arrow keys. Pressing **enter** confirms/stores the settings.

Return to measurement: Press **conf**. Press **enter** to confirm safety prompt. After 20 sec the Transmitter will be in measuring mode again.



Overview of configuration steps

Code	Menu	Selection / Default	
out1	Current output	(Factory setting bold print)	
o1.FCT	Select process medium: dissolved oxygen (DO) / gas	DO / GAS	
o1.SnSR	Select sensor type	Standard (Type A) / Traces (Type B)	
o1.UnIT	Select: Saturation (SAt) / Concentration (Conc) Only with GAS selected: Volume concentration	% / μg/l, mg/l, ppb, ppm % / ppm	
o1.4mA	Enter current start	xxxx (0000 %)	
o1.20mA	Enter current end	xxxx (0500 %)	
o1.FtME	Time constant of output filter	0000 0120 SEC (0000 SEC)	
o1.FAIL	22 mA signal for error messages	ON / OFF	
o1.HoLD	Signal behavior during HOLD	LAST / FIX	
o1.FIX	FIX: Enter FIX value	003.8 022.0 mA (021.0 mA)	
tc	Temperature compensation		
tc.UnIT	Select temperature unit	° C / °F	
tc.rTD	Select temperature probe	22 NTC / 30 NTC	
Corr	Correction		
Co.UPoL	Enter polarization voltage	0400 1000 mV (0675 mV)	
Co.UnIT	Select pressure unit	BAR / KPA / PSI	
Co.PrES	Enter process pressure correction	xxxx (1.013 BAR)	
Co.SAL	Enter salinity correction	00,00 45.00 g/kg (00.00 ppt)	
CAL	Calibration mode		
CA.MoD	Select calibration mode	SAt / Conc	
CA.tiME	Selection of cal timer interval	0000 9999 h / (0000 h)	
ALrt	Alarm settings		
AL.SnSO AL.dLY AL.LED	Select Sensocheck Enter alarm delay LED in HOLD mode	ON / OFF 0000 0600 SEC (0010 SEC) ON / OFF	


Individual settings

(Original for copy)


(Original	101 COP3/		
Code	Parameter	Default settings	Individual settings
o1.FCT	Process medium	DO	
o1.SnSR	Sensor type	Type A	
o1.UnIT	Variable/Unit	%	
o1.4mA	Current start	0000 %	
o1.20mA	Current end	0500 %	
o1.FtME	Filter time	<u>0000 SEC</u>	
o1.FAIL	22mA signal	<u>OFF</u>	
o1.HoLD	Hold behavior	LAST	
o1.FIX	FIX current	021.0 mA	
tc.UnIT	Unit °C / °F	<u>°C</u>	
tc.rTD	Temp probe	22 NTC	
Co.UPoL	Polarization voltage	675 mV	
Co.UnIT	Pressure unit	BAR	
Co.PrES	Process pressure corr.	1.013 BAR	
Co.SAL	Salinity correction	00.00 ppt	
CA.MoD	Calibration mode	SAt	
CA.tIME	Cal timer interval	<u>0000 h</u>	
AL.SnSO	Sensocheck	<u>OFF</u>	
AL.dLY	Alarm delay	<u>0010 SEC</u>	
AL.LED	LED in Hold mode	OFF	

Configuration

Current output Select process medium



- 1: Press conf kev.
- 2: Enter passcode 1200*.
- 3: Select **Current output** menu group using arrow keys. All items of this menu group are indicated by the code "o1."
- 4: Press **enter** to select menu, edit with arrow keys (see Pg 39). Confirm (and proceed) with **enter**.
- 5: To end, press conf, then enter




Code	Display	Action	Choices
о1.		Select configuration (Press conf key).	
	After correct input a wel- come text (CONF) is dis- played for approx. 3 sec	Enter passcode "1200"*. (Select position using arrow key ▶ and edit number using ▲. When the display reads "1200", press enter to confirm.)	
	A HOLD	The Transmitter is in HOLD mode (HOLD icon is on, red LED flashes when "HOLD ON" has been set.).	
	7.5	Select process medium:	DO (GAS)
	10 ♠ d. Fc <u>r</u>	Dissolved oxygen (DO)	
	685 ♠ al. Fc <u>r</u>	• Gas (GAS) Select with arrow key ► Proceed with: enter	


Note: Characters represented in gray are flashing and can be edited.

* Factory setting

Configuration Current output Select sensor type



- 1: Press conf key.
- 2: Enter passcode 1200*.
- 3: Select **Current output** menu group using arrow keys. All items of this menu group are indicated by the code "o1."
- 4: Press **enter** to select menu, edit with arrow keys (see Pg 41). Confirm (and proceed) with **enter**.
- 5: To end, press **conf**, then **enter**



* Factory setting

Code	Display	Action	Choices
o1.	と 	Select sensor type A / B (see table on left-hand side) Select with arrow key ► Proceed with: enter	Type A (InPro6800) Type B (InPro6900)

** Type A sensors (standard applications)


	•		
Sensor type	Screw cap	Sensor current in air (25 °C)	Detection limit
InPro6800	4-pole (T82) VP	50 110 nA typ. 60 nA	0.01 ppm 0.006 ppm

** Type B sensor (traces)


Sensor type	Screw cap	Sensor current in air (25 °C)	Detection limit
InPro6900	VP	typ. 350 nA	0.001 ppm

Note: Characters represented in gray are flashing and can be edited.

Configuration Current output Select process variable / unit



- 1: Press conf kev.
- 2. Enter passcode 1200*.
- 3. Select **Current output** menu group using arrow keys. All items of this menu group are indicated by the code "o1."
- Press enter to select menu, edit with arrow keys (see Pg 43). Confirm (and proceed) with enter.
- 5. To end, press **conf**, then **enter**




^{*} Factory setting

Code	Display	Action	Choices
o1.		Select process variable / unit (valid for all following settings): Select with arrow key ►. Proceed with enter • SAt: Percent saturation: 0.0 199,9 % 200 500 % • Conc: Concentration (µg/l, mg/l, ppb or ppm) 0.00 50.00 mg/l 0.00 50.00 ppm 0000 9999 µg/l	% (μg/l mg/l ppb ppm)
	GAS O/O A OF HOLITON B OF HOLITON	Only with measurement in gas (GAS) selected: Select process variable (valid for all following settings): Select with arrow key ► Proceed with: enter	% (ppm)


Note: Characters represented in gray are flashing and can be edited.

Configuration Current output Current start / end



Code	Display	Action	Choices
o1.		Current start Enter lower end of scale, depending on the measurement procedure selected (Saturation or Concentration) Select with ▶ key, edit number with ▲ key, proceed with enter	0000 % (corresponding to selected range: μg/l mg/l ppb ppm)
	O ∏ 1000/o A ol20 A⊡	Current end Enter upper end of scale, depending on the measurement procedure selected (Saturation or Concentration) Select with ▶ key, edit number with ▲ key, proceed with enter	0500 % (corresponding to selected range: μg/l mg/l ppb ppm)


Assignment of measured values: Current start and current end



Configuration

Current output

Output filter: Time constant



^{*} Factory setting


Code	Display	Action	Choices
o1.	O O O SEC A OLFERE	Default setting: 0 s (inactive).	0000 SEC (0000 - 0120 SEC)

Time constant of output filter (attenuation)

To smoothen the current output, a low-pass filter with adjustable filter time constant can be switched on. When there is a jump at the input (100 %), the output level is at 63 % after the time constant has been reached. The time constant can be set from 0 to 120 sec. If the time constant is set to 0 s, the current output follows the input.

Note:


The filter only acts on the current output, not on the display!



Time constant 0 to 120 sec


Configuration

Current output Output current during Error and HOLD




Code	Display	Action	Choices
о1.	A	22 mA signal for error message Select with arrow key ► Proceed with: enter	OFF (ON)
	O LFST	Output signal during HOLD LAST: During HOLD the last measured value is maintained at the output FIX: During HOLD a value (to be entered) is maintained at the output Select with arrow key Proceed with: enter	LAST (FIX)
	F IX A of Hot III O I	Only with FIX selected: Enter current which is to flow at the output during HOLD Select position with arrow key ▶ and edit number with ▲ key Proceed with: enter	021.0 mA (003.8 022.0 mA)

Output signal for HOLD:



Configuration

Temperature compensation Temperature unit, temperature probe




Factory setting

Code	Display	Action	Choices
tc.	♣ FCTIVITE	Specify temperature unit Select with arrow key ► Proceed with: enter	°C (°F)
	22NTE a od rije	Select temperature probe Select with arrow key ► Proceed with: enter	22NTC (30NTC)

Configuration

Correction: Polarization voltage, process pressure, salinity correction




" Factory setting

Code	Display	Action	Choices
Co.	1675	Enter polarization voltage Select with ►key, edit number with ▲key, proceed with enter	0675 mV (0400 1000 mV)
	JAR ♠ CoUnii <u>e</u>	Select pressure unit Select with arrow key ► Proceed with: enter	bar (kPa, PSI)
	10 13 BAR A Copressed	Process pressure correction Enter process pressure. This value is used to correct oxygen saturation. It has no influence on concentration measurement (Conc). Select position with arrow key and edit number with ▲ key. Proceed with: enter	1.013 bar (0.000 9.999 bar 0.000 999.0 kPa, 0.000 145.0 PSI)
	0000ppt ♠ Co. Sn <u>e</u>	Enter salinity correction Select position with arrow key and edit number with ▲ key. Proceed with: enter	00.00 ppt* (00.00 45.00 ppt)


^{*} ppt (parts per thousand) - corresponds to g/kg

Configuration Calibration mode




- 1: Press conf key.
- 2: Enter passcode 1200.
- 3: Select **Calibration mode** menu group using arrow keys.

 All items of this menu group are indicated by the code "CA."
- 4: Press **enter** to select menu, edit with arrow keys (see Pg 55). Confirm (and proceed) with **enter**.
- 5: To end, press **conf**, then **enter**




Code	Display	Action	Choices
CA.	SAL ACR MODE	Specify calibration mode (Calibration to saturation or concentration) Select with ►key, proceed with enter	SAt (Conc)
	∏∏∏h ≜ [AL, ME	Cal timer interval The cal timer reminds you to calibrate in time. Select with ▶ key, edit number with ▲ key, proceed with enter	0000 h (0000 9999 h)

Configuration Alarm settings



- 1: Press conf key.
- 2: Enter passcode 1200.
- 3: Select **Alarm settings** menu group using arrow keys.
 All items of this menu group are indicated by the code "AL."
- 4: Press **enter** to select menu, edit with arrow keys (see Pg 57). Confirm (and proceed) with **enter**.
- 5: To end, press conf, then enter



Code	Display	Action			Choices
AL.	CHECK A ALSoso	Select Sensoche (continuous mo Select with kenter With sensor type must be switch	onitoring ey, proce be B Sens	ed with	ON / OFF
		Alarm delay Select with ►key, edit number with ▲ key, proceed with enter		0010 SEC (0000 0600 SEC)	
	Kold <u>"</u>	LED in HOLD mode Select with ▶key, edit number with ▲key, proceed with enter LED state:		ON / OFF	
		Setting	Alarm	HOLD	
		ON	on	flashes	
		OFF	flashes	off	

Passcodes according to FDA 21 CFR Part 11

Access to the device functions can be protected with adjustable passcodes if required.

If such a protection is not required, you should use the preset passcodes.

To call up passcode editor: Press **conf** key and enter Administrator passcode (Factory setting: **1989**).

Display	Action	Remark
SPCL _{S7}	1: Press conf key. 2: Enter Administrator passcode (1989): Welcome text is displayed	This text is displayed for approx. 3 sec
	"Cal Info" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 0000
100 (cal ==	"Cal zero" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 1001
	"Calibration Sat/Conc" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 1100
1 10 5 13 13 13 13 13 13 13 13 13 13 13 13 13	"Product calibration" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 1105
-td <u>-</u>	"Temp probe adjustment" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 1015

Display	Action	Remark
	"Error Info" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 0000
	"Configuration" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 1200
6 666668 5666	"Sensor monitor" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 2222
	"Current source" Edit: Arrow keys Proceed with: enter Cancel: conf	Default setting: 5555
1989crs 580.	"Administrator passcode" Edit: Arrow keys Proceed with: enter. Cancel: conf	Default setting: 1989 Caution!
	New "Administrator passcode" Select "NO" / "YES" with arrow keys	If you have lost the Administrator pass- code, the Passcode Editor cannot be
1199 HD SPEL <u>=</u>	"NO" enter = old passcode Cancel: conf = old passcode	called up! Please consult our technical support!
6]	"YES" enter = take over new passcode Cancel: conf = old passcode	

Calibration

Calibration adjusts the device to the sensor.

Activate	cal	Activate with cal	
		Enter passcode • 1001: Zero calibration • 1100: Saturation/Concentration Volume concentration (GAS) • 1105: Product calibration Edit parameter with ▶ and ▲, confirm / continue with enter. (End with cal and enter.)	
Hold	HOLD	The loop current is frozen (at its last value or at a preset fixed value, depending on the configuration), Sensoface is off, mode indicator "Calibration" is on. Red LED flashes when "HOLD ON" has been set.	
During calibration the Transmitter remains in the Hold mode.	☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐		
Err in		The calibration parameters are checked during the input. In the case of an incorrect input "Err" is displayed for approx. 3 sec. The incorrect parameters cannot be stored. Input must be repeated.	
End	enter	End with cal . Safety prompt: The measured value and Hold are displayed alternately, "enter" flashes.	

60

Press **enter** to end the Hold mode. The measured value is displayed. The output current remains frozen for another 20 sec (HOLD icon on, "hourglass" flashes).

Calibration

It is always recommended to calibrate in air.

Compared to water, air is a calibration medium which is easy to handle, stable, and thus safe. In the most cases, however, the sensor must be dismounted for a calibration in air.

When dealing with biotechnological processes which require sterile conditions, the sensor cannot be removed for calibration. Here, calibration must be performed with aeration directly in the process medium (e.g. after sterilization).

In the field of biotechnology, for example, often saturation is measured and calibration is performed in the medium for reasons of sterility.

For other applications where concentration is measured (water control etc.), calibration in air has proved to be useful.

Common combination: process variable / calibration mode

Measurement	Calibration
Saturation	Water
Concentration	Air (synthetic air)
Volume concentration	Air

The calibration procedures for these two common applications are described on the following pages. Of course, other combinations of process variable and calibration mode are possible.

Note:

When a 2-point calibration is required, the zero calibration should be performed prior to saturation or concentration calibration, resp. (see Pg 70).

All calibration procedures must be performed by trained personnel.

Calibration to saturation (SAT)

Display	Action	Remark
	Activate calibration (Press cal .) Enter passcode 1100. Select with ▶ key, edit number with ▲ key, proceed with enter	SAT or Conc calibration is selected during configuration. If an invalid passcode is entered, the Transmitter returns to measuring mode.
	Place sensor in calibration medium Start with enter	Welcome (3 sec) The Transmitter is in the Hold mode.
	Enter relative humidity Select with ▶ key, edit number with ▲ key, proceed with enter	Default for relative humidity in aqueous media: rH = 100 % (in air approx. 50 %)
2.320 JAR	Enter calibration pressure Select with ►key, edit number with ▲key, proceed with enter	Default for calibration pressure is the process pressure configured
60.3 nA ≥265• c ⇒	Automatic drift check Display of sensor current (related to 25 °C and 1013 mbars normal pressure) and measuring temperature. The drift check might take some time.	Drift check can be stopped after > 10 sec by pressing cal (accuracy reduced).

Display	Action	Remark
	Enter desired value for saturation Select with ►key, edit number with ▲ key, proceed with enter	Default: last value entered
© 60.5 _{.4} <u>× ≜</u> 000s ₪	Display new slope and zero (related to 25°C and 1013 mbars). End calibration with enter.	
1002°/°° <u>*</u> 265°©	Place sensor in process. The percent saturation is shown in the main display alternately with "Hold"; "enter" flashes. Stop Hold with enter .	After end of calibration, the outputs remain in Hold mode for approx. 20 sec.

Information on saturation calibration (SAT)

- The calibration medium must be in equilibrium with air (percent saturation for water is 100 %). Oxygen exchange between water and air is very slow.
 To speed up the adjustment processes, make sure that there is a steady medium flow during calibration.
- If the percent saturation is known from a simultaneous measurement, it can be entered manually.
- For 2-point calibration, perform zero point calibration first, see Pg. 70

Calibration to concentration (Conc)

Display	Action	Remark
	Activate calibration (Press cal .) Enter passcode 1100. (Press ► key to select position, enter number using ▲ key, confirm with enter)	SAT or Conc calibration is selected during configuration. If an invalid passcode is entered, the Transmitter returns to measuring mode.
EAL	Place sensor in air Start with enter	The Transmitter is in the Hold mode
1050°/°	Enter relative humidity (Press ▶ key to select position, enter number using ▲ key, confirm with enter)	Default for relative humidity in air: rH = 50 %
RREEL O.	Enter calibration pressure (Press ► key to select position, enter number using ▲ key, confirm with enter)	Default for calibration pressure is normal pressure 1.013 bars.
60.3∩R □ 265•c ==	Automatic drift check Display of input current (related to 25 °C and 1013 mbars) and measuring temperature. The drift check might take some time.	Drift check can be stopped after > 10 sec by pressing cal (accuracy reduced).

Display	Action	Remark
[B.32]pm	Enter default for concentration (Press ► key to select position, enter number using ▲ key, confirm with enter)	Default value is calculated from rel. humidity, cal pressure and cal temperature.
		(The unit of measurement, ppm or mg/l, is preset during configuration.)
© 60.5 m	Display of new slope and zero (related to 25 °C and 1013 mbars)	
	Press enter to end concentration calibration.	
0875°ppm <u>*</u> 26.5°cr=	Place sensor in process The new value is shown in the main display alternately with "Hold"; "enter" flashes. End with enter .	After end of calibration, the outputs remain in Hold mode for approx. 20 sec.

Information on concentration calibration (Conc)

Calibration in air. This calibration method is recommended when the sensor can be removed for calibration. Air has a stable oxygen content. Therefore the adjustment processes during calibration run more quickly.

• For 2-point calibration, perform zero point calibration first, see Pg. 70

Calibration to volume concentration (GAS) Calibration medium: air

(Pr En (Pr	ctivate calibration Press cal .) Inter passcode 1100. Press ► key to select position,	GAS must have been selected
en	nter number using ▲ key, onfirm with enter)	during configura- tion. If an invalid passcode is entered, the Transmitter returns to measuring mode.
[AL Pla	lace sensor in air	Welcome (3 sec) the Transmitter is in the Hold mode.
[inter relative humidity Press ►key to select position, enter number using ▲ key, confirm with enter)	Default for relative humidity in air: rH = 50 %
(P	Enter calibration pressure Press ► key to select position, enter number using ▲ key, confirm with enter)	Default for calibration pressure is normal pressure 1.013 bars.
Fig. 265cm D. C. (ru ar	Automatic drift check Display of input current related to 25 °C and 1013 mbars) and measuring temperature. The drift check might take some time.	Drift check can be stopped after > 10 sec by pressing cal (accuracy reduced).

Display	Action	Remark
© 60.5 na	Display of new slope and zero (related to 25 °C and 1013 mbars)	
	Press enter to end concentration calibration.	
<u> </u>	Place sensor in process The new value is shown in the main display alternately with "Hold"; "enter" flashes. Stop Hold with enter .	After end of calibration, the outputs remain in Hold mode for approx. 20 sec.

Please note:

• For 2-point calibration, perform zero point calibration first, see Pg. 70

Product calibration

Calibration with sampling

Can be performed for all process variables: saturation, concentration, volume concentration.

During product calibration the sensor remains in the process.

The measurement process is only interrupted briefly.

Procedure: During sampling the currently measured value is stored in the Transmitter. The Transmitter immediately returns to measuring mode.

The calibration mode indicator flashes and reminds you that calibration has not been terminated. The comparison value is measured on the site, e.g. using a portable DO meter in a bypass. This value is then entered in the Transmitter. The new value for slope or zero is calculated from the stored value and the comparison value. From the measured value, the Transmitter automatically recognizes whether a new slope or zero must be calculated (above approx. 5 % saturation: slope, below: zero).

If the sample is invalid, you can take over the measured value stored during sampling instead of the comparison value. In that case the old calibration values remain stored. Afterwards, you can start a new product calibration.

The following describes a product calibration with slope correction – a product calibration with zero correction is performed correspondingly.

Display	Action	Remark
	Product calibration step 1: Activate calibration (Press cal key). Enter passcode 1105. (Select with ▶ key, edit number with ▲ key, proceed with enter)	The type of product calibration (SAT or Conc) is selected during configuration (Process variable). If an invalid passcode is entered, the Transmitter returns to measuring mode.
		Display for approx. 3 sec

METTLER TOLEDO

Display	Action	Remark
06.20 PPm SkorE ==	Take sample and store the currently measured value. Proceed with enter	Now the comparison value must be determined. The Transmitter goes to measuring mode.
	Measuring mode	From the flashing CAL mode indicator you see that product calibration has not been terminated.
EAL PR] <u>△</u> 2	Product calibration 2nd step: When a comparison value has been determined, call up the product calibration once more (cal key, passcode 1105).	Display (approx. 3 sec)
05.6 3 PPm 	Enter the comparison value. Confirm with enter .	Calculation of new slope.
© 60.5 m <u>∞ 4</u> 000s m	Display of new slope and zero point (related to 25 °C at 1013 mbars) End calibration with enter	New calibration: Press cal key.
05.6 3 ppm <u>*</u> 283 c =	The measured value is shown in the main display alternately with "Hold"; "enter" flashes. Stop Hold with enter .	After end of calibration, the outputs remain in Hold mode for approx. 20 sec.

Zero calibration

Zero calibration

The Series InPro6800/InPro6900 sensors have a very low zero point current. Therefore, a zero point calibration is only recommended for measurement of oxygen traces. If a zero calibration is performed, the DO sensor should remain for at least 10 to 30 minutes in the calibration medium in order to obtain stable, non-drifting values (InPro6900: at least 60 min). During zero point calibration, a drift check is not performed. Zero point current of a properly functioning sensor is notably less than 0.5 % of air current. The display (secondary: measured value, main: entered value) does not change until an input current is entered for the zero point. When measuring in an oxygen-free medium, the displayed current can be taken directly.

Display	Action	Remark
	Activate calibration (Press cal key). Enter passcode 1001. Select with ▶ key, edit number with ▲ key, proceed with enter	The Transmitter is in the Hold mode. If an invalid pass- code is entered, the Transmitter returns to measuring mode.
ZRO	Place sensor in oxygen-free medium	Welcome (3 sec)
0.005 **	Main display: Zero point current; store with enter or correct with arrow keys and then store with enter . Lower display: Sensor current measured	
© 6 0.5	Display of slope Display of new zero point current End calibration with enter key, place sensor in process	
7 14º/₀ <u>△</u> 249•c	The oxygen value is shown in the main display alternately with "Hold"; "enter" flashes. Stop Hold with enter .	The outputs remain in Hold mode for approx. 20 sec.

Temperature probe adjustment

Display	Action	Remark
	Activate calibration (press cal key) Enter passcode 1015. Select with ►key, edit number with ▲ key, proceed with enter	Wrong settings change the measurement properties! If an invalid passcode is entered, the Transmitter returns to measuring mode.
Total	Ready for calibration	The Transmitter is in the Hold mode. Display for approx. 3 s
1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Measure the temperature of the process medium using an external thermometer. Enter measured temperature value: Select with ▶ key, edit number with ▲ key, proceed with enter.	Default: Current value of secondary display
	End adjustment with enter . HOLD will be deactivated after 20 sec.	

Measurement

Display	Action / Remarks
© 978°/° 243°C	In the measuring mode the main display shows the configured process variable (%, mg/l or ppm), the lower display shows the temperature. During calibration you can return to measuring mode by pressing the cal key, during configuration by pressing the conf key. (Waiting time for measured value stabilization approx. 20 sec).

Diagnostics functions

Display	Action / Remarks
O 1794	Display of the output current Press enter while in measuring mode. The output current is shown in the main display. After 5 sec the Transmitter returns to measuring mode.
0000 * 060.5 nR 0005.4 mg	Display of calibration data (Cal Info) Press cal while in measuring mode and enter pass- code 0000*. The slope is shown in the main display, the zero point current in the secondary display. After 20 sec the Transmitter returns to measuring mode (immediate return at pressing enter).
2222 * 70.2 nR 332 °C=	Display of sensor current (Sensor monitoring for validation of sensor and complete measured-value processing) Press conf while in measuring mode and enter passcode 2222*. The (uncompensated) sensor current is shown in the main display, the measuring temp in the secondary display. Press enter to return to measurement.
0000 * LASE Errose	Display of last error message (Error info) Press conf while in measuring mode and enter passcode 0000*. The last error message is displayed for approx. 20 sec. After that the message will be deleted (immediate return to measurement at pressing enter).

^{*} Factory setting

Diagnostics functions

Display	Action / Remarks
5555 * O D D L 1 1 1 1 1 1 1 1 1 1 1 1	Specify output current for testing the connected peripherals Press conf while in measuring mode and enter passcode 5555*. The actually measured current is shown in the secondary display. The output current indicated in the main display can be modified. Select with ► key, edit number with ▲ key. Confirm with enter key. Then the entered value will be shown in the secondary display. The Transmitter is in Hold mode. Press conf, then enter to return to measurement (Hold remains active for another 20 sec).

Cleaning

To remove dust, dirt and spots, the external surfaces of the device may be wiped with a damp, lint-free cloth. A mild household cleaner may also be used if necessary.

Factory setting

Operating states

Operating state	Out	TED .	Time out
Measuring			
Cal Info (cal) 0000*			20 sec
Error Info (conf) 0000*			20 sec
Calibration (cal) 1100*			
Zero calibration (cal) 1001*			
Temp adjustment (cal) 1015*			
Product calibration step 1 (cal) 1105* step 2 (cal) 1105*			
Configuration (conf) 1200*			20 min
Sensor monitor (conf) 2222*			20 min
Current source (conf) 5555*			20 min

Explanation:

active

as configured (Last/Fix or Last/Off)

LED flashes during HOLD (configurable)

* Factory setting

Error messages (Error Codes)

Error	Display	Problem Possible causes	Red LED	Out 1 (22 mA)
ERR 01	Measured value flashes	SAT range Sensor defective Wrong sensor connected Measurement range exceeded	х	Х
ERR 02	Measured value flashes	Conc range Sensor defective Wrong sensor connected Measurement range exceeded	х	х
ERR 98	"Conf" flashes	System error Configuration or calibration data defective; completely reconfigure and recalibrate the device Memory error in device program (PROM defective)	х	х
ERR 99	"FAIL" flashes	Factory settings EEPROM or RAM defective This error message only occurs in the case of a total defect. The Transmitter must be repaired and recalibrated at the factory.	х	х

Error	Icon (flashes)	Problem Possible causes	Red LED	Out 1 (22 mA)
ERR 03		Temperature probe Open or short circuit Temperature range exceeded	Х	х
ERR 11	(IIA)	Current output Current below 0 (3.8) mA	х	Х
ERR 12	(1)A	Current output Current above 20.5 mA	х	х
ERR 13	(1)A	Current output Current span too small / too large	х	Х
ERR 33	\$	Sensocheck Sensor: Connecting cable defective	Х	х
	Zero ●	• Zero error, Sensoface active, see Po	g 80	
	Slope Zero	• Slope error, Sensoface active, see Pg 80		
	Ŀ	Response time exceeded, Sensoface active, see Pg 80		
	M	• Calibration interval expired, Sensof see Pg 80	ace acti	ve,

Calibration error messages

Symbol flashes:	Problem Possible causes
Slope Zero	Slope out of range Wrong calibration values specified (relative humidity, pressure, saturation, concentration) Wrong calibration medium
(L)	Calibration aborted after 12 minutes •Sensor defective or dirty
In addition "CAL Err" flashes.	No electrolyte in the sensor Sensor cable insufficiently shielded or defective Strong electric fields influence the measurement Temperature fluctuation of calibration solution

Sensoface

(Sensocheck must have been activated during configuration.)

The little smiley in the display (Sensoface) alerts to sensor problems (defective cable, maintenance required).

The permitted calibration ranges and the conditions for a friendly, neutral, or sad Sensoface are summarized in the following chart. Additional icons refer to the error cause. Replace membrane module or filling solution, if required.

Type A sensors (InPro6800)

	Slope	Zero point	Response time	Cal timer
Adm. Range	25 130 nA	-2 +2 nA	max. 720 sec	
<u></u>	> 35 < 90 nA	> - 0,5 < 0.5 nA	≤ 300 sec	≤ 80 % expired
<u></u>	Slope	Zero ●	(
	30 35 nA or 90 110 nA	-1.00.5 nA or +0.5 +1.0 nA	300 600 sec	80 ≤ 100 % expired
<u></u>	Slope	Zero •	<u>(</u>)
	< 30 nA or > 110 nA	< -1.0 nA or > +1.0 nA	> 600 sec	Timer expired

Note

The worsening of a Sensoface criterion leads to the devaluation of the Sensoface indicator (Smiley becomes "sad").

An improvement of the Sensoface indicator can only take place after calibration or removal of a sensor defect.

Type B sensor (InPro6900)

	Slope	Zero point	Response time	Cal timer
Adm. Range	200 550 nA	-2 +2 nA	max. 720 s	
<u></u>	> 250 < 500 nA	> - 0.5 < 0.5 nA	< 300 sec	< 80 % expired
<u> </u>	225 250 nA or 500 525 nA		300 600 sec	80 ≤ 100 % expired
:	SlopeZero< 225 nA or> 525 nA	Zero • < -1.0 nA or > +1.0 nA	> 600 s	Timer expired





Thermometer and Sensoface:

Temperature out of concentration or saturation range

Sensocheck

Continuously monitors the sensor and leads for short circuits or open circuits. Critical values make the Sensoface "sad" and the corresponding icon flashes:

The Sensocheck message is also output as error message Err 33. The alarm contact is active, the red LED is lighted, output current 1 is set to 22 mA (when configured correspondingly). Sensocheck can be switched off during configuration (then Sensoface is also disabled). Exception: After a calibration a Smiley is always displayed for confirmation.

Note

With sensor type B Sensocheck must be switched off!

Appendix

Product line and accessories

Devices	Order no.
O ₂ Transmitter 4100 e/2H	52 121 215
O ₂ Transmitter 4100 e/2XH	52 121 168
Mounting accessories	
Pipe-mount kit	52 120 741
Panel-mount kit	52 120 740
Protective hood	52 120 739

Sensors

Mettler-Toledo GmbH, Process Analytics offers a wide range of sensors for the following fields of applications:

- Chemical process industry
- Pharmaceutical industry
- Food and beverage industry
- Water/waste-water

For more information concerning our sensors and housings program, please refer to http://www.mt.com.

Specifications

DO input	Sensor Type A: Sensor Type B:	InPro6800 InPro6900		
Measuring current Measurement error ^{1,2,3)} Ranges*	0 1200 nA, 0.5%m.val. + 0.05 n/ Saturation (-10 80 0.0 199.9 % / 200 (autom. switchover in Concentration (-10	°C) 500 % display) .80 °C) 0.00 50.00 mg/l 0.00 50.00 ppm 0000 9999 µg/l 0000 9999 ppb		
	Volume concentration in gas (-10 80 °C) 0000 9999 ppm 0.0 120 % Display: (0.00 29.99 % / 30.0 120.0 %)			
Adm. guard current Polarization voltage* Process pressure* Salinity correction* Sensor standardization	≤ 20 μA 400 1000 mV 0.000 9.999 bars (00.00 45.00 g/kg	999.9 kPa / 145.0 PSI)		
Operating modes*	 O₂ saturation (automatic) O₂ concentration (automatic) Volume concentration (gas) Product calibration Zero calibration 			
Calibration range Sensor type A	Zero point Slope	± 2 nA 25 130 nA (at 25°C, 1013 mbars)		
Calibration range Sensor type B	Zero point Slope	± 2 nA 200 550 nA (at 25°C, 1013 mbars)		

0000 ... 9999 h

84

Cal timer*

Pressure correction*

0.000 ... 9.999 bars (... 999.9 kPa / ... 145.0 PSI)

Sensocheck Monitoring for short circuits / open circuits

(can be disabled), delay: 30 sec

Sensoface Provides information on the sensor condition

evaluation of zero point/slope, response time,

calibration interval, Sensocheck

Sensor monitor Direct display of measured values from sensor for

validation (uncompensated sensor current.

measuring temp)

Temperature input* NTC 22 kOhm / NTC 30 kOhm

2-wire connection, adjustable

Range -20.0 ... +150.0 °C / -4 ... +302 °F

Adjustment range 10 K

Resolution 0.1 °C / 1 °F

Measurement error $^{1,2,3)}$ < 0.5 K (< 1 K at T > 100°C)

Supply/Output

Loop current 4 ... 20 mA (22 mA), floating

(3.8 ... 20.5 mA)

Supply voltage 12 ... 30 V, $I_{max} = 100$ mA, $P_{max} = 0.8$ W (Ex)

Process variable* O₂ saturation / O₂ concentration

Characteristic Linear

Overrange*) 22 mA in the case of error messages

Output filter* Low-pass, PT₁, filter time constant 0 ... 120 sec

Measurement error 1) < 0.3 % current value + 0.05 mA

Start/end of scale As desired within range

Adm. span 2 % ... 500 % Gas: 500 ... 9999 ppm

200 ... 9999 μg/l 1 ... 120 %

200 ... 9999 ppb 0.5 ... 50 mg/l 0.5 ... 50 ppm

Current source function 3.8 mA ... 22 mA

Specifications

HART communication Digital communication by FSK modulation

of loop current, reading of device identification,

measured values, status and messages, reading and writing of parameters, start of product calibration, signaling of configuration changes according to

FDA 21 CFR Part 11

Display LC display, 7-segment with icons

Main display Character height 17 mm, unit symbols 10 mm
Secondary display Character height 10 mm, unit symbols 7 mm

Sensoface 3 status indicators (friendly, neutral, sad)

Mode indicator 4 indicators: "meas", "cal", "alarm", "config"

18 further icons for configuration and messages

Alarm indication Red LED in case of alarm or HOLD, user defined

Keypad 5 keys: [cal] [conf] [▶] [▲] [enter]

Service functions

Current source Loop current specifiable 3.8 ... 22.00 mA
Device self-test Automatic memory test (RAM, FLASH, EEPROM)

Display test Display of all segments
Last Error Display of last error occurred

Sensor monitor Display of direct, uncorrected sensor signal

(sensor current/temperature)

Passcodes Modifiable according to FDA 21 CFR Part 11

"Electronic Signatures"

* User-defined

1) To IEC 746 Part 1, at nominal operating conditions

2) ± 1 count

3) Plus sensor error

Data retention Parameters and cal data > 10 years (EEPROM)

EMC EN 61326

Emitted interference: Class B (residential area)

Class A

Immunity to interference: Industry

Lightning protection EN 61000-4-5, Installation Class 2

Explosion protection

4100 e/2XH: ATEX: TÜV 04 ATEX 2431

II 2(1) G EEx ib[ia] IIC T6

FM: FMRC 3023119

IS/I/1/ABCD/T4; Entity; Type 2 I/0/AEx ia IIC T4; Entity; Type 2

NI/I/2/ABCD/T4; Type 2 AIS/I, II, III/1/ABCDEFG

CSA: 1662790

CI I, Div 1, Gr ABC & D T4; Ex ib [ia] IIC T4

CI I, Div 2, Gr ABC & D, T4; Ex nAL[L] IIC T4

4100 e/2H: FM: FM 300580 / FM 3023119

NI/I/2/ABCD/T4

Nominal operating conditions

Ambient temperature -20 ... +55 °C Transport/Storage temp -20 ... +70 °C

Relative humidity 10 ... 95% not condensing

Supply voltage 12... 30 V

Specifications

Protection

Enclosure Molded enclosure made of PBT

(polybutylene terephtalate)

Color Bluish grav RAL 7031 Assembly · Wall mounting

· Pipe mounting:

Ø 40 ... 60 mm,

30 to 45 mm

• Panel mounting, cutout to DIN 43 700

Sealed against panel

H 144 mm, W 144 mm, D 105 mm Dimensions

IP 65/NFMA 4X

(USA, Canada: indoor use only) Cable glands


3 breakthroughs for M20x1.5 cable glands


2 breakthroughs for NPT 1/2" or

Rigid Metallic Conduit

Approx. 1 kg Weight

FM Control Drawing





Explosion protection



Certificate of Compliance

Certificate: 1662790 Master Contract: 220331

Project: 1662790 Date Issued: May 18, 2005

Issued to: Mettler-Toledo GmbH

Im Hackacker 15 Urdorf, 8902 SWITZERLAND

Attention: Mr. Michael Haas

The products listed below are eligible to bear the CSA Mark shown



Issued by: K. Atkins

Authorized by: Nick Alfano

erations Manager

PRODUCTS

CLASS 2258 04 - PROCESS CONTROL EQUIPMENT - Intrinsically Safe Entity - For Hazardous Locations

Class I, Division I, Groups A, B, C and D

Ex ib fial IIC

Transmitters Models 2100/2XH, 7100/2XH and Models 4100/2XH, input rated 30V, 4-20 mA, intrinsically safe devices provides intrinsically safe outputs to simple apparatus, ph, conductivity and oxygen probes when connected per control drawings 194.120-170, 194.220-190 and 194.320-190, 194.401-120. Maximum Ambient Temperature 55°C, Temperature Code 74.

For all models the input entity parameters are:

Terminals Ui, Vmax Ii, Imax Pi, Pmax Ci Li 10, 11 or 14.15 30V 100mA 0.8W 32.4nF 0.24mH

Output entity parameters are:

2100/2XH

DQD 507 Rev. 2003-01-31

Warnings and notes to ensure safe operation

Warning: Do not disconnect equipment unless power has

been switched off.

Warning: Clean only with antistatic moistened cloth.

Warning: Substitution of components may impair suitability


for hazardous locations.


 The equipment shall be installed and protected from mechanical impact and ultraviolet (UV) sources.

- Clean only with a moistened antistatic cloth as potential electrostatic hazard may exist. Service equipment only with conductive clothing, footwear and personal grounding devices to prevent electrostatic accumulation.
- Internal grounding provisions shall be provided for field wiring.
 Bonding between conduit shall be provided during installation, and all exposed non-current carrying metallic parts shall be bonded and grounded.
- Installation in a Class I, Division 2 or Class I, Zone 2 hazardous location shall be in accordance with the Canadian Electrical Code (CEC Part 1) Section 18 Division 2 wiring methods.

OBSERVE THE SPECIFICATIONS OF THE CONTROL DRAWING!

CSA Control Drawing





Index

2-point calibration		
A Alarm		 56 83 14 47
Calibration	Conc)	

Index

Configuration: Current output	
Current start / end	
Output current during Error 48	
Output signal during HOLD49	
Process medium	
Process variable / unit	
Sensor type41	
Temperature compensation	
Time constant of output filter	
Connection	
Contents	
Control Drawing	
CSA Control Drawing	
Current source	
Current start/end	
Current starvena	
D	
Diagnostics functions	
Display of calibration data	
Display of last error message	
Display of output current	
Display of sensor current	
Output current fixed	
Display	
_	
E	
EC Declaration of Conformity	
EC-Type-Examination Certificate	
Error messages	
Calibration error messages	
Display of last error message	
Explosion protection	
Cleaning in a hazardous location 6	

F FDA 21 CFR Part 11 9 FM Control Drawing 90
H HART communication 31 Hold mode 30 Configuration 48 LED in HOLD mode 56 Output signal for HOLD 49
lnstallation
K Keypad
Loop current
Mounting plan
Operating states .75 Output filter .46 Overview .2, 13
P Package contents

Index

Panel-mount kit 17 Parameter set - individual settings 37 Passcode editor 58 Pipe mounting 16 Pipe-mount kit 16 Polarization voltage 52 Power supply 6 Process pressure 52 Product calibration 68 Product line and accessories 83 Protective hood 16
Safety functions 29 Safety information 5 Installation 6 Salinity correction 52 Self test 29 Sensocheck 29, 81 ON / OFF 57 Sensoface 80 Calibration ranges 80 Specifications 84 Supply units 6
T50Temperature compensation50Temperature probe adjustment72Terminal assignments19Time constant of output filter47
U User interface

V, VV	
Volume concentration (Gas)	66
VP cable	24
Connection sensor/VP cable	24
Connection VP cable/meter	23
Connector assignment	25
Wiring examples	22
Sensors with connection via VP cable 23, 24,	25
Wiring examples	22
_	
Z	70
Zero calibration	70

BR Mettler-Toledo Ind. e Com. Ltda.,

Alameda Araguaia, 451 - Alphaville BR - 06455-000 Barueri / SP, Brazil Phone +55 11 4166 74 00 Fax +55 11 4166 74 01

CH Mettler-Toledo (Schweiz) AG,

Im Langacher, CH-8606 Greifensee, Switzerland Phone +41 44 944 45 45 Fax +41 44 944 45 10

D Mettler-Toledo GmbH, Prozeβanalytik,

Ockerweg 3, D-35396 Gießen, Germany Phone +49 641 507-333 Fax +49 641 507-397

F Mettler-Toledo Analyse Industrielle Sàrl,

30 Bld. de Douaumont, BP 949, F-75829 Paris Cedex 17, France Phone +33 1 47 37 06 00 Fax +33 1 47 37 46 26

USA Mettler-Toledo Ingold, Inc.,

36 Middlesex Turnpike, USA - Bedford, MA 01730, USA Phone +1 781 301-88 00 Fax +1 781 271-06 81





Subject to technical changes.

© Mettler-Toledo GmbH, Process Analytics
06/05 Printed in Switzerland. 52 121 169

Mettler-Toledo GmbH, Process Analytics Industrie Nord, CH-8902 Urdorf, Switzerland Phone + 41 44 736 22 11, Fax +41 44 736 26 36

www.mtpro.com