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Operating instructions for pH-electrodes InPro® 426X and InPro® 426X i

1. Introduction

The METTLER TOLEDO InPro® 426X pH-single rod measuring cell is a low maintenance and high-pressure resistant pH electrode with a solid polymer electrolyte (XEROLYT EXTRA) together with a silver ions barrier. The diaphragm is of an open aperture type, so that direct contact exists between the medium being measured and the electrolyte. Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. We recommend that you only operate the electrode in conjunction with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

2. Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn. Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

3. Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO InPro® 426X i /yyy/zzz	manufacturer of the electrode Type designation X = membrane glass, i = ISM Digital /SG = auxillary electrode / yyy = rod length in mm/zzz = type of temperature sensor/ISM (Analog)
combination pH pH 0 ...14 0...100 °C (130 °C) Order No. 30 04X XXX Observe the operating instructions	Type of electrode pH-measurement range temperature range for operation (for sterilisation) order number instruction (see also: instruction manual)

In addition each electrode is provided with a serial number on the red part of the plug-in head in order to enable identification.

4. Installation and commissioning

- While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
- Remove the protective cap and rinse the electrode for a short time with deionised water. If the electrode has already been used the sensitive part can be cleaned if necessary with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionised water. After rinsing the electrode should only be dabbed. Any rubbing of the pH-sensitive glass can increase the response time as a result of electrostatic charge accumulation.
- Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode lightly in a vertical direction.
- Place the electrode into the armature as described in the instructions for the armature.
- Connect the electrode and the pH transmitter with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.
- Intelligent ISM sensors such as InPro 426X i enable «plug and measure» and expanded diagnostics. Operation of the electrodes in the ISM mode requires the M700 transmitter (software) as well as an ISM pH-module (2700 i X). For installation, commissioning and operation of the ISM system please observe the instructions for the transmitter, module and cable.

5. Operation: calibration of the electrode and pH transmitter

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connecting it to the transmitter. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the pH transmitter.

6. Maintenance

- After each working cycle the electrode tip has to be carefully rinsed with deionised water. Any drying out of residues from the solution being measured on these parts must be avoided at all costs!
- If the electrode is not in use it is stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9848 FriscoIy or 9823/3M KCl).
- If the electrode is stored in its armature the rules described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the armature parts.
- If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
- The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionised water or alcohol and then dry it off carefully.

7. Removal of sources of malfunction

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used.

8. Disposal

It is the user’s responsibility to dispose of the sensor in a professional manner. The sensor can contain electronic components that require correct disposal such that no persons or the environment are endangered.

9. Garantie

On manufacturing defects 12 months after delivery.

10. Storage conditions

We recommend that the electrodes are stored at room temperature.

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Subject to technical changes. 05./2012. © Mettler-Toledo AG. Printed in Switzerland.
D Mettler-Toledo GmbH, Prozeßanalytik, Ockerweg 3, D - 35396 Gißen, Germany, Phone +49 641 507-333, Fax +49 641 507-397